

Well No.																																										
HLA-B allele ³	1	2	3	4	5	6	7	8	9	0	1	2	3	4	5	6	7	8	9	0	1	2	3	4	5	6	7	8	9	0	1	2	3	4	5	6	7	8	9	0		
*1542	+	+																																								
*1543	+	+	+																			+	+																			
*1544		+			+	+							+																			+										
*1545	+	+	+								+			+									+																			
*1546	+	+	+				+																+																			
*1547	+						+																																+			
*1548	+	+											+										+									+										
*1549	+						+					+																											+			
*1550	+	+	+																														+									
*1551					+		+						+										+													+						
*1552					+		+		+				+																													
*1553	+		+				+																+			+																
*1554	+		+				+																+																			
*1555	+	+			+				+		+			+																												
*1556		+	+									+											+																			
*1557	+	+	+																				+																			
*1558	+	+	+																				+																			
*1560	+	+	+																				+																		+	
*1561	+						+																+																+			
*1562	+					+	+																+																	+		
*1563	+	+	+										+										+																			
*1564	+				+		+																+																			
*1565	+	+																					+															+				
*1566		+	+									+											+																			
*1567							+											+		+	+	+	+																	+		
*1568			+								+		+										+																			
*1569	+						+																									+										
*1570	+	+	+																				+																+			
*1571	+	+	+																				+																		+	
*1572					+		+					+											+			+																
*1573	+	+	+																																	+	+					
*1574							+																+																			
*1575	+	+	+																				+																			
*1576		+	+		+																		+																+			
*1577	+	+	+			+																																				
*1578	+	+																					+		+																	
*1579N	+	+	+									+											+																			
*1580					+	+	+					+											+																			
*1581	+	+	+																				+																			
*1582	+	+	+																				+															+				
*1583	+	+										+		+																		+										
*1584	+	+							+														+															+				
*1585	+	+	+			+																	+		+																	
*1586	+	+																							+							+										
*1587	+	+	+																+				+																			
*1588	+	+			+	+																																				
*1589	+	+			+	+																														+						
*1590			+		+		+					+	+										+																			
*1591	+						+		+														+																			
*1592	+	+	+																				+																			
*1593					+		+					+																				+										
*1594N	+	+	+																				+																			
*1595							+												+		+	+	+	+															+	+		
*1596	+	+	+																				+																			
*1597		+	+																				+																			
*1598	+						+																+																			

[illegible]

[illegible]

HLA-B allele ³	1	2	3	4	5	6	7	8	9	0	1	2	3	4	5	6	7	8	9	0	1	2	3	4	5	6	7	8	9	0	1	2	3	4	5	6	7	8	9	0
---------------------------	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Well No.																																											
HLA-B allele ³	1	2	3	4	5	6	7	8	9	0	1	2	3	4	5	6	7	8	9	0	1	2	3	4	5	6	7	8	9	0	1	2	3	4	5	6	7	8	9	0			
*3505, 3522, 3551, 3572, 4071, 5314, 5818									+		+											+																					
*351401-351402, 3543, 3562, 3567, 3579			+																			+																					
*3529, 3569, 3576, 5403, 5601, 560501-5606, 5608, 5609, 5611, 5613-5617, 5620, 5621, 7801, 7803																						+																		+			
*3532, 356802, 4063, 4092, 5816									+													+																					
*3544			+																			+																					
*3545																				+		+																	+				
*3546, 570101-570102, 570301-570302, 5706-5710, 5714, 5715, 5717																						+																					
*3558											+											+																					
*3566											+											+																	+				
*3571, 4443, 5807																				+																							
*3903, 3914, 3924, 3929, 3937, 4454											+																																
*390601-390602, 3934, 4086, 4502, 4905, 5903, 7301								+																																			
*3908, 3911																							+																	W			
*3932, 5145																						+						+															
*3933								+														+																	+				
*3936			+																																			+					
*400106									+					+																													
*4012	+						+		+														+																				
*4026														+									+																				
*4052								+												+				+																			
*4059																			+																								
*4077, 4087								+																																			
*4209																																						+					
*4417																				+			+															+					
*4422																																											
*4509								+													+		+														+						
*460101, 4604, 4607N, 4610, 4614-4616	+	+																				+																			W		
*460102, 4602, 4605	+	+																				+																					
*4603	+																																			+		W					
*4606	+							+						+																										W			
*4608	+																					+																		W			
*4609	+	+																				+																		W			
*4611	+																			+		+														+				W			
*4612	+	+									+			+								+																		W			
*4613	+							+														+																	</				

[illegible]

[illegible]

																														Well No.									
																														HLA-B allele ³									
																														*0745, 0752									
																														*0748, 1310, 1311, 1810, 1811, 2707, 2711, 2733, 2734, 2743, 3531, 3930, 400101-400103, 400105, 400201-4003, 4007-4011, 4013-401403, 4018-4022N, 4024, 4025, 4027, 4029, 4031, 4033, 4036-4040, 4042, 4045-4050, 4054-4058, 4061, 4062, 4065-4067, 4069, 407201- 4074, 4076, 4078-4082, 4084, 4085, 4088 4091, 4431, 480101-480102, 480301, 4804, 4806, 4807, 4809-4811, 4813, 4816, 4818, 8101-8104N									
																														*0765, 5410, 5509									
																														*080101-080105, 0804, 0805, 0807-0812, 0814-0819N, 0822, 0824, 0827, 0829- 0834, 0836, 4506, 5503, 5512, 5516, 5525									
																														*0806, 0820									
																														*0813, 4043									
																														*0821									
																														*0825									
																														*0826, 3904, 5401, 5404, 5405N, 5407, 5408N, 5411-5413, 5415, 550101- 550202, 5505, 5507, 5511, 5514, 5515, 5517, 5519, 5520, 5522, 5524, 5526- 5530, 5612, 5619N, 8301									
																														*0828, 0835									
																														*1301, 1306, 1307N, 1312, 1313, 1317, 1320-1322, 4408, 4457, 4460, 570103, 5704, 5705									
																														*1303, 350101-350403, 3506, 3507, 350901-3513, 3516, 3517, 3519-3521, 3523-3528, 3530, 3534, 3536, 3538- 3542, 3547-3550, 3552-3557, 3559, 3563, 3564, 356801, 3570, 3574, 3575, 3577, 3578, 3581-3585, 3590-3592, 3594, 4005, 4028, 4410, 4414, 4445, 4451, 4462, 4802, 510101-5104, 5107- 510902, 5111N, 5112, 5114, 5115, 5119, 5120, 5122-512402, 5126-5128, 5130, 5132, 5133, 5135, 5138-5141N, 5143, 5144N, 5146, 5148-5153, 5155-5157, 520101-5207, 5209-5213, 530101-5306, 5308-5310, 5312, 5313, 5315, 5607, 5711, 5716, 580101-5802, 5804, 5806, 5808-5815, 5817N, 5819, 780201- 780202, 7804-7806, 8201, 8202									
																														*1304, 3537, 4904, 511301-511302, 5137, 5214									
																														*1405									
																														*140601, 1815, 4051									
																														*180101-180103, 1803-1807, 1809, 1812, 1817N, 1820, 1823N-1826, 1828, 1829, 3943									
																														*1813									
																														*1814									
																														*1819									
																														*2725, 3586, 4505									
																														*2732, 4035									
																														Well No.									

1	2	3	4	5	6	7	8	9	0	1	2	3	4	5	6	7	8	9	0	1	2	3	4	5	6	7	8	9	0	HLA-B allele ³
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---------------------------

Well No.																															
HLA-B allele ³																															
*3505, 3522, 3551, 3572, 4071, 5314, 5818																															
*351401-351402, 3543, 3562, 3567, 3579																															
*3529, 3569, 3576, 5403, 5601, 560501-5606, 5608, 5609, 5611, 5613-5617, 5620, 5621, 7801, 7803																															
*3532, 356802, 4063, 4092, 5816																															
*3544																															
*3545																															
*3546, 570101-570102, 570301-570302, 5706-5710, 5714, 5715, 5717																															
*3558																															
*3566																															
*3571, 4443, 5807																															
*3903, 3914, 3924, 3929, 3937, 4454																															
*390601-390602, 3934, 4086, 4502, 4905, 5903, 7301																															
*3908, 3911																															
*3932, 5145																															
*3933																															
*3936																															
*400106																															
*4012																															
*4026																															
*4052																															
*4059																															
*4077, 4087																															
*4209																															
*4417																															
*4422																															
*4509																															
*460101, 4604, 4607N, 4610, 4614-4616																															
*460102, 4602, 4605																															
*4603																															
*4606																															
*4608																															
*4609																															
*4611																															
*4612																															
*4613																															
*4812, 4814																															
*4819																															
*5106																															
*5406, 5603																															
*5409, 5504, 5610, 5618																															
*5414, 5513, 5523, 5532																															
*5508, 5602, 5604																															
*5521																															
*5622																															
*5702, 5712																															
*5713																															
HLA-B allele ³																															
Well No.																															
Cw*0203																															
Cw*0348																															
Cw*0746																															
Cw*1508																															
Cw*1515																															
Well No.																															