

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	
HLA-B allele ^{1,2,3}																																																
*27:01	1					6		8															23																									
*27:02:01:01-27:02:03, 27:77, 27:95, 27:126		2				6		8															23																									
*27:03, 27:111, 27:151			3	4		6		8															23																								w	
*27:04:01-27:04:06, 27:100, 27:103, 27:105, 27:109, 27:112-27:114, 27:147					5	6		8	9													21	23																									
*27:05:02:01-27:05:03, 27:05:05-27:05:19, 27:05:21-27:05:32, 27:84, 27:87-27:88, 27:90:01-27:90:04, 27:93, 27:96:01-27:96:02, 27:98-27:99, 27:117, 27:121-27:122, 27:124, 27:128-27:129, 27:132-27:133, 27:135, 27:140, 27:143-27:145, 27:148, 27:152, 27:155, 27:158-27:160				4		6		8															23																									
*27:05:02:04Q				4		6		8															23																									
*27:05:04, 27:116				4				8															23																									
*27:05:20				4		6		8																																								
*27:06					5	6			9	10												21	23																									
*27:07:01-27:07:05, 27:130				4		6					11																																					
*27:08					5		7	8															23																									
*27:09				4		6														20			23																									
*27:10				4		6		8	9														23																									
*27:11					5	6					11																																					
*27:12					5		7	8					13										23																									
*27:13				4		6		8				12											23																									
*27:14, 27:81				4		6											17																															
*27:15					5	6		8	9													22																										
*27:16				4		6		8					13										23																									
*27:17				4		6		8					14										23																									
*27:18	1				5		7	8	9					15									23	24																								
*27:19				4		6								16																																		
*27:20					5	6			9		11												23																									
*27:21					5	6																																										
*27:23, 27:131	1				5	6		8															23																									
*27:24					5	6					11							18				21																										
*27:25					5	6		8	w										19	20		21																										
*27:26					5		7	8															22	23																								
*27:27				4		6					11												23																									
*27:28				4		6		8															22																									
*27:29	1			4		6		8						15									23	24																								
*27:30		2				6								16																																	w	
*27:31					5	6		8														22	23																									
*27:32				4		6																																										
*27:33					5		7				11																																					
*27:34				4		6							14																																			
*27:35				4		6					10								19				23																									
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	

Negative Control

[illegible]

[illegible]

[illegible]

[illegible]

[illegible]

[illegible]