

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	
HLA-B allele ^{1,2,3}																																																
*27:01	1					6		8															23																			44						
*27:02:01:01-27:02:05, 27:77, 27:95, 27:126, 27:163, 27:171-27:172, 27:176N, 27:188, 27:197, 27:203		2				6		8															23																						w			
*27:03, 27:111, 27:151, 27:182			3	4		6		8															23																									
*27:04:01-27:04:06, 27:100, 27:103, 27:105, 27:109, 27:112-27:114, 27:147, 27:210-27:211, 27:214, 27:226					5	6		8	9												21	23																										
*27:05:02:01-27:05:02:03, 27:05:02:05-27:05:03, 27:05:05-27:05:19, 27:05:21-27:05:36, 27:05:38-27:05:47, 27:87-27:88, 27:90:01-27:90:04, 27:93, 27:96:01-27:96:02, 27:98-27:99, 27:117, 27:121-27:122, 27:124, 27:128-27:129, 27:132-27:133, 27:135, 27:140, 27:143-27:145, 27:148, 27:152, 27:155, 27:158-27:160, 27:162, 27:166-27:167, 27:170, 27:174-27:175, 27:177-27:179, 27:184-27:186, 27:189-27:191, 27:193-27:195, 27:198-27:201, 27:205-27:207, 27:209, 27:212N, 27:215N-27:217, 27:220-27:221, 27:223N-27:225N				4		6		8															23																									
*27:05:02:04Q				4		6		8															23																			43						
*27:05:04, 27:05:37, 27:116, 27:169, 27:196				4				8															23																									
*27:05:20				4		6		8																																								
*27:06:01:01-27:06:01:02, 27:208, 27:222					5	6			9	10												21	23																									
*27:07:01-27:07:06, 27:130, 27:164, 27:180, 27:187, 27:218-27:219				4		6					11																																					
*27:08, 27:204					5		7	8															23																									
*27:09				4		6														20			23																									
*27:10				4		6		8	9														23																									
*27:11					5	6					11																																					
*27:12:01:01-27:12:01:03					5		7	8					13										23																									
*27:13:01-27:13:02				4		6		8				12											23																									
*27:14, 27:81				4		6											17																															
*27:15					5	6		8	9													21	22																									
*27:16				4		6		8					13										23																									
*27:17				4		6		8					14										23																									
*27:18	1				5		7	8	9					15									23	24																								
*27:19:01:01-27:19:01:02				4		6								16																									39									
*27:20					5	6			9		11											21	23																									
*27:21					5	6																																										
*27:23, 27:131	1				5	6		8															23																									
*27:24					5	6					11							18				21																										
*27:25					5	6		8	w										19	20		21																										
*27:26					5		7	8															22	23																								
*27:27				4		6					11												23																									
*27:28				4		6		8															22																									
*27:29	1			4		6		8						15									23	24																								
*27:30		2				6									16																															w		
*27:31					5	6		8															22	23																								
*27:32				4		6																																										

Negative Control

[illegible]

[illegible]

[illegible]

[illegible]

[illegible]

[illegible]

[illegible]