

INTERPRETATION TABLE		
HLA-B*39 Add-on SSP typing		
	Well	
	1	2
Length of spec.	105	235
PCR product		
Length of int.	800	1070
pos. control ¹		
5'-primer(s) ²	246	-149
	5' -gAA 3'	5' -ggg 3'
	246	
	5' -gAA 3'	
3'-primer(s) ³	311	44
	5' -ggT 3'	5' -CCg 3'
Well No.	1	2
HLA-B allele ⁴		
*39:01:01:03, 39:01:03-39:01:08, 39:01:10-39:01:14, 39:01:16-39:02:01, 39:03-39:09, 39:11-39:12, 39:14-39:15, 39:18, 39:19:02, 39:22, 39:24:01-39:39:01, 39:40:01N-39:48, 39:50-39:62, 39:64-39:82, 07:02:32, 14:01:01-14:02:02, 14:02:04, 14:02:07-14:20, 14:22-14:33, 15:189, 35:26, 39:01:01:01, 67:01:01, 67:03	1	
*39:01:01:02L	1	2
HLA-B allele		
Well No.	1	2