

*08:62, 42:02:01:01-42:02:01:02, 42:09, 42:18	2	3		5				10	12															
*08:71, 42:10	2							10	12					18										
*08:96	2		4					10	12															
*08:107			4					10																
*08:116, 08:134, 08:175, 42:08, 42:14	2								12				16											
*08:165	2							10	12					17										
*13:02:01:01-13:04, 13:08-13:09, 13:14-13:16, 13:18-13:19, 13:27, 13:31-13:35, 13:37-13:38, 13:40-13:42, 13:44-13:45, 13:47-13:49N, 13:54, 13:56:01N-13:56:02N, 13:58-13:59, 13:62, 13:65-13:70, 13:72, 13:74-13:75, 13:81, 13:84-13:85, 13:88-13:89, 13:91, 13:93-13:94, 13:96-13:105, 15:04:01:01-15:04:03, 15:16:01:01-15:16:03, 15:67, 15:95, 15:137, 15:155, 15:222, 15:254, 15:293, 15:303, 15:310, 15:362, 27:14, 27:81, 27:130, 35:37, 35:235, 35:265, 35:273, 40:75, 40:103, 40:162, 46:11, 46:32, 51:01:01:01-51:03, 51:05, 51:07:01-51:07:02, 51:09:01-51:19, 51:22-51:24:05, 51:26-51:35, 51:37-51:41N, 51:43, 51:48-51:55, 51:57-51:58, 51:60-51:61:02, 51:63:01-51:63:02, 51:65-51:80, 51:82-51:95, 51:98N-51:100, 51:102, 51:104-51:115, 51:117-51:134, 51:136-51:138, 51:140, 51:142, 51:144-51:147, 51:150-51:152, 51:154-51:156, 51:158:01-51:193, 51:195-51:196, 51:198-51:202, 51:204-51:212, 51:214-51:219, 51:221-51:226, 51:228-51:229, 52:01:01:01-52:18, 52:20-52:36, 52:38-52:44, 52:46-52:50, 52:52-52:70, 52:72-52:73, 54:03, 54:26, 54:33, 55:09, 55:22, 55:24, 56:01:01:01-56:01:07, 56:01:09, 56:05:01-56:08, 56:14-56:17, 56:19N-56:22, 56:24-56:30, 56:33-56:34, 56:36-56:50, 56:52-56:53, 56:55, 58:08:01-58:08:02, 59:04, 59:09, 73:01-73:02, 78:01:01:01-78:06, 78:08-78:09																								
*13:23, 15:71, 15:175, 15:204, 15:225, 15:344, 15:381, 18:80, 18:98, 18:112, 18:120, 37:22, 37:41, 40:170, 40:185, 40:215, 40:310							8																	
*13:30						7					13			19										
*13:46, 15:83, 51:36, 51:153, 55:56					6	7			11	13														
*13:53, 15:42, 35:60, 35:276, 46:18, 51:21, 51:101, 51:143, 51:157, 51:197, 52:37, 52:71, 54:01:01-54:01:07, 54:04-54:05N, 54:07-54:08N, 54:10-54:25, 54:27-54:32, 54:34-54:37, 55:01:01-55:01:06, 55:01:08-55:03, 55:05, 55:07, 55:10-55:13, 55:15, 55:17-55:19, 55:21, 55:23, 55:25-55:48, 55:50, 55:52-55:55N, 55:57-55:66, 55:68-55:85, 55:87, 56:23, 56:54, 59:01:01:01-59:03, 59:05-59:08, 59:10N, 78:07, C*05:152, C*12:183						7			11	13														
*13:55, 40:83, 40:362, 49:03, 52:51						7	8				13													
*13:82											13													
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
*14:01:01:01-14:04, 14:06:01-14:07N, 14:10, 14:12, 14:31, 14:33-14:34, 14:41N, 14:46, 14:49, 14:59, 15:09:02, 35:01:25, 35:35, 38:01:01:01-38:02:07, 38:05-38:06, 38:09, 38:11-38:12, 38:14, 38:17-38:18, 38:24, 38:28, 38:35, 38:46-38:50, 38:55Q, 38:58, 38:60, 38:67-38:68Q, 39:01:01:01-39:01:01:07, 39:01:03:01-39:02:03, 39:04-39:05:02, 39:07-39:13:02, 39:15, 39:19:01, 39:20, 39:22-39:23, 39:25N, 39:28, 39:36, 39:38Q-39:39:01, 39:40:01N-39:42, 39:44-39:47, 39:49, 39:58-39:60, 39:63, 39:77, 39:85-39:88, 39:92, 39:97N-39:98, 39:103-39:104, 39:116N, 39:119, 39:121-39:123, 39:130-39:131, 58:20, 67:01:01-67:04, 67:06-67:07																						23		
*14:05, 14:53, 39:03:01:01-39:03:01:02, 39:14:01:01-39:14:01:02, 39:24:01-39:24:02, 39:37															16		18						23	
*14:08:01-14:09, 14:11, 14:14-14:30, 14:32, 14:35-14:36, 14:38-14:40, 14:42-14:45, 14:47-14:48, 14:50-14:52, 14:54-14:55, 14:57-14:58, 38:03-38:04, 38:07-38:08, 38:10, 38:13, 38:15-38:16, 38:20-38:21, 38:23, 38:25-38:27, 38:29, 38:31-38:34N, 38:36-38:45, 38:51-38:54, 38:56-38:57, 38:59, 38:61-38:66, 39:16-39:18, 39:19:02, 39:26-39:27, 39:30, 39:32, 39:35, 39:39:02, 39:48, 39:51-39:56, 39:61, 39:65-39:75, 39:78-39:82, 39:84, 39:89, 39:91, 39:93-39:96, 39:99-39:100, 39:105-39:108, 39:110-39:115, 39:117-39:118, 39:124-39:126																						?		
*14:13, 38:19, 39:29, 39:76, 39:120															16		18						?	
*14:37, 39:57, 39:83, 39:109, 39:127						7				11		13											?	pl

Negative Control

*15:04:04, 39:33, 39:50								7						13									23	
*15:46, 15:106, 35:47, 35:154, 40:10:01:01-40:10:02, 44:269, 53:15, 57:45				4																			Negative Control	
*15:51, 15:179:01-15:179:02, 15:199, 15:218Q, C*07:04:03, C*07:199:02													12											
*15:53, 15:336, 15:347, 18:107, 27:119, 40:49, 40:129, 40:191, 40:210, 40:328, 44:02:01:01-44:02:01:06, 44:02:03-44:02:04, 44:02:06-44:02:32, 44:02:34-44:02:47, 44:03:01:01-44:03:37, 44:03:39-44:05:04, 44:07-44:14, 44:16-44:17, 44:19N, 44:21-44:34:01, 44:35-44:40, 44:41:02-44:46, 44:48-44:53:02, 44:56N-44:59:02, 44:61N-44:79:02, 44:81-44:82, 44:85-44:89, 44:91-44:92, 44:94-44:96, 44:98-44:99, 44:101-44:102, 44:104-44:105, 44:107-44:130, 44:132-44:134, 44:136, 44:138Q, 44:140-44:155, 44:157, 44:159-44:165, 44:167-44:183, 44:185-44:187, 44:189-44:193, 44:195N-44:196, 44:198N-44:204, 44:206-44:211, 44:214-44:218, 44:220-44:224, 44:226-44:229, 44:231, 44:233-44:234, 44:236, 44:238, 44:240-44:241, 44:243-44:244, 44:247-44:262, 44:264-44:268, 44:270-44:277, 47:01:01:01-47:10, 57:51, 57:69, C*16:85																					22			
*15:159, 15:260, 35:302, 46:52														15										
*15:210, 18:01:01:01-18:01:12, 18:01:14-18:02, 18:05:01:01-18:06, 18:08-18:13, 18:15, 18:17N-18:24, 18:26-18:28, 18:30-18:32, 18:34-18:42, 18:44:01-18:55, 18:57-18:78, 18:81-18:97, 18:99-18:101, 18:103-18:106, 18:108-18:111, 18:114-18:119, 18:121-18:124, 18:126-18:135, 18:137-18:143, 37:01:01:01-37:08, 37:10-37:11, 37:13-37:20, 37:23-37:25, 37:27-37:36, 37:38-37:40, 37:42N-37:43, 37:45-37:54, 37:58-37:59, 37:61, 37:63-37:66, 37:68-37:70, 48:22, 48:30				3				8																
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
*15:212, 18:102, 27:18, 27:29, 35:63, 35:317, 40:01:01-40:01:38, 40:01:40-40:01:41, 40:01:43-40:01:49, 40:04, 40:11:01-40:11:02, 40:14:01-40:16, 40:20:01:01-40:20:01:02, 40:22N-40:23, 40:25-40:26, 40:28, 40:30-40:34, 40:36, 40:38, 40:42, 40:45-40:48, 40:51-40:52, 40:54-40:55, 40:59-40:60, 40:62-40:67, 40:69, 40:72:01-40:74, 40:76-40:77, 40:80-40:81, 40:84, 40:87:01-40:88, 40:92, 40:98-40:102, 40:108, 40:112-40:114:02, 40:116-40:118N, 40:120-40:121, 40:123-40:126, 40:128, 40:130:01-40:130:02, 40:132, 40:134-40:135, 40:138-40:141, 40:146-40:147, 40:150-40:156, 40:158, 40:160:01-40:160:02, 40:163, 40:168, 40:171-40:172, 40:175, 40:178-40:179, 40:182-40:184, 40:186:01-40:188, 40:192-40:199, 40:204, 40:207-40:208, 40:212-40:213, 40:216N-40:218, 40:221-40:223, 40:225, 40:227-40:228, 40:233-40:243, 40:245, 40:247, 40:249-40:251, 40:253, 40:257-40:265N, 40:272-40:273, 40:277-40:278, 40:280-40:281, 40:285-40:286N, 40:288, 40:298:01-40:301, 40:308, 40:312, 40:315-40:316, 40:319, 40:321, 40:323-40:324, 40:326-40:327, 40:329, 40:332-40:333, 40:336, 40:338N-40:339, 40:343-40:344, 40:346, 40:352-40:353, 40:355, 40:357-40:358, 40:360, 44:55, 44:103, 44:131, 44:188, 44:212, 49:47, 50:12, 50:17, 50:51								8														22		
*15:408, 52:45, 56:01:08								7																
*15:426, 18:07:01-18:07:02, 18:33, 18:79, 18:136, 27:75, 35:162, 37:21, 37:26, 37:57, 37:60, 37:67, C*06:77				3																				
*18:03, 37:37				3	5			8																
*18:14				3				8							16		18							
*18:56, 18:113, 37:12, 37:44, 37:56				3				8				12												
*27:04:01-27:04:06, 27:06, 27:08, 27:12:01:01-27:12:01:02, 27:15, 27:20-27:21, 27:25-27:26, 27:31, 27:36, 27:40, 27:61, 27:63, 27:66N, 27:68-27:69, 27:79, 27:86, 27:91, 27:100, 27:103, 27:105-27:109, 27:112-27:115, 27:120, 27:125, 27:147, 27:149, 27:154, 27:165, 40:166					5																			
*27:11, 27:24, 27:33					5											16								
*27:89, 46:14, C*03:278, C*16:86					w																			
*27:138					5											16		18						
*27:153, 56:35					5		7						13											
*37:09			2	3				8				12				16								

*38:22																				22	?	?		
*38:30											12										23			
*39:06:01-39:06:05, 39:34, 39:62, 39:64, 39:129							7				11		13								23			
*39:31:01:01-39:31:01:02			3		5																23			
*39:90											11		13								?			
*39:101			3		5		8														?			
*39:102							7						13								?			
*39:128			3		5																?			
*40:01:42								8													22	23		
*40:02:01:01-40:02:02, 40:02:05-40:03:01:02, 40:05:01:01-40:05:01:02, 40:09, 40:18-40:19, 40:24, 40:27:01, 40:29, 40:37, 40:39, 40:56-40:57, 40:78, 40:82, 40:85, 40:89-40:91, 40:97, 40:104-40:105, 40:107, 40:111, 40:115, 40:119, 40:122, 40:133Q, 40:142N-40:145, 40:157, 40:169, 40:173-40:174, 40:176, 40:180-40:181, 40:189, 40:201-40:203, 40:205- 40:206, 40:211, 40:214, 40:224, 40:246, 40:254-40:255, 40:266-40:267, 40:271, 40:274, 40:276, 40:283-40:284, 40:287, 40:289-40:292, 40:296-40:297, 40:302-40:305, 40:309, 40:314, 40:320, 40:331, 40:334-40:335, 40:337N, 40:345N, 40:347, 40:356, 40:359								8							16					22				
*40:02:04, 40:94, 40:209								8							16									
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
*40:06:01:01-40:06:02, 40:06:04:01-40:06:14, 40:06:16, 40:44, 40:53, 40:70:01-40:70:02, 40:86, 40:93, 40:95- 40:96, 40:109, 40:127, 40:131, 40:148, 40:159, 40:161, 40:165, 40:167, 40:177, 40:190, 40:230, 40:244, 40:256N, 40:268-40:270, 40:275, 40:294-40:295, 40:306- 40:307, 40:311, 40:313, 40:317, 40:330, 40:340-40:341, 40:348-40:350, 40:354, 40:361N, 40:363, 49:01:01:01- 49:02, 49:04:01-49:09, 49:11-49:13, 49:16-49:23, 49:27- 49:30, 49:32-49:46, 50:01:01:01-50:02, 50:04-50:10, 50:14-50:16, 50:18-50:20, 50:31-50:42, 50:44, 50:46- 50:48, 50:50, 50:53-50:54								7	8				13								22			
*40:06:03								7	8				13									22	23	
*40:06:15, 40:110, 40:318, 40:342, 44:47, 44:197, 49:10, 49:24-49:26, 49:31, 50:13, 50:43, 51:103, 51:203								7					13									22		
*40:13, 40:40, 40:58, 40:164, 40:219, 40:229															16							22		
*40:35:01-40:35:02, 40:71								8							16		18					22		
*40:43					5			8														22		
*40:50					5			8							16							22		
*40:79								8							15							22		
*40:149			3					8														22		
*40:248				4				8							16									