

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
HLA-B allele^{1,2,3}																
*42:01:01:01-42:01:04, 42:22-42:24, 42:28	1	2				6			9							
*42:02:01:01-42:02:02		2	3			6			9							
*42:04	1			4		6			9							
*42:05:01-42:05:02, 07:349 ⁵	1	2			5	6			9							
*42:06	1	2				6	7		9							
*42:07, 42:26, 08:94, 08:226 ⁵		2				6		8								
*42:08	1	2				6			9	10						
*42:09		2	3				7		9							
*42:10	1	2				6			9		11					
*42:11		2				6		8				12				
*42:12	1	2				6			9				13			
*42:13	1					6			9					14		
*42:14, 42:25 ⁴	1	2				6			9						15	
*42:15	1	2				6			9					14		
*42:16, 42:20	1	2					7		9							
*42:17		2	3			6			9		11					
*42:18		2	3			6			9					14		
*42:19, 07:388 ⁵	1								9		11					
*42:21	1	2				6										
*42:27		2	3			6										
*07:02:01:01-07:02:27, 07:02:29-07:02:81, 07:05:01:01-07:07:03, 07:09:01-07:09:02, 07:12, 07:14:01-07:15, 07:17-07:18:02, 07:20-07:24, 07:26, 07:28, 07:30-07:31, 07:34-07:36, 07:39-07:49N, 07:51-07:52, 07:54-07:55, 07:58-07:59, 07:61-07:64, 07:66-07:68:03, 07:70, 07:73-07:74, 07:76-07:82, 07:84, 07:86-07:99, 07:101-07:115, 07:117-07:124, 07:126-07:142, 07:144-07:145, 07:147-07:159, 07:161N-07:176, 07:178-07:179, 07:181N-07:182N, 07:184-07:185, 07:187-07:205, 07:208-07:218, 07:220-07:227, 07:229-01-07:239, 07:241-07:252, 07:254-07:258, 07:260-07:262, 07:264-07:291, 07:293-07:305, 07:307-07:312, 07:314-07:339, 07:341-07:348, 07:350-07:367, 07:369-07:382, 07:384-07:387, 07:389-07:399, 07:401, 14:21, 15:255, 15:509, 38:26, 54:21, 55:10, 56:16, 67:01:01-67:01:03, 67:03-67:07, 81:01:01-81:04N, 81:06-81:09, 83:01	1								9							
*07:03, 07:08:01-07:08:02, 07:16, 07:32, 07:38, 07:69, 07:85:01-07:85:02, 07:180, 07:207, 07:228:01-07:228:02, 07:306, 07:400, 13:14, 13:83, 14:01:01-14:01:05, 14:01:07-14:01:09, 14:07N-14:08:02, 14:10, 14:12, 14:14, 14:19, 14:26, 14:32, 14:40, 14:46-14:47, 14:49, 14:54-14:55, 14:65, 14:71, 14:77-14:78N, 14:80-14:83, 14:86-14:87, 15:71, 15:175, 15:381, 15:556, 18:01:01-18:01:12, 18:01:14-18:03:02, 18:05:01-18:13, 18:15, 18:17N-18:24, 18:26-18:36, 18:38-18:48, 18:50-18:55, 18:57:01-18:67, 18:69, 18:71-18:78, 18:80-18:106Q, 18:109-18:112, 18:114-18:124, 18:126-18:131:01:02, 18:133-18:141, 18:143-18:154N, 18:156:01-18:181, 18:183-18:192, 27:12:01-27:12:01:03, 27:16, 27:23, 27:92, 27:119, 27:153, 35:50, 35:84, 35:162, 35:197, 35:231, 35:267, 35:269, 35:411, 35:423, 37:02, 37:07, 37:25, 37:39-37:40, 37:66, 37:83, 38:01:01-38:02:09, 38:04, 38:06-38:18, 38:20-38:25, 38:27-38:29, 38:31-38:32, 38:35-38:53, 38:55Q-38:83N, 38:85-38:86, 38:156-38:167, 39:01:01-39:01:01:12, 39:01:03-39:01:32, 39:05:01-39:07, 39:09:01-39:11, 39:15-39:19:02, 39:22, 39:25N-39:28, 39:30-39:36, 39:38Q, 39:40:01N-39:48, 39:50-39:51, 39:53-39:62, 39:64-39:71, 39:73-39:75, 39:77-39:87N, 39:89-39:98, 39:100, 39:102-39:104, 39:106-39:109, 39:111-39:119, 39:121-39:123, 39:125-39:128, 39:130-39:143, 39:145-39:151, 39:153, 39:155, 39:157N-39:161N, 39:163, 40:04:01-40:04:02, 40:06:01-40:06:14, 40:06:16-40:06:22, 40:11:01-40:11:02, 40:20:01-40:20:01:02, 40:26, 40:27:03-40:28, 40:44, 40:64:01-40:64:02, 40:68, 40:70:01-40:70:02, 40:74-40:75, 40:83, 40:86, 40:95-40:96, 40:98-40:99, 40:109, 40:120, 40:127, 40:131, 40:148-40:149, 40:159, 40:161-40:162, 40:165, 40:167, 40:177, 40:184, 40:190, 40:210, 40:225:01-40:225:02, 40:230, 40:243-40:244, 40:256N, 40:268-40:270, 40:275, 40:294-40:295, 40:306-40:307, 40:311, 40:313, 40:317, 40:330, 40:340-40:341, 40:343, 40:348-40:351, 40:354, 40:361N-40:364, 40:366, 40:373-40:374, 40:376, 40:387, 40:389, 40:408, 40:411, 40:423, 40:425, 40:427, 40:429, 40:434, 40:445, 40:451, 44:103, 44:131, 44:212, 44:224, 48:06, 48:31, 48:41, 49:31, 50:43, 51:57, 51:103, 51:116, 51:218, 53:21								8								
*07:04:01-07:04:02, 07:25, 07:146	1	2			5		7		9							
*07:11, 07:57, 07:75:01:01-07:75:01:02, 07:160, 07:240, 07:253, 07:340	1								9					14		
*07:19	1	2					7		9	10						
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
*07:27, 07:50, 18:37:01-18:37:02, 27:18, 27:29, 38:03, 39:02:01-39:02:03, 39:08, 39:13:01-39:13:02, 39:23, 39:39:01-39:39:02, 39:49, 39:88, 39:101, 39:105, 39:110, 39:162, 39:164, 40:93, 44:130, 48:01:01-48:03:02, 48:05, 48:07-48:12, 48:14-48:21, 48:23, 48:25-48:30, 48:32-48:40, 48:42-48:50								8				12				

Negative Control

*07:29, 07:71, 07:116, 07:186, 07:313, 07:383, C*07:335										9								
*07:33:01-07:33:03, 07:53, 07:60, 07:100	1									9	10							
*07:56:01-07:56:02, 07:65, 07:72, 07:177, 07:206, 07:263, 15:76, 15:101, 35:76, 44:90, 54:01:01:01-54:01:08, 54:03-54:06, 54:08N-54:20, 54:22-54:23, 54:25-54:29, 54:31-54:41, 55:01:01:01-55:03, 55:05, 55:07-55:09, 55:11-55:15, 55:17, 55:19, 55:21-55:24, 55:26-55:33, 55:35-55:46, 55:48, 55:50, 55:52-55:55N, 55:57-55:63, 55:65-55:77, 55:79-55:85, 55:87-55:107, 56:01:01:01-56:07, 56:09-56:12, 56:15, 56:18-56:22, 56:24-56:32, 56:34, 56:36-56:41, 56:43, 56:45-56:52, 56:55:01:02-56:66, 56:68-56:73, C*14:62	1																	
*07:143			3							9								
*08:01:01:01-08:01:22, 08:01:24-08:03, 08:05, 08:08N, 08:10-08:11, 08:15, 08:18-08:19N, 08:21-08:24, 08:26:01-08:27, 08:29-08:36, 08:38:01-08:39, 08:41-08:46, 08:48, 08:50-08:52, 08:56:01-08:59:02, 08:61-08:68, 08:72N-08:73, 08:75, 08:78, 08:80-08:83, 08:85-08:86N, 08:88, 08:90-08:93, 08:95-08:102, 08:104-08:106, 08:109, 08:112-08:115, 08:117-08:122, 08:124-08:128, 08:130-08:131, 08:133, 08:135-08:142, 08:144, 08:146-08:155, 08:157-08:163, 08:165-08:170, 08:172-08:174, 08:176-08:180, 08:182-08:183, 08:185-08:186, 08:188-08:189, 08:191-08:199, 08:201-08:204, 08:206-08:225, 08:228-08:231, 08:233-08:244, 08:246-08:250, 08:252N-08:253, 35:87, 41:02:01:01-41:02:09, 41:11, 41:13, 41:19, 41:23, 41:27, 41:31, 41:36, 41:38-41:43, 41:45N-41:47, 41:55, 41:58, 44:166	2					6												
*08:01:23, 08:14, 08:28, 08:74, 08:77, 08:87, 08:108, 08:111, 08:164, 08:227, 08:245, 08:251, 41:10, 41:49, 41:52	2																	
*08:04:01-08:04:02, 08:54:01-08:54:02, 08:110	2					6						12						
*08:07, 08:37, 08:200, 41:04, 41:18, 41:51	2						7											
*08:09, 08:84, 08:254, 15:83:01:01-15:83:01:02, 41:01:01:01-41:01:07, 41:05, 41:07, 41:09, 41:12, 41:14, 41:16-41:17, 41:20, 41:22, 41:25-41:26, 41:28-41:30, 41:32-41:35, 41:37, 41:44, 41:53-41:54, 41:56Q-41:57, 41:59-41:61, 41:64, 51:08:01:01-51:08:03, 51:20, 51:36, 51:44N, 51:97, 51:141, 51:153, 51:263, 51:275, 52:19				4		6												
*08:12:01-08:12:03, 08:184, 08:232, 41:03:01-41:03:02, 44:14, 44:62, 57:09						6												
*08:16, 08:60, 08:129, 08:181, 41:48						6										14		
*08:17	2					6						12				14		
*08:20:01-08:20:02, 15:516, 27:83											11							
*08:47, 08:76, 08:145	2					6										14		
*08:49							7				11					14		
*08:53:01-08:53:02, 08:71, 41:24	2					6					11							
*08:69	2					6	7											
*08:89, 08:107, 15:51, 15:179:01-15:179:02, 15:199, 15:218Q, 15:518, 35:38, 35:115, 35:169, 35:260, 37:55, 37:60, 37:79N, 37:85, 41:63, 44:17, 44:51, 44:123, 44:140, 44:230, 53:22, 57:24, 57:119, C*05:62, C*05:146, C*05:163, C*07:626, C*08:142, C*15:25							7											
*08:103, 08:134, 08:175, 44:213	2					6					10							
*08:116, 41:15	2					6											15	
*08:123	2					6				9						14		
*08:132	2					6				9								
*08:187	2				5	6												
*08:190	2					6		w										
*13:46, 41:21, 44:15:01:01-44:15:01:02, 44:18, 45:01:01:01-45:02, 45:04-45:05, 45:07-45:13, 45:15-45:25				4			7											
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16		
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16		
*13:61, 52:92																15		
*13:71, 13:108, 14:05, 14:13, 15:07:01:01-15:07:03, 15:45, 15:55, 15:126, 15:197:02, 15:207, 15:324, 15:329, 15:331, 15:360, 15:383, 15:405, 15:421, 15:431, 15:450, 15:524, 27:07:01-27:07:06, 27:11, 27:24, 27:32-27:34, 27:43, 27:138, 27:180, 27:219, 35:05:01:01-35:05:03, 35:22:01:01-35:22:01:02, 35:31, 35:51, 35:58, 35:72, 35:89, 35:97, 35:114, 35:199, 35:232, 35:368, 35:389, 35:416, 35:427N, 35:440, 35:459N, 35:472, 39:76, 40:27:02, 40:57-40:58, 40:164, 40:200, 40:219, 40:226, 40:420, 41:62, 44:54, 44:135, 44:184, 46:12, 51:64, 51:81, 51:148, 53:14, 58:18, 58:27, C*02:60, C*03:05, C*03:25, C*03:27, C*03:167, C*03:296:01-03:296:02, C*03:335, C*03:386, C*03:407, C*03:482, C*03:494, C*06:264, C*07:294, C*07:526:01, C*07:530, C*07:681, C*07:850, C*08:08:01-08:08:02										10								
*13:97, 14:73, 15:03:01:01-15:03:09, 15:03:11, 15:47:01-15:47:02, 15:54, 15:61-15:62, 15:64:01-15:64:02, 15:69, 15:91, 15:98, 15:103, 15:123, 15:127, 15:131-15:132, 15:151, 15:156, 15:158, 15:173, 15:210, 15:220:01:01-15:220:02, 15:235, 15:242:01-15:242:02, 15:253, 15:259, 15:266, 15:274, 15:281-15:282, 15:369, 15:379, 15:389, 15:397, 15:444, 15:448, 15:475, 15:502, 15:511, 15:517, 15:527, 15:539, 15:560, 18:125, 35:361, 40:12:01-40:12:02, 40:46, 40:430, 44:40, 44:44, 44:156, 44:210:01-44:210:02, 48:22, 52:16, 55:18, 56:44												12						

Negative Control

Negative Control

[illegible]