

[illegible]

*07:78, 07:84, 13:02:01-13:02:07, 13:08Q-13:09, 13:14-13:16, 13:18-13:19, 13:27, 13:30-13:35, 13:37, 27:14, 39:50, 40:93, 40:96, 49:07, 51:10, 51:16, 51:31, 51:34, 51:82, 51:93, 52:08, 55:09, 55:22, 55:24, 73:01						6					11		*07:78, 07:84, 13:02:01-13:02:07, 13:08Q-13:09, 13:14-13:16, 13:18-13:19, 13:27, 13:30-13:35, 13:37, 27:14, 39:50, 40:93, 40:96, 49:07, 51:10, 51:16, 51:31, 51:34, 51:82, 51:93, 52:08, 55:09, 55:22, 55:24, 73:01
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	Well No.
Length of spec.	255	145	160	160	195	255	230	105	245	160	245	205	Length of spec.
PCR product													PCR product
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	Well No.
*08:09, 15:42, 15:83, 35:60, 39:06:01-39:06:02, 39:34, 39:57, 41:05, 42:04, 46:18, 51:21, 51:36, 54:01-54:02, 54:04-54:05N, 54:07-54:08N, 54:10-54:22, 55:01:01-55:03, 55:05, 55:07, 55:10-55:13, 55:15-55:21, 55:23, 55:25-55:40, 56:23, 59:01-59:03, 59:05, 78:07						6			9				*08:09, 15:42, 15:83, 35:60, 39:06:01-39:06:02, 39:34, 39:57, 41:05, 42:04, 46:18, 51:21, 51:36, 54:01-54:02, 54:04-54:05N, 54:07-54:08N, 54:10-54:22, 55:01:01-55:03, 55:05, 55:07, 55:10-55:13, 55:15-55:21, 55:23, 55:25-55:40, 56:23, 59:01-59:03, 59:05, 78:07
*08:28, 08:35, 08:37, 15:51, 15:179, 42:09, 44:17, 53:22, 57:24, C*07:68										10			*08:28, 08:35, 08:37, 15:51, 15:179, 42:09, 44:17, 53:22, 57:24, C*07:68

*13:03-13:04, 15:04, 15:16, 15:67, 15:95, 15:137, 15:155, 35:37, 39:33, 46:11, 49:01:01-49:06, 49:08-49:10, 50:06, 51:01:01-51:03, 51:05, 51:07:01-51:09:02, 51:11N-51:15, 51:17-51:20, 51:22, 51:24:01-51:24:03, 51:26-51:30, 51:32-51:33, 51:35, 51:37-51:41N, 51:43-51:44N, 51:48-51:55, 51:57-51:58, 51:60-51:61, 51:63, 51:65-51:69, 51:71-51:80, 51:83-51:92, 51:94-51:95, 52:01:01-52:07, 52:09-52:14, 52:16-52:20, 54:03, 56:01:01-56:01:04, 56:05:01-56:08, 56:14-56:17, 56:19N-56:22, 56:24-56:29, 58:08, 59:04, 78:01-78:06						6							*13:03-13:04, 15:04, 15:16, 15:67, 15:95, 15:137, 15:155, 35:37, 39:33, 46:11, 49:01:01-49:06, 49:08-49:10, 50:06, 51:01:01-51:03, 51:05, 51:07:01-51:09:02, 51:11N-51:15, 51:17-51:20, 51:22, 51:24:01-51:24:03, 51:26-51:30, 51:32-51:33, 51:35, 51:37-51:41N, 51:43-51:44N, 51:48-51:55, 51:57-51:58, 51:60-51:61, 51:63, 51:65-51:69, 51:71-51:80, 51:83-51:92, 51:94-51:95, 52:01:01-52:07, 52:09-52:14, 52:16-52:20, 54:03, 56:01:01-56:01:04, 56:05:01-56:08, 56:14-56:17, 56:19N-56:22, 56:24-56:29, 58:08, 59:04, 78:01-78:06
*15:46, 15:53, 15:106, 15:143, 18:48, 35:19, 35:47, 35:63, 40:01:01-40:05, 40:07-40:11, 40:14:01-40:16, 40:18, 40:20, 40:22N-40:40, 40:42-40:43, 40:45, 40:48-40:52, 40:54-40:69, 40:71-40:74, 40:77-40:82, 40:84-40:85, 40:87-40:92, 40:94, 40:97-40:99, 41:02:01-41:04, 41:08, 41:10-41:11, 44:09, 44:46, 44:75, 47:02		2											*15:46, 15:53, 15:106, 15:143, 18:48, 35:19, 35:47, 35:63, 40:01:01-40:05, 40:07-40:11, 40:14:01-40:16, 40:18, 40:20, 40:22N-40:40, 40:42-40:43, 40:45, 40:48-40:52, 40:54-40:69, 40:71-40:74, 40:77-40:82, 40:84-40:85, 40:87-40:92, 40:94, 40:97-40:99, 41:02:01-41:04, 41:08, 41:10-41:11, 44:09, 44:46, 44:75, 47:02
*40:06:01:01-40:06:02, 40:44, 40:53, 40:70, 40:75, 40:83, 40:86, 40:95		2				6					11		*40:06:01:01-40:06:02, 40:44, 40:53, 40:70, 40:75, 40:83, 40:86, 40:95
*41:01, 41:06-41:07, 41:09, 41:12		2				6			9				*41:01, 41:06-41:07, 41:09, 41:12
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	Well No.
Length of spec.	255	145	160	160	195	255	230	105	245	160	245	205	Length of spec.

[illegible]

[illegible]