

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
HLA-B allele ^{1,2,3}																				
*48:01:01:01-48:01:07, 48:32-48:33, 48:35-48:38, 48:41-48:48	1	2	3								11									
*48:02:01-48:02:03	1			4							11									
*48:03:01-48:03:02, 48:40	1	2			5						11									
*48:04:01-48:04:02	1		3							10	11									
*48:05	1					6			9		11		13			16				
*48:06		2	3				7											18		
*48:07	1	2	3					8			11									
*48:08	1								9		11					16				
*48:09, 48:11 ⁵	1	2	3								11			14						
*48:10, 39:13:02, 42:11 ⁴	1										11		13							
*48:12	1										11		13		15					
*48:13	1									10	11		13							
*48:14	1										11		w	14	15					
*48:15	1	2	3								11					16				
*48:16		2	3					8			11									
*48:17	1	2			5			8			11									
*48:18, 48:29 ⁵	1	2	3								11				15					
*48:19, 48:27 ⁵	1	2	3						9		11									
*48:20	1	2	3			6					11									
*48:21	1	2			6						11									
*48:22, 48:30	1	2	3								11	12								
*48:23	1	2			5						11	12								
*48:24	1		3						9	10	11									
*48:25	1			4							11								19	
*48:26	1	2			5	6					11									
*48:28	1	2	3													16				
*48:31	1	2	3															18		
*48:34	1	2	3								11						17			
*48:39		2	3								11									
*07:02:01:01-07:02:25, 07:02:27-07:02:38, 07:02:40-07:03, 07:08:01-07:08:02, 07:10, 07:15-07:16, 07:20-07:24, 07:28-07:31, 07:35, 07:37:01-07:37:02, 07:39, 07:41-07:42, 07:44N, 07:46-07:47, 07:49N, 07:55-07:59, 07:61-07:64, 07:67N, 07:70-07:76, 07:79, 07:82-07:83, 07:85:01-07:96:02, 07:98-07:99, 07:101-07:104, 07:106-07:111N, 07:113-07:114, 07:116-07:121, 07:124-07:130, 07:133-07:134, 07:136:01-07:136:02, 07:141-07:142, 07:144-07:145, 07:148-07:154, 07:156-07:161N, 07:166-07:168, 07:170-07:173, 07:175, 07:177, 07:179, 07:184-07:192, 07:194-07:195, 07:199-07:201N, 07:203-07:205, 07:208, 07:211-07:212, 07:215-07:218, 07:221, 07:224, 07:229-07:234, 07:236, 07:238, 07:240-07:247, 07:250-07:254, 07:256-07:257, 07:260, 07:263, 07:265-07:268, 07:272N-07:277, 07:280-07:282, 07:284-07:286, 07:288-07:289, 07:291-07:292, 07:294-07:300, 07:302-07:303:02, 07:305-07:316N, 07:318N-07:323, 07:325N-07:326, 07:328-07:331, 07:334-07:339, 07:341-07:342, 07:344-07:345, 07:347-07:348, 07:353-07:355, 08:156							7		9						16					
*07:02:26, 07:02:39, 07:04:01-07:04:02, 07:26, 07:43, 07:45, 07:52, 07:66, 07:84, 07:131, 07:135N, 07:146-07:147, 07:178, 07:196, 07:198, 07:227, 07:239, 07:350							7									16				
*07:05:01:01-07:06:03, 07:32, 07:69, 07:97, 07:105, 07:112, 07:138, 07:140, 07:176, 07:182N, 07:206, 07:210, 07:213, 07:249, 07:258, 07:264, 07:269-07:270, 07:278, 07:283, 07:293, 07:304, 07:317, 07:324, 07:333, 07:340, 07:352, 07:356-07:358, 08:79							7		9				13		16					
*07:07:01-07:07:02, 07:12, 07:14, 07:18:01-07:18:02, 07:77, 07:163, 07:193, 07:248					5		7								16					
*07:09:01-07:09:02, 07:11							7		9				14		16			w		
*07:13, 07:54, 07:65, 07:132, 07:181N, 07:259, 07:271, 07:290, 07:351N									9						16					
*07:17							7		9				14		16					
*07:19, 07:25, 07:48, 07:164, 07:223, 07:226, 08:01:12, 08:01:23, 08:07, 08:12:02, 08:14, 08:64, 08:74, 08:83, 08:87, 08:89, 08:107, 08:111, 08:200, 08:211, 14:01:01:01-14:03, 14:06:01-14:09, 14:11-14:12, 14:14-14:48, 14:50-14:52, 14:54-14:56, 14:58-14:68, 38:02:06, 38:29, 38:72, 39:01:19, 39:07, 39:15, 39:17, 39:25N, 39:43-39:44, 39:55, 39:91, 39:96, 39:141, 42:16, 42:20, 45:22, 50:15							7													
*07:27, 07:50							7		9		11				16					
*07:33:01, 07:60, 07:139:01							7		9	10										
*07:33:02-07:33:03, 07:51, 07:122, 07:139:02, 07:155, 07:165, 07:220, 07:261, 07:301, 07:343N, 07:346							7		9											
*07:34, 07:80, 07:225, 07:349							7						13		16					
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
*07:36, 07:38, 07:68:01-07:68:03, 07:81, 07:180							7		9						15	16				
*07:40, 07:332, 08:20:01-08:20:02, 08:40, 08:70, 08:77							7		9				13							
*07:53							7		9	10			13							
*07:78							7	8					13		16					

Negative Control

*07:100								7		9	10				15					
*07:115								7	8	9					16					
*07:123										9					13		16			
*07:137						5		7							13		16			
*07:143								7		9			12				16			
*07:162, 07:255						5		7								16			w	
*07:169								7		9						16	17			
*07:174, 07:222, 07:279	2	3						7		9						16				
*07:183									w	9						16				
*07:197						5	6	7								16				
*07:202	2	3						7		9										
*07:207								7		9					13		16		18	
*07:209, 08:01:01:01-08:01:11, 08:01:13-08:01:22, 08:01:24-08:02, 08:04:01-08:05, 08:08N-08:09, 08:11-08:12:01, 08:12:03-08:13, 08:15-08:16, 08:18-08:19N, 08:21-08:22, 08:24-08:27, 08:29-08:34, 08:36, 08:39, 08:41-08:48, 08:50-08:51, 08:53:01-08:53:02, 08:55-08:63, 08:65-08:69, 08:72N-08:73, 08:75-08:76, 08:80-08:82N, 08:85-08:86N, 08:88, 08:90-08:100, 08:102-08:106, 08:108-08:110, 08:112-08:115, 08:117-08:132, 08:135-08:146, 08:149-08:155, 08:157-08:166, 08:168-08:174, 08:176-08:199, 08:201-08:210, 08:212, 08:214N, 08:216-08:224, 38:02:01-38:02:05, 38:02:07-38:02:08, 38:08, 38:15, 38:17-38:18, 38:23, 38:35, 38:43-38:48, 38:50, 38:62, 38:64, 38:74-38:76, 38:79, 38:82, 39:01:01:01-39:01:01:09, 39:01:03:01-39:01:18, 39:01:20-39:01:26, 39:04-39:06:06, 39:09:01:01-39:10:01, 39:11-39:12:01:02, 39:16, 39:18-39:20, 39:26-39:28, 39:30-39:36, 39:38Q-39:41, 39:45-39:48, 39:50-39:54, 39:57-39:75, 39:77, 39:79-39:90, 39:92-39:95N, 39:97N, 39:99-39:100, 39:102-39:104, 39:106-39:108, 39:111-39:119, 39:121-39:140, 39:142N-39:143, 42:01:01-42:01:03, 42:04-42:07, 42:12-42:13, 42:15, 42:19, 42:21-42:25, 67:01:01-67:01:03, 67:03-67:07								7						13						
*07:214						5		7	8						13		16			
*07:219										9						15	16			
*07:228:01-07:228:02								7		9							16		18	
*07:235, 08:49, 14:10				4				7												
*07:237								7		9					14	15	16		w	
*07:262								7		9					13		15	16		
*07:287								7		9					13	14		16		
*07:327								7	8	9					13		16			
*08:03, 08:52, 08:71, 08:78, 38:01:01:01-38:01:01:02, 38:01:03-38:01:15, 38:05-38:07, 38:09-38:14, 38:16, 38:21, 38:24, 38:26-38:28, 38:30-38:34N, 38:36-38:37, 38:39-38:42, 38:51, 38:53-38:55Q, 38:57-38:61, 38:65-38:69, 38:71, 38:73, 38:77-38:78, 38:80N-38:81, 42:10								7							13		15			
*08:10, 08:147-08:148N, 08:215N, 15:09:01-15:09:02, 15:229, 15:243, 15:252, 15:303, 15:361, 15:435, 15:486, 27:20, 27:27, 27:130, 27:153, 27:164, 27:168, 27:187, 39:10:02, 39:78, 40:04:01-40:04:02, 40:06:01:01-40:06:09, 40:06:11-40:06:15, 40:06:17-40:06:18, 40:11:01-40:11:02, 40:26, 40:64:01:01-40:64:01:02, 40:70:01-40:70:02, 40:99, 40:244, 40:268-40:270, 40:373, 40:387, 40:389, 40:408, 40:411, 44:166, 46:11, 46:18, 46:61, 46:74, 51:01:01:01-51:01:18, 51:01:20-51:01:30, 51:01:32-51:01:49, 51:01:51-51:02:02, 51:02:04-51:03, 51:05-51:06:04, 51:08:01:01-51:16, 51:19-51:24:05, 51:26, 51:28-51:36, 51:38-51:41N, 51:43-51:44N, 51:48-51:55, 51:57-51:63:02, 51:65-51:74, 51:76-51:80, 51:82-51:83, 51:86-51:89, 51:91, 51:95-51:115, 51:117-51:121, 51:123-51:129, 51:131-51:134, 51:136-51:138, 51:140-51:145, 51:147, 51:149N-51:150, 51:152-51:162, 51:164-51:171, 51:173Q-51:176, 51:178N-51:184N, 51:186-51:194, 51:196, 51:198-51:205, 51:207-51:213, 51:215, 51:217-51:219, 51:221, 51:224-51:226, 51:228-51:246, 51:248-51:250, 51:253, 51:255-51:260, 51:262-51:267, 52:01:13, 52:06:01-52:06:03, 52:22-52:23, 52:70, 54:01:01-54:05N, 54:07-54:08N, 54:10, 54:12-54:13, 54:17-54:21, 54:23-54:26, 54:28-54:38, 55:01:01:01-55:03, 55:05, 55:07-55:13, 55:15-55:22, 55:24-55:26, 55:28-55:40, 55:42-55:45, 55:47-55:48, 55:52-55:57, 55:59-55:67, 55:69-55:85, 55:88-55:97N, 56:01:01:01-56:02:01, 56:04:01-56:08, 56:10, 56:13-56:17, 56:19N-56:30, 56:33-56:49, 56:51-56:59, 56:62-56:64, 57:82, 58:08:01-58:08:02, 58:20, 59:01:01:01-59:10N, 67:02, 78:01:01:01-78:01:02, 78:03, 78:07-78:08, 82:01:01:01-82:03, C*01:90, C*01:116, C*02:68, C*02:125, C*03:317, C*06:21, C*14:92, C*15:24, C*15:91, C*15:128														13						
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
*08:17, 08:38:01-08:38:02, 08:54:01-08:54:02, 08:101, 38:03-38:04, 39:02:01-39:02:03, 39:08, 39:13:01, 39:22, 39:101, 39:105, 39:110, 40:317, 41:01:01:01-41:03:02, 41:05-41:07, 41:09, 41:11-41:14, 41:16-41:17, 41:19-41:20, 41:22-41:23, 41:26-41:36, 41:38-41:45N, 41:47-41:48, 41:50, 41:54-41:56, 45:12, 45:14, 50:36	1														13					
*08:23, 08:133							7								13				18	
*08:28, 08:37							7								14					
*08:35							7								13	14				
*08:84, 08:167, 38:49, 39:56, 39:98, 39:109							7							w						

Negative Control

Negative Control

[illegible]

*15:07:01:01-15:07:03, 15:126, 15:207, 15:331, 15:405, 15:431, 15:450										9	10	11		14	15					
*15:08:01-15:08:02, 15:56, 15:191, 15:479, 27:23, 35:01:55, 35:32:01-35:32:02, 35:37, 35:43:01, 35:43:03, 35:53N, 35:64:01-35:64:02, 35:67-35:68:02, 35:75, 35:79, 35:86, 35:102, 35:117-35:119, 35:127, 35:151:01-35:151:02, 35:174, 35:185, 35:205, 35:230, 35:234-35:235, 35:260, 35:265, 35:331, 35:369, 35:409, C*06:147																		18		
*15:21:01:01-15:21:01:02, 15:44, 15:139, 15:394, 15:404, 15:420, 15:468, 27:19, 27:30, 27:127, 27:150, 35:03:12, 35:03:18, 35:26, 35:29:01-35:29:02, 35:63, 35:69, 35:76, 35:82, 35:85, 35:112, 35:136, 35:153, 35:186, 35:198, 35:203, 35:212, 35:229, 35:243, 35:249, 35:254, 35:294, 35:305, 35:317, 35:355, 37:01:07, 44:02:01:01-44:02:33, 44:02:35-44:02:57, 44:03:01:01-44:09, 44:11-44:13, 44:16-44:17, 44:19N, 44:21-44:30, 44:32-44:40, 44:42-44:46, 44:48-44:52N, 44:55-44:59:02, 44:61N, 44:63-44:64:02, 44:66-44:76, 44:78-44:81, 44:84:02-44:98, 44:101-44:105, 44:107-44:109, 44:111-44:134, 44:136-44:137, 44:139-44:149N, 44:151-44:157, 44:159-44:165, 44:167-44:181, 44:183, 44:185-44:191, 44:195N-44:196, 44:198N-44:204, 44:206-44:210:02, 44:212, 44:214-44:224, 44:226-44:229, 44:231-44:235, 44:237N-44:245, 44:247-44:268, 44:270:01-44:304, 44:306N-44:333N, 46:33, 46:67, 51:42, 51:227, 53:01:01:01-53:04, 53:07, 53:10-53:13, 53:15, 53:18, 53:20-53:27, 53:29-53:32, 53:34-53:35, 53:37, 53:39, 53:42-53:43, 53:45-53:46, 53:48N-53:49, 53:51-53:53, 56:09, 56:11, 56:60, 57:01:01:01-57:01:04, 57:01:06-57:01:36, 57:04:01-57:06, 57:08, 57:10, 57:14:01-57:16, 57:18-57:20, 57:23-57:27, 57:29-57:30, 57:32-57:38, 57:40-57:41, 57:43-57:45, 57:48-57:56, 57:58, 57:60-57:62, 57:64-57:65, 57:67:01-57:69, 57:71-57:72, 57:75-57:79N, 57:81, 57:85-57:93, 57:95, 57:97-57:100, 57:102-57:106, 57:108-57:114, 58:01:01:01-58:01:02, 58:01:04-58:01:15, 58:01:17-58:01:31, 58:04, 58:09-58:15, 58:17N, 58:21-58:24, 58:29, 58:31N-58:37, 58:39N-58:40, 58:42, 58:44-58:45:02, 58:47-58:59:02, 58:62-58:64, 58:66-58:72N, 58:74-58:90, 58:92-58:100, 83:01, C*03:102, C*03:263:01-03:263:03																			Negative Control	
*15:29																		18		19
*15:45										9	10	11			15					
*15:47:01-15:47:02, 15:49, 15:221, 15:379, 46:03, 46:62, 57:13, C*01:102, C*15:143																				w
*15:55, 15:360, 15:383, 35:72											10	11			14	15				
*15:62, 15:80, 53:06, 53:08:01-53:08:02, 55:14, 57:22, 58:19, 58:91, C*15:126				4															19	
*15:65												11							w	
*15:68	1									9	10	11			14	15				
*15:124, 18:01:01:01-18:01:12, 18:01:14-18:03:02, 18:05:01:01-18:06, 18:08-18:11, 18:13, 18:17N-18:20, 18:23N-18:24, 18:27-18:28, 18:31-18:34, 18:36-18:42, 18:44:01-18:56, 18:58-18:60, 18:62-18:68, 18:70-18:71, 18:73-18:78, 18:80-18:96, 18:98-18:100, 18:103-18:104, 18:106, 18:108-18:109, 18:111-18:124, 18:126-18:127, 18:129-18:148, 18:151, 18:153-18:161, 27:75, 46:14, 46:43, C*07:123, C*07:173													12							
*15:138					5							11				16			19	
*15:194, 15:223, 15:402, 53:28				4								11							19	
*15:197:02										10					14	15			19	
*15:230, C*02:23, C*03:260					5											16				
*15:241					5							11				16				
*15:297				4	5							11				16				
*15:324										9	10	11		13	15					
*15:329, 27:180, 44:184, 44:213, 51:64										10				13	15					
*15:421										10	11			14						
*18:07:01-18:07:02, 18:29, 18:79													12					18		
*18:12:01-18:12:02													12			16				
*18:14										10	12			14	15					
*18:15, 18:57:01-18:57:02, 18:101, 18:152													12					19		
*18:21, 18:30													12					w		
*18:22				4							10	12								
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
*18:26, 18:61, 18:97, 51:116, 51:151												12	13							
*18:35, 18:110					5							12								
*18:69, 18:105, 18:149-18:150, 37:21, 37:26, 37:57, 37:60, 37:79N, 57:59, 57:73, 58:05, 58:41				4								12								
*27:07:01-27:07:06, 27:11, 27:33, 27:43, 40:02:01:01-40:02:03, 40:02:05-40:02:09, 40:02:11-40:02:17, 40:02:19-40:02:26, 40:05:01:01-40:05:01:02, 40:13, 40:27:01, 40:37, 40:39-40:40, 40:50, 40:90, 40:94, 40:97, 40:111, 40:119, 40:122, 40:142N-40:144N, 40:214, 40:229, 40:254-40:255, 40:266, 40:271, 40:296, 40:302-40:304, 40:309, 40:314, 40:320, 40:331, 40:334, 40:345N, 40:347, 40:356, 40:369, 40:375, 40:378, 40:384, 40:390, 40:399N, 40:409-40:410, 40:412-40:413, 51:81, 51:148, 55:04, 55:49, 55:51, 58:27											10			13						
*27:24				3							10									
*27:32, 27:138, 44:54, 44:106, 44:135, 44:158										10					15					

*27:34, 40:02:10, 40:02:18, 40:89, 40:414, C*03:05, C*03:25, C*03:27, C*03:143, C*03:167, C*03:296:01-03:296:02, C*03:335, C*03:346, C*03:386, C*03:407, C*08:08:01-08:08:02, C*15:140									10											
*27:36	2				5	6														
*27:125, 40:10:01:01-40:10:01:02, 40:12, 40:69, 40:353, 40:357, 44:31	2				5															
*35:01:01:01-35:01:54, 35:03:01:01-35:03:11, 35:03:13-35:03:17, 35:03:19-35:03:25, 35:07-35:08:10, 35:10, 35:13-35:15:02, 35:19, 35:23-35:25, 35:27, 35:33-35:36, 35:38-35:42:02, 35:45-35:49, 35:52, 35:54-35:57, 35:61:01-35:62, 35:65Q, 35:70-35:71, 35:74, 35:77-35:78, 35:80, 35:90-35:94, 35:96, 35:98, 35:100-35:101:02, 35:103-35:111, 35:115-35:116, 35:120-35:124, 35:126, 35:128, 35:130N-35:134N, 35:137-35:145N, 35:147-35:148, 35:150:01-35:150:02, 35:152, 35:155-35:156, 35:158-35:161, 35:163, 35:165N-35:171, 35:173:01N-35:173:02N, 35:175-35:181, 35:187, 35:189-35:197, 35:200, 35:202, 35:204, 35:206-35:210, 35:213-35:219, 35:221-35:225, 35:227-35:228, 35:236-35:241, 35:244-35:248, 35:250, 35:252-35:253, 35:255-35:257, 35:259, 35:261-35:264, 35:267-35:269, 35:271-35:272, 35:274-35:275, 35:278-35:279, 35:281-35:284, 35:286-35:291, 35:293, 35:295, 35:297-35:304, 35:306-35:308, 35:310, 35:312-35:315, 35:319-35:320, 35:322, 35:324-35:330, 35:332-35:334, 35:336, 35:338, 35:340-35:345, 35:347-35:354, 35:356, 35:358-35:360, 35:362-35:365, 35:370-35:371, 35:373, 35:375-35:376, 35:380-35:383, 35:385-35:386, 35:388, 35:390N, 35:392-35:393, 35:395, 35:397-35:402, 35:404-35:406, 35:408, 35:411-35:412, 53:05:01-53:05:02, 53:16, 53:33	4															18				
*35:02:01:01-35:02:12, 35:04:01-35:04:03, 35:06, 35:12:01-35:12:04, 35:59:01-35:59:02, 35:81, 35:83, 35:129N, 35:146, 35:149, 35:154, 35:157, 35:172, 35:182-35:184, 35:201, 35:211, 35:220, 35:233, 35:258, 35:266, 35:270, 35:285, 35:309, 35:311, 35:335, 35:339, 35:357, 35:366, 35:374, 35:377-35:379, 35:384, 35:387, 35:391, 35:394, 35:407, 35:413				4									13				18			
*35:05:01:01-35:05:03, 35:51, 35:89, 35:97, 35:368, 35:389									10				14	15			18			
*35:09:01-35:09:03, 35:18, 35:44, 35:60, 35:87-35:88, 35:251, 35:273, 35:276, 35:292, 35:321, 35:367, 35:372, 51:185, 51:247, 78:02:01-78:02:02, 78:04, 78:09-78:10												13					18			
*35:11:01-35:11:02, 35:21, 35:188, 35:277				4													18 19			
*35:11:03, 35:337:01-35:337:02				4													18 w			
*35:16-35:17:02, 35:30, 35:113, 35:125, 35:164, 35:242, 35:318, 35:346				4					10								18			
*35:22:01:01-35:22:01:02									10			13		15			18			
*35:31, 40:08									10			13					18			
*35:50, 35:84, 35:231, 35:280, 37:67				4								12					18			
*35:58									10				14	15			18 19			
*35:66									9				14	15			18			
*35:99																	18 w			
*35:114									10				14				18			
*35:162, 35:323				4								12	13				18			
*35:199									10					15			18			
*35:232									9	10			14	15			18 w			
*35:316, 35:403, 44:10, 44:14, 44:62, 44:77, 44:82, 51:04:01-51:04:02, 51:46, 51:56:01-51:56:03, 51:139, 53:19, 53:36, 53:50, 56:12, 57:02:01-57:03:04, 57:07, 57:09, 57:12, 57:17, 57:28N, 57:46, 57:57, 57:63, 57:66, 57:70, 57:80, 57:84, 57:94, 57:96, 57:101, 57:107, 58:28:01-58:28:02, 58:65, 58:73				4									13							
*35:396, 35:410				4								w					18			
*37:07									9		11	12			16					
*37:09									10	11	12									
*37:12, 40:362, 52:51											11	12	13							
*38:01:02, 38:22, 38:38, 38:52, 38:70, 49:26						7								15						
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
*38:20:01-38:20:02				4			7						13	15						
*38:25, 41:46	1												13	15						
*38:56													13	15						
*38:63		2			5		7	8						15						
*39:23, 39:49, 40:51, 40:130:02, 40:398, 41:04, 41:10, 41:18, 41:21, 41:37, 41:52, 44:15:01:02, 45:01:01:01-45:05, 45:07-45:11, 45:13, 45:15-45:21, 45:23-45:24, 49:02, 50:01:01:01-50:02, 50:04-50:14, 50:16, 50:18-50:20, 50:31-50:35, 50:37-50:61	1																			
*39:42				4			7						13							
*40:01:01-40:01:17, 40:01:19-40:01:29, 40:01:31-40:01:45, 40:01:47-40:01:60, 40:42-40:43, 40:54, 40:60, 40:66, 40:73, 40:79, 40:113-40:114:01, 40:124:01, 40:125:02, 40:138-40:140, 40:146-40:147, 40:150, 40:152-40:153, 40:155:01N-40:155:02N, 40:170-40:172, 40:186:01, 40:245, 40:249-40:253, 40:258-40:265N, 40:272-40:273, 40:277-40:278, 40:285-40:286N, 40:301, 40:310, 40:321, 40:323, 40:332, 40:346, 40:367, 40:370, 40:372N, 40:379-40:383, 40:385-40:386, 40:391-40:393, 40:395, 40:397, 40:401-40:402, 40:404, 40:406-40:407	1	2			5															

*40:01:18, 40:30, 40:316, 40:403	1	2																	
*40:01:30, 40:01:46	1				5														
*40:02:04										10	11		13						
*40:03:01-01-40:03:01:02, 40:267, 40:284										10				14					
*40:06:16, 40:53, 40:75, 40:86, 40:93, 40:95-40:96, 40:103, 40:110, 40:127, 40:129, 40:148, 40:161, 40:165, 40:167, 40:177, 40:184, 40:190, 40:225, 40:230, 40:243, 40:256N, 40:275, 40:294-40:295, 40:306, 40:311, 40:318, 40:342-40:343, 40:348-40:351, 40:354, 40:361N, 40:363, 40:366, 40:374, 40:376, 40:388	?												13						
*40:07	1	2			5													18	
*40:09, 40:18, 40:24, 40:85, 40:91, 40:180	?									10									
*40:10:02		2			5						11								
*40:14:01-40:14:02, 40:121, 40:136	?				5								13						
*40:14:03, 40:237, 41:08	1				5								13						
*40:15	?									9			13				16		
*40:16	1									9			13				16		
*40:19, 40:292	?									10			13		?				
*40:21, 40:137, 40:365	?	2			5						11								
*40:22N, 40:33, 40:38, 40:46-40:47, 40:49, 40:55, 40:61, 40:63, 40:67, 40:74, 40:76, 40:81, 40:84, 40:87:01-40:88, 40:100-40:102, 40:112, 40:114:02, 40:116, 40:118N, 40:123, 40:124:02-40:125:01, 40:126, 40:128, 40:132, 40:134-40:135, 40:141, 40:149, 40:151, 40:154, 40:156, 40:160:01-40:160:02, 40:168, 40:175, 40:179, 40:182, 40:186:02-40:188, 40:191-40:197, 40:199, 40:204, 40:207-40:208, 40:210, 40:212-40:213, 40:215-40:218, 40:221-40:223, 40:227-40:228, 40:231, 40:234-40:236, 40:239-40:242, 40:247, 40:257, 40:279, 40:281-40:282, 40:288, 40:299-40:300, 40:308, 40:315, 40:319, 40:324-40:325, 40:327-40:329, 40:333, 40:338N-40:339, 40:344, 40:358, 40:360, 40:396, 40:400	?	2			5														
*40:23, 40:183, 40:355	1	2			5												16		
*40:25, 40:166		2			5		?												
*40:27:02, 40:29, 40:56-40:58, 40:78, 40:82, 40:104, 40:107, 40:115, 40:133Q, 40:145, 40:157, 40:164, 40:169, 40:173-40:174, 40:176, 40:181, 40:189, 40:200-40:203, 40:205-40:206, 40:211, 40:219, 40:224, 40:226, 40:246, 40:274, 40:276, 40:283, 40:287, 40:289, 40:291N, 40:293, 40:297, 40:322, 40:335, 40:337N, 40:359	?									10			13						
*40:28	?			4									13						
*40:31	1	2	3																
*40:32	1									9							16		
*40:34, 40:233, 40:298:01-40:298:02	?	2																	
*40:35:01-40:35:02, 40:71	?									10			13		15				
*40:36, 40:280	?	2			5			8											
*40:44, 40:92, 40:120, 40:130:01, 40:131, 40:159, 40:163, 40:178, 40:307, 40:330, 40:336, 40:341	?																		
*40:45, 40:80	?	2	3																

*42:17								7				12	13	15					Negative	
*44:18, 49:01:01:01-49:01:13, 49:05-49:14, 49:16-49:25, 49:27-49:59	1													15						
*44:110, 44:211				4							10									
*44:150				4	5											16				
*46:06		2			5							12								
*46:12										9	10			14	15					
*49:03	1											11	12		15					
*49:04:01-49:04:02	1														15			19		
*50:17	1			4																
*51:01:19, 51:75, 51:206, 51:252, 56:02:02, A*26:01:19													w							
*53:14, 58:18											10			14	15					
*57:39				4									12	13						
*57:42				4										w						
*73:01-73:02							w						12							
*81:01:01-81:01:03, 81:03-81:08		2	3																	
*81:02		2	3				7													
C*02:60										9	10			14						
C*04:77					5			8								16				
C*07:38:01-07:38:02					5															
C*07:294													12	14						
C*07:526:01, C*07:530, C*07:681														14						
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20