

INTERPRETATION TABLE		
HLA-B*5111N SSP typing		
	Well	
	1	2
Length of spec.	95	495
PCR product		
Length of int.	800	1070
pos. control ¹		
5'-primer(s) ²	206	3rd I
	5'-gAC ^{3'}	5'-CTT ^{3'}
3'-primer(s) ³	259	620
	5'-gTT ^{3'}	5'-ggg ^{3'}
Well No.	1	2
HLA-B allele		
B*5111N	1	2
*0832, 180101-1811, 1813-1815, 1817N-1829, 350101-350902, 3511, 3512, 351401-3515, 3517, 3518, 352001-3524, 3527, 3529-3545, 3548, 3550-3562, 3564-356802, 3570-3572, 3574-3579, 3581-3588, 3708, 3806, 3807, 391901-391902, 4406, 510101-5110, 5112-512402, 5126-5146, 5148-5155, 530101-5316, 7801-7804	1	
HLA-B allele		
Well No.	1	2