

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21
HLA-B allele ^{1,2,3}																					
*52:01:01:01-52:01:28, 52:22-52:24, 52:33-52:35, 52:37-52:42, 52:44-52:48, 52:50-52:52, 52:54-52:56	1		3									12	13								
*52:02:01-52:02:02	1			4								12	13								
*52:03			3		5	6						12	13								
*52:04, 52:09 ⁵	1	2	3									12	13								
*52:05	1		3		5	6						12	13								
*52:06:01	1		3						9		11		13								
*52:06:02-52:06:03	1		3								11		13								
*52:07	1		3									12		14							
*52:08	1		3							10		12	13								
*52:10:01-52:10:03			3		6							12	13								
*52:11			3		5							12	13								
*52:12, 52:26	1		3									12	13		15						
*52:13	1		3			7						12	13								
*52:14, 52:49N ⁵	1		3						9			12	13								
*52:15	1		3		6							12	13								
*52:16	1											12	13							20	
*52:17	1		3								11	12	13								
*52:18	1									10			13								
*52:19			3			7						12	13								
*52:20, 52:25 ⁵	1		3				8					12	13								
*52:21	1		3		5							12	13								
*52:27	1		3									12	13							20	
*52:28-52:29, 53:17:01-53:17:02 ⁴			3									12	13								
*52:30	1		3									12	13					18			
*52:31:01-52:31:02, 53:28 ⁵	1		3									12	13		16						
*52:32	1											12	13								
*52:36	1		3	4								12	13								
*52:43	1		3									12	13				17				
*52:53	1		3									12	13						19		
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21
*07:02:01-07:02:10, 07:02:12-07:05:02, 07:05:04-07:26, 07:28-07:37, 07:39-07:49N, 07:51-07:57, 07:59-07:77, 07:79-07:83, 07:85-07:217, 07:219-07:256, 08:01:01-08:01:31, 08:05, 08:07-08:08N, 08:10-08:16, 08:18-08:51, 08:53:01-08:53:02, 08:55-08:77, 08:79-08:83, 08:85-08:109, 08:111-08:116, 08:118-08:125, 08:127-08:156, 13:07N, 13:36, 13:39, 14:41N, 14:43, 14:46, 15:01:01:01-15:01:01:04, 15:01:03-15:01:04, 15:01:06-15:02:09, 15:05:01-15:08:02, 15:10:01-15:12, 15:14-15:15, 15:18:01-15:19, 15:21, 15:23, 15:25:01-15:35, 15:39:01-15:40, 15:42, 15:44-15:46, 15:48-15:53, 15:55-15:58, 15:60, 15:63, 15:65-15:66, 15:70-15:82, 15:84-15:86, 15:88-15:90, 15:92-15:94N, 15:96-15:97, 15:99, 15:101, 15:104-15:122, 15:124-15:126, 15:128-15:129, 15:133-15:136, 15:138-15:150, 15:152-15:154, 15:159-15:161, 15:163-15:167, 15:169-15:172, 15:174-15:176, 15:178-15:195, 15:197-15:201, 15:203-15:207, 15:209N, 15:211-15:215, 15:217-15:219, 15:221, 15:223-15:227, 15:229, 15:231-15:234, 15:236-15:238, 15:240-15:241, 15:244-15:250, 15:255, 15:257-15:265, 15:267,																					

Negative Control

[illegible]

[illegible]

*51:10										10	11	13									
*51:15, 53:11-53:13, 53:31, 53:36								8			11	13									
*51:16, 51:82	1									10	11	13									
*51:31						6				10	11	13									
*51:34						5				10	11	13									
*51:36						5		7			11	13									
*51:37, 51:90	1										11	13		15							
*51:40, 51:108						5	6				11	13									
*51:54						5	w		8		11	13									
*51:61:01-51:61:02	1					5					11	13									
*51:62, 51:78:01-51:78:02	1							8			11	13									
*51:65	1										11	13						17			
*51:69	1										11	13							18		
*51:92:01-51:92:02	1								9		11	13		15							
*51:93						5	6			10	11	13									
*51:104	1			4							11	13									
*51:139							6				11	13						16			
*51:157								8			11			14							
*51:164	1										11										
*51:185	1	2									11	13									
*53:05, 53:16, 53:33		2									11	13									
*53:38			3					8				12	13								
*57:01:11				4																	
*57:60		2		4											14						
*58:36				4										13							
*59:03									9		11			14							
*59:09										10	11			14							
*73:01-73:02									9	10											
*78:05-78:06	1		3											13							
*78:08						6				10				13	14						
*83:01														?	?						
C*01:30, C*08:51, C*08:114, C*12:87, C*15:39	1																				
C*02:02:23, C*12:03:31						5															
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21