

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21
HLA-B allele ^{1,2,3}																					
*54:01:01, 54:01:03-54:01:04, 54:25-54:30	1	2	3																		
*54:01:02			2	3																	
*54:02	1	2		4																	
*54:03	1	2	3		5																
*54:04	1	3				6															
*54:05N	1	2					7														
*54:06	1	3						8													
*54:07	1	2	3											14							
*54:08N	1	2	3		5	6															
*54:09	1	3						8	9												
*54:10	1	2	3							10											
*54:11	1	3			6						11										
*54:12	1	2	3									12									
*54:13	1	2	3										13								
*54:14	1	3						8							15						
*54:15	1	3			6										15						
*54:16	1	2	3												15						
*54:17-54:18 ⁴	1	2	3													16					
*54:19	1	2	3	4																	
*54:20	1	2	3							10								17			
*54:21	1	2	3																18		
*54:22	1	2	3								11										
*54:23	1	2	3																	19	
*54:24	1	2	3																		20
*07:02:19, 37:01:04	1																				

Negative Control

[illegible]