

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21
*13:65											11										
*13:86, 35:02:01:01-35:02:09, 35:04:01-35:04:03, 35:34, 35:81, 35:83, 35:95, 35:129N, 35:154, 35:162, 35:172, 35:182-35:184, 35:211, 35:220, 35:258, 35:270, 35:285, 35:309, 35:311, 35:316, 35:323, 35:335, 35:339, 35:357, 35:366, 40:28, 53:36, 53:42	1									10							17				
*13:97	1								9		11										
*14:54, 15:66, 15:189, 35:135, 38:28, 39:135, 40:95, 51:01:15, 51:14, 51:63:01, 51:85, 51:97, 51:182, 51:192, 51:214, 51:216, 52:26, 52:39, 52:62-52:63, 78:07-78:08, C*14:92							7														
*14:56, 15:368, 18:01:01:01-18:01:13, 18:01:15-18:01:17, 18:01:19, 18:01:21-18:09, 18:12:01-18:14, 18:17N-18:19, 18:23N-18:25, 18:27-18:29, 18:31-18:34, 18:38-18:40, 18:42-18:49, 18:51, 18:53-18:55, 18:59-18:60, 18:62-18:68, 18:70-18:72:01, 18:72:03-18:85, 18:87-18:96, 18:98-18:100, 18:102-18:104, 18:107-18:109, 18:111-18:121, 18:123-18:124, 18:126-18:147, 39:43									9												
*15:03:01:01-15:03:07, 15:18:01:01-15:18:02, 15:18:04-15:18:07, 15:29, 15:64:01-15:64:02, 15:72, 15:74, 15:98, 15:103, 15:108, 15:127, 15:131-15:132, 15:134, 15:153, 15:158, 15:161, 15:173, 15:176, 15:197:01-15:198, 15:210, 15:220, 15:226N, 15:238, 15:242:01, 15:263, 15:266, 15:274, 15:281-15:282, 15:290, 15:292-15:293, 15:306-15:307, 15:311, 15:314, 15:337, 15:351, 15:354, 15:369, 15:380N, 15:382, 15:388-15:389, 15:414, 15:426, 15:433, 15:444, 15:447-15:448						6															21
*15:09:01-15:10:04, 15:37, 15:90, 15:99, 15:243, 15:312-15:313, 15:329, 15:338, 15:427, 15:435, 15:439, 15:445																	17				21
*15:16:01:01-15:16:03, 15:67, 15:95, 15:254, 15:362, 15:408, 35:58, 49:04:01-49:04:02, 58:06	1					6															
*15:18:03				4		6							14								21
*15:21:01:01-15:21:01:02, 15:44, 15:394, 15:404, 15:420							7			10											
*15:47:01-15:47:02, 15:49, 15:221, 15:379						w															21
*15:58, 35:205, 38:30, 38:65, 39:32-39:33, 39:74, 46:61, 51:13:01, 51:92:01, 52:14													13								
*15:62, 48:25	1					6				10											21
*15:69, 15:93						6															21
*15:73			3	4							w		14								
*15:80						6				10											21
*15:86	1			4							12										
*15:115, 15:364, 18:57, 51:45						6			9												
*15:133, 15:229, 15:253													13								21
*15:186, 39:07																					21
*15:194, 35:11:01-35:11:02, 53:08:01-53:08:02, 58:19, 58:91	1					6				10											
*15:223, 15:402, 57:22						6				10											
*15:224				4								12	14								
*15:242:02	1					6															21
*15:303	1		3								11	w									
*15:335						6			9												21
*15:361, 51:01:04, 51:01:13, 51:02:04, 51:06:01-51:06:04, 51:24:01-51:24:02, 51:38, 51:59, 51:64, 51:148, 51:194, 51:199, 51:213, 52:01:05, 52:01:15, 57:82								8									17				
*18:01:20, 18:72:02									9												20
*18:15						6			9												
*18:22, 18:105, 35:190, 58:76	1								9	10											
*18:30						w			9												
*18:37:01-18:37:02									9												21
*18:52	1								9												
*18:69, 18:149									9	10											
*18:125						6			9												21
*27:77, 27:129		2																			21
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21

HLA-B56 2018 0H2

*49:40			3							11											
*49:47	1		3	4																	
*50:01:02			3							11											
*50:04			3							11		14									
*50:12, 50:16, 50:41-50:42	1		3																		
*50:14	1						7												17		
*50:51			3	4																	
*51:01:01:01-51:01:03, 51:01:05-51:01:07, 51:01:09-51:01:12, 51:01:14, 51:01:16-51:01:17, 51:01:19-51:01:31, 51:01:33-51:01:38, 51:01:41-51:01:42, 51:01:44-51:01:58, 51:01:60-51:02:03, 51:02:05-51:03, 51:05, 51:07:01-51:09:03, 51:11N-51:12, 51:18-51:20, 51:22-51:23, 51:24:03-51:24:05, 51:26-51:30, 51:32-51:33, 51:35, 51:39-51:41N, 51:43-51:44N, 51:48-51:49, 51:51-51:53, 51:55, 51:57-51:58, 51:60-51:61:02, 51:65-51:77, 51:79-51:80, 51:83-51:84, 51:86-51:89, 51:91, 51:94-51:96, 51:98N-51:100, 51:102-51:105, 51:107-51:134, 51:136-51:138, 51:140-51:142, 51:145-51:147, 51:149N-51:152, 51:154-51:156, 51:158:01, 51:159-51:171, 51:173Q-51:181, 51:183-51:191, 51:193, 51:195-51:196, 51:198, 51:200-51:212, 51:215, 51:217-51:219, 51:221-51:222, 51:224, 51:226, 51:228-51:236, 52:01:01:01-52:01:04, 52:01:06-52:01:14, 52:01:16-52:01:20, 52:01:22-52:01:28, 52:01:30-52:07, 52:09:01-52:11, 52:13, 52:15-52:19, 52:21:01-52:24, 52:27-52:36, 52:38, 52:40-52:61, 52:64-52:66, 52:69-52:70, 52:72-52:74, 52:76, 58:08:01-58:08:02							7	8										17			
*51:01:18								7	8										17	20	
*51:01:32, 51:01:39, 51:01:59, 51:16-51:17, 51:21, 51:31, 51:34, 51:36, 51:82, 51:93, 51:101, 51:143-51:144, 51:153, 52:01:21, 52:08, 52:37, 52:67, 52:71								7	8												
*51:04:01-51:04:02, 51:46, 51:139, 57:02:01-57:03:03, 57:07, 57:09, 57:17, 57:39, 57:42, 57:46, 57:57, 57:63, 57:66, 57:70, 57:84, 57:94, 57:96									8	10									17		
*51:13:02, 51:92:02, 51:172, 52:75								7				13									
*51:37							6	7	9												
*51:56:01-51:56:03, 53:19, 57:05, 58:28:01-58:28:02, 58:65	1							8	10										17		
*51:62			3	4							12	14									
*51:106:01, 52:25:02			3					7			12										
*51:106:02, 52:25:01			3							11	12										
*51:223			3					7													
*54:01:01, 54:01:03-54:01:05, 54:01:07-54:02, 54:05N, 54:08N, 54:10, 54:12-54:13, 54:17-54:20, 54:22-54:23, 54:25, 54:27-54:30, 54:32, 54:34-54:36, 55:02:07	1	2									11	12							20		
*54:01:02, 55:35	1	2				5					11	12									
*54:01:06, 54:37		2									11	12							20		
*54:03	1	2	3								11								20		
*54:04, 54:15	1	2																	20		
*54:06, 54:11		2																	20		
*54:07, 54:24	1										11	12							20		
*54:09		2																	20		
*54:14	1	2																	20		
*54:16, 54:26	1	2									11								20		
*54:21	1	2									11	12					18		20		
*54:31	1	2									12								20		
*54:33	1	2	3								11	w							20		
*54:38		2		4							12								20		
*55:01:01-55:01:03, 55:01:05, 55:01:08-55:02:03, 55:02:05-55:02:06, 55:02:08-55:03, 55:05, 55:07, 55:09, 55:11-55:12, 55:15-55:16, 55:19, 55:21, 55:24-55:26, 55:29-55:30, 55:33, 55:36-55:37, 55:40, 55:43, 55:45, 55:48, 55:50, 55:52, 55:54-55:55N, 55:57, 55:59-55:60, 55:62-55:63, 55:65-55:67, 55:69-55:71, 55:73-55:77, 55:79-55:80, 55:82-55:85, 55:87-55:89N	1	2									11	12									
*55:01:04	1	2						7				12									
*55:01:06, 55:02:04, 55:72	1	2									11	12									
*55:01:07, 55:39, 55:42, 55:44, 55:53, 55:61	1	2									12										
*55:04, 55:17, 55:20, 55:49, 55:86, 81:01:01-81:01:02		2																			
*55:08, 55:51		2																	17		
*55:10	1	2									11	12							18		
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21
*55:14		2				6				10											
*55:23	1	2																			
*55:38		2									11	12		14							
*55:56	1	2									11	12			15						
*55:68	1	2				6															
*57:58										10								18			
*57:80								8		10							17	18			
*58:01:08	1									10										20	
*59:08	1						7		9		11	12									
*67:06-67:07																			?		
*81:03-81:04N, 81:06-81:07		?																			
*81:08		?																			21
*82:01:01:01-82:03	1	2	3	4											15			18			
*83:01		2								10						16					
A*02:339, C*15:26														w							
A*02:510, C*03:307, C*03:318N														14							
C*01:90				4								12									
C*01:136, C*01:145N, C*02:06:01-02:06:02, C*02:47, C*05:132, C*06:168, C*07:289, C*12:15, C*12:208, C*15:37, C*15:58, C*15:72				4																	
C*03:19			w												w						
C*03:102			w							10					w						
C*07:536				4																	
C*08:113, C*15:09, C*15:19, C*15:27, C*15:61, C*15:107						6															
C*15:91																					
C*15:128												12									
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21

22		
22		
Negative Control		
22	23	24
22	23	24

