

INTERPRETATION TABLE								
HLA-B*5701 SSP typing								
	Well ⁴							
	1	2	3	4	5	6	7	8
Length of spec.	90	220	95	105	90	90	165	195
PCR product			170		250	210	215	
Length of int.	800	800	800	1070	800	1070	1070	1070
pos. control ¹								
5'-primer(s) ²	209	362	362	209	209	320	362	103
	5' -ggC 3'	5' -ggT 3'	5' -ggT 3'	5' -ggC 3'	5' -ggA 3'	5' -CCC 3'	5' -ggT 3'	5' -CCT 3'
					362	362		
					5' -ggT 3'	5' -ggT 3'		
3'-primer(s) ³	256	539	418	271	256	2 nd	486	256
	5' -CCC 3'	5' -TCA 3'	5' -gTC 3'	5' -CAC 3'	5' -CCC 3'	5' -TCg 3'	5' -gCg 3'	5' -CCC 3'
			481		572	412	538	
			5' -gTA 3'		5' -gCg 3'	5' -gTT 3'	5' -gTC 3'	
			500					
			5' -ggA 3'					
Well No.	1	2	3	4	5	6	7	8
HLA-B allele								
*570101-570103	1	2						
*5702, 5712	1					6		
*570301, 570302, 5708, 5717	1	2				6		
*5704	1		3					
*5705, 5711, 5713, 5719	1							
*5706, 5718	1	2	3					
*5707	1	2			5	6		
*5709	1					6	7	
*5710	1	2						8
*5714	1	2					7	
*5715	1	2		4				
*5716		2			5			
*4030, 4034		2				6		
*5514, 5814		2					7	
HLA-B allele								
Well No.	1	2	3	4	5	6	7	8