

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
HLA-B allele ^{1,2,3}																				
*57:01:01-57:01:19	1	2									11									
*57:02:01-57:02:02, 57:12, 57:42	1					6														
*57:03:01-57:03:02, 57:17, 57:39, 57:66, 57:70	1	2				6														
*57:04	1		3														17			
*57:05, 57:11, 57:28N, 58:36	1																			
*57:06, 57:18, 57:27	1	2	3								11									
*57:07, 57:46	1	2			5	6														
*57:08	1	2				6					11	12								
*57:09	1					6	7													
*57:10, 57:44	1	2						8			11									
*57:13	1						7			10	11									
*57:14, 57:25, 57:31, 57:50	1	2					7				11									
*57:15, 57:29	1	2		4							11									
*57:16		2			5						11									
*57:19, 57:30, 57:52, 57:61	1										11									
*57:20, 57:62		2		4							11									
*57:21, 57:33	1	2							9		11									
*57:22, 57:43	1	2								10	11									
*57:23, 57:26, 57:34, 57:64	1	2			5						11									
*57:24	1						7				11									
*57:32	1	2	3																	
*57:35, 57:38, 57:60	1	2									11	12								
*57:36		2									11	12								
*57:37	1	2									11							18		
*57:40	1	2							9		w									
*57:41, 57:68	1	2									11						17			
*57:45, 57:51, 57:69		2									11		13							
*57:47, 57:54, 57:59	1	2									11			14						
*57:48, 57:55-57:56	1	2									11				15					
*57:49, 57:58, 57:71	1	2									11		13							
*57:53	1	2					7													
*57:57	1	2				6				10										
*57:63	1					6				10										
*57:65	1	2				6					11									
*57:67	1	2									11					16				
*57:72	1	2									11								19	
*07:120, 15:214, 18:81, 35:250, 40:150, 51:165				4																
*07:219								8					13							
*07:227							7			10										
*14:20, 35:127									9											
*15:33, 15:116, 15:248, 40:63, 40:92, 44:169, 44:182, 46:43, 49:22, 51:126, C*03:87, C*03:129, C*03:234, C*05:27, C*05:39, C*17:07														14						
*15:87												12								
*35:208		2									11									
*40:30, 40:34		2				6	7													
*44:153		2	3		5													17		
*55:14		2			5					10	11									
*58:12, A*02:42, A*02:285, A*02:310													13							
*58:14		2									11					16				
*58:54					5															
A*33:12															15					
C*06:72					5		7			w	11					16	17			
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20

Negative Control