

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21
HLA-B allele ^{1,2,3}																					
*57:01:01-57:01:22	1	2									11										
*57:02:01-57:02:02, 57:12, 57:42	1					6															
*57:03:01-57:03:02, 57:17, 57:66, 57:70	1	2				6															
*57:04	1		3														17				
*57:05, 57:11, 57:28N, 58:36	1																				
*57:06, 57:18, 57:27	1	2	3								11										
*57:07, 57:46	1	2			5	6															
*57:08	1	2				6					11	12									
*57:09	1					6	7														
*57:10, 57:44, 57:73	1	2						8			11										
*57:13	1						7			10	11										
*57:14:01-57:14:02, 57:25, 57:31, 57:50	1	2					7				11										
*57:15, 57:29	1	2		4							11										
*57:16		2			5						11										
*57:19, 57:30, 57:52, 57:61	1										11										
*57:20, 57:62		2		4							11										
*57:21, 57:33, 57:74	1	2							9		11										
*57:22, 57:43	1	2								10	11										
*57:23, 57:26, 57:34, 57:64	1	2			5						11										
*57:24	1						7				11										
*57:32	1	2	3																		
*57:35, 57:38, 57:60	1	2									11	12									
*57:36		2									11	12									
*57:37, 57:79N	1	2									11						18				
*57:39	1	2				6		8													
*57:40	1	2							9		w										
*57:41, 57:68	1	2									11						17				
*57:45, 57:69		2									11		13								
*57:47, 57:54, 57:59, 57:77	1	2									11			14							
*57:48, 57:55-57:56, 57:78	1	2									11				15						
*57:49, 57:58, 57:71	1	2									11		13								
*57:51		2						8			11		13								
*57:53	1	2					7														
*57:57	1	2				6				10											
*57:63	1					6				10											
*57:65	1	2				6					11										
*57:67	1	2									11					16					
*57:72	1	2									11							19			
*57:75	1	2									11								20		
*57:76	1	2																	19		
*07:120, 15:214, 18:81, 35:250, 40:150, 51:165				4																	
*07:219								8					13								
*07:227							7			10											
*14:20, 35:127, 51:186									9												
*15:33, 15:116, 15:248, 40:63, 40:92, 44:36, 44:169, 44:182, 46:43, 49:22, 51:126, C*03:87:01-03:87:02, C*03:129, C*03:232, C*03:234, C*05:27, C*05:39, C*08:115, C*17:07														14							
*15:87												12									
*15:340, 40:218, 46:28, 51:58, A*33:12															15						
*35:208		2									11										
*40:30, 40:34		2				6	7														
*44:153		2	3		5												17				
*55:14		2			5					10	11										
*58:12, A*02:42, A*02:285, A*02:310													13								
*58:14		2									11					16					
*58:41							8														
*58:54					5																
C*06:72					5		7			w	11					16	17				
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21

Negative Control