

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23
HLA-B allele ^{1,2,3}																							
*57:01:01-57:01:22	1	2									11												
*57:02:01-57:02:02, 57:12, 57:42	1					6																	
*57:03:01:01-57:03:02, 57:17, 57:66, 57:70	1	2				6																	
*57:04:01	1		3														17						
*57:04:02			3														17						
*57:05, 58:36	1																					22	
*57:06, 57:18, 57:27	1	2	3								11												
*57:07, 57:46	1	2			5	6																	
*57:08	1	2				6					11	12											
*57:09	1					6	7																
*57:10, 57:44, 57:73	1	2						8			11												
*57:11, 57:28N, 57:82-57:83	1																						
*57:13	1						7			10	11												
*57:14:01-57:14:02, 57:25, 57:31, 57:50	1	2					7				11												
*57:15, 57:29	1	2		4							11												
*57:16		2			5						11												
*57:19, 57:30, 57:52, 57:61	1										11												
*57:20, 57:62		2		4							11												
*57:21, 57:33, 57:74	1	2							9		11												
*57:22, 57:43	1	2								10	11												
*57:23, 57:26, 57:34, 57:64	1	2			5						11												
*57:24	1						7				11												
*57:32	1	2	3																				
*57:35, 57:60	1	2									11	12											
*57:36		2									11										21		
*57:37, 57:79N, 57:81	1	2									11						18						
*57:38	1	2									11										21		
*57:39	1	2				6		8															
*57:40	1	2							9		w												
*57:41, 57:68	1	2									11						17						
*57:45, 57:69		2									11		13										
*57:47, 57:54, 57:59, 57:77	1	2									11			14									
*57:48, 57:55-57:56, 57:78	1	2									11				15								
*57:49, 57:58, 57:71	1	2									11		13										
*57:51		2						8			11		13										
*57:53	1	2					7																
*57:57	1	2				6				10													
*57:63	1					6				10													
*57:65	1	2				6					11												
*57:67:01	1	2									11					16						22	
*57:67:02	1	2									11											22	
*57:72	1	2									11							19					
*57:75	1	2									11								20				
*57:76	1	2																19					
*57:80	1	2				6							13										
*07:120, 15:214, 18:81, 40:150, 51:165				4																			
*07:219								8					13										
*07:227							7			10													
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23

Negative Control

[illegible]