

INTERPRETATION TABLE			
HLA-B*73 SSP subtyping			
Amplification patterns of the B*73:01 to 73:02 alleles			
	Well		
	1	2	3
Length of spec.	150	220	220
PCR product			
Length of int.	800	1070	1070
pos. control¹			
5'-primer²	165	379	379
	5' -ACC 3'	5' -ACA 3'	5' -ACA 3'
3'-primer³	272	556	556
	5' -TgC 3'	5' -gCC 3'	5' -gCT 3'
Well No.	1	2	3
HLA-B allele⁴			
*73:01	1	2	
*73:02	1		3
*27:75	1		
HLA-B allele			
Well No.	1	2	3