

Well No.	1	2	3	4	5	6	7	8
HLA-B allele ^{1,2,4}								
*81:01	1	2		4				Negative Control
*81:02, 07:174 ³	1		3	4				
*81:03	1	?	?	4	5			
*81:04N	1	?	?	4		6		
*81:05	1	?	?	4	5	6		
*81:06	1	?	?	4			7	
*07:02:01-07:02:38, 07:04-07:06, 07:09-07:11, 07:15, 07:17, 07:20-07:26, 07:28-07:31, 07:34-07:36, 07:39, 07:41-07:46, 07:48-07:49N, 07:51-07:52, 07:54-07:59, 07:61-07:64, 07:66-07:68:03, 07:70-07:76, 07:79-07:82, 07:86-07:99, 07:101-07:109, 07:111N-07:121, 07:123-07:136, 07:138, 07:140-07:161N, 07:164, 07:166-07:173, 07:175-07:179, 07:181N-07:192, 07:194-07:195, 42:05:01-42:05:02			3	4				
*07:02:39, 07:07, 07:12, 07:14, 07:18:01-07:19, 07:33:01-07:33:02, 07:40, 07:53, 07:60, 07:77-07:78, 07:84, 07:100, 07:122, 07:137, 07:139, 07:162-07:163, 07:165, 07:193, 07:196-07:198, 14:21, 38:26, 42:01:01-42:02, 42:04, 42:06, 42:08-42:10, 42:12-42:20, 45:06, 67:01:01-67:01:03, 67:03, 73:01-73:02			3					
*07:03, 07:08, 07:16, 07:27, 07:32, 07:37-07:38, 07:47, 07:50, 07:69, 07:85, 07:180, 08:79, 37:07, 40:15-40:16, 40:32, 40:98, 48:05, 48:08, 48:10, 48:12, 48:14				4				
*07:13, 07:110			3	4			7	
*07:65		2		4				
*07:83			3	4	5	6		
*15:20, 15:228, 35:01:01:01-35:01:27, 35:01:29-35:30, 35:32:01-35:42:02, 35:44-35:56, 35:58-35:59, 35:61:01-35:65Q, 35:67-35:72, 35:74-35:75, 35:77-35:78, 35:80-35:86, 35:88-35:184, 35:187-35:228, 35:231-35:233, 48:02:01-48:02:03, 51:01:01-51:24:04, 51:26-51:41N, 51:43-51:46, 51:48-51:154, 52:01:01:01-52:06:02, 52:08-52:31, 53:01:01-53:31, 58:01:01-58:02, 58:04-58:10N, 58:12-58:19, 58:21-58:29, 58:31N-58:46, 78:01:01-78:07					5			
*27:01-27:11, 27:13-27:15, 27:17, 27:19-27:21, 27:25-27:28, 27:30-27:38, 27:40-27:58, 27:60-27:76, 27:78-27:91, 27:93-27:106, 44:90, 44:97, 54:01:01-54:01:04, 54:03-54:29, 55:01:01-55:05, 55:07-55:17, 55:19-55:59, 56:01:01-56:04, 56:06-56:16, 56:18-56:20:02, 56:22, 56:24-56:32, 56:34-56:35, 56:37-56:39, 56:41, 82:01-82:03		2						
*27:24	1	2						
*35:66				4	5			
*35:76, 35:186, 35:229, 56:05:01-56:05:02, 56:21, 56:36		2			5			
*38:11, C*03:137, C*03:182, C*12:82							7	
*40:31, 40:45, 40:80, 48:01:01-48:01:04, 48:06-48:07, 48:09, 48:11, 48:15-48:16, 48:18-48:20, 48:22, 48:27-48:31	1			4				
*48:04, 48:24	1							
*54:02		2					7	
*67:02			3				7	
*83:01		2			?			
C*07:102							w	
Well No.	1	2	3	4	5	6	7	8