

INTERPRETATION TABLE					
HLA-B*82 SSP subtyping					
Amplification patterns of the B*82:01 to 82:03 alleles					
	Well				
	1	2	3	4	5
Length of spec.	195	140	230	195	155
PCR product					
Length of int.	800	800	1070	1070	1070
pos. control¹					
5'-primer²	420	557	368	105	105
	5' -TTA 3'	5' -ggA 3'	5' -gTT 3'	5' -gCT 3'	5' -gCT 3'
3'-primer³	572	3rd I	557	259	219
	5' -gCg 3'	5' -TAT 3'	5' -ggC 3'	5' -gTT 3'	5' -ggg 3'
				259	
				5' -gTT 3'	
Well No.	1	2	3	4	5
HLA-B allele					
*82:01	1	2		4	
*82:02	1		3	4	
*82:03	1	2		4	5
*15:06, 15:27:01-15:27:03, 15:84, 15:109, 15:195, C*04:08, C*04:34			3		
*44:10, 44:15, 44:18, 45:01, 45:05-45:07, 45:11-45:13, 50:02	1				
HLA-B allele					
Well No.	1	2	3	4	5