

HLA-B82 2019 4H2

	A	B	C	D	E	F	G	H
1	Well No.	1	2	3	4	5	6	
2	HLA-B allele^{1,2}							
3	*82:01:01:01-82:01:02	1	2		4			
4	*82:02:01:01-82:02:02	1		3	4			
5	*82:03	1	2			5		
6	*15:06, 15:27:01-15:27:04, 15:84, 15:109, 15:195, 15:327, 15:344, 15:398, C*03:89, C*03:271, C*03:338, C*04:08, C*04:34, C*04:147, C*04:212, C*18:08			3			Negative Control	
7	*44:10, 44:15:01:01-44:15:01:02, 44:18, 44:140, 45:01:01-45:01:04, 45:05-45:07, 45:11-45:23, 49:20, 50:02, 50:56	1						
8	Well No.	1	2	3	4	5	6	