

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	
HLA-C allele ^{1,2}																																	
*01:02:01:01-01:02:17, 01:02:19-01:02:33, 01:02:35-01:02:48, 01:02:50-01:02:67, 01:02:69-01:02:70, 01:02:72-01:02:76, 01:06-01:07:01, 01:08, 01:11-01:13, 01:15:01-01:20, 01:23, 01:25-01:33, 01:37N-01:42, 01:44-01:47, 01:51-01:54, 01:56N-01:58, 01:61-01:66, 01:68-01:77, 01:80-01:89N, 01:91-01:96, 01:98N-01:100, 01:103-01:112, 01:115, 01:117N, 01:119, 01:121Q, 01:123-01:124, 01:126-01:130, 01:132-01:135, 01:137N-01:144, 01:145:02N, 01:147-01:153, 01:155-01:157, 01:159, 01:161-01:167, 01:170-01:175, 01:178-01:185Q, 01:187-01:211N. 01:213	1																			20						26							
*01:02:18	1											12									20						26						
*01:02:34	1																				20						26	27					
*01:02:49, 01:136, 01:145:01N, 01:146:01-01:146:02	1							8													20						26						
*01:02:68, 01:02:71	1														15						20						26						
*01:03:01-01:03:02, 01:24, 01:34, 01:78	1																				20	21					26						
*01:04	1											12		14							20												
*01:05, 01:07:02, 01:22, 01:35-01:36:02, 01:79:01-01:79:02, 01:120, 01:158	1																				20												
*01:09:01-01:09:02	1													14							20												
*01:10	1	2																			20						26						
*01:14	1											11									20						26					31	
*01:21	1											12	13															27		30			
*01:43	1	2									10										20						26						
*01:48	1		3																		20						26						
*01:49:01-01:50, 01:55, 01:176	1																				20	21											
*01:59, 01:118	1											11									20						26						
*01:60	1																		19		20						26						
*01:67, 01:122	1			4												16					20						26						
*01:90, 01:101-01:102, 01:113	1							8									17				20						26						
*01:97, 01:114:01-01:114:02, 01:168-01:169:01	1																				20						26		28				
*01:116	1																17				20						26						
*01:125	1															16					20						26						
*01:131	1			4																	20	21											
*01:154	1																				20						26					31	
*01:160	1													14				18			20												
*01:169:02, 03:94 ⁴	1		3																		20		22				26		28				
*01:177	1															16					20	22				25	26						
*01:186	1											11									20					25	26						
*02:02:01-02:02:03, 02:02:06-02:02:12, 02:02:15-02:02:25, 02:02:27-02:02:28, 02:02:30-02:02:34, 02:02:36, 02:02:39-02:02:43, 02:02:45-02:02:52, 02:02:54-02:02:61, 02:04, 02:07-02:09, 02:11, 02:15, 02:19-02:20, 02:24-02:25Q, 02:28-02:31, 02:34-02:35, 02:37-02:38:01N, 02:42-02:43:01, 02:44, 02:46, 02:48, 02:50, 02:52N-02:54, 02:56-02:57, 02:59-02:61:02, 02:63-02:64, 02:66-02:67Q, 02:69-02:71, 02:74, 02:77-02:86, 02:88, 02:90-02:93, 02:96-02:100, 02:102-02:106, 02:108, 02:112-02:114, 02:116-02:117, 02:119-02:123, 02:127-02:130, 02:132, 02:136, 02:138-02:142, 02:144-02:151, 02:153-02:160, 02:162-02:163, 02:165N-02:166, 02:170-02:173, 02:175, 02:177-02:179, 02:183-02:186, 02:188, 02:190, 02:192N-02:196, 02:198-02:203		2	3			6					11				15									24								31	
*02:02:13		2	3									11			15							22										31	
*02:02:14, 02:02:37, 02:10:01:01-02:10:07, 02:32, 02:38:02N, 02:43:02, 02:51, 02:73, 02:76, 02:89, 02:95, 02:109-02:110, 02:124-02:125, 02:133, 02:143, 02:161, 02:168-02:169N, 02:174, 02:176, 02:181-02:182		2	3			6						11			15																		31
*02:02:26, 02:21, 02:26:01-02:26:03, 02:39, 02:45, 02:189, 02:191		2	3									11			15									24									31
*02:02:29		2	3									11			15							22		24									31
*02:02:35, 02:03, 02:16:02:01-02:16:02:02, 02:72, 02:135N, 02:137		2	3			6						11			15																		
*02:02:38		2	3			6						11			15							22		24									31
*02:02:44		2	3			6						11	12		15									24									31
*02:02:53		2	3			6						11			15							22											31
*02:05:01, 02:17		2	3			6						11		14	15																		31
*02:05:02-02:05:03		2	3			6						11		14	15									24									31
*02:06:01, 02:36:02		2	3			6		8				11			15		17																31
*02:06:02, 02:23, 02:36:01		2	3			6		8				11			15		17						24										31
*02:12		2	3			6						11		13	15												26	w			30		31
*02:13, 02:75		2	3			6						11			15			18						24									31
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	

Negative Control

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	
*02:14:01-02:14:02		2	3			6					11																					31	
*02:18		2	3			6					11				15			18															
*02:22, 02:62, 02:118		2	3			6					11				15						21			24								31	
*02:27:01, 02:131		2	3			6					11												24				27						
*02:27:02		2	3			6					11																27						
*02:33		2	3			6					11				15			18														31	
*02:40:01		2	3			6					11				15				19					24								31	
*02:40:02		2	3			6					11				15				19													31	
*02:47, 02:111		2	3			6		8			11				15									24								31	
*02:49		2	3			6					11		13		15			18						24		26				30	31		
*02:55:01		2	3			6					11		13		15									24		26				30	31		
*02:55:02		2	3			6					11		13		15											26				30	31		
*02:58			3			6					11				15																	31	
*02:65		2	3			6					11													24									
*02:68		2	3			6					11				15		17							24								31	
*02:87		2								10														24				27					
*02:94		2	3			6	7				11				15																	31	
*02:101		2	3			w					11				15									24								31	
*02:107		2									11													24								31	
*02:115		2	3			6					11		13					18						24		26	27			30			
*02:126		2	3			6					11				15									24			w					31	
*02:134, 02:152, 02:167		2									11				15									24								31	
*02:164, 02:187		2	3			6					11													24								31	
*02:180		2	3			6					11				15			18	19														
*02:197			3			6					11				15						21											31	
*03:02:01-03:02:09, 03:02:11-03:02:25, 03:05:01:01-03:05:01:02, 03:17:01-03:17:02, 03:25, 03:33, 03:35:01-03:36, 03:60, 03:71, 03:84, 03:89, 03:95, 03:108, 03:110, 03:121N, 03:135, 03:139, 03:146, 03:154, 03:169Q, 03:178, 03:190, 03:194, 03:197-03:198, 03:200-03:201N, 03:216, 03:221-03:222, 03:224N-03:226, 03:238, 03:249, 03:258, 03:292, 03:296:01-03:296:02, 03:298-03:301, 03:315, 03:326, 03:329-03:330, 03:332, 03:335, 03:338, 03:349-03:350, 03:365, 03:371, 03:373, 03:386, 03:391N, 03:402-03:403, 03:406-03:407, 03:411-03:412, 03:425, 03:429, 03:449N, 03:452, 03:464, 03:468, 03:473, 03:475, 03:482, 03:489, 03:493-03:494, 03:497, 03:505, 03:532, 03:544, 03:547, 03:549, 03:561			3	4																		22											
*03:02:10, 03:264			3	4																													
*03:03:01:01-03:03:20, 03:03:22-03:03:50, 03:03:53-03:03:61, 03:11:01-03:11:02, 03:18:02-03:18:03, 03:20N-03:22Q, 03:30, 03:43:01-03:43:02, 03:49-03:50, 03:52-03:53, 03:55-03:56, 03:59, 03:61-03:62, 03:66, 03:68, 03:75-03:76, 03:79, 03:83, 03:85, 03:88, 03:96, 03:102-03:103, 03:112, 03:116:01-03:116:02, 03:119:01-03:120, 03:122, 03:124, 03:126-03:127, 03:133, 03:141-03:142, 03:150, 03:152, 03:158, 03:160-03:161, 03:168, 03:171, 03:176, 03:185, 03:187-03:189N, 03:192, 03:196, 03:202-03:207, 03:214, 03:217, 03:220, 03:223, 03:227-03:230, 03:237, 03:241-03:242, 03:253-03:254, 03:272-03:273, 03:275-03:276, 03:284-03:285, 03:288-03:291, 03:304, 03:307, 03:312, 03:316N, 03:319-03:321, 03:325, 03:327, 03:336, 03:339, 03:341, 03:345-03:346, 03:352, 03:356-03:357, 03:360, 03:363N-03:364, 03:367, 03:370, 03:372:01-03:372:02, 03:374-03:375, 03:377N-03:378, 03:380N, 03:383, 03:389, 03:395, 03:398, 03:413-03:414, 03:418, 03:421N-03:422, 03:427-03:428, 03:430, 03:432N-03:433, 03:436, 03:438, 03:440, 03:444N, 03:447N, 03:457, 03:460, 03:471-03:472, 03:476-03:477, 03:481, 03:490, 03:495, 03:503, 03:506-03:507, 03:509N, 03:514-03:517, 03:519-03:520, 03:524, 03:527-03:529, 03:533, 03:537, 03:541-03:542, 03:545, 03:550-03:551, 03:556, 03:559-03:560N				4	5																21	22											
*03:03:21, 03:03:51, 03:18:01, 03:39, 03:64:02, 03:144, 03:182, 03:195, 03:344:01-03:344:01:02, 03:504				4																		21											
*03:03:52, 03:151, 03:165:01-03:165:02, 03:251, 03:267, 03:451, 03:534-03:535, 03:540					5																	21	22										
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	

Negative Control

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	
*03:04:01:01-03:04:06, 03:04:08-03:04:14, 03:04:16-03:04:24, 03:04:27-03:04:49, 03:06:01-03:06:02, 03:09, 03:14, 03:19, 03:23N-03:24, 03:26:01-03:26:02, 03:28, 03:32, 03:34, 03:37:01-03:37:02, 03:40:01-03:40:04, 03:42:01-03:42:02, 03:44, 03:46-03:48, 03:51, 03:54, 03:57:01-03:57:02, 03:63-03:64:01, 03:65, 03:70, 03:72-03:74, 03:77-03:78, 03:80:02, 03:82, 03:87:01-03:87:02, 03:90-03:93, 03:98, 03:100-03:101, 03:104-03:107, 03:109, 03:111, 03:114-03:115:02, 03:117, 03:123, 03:125, 03:128-03:129, 03:131, 03:134, 03:138, 03:143, 03:147-03:149, 03:153, 03:157, 03:159, 03:162, 03:164, 03:173-03:174, 03:179-03:181, 03:183-03:184:02, 03:186:01-03:186:02, 03:191, 03:193, 03:208N, 03:210-03:213, 03:215, 03:218-03:219, 03:232-03:236, 03:239-03:240, 03:244Q, 03:247, 03:250, 03:252, 03:255-03:257, 03:259, 03:261, 03:263:01-03:263:03, 03:265N-03:266, 03:269-03:270, 03:277N, 03:279, 03:281-03:283, 03:287:02, 03:294, 03:303, 03:305-03:306, 03:309-03:311, 03:313-03:314, 03:318N, 03:322-03:323N, 03:328, 03:331, 03:333-03:334, 03:337, 03:340, 03:343, 03:347-03:348, 03:353-03:355, 03:358-03:359, 03:361-03:362, 03:366N, 03:368-03:369, 03:376, 03:379, 03:381-03:382, 03:385, 03:387, 03:390, 03:392N-03:394, 03:396:01N-03:397, 03:399-03:401, 03:404:01:01-03:405, 03:408-03:410, 03:415, 03:417, 03:419-03:420, 03:423-03:424N, 03:426, 03:434, 03:437, 03:439, 03:441-03:443, 03:445N-03:446N, 03:448Q, 03:453, 03:455-03:456, 03:458, 03:462N-03:463N, 03:465-03:467, 03:469-03:470N, 03:478, 03:480, 03:483-03:488, 03:491:01-03:492, 03:496, 03:498-03:502, 03:508N, 03:510N-03:511N, 03:513, 03:518, 03:521-03:523, 03:525-03:526, 03:531N, 03:538-03:539, 03:543, 03:546, 03:548, 03:552-03:554, 03:558																																	
*03:04:07, 03:41:01-03:41:02, 03:80:01, 03:113:01-03:113:02, 03:260, 03:286, 03:293, 03:384, 03:435, 03:512				3																		21	22										
*03:04:15, 03:04:26, 03:118, 03:137, 03:145, 03:156, 03:166, 03:170, 03:172, 03:342, 03:454, 03:459				3	4																	21											
*03:04:25				3	4																	21						27					
*03:07:01:01-03:07:02, 03:45				3	4						11				15							21	22										
*03:08				3	4																	21	22								29		
*03:10				3	4										15							21	22										
*03:13:01:01-03:13:02, 03:97, 03:132, 03:167, 03:177, 03:271, 03:474					4	5																22											
*03:15				3	4						11				15							22				26							
*03:16, 03:536, B*46:77 ³				3																		22											
*03:27, 03:431				3	4																	22				26							
*03:29				3	4										15							21	22								29		
*03:31				4	5																	21	22								29		
*03:38:01-03:38:02, 03:136				3	4																	21	22			26							
*03:58	1				5															20	21	22			26		28						
*03:67, 03:295				4	w																	21	w										
*03:69, 03:274				4	5																	21	22			26							
*03:81, 03:175				4	5																	21	22										
*03:86	1			5														17			20	22			26		28						
*03:130				3	4						11											21	22			26							
*03:140:01:01-03:140:01:02				3	4						11											21	22										
*03:155				3	4											16						21	22										
*03:163, 03:297				3	4						11				15							21	22			26							
*03:199, 03:245				3	4													17				22											
*03:209				3	4																	21	22						28				
*03:231				4		6					11				16							21			25								
*03:243				4	5						11											21	22										
*03:246				3	4																	21	22			26			29				
*03:248				3	4						11														25								
*03:262				4	5																	?	22										
*03:268				4	5						11				15							21	22										
*03:278				3	4																	21	w										
*03:280				3	4					10												21	22										
*03:287:01				3	4																	21	22		24								
*03:302				3	4										15							22						28					
*03:308				4	5						11											21	22		25								
*03:317, 03:388				3	4													17				21	22										
*03:324, 03:479				4	5																	21											
*03:351				4	5										16							21	22										
*03:416	1			4	5																	21	22										
*03:450				4	5										15							21	22										
*03:461				4	5						11				15							21	22			26							
*03:530				4						10												21											
*03:555				3	4																	22							29				
*03:557				4											16							21	22		25								
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32				
*04:01:01:01-04:01:01:29, 04:01:01:31-04:01:15, 04:01:17-04:01:22, 04:01:24-04:01:27, 04:01:29-04:01:54, 04:01:56-04:01:60, 04:01:62, 04:01:64, 04:01:66-04:01:73, 04:01:75-04:01:76, 04:01:78-04:01:83, 04:01:85-04:01:139, 04:05-04:07:01, 04:08, 04:12, 04:15:02, 04:17-04:20, 04:24-04:28, 04:30-04:31, 04:35:01-04:35:02, 04:38-04:41, 04:43-04:53, 04:56-04:57, 04:59Q-04:67, 04:69-04:70, 04:73-04:76, 04:78-04:79, 04:81-04:99, 04:101-04:102, 04:104-04:106, 04:109-04:110, 04:113, 04:115N, 04:117-04:121, 04:123N-04:128, 04:131-04:139, 04:141-04:146, 04:148-04:152, 04:154-04:159, 04:161-04:166:01, 04:166:03-04:168, 04:170N, 04:173N-04:177, 04:179-04:187, 04:189, 04:191N, 04:193, 04:195, 04:197-04:201:01, 04:202-04:209, 04:211, 04:216-04:219, 04:221-04:225N, 04:228-04:229, 04:232-04:235, 04:237-04:248, 04:250-04:255N, 04:257-04:262, 04:264, 04:266-04:277, 04:280-04:281, 04:284-04:285, 04:287-04:290, 04:292-04:293, 04:295-04:298, 04:300N-04:303, 04:305N-04:310, 04:312-04:334, 04:336, 04:338Q-04:340, 04:342-04:346, 04:348-04:350N, 04:352-04:356, 04:358-04:359, 04:361-04:362N, 04:364N-04:366, 04:369N-04:380, 04:382Q, 04:384-04:392, 04:394-04:396N, 04:398-04:399, 04:401, 04:403-04:405, 04:407-04:413, 04:415-04:421, 04:423-04:431, 04:433-04:435, 04:437-04:438, 04:440-04:441																																				
*04:01:16, 04:54:01-04:54:02, 04:72, 04:341, 04:439						6					11	13					18								25							31				
*04:01:23, 04:33, 04:116, 04:130, 04:166:02, 04:231, 04:283						6					11	13								21				25								31				
*04:01:28, 04:01:61, 04:01:63, 04:01:65, 04:140, 04:196, 04:236N, 04:304:01-04:304:02:02, 04:311											11	13					18		21					25								31				
*04:01:55, 04:15:01, 04:15:03						6					11	13					18		21				25							30	31					
*04:01:74						6					11	13					w		21				25									31				
*04:01:77, 04:04:01:01-04:04:02, 04:13:01:01-04:13:01:02, 04:34, 04:111, 04:122, 04:188, 04:282, 04:291						6					11	13					18		21					25												
*04:01:84, 04:09N						6					11	13					18		21					25			28					31				
*04:03:01:01-04:03:01:02, 04:03:03-04:03:07, 04:80, 04:147, 04:171, 04:190, 04:256, 04:286, 04:337, 04:363, 04:393, 04:400, 04:402			3			6					11	13		15			18		21					25	26							31				
*04:03:02			3			6					11	13		15			18		21					25									31			
*04:06:01-04:06:03, 04:381			3			6					11	13		15			18		21					25	26											
*04:07:02											11									21	22			25									31			
*04:10, 04:153, 04:192, 04:194, 04:210, 04:215N, 04:227, 04:278, 04:347						6					11						18		21					25									31			
*04:11						6					11				16					21				25		27										
*04:14, 04:68						6					11	13					18		21					25						30						
*04:16			3								11	13		15			18		21					25	26								31			
*04:23						6		8			11	13					18							25										31		
*04:29						6					11	13			16		18			21					25		27									
*04:32, 04:77			2			6					11	13					18		21	22				25									31			
*04:36, 04:55						6					11				16		18		21					25		27										
*04:37			1			6					11	13		15			18		20	21				25	26								31			
*04:42:01-04:42:02, 04:220											11	13	14				18		21		23			25									31			
*04:58, 04:368						6					11	13					18	19		21				25												
*04:71, 04:436N			1			6					11	13					18		21					25									31			
*04:100						6		8			11	13					18		21					25										31		
*04:103:01-04:103:02											11	13		15			18		21					25	26									31		
*04:107			3			6					11	13		15						21				25	26									31		
*04:108						6		8			11	13					17	18						25										31		
*04:112						6					11	13					17	18		21				25										31		
*04:114						6					11	13			16		18		21					25												
*04:129						6	7				11	13					18		21					25	26									31		
*04:160:01-04:160:02			3			6					11	13		15			18	19	21					25	26											
*04:169						6					11						17			21				25										31		
*04:172						6					11	13			16					21				25		27										
*04:178, 04:422						6		8			11	13					17	18		21				25												
*04:201:02				4		6					11	13					18		21					25										31		
*04:212				4		6					11	13					18		21					25												
*04:213											11	13					18		21															31		
*04:214						6					11						18		21					25			w							31		
*04:226						6					11	13					18		21																31	
*04:230						6					11	13					18		20	21				25	26									31		
*04:249						6					11						18			21				25												
*04:263						6					11	13					18		20	21				25	26		28							31		
*04:265			2			6					11	13					18			21	22			25												
*04:279N						6					11	w					18		21					25										31		
*04:294			3			6					11	13		15			18		21					25	26					30						
*04:299			3			6					11			15			18		21					25												
*04:335			3			6		8			11	13		15			18							25	26									31		
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32				

Negative Control

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
*04:351:01-04:351:02			3								11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*04:357		2	3			6					11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*04:360						6					11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*04:367						6					11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*04:383			3			6					11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*04:397, 04:442						6					11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*04:406		2				6					11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*04:414						6					11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*04:432						6					11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*05:01:01:01-05:01:50, 05:01:52-05:01:55, 05:01:57, 05:01:59, 05:01:61-05:01:64, 05:03, 05:05:01-05:10, 05:12-05:15, 05:18:01-05:19, 05:21, 05:24, 05:28-05:31, 05:33-05:35, 05:37-05:41, 05:44:01-05:45, 05:47-05:50, 05:52-05:54, 05:56-05:57, 05:59-05:61, 05:63-05:67, 05:69, 05:71-05:75, 05:77, 05:80-05:84, 05:86-05:102, 05:104, 05:108-05:111, 05:114, 05:116-05:122, 05:124-05:128N, 05:130-05:131, 05:133, 05:136-05:140, 05:142, 05:144-05:145, 05:148-05:150, 05:152-05:155, 05:157-05:162, 05:164-05:183, 05:185-05:190, 05:192-05:199, 05:201, 05:203-05:205, 05:208N-05:214, 05:216:01:01-05:237, 05:239N-05:240, 05:242-05:251N, 05:253N-05:255							7				11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*05:01:51, 05:01:58, 05:04:01-05:04:02, 05:103:01-05:103:02, 05:123, 05:135, 05:147, 05:191, 05:207							7				11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*05:01:56							7	8			11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*05:01:60, 05:106:01-05:106:02							7				11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*05:11, 05:27, 05:32, 05:51Q, 05:70, 05:184							7				11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*05:16							7				11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*05:17, 05:68, 05:79, 05:115							7				11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*05:20							7				11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*05:22:01-05:22:02, 05:252							7	8			11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*05:23							7				11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*05:25, 05:55, 05:76							7				11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*05:26						6	7				11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*05:36							7				11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*05:42							7				11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*05:43, 05:113N, 05:156, 05:202Q							7				11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*05:46:01:01-05:46:01:02							7				11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*05:58:01-05:58:04			3				7				11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*05:62, 05:146, 05:151, 05:163							7				11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*05:78:01-05:78:02						6	7				11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*05:85	1						7				11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*05:105		2					7				11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*05:107	1						7				11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*05:112							7				11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*05:129							7				11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*05:132, 05:215							7	8			11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*05:134							7				11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*05:141							7				11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*05:143							7				11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*05:200	1	2					7				11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*05:206		2					7				11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*05:238							7				11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*05:241	1						7	8			11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*06:02:01:01-06:02:01:60, 06:02:03-06:02:11, 06:02:13-06:02:15, 06:02:17-06:02:36, 06:02:40-06:02:42, 06:02:44-06:02:64, 06:02:66-06:02:71, 06:02:73-06:02:90, 06:02:92, 06:07, 06:10, 06:12-06:13, 06:15-06:16N, 06:18-06:30, 06:33-06:34:02, 06:36-06:39, 06:41-06:51, 06:53:01-06:66, 06:68-06:69, 06:71, 06:73-06:75, 06:78, 06:80, 06:83-06:88, 06:90-06:93, 06:96-06:98, 06:100, 06:102:01-06:117, 06:119, 06:121, 06:125-06:126, 06:128N-06:131, 06:133-06:135, 06:137-06:142, 06:145-06:146, 06:148-06:152N, 06:154-06:159, 06:161-06:179, 06:181-06:182, 06:185-06:196, 06:198-06:202, 06:205-06:209, 06:211:01:01N-06:216, 06:218-06:223, 06:225-06:235, 06:237-06:243, 06:245-06:247, 06:249-06:251, 06:253-06:263N, 06:265-06:270, 06:272-06:279, 06:281N-06:303, 06:306-06:316N, 06:318-06:325							8				11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
*06:02:12, 06:122, 06:124, 06:280							8				11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32

Negative Control

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	
*06:02:16, 06:02:43, 06:81, 06:101, 06:123, 06:127:01-01-06:127:01:02, 06:136, 06:203, 06:304, 06:317								8			11												23								31		
*06:02:37								8			11			14									23							30	31		
*06:02:38								8			11	12		14									23								31		
*06:02:39, 06:32, 06:94-06:95, 06:99, 06:180, 06:236, 06:244											11			14									23								31		
*06:02:65, 06:11, 06:52, 06:147, 06:217, 06:248, 06:252								8						14									23										
*06:02:72											11		13	14			18						23								31		
*06:02:91, 06:82, 06:160, 06:183-06:184, 06:210								8			11			14									23										
*06:03:01, 06:132:01			3					8			11			14	15							22				26					31		
*06:03:02, 06:76:01								8			11			14	15											26					31		
*06:04:01-06:04:02:02, 06:40, 06:153:01-06:153:02, 06:197								8			11												23										
*06:05							7	8			11						18		20				23	24		26					31		
*06:06:01:01-06:06:01:02, 06:271								8			11								20				23			26					31		
*06:08		2						8			11			14									23								31		
*06:09:01-06:09:02								8			11			14							21		23								31		
*06:14, 06:35, 06:72, 06:143-06:144								8			11										21		23								31		
*06:17											11			14					20				23								31		
*06:31								8			11			14			18		20				23								31		
*06:67							7	8			11			14									23	24		26					31		
*06:70:01-06:70:02, 06:305			3					8			11			14									23								31		
*06:76:02											11	12	13	14	15		18									26					31		
*06:77								8			11			14																	31		
*06:79N								8			11								20				23								31		
*06:89, 06:224								8			11			14			17			20			23								31		
*06:118								8			11			14			18	19					23									31	
*06:120								8			11			14									23		25						31		
*06:127:02, 06:264								8			11												23							30	31		
*06:132:02			3								11			14	15							22				26					31		
*06:204	1							8			11									20			23			26					31		
*07:01:01-07:01:02:10, 07:01:04-07:01:06, 07:01:08-07:01:40, 07:01:42-07:01:73, 07:01:75-07:01:92, 07:01:94-07:02:07, 07:02:10-07:02:28, 07:02:30-07:02:40, 07:02:42-07:02:90, 07:02:92-07:04:02:02, 07:04:04-07:06:04, 07:08, 07:11-07:13:01:02, 07:15-07:19, 07:21-07:25, 07:27:01, 07:29:01-07:30, 07:32N, 07:35-07:38:02, 07:42, 07:44-07:48, 07:50-07:52, 07:54-07:63, 07:65-07:75, 07:77-07:91, 07:93:01-07:95, 07:97-07:100, 07:102-07:122, 07:124, 07:126, 07:128-07:132, 07:133:02-07:138, 07:140:01-07:147, 07:149-07:151, 07:153, 07:155, 07:157-07:160, 07:162-07:172:02, 07:174-07:176, 07:178-07:183, 07:185-07:195, 07:197-07:199:01, 07:200-07:209, 07:211-07:215, 07:217-07:237, 07:239-07:241, 07:243-07:245, 07:248-07:262, 07:264N-07:266, 07:268-07:292, 07:297-07:313, 07:318-07:321, 07:323-07:327, 07:329N, 07:330:02-07:334, 07:336-07:337, 07:339-07:350N, 07:352-07:355, 07:357-07:366, 07:368:01-07:389, 07:391-07:401, 07:404, 07:405:02, 07:407-07:414, 07:417-07:448, 07:450-07:492, 07:494Q-07:513Q, 07:515-07:516, 07:518, 07:520-07:525, 07:526:02-07:530, 07:532-07:535, 07:537-07:561, 07:564-07:567, 07:569-07:577, 07:580-07:582Q, 07:584-07:597, 07:599-07:625, 07:627-07:680, 07:682:01-01-07:713Q, 07:715:01-07:722, 07:726N-07:747N, 07:749N-07:755, 07:757-07:774, 07:776N-07:820N, 07:822-07:844, 07:846:01-07:849N, 07:851-07:856N, 07:858-07:869, 07:871-07:890, 07:893, 07:895-07:931, 07:933-07:935, 07:937-07:939, 07:941-07:944									9														23										
*07:01:07, 07:04:03, 07:14, 07:27:02, 07:199:02, 07:526:01, 07:681									9														23								30		
*07:01:41, 07:02:41, 07:415-07:416								9															w										
*07:01:74					6			9					13				17	18					23				27						
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	

Negative Control

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	
*07:01:93, 07:02:08, 07:02:29, 07:02:91, 07:33N, 07:139, 07:152N, 07:156, 07:246:01-07:246:02, 07:267, 07:296, 07:322, 07:405:01, 07:493, 07:517, 07:714, 07:724, 07:936N									9																								
*07:02:09					6				9				13				18					23					27						
*07:07, 07:09, 07:49, 07:76:01-07:76:02, 07:247, 07:315, 07:328, 07:403, 07:406, 07:598									9	11									20			23											
*07:10, 07:28, 07:41, 07:43:01, 07:184, 07:196, 07:367, 07:536, 07:568, 07:857, 07:892									9												21	23											
*07:20:01-07:20:02, 07:263			3						9													23							29				
*07:26:01:01-07:26:03, 07:92, 07:314:01, 07:314:03, 07:317, 07:351, 07:821N, 07:845, 07:870, 07:891, 07:894									9																	26							
*07:31:01, 07:514									9				14				18					23											
*07:31:02									9				14				18					23								30			
*07:39-07:40, 07:335									9											20													
*07:43:02									9												21	23								30			
*07:53, 07:216, 07:932									9	11												23									30		
*07:64, 07:402									9						16							22											
*07:96:01-07:96:02			3						9													22				26			29				
*07:101, 07:148, 07:161		2							9	10												23			26								
*07:123:01-07:123:02, 07:173, 07:294, 07:626, 07:725, 07:775, 07:850									9								17					23											
*07:125:01, 07:356, 07:531						6			9	11			13	14								23											
*07:125:02						6			9	11			13	14				18				23				27							
*07:127:01-07:127:02, 07:390, 07:562			3						9													23											
*07:133:01, 07:242, 07:330:01, 07:748				4					9													23											
*07:154									9									18				23											
*07:177									9				14					18		20			23										
*07:210, 07:238									9	11										20													
*07:293, 07:519									9						16							23											
*07:314:02, 07:578, 07:756			3						9													22				26							
*07:316, 07:338, 07:579	1								9														23										
*07:449									9																		27						
*07:563									9											20			23										
*07:583		2				6			9	10			13					18								26	27						
*07:723						6			9				13													26	27						
*07:940									9			12											23										
*08:01:01:01-08:01:17, 08:01:20-08:01:29, 08:01:31-08:01:32, 08:03:01:01-08:03:04, 08:06, 08:08:01-08:08:02, 08:14, 08:16:01-08:16:02, 08:20, 08:22:01:01-08:22:01:02, 08:24, 08:26N, 08:28, 08:36N, 08:38, 08:40-08:42, 08:44, 08:46, 08:50, 08:56-08:61, 08:65, 08:72:01-08:72:02, 08:78-08:82, 08:84-08:89N, 08:91, 08:96-08:99, 08:101-08:102, 08:106, 08:109, 08:117-08:119, 08:121N-08:122, 08:124, 08:127N-08:131, 08:133, 08:135-08:136, 08:141Q, 08:143-08:145, 08:147-08:149, 08:153-08:155, 08:157, 08:162, 08:164, 08:168, 08:173N-08:176, 08:178, 08:186-08:187, 08:189-08:190, 08:192-08:194, 08:196-08:197, 08:199, 08:203-08:204, 08:208N-08:210, 08:212, 08:217, 08:219-08:221									10								18				21						27						
*08:01:18, 08:02:01:01-08:02:01:34, 08:02:03-08:02:10, 08:02:12-08:02:33, 08:04:01:01-08:04:03, 08:13, 08:15:01-08:15:02, 08:17-08:19:02, 08:23, 08:27, 08:29-08:30, 08:32-08:35, 08:37, 08:39, 08:43, 08:45, 08:49, 08:51-08:55N, 08:62:01-08:63, 08:66-08:71, 08:73, 08:75-08:77, 08:90, 08:92-08:94, 08:100, 08:103, 08:107, 08:110-08:112, 08:114-08:116, 08:120, 08:123, 08:125, 08:132, 08:134, 08:138-08:140, 08:146, 08:150-08:151, 08:156, 08:158-08:161N, 08:165-08:167, 08:169-08:172, 08:177, 08:179-08:182, 08:184, 08:191, 08:195, 08:198, 08:200-08:202, 08:206, 08:211, 08:213-08:216, 08:218, 08:222-08:225Q									10												21						27						
*08:01:19									10								18				21		24				27						
*08:01:30									10								18				21												
*08:02:02, 08:126									10												21		24				27						
*08:02:11:01-08:02:11:02									10												21												
*08:05, 08:25									10				13								21					26	27						
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	

Negative Control

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
*08:07, 08:104, 08:142										10									19		21						27					
*08:09, 08:11, 08:205										10								18									27					
*08:10							7			10	11				15			18			21											
*08:12:01:01-08:12:01:02										10										20	21			24		26	27					
*08:21, 08:137										10			13					18			21					26	27					
*08:31		2								10											21						27					
*08:47, 08:188										10								18	19		21						27					
*08:48, 08:152										10																	27					
*08:74										10						16					21						27					
*08:83								8		10								18									27					
*08:95										10						16		18			21						w					
*08:105			3							10								18			21						27					
*08:108								8		10																	27					
*08:113								8		10							17										27					
*08:163						6				10						16		18			21				25	27						
*08:183						6				10						16					21				25	27						
*08:185		2								10									19		21					27						
*08:207								8		10											21					27						
*12:02:01-12:02:02:15, 12:02:04-12:02:07, 12:02:09-12:02:10, 12:02:13, 12:02:16, 12:02:18, 12:02:20-12:02:22, 12:02:24-12:02:27, 12:02:31, 12:02:33, 12:02:36, 12:02:38-12:02:39, 12:02:41-12:02:42, 12:17-12:18:01, 12:22, 12:30, 12:36, 12:40, 12:56, 12:64, 12:67-12:69, 12:73-12:74, 12:80N, 12:83-12:86, 12:103-12:106, 12:112, 12:114, 12:117, 12:124, 12:126, 12:128, 12:130, 12:134, 12:136-12:137, 12:142, 12:145, 12:148N, 12:151, 12:161, 12:164, 12:166, 12:168, 12:177, 12:179, 12:183, 12:193, 12:196, 12:200, 12:204, 12:207, 12:212, 12:214, 12:219N, 12:226, 12:228, 12:231, 12:234, 12:236N, 12:240-12:241, 12:243, 12:247, 12:250, 12:255, 12:261, 12:263, 12:268, 12:270N, 12:275, 12:279-12:281, 12:285, 12:294, 12:303-12:304, 12:307-12:308, 12:315, 12:327N, 12:329N, 12:335-12:336, 12:338												12	13													26	27			30		
*12:02:03, 12:02:08, 12:02:12, 12:03:47, 12:10:01-12:10:02, 12:20, 12:96, 12:118, 12:153, 12:155Q, 12:162, 12:175, 12:221, 12:251, 12:258, 12:296												12	13														26	27				
*12:02:11								8				12															26				30	
*12:02:14-12:02:15, 12:02:19, 12:02:30, 12:02:37, 12:02:40, 12:18:02, 12:224													13														26	27			30	
*12:02:17, 12:239			3									12	13														26	27			30	
*12:02:23												12	13														26	w			30	
*12:02:28												12	w														26	w			30	
*12:02:29, 12:123												12															26				30	
*12:02:32, 12:02:34													13											24		26	27			30		
*12:02:35								8				12	13														26	27			30	
*12:03:01:01-12:03:01:42, 12:03:03, 12:03:05-12:03:07, 12:03:10-12:03:12, 12:03:14-12:03:15, 12:03:18-12:03:19, 12:03:21-12:03:22, 12:03:25, 12:03:27-12:03:28, 12:03:31, 12:03:33, 12:03:35-12:03:46, 12:03:48-12:03:65, 12:03:67-12:03:69, 12:03:71-12:03:72, 12:06-12:07, 12:11-12:13:01:02, 12:19, 12:23, 12:25, 12:29, 12:32, 12:34-12:35, 12:37-12:39N, 12:42Q-12:43, 12:45-12:48, 12:50-12:53, 12:57:01-12:57:02, 12:61-12:63, 12:65-12:66, 12:70-12:71, 12:75-12:79, 12:82, 12:87-12:91, 12:93-12:95, 12:97, 12:99:01-12:102, 12:107-12:111, 12:115-12:116, 12:120-12:122, 12:125, 12:129, 12:131, 12:133, 12:138-12:139Q, 12:141, 12:143, 12:149-12:150, 12:152, 12:156-12:160, 12:163, 12:165, 12:167, 12:170-12:174, 12:180, 12:182, 12:184, 12:186-12:187, 12:189-12:192, 12:197, 12:199, 12:201-12:202, 12:205-12:206, 12:209-12:211, 12:213, 12:216, 12:218, 12:220, 12:223, 12:225, 12:227, 12:229-12:230, 12:232N, 12:235, 12:237-12:238, 12:242, 12:244-12:246, 12:248-12:249, 12:253-12:254, 12:256-12:257, 12:259-12:260:02, 12:262, 12:264-12:267, 12:271-12:274:01N, 12:276-12:278, 12:283-12:284, 12:286, 12:288-12:293, 12:295N, 12:297, 12:300, 12:302, 12:305-12:306, 12:310-12:314, 12:316-12:324N, 12:326, 12:328, 12:330N-12:332, 12:334, 12:337, 12:339-12:342Q												12	13	14													26	27				
*12:03:02, 12:03:08, 12:03:16, 12:03:20, 12:03:30, 12:26, 12:59, 12:185													13	14													26	27				
*12:03:04												12		14													26					
*12:03:09								8				12		14													26					
*12:03:13												12	w	14													w	27				
*12:03:17												12	13	14										24		26	27					
*12:03:23, 12:140:01-12:140:02			3									12	13	14													26	27				
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
Negative Control																																

Negative Control

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
*12:03:24												12	w	14												26	w					
*12:03:26, 12:03:66								8				12	13	14												26	27					
*12:03:29, 12:03:32												12	13	14												26	27			30		
*12:03:34:01-12:03:34:02													14													26						
*12:03:70												12	13	14												26	w					
*12:04:01											11			14	15											26					31	
*12:04:02:01-12:04:02:02, 12:60											11	12	13	14	15											26					31	
*12:05:01-12:05:02											11	12		14	15																31	
*12:08												12	13				17									26	27			30		
*12:09											11				15					20						26					31	
*12:14:01													13				18	19								26	27			30		
*12:14:02												12	13				18	19								26	27			30		
*12:15, 12:113								8				12	13	14			17									26	27					
*12:16:01, 12:217												12	13									23					27			30		
*12:16:02												12										23								30		
*12:21							7				11	12			15									24							31	
*12:24												12	13							20						26	27					
*12:27												12	13													26				30		
*12:28, 12:58, 12:325												12	13	14			18									26	27					
*12:31:01-12:31:02												12	13	14							21					26	27					
*12:33							7				11			14	15									24							31	
*12:41											11	12	13		15											26				30	31	
*12:44												12	13								21					26	27					
*12:49												12														27						
*12:54											11	12	13	14												26						
*12:55												12	13	14		16											27					
*12:72												12	13		15											26	w			30	31	
*12:81, 12:282												12	13	14			17									26	27					
*12:92												12	13	14				19								26	27					
*12:98, 12:301												12		14													27					
*12:119		2										12	13	14												26	27					
*12:127										10		12	13													26	27			30		
*12:132												12	13					18							25	26	27		30			
*12:135												12	13	14	15		18									26	w				31	
*12:144													13	14							21					26	27					
*12:146											11	12	13		15		18							25	26				30	31		
*12:147, 12:195:02												12	13	14								23					27					
*12:154												12	13	14	15											26	w				31	
*12:169, 12:252												12	13				18									26	27			30		
*12:176												12	13	14			18	19	20							26	27					
*12:178												12	13	14												26	27	28				
*12:181									9				13				18									26	27					
*12:188											11	12					17															
*12:194											11	12	13	14											25	26	27					
*12:195:01												12		14									23									
*12:195:03														14									23									
*12:198		2										12	13													26	27			30		
*12:203										10		12	13	14										24		26	27					
*12:208								8					13				17									26	27					
*12:215	1											12	13	14												26	27					
*12:222, B*67:02:01-01-67:02:01:02 ³													13													26	27					
*12:233													13								21					26	27					
*12:269						6						12	13	14		16										25		27				
*12:287												12	13				17	18							25	26	27		30			
*12:298								8					13													26	27			30		
*12:309								8					13				17				21					26	27					
*12:333												12	13	14			18	20								26	27					
*14:02:01:01-14:02:02, 14:02:04-14:02:07, 14:02:09-14:02:20, 14:02:22-14:02:31, 14:02:33-14:02:39, 14:07N, 14:11, 14:13-14:14, 14:18-14:21N, 14:23-14:24:02, 14:27, 14:29-14:34, 14:36-14:37, 14:39-14:40, 14:42-14:44, 14:46:01-14:48, 14:50-14:52, 14:55-14:57, 14:59-14:60, 14:62, 14:65-14:69, 14:71-14:76, 14:80-14:85, 14:88-14:89, 14:91, 14:93N-14:104, 14:106-14:107, 14:109, 14:111, 14:114-14:115, 14:117N, 14:119-14:120, 14:122-14:125, 14:127-14:128																16				20						26		28				
*14:02:03, 14:78																16				20						26	27	28				
*14:02:08, 14:02:21												12				16				20						26		28				
*14:02:32																16				20												
*14:03:01:01-14:03:01:02, 14:03:03-14:03:04, 14:22, 14:35N, 14:38, 14:41, 14:61, 14:70, 14:79, 14:86, 14:112-14:113, 14:118, 14:126			3													16				20						26	27	28				
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32

Negative Control

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
*14:03:02			3													16				20						26	28					
*14:03:05			3									12				16				20						26	27	28				
*14:04, 14:64											11					16				20						26	28					
*14:05																16				20						26						
*14:06, 14:09, 14:28:01-14:28:02, 14:63, 14:105Q, A*29:86, A*30:62, A*31:187, A*33:63, B*08:249, B*18:64, B*35:183, B*35:252, B*35:502, B*39:114, B*39:140, B*40:138, B*41:59 ³																16																
*14:08																16											27					
*14:10																16		18		20						26	27	28				
*14:12, 14:49, 14:108											11									20						26	28			31		
*14:15, 14:87														14		16		18														
*14:16														14		16																
*14:17																16		18	19	20						26	28					
*14:25		2														16				20		22				26	28					
*14:26								8								16				20						26	28					
*14:45		1													15	16				20						26						
*14:53			3													16											27					
*14:54			3													16					21						27					
*14:58		1														16				20						26	28					
*14:77											11					16																
*14:90			3								11		13			16											27			30		
*14:92																16	17															
*14:110			3													16		18		20						26	28					
*14:116			3													16		18		20						26	27	28				
*14:121									9							16				20						26	28					
*15:02:01:01-15:02:09, 15:02:11-15:02:16, 15:02:18-15:02:47, 15:05:01:01-15:05:07, 15:05:09-15:05:10, 15:05:13-15:05:17, 15:06:01-15:06:03, 15:09:01:01-15:09:02, 15:13:01:01-15:13:02, 15:15, 15:18-15:19, 15:26, 15:28-15:29, 15:32Q-15:35, 15:38-15:39, 15:41, 15:44:01-15:54, 15:56-15:57, 15:59-15:61, 15:64, 15:67-15:70, 15:73, 15:78:01-15:84Q, 15:86-15:101, 15:103, 15:105Q-15:106, 15:108-15:110, 15:112, 15:114-15:115N, 15:117, 15:119-15:123, 15:126-15:137, 15:139-15:142, 15:145N-15:147, 15:150-15:157, 15:159-15:161, 15:163-15:168, 15:170-15:177N, 15:179-15:180, 15:182-15:189N, 15:191, 15:193-15:200, 15:202-15:207, 15:209-15:210, 15:212-15:213N, 15:215-15:217, 15:219-15:221, 15:223-15:224, 15:226-15:233			3					8			11				15		17				21								29			
*15:02:10, 15:02:17			3					8			11				15		17				21	22							29			
*15:03			3					8			11		13		15		17				21					26			29			
*15:04:01:01-15:04:05, 15:24, 15:30, 15:40, 15:66, 15:76-15:77, 15:211, 15:222			3								11				15		17				21								29			
*15:05:11-15:05:12, 15:71			3					8							15		17				21								29			
*15:05:18, 15:22, 15:55, 15:58, 15:72, 15:125, 15:149, 15:225			3					8			11				15		17												29			
*15:07:01:01-15:07:01:02			3					8									17				21						27		29			
*15:08:01-15:08:02, 15:74, 15:143, 15:148, 15:214			3					8			11				15		17				21								29	31		
*15:10:01, 15:10:04, 15:208								8			11				15		17				21											
*15:10:02-15:10:03, 15:17, 15:27, 15:113, 15:124, 15:192			3					8			11				15		17				21											
*15:11			3			6		8			11				15		17				21											
*15:12					5			8			11				15		17	18			21								29			
*15:16			3					8			11		13		15		17				21					26						
*15:21, 15:116			3					8							15		17				21						w		29			
*15:23:01-15:23:02, 15:85, 15:138, 15:158, 15:162, 15:181			3					8			11						17				21								29			
*15:25			3										13					18	19		21					26	27					
*15:31, 15:62, 15:169			3					8			11				15		17	18			21								29			
*15:36			3			6		8			11						17				21								29			
*15:37			3					8			11				15		17											28	29			
*15:42, 15:111		2	3					8			11				15		17				21								29			
*15:43			3					8									17				21	22										
*15:63								8			11						17				21											
*15:65			3								11				15		17	18	19										29			
*15:75			3					8			11		13		15		17	18			21								29	30		
*15:102		1	3								11				15		17			20						26	28	29		31		
*15:104		1	3								11				15		17				21								29			
*15:107							7	8			11				15		17				21											
*15:118			3			6		8			11				15		17				21								29			
*15:144					5			8									17	18			21						27		29			
*15:178			3				7	8			11				15		17				21								29			
*15:190			3								11				15		17												29			
*15:201		2	3					8			11				15		17		19		21								29			
*15:218			3	4				8			11				15		17				21								29			
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32

Negative Control

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32		
*16:01:01:01-16:01:01:23, 16:01:03-16:01:15, 16:01:18-16:01:25, 16:01:27-16:01:41, 16:06-16:08, 16:11, 16:13-16:14, 16:16Q-16:17, 16:21-16:24, 16:26-16:28, 16:30N-16:32, 16:36, 16:38-16:39:02, 16:41, 16:44, 16:49, 16:51-16:52, 16:54, 16:56, 16:58-16:59, 16:65, 16:72-16:73, 16:75-16:76, 16:79-16:81, 16:83, 16:86, 16:92-16:93, 16:95-16:98, 16:100, 16:105, 16:111-16:112, 16:114, 16:116-16:119, 16:122, 16:125, 16:127-16:131, 16:134-16:135, 16:137-16:138, 16:141-16:142, 16:146-16:148, 16:151-16:152, 16:154, 16:157-16:162, 16:164-16:165, 16:168, 16:171-16:175, 16:177-16:178, 16:180, 16:183											11		13				18										27			30				
*16:01:02, 16:01:26, 16:57, 16:85, 16:94, 16:106, 16:126											11		13				18														30			
*16:01:16								8			11		13				18														30			
*16:01:17											11	12	13				18										27				30			
*16:02:01:01-16:02:12, 16:02:14-16:02:19, 16:09, 16:12, 16:19, 16:46-16:47, 16:60, 16:63, 16:69, 16:74, 16:77N, 16:84, 16:89N-16:90, 16:99, 16:101-16:102, 16:104, 16:107-16:108, 16:115, 16:120, 16:123N, 16:132N-16:133, 16:136, 16:143-16:145, 16:153, 16:155-16:156, 16:163, 16:166-16:167, 16:176, 16:179, 16:181, 16:184											11		13		15		18														30			
*16:02:13											11	12	13		15		18														30			
*16:04:01:01-16:04:01:03, 16:04:03-16:04:05, 16:33, 16:42, 16:55, 16:61, 16:66, 16:78, 16:82, 16:124, 16:149, 16:182														14			18		20								27							
*16:10													13				18										27				30			
*16:15:01											11		13				18									26	27				30			
*16:15:02											11	12	13				18									26	27				30			
*16:18											11		13				18										27	28			30			
*16:20								8			11		13				17	18									27				30			
*16:25											11		13		15		18									26					30			
*16:29		2												14			18		20								27							
*16:34			3			6					11		13				18								25		27				30			
*16:35														14			18										27							
*16:37											11		13		15		18										w				30			
*16:40, 16:110							7								16												27							
*16:43											11		13														27				30			
*16:45													13				18			21							27				30			
*16:48											11			14	15		18																	
*16:50		2									11		13				18										27				30			
*16:53							7								16				20							26	27							
*16:62											11		13				18								25		27				30			
*16:64											11		13				17	18								26	27				30			
*16:67													13				18										27							
*16:68																	18		20							26	27							
*16:70			3								11		13		15		17	18												29	30			
*16:71											11		13				18										w				30			
*16:87											11		13				17	18									27				30			
*16:88											11				15		18		20							26						31		
*16:91											11			14	15		18		20													31		
*16:103											11				15		18																	
*16:109								8						14			17	18		20							27							
*16:113							7				11		13		16												27				30			
*16:121			3			6					11		13		15		18															30		
*16:139										10	11		13				18							24			27				30			
*16:140											11				15																			
*16:150						6								14	16		18		20						25		27							
*16:169								8					13				17	18									27				30			
*16:170											11		13		16		18										w				30			
*17:01:01:02-17:01:03, 17:01:05-17:01:09, 17:01:11-17:04, 17:06, 17:08-17:16:02, 17:18-17:19, 17:21, 17:24-17:26, 17:28-17:32, 17:34-17:39, 17:41-17:59		2									11				15			19		21						26								
*17:01:04		2									11		13		15			19		21						26								
*17:01:10		2									11				15			19		21						26				30				
*17:05		2					7				11				15			19		21						26								
*17:07, 17:23, 17:33											11				15			19		21						26								
*17:17											11				15			19		21						26						31		
*17:20		2									11							19		21		23												
*17:22		2																19		21						26								
*17:27N		2									11				15			19		21						w								
*17:40		2									11				15			19		21														
*18:01:01:01-18:02:05, 18:04-18:08, 18:10-18:13											11								20	21		23										31		
*18:03		2									11								20	21	22	23										31		
*18:09											11				15				20	21						26						31		
*18:14											11								20			23										31		
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32		

Negative Control

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	
A*01:110, A*02:636, A*03:210, A*11:135, A*24:137, A*32:96, A*68:160, B*07:234, B*15:125:02, B*15:315, B*18:126, B*39:138, B*44:02:44, B*44:85, B*46:01:33, B*50:56, B*51:226, B*55:30, B*58:74, B*58:107			3																														
A*24:96, A*24:106																																30	
A*24:174, B*57:129											11																					30	
A*25:55, B*07:13, B*07:15, B*07:160																										26							
B*07:02:77, B*08:01:36, B*15:01:50, B*15:436, B*18:03:02, B*35:08:02, B*35:08:05, B*55:01:25																											27						
B*07:239, B*14:03																				20						26							
B*08:145																								23									
B*13:13:02, B*40:01:49-40:01:50																							22										
B*15:96, B*44:274				4																													
B*15:277, B*15:430, B*35:222, B*42:22																15																	
B*15:510, B*35:475																																	
B*15:539, B*39:76, B*54:18	1																																
B*18:01:42, B*57:01:24												12																					
B*46:84														14																	28		
B*54:38																															28		
B*58:02:01-01:58:02:01:03								8																									
B*58:02:02								8													20												
B*73:01:01-01:73:02																						w											
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	

Negative Control