

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42
HLA-C allele^{1,2,3}																																										
*08:01:01:01-08:01:06, 08:01:09 -08:01:20, 08:01:22-08:01:32, 08:22:02, 08:46, 08:50, 08:58, 08:60, 08:79, 08:81, 08:83-08:85, 08:91, 08:96-08:98, 08:105-08:106, 08:109, 08:119, 08:122, 08:124, 08:131, 08:133, 08:136, 08:145, 08:147-08:148, 08:153, 08:155, 08:157, 08:162, 08:164, 08:173N, 08:175-08:176, 08:186-08:187, 08:190, 08:192-08:194, 08:196-08:197, 08:199, 08:203-08:204, 08:208N, 08:210, 08:212, 08:217, 08:219-08:220, 08:228, 08:232, 08:234-08:236N, 08:241, 08:244, 08:249-08:250	1	2										12									21												34								41	
*08:01:07	1	2										12									21							28					34								41	
*08:01:21, 08:163, 08:221	1	2										12									21																				41	
*08:02:01:01 -08:02:05, 08:02:08 -08:02:24, 08:02:26-08:02:29, 08:02:31-08:02:38, 08:45, 08:48, 08:67-08:68, 08:71, 08:76, 08:100, 08:108, 08:111, 08:114, 08:116, 08:120, 08:123, 08:125, 08:132, 08:134, 08:140, 08:146, 08:150-08:151, 08:156, 08:158-08:159, 08:161N, 08:166-08:167, 08:169-08:172, 08:179-08:181N, 08:184, 08:191, 08:195, 08:198, 08:200-08:201, 08:206-08:207, 08:211, 08:214N-08:216, 08:218, 08:222-08:227, 08:229-08:231, 08:237-08:240, 08:242, 08:245, 08:247-08:248			3									12	13																				34									
*08:02:06			3									12	13																	31			34									
*08:02:07			3									12	13															28					34									
*08:02:25, 08:02:30, 08:183			3									12	13																													
*08:03:01:01 -08:03:05, 08:117	1	2		4								12									21												34									
*08:04:01:01-08:04:03, 08:93, 08:139	1		3									12																					34									
*08:05			3		5							12	13				17																34									
*08:06	1	2		4		6						12									21												34									
*08:07			3				7					12	13																				34									
*08:08:01-08:08:02	1	2										12		14							21												34								41	
*08:09	1	2								10	11	12									21				25								34								41	
*08:10	1	2							9												21												34								41	
*08:11 , 08:59	1	2										11	12								21				25								34								41	
*08:12:01:01 -08:12:01:02			3									11	12	13																			34									
*08:13	1		3												15																		34									
*08:14	1	2		4		6						12									21				26								34									
*08:15:01-08:15:02 , 08:51, 08:246					5							12																					33	34								
*08:16:01-08:16:02	1	2												15							21												34								41	
*08:17			3							10		12	13																					34								
*08:18			3										13	14																				34								
*08:19:01			3									12	13														27				31			34								
*08:19:02			3									12	13														27							34								
*08:20 , 08:99 ⁴	1	2										12									21									30			34								41	
*08:21	1	2			5							12					17				21												34								41	
*08:22:01:01 -08:22:01:02, 08:154, 08:178	1	2										12			16						21												34								41	
*08:23			3			6						12	13																					34								
*08:24, 08:88N ⁴	1	2							9			12									21												34								41	
*08:25			3										13	15		17																	34									
*08:26N	1	2										12					18				21												34								41	

Negative Control

[illegible]

[illegible]

[illegible]

[illegible]

[illegible]

[illegible]

[illegible]

[illegible]

[illegible]