

Well No.	1	2	3	4	5	6	7	8
HLA-C allele ^{1,2,3}								
*18:01:01:01-18:01:02	1	2		4				Negative Control
*18:02:01:01-18:02:05, 18:10-18:15	1		3	4				
*18:03	1		3	4	5			
*18:04	1		3	4		6		
*18:05	1	?	?				7	
*18:06	1	?	?	4	5	6		
*18:07N	1	?	?	4			7	
*18:08-18:09	1	?	?	4				
*02:22, 02:118, 05:08, 05:52, 05:59, 05:106:01-05:106:02, 06:09:01-06:09:02, 06:144, 08:27, 08:29, 08:31, 08:160, 12:31:01-12:31:02, 12:144, B*15:137, B*15:488, B*37:66, B*50:47					5			
*02:197, 03:540, 04:01:01:01-04:01:01:29, 04:01:01:31-04:01:143, 04:03:01:01-04:20, 04:23-04:93N, 04:95N-04:111, 04:113-04:122, 04:124-04:168, 04:170N-04:257, 04:268, 04:274-04:359, 04:361-04:396N, 04:398-04:421, 04:423-04:430, 04:433-04:441, 04:443-04:459, 04:461-04:468, 05:192, 12:233, 12:348			3					
*03:08, 03:29, 03:31, 03:246, 03:555, 14:73, 15:77							7	
*04:94:01-04:94:02			3		5			
*04:112, 04:123N, 04:169			3				7	
*04:397, 04:432, 04:442	1		3			6		
*06:02:08, 06:34:02				4				
*06:17	1							
*07:01:01:01-07:01:02:14, 07:01:04-07:06:05, 07:08, 07:10-07:19, 07:21-07:33N, 07:35-07:48, 07:50-07:63, 07:65-07:72, 07:74-07:75, 07:77-07:91, 07:93:01-07:95, 07:97-07:141:02, 07:142-07:171, 07:173-07:209, 07:211-07:237, 07:239-07:246:02, 07:248-07:294, 07:296-07:314:03, 07:316-07:327, 07:329N-07:347N, 07:349-07:389, 07:391-07:402, 07:404-07:405:02, 07:407-07:577, 07:579-07:582Q, 07:584-07:597, 07:599-07:605, 07:607-07:617, 07:619-07:722, 07:724-07:856N, 07:858-07:890, 07:892-07:952, 07:954N-07:1001N, 12:181		2						
*07:07, 07:09, 07:49, 07:76:01-07:76:02, 07:210, 07:238, 07:247, 07:315, 07:328, 07:403, 07:406, 07:598	1	2						
*07:20:01-07:20:02, 07:64, 07:73:01-07:73:02, 07:92, 07:96:01-07:96:02, 07:172:01-07:172:02, 07:390, 07:578, 07:583, 07:606, 07:723, 07:891		2				6		
*07:141:03, 07:857		2					7	
*07:953		2		4				
*12:349, 14:121						6		
Well No.	1	2	3	4	5	6	7	8