

KIR Genotyping 2021 4N4 R01 Interpretation Table.xlsx

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30
KIR allele ^{1, 2, 3}																														
2DL1*0010101-068	1																													
2DL2*0010101-00304, 005-015		2																												
2DL2*004		2																									27			
2DL3*0010101-0020103, 004 - 011, 013-015, 017-026, 028-036			3																									28		
2DL3*003, 016			3																											
2DL3*01201-01202, 027																												28		
2DL4*00101-00604, 010, 01201-01202, 014-016, 018, 021-026, 028-034, 036-042, 045, 049, 050, 054, 056, 058				4																				24						
2DL4*007-00902, 0110101-01102, 013, 017, 019, 020, 027, 035, 043-044, 046-048, 051-053, 055, 057				4																										
2DL5A*0010101-0010902, 0050101, 0050103-0050105, 01201-01202, 0210101-0210102, 0230101-0230102, 025-0340102					5	6																				26				
2DL5A*0050102					5	6																								
2DL5A*014-015, 022, 024, 040					5	6																				?				
2DL5B*0020101, 0020103, 0020107-0020111, 00205, 004, 00601, 0070101, 0080101, 00802, 010-011, 035-0390102					5		7																		25					
2DL5B*0020102, 0020104, 00203-00204, 00603, 0080102, 00803-009, 001301-01302, 01304, 016-020, 041					5		7																		?					
2DL5B*0020106, 0020201-003, 00602, 01303					5		7																							
2DS1*001								8													21								29	
2DS1*0020101-00202, 004, 006, 008-0120102								8													22								29	
2DS1*00301-00302, 00501-00502								8													22									
2DS1*013	1																												29	
2DS2*0010101-0210103									9																					
2DS3*00101-024										10																				
2DS4*0010101-00106, 01101-01102, 014- 017, 019-022											11																			
2DS4*0030101-0040102, 0060101-010, 012, 013, 018												12																		
2DS5*001, 00202, 003-00501, 00801-009, 011-017, 038													13													?				
2DS5*0020102-0020133, 00203-0020702, 018-037													13													26				
2DS5*00502-0070102, 010													13																	
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30

Negative Control

KIR Genotyping 2021 4N4 R01 Interpretation Table.xlsx

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	
3DL1*0010101-0020105, 0050101-0090104, 0150101-018, 0200101-0200103, 022-03502, 038, 040-044, 051-054, 057, 05901-071, 073-077, 079-081N, 086-090, 092-103, 109, 111-116, 118-1190102														14																	
3DL1*0040101-00404, 019, 021, 036-037, 039, 056, 072, 091, 110, 117														14									23								
3DL2*0010101-116															15																
3DL3*0010101-023, 02501-114																16															
3DS1*010-014, 045-050, 055, 058, 082-085, 104-108																	17														
3DS1*078														14			17														
2DP1*00101-004, 007-008, 011-013, 015-024																	18														
2DP1*006, 009, 010																	18													29	
3DP1*001, 0090101-0090103																		19				22			25						
3DP1*0030101-0030402, 005-00602, 008, 0100101-0100102, 01003-01005, 0140101-052																			19	20											
3DP1*004																		19			22				26						
3DP1*002, 007, 00902, 01002, 011-012																		19			22				?						
3DP1*013																			19	20					?						
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	

Negative Control