

KIR Genotyping 2021 9L3 R01.xlsx

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26
KIR allele ^{1, 2, 3}																										
2DL1*0010101-038	1																									
2DL2*0010101-00304, 005-015		2																								
2DL2*004		2																								
2DL3*0010101-011, 013-026, 028-036			3																							
2DL3*01201-01202, 027		w																								
2DL4*00101-00604, 010, 01201-01202, 014-016, 018, 021-026, 028-034, 036-042, 045, 049-050, 054				4																				24		
2DL4*007-00902, 0110101-01102, 013, 017, 019, 020, 027, 035, 043-044, 046-048, 051-053				4																						
2DL5A*0010101-00105, 0050101, 0050103-0050104, 01201-01202					5	6																				26
2DL5A*0050102					5	6																				
2DL5A*014-015, 021-024					5	6																				?
2DL5B*0020101, 0020103, 0020107-0020108, 004, 00601, 0070101, 0080101, 00802, 010-011					5		7																		25	
2DL5B*0020102, 0020104, 0020106, 00203-00204, 00603, 0080102, 00803-009, 01301-01302, 01304, 016-020					5		7																		?	
2DL5B*0020106, 0020201-0020202, 003, 00602, 0070102, 01303					5		7																			
2DS1*001								8													21					
2DS1*0020101-009								8														22				
2DS2*0010101-009									9																	
2DS3*00101-008										10																
2DS4*0010101-00106, 01101-01102, 014-017, 019-020											11															
2DS4*0030101-0030104, 0040101-0040102, 0060101-0060102, 007-010, 012-013, 018												12														
2DS5*001, 00202-00502, 00801-009, 011-017													13													?
2DS5*0020101-0020104													13													26
2DS5*006, 007, 010													13													
3DL1*0010101-0020105, 0050101-0090104, 0150101-018, 020010-0200103, 022-03502, 038, 040-044, 051-054, 057, 05901-077, 079-081N, 086-090, 092-103, 109, 111-116, 118-1190102														14												
3DL1*00401-00404, 019, 021, 036-037, 039, 056, 072, 091, 110, 117														14									23			
3DL2*0010101-115															15											
3DL3*00101-094																16										
3DS1*010-014, 045-049N, 050, 055, 058, 082-085, 104-108																	17									
3DS1*078														14			17									
2DP1*00101-004, 006-024																		18								
3DP1*001, 007, 0090101-00902																			19			22			25	
3DP1*0030101-0030402, 005-006, 008, 01001																				19	20					
3DP1*004																				19		22				26
3DP1*002, 01002-012																				19		22			?	
3DP1*013-015																				19	20				?	

KIR Genotyping 2021 9L3 R01.xlsx

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26
----------	---	---	---	---	---	---	---	---	---	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----

