

INTERPRETATION TABLE					
KIR HLA Ligand					
	Well				
	1	2	3	4	5
Length of spec.	340	340	350	350	370
PCR product					
Length of int.	800	800	800	1070	1070
pos. control ¹					
5'-primer(s) ²	1 st I	1 st I	1 st I	1 st I	1 st I
	5'-CgA ^{3'}	5'-CgA ^{3'}	5'-CAG ^{3'}	5'-CAG ^{3'}	5'-gCA ^{3'}
3'-primer ³	302	302	309	309	317
	5'-ggC ^{3'}	5'-ggT ^{3'}	5'-gTg ^{3'}	5'-ATC ^{3'}	5'-ggA ^{3'}
Well No.	1	2	3	4	5
HLA allele					
Cw*010201-0113, 0115-0133, 0227, 030201-0306, 0308, 0309, 031101-0314, 0316-0328, 0330-0344, 0346-0376, 0411, 0429, 0436, 0611, 070101-0706, 0708, 0710-0733N, 0735-0748, 0750-0775, 0777-0791, 080101-0809, 0811-0829, 120201-120308, 1206-1208, 1210-1220, 1222-1226, 1228, 140201-1403, 1405-1411, 1413, 1414, 1507, 160101-160103, 160401, 1606-1608, 1610, 1611, 1613-1617	1				
Cw*0212, 0310, 0329, 1521	W				
Cw*0114, 020201-020203, 020205-0211, 0213-0226, 0228-0231, 0307, 0315, 0345, 04010101-040112, 0403-0410, 0412-0420, 0423-0428, 0430-0435, 0437-0448, 05010101-050110, 0503-0530, 06020101, 06020102, 060203-0610, 0612-0628, 0707, 0709, 0749, 0776, 0810, 120401-1205, 1209, 1221, 1404, 1412, 150201-1506, 1508-1513, 1515-1520, 1522, 1523, 1602, 1609, 1612, 17010101-1705, 1801-1803	2				
B*0802, 130101-1304, 1306-1308Q, 1310-1323, 1325-1330, 1536, 1589, 9515, 1809, 2701, 3710, 380201-3804, 3808, 3815, 3818, 4047, 4096, 44020101-4405, 4407, 4408, 4410, 4412-4417, 4419N-4424, 4426-4445, 4447-4449, 4451-4474, 4476-4488, 4704, 4902, 5154, 5178, 5309, 5311-5313, 5607			3		
Well No.	1	2	3	4	5
Length of spec.	340	340	350	350	370
PCR product					
Well No.	1	2	3	4	5

B*0736, 0738, 0781, 0803, 1513, 1516-151702, 1523, 1524, 1567, 1587, 1595, 9557, 9562, 9568, 9577, 2702, 2730, 2753, 2757, 380101-380104, 3805-3807, 3809-3814, 3816, 3819, 3820, 4013, 4019, 4406, 4418, 4425, 4450, 4818, 490101, 490102, 4903-4907, 510101-512402, 5126-5146, 5148-5153, 5155-5177, 5179-5184, 520101-5216, 530101-5302, 5304-5308, 5310, 5314-5319, 5412, 5621, 570101-5711, 5713-5726, 580101-5802, 5804-5816, 5818-5824, 5901-5905	4				
A*0281, 0287, 9212, 9224, 9229, 9236, 230101-2322, 24020101-240302, 2405-2411N, 241301-2415, 2417, 2418, 2420-2425, 2427, 2429-2443, 2445N-2464, 2466-2488, 2490N-2499, 250101-2510, 2913, 3107, 3108, 3110, 320101-3220, 6836					5
HLA allele					
Well No.	1	2	3	4	5