

Well No.	1	2	3	4	5	6	7
HLA allele ^{1,2}							
C*01:02:01-01:13, 01:15-01:45, 01:47-01:107, 02:27:01-02:27:02, 02:65, 02:87, 03:02:01-03:06:02, 03:08-03:09, 03:11:01-03:11:02, 03:13:01-03:14, 03:16-03:28, 03:30-03:44, 03:46-03:162, 03:164-03:267, 03:269-03:277N, 04:11, 04:29, 04:36, 04:55, 04:114, 04:172, 05:20, 06:11, 06:82, 06:147, 07:01:01-01-07:06, 07:08, 07:10-07:33N, 07:35-07:75, 07:77-07:294, 07:296-07:314, 07:316-07:327, 07:329N-07:405, 07:407-07:409, 08:01:01-08:09, 08:11-08:63, 08:65-08:115, 12:02:01-12:03:34, 12:06-12:08, 12:10:01-12:20, 12:22-12:32, 12:34-12:40, 12:42Q-12:59, 12:61-12:71, 12:73-12:134, 12:136-12:145, 12:147-12:148N, 14:02:01-14:11, 14:13-14:48, 14:50-14:69, 15:07, 15:25, 15:43, 15:85, 16:01:01-16:01:19, 16:04:01, 16:04:03, 16:06-16:08, 16:10-16:11, 16:13-16:18, 16:20-16:24, 16:26-16:36, 16:38-16:45, 16:49-16:59, 16:61-16:62, 16:64-16:68, 16:71-16:73, 16:75-16:76, 16:78-16:79, 17:22	1						
C*01:14, 02:02:01-02:02:03, 02:02:05-02:26:03, 02:28-02:40:02, 02:42-02:64, 02:66-02:86, 02:88-02:95, 03:07, 03:10, 03:15, 03:29, 03:45, 03:163, 03:268, 04:01:01-01-04:01:66, 04:03:01-04:10, 04:12-04:20, 04:23-04:28, 04:30-04:35, 04:37-04:54, 04:56-04:113, 04:115N-04:171, 04:173N-04:194, 05:01:01-01-05:01:31, 05:03-05:19, 05:21-05:114, 06:02:01-01-06:02:01:03, 06:02:03-06:10, 06:12-06:81, 06:83-06:146, 06:148-06:149, 07:07, 07:09, 07:76:01-07:76:02, 07:315, 07:328, 07:406, 08:10, 12:04:01-12:05, 12:09, 12:21, 12:33, 12:41, 12:60, 12:72, 12:135, 12:146, 14:12, 14:49, 15:02:01-01-15:06:03, 15:08-15:13, 15:15-15:19, 15:21-15:24, 15:26-15:42, 15:44-15:84Q, 15:86-15:106, 16:02:01-16:02:13, 16:09, 16:12, 16:19, 16:25, 16:37, 16:46-16:48, 16:60, 16:63, 16:69-16:70, 16:74, 16:77N, 17:01:01-01-17:21, 17:23-17:28, 18:01-18:09	2						
B*07:27, 07:236, 08:126, 15:43, 18:54, 27:03-27:07:04, 27:09-27:11, 27:13-27:17, 27:19-27:21, 27:23-27:25, 27:27-27:29, 27:31-27:32, 27:34-27:39, 27:41, 27:43, 27:45-27:48, 27:50:01-27:52, 27:54-27:56, 27:58-27:61, 27:63-27:64N, 27:66N-27:74, 27:76, 27:78-27:82, 27:84-27:88, 27:90:01-27:94N, 27:96:01-27:101, 27:103-27:118, 27:120-27:125, 27:127-27:133, 27:135, 37:01:01, 37:01:03-37:04:02, 37:06:01-37:09, 37:12-37:13, 37:15-37:33N, 37:35-37:36, 37:38-37:54, 38:17, 40:188, 47:01:01-01-47:01:02, 47:05-47:09, 53:03						6	Negative Control
B*07:36, 07:38, 07:81, 07:180, 07:219, 08:03, 08:52, 08:78, 15:13:01-15:13:02, 15:16:01-15:17:03, 15:23-15:24:02, 15:67, 15:87, 15:95, 15:157, 15:162, 15:168, 15:177, 15:196, 15:208, 15:216, 15:222, 15:230, 15:254, 15:268, 15:273, 18:67, 27:02:01-27:02:02, 27:30, 27:53, 27:57, 27:62, 27:65N, 27:75, 27:77, 27:83, 27:95, 27:102, 27:119, 27:126, 27:134, 37:34, 38:01:01-38:01:07, 38:05-38:07, 38:09-38:14, 38:16, 38:19-38:22, 38:24-38:28, 38:30-38:34N, 38:36-38:42, 38:51-38:56, 40:13, 40:19, 40:109, 40:117, 40:292, 44:06, 44:18, 44:25, 44:50:01, 44:95, 48:18, 49:01:01-49:01:07, 49:03-49:14, 49:16-49:35, 51:01:01-01-51:01:39, 51:01:41-51:24:05, 51:26-51:46, 51:48-51:53, 51:55-51:77, 51:79-51:185, 52:01:01-01-52:19, 52:21-52:44, 53:01:01-53:02, 53:04-53:08:02, 53:10, 53:14-53:30, 53:32-53:35, 53:37, 53:39, 54:12, 56:21, 57:01:01-57:11, 57:13-57:74, 58:01:01-58:02, 58:04-58:16:02, 58:18-58:29, 58:31N-58:63, 58:65-58:69, 59:01:01-01-59:07				4			
B*07:149, 08:02, 08:117, 13:01:01-13:04, 13:06-13:08, 13:10-13:12:01, 13:13-13:23, 13:25-13:38, 13:40-13:66, 13:68-13:83, 15:36, 15:89, 15:115, 15:256, 15:339, 18:09, 27:01, 37:10, 38:02:01-38:02:06, 38:03-38:04, 38:08, 38:15, 38:18, 38:23, 38:29, 38:35, 38:43-38:50, 40:47, 40:96, 40:110, 40:157, 40:201, 44:02:01-01-44:02:18, 44:02:20-44:05:02, 44:05:04, 44:07-44:08, 44:10, 44:12-44:17, 44:19N-44:24, 44:26-44:45, 44:47-44:49, 44:51-44:74, 44:76-44:89, 44:91-44:94, 44:96-44:128, 44:130, 44:132-44:151, 44:153-44:169, 44:171N-44:195N, 44:197-44:214, 47:04, 49:02, 51:54, 51:78:01-51:78:02, 52:20, 53:09, 53:11-53:13, 53:31, 53:36, 53:38, 56:07, 59:08			3				
A*01:95, 02:81, 02:87, 02:112, 02:124, 02:129, 02:136, 03:152, 11:178, 11:190, 23:01:01-23:69, 24:02:01-01-24:03:02, 24:05-24:11N, 24:13:01-24:15, 24:17-24:18, 24:20-24:27, 24:29-24:43, 24:45N-24:64, 24:66-24:88, 24:90N-24:99, 24:101-24:108, 24:110-24:128, 24:130-24:210, 24:212-24:240N, 24:242-24:289, 24:291-24:293, 25:01:01-25:30, 29:13, 31:07-31:08, 31:10, 32:01:01-32:69, 68:36					5		
Well No.	1	2	3	4	5	6	7