

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

Olerup SSP® HLA-A-B-DR-DQ SSP Combi Tray

Product number:	101.708-24 – including <i>Taq</i> pol. 101.708-24u – without <i>Taq</i> pol.
Lot number:	43R
Expiry date:	2015-March-01
Number of tests:	24 tests
Number of wells per test:	95 +1
Storage - pre-aliquoted primers:	dark at -20°C
- PCR Master Mix:	-20°C
- Adhesive PCR seals	RT
- Product Insert	RT

This Product Description is only valid for Lot No. 43R.

**CHANGES COMPARED TO THE PREVIOUS OLERUP SSP®
HLA-A-B-DR-DQ SSP COMBI TRAY LOT (96N)**

The Lot-specific information for HLA-A-B-DR-DQ Combi Tray including and without *Taq* polymerase is now described in one common Product Insert.

The **HLA-A low resolution** specificity and interpretation tables have been updated for the HLA-A alleles described since the previous *Olerup SSP®* HLA-A-B-DR-DQ SSP Combi Tray lot was made (**Lot No. 96N**).

The **A low resolution** primer set is unchanged compared to the previous lot.

The **HLA-B low resolution** specificity and interpretation tables have been updated for the HLA-B alleles described since the previous *Olerup SSP®* HLA-A-B-DR-DQ SSP Combi Tray lot was made (**Lot No. 96N**).

The primers of the wells detailed below have been exchanged, modified or added compared to the previous lot.

Well	5'-primer	3'-primer	rationale
47	-	-	Exchanged positive control primer pair.
49	Modified	-	Exchanged positive control primer pair.

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: 43R**Lot-specific information**

The **HLA-DR low resolution** specificity and interpretation tables have been updated for the HLA-DRB alleles described since the previous *Olerup* SSP® HLA-A-B-DR-DQ SSP Combi Tray lot was made (**Lot No. 96N**).

The primers of the wells detailed below have been exchanged, modified or added compared to the previous lot.

Well	5'-primer	3'-primer	rationale
79	-	Added	Improved yield of HLA-specific PCR product.

The **DQ low resolution** specificity and interpretation tables have been updated for the HLA-DQB1 alleles described since the previous *Olerup* SSP® HLA-A-B-DR-DQ SSP Combi Tray lot was made (**Lot No. 96N**).

The primers of the wells detailed below have been exchanged, modified or added compared to the previous lot.

Well	5'-primer	3'-primer	rationale
95	Exchanged	-	Improved yield of HLA-specific PCR product.

Changes in revision R01 compared to R00:

1. Primer mix 39 does not amplify the B*15:101 allele. This has been corrected in the Specificity and Interpretation tables.

Change in revision R02 compared to R01:

1. The HLA-A*03:01:03, 03:09, 03:23:01, 03:89, 03:108, 11:06, 11:18, 11:121, 24:28, 24:89, 26:03:01-26:03:02, 26:06, 26:07:02, 26:21, 26:30, 30:13, 30:16, 30:44, 30:46, 68:05, 68:15, 68:20 and 68:30 alleles are weakly amplified by primer mix 2. This has been corrected in the Specificity and Interpretation Tables.

Changes in revision R03 compared to R02:

1. Primer mix 39 does not amplify the B*44:10 allele. This has been corrected in the Specificity and Interpretation Tables.

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

Well 96 contains Negative Control primer pairs, that will amplify more than 95% of the *Olerup* SSP® HLA Class I, DRB, DQB1 and DPB1 amplicons as well as amplicons generated by a control primer pair.

PCR product sizes range from 75 to 430 base pairs.
The PCR product generated by the control primer pair is 430 base pairs.

Length of PCR product	105	200	105	80	75	80
5'-primer¹	164	340	440	45	45	43
	5'-CAC ³	5'-Agg ³	5'-TTA ³	5'-Tgg ³	5'-Tgg ³	5'-Tgg ³
3'-primer²	231	2nd I	507	59	58	57
	5'-TgC ³	5'-AAA ³	5'-TTg ³	5'-CTC ³	5'-ggC ³	5'-CTC ³
A*	+	+	+			
B*	+	+	+			
C*	+	+	+			
DRB1				+	+	
DRB3				+	+	
DRB5				+		
DQB1					+	
DPB1						+

¹The nucleotide position for HLA class I genes and the codon for HLA class II genes, in the 2nd or 3rd exon, matching the specificity-determining 3'-end of the primer is given. Nucleotide and codon numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

²The nucleotide position for HLA class I genes and the codon for HLA class II genes, in the 2nd or 3rd exon or the 2nd intron, matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Nucleotide and codon numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

101.708-24 – including *Taq* polymerase, IFU-01
101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

PRODUCT DESCRIPTION**HLA-A-B-DR-DQ SSP Combi Tray****CONTENT**

The primer set contains 5'- and 3'-primers for grouping the HLA-A*01:01 to A*80:02 alleles into the corresponding serological groups A1 to A80.

The primer set contains 5'- and 3'-primers for grouping the B*07:02 to B*83:01 alleles into the corresponding serological groups B7 to B81 as well as primer pairs for recognizing the Bw4 and Bw6 sequence motifs.

The primer set contains 5'- and 3'-primers for grouping the DRB1*01:01 to DRB1*10:04 alleles into the corresponding serological groups DR1 to DR18 as well as primer pairs for recognizing the DRB3, DRB4 and DRB5 groups of alleles.

The primer set contains 5'- and 3'-primers for grouping the DQB1 alleles into the serological groups DQ2 to DQ9.

Please note that DQB1 amplifications usually are somewhat less pronounced than e.g. DRB and DQA1 amplifications even when using the same DNA preparation and exactly the same experimental procedures.

PLATE LAYOUT

Each test consists of 96 PCR reactions in a 96 well PCR plate.

1	2	3	4	5	6	7	8
9	10	11	12	13	14	15	16
17	18	19	20	21	22	23	24
25	26	27	28	29	30	31	32
33	34	35	36	37	38	39	40
41	42	43	44	45	46	47	48
49	50	51	52	53	54	55	56
57	58	59	60	61	62	63	64
65	66	67	68	69	70	71	72
73	74	75	76	77	78	79	80
81	82	83	84	85	86	87	88
89	90	91	92	93	94	95	ctrl

Wells 1 to 21 – HLA-A low resolution primers.

Wells 22 to 64 – HLA-B low resolution primers.

Wells 65 to 87 – HLA-DR low resolution primers.

Wells 88 to 95 – HLA-DQ low resolution primers.

Well 96 – Negative Control.

The 96 well PCR plate is marked with 'A-B-DR-DQ' in silver/gray ink.

Well No. 1 is marked with the Lot No. '43R'.

A faint row of numbers is seen between wells 1 and 2 or wells 7 and 8 of the PCR trays. These stem from the manufacture of the trays, and should be disregarded.

The PCR plates are covered with a PCR-compatible foil.

101.708-24 – including *Taq* polymerase, IFU-01
101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

INTERPRETATION

Only HLA-A alleles will be amplified by the 21 wells of the HLA-A low resolution primer set, **wells 1 to 21**, except that primer mix 5 will amplify the B*18:27 allele. Thus, the interpretation of HLA-A low resolution is only influenced by this HLA-B allele and not by other HLA Class I genes.

Only HLA-B alleles will be amplified by the 43 wells of the HLA-B low resolution, primer set, **wells 22 to 64**, except that the HLA-C*03:05, 03:25, 03:27 and 03:143 alleles will be amplified by primer mix 24, the C*01:30 allele will be amplified by primer mix 25, the A*23:31, A*24:106, C*07:231 and C*16:10 alleles will be amplified by primer mix 27, the C*07:02:30 and C*08:16:02 alleles will be amplified by primer mix 28, the C*07:46 allele will be amplified by primer mix 29, the A*24:174 allele will be amplified by primer mix 33, the C*03:102 allele will be amplified by primer mix 37, the C*15:51 allele will be amplified by primer mix 41, the C*03:129 allele will be amplified by primer mix 45, the C*15:25 allele will be amplified by primer mix 48, the C*15:39 allele will be amplified by primer mix 49, the C*15:02:04 allele will be amplified by primer mix 58, the C*03:12, 03:19 and 03:102 alleles will be amplified by primer mix 59, the C*06:72 allele will be amplified by primer mix 60 and the A*26:68, A*68:56, C*02:56, C*06:20 and C*12:50 alleles will be amplified by primer mix 62. Thus, the interpretation of HLA-B low resolution typings is only influenced by these HLA-A and HLA-C alleles and not by other HLA class I genes.

Only HLA-DRB alleles will be amplified by the 23 wells the DR low resolution primer set, **wells 65 to 87**. Thus, the interpretation of DR low resolution typings is not influenced by other HLA class II genes.

Only the DQB1 alleles will be amplified by the 8 wells the DQ low resolution primer set, **wells 88 to 95**. Thus, the interpretation of DQ low resolution typings is not influenced the DQB2 and DQB3 genes.

UNIQUELY IDENTIFIED ALLELES

All the HLA-A alleles, i.e. **A*01:01 to A*80:02**, recognized by the HLA Nomenclature Committee in April 2012¹ will be amplified by the primers in the HLA-A low resolution primer set, **wells 1 to 21**. The HLA-A alleles will be grouped into their corresponding serological specificities².

All the HLA-B alleles, i.e. **B*07:02 to B*83:01**, recognized by the HLA Nomenclature Committee in April 2012¹ will be amplified by the primers in the HLA-B low resolution primer set, **wells 22 to 64**. The HLA-B alleles will be grouped into their corresponding serological specificities³.

All the HLA-DRB1, -DRB3, -DRB4⁴ and -DRB5 alleles, i.e. **DRB1*01:01:01 to 10:04, DRB3*01:01:02:01 to DRB3*03:03, DRB4*01:01:01:01 to DRB4*01:08, DRB5*01:01:01 to DRB5*02:05**, recognized by the HLA Nomenclature Committee in April 2012¹ will be amplified by the primers in the DR low resolution primer set,

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

wells 65 to 87. The HLA-DRB alleles will be grouped into their corresponding serological specificities⁵.

All the DQB1 alleles, i.e. **DQB1*05:01:01:01 to 05:14**, **DQB1*06:01:01 to 06:49**, **DQB1*02:01:01 to 02:06**, **DQB1*03:01:01:01 to 03:40** and **DQB1*04:01:01 to 04:08**, recognized by the HLA Nomenclature Committee in April 2012¹ will be amplified by the primers in the DQ low resolution primer set. The DQB1 alleles will be grouped into their corresponding serological specificities⁶, i.e.:

DQ5(1) =	DQB1*05:01:01-05:05
DQ6(1) =	DQB1*06:01:01-06:44
DQ2 =	DQB1*02:01:01-02:05
DQ3 =	DQB1*03:06, 03:10, 03:14
DQ7(3) =	DQB1*03:01:01-03:01:06, 03:04, 03:09, 03:13, 03:16, 03:19
DQ8(3) =	DQB1*03:02:01-03:02:05, 03:05:01-03:05:04, 03:07-03:08, 03:11, 03:18
DQ9(3) =	DQB1*03:03:02:01-03:03:04, 03:12, 03:15, 03:17, 03:20
DQ4 =	DQB1*04:01-04:02

¹HLA-A, HLA-B, HLA-DRB and HLA-DQB1 alleles listed on the IMGT/HLA web page 2012-April-12, release 3.8.0, www.ebi.ac.uk/imgt/hla.

²The A*36 alleles will give rise to identical amplification patterns as some A*01 alleles. These alleles can be separated by the A*01 and A*36 high resolution SSP primer sets.

The A*03:01:23, 03:08, 03:32, 03:36N, 03:57, 03:59, 03:72, 03:89, 03:107-03:108, 03:111 and 03:142 and the A*24:92 alleles will give rise to identical amplification patterns. These alleles can be separated by the respective high resolution SSP primer sets.

The A*23:14 and the A*24:05, 24:13:02, 24:24 and 24:199 alleles will give rise to identical amplification patterns. These alleles can be separated by the respective high resolution SSP primer sets.

³The B*08:26, 08:50, 08:62 and 08:85 and the B*42:07 alleles give rise to identical amplification patterns with the HLA-B low resolution primer set. These alleles can be separated by the respective high resolution primer sets.

The B*13:04, 13:10, 13:21 and 13:35 and the B*44:135 alleles give rise to identical amplification patterns with the HLA-B low resolution primer set. These alleles can be separated by the respective high resolution primer sets.

The B*14:08 and the B*39:25N, 39:30, 39:32-39:34, 39:43, 39:47 and 39:50 alleles give rise to identical amplification patterns with the HLA-B low resolution primer set. These alleles can be separated by the respective high resolution primer sets.

The B*18:29 and the B*35:09:01-35:09:03, 35:18, 35:31-35:32:02, 35:37, 35:53N, 35:64, 35:68:01-35:68:02, 35:75, 35:88, 35:99, 35:118-35:119, 35:127, 35:151 and 35:174 alleles give rise to identical amplification patterns with the HLA-B low resolution primer set. These alleles can be separated by the respective high resolution primer sets.

The B*41:09 and the B*45:02-45:03 alleles give rise to identical amplification patterns with the HLA-B low resolution primer set. These alleles can be separated by the respective high resolution primer sets.

The B*51:104 and 51:118N and the B*58:08:01-58:08:02 alleles give rise to identical amplification patterns with the HLA-B low resolution primer set. These alleles can be separated by the respective high resolution primer sets.

The B*54:01:02 and the B*55:01:07, 55:02:01-55:02:06, 55:07, 55:10, 55:12, 55:16, 55:19, 55:26, 55:30, 55:35, 55:37, 55:39, 55:41-55:43, 55:47-55:48, 55:50 and the B*56:10 alleles give rise to identical amplification patterns with the HLA-B low resolution primer set. These alleles can be separated by the respective high resolution primer sets.

The B*55:04, 55:08, 55:13, 55:23, 55:27, 55:32, 55:46 and 55:49 and the B*56:15, 56:18-56:19N, 56:22 and 56:31-56:32 alleles give rise to identical amplification patterns with the HLA-B low resolution primer set. These alleles can be separated by the respective high resolution primer sets.

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: 43R**Lot-specific information**

The B*57:01:01-57:01:04, 57:01:06-57:03:02, 57:05-57:06, 57:08, 57:10, 57:15-57:20, 57:22-57:23, 57:25-57:30, 57:32-57:44, 57:46, 57:48-57:50 and 57:52-57:57 and the B*58:36 alleles give rise to identical amplification patterns with the HLA-B low resolution primer set. These alleles can be separated by the respective high resolution primer sets.

⁴The DRB4*02:01N and DRB4*03:01N null alleles will not be amplified by the DR low resolution primer set.

⁵The DRB1*08:09 and the DRB1*14:15 alleles yield identical amplification patterns with the DR low resolution primer set. These alleles can be separated by the respective high resolution primer sets. The DRB1*08:20 and the DRB1*13:18, 13:47 and 13:55 alleles yield identical amplification patterns with the DR low resolution primer set. These alleles can be separated by the respective high resolution primer sets.

The DRB1*08:31, 08:41 and DRB1*11:67 alleles yield identical amplification patterns with the DR low resolution primer set. These alleles can be separated by the respective high resolution primer sets. The DRB1*13:13 and 13:119 and the DRB1*14:84 and 14:116 alleles yield identical amplification patterns with the DR low resolution primer set. These alleles can be separated by the respective high resolution primer sets.

⁶The serological split of the DQB1*05:05-05:14, DQB1*06:06 to 06:07 alleles, the DQB1*06:10, 06:13, 06:15-06:24 and 06:27 to 06:49, the DQB1*02:04-02:06, the DQB1*03:07-03:09 and 03:11-03:40 alleles and the DQB1:04:0301-04:08 alleles is not known. In this table we have used the expert-assigned serological grouping in Tissue Antigens (2009) 73:95-170, and also inferred the serological grouping from the naming of the sequence-defined allele.

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

SPECIFICITY TABLE

HLA-A low resolution primer set

Specificities and sizes of the PCR products of the 21 primer mixes used for
 HLA-A low resolution SSP typing

Primer Mix	Size of spec. PCR product ¹	Size of control band ²	HLA-A serology ³	Amplified HLA-A alleles ^{4,5,6}
1^{7,8}	120, 140, 225 bp	800 bp	A1, A36	*01:01:01:01-01:04N, 01:06-01:33, 01:35-01:109, 03:18, 03:135, 11:94, 11:112, 36:01-36:05
2⁸	210, 255, 365, 415 bp	800 bp	A2, A19, A28, A203, A210, A3, A11, A9, A26, A30, A68	*02:01:01:01-02:01:15, 02:01:17-02:01:19, 02:01:21-02:22:02, 02:24:01-02:35:01, 02:35:03-02:47, 02:49-02:77, 02:78 ^w , 02:79:01-02:97:02, 02:99, 02:101:01-02:128, 02:130-02:358, 03:01:03 ^w , 03:09 ^w , 03:23:01 ^w , 03:89 ^w , 03:108 ^w , 11:06 ^w , 11:18 ^w , 11:121 ^w , 24:28 ^w , 24:89 ^w , 26:03:01-26:03:02 ^w , 26:06 ^w , 26:07:02 ^w , 26:21 ^w , 26:30 ^w , 30:13 ^w , 30:16 ^w , 30:44 ^w , 30:46 ^w , 68:05 ^w , 68:15 ^w , 68:20 ^w , 68:30 ^w
3¹⁰	205, 235 bp	1070 bp	A1, A3, A11, A32, A34, A36	*01:12, 01:19, 01:21, 02:338, 03:01:01:01-03:17, 03:19-03:74, 03:76-03:94, 03:96-03:134, 03:136-03:143, 11:25, 11:60, 24:92, 32:04, 34:02:01, 34:02:02 ^w , 34:02:03-34:04, 34:07-34:09, 36:02
4	190 bp	800 bp	A1, A2, A3, A11, A26, A30, A36, A68	*01:01:01:01-01:01:22, 01:01:24-01:01:43, 01:01:44 ^w -01:01:45 ^w , 01:02-01:04N, 01:06-01:07, 01:09-01:11N, 01:13, 01:16N-01:18N, 01:20-01:29, 01:31N-01:33, 01:35-01:78, 01:80-01:98, 01:100-01:104, 01:105 ^w -01:109 ^w , 02:78, 02:169, 03:12, 03:18, 03:88, 03:135, 11:01:01-11:27, 11:29-11:52Q, 11:54-11:122, 26:19, 26:72, 30:08, 36:04, 68:13, 68:66
5¹¹	135, 175, 210 bp	800 bp	A2, A23, A24, A29, A80	*23:01:01-23:51, 24:05, 24:13:02, 24:24, 24:199, 29:07, 31:29, 80:01-80:02, B*18:27
6⁹	175, 200 bp	1070 bp	A2, A23, A24, A26, A33	*02:17:01 ^w -02:17:02 ^w , 23:14, 24:02:01:01-24:11N, 24:13:01-24:13:02, 24:17-24:50, 24:54-24:56, 24:58-24:63, 24:66-24:91, 24:93, 24:95-24:113, 24:115-24:137, 24:139-24:187, 24:189-24:199, 26:16, 33:19, 68:45
7	165, 200 bp	800 bp	A2/A28, A3, A10, A11, A25, A26, A32, A34, A66, A68, A69	*01:51, 02:55, 03:24, 03:50, 11:10, 25:01:01-25:16, 26:01:01-26:06, 26:08-26:15, 26:17-26:18, 26:20-26:43:02, 26:45-26:63, 26:65-26:71N, 26:73-26:74, 29:28, 32:15, 33:51, 34:01:01-34:09, 66:01-66:16, 68:01:01:01-68:88, 69:01

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

8^{7,13}	75 bp	800 bp	A25, A32	*25:01:01-25:16, 32:01:01-32:02, 32:04, 32:06-32:37, 32:39-32:41
9⁷	85 bp	1070 bp	A2, A26, A33	*01:51, 01:83, 02:146, 26:01:01-26:02, 26:04, 26:07:01-26:18, 26:20, 26:22-26:29, 26:31-26:43:02, 26:45-26:71N, 26:73-26:74, 33:13, 33:48, 68:84
10^{7,8}	80, 175, 500 bp	1070 bp	A1, A11, A24, A26, A31, A34, A66	*01:13, 01:28, 03:63, 03:88, 11:01:01-11:27, 11:29-11:52Q, 11:54-11:122, 24:19, 24:44, 26:03:01-26:03:02, 26:06, 26:21, 31:03, 34:01:01-34:08, 66:01, 66:04-66:11, 66:13-66:15, 80:02
11⁷	125, 185 bp	800 bp	A3, A10, A25, A26, A31, A34, A43, A66	*02:309, 03:01:19, 11:11, 25:05-25:06, 26:09, 26:54, 31:03-31:04, 34:01:01-34:09, 43:01, 66:02-66:03, 66:16
12	175, 225 bp	1070 bp	A1, A2, A3, A25, A26, A34, A43, A66	*01:13, 02:34-02:35:03, 02:56:01-02:56:02, 02:62, 02:103, 02:135, 03:01:01-03:01:22, 03:01:24-03:07, 03:09-03:11N, 03:13-03:31, 03:33-03:35, 03:37-03:40, 03:42-03:56, 03:58, 03:60-03:71, 03:73-03:87, 03:90-03:106, 03:109-03:110, 03:112-03:141, 03:143, 25:01:01-25:05, 25:07-25:16, 26:01:01-26:01:20, 26:01:22-26:01:24, 26:02 ^w , 26:03:01-26:03:02, 26:05-26:08, 26:10-26:33, 26:35-26:43:02, 26:45-26:72, 26:74, 30:55, 34:08, 43:01, 66:01, 66:04-66:15, 68:71, 74:13
13^{7,11}	100, 200, 240 bp	800 bp	A29, A33	*02:237, 02:309, 03:95, 26:19, 26:22, 29:01:01-29:35, 30:57, 33:13, 33:48, 66:09
14^{7,8,9,12}	90, 135, 205 bp	1070 bp	A1, A30	*01:13, 01:28, 03:43, 03:82, 11:113, 30:01:01-30:04:02, 30:06-30:20, 30:22-30:63, 31:35
15	240, 370, 395 bp	1070 bp	A31, A32	*02:237, 03:95, 29:14, 31:01:02-31:62, 32:05, 33:53
16	140, 180 bp	1070 bp	A29, A32	*01:95, 03:43, 03:82, 29:13, 31:35, 32:01:01-32:03, 32:05-32:41, 74:07
17	200 bp	1070 bp	A33, A68	*02:243, 33:01:01-33:01:06, 33:03:01-33:50, 33:52-33:59, 68:29
18	160, 200 bp	800 bp	A74	*29:19, 74:01-74:15
19¹¹	220, 245 bp	800 bp	A2, A210, A25, A30, A68	*02:34-02:35:03, 02:46, 02:48, 02:56:01-02:56:02, 02:62, 02:70, 02:78, 02:103, 02:129, 25:05, 26:54, 68:01:01-68:88
20	240, 375 bp	800 bp	A2, A24, A26, A68, A69	*02:55, 02:243, 24:82, 26:22, 33:22, 66:09, 68:29, 69:01
21^{7,11}	75, 160, 240, 495 bp	800 bp	A2, A26, A68, A80	*02:55, 03:41, 03:63, 03:75, 03:88, 24:18, 26:03:01-26:03:02, 26:05-26:06, 26:21, 26:30, 33:24, 36:02, 68:05, 68:15, 68:20, 80:01

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: 43R**Lot-specific information**

¹Alleles are assigned by the presence of specific PCR product(s). However, the sizes of the specific PCR products may be helpful in the interpretation of HLA-A low resolution SSP typings.

When the primers in a primer mix can give rise to HLA-specific PCR products of more than one length this is indicated if the size difference is more than 20 base pairs. Size differences of 20 base pairs or less are not given. For high resolution SSP kits the respective lengths of the HLA-specific PCR product(s) are given for the alleles amplified by these primer mixes.

Nonspecific amplifications, i.e. a ladder or a smear of bands, may sometimes be seen. GC-rich primers have a higher tendency of giving rise to nonspecific amplifications than other primers.

PCR fragments longer than the control bands may sometimes be observed. Such bands should be disregarded and do not influence the interpretation of the SSP typings.

PCR fragments migrating faster than the control bands, but slower than a 400 bp fragment may be seen in some gel read-outs. Such bands can be disregarded and do not influence the interpretation of the SSP typings.

Some primers may give rise to primer oligomer artifacts. Sometimes this phenomenon is an inherent feature of the primer pair(s) of a primer mix. More often it is due to other factors such as too low amount of DNA in the PCR reactions, taking too long time in setting up the PCR reactions, working at elevated room temperature or using thermal cyclers that are not pre-heated.

²The internal positive control primer pairs amplify segments of the human growth hormone gene. The two different control primer pairs give rise to either an internal positive control band of 1070 base pairs, for most wells, or a band of 800 base pairs, for some wells.

Well number 1 contains the primer pair giving rise to the shorter, 800 bp, internal positive control band in order to help in the correct orientation of the HLA-A low resolution typing. In addition, wells number 2, 4, 5, 7, 8, 9, 11, 13 and 18 to 21 contain the primer pair giving rise to the shorter, 800 bp, internal positive control band in order to allow kit identification.

In the presence of a specific amplification the intensity of the control band often decreases.

³The serological reactivity of all HLA-A alleles is not known. The grouping of not serologically defined alleles is taken from Tissue Antigens 73, 95-170, 2009.

⁴For several HLA Class I alleles 1st and/or 4th exon(s) and beyond, as well as intron nucleotide sequences, are not available, and for HLA Class II alleles 1st and/or 3th exon(s) and beyond, as well as intron nucleotide sequences, are not available. In these instances it is not known whether some of the primers of the SSP sets are completely matched with the target sequences or not. We assume that unknown sequences in these regions are conserved within allelic groups.

⁵The A*36 alleles will give rise to identical amplification patterns as a number of A*01 alleles. These alleles can be separated by the A*01 and A*36 high resolution SSP primer sets.

The A*03:01:23, 03:08, 03:32, 03:36N, 03:57, 03:59, 03:72, 03:89, 03:107-03:108, 03:111 and 03:142 and the A*24:92 alleles will give rise to identical amplification patterns. These alleles can be separated by the respective high resolution SSP primer sets.

The A*23:14 and the A*24:05, 24:13:02, 24:24 and 24:199 alleles will give rise to identical amplification patterns. These alleles can be separated by the respective high resolution SSP primer sets.

⁶Primer mix 5 will amplify the B*18:27 allele.

⁷HLA-specific PCR products shorter than 125 base pairs have a lower intensity and are less sharp than longer PCR products.

⁸The primer pairs in vials 1, 2, 10 and 14 will in many samples give rise to two or three HLA-specific PCR fragments.

⁹Primer mixes 6 and 14 may give rise to a lower yield of HLA-specific PCR product than the other HLA-A low primer mixes.

¹⁰Primer mix 3 may faintly amplify the A*30:04, 30:06, 30:17 and 30:29 alleles.

¹¹Primer mix 5, 13, 19 and 21 may have a tendency to giving rise to primer oligomer formation..

¹²Primer mix 14 may have tendencies of unspecific amplifications.

¹³Primer mix 8 may weakly amplify the A*34 alleles.

‘w’, might be weakly amplified.

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

SPECIFICITY TABLE**HLA-B low resolution primer set**

Specificities and sizes of the PCR products of the 43 primer mixes used for HLA-B low resolution SSP typing

Primer Mix	Size of spec. PCR product ¹	Size of control band ²	HLA-B serol. ³	Amplified HLA-B alleles ^{4,5}
22¹⁰	285 bp	800 bp	7, 15, 37, 40, 41, 42, 48	*07:02:01-07:06, 07:08-07:18:02, 07:20-07:32, 07:34-07:39, 07:41-07:52, 07:54-07:59, 07:61-07:76, 07:79-07:83, 07:85-07:99, 07:101-07:121, 07:123-07:138, 07:140-07:154, 08:01:01-08:05, 08:07-08:08N, 08:10-08:11, 08:13-08:15, 08:17-08:83, 08:85-08:88, 35:66, 35:87, 37:07, 40:15-40:16, 40:30-40:32, 40:34, 40:45, 40:59, 40:80, 40:98, 40:137, 40:160, 41:02:01-41:02:04, 41:04, 41:10-41:11, 41:13, 41:18-41:19, 42:01:01-42:02, 42:05:01-42:07, 42:09-42:13, 42:15-42:16, 44:150, 48:01:01-48:01:03, 48:05-48:12, 48:14-48:20, 48:22, 48:27, 53:15, 81:01-81:05
23	215 bp	1070 bp	8, 15, 44	*08:01:01-08:05, 08:07-08:25, 08:27-08:49, 08:51-08:61, 08:63-08:64, 08:66-08:84, 08:86N-08:88, 15:142, 15:180, 44:49, 51:68
24	140, 235 bp	1070 bp	7, 8, 13, 15, 35, 40, 4005, 44, 46, 53, 57, 58	*07:20, 07:24, 07:60, 07:100, 07:131, 08:21, 08:25, 13:01:01-13:04, 13:06-13:08Q, 13:10-13:23, 13:25-13:38, 13:40-13:54, 15:07:01-15:07:02, 15:36 ^w , 15:45, 15:55, 15:68, 15:89 ^w , 15:126, 15:207, 35:05:01-35:05:02, 35:16-35:17, 35:22, 35:30, 35:51, 35:58, 35:72, 35:89, 35:97, 35:113-35:114, 35:125, 35:164, 35:199, 40:05, 40:71, 40:174, 44:08 ^w , 44:54, 44:57 ^w , 44:60, 44:106, 44:110, 44:135, 46:12, 51:64, 53:14, 55:51, 58:18, C*03:05, C*03:25, C*03:27, C*03:143
25^{7,8}	130, 265 bp	800 bp	13, 14, 15, 35, 40, 41, 44, 45, 47, 49, 50	*13:01:01-13:04, 13:06-13:13, 13:15-13:23, 13:25-13:54, 14:01:01-14:04, 14:07N, 14:09, 14:11-14:12, 14:14-14:32, 15:46, 15:53, 15:106, 15:212, 18:44, 35:46-35:47, 35:63, 35:154, 40:01:01-40:01:24, 40:07, 40:10:01-40:10:02, 40:14:01-40:16, 40:22N-40:23, 40:25, 40:30-40:34, 40:36, 40:38, 40:42-40:43, 40:45, 40:47-40:49, 40:51-40:55, 40:58-40:63, 40:65-40:67, 40:69, 40:72:01-40:73, 40:76-40:77, 40:79-40:81, 40:84, 40:87:01-40:88, 40:92, 40:100-40:102, 40:106, 40:108, 40:110, 40:112-40:114, 40:116-40:118N, 40:121, 40:123-40:126, 40:128-40:130, 40:132, 40:134-40:136, 40:139-40:141, 40:146-40:147, 40:150-40:156, 40:158, 40:160, 40:163-40:164, 40:166, 40:168, 40:171-40:172, 40:175, 40:178-40:179, 40:182-40:183,

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

				40:186-40:188, 40:191-40:193, 41:01-41:20, 44:02:01:01-44:03:06, 44:03:08-44:05:03, 44:09- 44:39, 44:41:01-44:43:02, 44:45-44:56N, 44:58N- 44:80, 44:82, 44:84-44:102, 44:104-44:110, 44:112- 44:129, 44:132-44:133, 44:135-44:144, 44:146- 44:152, 45:01-45:14, 47:01:01:01-47:08, 49:01:01- 49:02, 49:04-49:20, 50:01:01-50:02, 50:04-50:06, 50:08-50:16, 57:45, 57:51, C*01:30
26⁸	185, 235 bp	800 bp	14, 16, 27, 37, 38, 39, 47, 64, 67	*07:02:32, 14:01:01-14:01:02, 14:07N-14:08, 14:10, 14:12, 14:14, 14:19, 14:26, 14:32, 27:01-27:05:15, 27:05:17-27:06, 27:08-27:10, 27:12-27:13, 27:16- 27:18, 27:20, 27:23, 27:26-27:27, 27:29, 27:31, 27:35-27:37, 27:39-27:42, 27:44-27:46, 27:48-27:61, 27:64N-27:69, 27:72-27:75, 27:77-27:80, 27:82- 27:89, 37:02, 38:01:01-38:02:02, 38:03, 38:07-38:24, 38:26-38:32, 38:34N-38:40, 39:01:01:01- 39:01:01:02L, 39:01:03-39:01:08, 39:01:10-39:02:01, 39:03, 39:05:01-39:09, 39:11, 39:14-39:15, 39:18, 39:19:02, 39:22, 39:24:01-39:39:01, 39:40N-39:44, 39:46-39:48, 39:50-39:62, 39:64-39:71, 39:73, 47:04- 47:05, 48:21, 48:26, 67:01:01, 67:03
27	190 bp	800 bp	14, 35, 38, 39, 65	*07:28, 08:87, 14:02:01-14:02:02, 14:02:04-14:02:05, 14:03-14:06:02, 14:09, 14:11, 14:13, 14:15-14:18, 14:20, 14:22-14:23, 14:25, 14:27, 14:29-14:31, 15:77, 15:189, 15:233, 35:26, 38:05, 38:33, 39:04, 44:16, 44:37:01-44:37:02, 44:64:01-44:64:02, 44:91, 44:132, 44:150, 52:26, 57:04, A*23:31, A*24:106, C*07:231, C*16:10
28	290 bp	1070 bp	15, 62, 63, 75, 76, 77, 46	*15:01:01:01-15:01:03, 15:01:06-15:02:05, 15:03:03- 15:08, 15:10:02-15:11:03, 15:11:05-15:15, 15:17:01:01-15:17:02, 15:19-15:21, 15:24-15:28, 15:30-15:36, 15:38:01-15:40, 15:42-15:46, 15:48, 15:50, 15:55-15:58, 15:60, 15:63, 15:65-15:66, 15:70-15:71, 15:73, 15:75-15:79N, 15:81-15:89, 15:92, 15:94N, 15:96-15:97, 15:101-15:102, 15:104- 15:107, 15:109-15:113, 15:116-15:118, 15:120- 15:122, 15:125-15:126, 15:128-15:129, 15:135- 15:150, 15:152, 15:154-15:155, 15:157, 15:159- 15:160, 15:162-15:172, 15:174-15:175, 15:177- 15:179, 15:181N-15:185, 15:187-15:196, 15:199, 15:201-15:209N, 15:211, 15:213-15:219, 15:223- 15:225, 15:227-15:228, 15:230-15:234, 15:236- 15:237, 15:239-15:241, 15:244, 46:01:01-46:30, C*07:02:30, C*08:16:02
29	165, 220, 330 bp	1070 bp	8, 15, 35, 37, 40, 41, 42, 44, 45, 51, 56, 57, 70, 71, 72	*08:01:01-08:05, 08:08N-08:12:03, 08:15-08:19N, 08:21-08:24, 08:26-08:27, 08:29-08:36, 08:38-08:39, 08:41-08:48, 08:50-08:54, 08:56-08:69, 08:71-08:73, 08:75-08:76, 08:78, 08:80-08:86N, 08:88, 13:46, 15:03:01-15:03:03, 15:09-15:10:02, 15:18:01- 15:18:04, 15:23, 15:29, 15:37, 15:46-15:47:02, 15:49, 15:51-15:54, 15:61-15:62, 15:64, 15:69,

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

				15:72, 15:74, 15:80, 15:83, 15:90-15:91, 15:93, 15:98-15:99, 15:103, 15:106, 15:108, 15:114-15:115, 15:119, 15:123-15:124, 15:127, 15:131-15:134, 15:143, 15:151, 15:153, 15:156, 15:158, 15:161, 15:173, 15:176, 15:186, 15:197-15:198, 15:200, 15:210, 15:212, 15:220-15:221, 15:226N, 15:229, 15:235, 15:238, 15:242-15:243, 35:87, 37:12, 38:30, 40:12, 40:136, 41:01-41:03:02, 41:05-41:09, 41:11-41:17, 41:19-41:20, 42:01:01-42:02, 42:04-42:08, 42:10-42:15, 44:14-44:15, 44:18, 44:20, 44:62, 44:100, 45:01-45:14, 51:08, 51:20, 51:36, 51:44N, 51:97, 52:19, 55:20, 55:56, 56:13, 57:09, 82:01-82:03, C*07:46
30^{8,11}	165, 190, 390 bp	1070 bp	15, 35, 40, 44, 51, 52, 56, 58, 62, 63, 70, 78	*15:01:02, 15:09, 15:16:01-15:17:02, 15:67, 15:95, 15:162, 15:168, 15:177, 15:196, 15:208, 15:216, 15:222, 15:230, 15:243, 35:01:10, 35:04:02, 40:26, 40:28, 44:62, 49:18, 50:14, 51:01:01-51:01:03, 51:01:05-51:02:03, 51:02:05-51:09:02, 51:11N-51:12, 51:13:02-51:14, 51:16-51:23, 51:24:03-51:24:04, 51:26-51:33, 51:37-51:41N, 51:43-51:44N, 51:46, 51:48-51:61, 51:63-51:80, 51:82-51:91, 51:94-51:98N, 51:100-51:130, 52:01:01:01-52:01:04, 52:01:06-52:13, 52:15-52:24, 52:26-52:29, 56:05:01-56:06, 56:21, 58:08:01-58:08:02, 78:01:01-78:07
31	180 bp	800 bp	13, 15, 40, 55, 62, 70, 71, 72, 75	*07:78, 13:01:01-13:02:12, 13:07N-13:09, 13:11, 13:14-13:20, 13:22:01-13:23, 13:25, 13:27-13:34, 13:36-13:45, 13:47, 13:49N-13:50, 13:52, 13:54, 40:48, 45:10, 49:07, 54:26, 55:09, 55:22, 55:24
32⁶	105, 195 bp	1070 bp	7, 8, 15, 27, 38, 39, 40, 44, 45, 48, 55, 72	*07:27, 07:50, 08:04, 08:17, 08:54, 15:03:01-15:03:03, 15:47:01-15:47:02, 15:49, 15:54, 15:61-15:62, 15:64, 15:68-15:69, 15:91, 15:98, 15:103, 15:123, 15:127, 15:131-15:132, 15:151, 15:156, 15:158, 15:173, 15:210, 15:220, 15:235, 15:242, 18:37, 27:18, 27:29, 37:28, 38:03, 39:02:01-39:02:02, 39:08, 39:13:01-39:13:02, 39:23, 39:39:01-39:39:02, 39:49, 40:12, 40:46, 40:93, 42:11, 44:10, 44:15, 44:18, 44:40, 44:44, 44:130, 44:140, 45:01, 45:05-45:07, 45:11-45:14, 48:01:01-48:05, 48:07-48:27, 49:20, 50:02, 52:16, 55:18, 82:01-82:03
33	195 bp	1070 bp	15, 35, 37, 38, 44, 45, 58, 76	*15:12, 15:14, 15:19, 15:91, 15:131, 15:161, 18:54, 18:56, 35:45, 35:71, 37:01:01-37:09, 37:12-37:13, 37:15-37:21, 37:23-37:28, 37:30N-37:33N, 38:17, 44:17, 44:43:01-44:43:02, 44:144, 45:09, 46:17, 53:22, 58:07, A*24:174
34^{6,10}	105, 395, 435 bp	1070 bp	7, 15, 41, 42, 59, 62, 63, 77	*07:04, 07:25, 07:146, 15:09-15:10:02, 15:13:01-15:13:02, 15:16:01-15:18:04, 15:21, 15:23-15:24, 15:37, 15:44, 15:51-15:52, 15:66-15:67, 15:72, 15:80, 15:87, 15:90, 15:93, 15:95, 15:99, 15:108, 15:114-15:115, 15:119, 15:124, 15:133-15:134, 15:153, 15:157, 15:161-15:162, 15:168,

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

				15:176-15:177, 15:186, 15:189, 15:196-15:198, 15:200, 15:208, 15:216, 15:221-15:222, 15:226N, 15:229-15:230, 15:238, 15:243, 40:136, 41:08, 42:05:01-42:05:02
35^{6,10}	115, 150 bp	1070 bp	18	*18:01:01:01-18:15, 18:17N-18:28, 18:30-18:71
36⁶	80 bp	1070 bp	27, 40, 44	*07:73, 08:88, 27:01-27:21, 27:23-27:51, 27:53-27:66N, 27:68-27:74, 27:76-27:89, 38:22, 40:46, 40:93, 44:40, 44:44, 44:130
37	150 bp	800 bp	13, 15, 18, 27, 35, 37, 39, 40, 44, 48, 51, 53, 55, 56, 57, 58, 62, 75, 77	*08:49, 13:01:01-13:01:06, 13:06-13:07N, 13:12-13:13, 13:17, 13:20, 13:22:01-13:23, 13:25-13:26, 13:28-13:29, 13:36, 13:39, 13:43, 13:50-13:52, 14:10, 15:02:01-15:02:05, 15:13:01-15:13:02, 15:20-15:21, 15:25:01-15:25:03, 15:36, 15:44, 15:62, 15:77, 15:80, 15:85, 15:88-15:89, 15:106, 15:112, 15:121, 15:139, 15:144, 15:154, 15:165, 15:170, 15:194, 15:204, 15:213-15:214, 15:223, 15:240, 18:22, 18:69, 27:19, 27:30, 35:01:01:01-35:04:03, 35:06-35:08:05, 35:10-35:17, 35:19-35:21, 35:23-35:30, 35:33-35:36, 35:38-35:42:02, 35:45-35:50, 35:52, 35:54-35:57, 35:59, 35:61-35:63, 35:65Q, 35:69-35:71, 35:74, 35:76-35:78, 35:80-35:85, 35:90-35:96, 35:98, 35:100-35:101:02, 35:103-35:113, 35:115-35:116, 35:120-35:126, 35:128-35:134N, 35:136-35:150, 35:152-35:173N, 35:175-35:184, 35:186-35:198, 35:200, 37:01:01-37:01:07, 37:03N-37:06, 37:08, 37:10-37:11, 37:13-37:33N, 38:20, 39:42, 40:28, 44:02:01:01-44:14, 44:16-44:17, 44:19N, 44:21-44:30, 44:32-44:40, 44:42-44:46, 44:48-44:52N, 44:55-44:64:02, 44:66-44:98, 44:101-44:105, 44:107-44:134, 44:136-44:137, 44:139-44:152, 48:02:01-48:02:02, 48:25, 51:04, 51:42, 51:46, 51:56:01-51:56:02, 53:01:01-53:13, 53:15-53:27, 55:14, 56:09, 56:11-56:12, 57:01:01-57:01:04, 57:01:06-57:10, 57:12, 57:14-57:20, 57:22-57:30, 57:32-57:46, 57:48-57:57, 58:01:01-58:01:02, 58:01:04-58:01:11, 58:04-58:05, 58:09-58:15, 58:17N, 58:19, 58:21-58:24, 58:28-58:36, 83:01, C*03:102
38	135 bp	1070 bp	7, 8, 18, 35, 37, 39, 56, 78	*07:65, 07:134, 08:32, 15:202, 15:239, 18:01:01:01-18:08, 18:10-18:15, 18:17N-18:36, 18:38-18:47, 18:50-18:53, 18:55-18:65, 18:68-18:71, 35:01:01:01-35:03:02, 35:03:04-35:08:04, 35:09:01-35:18, 35:20:01-35:24:02, 35:28-35:45, 35:48, 35:50-35:55, 35:57-35:62, 35:64-35:72, 35:74-35:153, 35:155-35:197, 35:199-35:200, 37:11, 39:19:01-39:19:02, 40:185, 56:06, 78:01:01-78:05, 78:07
39^{6,10}	60, 245, 400 bp	1070 bp	15, 62, 63, 70, 71, 72, 75, 76, 77, 40, 44	*15:01:01:01-15:01:04, 15:01:06-15:21, 15:23-15:40, 15:42-15:58, 15:60-15:67, 15:69-15:99, 15:102-15:129, 15:131-15:179, 15:181N-15:235, 15:237-15:244, 40:12, 82:01-82:03

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

40^{9,10}	210 bp	800 bp	7, 8, 15, 18, 37, 38	*07:27, 08:49, 18:54, 37:01:01-37:33N, 38:17
41⁸	170 bp	800 bp	15, 35, 38, 39, 67	*08:55, 15:69, 15:186, 35:35, 38:01:01-38:09, 38:11-38:21, 38:23-38:29, 38:31-38:35, 38:38-38:40, 39:01:01:01-39:01:01:02L, 39:01:03-39:20, 39:22-39:24:02, 39:26-39:29, 39:31, 39:35, 39:37-39:42, 39:44-39:46, 39:48-39:49, 39:51-39:73, 40:184, 51:101, 58:20, 67:01:01-67:03, C*15:51
42^{6,10}	110 bp	1070 bp	38	*38:01:01-38:02:02, 38:03, 38:05, 38:08-38:24, 38:26-38:34N, 38:36-38:40, 52:01:10
43	395 bp	1070 bp	14, 38, 39, 64, 65	*07:69, 07:85, 08:65, 14:01:01-14:20, 14:22-14:32, 38:01:01-38:02:03, 38:04-38:05, 38:08-38:25, 38:27-38:40, 39:01:01:01-39:01:01:02L, 39:01:03-39:01:12, 39:03-39:07, 39:09, 39:11-39:12, 39:14-39:15, 39:18-39:19:02, 39:22, 39:24:01-39:38Q, 39:40N-39:48, 39:50-39:57, 39:59-39:62, 39:64-39:73
44	160, 425 bp	1070 bp	7, 8, 14, 39, 40, 41, 42, 45, 48, 50, 60, 64, 65, 67, 73	*07:02:01-07:26, 07:28-07:35, 07:37, 07:39-07:64, 07:66-07:80, 07:82-07:148, 07:150-07:154, 08:01:01-08:01:17, 08:04-08:05, 08:07-08:35, 08:37-08:51, 08:53-08:77, 08:79-08:88, 14:01:01-14:02:04, 14:02:06-14:29, 14:31-14:32, 27:14, 27:81, 39:01:01:01-39:01:01:02L, 39:01:03-39:20, 39:22-39:24:01, 39:25N-39:73, 40:01:01-40:01:24, 40:06:01:01-40:06:04, 40:07, 40:09?, 40:11:01?-40:11:02?, 40:14:01?-40:15?, 40:16, 40:18?, 40:21?-40:22N?, 40:23, 40:24?-40:38?, 40:42?-40:43?, 40:44, 40:45?-40:46?, 40:48?-40:50?, 40:51, 40:52?, 40:53, 40:54?-40:58?, 40:59-40:60, 40:61?-40:65?, 40:66, 40:67?-40:69?, 40:70, 40:71?, 40:72:01-40:73, 40:74?, 40:75, 40:77, 40:78?, 40:79, 40:80?-40:82?, 40:83, 40:84?-40:85?, 40:86, 40:87:01?-40:88?, 40:90?-40:92?, 40:93, 40:95?, 40:96, 40:98?-40:102?, 40:103, 40:104?-40:108?, 40:109-40:110, 40:111?-40:116?, 40:118N?-40:121?, 40:123?, 40:124:01, 40:124:02?-40:126?, 40:127, 40:128?-40:130?, 40:131, 40:132?-40:134?, 40:136?-40:137?, 40:138-40:140, 40:141?, 40:145?, 40:146-40:148, 40:149?, 40:150, 40:151?, 40:152-40:153, 40:154?, 40:155N, 40:156?, 40:158?, 40:159, 40:160?, 40:161-40:162, 40:163?-40:164?, 40:165, 40:166?, 40:167, 40:168?-40:171?, 40:172, 40:173?-40:176?, 40:177, 40:178?-40:185?, 40:186, 40:187?, 40:189?, 40:190, 40:191?-40:193?, 41:01-41:20, 42:01:01-42:02, 42:04-42:16, 45:01-45:14, 48:01:01-48:17, 48:19-48:27, 50:01:01-50:02, 50:04-50:16, 51:01:09, 51:01:21, 51:01:25, 51:10, 51:24:02-51:24:04, 54:11, 56:20:02, 67:01:01-67:03, 73:01-73:02, 81:02
45¹³	180 bp	1070 bp	15, 40, 44, 46, 48, 60	*15:116, 15:124, 40:01:01-40:01:06, 40:01:08-40:01:17, 40:01:19-40:01:24, 40:07, 40:10:01-40:10:02, 40:12, 40:21-40:23, 40:25, 40:30, 40:33-

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

				40:34, 40:36, 40:38, 40:42-40:43, 40:46-40:49, 40:51-40:52, 40:54-40:55, 40:59-40:63, 40:65-40:67, 40:69, 40:73-40:74, 40:76, 40:79, 40:81, 40:84, 40:87:01-40:88, 40:92, 40:100-40:102, 40:106, 40:108, 40:112-40:114, 40:116-40:118N, 40:123-40:126, 40:128, 40:130, 40:132, 40:134-40:135, 40:137-40:141, 40:146-40:147, 40:149-40:156, 40:160, 40:163, 40:166, 40:168, 40:170-40:172, 40:175, 40:178-40:179, 40:182-40:183, 40:185-40:188, 40:191, 40:193, 44:31, 46:06, 48:03:01-48:03:02, 48:17, 48:23, C*03:129
46⁶	105 bp	1070 bp	7, 8, 35, 41, 42, 55	*07:04, 07:25, 07:146, 08:01:01-08:05, 08:07-08:12:03, 08:14-08:19N, 08:21-08:24, 08:26-08:39, 08:41-08:54, 08:56-08:69, 08:71-08:78, 08:80-08:88, 35:87, 40:136, 41:01-41:08, 41:10-41:20, 42:01:01-42:02, 42:04-42:07, 42:09-42:16, 55:20
47	325 bp	800 bp	7, 15, 27, 35, 42, 45, 46, 54, 55, 56, 67, 73, 81	*07:02:01-07:02:32, 07:04-07:07, 07:09-07:15, 07:17-07:26, 07:28-07:31, 07:33-07:36, 07:39-07:46, 07:47 ^w , 07:48-07:49N, 07:51-07:68:03, 07:70-07:84, 07:86-07:154, 14:21, 15:76, 15:101, 27:01-27:11, 27:13-27:15, 27:17, 27:19-27:21, 27:24-27:28, 27:30-27:38, 27:40-27:58, 27:60-27:76, 27:78-27:89, 35:76, 35:186, 38:26, 40:166, 42:01:01-42:02, 42:04-42:06, 42:08-42:10, 42:12-42:16, 44:90, 44:97, 45:06, 46:01:01-46:01:07, 46:02 ^w , 46:03-46:30, 54:01:01-54:26, 55:01:01-55:05, 55:07-55:17, 55:19-55:56, 56:01:01-56:16, 56:18-56:22, 56:24-56:32, 56:34-56:35, 67:01:01-67:03, 73:01-73:02, 81:01-81:05, 82:01-82:03, 83:01
48⁶	115, 195, 225, 260 bp	1070 bp	35, 37, 44, 47, 51, 57	*08:49, 08:60, 08:76, 08:87, 35:38, 35:45, 35:71, 35:115, 35:169, 37:01:01-37:01:07, 37:03N-37:06, 37:08, 37:10-37:33N, 40:132, 40:178, 42:13, 44:02:01:01-44:14, 44:16-44:17, 44:19N-44:36, 44:38-44:63, 44:65-44:134, 44:136-44:149N, 44:151-44:152, 51:42, 53:22, 57:07, 57:09, 57:24, 83:01, C*15:25
49	130, 270 bp	800 bp	18, 35, 37, 40, 46, 51, 53, 57, 58	*14:10, 15:57 ^w , 18:22, 18:69, 35:21, 35:24:01-35:24:02, 35:81, 35:96, 35:109, 35:157, 35:188, 35:190, 37:04:01-37:04:02, 40:28, 46:01:01-46:30, 51:04, 51:46, 51:56:01-51:56:02, 53:02, 53:06, 57:14, 58:09, C*15:39
50^{6,14}	90, 410 bp	1070 bp	7, 27, 35, 47, 54, 55, 56, 81	*07:65, 27:01 ^w , 27:02:01-27:11, 27:13-27:15, 27:17, 27:19-27:21, 27:24-27:28, 27:30-27:38, 27:40-27:58, 27:60-27:76, 27:78-27:89, 35:76, 35:186, 44:90, 44:97, 47:01:01:01-47:03, 47:06-47:08, 54:01:01-54:26, 55:01:01-55:05, 55:07-55:17, 55:19-55:56, 56:01:01-56:16, 56:18-56:22, 56:24-56:32, 56:34-56:35, 81:01, 82:01-82:03, 83:01
51⁶	90, 175 bp	1070 bp	15, 27, 35, 40, 48, 62,	*13:36, 15:02:01-15:02:05, 15:08, 15:11:01-15:11:05, 15:13:01-15:13:02, 15:15, 15:21, 15:31, 15:44,

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

			75, 77, 81	15:55, 15:76 ^w , 15:88-15:89, 15:112, 15:121, 15:139, 15:144, 15:148, 15:170, 15:189, 15:191, 15:194, 15:209N, 15:213-15:215, 15:223, 27:24, 35:46, 40:31, 40:45, 40:80, 48:01:01-48:01:03, 48:04, 48:06-48:07, 48:09, 48:11, 48:15-48:16, 48:18-48:20, 48:22, 48:24, 48:27, 81:01-81:05
52	145, 430 bp	1070 bp	40, 44, 49	*40:13, 40:19, 40:109, 40:117, 44:18, 44:25, 44:50, 44:95, 49:01:01-49:01:03, 49:04-49:20, 51:112, 54:12, 56:21, 57:45, 57:51, 59:01:01:01-59:05
53⁶	120, 210 bp	800 bp	13, 15, 45, 49, 50, 51, 52, 55, 56, 58, 62, 63, 78	*07:78, 07:84, 13:16, 13:31 ^w , 13:48, 15:04 ^w , 15:16:01-15:16:03, 15:42, 15:67, 15:95, 15:137 ^w , 40:95, 40:148, 40:161, 49:01:01, 49:01:03-49:10, 49:12-49:20, 50:01:01-50:02, 50:04-50:11, 50:13-50:15, 51:01:01-51:03, 51:05, 51:07:01-51:09:02, 51:11N-51:14, 51:16-51:24:04, 51:26-51:41N, 51:43-51:44N, 51:48-51:55, 51:57-51:58, 51:60, 51:61 ^w , 51:63, 51:65-51:80, 51:82-51:126, 52:01:01:01-52:27, 54:20, 55:01:01-55:01:06, 55:01:08, 55:03, 55:05, 55:09, 55:11, 55:15, 55:17, 55:21 ^w , 55:24-55:25, 55:28-55:29, 55:31, 55:33, 55:36, 55:38, 55:40, 55:44-55:45, 55:52-55:54, 56:05:01-56:06, 56:21, 56:25, 58:08:01-58:08:02, 78:01:01-78:07
54	430 bp	1070 bp	27, 40, 44, 45, 51, 52, 53, 57, 58, 59	*18:67, 27:02:01-27:02:02, 27:30, 27:53, 27:57, 27:62, 27:65N, 27:75, 27:77, 27:83, 40:13, 40:19, 40:109, 40:117, 44:06, 44:25, 44:50, 44:95, 51:01:01-51:24:04, 51:26-51:46, 51:48-51:53, 51:55-51:77, 51:79-51:131, 52:01:01:01-52:19, 52:21-52:29, 53:01:01-53:02, 53:04-53:08:02, 53:10, 53:14-53:27, 57:01:01-57:11, 57:13-57:57, 58:01:01-58:02, 58:04-58:16, 58:18-58:36
55⁸	145 bp	1070 bp	15, 35, 40, 41, 44, 45, 47, 50, 72	*07:133, 15:46, 15:53, 15:106, 15:143, 15:212, 18:48, 35:19, 35:47, 35:63, 35:154, 39:69, 40:01:01-40:11:02, 40:14:01-40:16, 40:18, 40:20, 40:22N-40:40, 40:42-40:45, 40:48-40:75, 40:77-40:92, 40:94-40:95, 40:97-40:108, 40:111-40:116, 40:118N-40:136, 40:138-40:156, 40:158-40:184, 40:186-40:187, 40:189-40:193, 41:01-41:04, 41:05?, 41:06-41:20, 44:09, 44:46, 44:75, 44:90, 44:129, 44:131, 45:01-45:14, 47:02, 47:03 ^w , 50:01:01-50:02, 50:04-50:05, 50:07-50:16
56	300 bp	1070 bp	7, 8, 13, 15, 18, 27, 35, 37, 38, 39, 40, 41, 44, 45, 46, 47, 48, 49, 50, 52, 55, 57, 58, 62,	*07:54, 07:123, 08:17, 08:38, 08:54, 13:01:01-13:04, 13:06, 13:08Q-13:23, 13:25-13:35, 13:37-13:54, 15:01:01:01-15:01:04, 15:01:06-15:01:18, 15:01:20-15:01:28, 15:03:01-15:07:02, 15:12, 15:14, 15:19-15:20, 15:24-15:27:03, 15:30, 15:32-15:36, 15:38:01-15:40, 15:42-15:43, 15:45-15:50, 15:53-15:54, 15:56-15:58, 15:60-15:63, 15:65-15:66, 15:68-15:71, 15:73-15:75, 15:77-15:79N, 15:81-15:87, 15:91-15:92,

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

			72, 76, 78	15:94N, 15:96-15:98, 15:101-15:107, 15:109-15:111N, 15:113, 15:116-15:118, 15:122-15:123, 15:125-15:129, 15:131-15:132, 15:135-15:138, 15:140-15:142, 15:145-15:147, 15:150-15:152, 15:154-15:160, 15:163-15:167, 15:169, 15:171-15:175, 15:178-15:179, 15:181N-15:185, 15:187-15:188, 15:190N, 15:192-15:193, 15:195, 15:199, 15:201-15:207, 15:210-15:212, 15:217-15:220, 15:224-15:225, 15:227-15:228, 15:231-15:237, 15:239-15:244, 18:12, 27:01-27:21, 27:24-27:89, 35:10, 35:13, 35:16, 35:28, 35:63, 35:69, 35:80, 37:01:01-37:07, 37:09-37:33N, 38:03-38:04, 38:25, 39:02:01-39:02:02, 39:08, 39:13:01-39:13:02, 39:22-39:23, 39:49, 40:01:01-40:07, 40:09-40:16, 40:18-40:24, 40:26-40:40, 40:42-40:67, 40:70-40:105, 40:107-40:165, 40:167-40:193, 41:01-41:20, 42:11, 44:02:01:01-44:02:21, 44:02:23-44:05:03, 44:07-44:11, 44:13-44:152, 45:01-45:14, 47:01:01:01-47:08, 48:01:01-48:05, 48:07-48:15, 48:17-48:27, 49:01:01-49:20, 50:01:01-50:02, 50:04-50:14, 50:16, 52:01:01:01-52:29, 53:17:01-53:17:02, 55:18, 55:34, 55:40, 78:05-78:06
57^{8,9,12}	160, 330 bp	1070 bp	15, 37, 44, 54	*08:87, 44:02:01:01-44:02:24, 44:06, 44:08-44:09, 44:11-44:12, 44:16, 44:19N-44:24, 44:27:01-44:27:02, 44:33-44:34, 44:41:01-44:41:02, 44:44, 44:48-44:49, 44:52N-44:53, 44:55, 44:59, 44:63, 44:66-44:68, 44:71-44:74, 44:80, 44:83-44:84, 44:86-44:87, 44:89-44:91, 44:93, 44:95, 44:97, 44:99-44:102, 44:104, 44:106, 44:112-44:113, 44:116, 44:118-44:119, 44:121, 44:126-44:127, 44:131-44:132, 44:137-44:138Q, 44:142Q, 44:145, 44:148-44:149N, 44:151-44:152, 54:01:01, 54:02-54:26, 83:01
58	180, 210 bp	1070 bp	7, 13, 15, 18, 35, 39, 40, 45, 54, 55, 56, 59, 78	*07:65, 13:06, 13:53, 15:42, 15:86, 15:224, 18:04, 35:08:05, 35:42:02, 35:60, 35:185, 39:17, 39:63, 40:01:11, 40:58, 45:08, 46:18, 54:01:01-54:02, 54:05N, 54:07-54:08N, 54:10, 54:12-54:13, 54:16-54:25, 55:01:01-55:03, 55:05, 55:07, 55:10-55:12, 55:15-55:16, 55:18-55:19, 55:21, 55:25-55:26, 55:29-55:31, 55:33-55:45, 55:47-55:48, 55:50, 55:52-55:56, 56:05:01, 56:06, 56:10, 56:23, 59:01:01:01-59:01:01:02, 59:05, 78:01:01-78:01:02, 78:02:02-78:03, 78:07, C*15:02:04
59	180 bp	1070 bp	13, 15, 40, 44, 45, 46, 49, 50, 51, 54, 56	*13:03, 13:48, 15:73, 40:71, 44:10, 44:15, 44:18, 44:140, 45:01, 45:04-45:07, 45:11-45:14, 46:11, 49:01:01-49:03, 49:06, 49:08-49:17, 49:19N-49:20, 50:01:01-50:02, 50:04-50:08, 50:10-50:13, 50:15-50:16, 51:15, 51:62, 51:106, 52:25, 54:03, 56:01:01-56:02, 56:04, 56:07-56:08, 56:13-56:14, 56:16-56:17, 56:20:01-56:20:02, 56:24-56:30, 56:33-56:35, 59:04, 82:01-82:03, C*03:12, C*03:19, C*03:102

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

60^{6,10}	90, 240 bp	800 bp	55, 57	*55:14, 57:01:01-57:44, 57:46-57:50, 57:52-57:57, 58:36, C*06:72
61⁶	95 bp	1070 bp	7, 8, 18, 35, 37, 38, 39, 44, 51, 52, 53, 78	*07:65 ^w , 07:134 ^w , 08:32, 18:01:01:01-18:11, 18:13-18:15, 18:17N-18:36, 18:38-18:47, 18:49-18:65, 18:67-18:71, 35:01:01:01-35:08:04, 35:09:01-35:09:03, 35:11:01-35:12:03, 35:14:01-35:15, 35:17-35:18, 35:20:01-35:24:02, 35:27, 35:29:01-35:45, 35:48, 35:50-35:62, 35:64-35:68:02, 35:70-35:72, 35:74-35:75, 35:76 ^w , 35:77-35:79, 35:81-35:153, 35:155-35:185, 35:187-35:190, 35:192-35:197, 35:199-35:200, 37:08, 38:06-38:07, 39:19:01-39:19:02, 44:06, 51:01:01-51:24:04, 51:26-51:46, 51:48-51:103, 51:105-51:111, 51:113-51:117, 51:119-51:131, 53:01:01-53:16, 53:18-53:27, 56:06 ^w , 78:01:01-78:04, 78:07
62^{6,8,10}	115, 150 bp	1070 bp	15, 18, 27, 35, 39, 46, 54, 55, 56, 62, 73, 75, 76	*07:100, 13:31, 13:41, 15:01:01:01-15:01:04, 15:01:06-15:01:16, 15:01:18-15:01:28, 15:04, 15:07:01-15:08, 15:11:01-15:12, 15:14-15:15, 15:19, 15:24, 15:26N-15:28, 15:30, 15:32, 15:34-15:35, 15:38:01-15:38:02, 15:43, 15:45-15:46, 15:50, 15:53-15:54, 15:56-15:58, 15:60, 15:63, 15:66, 15:68, 15:70-15:71, 15:73, 15:75-15:77, 15:79N, 15:81-15:82, 15:85, 15:87, 15:92, 15:94N, 15:96-15:97, 15:101-15:102, 15:104-15:105, 15:109-15:111N, 15:113, 15:117-15:118, 15:120, 15:122, 15:125-15:126, 15:128-15:129, 15:135, 15:137, 15:140, 15:142-15:149N, 15:152, 15:154, 15:157, 15:159-15:160, 15:163-15:167, 15:169, 15:171-15:172, 15:174-15:175, 15:178, 15:180-15:184, 15:187, 15:189-15:193, 15:201-15:203, 15:205-15:207, 15:209N, 15:211-15:212, 15:215, 15:217, 15:225, 15:227-15:228, 15:231-15:234, 15:236, 15:239, 15:241, 15:244, 18:19, 27:25, 27:75, 35:14:01-35:14:02, 35:43-35:44, 35:62, 35:67, 35:79, 35:86, 35:102, 35:117, 35:135, 35:163, 35:185, 39:18, 39:36, 44:146, 46:01:01-46:02, 46:04-46:05, 46:07N, 46:09-46:10, 46:12, 46:14-46:17, 46:20, 46:22-46:24, 46:27-46:30, 51:61, 52:21, 54:06, 55:21, 56:03, 73:01-73:02, A*26:68, A*68:56, C*02:56, C*06:20, C*12:50
63¹⁵	360 bp	1070 bp	Bw4	
64	350 bp	1070 bp	Bw6	

¹Alleles are assigned by the presence of specific PCR product(s). However, the sizes of the specific PCR products may be helpful in the interpretation of HLA-B low resolution SSP typings.

When the primers in a primer mix can give rise to HLA-specific PCR products of more than one length this is indicated if the size difference is more than 20 base pairs. Size differences of 20 base pairs or less are not given. For high resolution SSP kits the respective lengths of the HLA-specific PCR product(s) are given for the alleles amplified by these primer mixes.

Nonspecific amplifications, i.e. a ladder or a smear of bands, may sometimes be seen. GC-rich primers have a higher tendency of giving rise to nonspecific amplifications than other primers.

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: 43R**Lot-specific information**

PCR fragments longer than the control bands may sometimes be observed. Such bands should be disregarded and do not influence the interpretation of the SSP typings.

PCR fragments migrating faster than the control bands, but slower than a 400 bp fragment may be seen in some gel read-outs. Such bands can be disregarded and do not influence the interpretation of the SSP typings.

Some primers may give rise to primer oligomer artifacts. Sometimes this phenomenon is an inherent feature of the primer pair(s) of a primer mix. More often it is due to other factors such as too low amount of DNA in the PCR reactions, taking too long time in setting up the PCR reactions, working at elevated room temperature or using thermal cyclers that are not pre-heated.

²The internal positive control primer pairs amplify segments of the human growth hormone gene. The two different control primer pairs give rise to either an internal positive control band of 1070 base pairs, for most wells, or a band of 800 base pairs, for some wells.

Well number 22 contains the primer pair giving rise to the shorter, 800 bp, internal positive control band.

In addition, wells number 25, 26, 27, 31, 37, 40, 41, 47, 49, 53 and 60 contain the primer pair giving rise to the shorter, 800 bp, internal positive control band.

In the presence of a specific amplification the intensity of the control band often decreases.

³The serological reactivity of all HLA-B alleles is not known. The grouping of not serologically defined alleles is taken from Tissue Antigens 73, 95-170, 2009.

⁴For several HLA Class I alleles 1st and/or 4th exon(s) and beyond, as well as intron nucleotide sequences, are not available, and for HLA Class II alleles 1st and/or 3th exon(s) and beyond, as well as intron nucleotide sequences, are not available. In these instances it is not known whether some of the primers of the SSP sets are completely matched with the target sequences or not. We assume that unknown sequences in these regions are conserved within allelic groups.

The B*08:26, 08:50, 08:62 and 08:85 and the B*42:07 alleles give rise to identical amplification patterns with the HLA-B low resolution primer set. These alleles can be separated by the respective high resolution primer sets.

The B*13:04, 13:10, 13:21 and 13:35 and the B*44:135 alleles give rise to identical amplification patterns with the HLA-B low resolution primer set. These alleles can be separated by the respective high resolution primer sets.

The B*14:08 and the B*39:25N, 39:30, 39:32-39:34, 39:43, 39:47 and 39:50 alleles give rise to identical amplification patterns with the HLA-B low resolution primer set. These alleles can be separated by the respective high resolution primer sets.

The B*18:29 and the B*35:09:01-35:09:03, 35:18, 35:31-35:32:02, 35:37, 35:53N, 35:64, 35:68:01-35:68:02, 35:75, 35:88, 35:99, 35:118-35:119, 35:127, 35:151 and 35:174 alleles give rise to identical amplification patterns with the HLA-B low resolution primer set. These alleles can be separated by the respective high resolution primer sets.

The B*41:09 and the B*45:02-45:03 alleles give rise to identical amplification patterns with the HLA-B low resolution primer set. These alleles can be separated by the respective high resolution primer sets.

The B*51:104 and 51:118N and the B*58:08:01-58:08:02 alleles give rise to identical amplification patterns with the HLA-B low resolution primer set. These alleles can be separated by the respective high resolution primer sets.

The*54:01:02 and the B*55:01:07, 55:02:01-55:02:06, 55:07, 55:10, 55:12, 55:16, 55:19, 55:26, 55:30, 55:35, 55:37, 55:39, 55:41-55:43, 55:47-55:48, 55:50 and the B*56:10 alleles give rise to identical amplification patterns with the HLA-B low resolution primer set. These alleles can be separated by the respective high resolution primer sets.

The B*55:04, 55:08, 55:13, 55:23, 55:27, 55:32, 55:46 and 55:49 and the B*56:15, 56:18-56:19N, 56:22 and 56:31-56:32 alleles give rise to identical amplification patterns with the HLA-B low resolution primer set. These alleles can be separated by the respective high resolution primer sets.

The B*57:01:01-57:01:04, 57:01:06-57:03:02, 57:05-57:06, 57:08, 57:10, 57:15-57:20, 57:22-57:23, 57:25-57:30, 57:32-57:44, 57:46, 57:48-57:50 and 57:52-57:57 and the B*58:36 alleles give rise to identical amplification patterns with the HLA-B low resolution primer set. These alleles can be separated by the respective high resolution primer sets.

⁵The HLA-C*03:05, 03:25, 03:27 and 03:143 alleles will be amplified by primer mix 24, the C*01:30 allele will be amplified by primer mix 25, the A*23:31, A*24:106, C*07:231 and C*16:10 alleles will be amplified by primer mix 27, the C*07:02:30 and C*08:16:02 alleles will be amplified by primer mix 28, the C*07:46 allele will be amplified by primer mix 29, the A*24:174 allele will be amplified by primer mix 33, the C*03:102 allele will be amplified by primer mix 37, the C*15:51 allele will be amplified by primer mix 41, the C*03:129 allele will be

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: 43R**Lot-specific information**

amplified by primer mix 45, the C*15:25 allele will be amplified by primer mix 48, the C*15:39 allele will be amplified by primer mix 49, the C*15:02:04 allele will be amplified by primer mix 58, the C*03:12, 03:19 and 03:102 alleles will be amplified by primer mix 59, the C*06:72 allele will be amplified by primer mix 60 and the A*26:68, A*68:56, C*02:56, C*06:20 and C*12:50 alleles will be amplified by primer mix 62.

⁶HLA-specific PCR products shorter than 125 base pairs have a lower intensity and are less sharp than longer PCR products.

⁷Primer mix 25 may give rise to a lower yield of HLA-specific PCR product than the other HLA-B low resolution primer mixes in B*40, B*41, B*45, B*49 and B*50 alleles.

⁸Primer mixes 25, 26, 30, 41, 55, 57 and 62 may give rise to a lower yield of HLA-specific PCR product than the other HLA-B low resolution primer mixes.

⁹Primer mixes 27, 40 and 57 may have tendencies of unspecific amplifications.

¹⁰Primer mixes 22, 34, 35, 39, 40, 42, 60 and 62 have a tendency to giving rise to primer oligomer formation.

¹¹The B*15, B*46, B*57, B*58 and C*03 alleles may be faintly amplified by primer mix 30.

¹²Primer mix 57 may give rise to a lower yield of HLA-specific PCR product than the other HLA-B low resolution primer mixes in B*54 alleles.

¹³The C*17:01 to C*17:04 alleles might be faintly amplified by primer mix 45.

¹⁴Primer mix 50 might generate a false band of about 800 base pairs. This band should be disregarded when interpreting HLA-B low resolution typings.

¹⁵The Bw4-associated HLA-A specificities A23, A24, A2403, A25 and A32 are not amplified by the primer pair in primer mix 63.

‘w’, might be weakly amplified.

‘?’, nucleotide sequence information not available for the primer matching sequence.

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

SPECIFICITY TABLE

Specificities and sizes of the PCR products of the 23 primer mixes used for
 DR low resolution SSP typing

Primer Mix	Size of spec. PCR product ¹	Size of control band ²	DR serology ³	Amplified HLA-DRB alleles ⁴
65 ^{6,7}	205, 255 bp	515 bp	1	*01:01:01-01:02:08, 01:04-01:38, 01:40N-01:45
66	200 bp	430 bp	103	*01:03, 01:39N, 01:42
67	200, 215 bp	430 bp	15	*15:01:01-01-15:71
68	210 bp	430 bp	16	*16:01:01-16:05:02, 16:07-16:19
69 ^{5,6,11}	120, 220 bp	430 bp	3, 17, 18, 11	*03:01:01-01-03:75, 03:77-03:80, 03:11:02, 11:07, 11:53, 11:103, 11:105, 11:107, 11:125, 15:25
70 ^{5,11}	80, 210 bp	430 bp	3, 17, 11, 13, 14	*03:01:01-01-03:01:19, 03:04-01-03:06, 03:08-03:16, 03:18-03:20, 03:22-03:23, 03:25-03:26, 03:28, 03:30-03:31, 03:33-03:34, 03:36-03:37, 03:43-03:48, 03:50-03:52, 03:54-03:68N, 03:70-03:73, 03:75-03:80, 03:11:02, 08:40, 11:02-01-11:03, 11:11-01-11:11:02, 11:14-01-11:14:02, 11:16, 11:20-11:21, 11:36, 11:40-11:41, 11:48, 11:59, 11:63, 11:65-01-11:65:02, 11:68, 11:70, 11:73, 11:76, 11:79-11:80, 11:83, 11:85-11:87, 11:93, 11:118, 11:122, 11:124, 13:01-01-13:04, 13:08, 13:10, 13:15-13:17, 13:19-13:20, 13:22-13:24, 13:27-13:29, 13:31-13:41, 13:43, 13:45, 13:48, 13:51-13:54, 13:57, 13:59, 13:61-01-13:61:02, 13:63-13:66:02, 13:68-13:76, 13:78-13:81, 13:83-13:85, 13:87-13:99, 13:101-13:102, 13:104-13:107, 13:109, 13:111-13:117, 13:120-13:131, 13:133, 13:135, 13:137N-13:139, 14:16, 14:19, 14:21, 14:82, 14:95, 14:109
71 ⁵	85, 210 bp	430 bp	3, 18, 11, 13, 14	*03:02:01-03:03, 03:27, 03:29, 03:38, 03:53, 03:74, 11:13-01 ^w -11:13:02 ^w , 11:26, 11:34, 13:15, 13:19, 13:26, 13:44, 13:53, 13:57, 13:85-13:86, 13:104, 14:02-14:03:02, 14:06-01-14:06:02, 14:09, 14:12-01-14:13, 14:17-14:21, 14:24, 14:27, 14:29-14:30, 14:32-01 ^w -14:32:02 ^w , 14:33, 14:40-14:41, 14:47-14:49, 14:51, 14:63, 14:65 ^w , 14:67, 14:77-14:78, 14:80-14:81, 14:83, 14:85, 14:89, 14:94, 14:98, 14:102, 14:106, 14:108-14:109, 14:115, 14:119, 14:121
72 ^{5,6,7}	100, 175 bp	430 bp	3, 4	*04:01:01-04:107
73	210, 235 bp	430 bp	7, 13, 14	*07:01:01-01-07:01:04, 07:03-07:22, 12:22, 13:17, 13:116, 14:50
74 ⁶	170, 215, 250 bp	515 bp	8, 12, 14	*08:01:01-08:19, 08:21-08:49, 11:67, 12:04, 12:16, 12:22, 14:11, 14:15, 14:68, 14:93

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: 43R		Lot-specific information		
75^{5,6}	90, 135, 180 bp	430 bp	3, 9, 11	*03:08, 03:65, 09:01:02-09:17, 11:07, 11:53, 11:103, 11:105, 11:107, 11:125
76⁷	175 bp	430 bp	10	*03:76, 10:01:01-10:04, 11:59, 11:80, 11:83, 11:87, 13:27, 13:41, 13:71, 13:129
77^{5,6}	100, 170 bp	430 bp	11, 3, 8	*03:08, 03:65, 08:31, 08:41, 11:01:01-11:70, 11:72-11:126
78^{5,6,9}	85, 105 bp	430 bp	12, 8	*08:32, 12:01:01-12:35
79^{7,8}	215 bp	430 bp	6, 8, 11, 13, 14, 1403	*03:76, 08:20-08:21, 11:01:01-11:04:08, 11:06:01-11:06:02, 11:08:01-11:12:02, 11:14:01-11:16, 11:18-11:21, 11:23-11:25, 11:27:01-11:33, 11:35-11:51, 11:54:01-11:54:02, 11:56-11:66, 11:68, 11:70, 11:72-11:81, 11:83-11:88, 11:90-11:97, 11:99-11:102:02, 11:106, 11:108-11:124, 11:126, 13:01:01-13:08, 13:10-13:16, 13:18-13:43, 13:45-13:85, 13:87-13:115, 13:117-13:128, 13:130-13:139, 14:03:01-14:03:02, 14:12:01-14:12:02, 14:16, 14:19, 14:21-14:22, 14:25, 14:27, 14:40, 14:53, 14:63, 14:67, 14:69, 14:74, 14:77-14:78, 14:84-14:85, 14:98, 14:102, 14:105, 14:109, 14:115-14:116, DRB3*02:27
80^{6,7,8}	195, 225 bp	430 bp	13, 8, 11, 12, 14	*08:01:01-08:02:04, 08:04:01-08:09, 08:11, 08:16-08:17, 08:20-08:22, 08:24, 08:26, 08:28, 08:31, 08:39, 08:41-08:44, 11:01:01-11:06:02, 11:09-11:12:02, 11:14:01-11:16, 11:20-11:21, 11:23-11:25, 11:27:01-11:30, 11:32-11:33, 11:35-11:41, 11:43-11:44, 11:46:01-11:51, 11:54:01-11:56, 11:58:01-11:63, 11:65:01-11:70, 11:72, 11:74:01-11:78, 11:80-11:88, 11:90-11:97, 11:99-11:102:02, 11:106, 11:108-11:118, 11:120-11:124, 11:126, 12:02:01-12:02:05, 12:13, 12:15-12:16, 12:18-12:21, 12:23, 12:26-12:27, 12:31N-12:33, 13:01:01-13:02:01, 13:02:03-13:02:05, 13:04-13:05:02, 13:07:01-13:09, 13:11:01-13:11:02, 13:14:01-13:24, 13:26-13:29, 13:31-13:32, 13:34-13:36, 13:38-13:43, 13:45-13:55, 13:57, 13:59, 13:61:01-13:65, 13:67-13:76, 13:78-13:80, 13:83-13:84, 13:87, 13:91-13:93, 13:96:01-13:100, 13:102-13:109, 13:111-13:114, 13:116-13:117, 13:121, 13:123-13:132, 13:135-13:136, 13:138-13:139, 14:15-14:16, 14:22, 14:24-14:25, 14:27, 14:37, 14:53, 14:73, 14:105
81¹¹	175 bp	430 bp	3, 13, 14, 8	*03:01:01:01-03:07, 03:09, 03:11:01-03:41, 03:43-03:45, 03:47-03:63, 03:66-03:68N, 03:70-03:80, 03:11:02, 08:20, 13:01:01-13:16, 13:18-13:42, 13:44, 13:46-13:66:02, 13:68-13:102, 13:104-13:115, 13:117-13:121, 13:123-13:139, 14:02-14:03:02, 14:05:01-14:06:02, 14:09, 14:12:01-14:14, 14:17-14:21, 14:23:01,

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

				14:23:03-14:24, 14:27, 14:29-14:30, 14:33, 14:36-14:37, 14:40-14:45, 14:47-14:48, 14:51, 14:56, 14:59, 14:63-14:65, 14:67, 14:77-14:78, 14:80-14:81, 14:83-14:85, 14:89, 14:91, 14:94-14:96, 14:98, 14:100, 14:102-14:103, 14:106, 14:108-14:109, 14:115-14:116, 14:121, 14:123
82⁵	100, 140, 155 bp	430 bp	14, 4, 8, 13	*04:62, 04:69, 04:73, 04:105, 08:08, 11:69, 11:82, 13:45, 14:01:01-14:01:02, 14:04, 14:07:01-14:07:02, 14:10, 14:16, 14:22, 14:25-14:26, 14:28, 14:31-14:32:02, 14:35, 14:37-14:39, 14:49-14:50, 14:53-14:54:01, 14:55, 14:57-14:58, 14:60-14:62, 14:68-14:71, 14:73-14:76, 14:79, 14:82, 14:86-14:88, 14:90, 14:93, 14:99, 14:101, 14:104-14:105, 14:107, 14:110-14:114, 14:117-14:120, 14:122, 14:124-14:125, DRB4*01:03:01:02N
83^{5,6,10}	110, 135, 170 bp	430 bp	14, 3, 9, 11, 12, 13, 15	*03:10, 09:01:02-09:01:05, 09:01:07-09:02:02, 09:04-09:17, 11:13:01-11:13:02, 11:17, 11:52, 13:43, 14:01:01-14:02, 14:04-14:11, 14:13-14:14, 14:16-14:18, 14:19 ^w , 14:20, 14:21 ^w , 14:22-14:23:03, 14:26, 14:28-14:36, 14:38-14:39, 14:41, 14:43-14:52, 14:54:01-14:57, 14:59-14:62, 14:64-14:65, 14:68, 14:70-14:76, 14:79-14:83, 14:86-14:88, 14:90-14:97, 14:99-14:101, 14:103-14:108, 14:109 ^w , 14:110-14:114, 14:117-14:125, 15:27, 15:34, 15:66
84^{5,6}	110, 175, 225 bp	430 bp	14, 3, 8, 11, 13, 15, 16	*03:10, 08:09, 08:20-08:21, 08:32, 08:35, 11:13:01-11:13:02, 11:17, 11:23, 11:25, 11:31, 11:45, 11:52, 11:55, 11:64, 11:89, 11:96, 11:119, 13:13, 13:18, 13:43, 13:45, 13:47, 13:55, 13:119, 14:01:01-14:01:03, 14:03:01-14:05:03, 14:07:01-14:08, 14:10-14:12:02, 14:14-14:16, 14:18, 14:22-14:23:03, 14:25-14:28, 14:31-14:32:02, 14:34-14:36, 14:38-14:40, 14:42-14:45, 14:49-14:50, 14:53-14:65, 14:67-14:79, 14:81-14:82, 14:84-14:93, 14:95-14:97, 14:99-14:105, 14:107, 14:110-14:120, 14:122-14:125, 15:21 ^w , 16:04 ^w , 16:18 ^w
85^{6,11}	160, 240 bp	430 bp	52	DRB3*01:01:02:01-01:15, DRB3*02:01:02:29N, DRB3*03:01:01-03:03
86^{7,12}	215 bp	430 bp	53	DRB4*01:01-01:08
87	175 bp	430 bp	51	DRB5*01:01:01-01:14, DRB5*02:02-02:05

¹Alleles are assigned by the presence of specific PCR product(s). However, the sizes of the specific PCR products may be helpful in the interpretation of DR low resolution SSP typings.

When the primers in a primer mix can give rise to specific PCR products of more than one length this is indicated if the size difference is 20 base pairs or more. Size differences shorter than 20 base pairs are not given. For high resolution SSP kits the respective lengths of the specific PCR product(s) of the alleles amplified by these primer mixes are given.

Nonspecific amplifications, i.e. a ladder or a smear of bands, may sometimes be seen. GC-rich

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: 43R**Lot-specific information**

primers have a higher tendency of giving rise to nonspecific amplifications than other primers, e.g. the primers in wells 67, 72, 82, 83 and 84.

PCR fragments longer than the control bands may sometimes be observed. Such bands should be disregarded and do not influence the interpretation of the SSP typings.

PCR fragments migrating faster than the control bands, but slower than a 400 bp fragment may be seen in some gel read-outs. Such bands can be disregarded and do not influence the interpretation of the SSP typings.

Some primers may give rise to primer oligomer artifacts. Sometimes this phenomenon is an inherent feature of the primer pair(s) of a primer mix. More often it is due to other factors such as too low amount of DNA in the PCR reactions, taking too long time in setting up the PCR reactions, working at elevated room temperature or using thermal cyclers that are not pre-heated.

²The internal positive control primer pairs amplify segments of the human growth hormone gene. The two different control primer pairs give rise to either an internal positive control band of 430 base pairs, for most wells, or a band of 515 base pairs, for some wells.

Well number 65 contains the primer pair giving rise to the longer, 515 bp, internal positive control band. In addition, well number 74 contains the primer pair giving rise to the longer, 515 bp, internal positive control band.

In the presence of a specific amplification the intensity of the control band often decreases.

³The serological reactivity of all DRB alleles is not known. The grouping of not serologically defined alleles is taken from Tissue Antigens 73, 95-170, 2009.

⁴For several DRB1 alleles 1st and/or 3rd exon(s) and beyond, as well as intron nucleotide sequences, are not available. In these instances it is not known whether some of the primers of the SSP sets are completely matched with the target sequences or not. We assume that unknown sequences in these regions are conserved within allelic groups.

The DRB1*08:09 and the DRB1*14:15 alleles yield identical amplification patterns with the DR low resolution primer set. These alleles can be separated by the respective high resolution primer sets.

The DRB1*08:20 and the DRB1*13:18, 13:47 and 13:55 alleles yield identical amplification patterns with the DR low resolution primer set. These alleles can be separated by the respective high resolution primer sets.

The DRB1*08:31, 08:41 and DRB1*11:67 alleles yield identical amplification patterns with the DR low resolution primer set. These alleles can be separated by the respective high resolution primer sets.

The DRB1*13:13 and 13:119 and the DRB1*14:84 and 14:116 alleles yield identical amplification patterns with the DR low resolution primer set. These alleles can be separated by the respective high resolution primer sets.

⁵Specific PCR fragments shorter than 125 base pairs have a lower intensity and are less sharp than longer PCR bands.

⁶Individual alleles can give rise to two differently sized specific PCR fragments in primer mixes 65 and 69, 72, 74, 75, 77, 78, 80 and 83 to 85.

⁷Primer mixes 65, 72, 76, 79, 80 and 86 have a tendency to giving rise to primer oligomer formation.

⁸Primer mixes 79 and 80 may give rise to a lower yield of HLA-specific PCR product than the other DR low resolution primer mixes.

⁹Primer mix 78 may have tendencies of unspecific amplifications.

¹⁰Primer mix 83 has a tendency of primer oligomer formation and also has an intense primer cloud due to the high number of primers present in the primer mix.

¹¹Due to sharing of sequence motifs in codon 38, DRB3*01:14 will also be amplified in primer mixes 69, 70 and 81 in addition to primer mix 85.

¹²The DRB4*01:03:01:02N allele is amplified by primer mix 86, whereas the DRB4*02:01N and DRB4*03:01N null alleles are not amplified by this primer mix.

'w', might be weakly amplified.

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

SPECIFICITY TABLE**DQ low resolution SSP typing**

Specificities and sizes of the PCR products of the 8 primer mixes used for DQ low resolution SSP typing

Primer Mix	Size of spec. PCR product ¹	Size of control band ²	DQ serology ³	Amplified DQB1 alleles ⁴
88	225 bp	515 bp	5	*05:01:01:01-05:14
89⁵	220, 270 bp	430 bp	1, 5, 6	*06:01:01-06:49
90	210 bp	430 bp	2	*02:01:01-02:06
91⁶	220 bp	515 bp	3, 7	*03:01:01:01-03:01:06, 03:04, 03:09-03:10, 03:13-03:14, 03:16, 03:19, 03:21-03:22, 03:24, 03:27-03:29, 03:35-03:36
92⁶	130 bp	515 bp	6, 8	*03:02:01-03:02:05, 03:05:01-03:05:04, 03:07-03:08, 03:11, 03:18, 03:32, 03:37, 06:29
93^{5,6,7}	135 bp	515 bp	2, 3, 9	*02:03, 03:03:02:01-03:03:04, 03:06, 03:12, 03:15, 03:20, 03:25-03:26, 03:30-03:31, 03:33-03:34, 03:38-03:40, 04:03:01-04:03:02
94^{6,7}	145, 185 bp	515 bp	3, 7, 8, 9	*03:01:01:01-03:40
95⁶	160 bp	430 bp	4	*04:01:01-04:08
96⁸				Negative control

¹ Alleles are assigned by the presence of specific PCR product(s). However, the sizes of the specific PCR products may be helpful in the interpretation of DQ low resolution SSP typings.

When the primers in a primer mix can give rise to specific PCR products of more than one length this is indicated if the size difference is 20 base pairs or more. Size differences shorter than 20 base pairs are not given. For high resolution SSP kits the respective lengths of the specific PCR product(s) of the alleles amplified by these primer mixes are given.

Nonspecific amplifications, i.e. a ladder or a smear of bands, may sometimes be seen. GC-rich primers have a higher tendency of giving rise to nonspecific amplifications than other primers.

PCR fragments longer than the control bands may sometimes be observed. Such bands should be disregarded and do not influence the interpretation of the SSP typings.

PCR fragments migrating faster than the control bands, but slower than a 400 bp fragment may be seen in some gel read-outs. Such bands can be disregarded and do not influence the interpretation of the SSP typings.

Some primers may give rise to primer oligomer artifacts. Sometimes this phenomenon is an inherent feature of the primer pair(s) of a primer mix. More often it is due to other factors such as too low amount of DNA in the PCR reactions, taking too long time in setting up the PCR reactions, working at elevated room temperature or using thermal cyclers that are not pre-heated.

² The internal positive control primer pairs amplify segments of the human growth hormone gene. The two different control primer pairs give rise to either an internal positive control band of 430 base pairs, for most wells, or a band of 515 base pairs, for some wells.

Well number 88 contains the primer pair giving rise to the longer, 515 bp, internal positive control band. In addition, wells number 91 to 94 contain the primer pair giving rise to the longer, 515 bp, internal positive control band.

In the presence of a specific amplification the intensity of the control band often decreases.

101.708-24 – including *Taq* polymerase, IFU-01
101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: 43R**Lot-specific information**

³The serological split of the DQB1*05:05, DQB1*06:06 to 06:07 alleles, the DQB1*06:10, 06:13, 06:15-06:24 and 06:27 to 06:44, the DQB1*02:04-02:06, the DQB1*03:07-03:09 and 03:11- 03:38 alleles and the DQB1*04:0301-04:08 alleles is not known. In this table we have used the expert-assigned serological grouping in Tissue Antigens (2009) 73:95-170, and also inferred the serological grouping from the naming of the sequence-defined allele.

⁴For several DQB1 alleles only partial second exon nucleotide sequences are available. In these instances it is not known whether some of the primers of the SSP sets are completely matched with the target sequences or not. We assume that unknown sequences in the 5'- and 3'-ends of the second exon of the DQB1 gene are conserved within allelic groups.

⁵Primer mixes 89 and 93 have a tendency to giving rise to primer oligomer formation.

⁶Primer mixes 91, 92, 93, 94 and 95 may give rise to a lower yield of HLA-specific PCR product than the other DQ low resolution primer mixes.

⁷Primer mixes 93 and 94 may have tendencies of unspecific amplifications.

⁸Primer mix 96 contains a negative control, which will amplify more than 95% of HLA amplicons as well as the amplicons generated by control primer pairs. PCR product sizes range from 75 to 200 base pairs. The PCR product generated by the control primer pair is 430 base pairs.

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

INTERPRETATION TABLE																						
HLA-A low resolution SSP typing																						
Amplification patterns of the A*01:01 to A*80:02 alleles																						
		Well ¹⁰																				
		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21
Length of spec.		120	210	205	190	135	175	165	75	85	80	125	175	100	90	240	140	200	160	220	240	75
PCR product(s)		145	255	235		175	205	200			175	185	225	200	135	370	180		200	245	375	160
		225	365			210					500			240	205	395						240
			415																			495
Length of int.		800	800	1070	800	800	1070	800	800	1070	1070	800	1070	800	1070	1070	1070	1070	800	800	800	800
pos. control ¹																						
5'-primer(s) ²		98	48	363	98	176	98	98	266	257	301	103	98	98	203	41	180	98	180	78	28	176
		5'-CTT ³	5'-gCT ³	5'-ATA ³	5'-CTA ³	5'-gCA ³	5'-CTC ³	5'-CTA ³	5'-ACg ³	5'-Cgg ³	5'-Cgg ³	5'-CCT ³	5'-CTT ³	5'-CAC ³	5'-gAA ³	5'-CTT ³	5'-TTT ³	5'-CAC ³	5'-TTT ³	5'-TCT ³	5'-TCg ³	5'-gCA ³
		103	78		413	368	368	102	266	259	302	423	423	238	362	355	203			106	261	261
		5'-CCT ³	5'-TCT ³		5'-CCg ³	5'-gTT ³	5'-gTT ³	5'-ACA ³	5'-ACg ³	5'-CgA ³	5'-ggA ³	5'-gCT ³	5'-gCT ³	5'-AgA ³	5'-ggT ³	5'-CCg ³	5'-gAA ³			5'-CCA ³	5'-AAC ³	5'-AAC ³
		123	106					413	266	261	385			355	363		418			2 nd		341
		5'-AgT ³	5'-CCA ³					5'-CCg ³	5'-Agg ³	5'-AAC ³	5'-ggC ³			5'-CCg ³	5'-ATA ³		5'-AgC ³			5'-CCT ³		5'-ggA ³
		363													363							355
		5'-ATA ³													5'-ATA ³							5'-CCC ³
															5'-TAC ³							
3'-primer(s) ³		203	240	527	256	270	259	259	302	299	341	257	282	257	299	238	290	256	299	265	97	292
		5'-TCT ³	5'-ggA ³	5'-CCA ³	5'-CTg ³	5'-ACA ³	5'-gTT ³	5'-gTT ³	5'-ggC ³	5'-TCg ³	5'-CgT ³	5'-gCA ³	5'-gAC ³	5'-gCA ³	5'-CCA ³	5'-CCT ³	5'-CAA ³	5'-CCC ³	5'-CCA ³	5'-CCC ³	5'-ggT ³	5'-gTg ³
		545	292	527	559	502	502	259			521	506	282	299	411	243	317	256	299	282	355	292
		5'-AgA ³	5'-gTg ³	5'-CCT ³	5'-CCg ³	5'-CTg ³	5'-CTT ³	5'-gTT ³			5'-ggg ³	5'-TgT ³	5'-gAC ³	5'-TCg ³	5'-TCA ³	5'-TCA ³	5'-ggA ³	5'-CTC ³	5'-CCA ³	5'-gAC ³	5'-gAC ³	5'-gTT ³
				555		538	539	538				559	559	418	526	265	555	259	341	282		299
				5'-gCA ³		5'-Cag ³	5'-TCT ³	5'-CCA ³				5'-CTC ³	5'-CCg ³	5'-gTC ³	5'-CCA ³	5'-CCC ³	5'-CCA ³	5'-gTT ³	5'-CgT ³	5'-gAC ³		299
				555								559		555		555				502		555
				5'-CCA ³								5'-CgT ³		5'-CCA ³		5'-CCA ³				5'-CTT ³		5'-CCA ³
																				506		
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21
HLA-A allele ⁴	ser ⁵																					
*01:01:01:01-01:01:22, 01:01:24-01:01:43, 01:02- 01:04N, 01:06-01:07, 01:09- 01:11N, 01:16N-01:18N, 01:20, 01:22N-01:27N, 01:29, 01:31N-01:33, 01:35- 01:50, 01:52N-01:78, 01:80- 01:82, 01:84-01:94, 01:96- 01:98, 01:100-01:104, 36:04 ⁶	A1, A36, Null, -	1			4																	
*01:01:23, 01:08, 01:14- 01:15N, 01:30, 01:79, 01:99, 36:01, 36:03, 36:05 ⁶	A1, A36, Nu ll, -	1																				
*01:01:44-01:01:45, 01:105- 01:109	-	1			w																	
*01:12, 01:19	A1, -	1		3																		
*01:13	-	1			4						10		12		14							
*01:21	-	1		3	4																	
*01:28	-	1			4						10				14							
*01:51	-	1			4			7		9												
*01:83	-	1			4					9												
*01:95	-	1			4											16						
*02:01:01:01-02:01:15, 02:01:17-02:01:19, 02:01:21- 02:16, 02:18-02:22:02, 02:24:01-02:33, 02:36- 02:45, 02:47, 02:49-02:54, 02:57-02:61, 02:63-02:69, 02:71-02:77, 02:79:01- 02:97:02, 02:99, 02:101:01- 02:102, 02:104-02:128, 02:130-02:134, 02:136- 02:145, 02:147-02:168, 02:170-02:236, 02:238- 02:242, 02:244-02:308, 02:310-02:337, 02:339- 02:358	A2, 19, A203, A210, Low A2, Null, -		2																			
*02:17:01-02:17:02	A2		2			w																
*02:34-02:35:01, 02:35:03, 02:56:01-02:56:02, 02:62, 02:103	A2, -		2										12							19		
*02:35:02	A2												12							19		
*02:46, 02:70	A2		2																	19		
*02:48, 02:129	-																			19		
*02:55	-		2					7													20	21
*02:78	A2		w		4															19		
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21
HLA-A allele ⁴	ser ⁵																					
*02:135	-		2										12									
*02:146	-		2							9												
*02:169	-		2		4																	
*02:237	-		2											13		15						
*02:243	-		2															17			20	
*02:309	-		2									11		13								
*02:338	-		2	3																		
*03:01:01:01-03:01:02, 03:01:04-03:01:18, 03:01:20- 03:01:22, 03:01:24-03:07, 03:10-03:11N, 03:13-03:17, 03:19-03:22, 03:23:02, 03:25- 03:31, 03:33-03:35, 03:37- 03:40, 03:42, 03:44-03:49, 03:51-03:56, 03:58, 03:60- 03:62, 03:64-03:71, 03:73- 03:74, 03:76-03:81, 03:83- 03:87, 03:90-03:94, 03:96- 03:106, 03:109-03:110, 03:112-03:134, 03:136- 03:141, 03:143	A3, Null, -			3									12									
*03:01:03, 03:09, 03:23:01	A3, -		w	3									12									
*03:01:19	-			3								11	12									
*03:01:23, 03:08, 03:32, 03:36N, 03:57, 03:59, 03:72, 03:107, 03:111, 03:142, 24:92 ⁷	Null, -			3																		
*03:12	-			3	4																	
*03:18, 03:135	-	1			4								12									
*03:24, 03:50	A3, -			3				7					12									
*03:41	-			3																		21
*03:43, 03:82	-			3									12		14		16					
*03:63	-			3							10		12									21
*03:75	-												12									21
*03:88	-			3	4						10											21
*03:89, 03:108	-		w	3																		
*03:95	-												12	13		15						
*11:01:01-11:05, 11:07- 11:09, 11:12-11:17, 11:19- 11:24:02, 11:26-11:27, 11:29- 11:52Q, 11:54-11:59, 11:61- 11:93, 11:95-11:111, 11:114- 11:120, 11:122	A11, -, Null				4						10											
*11:06, 11:18, 11:121	A11, -		w		4						10											
*11:10	A11				4			7			10											
*11:11	-				4						10	11										
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21
HLA-A allele ⁴	ser ⁵																					
*11:25, 11:60	A11, -			3	4						10											
*11:94, 11:112	-	1			4						10											
*11:113	-				4						10				14							
*23:01:01-23:13, 23:15-23:51, B*18:27 ⁸	A23(9), Null, -					5																
*23:14, 24:05, 24:13:02, 24:24, 24:199 ⁹	A9, A24(9), -					5	6															
*24:02:01:01-24:04, 24:06-24:11N, 24:13:01, 24:17, 24:20-24:23, 24:25-24:27, 24:29-24:43, 24:45N-24:50, 24:54-24:56, 24:58-24:63, 24:66-24:81, 24:83N-24:88, 24:90-24:91, 24:93, 24:95-24:113, 24:115-24:137, 24:139-24:187, 24:189-24:198	A24(9), Low A24(9), A2403, A9, Null, -						6															
*24:18	A24(9), A3						6															21
*24:19, 24:44	A9, -						6				10											
*24:28, 24:89	A9, -		w				6															
*24:82	-						6														20	
*25:01:01-25:04, 25:07-25:16	A25 (10), Null, -							7	8				12									
*25:05	-							7	8			11	12						19			
*25:06	-							7	8			11										
*26:01:01-26:01:20, 26:01:22-26:01:24, 26:08, 26:10-26:15, 26:17-26:18, 26:20, 26:23-26:29, 26:31-26:33, 26:35-26:43:02, 26:45-26:53, 26:55-26:63, 26:65-26:71N, 26:74	A26 (10), A10, Null, -							7		9			12									
*26:01:21, 26:04, 26:34, 26:73	A26 (10), -							7		9												
*26:02	A26 (10)							7		9			w									
*26:03:01-26:03:02, 26:06, 26:21	A26 (10), -		w					7			10		12								21	
*26:05	A26 (10)							7					12								21	
*26:07:01, 26:64	A26 (10), -									9			12									
*26:07:02	-		w							9			12									
*26:09	A26 (10)							7		9		11										
*26:16	-						6			9			12									
*26:19	-				4								12	13								
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21
HLA-A allele ⁴	ser ⁵																					
*26:22	A26 (10)							7		9			12	13							20	
*26:30	-		w					7					12									21
*26:54	-							7		9		11	12							19		
*26:72	-				4								12									
*29:01:01-29:06, 29:08N- 29:12, 29:15-29:18, 29:20- 29:27, 29:29-29:35	A29 (19), Null, -													13								
*29:07	-					5								13								
*29:13	-													13			16					
*29:14	-													13		15						
*29:19	-													13					18			
*29:28	-							7						13								
*30:01:01-30:04:02, 30:06- 30:07, 30:09-30:12, 30:14- 30:15, 30:17-30:20, 30:22- 30:43, 30:45, 30:47-30:54, 30:56, 30:58-30:63	A30 (19), Null, -														14							
*30:08	-				4										14							
*30:13, 30:16, 30:44, 30:46	A30 (19), -		w												14							
*30:55	-												12		14							
*30:57	-													13	14							
*31:01:02-31:02, 31:05- 31:28, 31:30-31:34, 31:36- 31:62	A31 (19), Null, -															15						
*31:03	-										10	11				15						
*31:04	A31 (19)											11				15						
*31:29	-					5										15						
*31:35	-														14	15	16					
*32:01:01-32:02, 32:06- 32:14, 32:16-32:37, 32:39- 32:41	A32 (19), Null, -								8								16					
*32:03, 32:38	-																16					
*32:04	-			3					8													
*32:05	-															15	16					
*32:15	-							7	8								16					
*33:01:01-33:01:06, 33:03:01-33:12, 33:14- 33:18, 33:20-33:21, 33:23, 33:25-33:47, 33:49-33:50, 33:52, 33:54-33:59	A33 (19, A19), -																	17				
*33:13, 33:48	-									9				13				17				
*33:19	-					6												17				
*33:22	-																	17			20	
*33:24	-																	17				21
*33:51	-							7														
*33:53	-															15		17				
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21
HLA-A allele ⁴	ser ⁵																					
*34:01:01-34:01:02, 34:05-34:06	A34 (10), -							7			10	11										
*34:02:01, 34:02:03-34:04, 34:07	A34 (10), -			3				7			10	11										
*34:02:02	-			w				7			10	11										
*34:08	-			3				7			10	11	12									
*34:09	-			3				7				11										
*36:02	-	1		3																		21
*43:01	A43											11	12									
*66:01, 66:04-66:08, 66:10-66:11, 66:13-66:15	A66 (10), -							7			10		12									
*66:02-66:03, 66:16	A66 (10), A10, -							7				11										
*66:09	-							7			10		12	13							20	
*66:12	-							7					12									
*68:01:01:01-68:04, 68:06-68:12, 68:14, 68:16-68:19, 68:21:01-68:28, 68:31-68:44, 68:46-68:65, 68:67-68:70, 68:72-68:83, 68:85-68:88	A68 (28), A28, Null, -							7												19		
*68:05, 68:15, 68:20	A68 (28), -		w					7												19		21
*68:13, 68:66	-				4			7												19		
*68:29	-							7										17		19	20	
*68:30	A68 (28)		w					7												19		
*68:45	-						6	7												19		
*68:71	-							7					12							19		
*68:84	-							7		9										19		
*69:01	A69 (28)							7													20	
*74:01-74:06, 74:08-74:12N, 74:14N-74:15	A74 (19), Null, -																			18		
*74:07	A74 (19)															16				18		
*74:13	-												12							18		
*80:01	A80					5																21
*80:02	-					5					10											
HLA-A allele ⁴	ser ⁵																					
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: **43R****Lot-specific information**

¹The internal positive control primer pairs amplify segments of the human growth hormone gene. The two different control primer pairs give rise to either an internal positive control band of 1070 base pairs, for most wells, or a band of 800 base pairs, for some wells.

Well number 1 contains the primer pair giving rise to the shorter, 800 bp, internal positive control band in order to help in the correct orientation of the HLA-A low resolution SSP typing.

In addition, wells number 2, 4, 5, 7, 8, 9, 11, 13 and 18 to 21 contain the primer pair giving rise to the shorter, 800 bp, internal positive control band in order to allow kit identification.

²The nucleotide position, in the 1st, 2nd or 3rd exon or the 2nd intron, matching the specificity-determining 3'-end of the primer is given. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

³The nucleotide position, in the 2nd or 3rd exon, matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

⁴The sequence of the A*0105N has been shown to be identical to A*01:04N.

The A*01:34N allele has been renamed A*01:01:38L

The A*020116 allele has been renamed to A*02:134.

The A*020120 has been renamed to A*02:01:18.

The sequence of the A*0223 allele has been shown to be identical to A*02:22:01.

The sequence of the A*0298 allele has been shown to be identical to A*02:96.

The A*1128 allele has been renamed to A*11:15:02.

The A*11:53 allele has been shown to be identical to the corrected A*11:02:01.

The sequence of the A*2401 allele has been shown to be in error.

The sequence of the A*2412 allele has been shown to be identical to A*24:08.

The A*2416 allele has been renamed to A*31:08.

The A*2465 allele has been renamed to A*24:13:02.

The A*26:44 allele has been renamed to A*26:43:02.

The sequence of the A*3005 allele has been shown to be identical to A*30:04.

The A*3021 allele has been renamed to A*30:11:02.

The sequence of the A*31011 allele has been shown to be identical to A*31:01:02.

The sequence of the A*3302 allele has been shown to be identical to A*33:03:01.

⁵The serological reactivity of all HLA-A alleles is not known. The grouping of not serologically defined alleles is taken from Tissue Antigens 73, 95-170, 2009.

⁶The A*36 alleles will give rise to identical amplification patterns as a number of A*01 alleles. These alleles can be separated by the A*01 and A*36 high resolution SSP primer sets.

⁷The A*03:01:23, 03:08, 03:32, 03:36N, 03:57, 03:59, 03:72, 03:89, 03:107-03:108, 03:111 and 03:142 and the A*24:92 alleles will give rise to identical amplification patterns. These alleles can be separated by the respective high resolution SSP primer sets.

⁸Primer mix 5 will amplify the B*18:27 allele.

⁹The A*23:14 and the A*24:05, 24:13:02, 24:24 and 24:199 alleles will give rise to identical amplification patterns. These alleles can be separated by the respective high resolution SSP primer sets.

¹⁰The primer pairs in vials 1, 2, 10 and 14 will in many samples give rise to two or three HLA-specific PCR fragments.

‘w’, may be weakly amplified.

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

INTERPRETATION TABLE		HLA-B low resolution SSP																			
Amplification patterns of the B*07:02 to B*83:01 alleles		Well																			
		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41
Length of spec. PCR product(s)		285	215	140	130	185	190	290	165	165	180	105	195	105	115	80	150	135	60	210	170
				235	265	235			220	190		195		395	150				245		
									330	390				435					400		
Length of int. pos. control ¹		800	1070	1070	800	800	800	1070	1070	1070	800	1070	1070	1070	1070	1070	1070	800	1070	1070	800
5'-primer(s) ²		355	97	209	103	103	103	45	45	45	420	206	142	45	161	167	355	206	45	142	409
		5'-TCA ^{3'}	5'-TCg ^{3'}	5'-ggC ^{3'}	5'-CCg ^{3'}	5'-CCT ^{3'}	5'-CCg ^{3'}	5'-ggA ^{3'}	5'-ggA ^{3'}	5'-ggA ^{3'}	5'-TTA ^{3'}	5'-AgA ^{3'}	5'-TCT ^{3'}	5'-ggA ^{3'}	5'-Cgg ^{3'}	5'-gCT ^{3'}	5'-TCA ^{3'}	5'-gAC ^{3'}	5'-ggA ^{3'}	5'-TCT ^{3'}	5'-ggC ^{3'}
		363		363	103	103	418		357	435		420	419	540					368	368	420
		5'-AgC ^{3'}		5'-AgC ^{3'}	5'-CCg ^{3'}	5'-CCT ^{3'}	5'-Agg ^{3'}		5'-Tgg ^{3'}	5'-AAA ^{3'}		5'-TTA ^{3'}	5'-gTC ^{3'}	5'-gAC ^{3'}					5'-gTT ^{3'}	5'-gTC ^{3'}	5'-TTA ^{3'}
				363	361	363			412										557		
				5'-AgC ^{3'}	5'-AgT ^{3'}	5'-AAT ^{3'}			5'-ATA ^{3'}										5'-ggA ^{3'}		
3'-primer(s) ³		603	272	309	193	246	246	165	206	266	559	272	301	272	234	204	463	302	259	311	544
		5'-gTg ^{3'}	5'-Tgg ^{3'}	5'-gTg ^{3'}	5'-CgT ^{3'}	5'-TAT ^{3'}	5'-TAT ^{3'}	5'-Tgg ^{3'}	5'-CCT ^{3'}	5'-TCC ^{3'}	5'-CTC ^{3'}	5'-Tgg ^{3'}	5'-gTC ^{3'}	5'-TgC ^{3'}	5'-TCT ^{3'}	5'-TCT ^{3'}	5'-gCT ^{3'}	5'-ggC ^{3'}	5'-gTT ^{3'}	5'-ggg ^{3'}	5'-ggT ^{3'}
		605	272	559	193	246	572		538	559		572	301	309	272				259	311	
		5'-gCT ^{3'}	5'-TgA ^{3'}	5'-Cag ^{3'}	5'-CgT ^{3'}	5'-TAT ^{3'}	5'-gCC ^{3'}		5'-gTC ^{3'}	5'-Cag ^{3'}		5'-gCg ^{3'}	5'-gTC ^{3'}	5'-ATC ^{3'}	5'-TgA ^{3'}				5'-gTT ^{3'}	5'-ggg ^{3'}	
					583	559				583			570	605	272				266	538	
					5'-gTg ^{3'}	5'-CTC ^{3'}				5'-gTg ^{3'}			5'-CCg ^{3'}	5'-gCT ^{3'}	5'-Tgg ^{3'}				5'-TCC ^{3'}	5'-gTC ^{3'}	
													572						272		
													5'-gCg ^{3'}						5'-Tgg ^{3'}		
																			292		
																			5'-gTA ^{3'}		
																			572		
																			5'-gCg ^{3'}		
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41

101.708-24 – including *Tag* polymerase, IFU-01

101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: 43R

Lot-specific information

[illegible]

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41
HLA-B allele ⁴	ser. ¹⁴																				
*07:02:01-07:02:31, 07:05:01-07:06, 07:09-07:15, 07:17-07:18:02, 07:21-07:23, 07:26, 07:29-07:31, 07:34-07:35, 07:39, 07:41-07:46, 07:48-07:49N, 07:51-07:52, 07:55-07:59, 07:61-07:64, 07:66-07:68:03, 07:70-07:72, 07:74-07:76, 07:79-07:80, 07:82-07:83, 07:86-07:99, 07:101-07:121, 07:124-07:130, 07:132, 07:135N-07:138, 07:140-07:145, 07:147-07:148, 07:150-07:154	B7, B42, Null, -	22																			
*07:02:32	-	22				26															
*07:03, 07:08, 07:16, 07:32, 07:37	B703, -	22																			
*07:04, 07:25, 07:146	B7, -	22												34							
*07:07, 07:19, 07:33, 07:40, 07:53, 07:77, 07:122, 07:139	B7, -																				
*07:20, 07:24, 07:131	B7, -	22		24																	
*07:27	-	22										32								40	
*07:28	-	22				27															
*07:36, 07:81, 07:149	-	22																			
*07:38	-	22																			
*07:47	-	22																			
*07:50	-	22										32									
*07:54, 07:123	-	22																			
*07:60	-			24																	
*07:65	-	22																38			
*07:69, 07:85	-	22																			
*07:73	-	22														36					
*07:78	-										31										
*07:84	-																				
*07:100	-			24																	
*07:133	-	22																			
*07:134	-	22																38			
*08:01:01-08:01:17, 08:05, 08:08N, 08:10-08:11, 08:15, 08:18-08:19N, 08:22-08:24, 08:27, 08:29-08:31, 08:33-08:35, 08:39, 08:41-08:48, 08:51, 08:53, 08:56-08:59, 08:61, 08:63-08:64, 08:66-08:69, 08:71-08:73, 08:75, 08:80-08:83, 08:86N	B8, -, Null	22	23						29												
*08:02-08:03, 08:52, 08:78	B8, -	22	23						29												
*08:04	-	22	23						29			32									
*08:07, 08:14, 08:28, 08:37, 08:74, 08:77	B8, -	22	23																		
*08:09, 08:12:01-08:12:03, 08:16, 08:84	B8, -		23						29												
*08:13, 08:20, 08:40, 08:70, 08:79	-	22	23																		
*08:17, 08:54	-	22	23						29			32									
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.
																							ser. ¹⁴	HLA-B allele ⁴
		44			47																	64	B7, B42, Null, -	*07:02:01-07:02:31, 07:05:01-07:06, 07:09-07:15, 07:17-07:18:02, 07:21-07:23, 07:26, 07:29-07:31, 07:34-07:35, 07:39, 07:41-07:46, 07:48-07:49N, 07:51-07:52, 07:55-07:59, 07:61-07:64, 07:66-07:68:03, 07:70-07:72, 07:74-07:76, 07:79-07:80, 07:82-07:83, 07:86-07:99, 07:101-07:121, 07:124-07:130, 07:132, 07:135N-07:138, 07:140-07:145, 07:147-07:148, 07:150-07:154
		44			47																	64	-	*07:02:32
		44																				64	B703, -	*07:03, 07:08, 07:16, 07:32, 07:37
		44		46	47																	64	B7, -	*07:04, 07:25, 07:146
		44			47																	64	B7, -	*07:07, 07:19, 07:33, 07:40, 07:53, 07:77, 07:122, 07:139
		44			47																	64	B7, -	*07:20, 07:24, 07:131
																					63	-	*07:27	
		44			47																	64	-	*07:28
					47																	63	-	*07:36, 07:81, 07:149
																						63	-	*07:38
		44			w																	64	-	*07:47
		44																				64	-	*07:50
		44			47									56								64	-	*07:54, 07:123
		44			47																	64	-	*07:60
					47			50								58		w				64	-	*07:65
43	44																					64	-	*07:69, 07:85
	44				47																	64	-	*07:73
	44				47						53											64	-	*07:78
	44				47						53											64	-	*07:84
	44				47															62		64	-	*07:100
	44				47							55										64	-	*07:133
	44				47													w				64	-	*07:134
		44			46																	64	B8, -, Null	*08:01:01-08:01:17, 08:05, 08:08N, 08:10-08:11, 08:15, 08:18-08:19N, 08:22-08:24, 08:27, 08:29-08:31, 08:33-08:35, 08:39, 08:41-08:48, 08:51, 08:53, 08:56-08:59, 08:61, 08:63-08:64, 08:66-08:69, 08:71-08:73, 08:75, 08:80-08:83, 08:86N
					46																63		B8, -	*08:02-08:03, 08:52, 08:78
		44			46																	64	-	*08:04
		44			46																	64	B8, -	*08:07, 08:14, 08:28, 08:37, 08:74, 08:77
		44			46																	64	B8, -	*08:09, 08:12:01-08:12:03, 08:16, 08:84
		44																				64	-	*08:13, 08:20, 08:40, 08:70, 08:79
		44			46									56								64	-	*08:17, 08:54
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41
*08:21	-	22	23	24					29												
*08:25	-	22	23	24																	
*08:26, 08:50, 08:62, 08:85, 42:07 ⁵	-	22							29												
*08:32	-	22	23						29									38			
*08:36	-	22	23						29												
*08:38	-	22	23						29												
*08:49	-	22	23														37		40		
*08:55	-	22	23																	41	
*08:60, 08:76	-	22	23						29												
*08:65	-	22							29												
*08:87	-	22	23				27														
*08:88	-	22	23						29							36					
*13:01:01-13:01:06, 13:17, 13:20, 13:22:01-13:23, 13:25, 13:28-13:29, 13:43, 13:50, 13:52	B13, -			24	25						31						37				
*13:02:01-13:02:12, 13:08Q, 13:11, 13:15, 13:18-13:19, 13:27, 13:30, 13:32-13:34, 13:37-13:38, 13:40, 13:42, 13:44-13:45, 13:47, 13:49N, 13:54	B13, Null, -			24	25						31										
*13:03	-			24	25																
*13:04, 13:10, 13:21, 13:35, 44:135 ⁶	B15, B21, -			24	25																
*13:06	-			24	25												37				
*13:07N	Null			24	25						31						37				
*13:09	-				25						31										
*13:12-13:13, 13:26, 13:51	B13, -			24	25												37				
*13:14	-			24							31										
*13:16	-			24	25						31										
*13:31	-			24	25						31										
*13:36	-			24	25						31						37				
*13:39	-				25						31						37				
*13:41	-			24	25						31										
*13:46	-			24	25				29												
*13:48	-			24	25																
*13:53	-			24	25																
*14:01:01-14:01:02, 14:07N, 14:12, 14:14, 14:19, 14:26, 14:32	B64(14), Null, -				25	26															
*14:02:01-14:02:02, 14:02:04, 14:03-14:04, 14:09, 14:11, 14:15-14:18, 14:20, 14:22-14:23, 14:25, 14:27, 14:29, 14:31	B65(14), B14, -				25		27														
*14:02:03, 14:02:06, 14:24, 14:28	B65(14), -				25																
*14:02:05, 14:30	B65(14), -				25		27														
*14:05-14:06:02, 14:13	B14, -						27														
*14:08, 39:25N, 39:30, 39:32-39:34, 39:43, 39:47, 39:50 ⁷	B39(16), Null, -					26															
*14:10	-					26											37				
*14:21	-				25																
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.
		44		46																		64	-	*08:21
		44																				64	-	*08:25
		44		46																		64	-	*08:26, 08:50, 08:62, 08:85, 42:07 ⁵
		44		46															61			64	-	*08:32
				46																		64	-	*08:36
		44		46										56								64	-	*08:38
		44		46		48																64	-	*08:49
		44																				64	-	*08:55
		44		46		48																64	-	*08:60, 08:76
43	44			46																		64	-	*08:65
	44			46		48									57							64	-	*08:87
	44			46																		64	-	*08:88
														56							63	B13, -		*13:01:01-13:01:06, 13:17, 13:20, 13:22:01-13:23, 13:25, 13:28-13:29, 13:43, 13:50, 13:52
														56							63	B13, Null, -		*13:02:01-13:02:12, 13:08Q, 13:11, 13:15, 13:18-13:19, 13:27, 13:30, 13:32-13:34, 13:37-13:38, 13:40, 13:42, 13:44-13:45, 13:47, 13:49N, 13:54
														56		59					63	-		*13:03
														56							63	B15, B21, -		*13:04, 13:10, 13:21, 13:35, 44:135 ⁶
														56		58					63	-		*13:06
														56							63	Null		*13:07N
														56							64	-		*13:09
														56							63	B13, -		*13:12-13:13, 13:26, 13:51
														56							63	-		*13:14
											53			56							63	-		*13:16
											w			56						62	63	-		*13:31
								51													63	-		*13:36
														56							64	-		*13:39
														56						62	63	-		*13:41
														56							63	-		*13:46
											53			56		59					63	-		*13:48
														56		58					63	-		*13:53
43	44																				64	B64(14), Null, -		*14:01:01-14:01:02, 14:07N, 14:12, 14:14, 14:19, 14:26, 14:32
43	44																				64	B65(14), B14, -		*14:02:01-14:02:02, 14:02:04, 14:03-14:04, 14:09, 14:11, 14:15-14:18, 14:20, 14:22- 14:23, 14:25, 14:27, 14:29, 14:31
43	44																				64	B65(14), -		*14:02:03, 14:02:06, 14:24, 14:28
43																					64	B65(14), -		*14:02:05, 14:30
43	44																				64	B14, -		*14:05-14:06:02, 14:13
43	44																				64	B39(16), Null, -		*14:08, 39:25N, 39:30, 39:32- 39:34, 39:43, 39:47, 39:50 ⁷
43	44						49															64	-	*14:10
	44				47																	64	-	*14:21
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	
*15:01:01:01-15:01:01:02N, 15:01:03, 15:01:06-15:01:16, 15:01:18, 15:01:20-15:01:28, 15:26N-15:27:03, 15:30, 15:32, 15:34-15:35, 15:38:01- 15:38:02, 15:50, 15:56, 15:58, 15:60, 15:63, 15:70-15:71, 15:75, 15:79N, 15:81-15:82, 15:92, 15:94N, 15:96-15:97, 15:102, 15:104-15:105, 15:109- 15:111N, 15:113, 15:117- 15:118, 15:122, 15:125, 15:128-15:129, 15:135, 15:140, 15:145-15:147, 15:152, 15:159-15:160, 15:163- 15:164, 15:166-15:167, 15:169, 15:171-15:172, 15:174- 15:175, 15:178, 15:181N- 15:184, 15:187, 15:190N, 15:192-15:193, 15:201, 15:203, 15:205-15:206, 15:211, 15:217, 15:225, 15:227-15:228, 15:231-15:232, 15:234, 15:241, 15:244	B62(15), Null, -							28										39				
*15:01:02	B62(15)							28		30									39			
*15:01:04	B62(15)																		39			
*15:01:17, 15:05:01-15:06, 15:33, 15:39:01-15:40, 15:48, 15:65, 15:78:01-15:78:03, 15:84, 15:107, 15:136, 15:138, 15:141, 15:150, 15:155, 15:179, 15:185, 15:188, 15:195, 15:199, 15:218Q-15:219, 15:237	B62(15), B15, -							28											39			
*15:01:19, 15:28, 15:120, 15:149N	B62(15), Null, -							28											39			
*15:02:01-15:02:05, 15:88, 15:112, 15:121, 15:139, 15:170, 15:194, 15:213- 15:214, 15:223	B75(15), B15, -							28									37		39			
*15:03:01-15:03:02, 15:47:01- 15:47:02, 15:49, 15:61, 15:98, 15:103, 15:123, 15:127, 15:132, 15:151, 15:156, 15:158, 15:173, 15:210, 15:220, 15:235, 15:242	B72(70), B70, -								29			32							39			
*15:03:03	B72(70)							28	29			32							39			
*15:04, 15:137	B62(15), -							28											39			
*15:07:01-15:07:02, 15:45, 15:126, 15:207	B62(15), -			24				28											39			
*15:08, 15:11:01-15:11:03, 15:11:05, 15:15, 15:148, 15:191, 15:209N, 15:215	B75(15), B62(15), Null, -							28											39			
*15:09	B70								29	30				34					39			
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.
														56						62		64	B62(15), Null, -	*15:01:01:01-15:01:01:02N, 15:01:03, 15:01:06-15:01:16, 15:01:18, 15:01:20-15:01:28, 15:26N-15:27:03, 15:30, 15:32, 15:34-15:35, 15:38:01- 15:38:02, 15:50, 15:56, 15:58, 15:60, 15:63, 15:70-15:71, 15:75, 15:79N, 15:81-15:82, 15:92, 15:94N, 15:96-15:97, 15:102, 15:104-15:105, 15:109 15:111N, 15:113, 15:117- 15:118, 15:122, 15:125, 15:128-15:129, 15:135, 15:140, 15:145-15:147, 15:152, 15:159-15:160, 15:163 15:164, 15:166-15:167, 15:169, 15:171-15:172, 15:174 15:175, 15:178, 15:181N- 15:184, 15:187, 15:190N, 15:192-15:193, 15:201, 15:203, 15:205-15:206, 15:211, 15:217, 15:225, 15:227-15:228, 15:231-15:232, 15:234, 15:241, 15:244
														56						62		64	B62(15)	*15:01:02
														56						62		64	B62(15)	*15:01:04
														56								64	B62(15), B15, -	*15:01:17, 15:05:01-15:06, 15:33, 15:39:01-15:40, 15:48, 15:65, 15:78:01-15:78:03, 15:84, 15:107, 15:136, 15:138, 15:141, 15:150, 15:155, 15:179, 15:185, 15:188, 15:195, 15:199, 15:218Q-15:219, 15:237
																				62		64	B62(15), Null, -	*15:01:19, 15:28, 15:120, 15:149N
									51													64	B75(15), B15, -	*15:02:01-15:02:05, 15:88, 15:112, 15:121, 15:139, 15:170, 15:194, 15:213- 15:214, 15:223
														56								64	B72(70), B70, -	*15:03:01-15:03:02, 15:47:01- 15:47:02, 15:49, 15:61, 15:98, 15:103, 15:123, 15:127, 15:132, 15:151, 15:156, 15:158, 15:173, 15:210, 15:220, 15:235, 15:242
														56								64	B72(70)	*15:03:03
											w			56						62		64	B62(15), -	*15:04, 15:137
														56						62		64	B62(15), -	*15:07:01-15:07:02, 15:45, 15:126, 15:207
									51											62		64	B75(15), B62(15), Null, -	*15:08, 15:11:01-15:11:03, 15:11:05, 15:15, 15:148, 15:191, 15:209N, 15:215
																						64	B70	*15:09
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41
*15:10:01, 15:18:01-15:18:04, 15:37, 15:51-15:52, 15:72, 15:90, 15:93, 15:99, 15:108, 15:114, 15:119, 15:133-15:134, 15:153, 15:176, 15:197-15:198, 15:200, 15:221, 15:226N, 15:229, 15:238	B71(70), B70, B15, Null, -								29					34					39		
*15:10:02	B71(70)							28	29					34					39		
*15:11:04	B75(15)																		39		
*15:12, 15:14, 15:19	B76(15)							28					33						39		
*15:13:01-15:13:02	B77(15), -							28						34			37		39		
*15:16:01-15:16:03, 15:67, 15:95	B63(15), -									30				34					39		
*15:17:01-15:17:02, 15:162, 15:168, 15:177, 15:196, 15:208, 15:216, 15:230	B63(15), -							28		30				34					39		
*15:20, 15:25:01-15:25:03, 15:204, 15:240	B62(15), -							28									37		39		
*15:21, 15:44	B75(15), -							28						34			37		39		
*15:23, 15:115	-								29					34					39		
*15:24, 15:87, 15:157	B62(15), -							28						34					39		
*15:29	B15								29										39		
*15:31	B75(15)							28											39		
*15:36	-			w				28									37		39		
*15:42	-							28											39		
*15:43	-							28											39		
*15:46	B72(70)				25			28	29										39		
*15:53, 15:212	-				25				29										39		
*15:54	-								29			32							39		
*15:55	B15			24				28											39		
*15:57	-							28											39		
*15:62	-								29			32					37		39		
*15:64	-								29			32							39		
*15:66	-							28						34					39		
*15:68	B35			24								32									
*15:69	-								29			32							39		41
*15:73	B62(15)							28											39		
*15:74	-								29										39		
*15:76	-							28											39		
*15:77	-						27	28									37		39		
*15:80	B70								29					34			37		39		
*15:83	-							28	29										39		
*15:85, 15:154, 15:165	-							28									37		39		
*15:86, 15:224	-							28											39		
*15:89	-			w				28									37		39		
*15:91, 15:131	-								29			32	33						39		
*15:101	-							28													
*15:106	-				25			28	29								37		39		
*15:116	-							28											39		
*15:124	-								29					34					39		
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.
																						64	B71(70), B70, B15, Null, -	*15:10:01, 15:18:01-15:18:04, 15:37, 15:51-15:52, 15:72, 15:90, 15:93, 15:99, 15:108, 15:114, 15:119, 15:133- 15:134, 15:153, 15:176, 15:197-15:198, 15:200, 15:221, 15:226N, 15:229, 15:238
																						64	B71(70)	*15:10:02
									51											62		64	B75(15)	*15:11:04
														56						62		64	B76(15)	*15:12, 15:14, 15:19
									51												63		B77(15), -	*15:13:01-15:13:02
											53										63		B63(15), -	*15:16:01-15:16:03, 15:67, 15:95
																					63		B63(15), -	*15:17:01-15:17:02, 15:162, 15:168, 15:177, 15:196, 15:208, 15:216, 15:230
														56								64	B62(15), -	*15:20, 15:25:01-15:25:03, 15:204, 15:240
									51													64	B75(15), -	*15:21, 15:44
																					63		-	*15:23, 15:115
														56						62	63		B62(15), -	*15:24, 15:87, 15:157
																						64	B15	*15:29
									51													64	B75(15)	*15:31
														56							63		-	*15:36
											53			56		58						64	-	*15:42
														56						62	63		-	*15:43
													55	56						62		64	B72(70)	*15:46
													55	56						62		64	-	*15:53, 15:212
														56						62		64	-	*15:54
									51													64	B15	*15:55
							w							56						62		64	-	*15:57
														56								64	-	*15:62
																						64	-	*15:64
														56						62		64	-	*15:66
														56						62		64	B35	*15:68
														56								64	-	*15:69
														56		59				62		64	B62(15)	*15:73
														56								64	-	*15:74
					47				w											62		64	-	*15:76
														56						62		64	-	*15:77
																						64	B70	*15:80
														56								64	-	*15:83
														56						62		64	-	*15:85, 15:154, 15:165
														56		58						64	-	*15:86, 15:224
									51												63		-	*15:89
														56								64	-	*15:91, 15:131
						47								56						62		64	-	*15:101
													55	56								64	-	*15:106
			45											56								64	-	*15:116
			45																			64	-	*15:124
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41
*15:142	-		23					28										39			
*15:143	-							28	29									39			
*15:144	-							28									37	39			
*15:161	-								29				33	34				39			
*15:180	-		23																		
*15:186	-								29					34				39		41	
*15:189	-						27	28						34				39			
*15:202, 15:239	-							28										38	39		
*15:222	-									30				34				39			
*15:233	-						27	28										39			
*15:236	-							28													
*15:243	-								29	30				34				39			
*18:01:01:01-18:03, 18:05-18:08, 18:10-18:11, 18:13-18:15, 18:17N-18:18, 18:20-18:21, 18:23N-18:28, 18:30-18:36, 18:38-18:43, 18:45-18:47, 18:50-18:53, 18:55, 18:57-18:65, 18:68, 18:70-18:71	B18, Null, -														35			38			
*18:04	B18														35			38			
*18:09	B18														35						
*18:12	-														35			38			
*18:19	-														35			38			
*18:22, 18:69	-														35		37	38			
*18:29, 35:09:01-35:09:03, 35:18, 35:31-35:32:02, 35:37, 35:53N, 35:64, 35:68:01-35:68:02, 35:75, 35:88, 35:99, 35:118-35:119, 35:127, 35:151, 35:174 ⁸	B35, Null, -																	38			
*18:37	-											32			35						
*18:44	-				25										35			38			
*18:48	-														35						
*18:49	-														35						
*18:54	-												33		35					40	
*18:56	-												33		35			38			
*18:66	-														35						
*18:67	-														35						
*27:01	B27					26										36					
*27:02:01-27:02:02, 27:53, 27:57, 27:65N, 27:83	B27, Null, -					26										36					
*27:03-27:05:15, 27:05:17-27:06, 27:09-27:10, 27:13, 27:17, 27:20, 27:27, 27:31, 27:35-27:37, 27:41, 27:45-27:46, 27:48-27:51, 27:54-27:56, 27:58, 27:60-27:61, 27:64N, 27:66N, 27:68-27:69, 27:72-27:74, 27:78-27:80, 27:82, 27:84-27:88	B27, Null, -					26										36					
*27:05:16, 27:07:01-27:07:02, 27:11, 27:15, 27:21, 27:28, 27:32, 27:34, 27:38, 27:43, 27:47, 27:63, 27:70-27:71, 27:76	B27, -															36					
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.
													56							62	64	-		*15:142
													55							62	64	-		*15:143
									51											62	64	-		*15:144
																					64	-		*15:161
																				62	64	-		*15:180
																					64	-		*15:186
									51											62	64	-		*15:189
													56							62	64	-		*15:202, 15:239
																				63		-		*15:222
													56							62	64	-		*15:233
													56							62	64	-		*15:236
													56								64	-		*15:243
																			61		64	B18, Null, -		*18:01:01:01-18:03, 18:05-18:08, 18:10-18:11, 18:13-18:15, 18:17N-18:18, 18:20-18:21, 18:23N-18:28, 18:30-18:36, 18:38-18:43, 18:45-18:47, 18:50-18:53, 18:55, 18:57-18:65, 18:68, 18:70-18:71
																58			61		64	B18		*18:04
																			61	63		B18		*18:09
													56								64	-		*18:12
																			61	62	64	-		*18:19
							49												61		64	-		*18:22, 18:69
																			61		64	B35, Null, -		*18:29, 35:09:01-35:09:03, 35:18, 35:31-35:32:02, 35:37, 35:53N, 35:64, 35:68:01-35:68:02, 35:75, 35:88, 35:99, 35:118-35:119, 35:127, 35:151, 35:174 ⁸
																					64	-		*18:37
																			61		64	-		*18:44
													55								64	-		*18:48
																			61		64	-		*18:49
																			61	63		-		*18:54
																			61		64	-		*18:56
																					64	-		*18:66
													54						61	63		-		*18:67
					47			w						56						63		B27		*27:01
					47			50				54		56						63		B27, Null, -		*27:02:01-27:02:02, 27:53, 27:57, 27:65N, 27:83
																								*27:03-27:05:15, 27:05:17-27:06, 27:09-27:10, 27:13, 27:17, 27:20, 27:27, 27:31, 27:35-27:37, 27:41, 27:45-27:46, 27:48-27:51, 27:54-27:56, 27:58, 27:60-27:61, 27:64N, 27:66N, 27:68-27:69, 27:72-27:74, 27:78-27:80, 27:82, 27:84-27:88
					47			50						56						63		B27, Null, -		*27:05:16, 27:07:01-27:07:02, 27:11, 27:15, 27:21, 27:28, 27:32, 27:34, 27:38, 27:43, 27:47, 27:63, 27:70-27:71, 27:76
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41
*27:08, 27:26, 27:40, 27:42, 27:44, 27:89	B2708, -					26										36					
*27:12	B27					26										36					
*27:14, 27:81	B27, -															36					
*27:16, 27:39, 27:59N	Null, -					26										36					
*27:18	-					26						32				36					
*27:19	B27															36	37				
*27:23	-					26										36					
*27:24	-															36					
*27:25	-															36					
*27:29	B27					26						32				36					
*27:30	B27															36	37				
*27:33	-															36					
*27:52, 27:67	-					26															
*27:62	-															36					
*27:75	-					26															
*27:77	-					26										36					
*35:01:01-01-35:01:09, 35:01:11-35:03:02, 35:03:04-35:04:01, 35:04:03, 35:06-35:08:04, 35:11:01-35:12:03, 35:15, 35:20:01-35:20:02, 35:23, 35:29:01-35:29:02, 35:33-35:34, 35:36, 35:39-35:42:01, 35:48, 35:50, 35:52, 35:54-35:55, 35:57, 35:59, 35:61, 35:65Q, 35:70, 35:74, 35:77-35:78, 35:82-35:85, 35:90-35:95, 35:98, 35:100-35:101:02, 35:103-35:108:02, 35:110-35:112, 35:116, 35:120-35:124, 35:126, 35:128-35:134N, 35:136-35:150, 35:152-35:153, 35:155-35:156, 35:158-35:162, 35:165N-35:168, 35:170-35:173N, 35:175-35:184, 35:187, 35:189, 35:192-35:197, 35:200	B35, Null, -																37	38			
*35:01:10, 35:04:02	B35									30							37	38			
*35:03:03, 35:27, 35:56	B35, -																37				
*35:05:01-35:05:02, 35:22, 35:51, 35:58, 35:72, 35:89, 35:97, 35:114, 35:199	B35, -			24														38			
*35:08:05	-																37				
*35:10, 35:13, 35:28, 35:69, 35:80	B35, -																37	38			
*35:14:01-35:14:02, 35:62, 35:163	B35, -																37	38			
*35:16	B35			24													37	38			
*35:17, 35:30, 35:113, 35:125, 35:164	B35, -			24													37	38			
*35:19	B35																37				
*35:21, 35:24:01-35:24:02, 35:81, 35:96, 35:109, 35:157, 35:188, 35:190	-																37	38			
*35:25, 35:49, 35:198	B35, -																37				
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.
					47			50						56								64	B2708, -	*27:08, 27:26, 27:40, 27:42, 27:44, 27:89
														56								64	B27	*27:12
		44			47			50						56								63	B27, -	*27:14, 27:81
														56								63	Null, -	*27:16, 27:39, 27:59N
														56								64	-	*27:18
					47			50						56								63	B27	*27:19
																						63	-	*27:23
					47			50	51					56								63	-	*27:24
					47			50						56							62	63	-	*27:25
														56								63	B27	*27:29
					47			50				54		56								63	B27	*27:30
					47			50						56								64	-	*27:33
					47			50						56								63	-	*27:52, 27:67
					47			50				54		56								63	-	*27:62
					47			50				54		56							62	63	-	*27:75
												54		56								63	-	*27:77
																								*35:01:01:01-35:01:09, 35:01:11-35:03:02, 35:03:04-35:04:01, 35:04:03, 35:06-35:08:04, 35:11:01-35:12:03, 35:15, 35:20:01-35:20:02, 35:23, 35:29:01-35:29:02, 35:33-35:34, 35:36, 35:39-35:42:01, 35:48, 35:50, 35:52, 35:54-35:55, 35:57, 35:59, 35:61, 35:65Q, 35:70, 35:74, 35:77-35:78, 35:82-35:85, 35:90-35:95, 35:98, 35:100-35:101:02, 35:103-35:108:02, 35:110-35:112, 35:116, 35:120-35:124, 35:126, 35:128-35:134N, 35:136-35:150, 35:152-35:153, 35:155-35:156, 35:158-35:162, 35:165N-35:168, 35:170-35:173N, 35:175-35:184, 35:187, 35:189, 35:192-35:197, 35:200
																			61			64	B35, Null, -	*35:01:10, 35:04:02
																			61			64	B35, -	*35:03:03, 35:27, 35:56
																			61			64	B35, -	*35:05:01-35:05:02, 35:22, 35:51, 35:58, 35:72, 35:89, 35:97, 35:114, 35:199
																58						64	-	*35:08:05
														56								64	B35, -	*35:10, 35:13, 35:28, 35:69, 35:80
																			61	62		64	B35, -	*35:14:01-35:14:02, 35:62, 35:163
														56								64	B35	*35:16
																			61			64	B35, -	*35:17, 35:30, 35:113, 35:125, 35:164
												55										64	B35	*35:19
							49													61		64	-	*35:21, 35:24:01-35:24:02, 35:81, 35:96, 35:109, 35:157, 35:188, 35:190
																						64	B35, -	*35:25, 35:49, 35:198
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41
*35:26	-						27										37				
*35:35	B35																37	38			41
*35:38, 35:115, 35:169	-																37	38			
*35:42:02	-																37	38			
*35:43-35:44, 35:67, 35:79, 35:86, 35:102, 35:117, 35:135	B35, -																	38			
*35:45, 35:71	-											33					37	38			
*35:46	B35				25												37				
*35:47, 35:154	-				25												37				
*35:60	B35																	38			
*35:63	-				25												37				
*35:66	-	22																38			
*35:76	B35, 22																37	38			
*35:87	-	22						29										38			
*35:185	-																	38			
*35:186	-																37	38			
*35:191	-																37	38			
*37:01:01-37:01:07, 37:03N, 37:06, 37:13, 37:15-37:21, 37:23-37:27, 37:30N-37:33N	B37, Null, -											33					37			40	
*37:01:08, 37:09	B37, -											33							40		
*37:02	-					26						33							40		
*37:04:01-37:04:02	B37, -											33					37		40		
*37:05	-											33					37		40		
*37:07	-	22										33							40		
*37:08	B37											33					37		40		
*37:10, 37:22, 37:29	B37, -																37		40		
*37:11	-																37	38	40		
*37:12	-							29				33							40		
*37:14	-																37		40		
*37:28	-											32	33				37		40		
*38:01:01-38:02:02, 38:08-38:09, 38:11-38:16, 38:18-38:19, 38:21, 38:23-38:24, 38:27-38:29, 38:31-38:32, 38:34N, 38:38-38:40	B38(16), B16, Null, -					26														41	
*38:02:03	-																			41	
*38:03	B16					26						32								41	
*38:04, 38:25	-																			41	
*38:05, 38:33	B38(16), -						27													41	
*38:06	-																			41	
*38:07	-					26														41	
*38:10, 38:36-38:37	-					26														41	
*38:17	-					26						33							40	41	
*38:20	-					26											37			41	
*38:22	-					26										36					
*38:26	-					26														41	
*38:30	-					26		29													
*38:35	-					26														41	
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.
																						64	-	*35:26
																			61			64	B35	*35:35
						48													61			64	-	*35:38, 35:115, 35:169
																58			61			64	-	*35:42:02
																			61	62		64	B35, -	*35:43-35:44, 35:67, 35:79, 35:86, 35:102, 35:117, 35:135
						48													61			64	-	*35:45, 35:71
								51														64	B35	*35:46
													55									64	-	*35:47, 35:154
																58			61			64	B35	*35:60
													55	56								64	-	*35:63
																			61			64	-	*35:66
						47		50											w			64	B35, 22	*35:76
				46															61			64	-	*35:87
																58			61	62		64	-	*35:185
						47		50														64	-	*35:186
																						64	-	*35:191
						48								56							63	B37, Null, -	*37:01:01-37:01:07, 37:03N, 37:06, 37:13, 37:15-37:21, 37:23-37:27, 37:30N-37:33N	
														56							63	B37, -	*37:01:08, 37:09	
														56							63	-	*37:02	
						48	49							56							63	B37, -	*37:04:01-37:04:02	
						48								56								64	-	*37:05
														56							63	-	*37:07	
						48													61		63	B37	*37:08	
						48								56							63	B37, -	*37:10, 37:22, 37:29	
						48								56								64	-	*37:11
						48								56							63	-	*37:12	
						48								56							63	w	-	*37:14
						48								56							63	-	*37:28	
42	43																				63	B38(16), B16, Null, -	*38:01:01-38:02:02, 38:08-38:09, 38:11-38:16, 38:18-38:19, 38:21, 38:23-38:24, 38:27-38:29, 38:31-38:32, 38:34N, 38:38-38:40	
	43																				63	-	*38:02:03	
42														56							63	B16	*38:03	
	43													56							63	-	*38:04, 38:25	
42	43																				63	B38(16), -	*38:05, 38:33	
																			61		63	-	*38:06	
																			61		63	-	*38:07	
42	43																				63	-	*38:10, 38:36-38:37	
42	43																				63	-	*38:17	
42	43																				63	-	*38:20	
42	43																				63	-	*38:22	
42						47															63	-	*38:26	
42	43																				63	-	*38:30	
	43																				63	-	*38:35	
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41
*39:01:01-39:01:01:02L, 39:01:03-39:01:08, 39:01:10- 39:01:12, 39:03, 39:05:01- 39:07, 39:09, 39:11, 39:14- 39:15, 39:24:01, 39:26-39:29, 39:31, 39:35, 39:37-39:38Q, 39:40N-39:41, 39:44, 39:46, 39:48, 39:51-39:57, 39:59- 39:62, 39:64-39:68, 39:70- 39:71, 39:73	B39(16), B3901, Low B39, Null, -					26															41
*39:01:09, 39:12, 39:45, 39:72	B39(16), B3901, -																				41
*39:02:01, 39:08	B39(16), B3902					26						32									41
*39:02:02, 39:13:01-39:13:02, 39:23, 39:49	B39(16), B3902, -											32									41
*39:04	B39(16)						27														41
*39:10:01-39:10:02, 39:16, 39:20	B39(16), -																				41
*39:17, 39:63	-																				41
*39:18	-					26															41
*39:19:01	B39(16)																	38			41
*39:19:02	B39(16)					26												38			41
*39:22	-					26															41
*39:24:02	-					26															41
*39:36	-					26															
*39:39:01	-					26						32									41
*39:39:02	-											32									41
*39:42	-					26											37				41
*39:58	-					26															41
*39:69	-					26															41
*40:01:01-40:01:06, 40:01:08- 40:01:10, 40:01:12-40:01:17, 40:01:19-40:01:24, 40:07, 40:23, 40:51, 40:60, 40:66, 40:73, 40:79, 40:124:01, 40:139-40:140, 40:146-40:147, 40:150, 40:152-40:153, 40:155N, 40:172, 40:186	B60(40), -, Null				25																
*40:01:07, 40:01:18, 40:53, 40:72:01-40:72:02, 40:77	B60(40), B61(40), -				25																
*40:01:11	-				25																
*40:02:01-40:04, 40:06:05, 40:20, 40:39-40:40, 40:89, 40:94, 40:97, 40:122, 40:142N- 40:144N	B61(40), B41, Null, -																				
*40:05	B4005			24																	
*40:06:01-40:06:04, 40:44, 40:70, 40:75, 40:83, 40:86, 40:103, 40:127, 40:131, 40:159, 40:162, 40:165, 40:167, 40:177, 40:190	B61(40), -																				
*40:08	-																				
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.
	43	44																				64	B39(16), B3901, Low B39, Null, -	*39:01:01:01-39:01:01:02L, 39:01:03-39:01:08, 39:01:10- 39:01:12, 39:03, 39:05:01- 39:07, 39:09, 39:11, 39:14- 39:15, 39:24:01, 39:26-39:29, 39:31, 39:35, 39:37-39:38Q, 39:40N-39:41, 39:44, 39:46, 39:48, 39:51-39:57, 39:59- 39:62, 39:64-39:68, 39:70- 39:71, 39:73
	43	44																				64	B39(16), B3901, -	*39:01:09, 39:12, 39:45, 39:72
		44												56								64	B39(16), B3902	*39:02:01, 39:08
		44												56								64	B39(16), B3902, -	*39:02:02, 39:13:01-39:13:02, 39:23, 39:49
	43	44																				64	B39(16)	*39:04
		44																				64	B39(16), -	*39:10:01-39:10:02, 39:16, 39:20
		44														58						64	-	*39:17, 39:63
	43	44																		62		64	-	*39:18
	43	44																	61			64	B39(16)	*39:19:01
	43	44																	61			64	B39(16)	*39:19:02
	43	44												56								64	-	*39:22
	43																					64	-	*39:24:02
	43	44																		62		64	-	*39:36
		44																				64	-	*39:39:01
		44																				64	-	*39:39:02
	43	44																				64	-	*39:42
		44																				64	-	*39:58
	43	44												55								64	-	*39:69
		44	45											55	56							64	B60(40), -, Null	*40:01:01-40:01:06, 40:01:08- 40:01:10, 40:01:12-40:01:17, 40:01:19-40:01:24, 40:07, 40:23, 40:51, 40:60, 40:66, 40:73, 40:79, 40:124:01, 40:139-40:140, 40:146-40:147, 40:150, 40:152-40:153, 40:155N, 40:172, 40:186
		44												55	56							64	B60(40), B61(40), -	*40:01:07, 40:01:18, 40:53, 40:72:01-40:72:02, 40:77
		44	45											55	56	58						64	-	*40:01:11
														55	56							64	B61(40), B41, Null, -	*40:02:01-40:04, 40:06:05, 40:20, 40:39-40:40, 40:89, 40:94, 40:97, 40:122, 40:142N 40:144N
														55	56							64	B4005	*40:05
		44												55	56							64	B61(40), -	*40:06:01:01-40:06:04, 40:44, 40:70, 40:75, 40:83, 40:86, 40:103, 40:127, 40:131, 40:159, 40:162, 40:165, 40:167, 40:177, 40:190
														55								64	-	*40:08
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41
*40:09, 40:11:01-40:11:02, 40:18, 40:24, 40:27:01- 40:27:02, 40:29, 40:35, 40:37, 40:50, 40:56-40:57, 40:64, 40:78, 40:82, 40:85, 40:90- 40:91, 40:99, 40:104-40:105, 40:107, 40:111, 40:115, 40:119-40:120, 40:133Q, 40:145, 40:169, 40:173, 40:176, 40:180-40:181, 40:189	B61(40), B40, -																				
*40:10:01-40:10:02, 40:135	B60(40), -				25																
*40:12	-								29			32							39		
*40:13, 40:19	-																				
*40:14:01-40:14:03, 40:121, 40:129, 40:158, 40:164, 40:192	B60(40), -				25																
*40:15, 40:32	-	22			25																
*40:16	B61(40)	22			25																
*40:21	-																				
*40:22N, 40:33, 40:36, 40:38, 40:42-40:43, 40:49, 40:52, 40:54-40:55, 40:61-40:63, 40:65, 40:67, 40:81, 40:84, 40:87:01-40:88, 40:92, 40:100- 40:102, 40:108, 40:112- 40:114, 40:116, 40:118N, 40:123, 40:124:02-40:126, 40:128, 40:130, 40:134, 40:141, 40:151, 40:154, 40:156, 40:163, 40:168, 40:171, 40:175, 40:179, 40:182-40:183, 40:187, 40:191, 40:193	B60(40), Null, -				25																
*40:25, 40:69, 40:106	-				25																
*40:26	B21									30											
*40:28	-									30							37				
*40:30, 40:34, 40:160	B60(40), -	22			25																
*40:31, 40:45, 40:80	B60(40), -	22			25																
*40:46	-											32				36					
*40:47, 40:188	B40, -				25																
*40:48	B60(40)				25					31											
*40:58	-				25																
*40:59	-	22			25																
*40:68	-																				
*40:71	-			24																	
*40:74, 40:149, 40:170	-																				
*40:76	-				25																
*40:93	-											32				36					
*40:95	-																				
*40:96	-																				
*40:98	-	22																			
*40:109	-																				
*40:110	-				25																
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.	
		?											55	56								64	B61(40), B40, -	*40:09, 40:11:01-40:11:02, 40:18, 40:24, 40:27:01- 40:27:02, 40:29, 40:35, 40:37, 40:50, 40:56-40:57, 40:64, 40:78, 40:82, 40:85, 40:90- 40:91, 40:99, 40:104-40:105, 40:107, 40:111, 40:115, 40:119-40:120, 40:133Q, 40:145, 40:169, 40:173, 40:176, 40:180-40:181, 40:189	
			45										55	56								64	B60(40), -	*40:10:01-40:10:02, 40:135	
			45											56								64	-	*40:12	
										52		54		56							63		-	*40:13, 40:19	
		?											55	56								64	B60(40), -	*40:14:01-40:14:03, 40:121, 40:129, 40:158, 40:164, 40:192	
		?											55	56								64	-	*40:15, 40:32	
		44											55	56								64	B61(40)	*40:16	
		?	45											56								64	-	*40:21	
		?	45										55	56								64	B60(40), Null, -	*40:22N, 40:33, 40:36, 40:38, 40:42-40:43, 40:49, 40:52, 40:54-40:55, 40:61-40:63, 40:65, 40:67, 40:81, 40:84, 40:87:01-40:88, 40:92, 40:100 40:102, 40:108, 40:112- 40:114, 40:116, 40:118N, 40:123, 40:124:02-40:126, 40:128, 40:130, 40:134, 40:141, 40:151, 40:154, 40:156, 40:163, 40:168, 40:171, 40:175, 40:179, 40:182-40:183, 40:187, 40:191, 40:193	
		?	45										55									64	-	*40:25, 40:69, 40:106	
		?											55	56								64	B21	*40:26	
		?					49						55	56								64	-	*40:28	
		?	45										55	56								64	B60(40), -	*40:30, 40:34, 40:160	
		?							51				55	56								64	B60(40), -	*40:31, 40:45, 40:80	
		?	45											56								64	-	*40:46	
			45											56							63	B40, -	*40:47, 40:188		
		?	45										55	56								64	B60(40)	*40:48	
		?											55	56		58						64	-	*40:58	
		44	45										55	56								64	-	*40:59	
		?											55									64	-	*40:68	
		?											55	56		59						64	-	*40:71	
		?	45										55	56								64	-	*40:74, 40:149, 40:170	
			45											56									-	*40:76	
		44												56								64	-	*40:93	
		?									53		55	56								64	-	*40:95	
		44												56							63		-	*40:96	
		?											55	56								64	-	*40:98	
		44								52		54		56								63		-	*40:109
		44												56								63		-	*40:110
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.	

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41
*40:117	-				25																
*40:132, 40:178	-				25																
*40:136	-				25				29					34							
*40:137	-	22																			
*40:138	-																				
*40:148, 40:161	-																				
*40:157	-																				
*40:166	-				25																
*40:174	-			24																	
*40:184	-																				41
*40:185	-																	38			
*41:01, 41:03:01-41:03:02, 41:06-41:07, 41:12, 41:14- 41:17, 41:20	B41, -				25				29												
*41:02:01-41:02:04, 41:11, 41:13, 41:19	B41, -	22			25				29												
*41:04, 41:10, 41:18	-	22			25																
*41:05	-				25				29												
*41:08	-				25				29					34							
*41:09, 45:02-45:03 ⁹	-				25				29												
*42:01:01-42:02, 42:06, 42:10, 42:12, 42:15	B42, -	22							29												
*42:04, 42:14	-								29												
*42:05:01-42:05:02	-	22							29					34							
*42:08	-								29												
*42:09, 42:16	-	22																			
*42:11	-	22							29			32									
*42:13	-	22							29												
*44:02:01:01-44:02:21, 44:02:23-44:02:24, 44:11, 44:19N, 44:21-44:24, 44:27:01- 44:27:02, 44:33-44:34, 44:48, 44:52N, 44:55, 44:59, 44:63, 44:66-44:68, 44:71-44:74, 44:80, 44:84, 44:86-44:87, 44:89, 44:93, 44:101-44:102, 44:104, 44:112-44:113, 44:116, 44:118-44:119, 44:121, 44:126, 44:137, 44:142Q, 44:148-44:149N, 44:151-44:152	B44(12), Null, -				25												37				
*44:02:22, 44:12	B44(12), -				25												37				
*44:03:01-44:03:06, 44:03:08- 44:05:03, 44:13, 44:26, 44:28:01-44:30, 44:32, 44:35- 44:36, 44:38-44:39, 44:42, 44:45, 44:51, 44:56N, 44:58N, 44:61N, 44:69-44:70, 44:76- 44:79, 44:82, 44:85, 44:88, 44:92, 44:94, 44:96, 44:98, 44:105, 44:107-44:109, 44:114- 44:115, 44:117, 44:120, 44:122-44:125, 44:128, 44:133, 44:136, 44:139, 44:141, 44:143, 44:147	B44(12), B21, Null, -				25												37				
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.
			45							52		54		56							63		-	*40:117
		?	45			48							55	56								64	-	*40:132, 40:178
		?		46									55	56								64	-	*40:136
		?	45											56								64	-	*40:137
		44	45										55	56								64	-	*40:138
		44									53		55	56								64	-	*40:148, 40:161
														56							63		-	*40:157
		?	45			47							55									64	-	*40:166
		?											55	56								64	-	*40:174
		?											55	56								64	-	*40:184
		?	45											56								64	-	*40:185
		44		46									55	56								64	B41, -	*41:01, 41:03:01-41:03:02, 41:06-41:07, 41:12, 41:14-41:17, 41:20
		44		46									55	56								64	B41, -	*41:02:01-41:02:04, 41:11, 41:13, 41:19
		44		46									55	56								64	-	*41:04, 41:10, 41:18
		44		46									?	56								64	-	*41:05
		44		46									55	56								64	-	*41:08
		44											55	56								64	-	*41:09, 45:02-45:03 ⁹
		44		46	47																	64	B42, -	*42:01:01-42:02, 42:06, 42:10, 42:12, 42:15
		44		46	47																	64	-	*42:04, 42:14
		44		46	47																	64	-	*42:05:01-42:05:02
		44			47																	64	-	*42:08
		44		46	47																	64	-	*42:09, 42:16
		44		46										56								64	-	*42:11
		44		46	47	48																64	-	*42:13
						48								56	57							63	B44(12), Null, -	*44:02:01:01-44:02:21, 44:02:23-44:02:24, 44:11, 44:19N, 44:21-44:24, 44:27:01-44:27:02, 44:33-44:34, 44:48, 44:52N, 44:55, 44:59, 44:63, 44:66-44:68, 44:71-44:74, 44:80, 44:84, 44:86-44:87, 44:89, 44:93, 44:101-44:102, 44:104, 44:112-44:113, 44:116, 44:118-44:119, 44:121, 44:126, 44:137, 44:142Q, 44:148-44:149N, 44:151-44:152
						48									57							63	B44(12), -	*44:02:22, 44:12
						48								56								63	B44(12), B21, Null, -	*44:03:01-44:03:06, 44:03:08-44:05:03, 44:13, 44:26, 44:28:01-44:30, 44:32, 44:35-44:36, 44:38-44:39, 44:42, 44:45, 44:51, 44:56N, 44:58N, 44:61N, 44:69-44:70, 44:76-44:79, 44:82, 44:85, 44:88, 44:92, 44:94, 44:96, 44:98, 44:105, 44:107-44:109, 44:114-44:115, 44:117, 44:120, 44:122-44:125, 44:128, 44:133, 44:136, 44:139, 44:141, 44:143, 44:147
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45
*44:03:07, 44:07, 44:81, 44:103, 44:111, 44:134	B44(12), -																37								
*44:06	B44(12)																37								
*44:08	B44(12)			w													37								
*44:09	B45(12)				25												37								
*44:10	B44(12)				25							32					37								
*44:14	B44(12)				25				29								37								
*44:15	B12				25				29			32													
*44:16, 44:91, 44:132	B47, -				25		27										37								
*44:17, 44:43:01-44:43:02, 44:144	B44(12), -				25								33				37								
*44:18	-				25				29			32													
*44:20, 44:100	-				25				29																
*44:25, 44:50	-				25												37								
*44:31	B44(12)				25																			45	
*44:37:01-44:37:02, 44:64:01-44:64:02	-				25		27										37								
*44:40, 44:130	B44(12), -											32				36	37								
*44:41:01-44:41:02, 44:53, 44:99, 44:138Q	-, B44(12)				25																				
*44:44	-											32				36	37								
*44:46, 44:75, 44:129	-				25												37								
*44:47, 44:65	-				25																				
*44:49	-		23		25												37								
*44:54	B44(12)			24	25																				
*44:57	-			w													37								
*44:60, 44:110	-			24	25												37								
*44:62	-				25				29	30							37								
*44:83, 44:145	-																37								
*44:90	-				25												37								
*44:95	-				25												37								
*44:97	-				25												37								
*44:106	-			24	25																				
*44:127	-				25												37								
*44:131	-																37								
*44:140	-				25							32					37								
*44:146	-				25												37								
*44:150	-	22			25		27										37								
*45:01, 45:05, 45:07, 45:11-45:14	B45(12), -				25				29			32												44	
*45:04	-				25				29															44	
*45:06	-				25				29			32												44	
*45:08	-				25				29															44	
*45:09	-				25				29				33											44	
*45:10	-				25				29		31													44	
*46:01:01-46:01:07, 46:04-46:05, 46:07N, 46:09-46:10, 46:14-46:16, 46:20, 46:22-46:24, 46:27-46:30	B46, Null, -							28																	
*46:02	B46							28																	
*46:03, 46:08, 46:13:01-46:13:03, 46:19, 46:21:01-46:21:02, 46:25-46:26	B46, -							28																	
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.
						48								56							63		B44(12), -	*44:03:07, 44:07, 44:81, 44:103, 44:111, 44:134
						48						54			57				61		63		B44(12)	*44:06
						48								56	57						63		B44(12)	*44:08
						48							55	56	57							64	B45(12)	*44:09
						48								56			59				63		B44(12)	*44:10
						48								56							63		B44(12)	*44:14
														56			59				63		B12	*44:15
						48								56	57						63		B47, -	*44:16, 44:91, 44:132
						48								56							63		B44(12), -	*44:17, 44:43:01-44:43:02, 44:144
										52				56			59				63		-	*44:18
						48								56	57						63		-	*44:20, 44:100
						48				52		54		56							63		-	*44:25, 44:50
			45			48								56							63		B44(12)	*44:31
														56							63		-	*44:37:01-44:37:02, 44:64:01-44:64:02
						48								56							63		B44(12), -	*44:40, 44:130
						48								56	57						63		-, B44(12)	*44:41:01-44:41:02, 44:53, 44:99, 44:138Q
						48								56	57						63		-	*44:44
						48							55	56								64	-	*44:46, 44:75, 44:129
						48								56							63		-	*44:47, 44:65
						48								56	57						63		-	*44:49
						48								56							63		B44(12)	*44:54
						48								56							63		-	*44:57
						48								56							63		-	*44:60, 44:110
						48								56							63		-	*44:62
						48								56	57						63		-	*44:83, 44:145
					47	48		50					55	56	57							64	-	*44:90
						48				52		54		56	57						63		-	*44:95
					47	48		50						56	57						63		-	*44:97
						48								56	57						63		-	*44:106
						48								56	57								-	*44:127
						48							55	56	57							64	-	*44:131
						48								56			59				63		-	*44:140
						48								56						62	63		-	*44:146
														56							63		-	*44:150
		44											55	56			59					64	B45(12), -	*45:01, 45:05, 45:07, 45:11-45:14
		44											55	56			59					64	-	*45:04
		44			47								55	56			59					64	-	*45:06
		44											55	56		58						64	-	*45:08
		44											55	56								64	-	*45:09
		44											55	56								64	-	*45:10
					47		49													62		64	B46, Null, -	*46:01:01-46:01:07, 46:04-46:05, 46:07N, 46:09-46:10, 46:14-46:16, 46:20, 46:22-46:24, 46:27-46:30
					w		49													62		64	B46	*46:02
					47		49															64	B46, -	*46:03, 46:08, 46:13:01-46:13:03, 46:19, 46:21:01-46:21:02, 46:25-46:26
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41
*46:06	-							28													
*46:11	-							28													
*46:12	-			24				28													
*46:17	-							28				33									
*46:18	-							28													
*47:01:01:01-47:01:02, 47:06-47:08	B47, -				25																
*47:02	B47				25																
*47:03	B47				25																
*47:04-47:05	-				25	26															
*48:01:01-48:01:03, 48:07, 48:09, 48:11, 48:15, 48:19-48:20, 48:22, 48:27	B48, -	22										32									
*48:02:01-48:02:02, 48:25	B48, -											32					37				
*48:03:01-48:03:02, 48:23	B48, -											32									
*48:04, 48:24	B48, -											32									
*48:05, 48:08, 48:10, 48:12, 48:14	B48, -	22										32									
*48:06	-	22																			
*48:13	-											32									
*48:16	-	22										32									
*48:17	B48, 40	22										32									
*48:18	-	22										32									
*48:21, 48:26	-					26						32									
*49:01:01, 49:01:03, 49:06, 49:08-49:10, 49:12-49:17, 49:19N	B49(21), Null, -				25																
*49:01:02, 49:11	B49(21), -				25																
*49:02	B49(21)				25																
*49:03	-																				
*49:04-49:05	B49(21), -				25																
*49:07	-				25							31									
*49:18	-				25					30											
*49:20	-				25							32									
*50:01:01-50:01:03, 50:04-50:05, 50:08, 50:10-50:11, 50:13, 50:16	B50(21), -				25																
*50:02	B45(12)				25							32									
*50:06	-				25																
*50:07	-																				
*50:09	-				25																
*50:12	-				25																
*50:14	-				25					30											
*50:15	-				25																
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.
			45		47		49															64	-	*46:06
					47		49										59					64	-	*46:11
					47		49													62		64	-	*46:12
					47		49													62		64	-	*46:17
					47		49										58					64	-	*46:18
								50						56							63		B47, -	*47:01:01-47:01:02, 47:06-47:08
								50					55	56								64	B47	*47:02
								50					w	56							63	w	B47	*47:03
														56							63		-	*47:04-47:05
		44							51					56								64	B48, -	*48:01:01-48:01:03, 48:07, 48:09, 48:11, 48:15, 48:19-48:20, 48:22, 48:27
		44												56								64	B48, -	*48:02:01-48:02:02, 48:25
		44	45											56								64	B48, -	*48:03:01-48:03:02, 48:23
		44							51					56								64	B48, -	*48:04, 48:24
		44												56								64	B48, -	*48:05, 48:08, 48:10, 48:12, 48:14
		44							51													64	-	*48:06
		44												56								64	-	*48:13
		44							51													64	-	*48:16
		44	45											56								64	B48, 40	*48:17
									51					56							63		-	*48:18
		44												56								64	-	*48:21, 48:26
										52	53			56			59				63		B49(21), Null, -	*49:01:01, 49:01:03, 49:06, 49:08-49:10, 49:12-49:17, 49:19N
										52				56			59				63		B49(21), -	*49:01:02, 49:11
											53			56			59				63		B49(21)	*49:02
											53			56			59				63		-	*49:03
										52	53			56							63		B49(21), -	*49:04-49:05
										52	53			56							63		-	*49:07
										52	53			56							63		-	*49:18
										52	53			56			59				63		-	*49:20
		44									53		55	56			59					64	B50(21), -	*50:01:01-50:01:03, 50:04-50:05, 50:08, 50:10-50:11, 50:13, 50:16
		44									53		55	56			59					64	B45(12)	*50:02
		44									53			56			59					64	-	*50:06
		44									53		55	56			59					64	-	*50:07
		44									53		55	56								64	-	*50:09
		44											55	56			59					64	-	*50:12
		44									53		55	56								64	-	*50:14
		44									53		55				59					64	-	*50:15
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41
*51:01:01-51:01:03, 51:01:05-51:01:08, 51:01:10-51:01:20, 51:01:22-51:01:24, 51:01:26-51:02:03, 51:02:05-51:03, 51:05, 51:07:01-51:07:02, 51:09:01-51:09:02, 51:11N-51:12, 51:13:02-51:14, 51:16-51:19, 51:21-51:23, 51:26-51:33, 51:37-51:41N, 51:43, 51:48-51:53, 51:55, 51:57-51:58, 51:60, 51:63, 51:65-51:67, 51:69-51:77, 51:79-51:80, 51:82-51:91, 51:94-51:96, 51:98N, 51:100, 51:102-51:103, 51:105, 51:107-51:111, 51:113-51:117, 51:119-51:130	B51(5), B5102, B5103, Null, -									30											
*51:01:04, 51:02:04, 51:13:01, 51:24:01, 51:34-51:35, 51:92-51:93, 51:99, 51:131	B51(5), B5102, -																				
*51:01:09, 51:01:21, 51:01:25, 51:24:03-51:24:04	B51(5), -									30											
*51:04, 51:46, 51:56:01-51:56:02	B51(5), -									30							37				
*51:06:01-51:06:02, 51:59	B51(5), -									30											
*51:08, 51:20, 51:44N, 51:97	B51(5), Null, -								29	30											
*51:10	-																				
*51:15, 51:62	-																				
*51:24:02	B51(5)																				
*51:36	-								29												
*51:42	B44(12)																37				
*51:45, 51:81	-																				
*51:54, 51:78:01-51:78:02	-									30											
*51:61	-									30											
*51:64	-			24						30											
*51:68	-		23							30											
*51:101	-									30											41
*51:104, 51:118N, 58:08:01-58:08:02 ¹⁰	B17, Null, -									30											
*51:106	-									30											
*51:112	-									30											
*52:01:01:01-52:01:04, 52:01:06-52:01:09, 52:01:11-52:13, 52:15, 52:17-52:18, 52:22-52:24, 52:27-52:29	B52(5), -									30											
*52:01:05, 52:14	B52(5), -																				
*52:01:10	-									30											
*52:16	-									30		32									
*52:19	-								29	30											
*52:20	-									30											
*52:21	-									30											
*52:25	-																				
*52:26	-							27		30											
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	Well No.
											53	54							61		63		B51(5), B5102, B5103, Null, - *51:01:01-51:01:03, 51:01:05- 51:01:08, 51:01:10-51:01:20, 51:01:22-51:01:24, 51:01:26- 51:02:03, 51:02:05-51:03, 51:05, 51:07:01-51:07:02, 51:09:01-51:09:02, 51:11N- 51:12, 51:13:02-51:14, 51:16- 51:19, 51:21-51:23, 51:26- 51:33, 51:37-51:41N, 51:43, 51:48-51:53, 51:55, 51:57- 51:58, 51:60, 51:63, 51:65- 51:67, 51:69-51:77, 51:79- 51:80, 51:82-51:91, 51:94- 51:96, 51:98N, 51:100, 51:102 51:103, 51:105, 51:107- 51:111, 51:113-51:117, 51:119 51:130
											53	54							61		63		B51(5), B5102, - *51:01:04, 51:02:04, 51:13:01, 51:24:01, 51:34-51:35, 51:92- 51:93, 51:99, 51:131
		44									53	54							61		63		B51(5), - *51:01:09, 51:01:21, 51:01:25, 51:24:03-51:24:04
							49					54							61		63		B51(5), - *51:04, 51:46, 51:56:01- 51:56:02
												54							61		63		B51(5), - *51:06:01-51:06:02, 51:59
											53	54							61		63		B51(5), Null, - *51:08, 51:20, 51:44N, 51:97
		44										54							61		63		- *51:10
												54				59			61		63		- *51:15, 51:62
		44									53	54							61		63		B51(5) *51:24:02
											53	54							61		63		- *51:36
						48						54							61		63		B44(12) *51:42
												54							61		63		- *51:45, 51:81
											53								61		63		- *51:54, 51:78:01-51:78:02
											w	54							61	62	63		- *51:61
												54							61		63		- *51:64
											53	54							61		63		- *51:68
											53	54							61		63		- *51:101
											53	54								63			B17, Null, - *51:104, 51:118N, 58:08:01- 58:08:02 ¹⁰
											53	54				59			61		63		- *51:106
										52	53	54								63			- *51:112
											53	54		56						63			B52(5), - *52:01:01:01-52:01:04, 52:01:06-52:01:09, 52:01:11- 52:13, 52:15, 52:17-52:18, 52:22-52:24, 52:27-52:29
											53	54		56						63			B52(5), - *52:01:05, 52:14
42											53	54		56						63			- *52:01:10
											53	54		56						63			- *52:16
											53	54		56						63			- *52:19
											53			56						63			- *52:20
											53	54		56					62	63			- *52:21
											53	54		56		59				63			- *52:25
											53	54		56						63			- *52:26
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	Well No.

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41
*53:01:01-53:01:07, 53:04-53:05, 53:07-53:08:02, 53:10, 53:16, 53:18-53:21, 53:23-53:27	B53, -																37				
*53:02, 53:06	-																37				
*53:03, 53:09, 53:11-53:13	B53, -																37				
*53:14	-			24																	
*53:15	-	22															37				
*53:17:01-53:17:02	-																37				
*53:22	-											33					37				
*54:01:01, 54:02, 54:05N, 54:07-54:08N, 54:10, 54:13, 54:16-54:19, 54:21-54:25	B54(22), Null, -																				
*54:01:02, 55:01:07, 55:02:01-55:02:06, 55:07, 55:10, 55:12, 55:16, 55:19, 55:26, 55:30, 55:35, 55:37, 55:39, 55:41-55:43, 55:47-55:48, 55:50, 56:10 ¹¹	B55(22), B22, -																				
*54:03	-																				
*54:04, 54:09, 54:14-54:15	-																				
*54:06	-																				
*54:11	B54(22)																				
*54:12	-																				
*54:20	-																				
*54:26	-										31										
*55:01:01-55:01:06, 55:01:08-55:01:09, 55:03, 55:05, 55:11, 55:15, 55:25, 55:29, 55:31, 55:33, 55:36, 55:38, 55:44-55:45, 55:52-55:55N	B55(22), -, B22, Null																				
*55:04, 55:08, 55:13, 55:23, 55:27, 55:32, 55:46, 55:49, 56:15, 56:18-56:19N, 56:22, 56:31-56:32 ¹²	B55(22), B56(22), Null, -																				
*55:09, 55:24	B22, -										31										
*55:14	-																37				
*55:17, 55:28	-																				
*55:18	-											32									
*55:20	-								29												
*55:21	-																				
*55:22	B55(22)										31										
*55:34	-																				
*55:40	-																				
*55:51	-			24																	
*55:56	-								29												
*56:01:01-56:02, 56:04, 56:08, 56:14, 56:16, 56:20:01, 56:24, 56:26-56:30, 56:34-56:35	B56(22), Null, -																				
*56:03	B22																				
*56:05:01	B56(22)									30											
*56:05:02	B56(22)									30											
*56:06	B78									30								38			
*56:07	B56(22)																				
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.
												54							61		63		B53, -	*53:01:01-53:01:07, 53:04-53:05, 53:07-53:08:02, 53:10, 53:16, 53:18-53:21, 53:23-53:27
							49					54							61		63		-	*53:02, 53:06
																			61		63		B53, -	*53:03, 53:09, 53:11-53:13
												54							61		63		-	*53:14
												54							61		63		-	*53:15
												54		56							63		-	*53:17:01-53:17:02
							48					54							61		63		-	*53:22
					47			50							57	58						64	B54(22), Null, -	*54:01:01, 54:02, 54:05N, 54:07-54:08N, 54:10, 54:13, 54:16-54:19, 54:21-54:25
					47			50								58						64	B55(22), B22, -	*54:01:02, 55:01:07, 55:02:01-55:02:06, 55:07, 55:10, 55:12, 55:16, 55:19, 55:26, 55:30, 55:35, 55:37, 55:39, 55:41-55:43, 55:47-55:48, 55:50, 56:10 ¹¹
					47			50							57	59						64	-	*54:03
					47			50							57							64	-	*54:04, 54:09, 54:14-54:15
					47			50							57						62	64	-	*54:06
	44				47			50							57							64	B54(22)	*54:11
					47			50		52					57	58						63	-	*54:12
					47			50			53				57	58						64	-	*54:20
					47			50							57							64	-	*54:26
					47			50			53					58						64	B55(22), -, B22, Null	*55:01:01-55:01:06, 55:01:08-55:01:09, 55:03, 55:05, 55:11, 55:15, 55:25, 55:29, 55:31, 55:33, 55:36, 55:38, 55:44-55:45, 55:52-55:55N
					47			50														64	B55(22), B56(22), Null, -	*55:04, 55:08, 55:13, 55:23, 55:27, 55:32, 55:46, 55:49, 56:15, 56:18-56:19N, 56:22, 56:31-56:32 ¹²
					47			50			53											64	B22, -	*55:09, 55:24
					47			50									60					64	-	*55:14
					47			50			53											64	-	*55:17, 55:28
														56		58						64	-	*55:18
				46	47			50														64	-	*55:20
					47			50			w					58					62	64	-	*55:21
					47			50														64	B55(22)	*55:22
					47			50						56		58						64	-	*55:34
					47			50			53			56		58						64	-	*55:40
					47			50														64	-	*55:51
					47			50								58						64	-	*55:56
					47			50								59						64	B56(22), Null, -	*56:01:01-56:02, 56:04, 56:08, 56:14, 56:16, 56:20:01, 56:24, 56:26-56:30, 56:34-56:35
					47			50													62	64	B22	*56:03
					47			50			53					58						64	B56(22)	*56:05:01
					47			50			53											64	B56(22)	*56:05:02
					47			50			53					58			w			64	B78	*56:06
					47			50								59					63		B56(22)	*56:07
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41
*56:09, 56:11-56:12	B56(22), B55(22), -																37				
*56:13	B56(22)								29												
*56:17, 56:33	-																				
*56:20:02	-																				
*56:21	-									30											
*56:23	-																				
*56:25	-																				
*57:01:01-57:01:04, 57:01:06- 57:03:02, 57:05-57:06, 57:08, 57:10, 57:15-57:20, 57:22- 57:23, 57:25-57:30, 57:32- 57:44, 57:46, 57:48-57:50, 57:52-57:57, 58:36 ¹³	B57(17), Null, -																37				
*57:01:05, 57:11, 57:13, 57:21, 57:31, 57:47	B57(17), -																				
*57:04	B57(17)						27										37				
*57:07, 57:24	-																37				
*57:09	-								29								37				
*57:12	-																37				
*57:14	-																37				
*57:45, 57:51	-				25												37				
*58:01:01-58:01:02, 58:01:04- 58:01:11, 58:04-58:05, 58:10N- 58:15, 58:19, 58:21-58:24, 58:28-58:35	B58(17), Null, -																37				
*58:01:03, 58:02, 58:06, 58:16, 58:25-58:27	B58(17), -																				
*58:07	-											33									
*58:09	-																37				
*58:17N	Null																37				
*58:18	-			24																	
*58:20	-																				41
*59:01:01-59:01:01:02, 59:05	B59, -																				
*59:02-59:03	B59, -																				
*59:04	-																				
*67:01:01, 67:03	B67, -					26															41
*67:01:02-67:02	B67, -																				41
*73:01-73:02	B73, -																				
*78:01:01-78:01:02, 78:02:02- 78:03, 78:07	B78, -									30								38			
*78:02:01, 78:04	B78, -									30								38			
*78:05	-									30								38			
*78:06	-									30											
*81:01	B81	22																			
*81:02	B81	22																			
*81:03-81:05	Null, -	22																			
*82:01-82:03	B82, -							29			32							39			
*83:01	-																37				
A*23:31, A*24:106, C*07:231, C*16:10						27															
A*24:174												33									
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	Well No.
					47			50														64	B56(22), B55(22), - *56:09, 56:11-56:12
					47			50									59					64	B56(22) *56:13
																	59					64	- *56:17, 56:33
		44			47			50									59					64	- *56:20:02
					47			50		52	53											63	- *56:21
																58						64	- *56:23
					47			50			53						59					64	- *56:25
												54						60				63	B57(17), Null, - *57:01:01-57:01:04, 57:01:06- 57:03:02, 57:05-57:06, 57:08, 57:10, 57:15-57:20, 57:22- 57:23, 57:25-57:30, 57:32- 57:44, 57:46, 57:48-57:50, 57:52-57:57, 58:36 ¹³
												54						60				63	B57(17), - *57:01:05, 57:11, 57:13, 57:21, 57:31, 57:47
												54						60				63	B57(17) *57:04
						48						54						60				63	- *57:07, 57:24
						48						54						60				63	- *57:09
																		60				64	- *57:12
							49					54						60				63	- *57:14
										52		54										63	- *57:45, 57:51
												54										63	B58(17), Null, - *58:01:01-58:01:02, 58:01:04- 58:01:11, 58:04-58:05, 58:10N 58:15, 58:19, 58:21-58:24, 58:28-58:35
												54										63	B58(17), - *58:01:03, 58:02, 58:06, 58:16, 58:25-58:27
												54										63	- *58:07
							49					54										63	- *58:09
																						63	Null *58:17N
												54										63	- *58:18
												54										63	- *58:20
										52						58						63	B59, - *59:01:01:01-59:01:01:02, 59:05
										52												63	B59, - *59:02-59:03
										52							59					63	- *59:04
		44			47																	64	B67, - *67:01:01, 67:03
		44			47																	64	B67, - *67:01:02-67:02
		44			47																62	64	B73, - *73:01-73:02
											53					58			61			64	B78, - *78:01:01-78:01:02, 78:02:02- 78:03, 78:07
											53								61			64	B78, - *78:02:01, 78:04
											53			56								64	- *78:05
											53			56								64	- *78:06
					47			50	51													64	B81 *81:01
		44			47				51													64	B81 *81:02
					47				51													64	Null, - *81:03-81:05
					47			50									59					64	B82, - *82:01-82:03
					47	48		50							57							64	- *83:01
																							A*23:31, A*24:106, C*07:231, C*16:10 A*24:174
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	Well No.

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41
A*26:68, A*68:56, C*02:56, C*06:20, C*12:50																					
C*01:30					25																
C*03:05, C*03:25, C*03:27, C*03:143				24																	
C*03:12, C*03:19																					
C*03:102																	37				
C*03:129																					
C*06:72																					
C*07:02:30, C*08:16:02								28													
C*07:46									29												
C*15:02:04																					
C*15:25																					
C*15:39																					
C*15:51																					41
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41

¹The internal positive control primer pairs amplify segments of the human growth hormone gene. The two different control primer pairs give rise to either an internal positive control band of 1070 base pairs, for most wells, or a band of 800 base pairs, for some wells. Well number 22 contains the primer pair giving rise to the shorter, 800 bp, internal positive control band.

In addition, wells number 25, 26, 27, 31, 37, 40, 41, 47, 49, 53 and 60 contain the primer pair giving rise to the shorter, 800 bp, internal positive control band.

²The nucleotide position, in the 1st, 2nd or 3rd exon or in the 1st intron, matching the specificity-determining 3'-end of the primer is given. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

³The nucleotide position, in the 2nd or 3rd exon, matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

⁴The sequence of the B*0701 allele has been shown to be in error.

The B*08:06 allele has been shown to be identical to B*08:20.

The sequence of the B*13:05 allele has been shown to be identical to B*13:04.

The B*13:24 allele has been renamed B*13:22:02

The B*15:01:05 allele has been corrected and renamed B*15:120

The B*15:22 allele has been renamed B*35:43.

The sequence of the B*15:41 allele has been shown to be identical to B*15:39.

The B*15:59 allele has been renamed B*35:44.

The sequence of the B*18:16 allele has been shown to be identical to B*18:14.

The sequence of the B*27:051 allele has been shown to be identical to B*27:05:02.

The sequence of the B*27:22 allele has been shown to be identical to the corrected B*27:06 sequence.

The B*35:43:02 allele has been renamed B*35:185.

The B*35:73 allele has been renamed B*35:08:03.

The sequence of the B*39012 allele has been shown to be identical to B*39:01:01:01.

The sequence of the B*3921 allele has been shown to be identical to B*39:24.

The sequence of the B*4017 allele has been shown to be identical to B*40:16.

The sequence of the B*4041 allele has been shown to be identical to B*40:40.

The sequence of the B*4203 allele has never been assigned.

The sequence of the B*4401 allele has been shown to be identical to B*44:02:01:01.

The sequence of the B*5003 allele has been shown to be identical to B*50:02.

The sequence of the B*5125 allele has been shown to be identical to B*51:22.

The B*5147 allele has been renamed B*51:09:02.

The sequence of the B*5506 allele has been shown to be identical to B*55:04. The sequence of the B*5803 allele has never been assigned.

101.708-24 – including *Taq* polymerase, IFU-01
101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

Lot No. 101								Lot Specific Information																		Well No.
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		A*26:68, A*68:56, C*02:56, C*06:20, C*12:50 C*01:30 C*03:05, C*03:25, C*03:27, C*03:143 C*03:12, C*03:19 C*03:102 C*03:129 C*06:72 C*07:02:30, C*08:16:02 C*07:46 C*15:02:04 C*15:25 C*15:39 C*15:51		
																				62						
																	59									
																	59									
			45																							
																		60								
																	58									
						48																				
							49																			
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.		

The B*7901 allele has been renamed B*15:18:01.

The B*9530 allele has been renamed B*15:27:02.

⁵The B*08:26, 08:50, 08:62 and 08:85 and the B*42:07 alleles give rise to identical amplification patterns with the HLA-B low resolution primer set. These alleles can be separated by the respective high resolution primer sets.

⁶The B*13:04, 13:10, 13:21 and 13:35 and the B*44:135 alleles give rise to identical amplification patterns with the HLA-B low resolution primer set. These alleles can be separated by the respective high resolution primer sets.

⁷The B*14:08 and the B*39:25N, 39:30, 39:32-39:34, 39:43, 39:47 and 39:50 alleles give rise to identical amplification patterns with the HLA-B low resolution primer set. These alleles can be separated by the respective high resolution primer sets.

⁸The B*18:29 and the B*35:09:01-35:09:03, 35:18, 35:31-35:32:02, 35:37, 35:53N, 35:64, 35:68:01-35:68:02, 35:75, 35:88, 35:99, 35:118-35:119, 35:127, 35:151 and 35:174 alleles give rise to identical amplification patterns with the HLA-B low resolution primer set. These alleles can be separated by the respective high resolution primer sets.

⁹The B*41:09 and the B*45:02-45:03 alleles give rise to identical amplification patterns with the HLA-B low resolution primer set. These alleles can be separated by the respective high resolution primer sets.

¹⁰The B*51:104 and 51:118N and the B*58:08:01-58:08:02 alleles give rise to identical amplification patterns with the HLA-B low resolution primer set. These alleles can be separated by the respective high resolution primer sets.

¹¹The B*54:01:02 and the B*55:01:07, 55:02:01-55:02:06, 55:07, 55:10, 55:12, 55:16, 55:19, 55:26, 55:30, 55:35, 55:37, 55:39, 55:41-55:43, 55:47-55:48, 55:50 and the B*56:10 alleles give rise to identical amplification patterns with the HLA-B low resolution primer set. These alleles can be separated by the respective high resolution primer sets.

¹²The B*55:04, 55:08, 55:13, 55:23, 55:27, 55:32, 55:46 and 55:49 and the B*56:15, 56:18-56:19N, 56:22 and 56:31-56:32 alleles give rise to identical amplification patterns with the HLA-B low resolution primer set. These alleles can be separated by the respective high resolution primer sets.

¹³The B*57:01:01-57:01:04, 57:01:06-57:03:02, 57:05-57:06, 57:08, 57:10, 57:15-57:20, 57:22-57:23, 57:25-57:30, 57:32-57:44, 57:46, 57:48-57:50 and 57:52-57:57 and the B*58:36 alleles give rise to identical amplification patterns with the HLA-B low resolution primer set. These alleles can be separated by the respective high resolution primer sets.

¹⁴The grouping of not serologically defined alleles is taken from Tissue Antigens 73, 95-170, 2009.

‘ser’, serological HLA specificity.

‘w’, might be weakly amplified.

‘?’, nucleotide sequence information not available for the primer matching sequence.

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

INTERPRETATION TABLE																									
DR low resolution SSP typing																									
Amplification patterns of the DRB1*01:01 to DRB1*10:04 alleles																									
		Well																							
		65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	
Length of spec.		255	205											170	100			175							
PCR product(s)			200											105	85										
		215	200											215				225	195						
		210																							
		220	120																						
		210	80																						
		210	85																						
		175	100																						
		235	210											235	210										
		250	215	170										250	215	170									
		430	180	135	90									430	180	135	90								
		430											430												
Length of int. pos. control ¹		515											515												
5'-primer(s) ²		5'-gAA 3'	14(129)	5'-A-T 3'	13(124)									5'-gTC 3'	13(125)			5'-gTC 3'	13(125)						
				5'-gAA 3'	14(129)													5'-gTC 3'	13(125)						
						5'-gTC 3'	13(125)							5'-gTC 3'	13(125)										
								5'-gTC 3'	13(125)					5'-gTC 3'	13(125)										
										5'-gTC 3'	13(125)			5'-gTC 3'	13(125)										
												5'-gTC 3'	13(125)												
														5'-gTC 3'	13(125)										
																5'-gTC 3'	13(125)								
																		5'-gTC 3'	13(125)						
																				5'-CAA 3'	1 ⁸¹ I				
																				5'-gTA 3'	26(164)				
																				5'-gTA 3'	34(189)				
																				5'-gTC 3'	13(125)				
																				5'-gTC 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				
																				5'-gCT 3'	10(116)				

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

Well No.	DR	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87
DRB1 allele ⁴	ser ⁹																							
*01:01:01-01:02:08, 01:04-01:38, 01:40N-01:41, 01:43-01:45	DR1, Null, -	65																						
*01:03, 01:39N	DR103, Null		66																					
*01:42	-	65	66																					
*03:01:01-03:01:19, 03:04:01-03:06, 03:09, 03:11:01-03:16, 03:18-03:20, 03:22-03:23, 03:25-03:26, 03:28, 03:30-03:31, 03:33-03:34, 03:36-03:37, 03:43-03:45, 03:47-03:48, 03:50-03:52, 03:54-03:63, 03:66-03:68N, 03:70-03:73, 03:75, 03:77-03:80, 03:11:02	DR17(3), DR3, Null, -					69	70											81						
*03:02:01-03:03, 03:27, 03:29, 03:38, 03:53, 03:74	DR18(3), -					69	71											81						
*03:07, 03:17, 03:21, 03:24, 03:32, 03:35, 03:39-03:41, 03:49	DR3, -					69												81						
*03:08, 03:65	-					69	70				75		77											
*03:10	DR17(3)					69	70												83	84				
*03:42, 03:69	-					69																		
*03:46, 03:64	-					69	70																	
*03:76	-					70						76			79			81						
*04:01:01-04:61, 04:63-04:68, 04:70-04:72:02, 04:74-04:104, 04:106-04:107	DR4, Null, -							72																
*04:62, 04:69, 04:73, 04:105	-							72										82						
*07:01:01-07:01:04, 07:03-07:22	DR7, Null, -									73														
*08:01:01-08:02:04, 08:04:01-08:07, 08:11, 08:16-08:17, 08:22, 08:24, 08:26, 08:28, 08:39, 08:42-08:44	DR8, -										74						80							
*08:03:02-08:03:03, 08:10, 08:12-08:15, 08:18-08:19, 08:23, 08:25, 08:27, 08:29-08:30:03, 08:33-08:34, 08:36-08:38, 08:45-08:49	DR8, -										74													
*08:08	-										74						80	82						
*08:09, 14:15 ⁵	DR8										74						80			84				
*08:20, 13:18, 13:47, 13:55 ⁶	DR13(6), -															79	80	81		84				
*08:21	-										74						79	80		84				
Well No.	DR	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

Well No.	DR	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87
DRB1 allele ⁴	ser ⁹																							
*08:31, 08:41, 11:67 ⁷	DR11(5), -									74				77			80							
*08:32	-									74					78					84				
*08:35	-									74										84				
*08:40	-						70			74														
*09:01:02-09:01:05, 09:01:07-09:02:02, 09:04-09:17	DR9, -										75								83					
*09:01:06, 09:03	DR9, -										75													
*10:01:01-10:04	DR10, -											76												
*11:01:01-11:01:16, 11:04:01-11:04:08, 11:06:01-11:06:02, 11:09-11:10:02, 11:12:01-11:12:02, 11:15, 11:24, 11:27:01-11:30, 11:32-11:33, 11:35, 11:37:01-11:39, 11:43-11:44, 11:46:01-11:47, 11:49:01-11:51, 11:54:01-11:54:02, 11:56, 11:58:01-11:58:02, 11:60-11:62, 11:66, 11:72, 11:74:01-11:75, 11:77-11:78, 11:81, 11:84:01-11:84:02, 11:88, 11:90-11:92, 11:94-11:95, 11:97, 11:99-11:102:02, 11:106, 11:108-11:117, 11:120-11:121, 11:123, 11:126	DR11(5), -													77		79	80							
*11:02:01-11:03, 11:11:01-11:11:02, 11:14:01-11:14:02, 11:16, 11:20-11:21, 11:36, 11:40-11:41, 11:48, 11:63, 11:65:01-11:65:02, 11:68, 11:70, 11:76, 11:85-11:86, 11:93, 11:118, 11:122, 11:124	DR11(5), DR13 (6), -						70							77		79	80							
*11:05	DR11(5)													77			80							
*11:07, 11:53, 11:103, 11:105, 11:107, 11:125	DR11(5), -					69					75			77										
*11:08:01-11:08:02, 11:18-11:19:03, 11:42, 11:57	DR11(5), -													77		79								
*11:13:01-11:13:02	DR11(5)							w						77						83	84			
*11:17, 11:52	DR11(5), -													77						83	84			
*11:22, 11:98, 11:104	-													77										
*11:23, 11:25, 11:96	DR11(5), -													77		79	80			84				
*11:26, 11:34	DR11(5), -						71							77										
*11:31, 11:45, 11:64, 11:119	-													77		79				84				
*11:55	-													77			80			84				
Well No.	DR	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

Well No.	DR	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87
DRB1 allele ⁴	ser ⁹																							
*11:59, 11:80, 11:83, 11:87	-						70						76	77		79	80							
*11:69, 11:82	-													77			80		82					
*11:73, 11:79	-						70							77		79								
*11:89	-													77						84				
*12:01:01-12:01:04, 12:03, 12:05-12:12, 12:14, 12:17, 12:24N-12:25, 12:28-12:30, 12:34-12:35	DR12(5), Null, -															78								
*12:02:01-12:02:05, 12:13, 12:15, 12:18-12:21, 12:23, 12:26-12:27, 12:31N-12:33	DR12(5), Null, -															78		80						
*12:04	-											74				78								
*12:16	-											74				78		80						
*12:22	-											73	74			78								
*13:01:01-13:02:01, 13:02:03-13:02:05, 13:04, 13:08, 13:16, 13:20, 13:22-13:24, 13:28-13:29, 13:31-13:32, 13:34-13:36, 13:38-13:40, 13:48, 13:51-13:52, 13:54, 13:59, 13:61:01-13:61:02, 13:63-13:65, 13:68-13:70, 13:72-13:76, 13:78-13:80, 13:83-13:84, 13:87, 13:91-13:93, 13:96:01-13:99, 13:102, 13:105-13:107, 13:109, 13:111-13:114, 13:117, 13:121, 13:123-13:128, 13:130-13:131, 13:135, 13:138-13:139	DR13(6), DR14(6), Null, -						70																	
*13:02:02, 13:03:01-13:03:06, 13:10, 13:33:01-13:33:03, 13:37, 13:66:01-13:66:02, 13:81, 13:88-13:90, 13:94-13:95, 13:101, 13:115, 13:120, 13:133, 13:137N	DR13(6), Null-,						70									79		81						
*13:05:01-13:05:02, 13:07:01-13:07:02, 13:11:01-13:11:02, 13:14:01-13:14:03, 13:21:01-13:21:02, 13:42, 13:46, 13:49-13:50:02, 13:62, 13:100, 13:108, 13:132, 13:136	DR13(6), -																79	80	81					
Well No.	DR	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

Well No.	DR	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87
DRB1 allele ⁴	ser ⁹																							
*13:06, 13:12:01-13:12:02, 13:25, 13:30, 13:56, 13:58, 13:60, 13:77, 13:82, 13:110, 13:118, 13:134	DR13(6), -															79	81							
*13:09	-																80	81						
*13:13, 13:119, 14:84, 14:116 ⁸	DR13(6), -															79	81			84				
*13:15, 13:19, 13:53, 13:57, 13:104	DR13(6), -						70	71								79	80	81						
*13:17, 13:116	DR13(6), -						70			73							80							
*13:26	-							71								79	80	81						
*13:27, 13:41, 13:71	DR13(6), -						70					76				79	80	81						
*13:43	-						70									79	80			83	84			
*13:44, 13:86	-							71										81						
*13:45	-						70									79	80		82	84				
*13:67, 13:103	-															79	80							
*13:85	-						70	71								79		81						
*13:122	-						70									79								
*13:129	-						70					76					80	81						
*14:01:01-14:01:02, 14:04, 14:07:01-14:07:02, 14:10, 14:26, 14:28, 14:31, 14:35, 14:38-14:39, 14:54:01, 14:55, 14:57, 14:60-14:62, 14:70- 14:71, 14:75-14:76, 14:79, 14:86-14:88, 14:90, 14:99, 14:101, 14:104, 14:107, 14:110-14:114, 14:117- 14:118, 14:120, 14:122, 14:124-14:125	DR14(6), DR1404, -																		82	83	84			
*14:01:03, 14:08, 14:23:02, 14:34, 14:54:02, 14:72, 14:92N, 14:97	DR14(6), Null, -																			83	84			
*14:02, 14:06:01-14:06:02, 14:09, 14:13, 14:17, 14:20, 14:29-14:30, 14:33, 14:41, 14:47-14:48, 14:51, 14:80, 14:83, 14:94, 14:106, 14:108, 14:121	DR14(6), DR6-,							71										81		83				
*14:03:01-14:03:02, 14:12:01- 14:12:02, 14:40, 14:63, 14:67, 14:77-14:78, 14:85, 14:102, 14:115	DR14(6), DR1403, -							71								79		81			84			
Well No.	DR	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

Well No.	DR	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87
DRB1 allele ⁴	ser ⁹																							
*14:05:01-14:05:03, 14:14, 14:23:01, 14:23:03, 14:36, 14:43-14:45, 14:56, 14:59, 14:64, 14:91, 14:96, 14:100, 14:103, 14:123	DR14(6), -																	81		83	84			
*14:11	DR14(6)									74										83	84			
*14:16	DR6						70									79	80		82	83	84			
*14:18, 14:81	DR6, -							71										81		83	84			
*14:19, 14:21, 14:109	DR14(6), -						70	71								79		81		w				
*14:22, 14:105	DR14(6), -															79	80		82	83	84			
*14:24	-							71									80	81						
*14:25, 14:53	DR13(6), -															79	80		82		84			
*14:27	DR14(6)							71								79	80	81			84			
*14:32:01-14:32:02	-							w											82	83	84			
*14:37	-																80	81	82					
*14:42	-																	81			84			
*14:46, 14:52	-																			83				
*14:49, 14:119	DR14(6), -							71											82	83	84			
*14:50	DR14(6)									73									82	83	84			
*14:58	-																		82		84			
*14:65	-							w										81		83	84			
*14:68, 14:93	-									74									82	83	84			
*14:69	-															79			82		84			
*14:73	-																80		82	83	84			
*14:74	-															79			82	83	84			
*14:82	-						70												82	83	84			
*14:89	-							71										81			84			
*14:95	-						70											81		83	84			
*14:98	-							71								79		81						
*15:01:01-15:20, 15:22-15:24, 15:26, 15:28-15:33, 15:35-15:65, 15:67-15:71	DR15(2), DR2, Null, -			67																				
*15:21	-			67																	w			
*15:25	-			67		69																		
*15:27, 15:34, 15:66	-			67																83				
*16:01:01-16:03, 16:05:01-16:05:02, 16:07-16:17, 16:19	DR16(2), DR2, Null, -				68																			
*16:04, 16:18	DR16(2), -				68																w			
Well No.	DR	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

Well No.	DR	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87
DRB1 allele ⁴	ser ⁹																							
DRB3*01:01:02:01-01:15, DRB3*02:01-02:26, DRB3*02:28-02:29N, DRB3*03:01:01-03:03	DR52, Null, -																					85		
DRB3*02:27	-														79							85		
DRB4*01:01-01:03:01:01, DRB4*01:03:01-03:01:08	DR53, -																						86	
DRB4*01:03:01:02N	Null																	82					86	
DRB5*01:01-01:01:14, DRB5*02:02-02:05	DR51, Null, -																							87
Well No.	ser ⁵	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87

¹The internal positive control primer pairs amplify segments of the human growth hormone gene. The two different control primer pairs give rise to either an internal positive control band of 430 base pairs, for most wells, or a band of 515 base pairs, for some wells. Well number 65 contains the primer pair giving rise to the longer, 515 bp, internal positive control band. In addition, well number 74 contains the primer pair giving rise to the longer, 515 bp, internal positive control band.

²The codon, and in parenthesis the nucleotide, in the 2nd exon or the 1st intron, matching the specificity-determining 3'-end of the primer is given. Codon and nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

³The codon, and in parenthesis the nucleotide, in the 2nd exon, matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Codon and nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

⁴The DRB1*0702 allele has been shown to be identical to DRB1*07:01:01:01.

The DRB1*08031 allele has been shown to be identical to DRB1*08:03:02.

The DRB1*09011 allele has been shown to be identical to DRB1*09:01:02.

The DRB1*03:11:02 allele has been renamed DRB1*03:81.

The DRB1*1171 allele has been shown to be identical to DRB1*11:02:01.

The DRB1*12031 allele has been shown to be identical to DRB1*12:01:01.

The DRB1*1466 allele has been shown to be identical to DRB1*14:23:02.

The DRB1*1606 allele has been shown to be identical to DRB1*16:05:01.

The DRB3*010101 allele has been shown to be identical to DRB3*01:01:02:01.

The DRB4*0101102N allele has been renamed DRB4*01:03:01:02N.

The DRB5*0201 allele has been shown to be identical to DRB5*02:02.

The DRB4*02:01N and DRB4*03:01N null alleles will not be amplified by the DR low resolution primer set.

⁵The DRB1*08:09 and the DRB1*14:15 alleles yield identical amplification patterns with the DR low resolution primer set. These alleles can be separated by the respective high resolution primer sets.

⁶The DRB1*08:20 and the DRB1*13:18, 13:47 and 13:55 alleles yield identical amplification patterns with the DR low resolution primer set. These alleles can be separated by the respective high resolution primer sets.

⁷The DRB1*08:31, 08:41 and DRB1*11:67 alleles yield identical amplification patterns with the DR low resolution primer set. These alleles can be separated by the respective high resolution primer sets.

⁸The DRB1*13:13 and 13:119 and the DRB1*14:84 and 14:116 alleles yield identical amplification patterns with the DR low resolution primer set. These alleles can be separated by the respective high resolution primer sets.

⁹The grouping of not serologically defined alleles is taken from Tissue Antigens 73, 95-170, 2009.

'ser', serological HLA specificity.

'w', may be weakly amplified.

'?', nucleotide sequence information for the primer matching region(s) is not known.

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

INTERPRETATION TABLE										
DQ low resolution SSP typing										
Amplification patterns of the DQB1*05:01 to DQB1*04:08 alleles										
		Well ⁶								
		88	89	90	91	92	93	94	95	96
Length of spec.		225	220	210	220	130	135	145	160	Negative control
PCR product(s)			270					185		
Length of int.		515	430	430	515	515	515	515	430	
pos. control ¹										
5'-primer(s) ²		26(173) 5'-ggg 3'	9(122) 5'-gTT 3'	30(185) 5'-AAg 3'	26(173) 5'-TTA 3'	28(179) 5'-gAC 3'	26(173) 5'-TCT 3'	38(210) 5'-gCA 3'	38(210) 5'-gCg 3'	
			26(173) 5'-TTA 3'			28(179) 5'-gAC 3'		48(240) 5'-CgC 3'		
			26(173) 5'-TCT 3'					55(260) 5'-gCC 3'		
								55(260) 5'-gCA 3'		
3'-primer(s) ³		87(356) 5'-ggT 3'	86(353) 5'-ACg 3'	86(353) 5'-gCT 3'	86(353) 5'-gCT 3'	57(266) 5'-Cgg 3'	57(266) 5'-CgT 3'	86(353) 5'-gCT 3'	77(327) 5'-ACg 3'	
			86(353) 5'-ACC 3'							
Well No.		88	89	90	91	92	93	94	95	96
DQB1 allele ⁴	ser. ⁵									
*05:01:01:01-05:14	DQ5(1), -	88								Negative control
*06:01:01:06:28, 06:30-06:49	DQ6(1), DQ1, Null, -		89							
*06:29	-		89			92				
*02:01:01:02:02, 02:04-02:06	DQ2, -			90						
*02:03	DQ2			90			93			
*03:01:01:01-03:01:06, 03:04, 03:09-03:10, 03:13-03:14, 03:16, 03:19, 03:21-03:22, 03:24, 03:27-03:29, 03:35- 03:36	DQ7(3), DQ8(3), -				91			94		
*03:02:01-03:02:05, 03:05:01- 03:05:04, 03:07-03:08, 03:11, 03:18, 03:32, 03:37	DQ8(3), -					92		94		
*03:03:02:01-03:03:04, 03:06, 03:12, 03:15, 03:20, 03:25- 03:26, 03:30-03:31, 03:33- 03:34, 03:38-03:40	DQ9(3), DQ3, -						93	94		
*03:17, 03:23	-							94		
*04:01:01-04:02:02, 04:04- 04:08	DQ4, -								95	
*04:03:01-04:03:02	-						93		95	
DQB1 allele ⁴	ser. ⁵									
Well No.		88	89	90	91	92	93	94	95	96

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: 43R**Lot-specific information**

¹The internal positive control primer pairs amplify segments of the human growth hormone gene. The two different control primer pairs give rise to either an internal positive control band of 430 base pairs, for most wells, or a band of 515 base pairs, for some wells.

Well number 88 contains the primer pair giving rise to the longer, 515 bp, internal positive control band in order to help in the correct orientation of the DQ low resolution typing.

In addition, wells number 91, 92, 93 and 94 contain the primer pair giving rise to the longer, 515 bp, internal positive control band in order to allow kit identification.

²The codon, and in parenthesis the nucleotide, in the 2nd exon, matching the specificity-determining 3'-end of the primer is given. Codon and nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

³The codon, and in parenthesis the nucleotide, in the 2nd exon, matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Codon and nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

⁴The sequence of the DQB1*03031 allele has been shown to be identical to DQB1*03:03:02.

⁵The serological split of the DQB1*05:05-05:14, DQB1*06:06 to 06:07 alleles, the DQB1*06:10, 06:13, 06:15-06:24 and 06:27 to 06:49, the DQB1*02:04-02:06, the DQB1*03:07-03:09 and 03:11-03:40 alleles and the DQB1:04:0301-04:08 alleles is not known. In this table we have used the expert-assigned serological grouping in Tissue Antigens (2009) 73:95-170, and also inferred the serological grouping from the naming of the sequence-defined allele.

⁶Primer mix 96 contains a negative control, which will amplify more than 95% of HLA amplicons as well as the amplicons generated by control primer pairs. PCR product sizes range from 75 to 200 base pairs. The PCR product generated by the control primer pair is 430 base pairs.

Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot-specific information

[illegible]

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)Lot No.: **43R**

Lot-specific information

CELL LINE VALIDATION SHEET										
HLA-A low resolution primer set										
					Well					
					16	17	18	19	20	21
				Lot No.:	201201917	201201918	201201919	201201920	201201921	201201923
	IHWC cell line		A*	A*						
1	9001	SA	*24:02		-	-	-	-	-	-
2	9280	LK707	*02:01		-	-	-	-	-	-
3	9011	E4181324	*01:01		-	-	-	-	-	-
4	9275	GU373	*30:01		-	-	-	-	-	-
5	9009	KAS011	*01:01		-	-	-	-	-	-
6	9353	SM	*02:01	*26:03	-	-	-	-	-	+
7	9020	QBL	*26:01		-	-	-	-	-	-
8	9025	DEU	*31:01		-	-	-	-	-	-
9	9026	YAR	*26:01		-	-	-	-	-	-
10	9107	LKT3	*24:02		-	-	-	-	-	-
11	9051	PITOUT	*29:02		-	-	-	-	-	-
12	9052	DBB	*02:01		-	-	-	-	-	-
13	9004	JESTHOM	*02:01		-	-	-	-	-	-
14	9071	OLGA	*31:01		-	-	-	-	-	-
15	9075	DKB	*24:02		-	-	-	-	-	-
16	9037	SWEIG007	*29:02		-	-	-	-	-	-
17	9282	CTM3953540	*03:01	*80:01	-	-	-	-	-	+
18	9257	32367	*33:03	*74:01	-	+	+	-	-	-
19	9038	BM16	*02:01		-	-	-	-	-	-
20	9059	SLE005	*02:01		-	-	-	-	-	-
21	9064	AMALA	*02:17		-	-	-	-	-	-
22	9056	KOSE	*02:01		-	-	-	-	-	-
23	9124	IHL	*02:01	*34:01	-	-	-	-	-	-
24	9035	JBUSH	*32:01		+	-	-	-	-	-
25	9049	IBW9	*33:01		-	+	-	-	-	-
26	9285	WT49	*02:05		-	-	-	-	-	-
27	9191	CH1007	*24:10	*29:01	-	-	-	-	-	-
28	9320	BEL5GB	*02:01	*29:02	-	-	-	-	-	-
29	9050	MOU	*29:02		-	-	-	-	-	-
30	9021	RSH	*30:01	*68:02	-	-	-	+	-	-
31	9019	DUCAF	*30:02		-	-	-	-	-	-
32	9297	HAG	*02:01		-	-	-	-	-	-
33	9098	MT14B	*31:01		-	-	-	-	-	-
34	9104	DHIF	*31:01		-	-	-	-	-	-
35	9302	SSTO	*32:01		+	-	-	-	-	-
36	9024	KT17	*02:06	*11:01	-	-	-	-	-	-
37	9065	HHKB	*03:01		-	-	-	-	-	-
38	9099	LZL	*02:17		-	-	-	-	-	-
39	9315	CML	*01:01	*03:01	-	-	-	-	-	-
40	9134	WHONP199	*02:07	*30:01	-	-	-	-	-	-
41	9055	H0301	*03:01		-	-	-	-	-	-
42	9066	TAB089	*02:07		-	-	-	-	-	-
43	9076	T7526	*02:06	*02:07	-	-	-	-	-	-
44	9057	TEM	*66:01		-	-	-	-	-	-
45	9239	SHJO	*23:01	*24:02	-	-	-	-	-	-
46	9013	SCHU	*03:01		-	-	-	-	-	-
47	9045	TUBO	*02:16	*03:01	-	-	-	-	-	-
48	9303	TER-ND	*02:01	*11:01	-	-	-	-	-	-

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

CELL LINE VALIDATION SHEET																			
HLA-B low resolution SSP kit																			
					Well														
					22	23	24	25	26	27	28	29	30	31	32	33	34	35	36
				Prod. No.:	201204248	201204202	201204203	201204204	201204205	201204206	201204207	201204208	201207209	201204210	201204211	201204213	201204214	201204215	201204216
	IHC cell line		B*																
1	9001	SA	*07:02		+	-	-	-	-	-	-	-	-	-	-	-	-	-	-
2	9280	LK707	*52:01	*73:01	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-
3	9011	E4181324	*52:01		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-
4	9275	GU373	*15:10	*53:01	-	-	-	-	-	-	-	+	-	+	-	-	+	-	-
5	9009	KAS011	*37:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	-
6	9353	SM	*39:01	*51:01	-	-	-	-	+	-	-	-	+	-	-	-	-	-	-
7	9020	QBL	*18:01		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-
8	9025	DEU	*35:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
9	9026	YAR	*38:01		-	-	-	-	+	-	-	-	-	-	-	-	-	-	-
10	9107	LKT3	*54:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
11	9051	PITOUT	*44:03		-	-	-	+	-	-	-	-	-	-	-	-	-	-	-
12	9052	DBB	*57:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
13	9004	JESTHOM	*27:05		-	-	-	-	+	-	-	-	-	-	-	-	-	-	+
14	9071	OLGA	*15:01	*15:20	-	-	-	-	-	-	+	-	-	-	-	-	-	-	-
15	9075	DKB	*40:01		-	-	-	+	-	-	-	-	-	-	-	-	-	-	-
16	9037	SWEIG007	*40:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
17	9282	CTM3953540	*08:01	*55:01	+	+	-	-	-	-	-	+	-	-	-	-	-	-	-
18	9257	32367	*14:01	*56:01	-	-	-	+	+	-	-	-	-	-	-	-	-	-	-
19	9038	BM16	*18:01		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-
20	9059	SLE005	*40:01		-	-	-	+	-	-	-	-	-	-	-	-	-	-	-
21	9064	AMALA	*15:01		-	-	-	-	-	-	+	-	-	-	-	-	-	-	-
22	9056	KOSE	*35:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
23	9124	IHL	*40:02	*56:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
24	9035	JBUSH	*38:01		-	-	-	-	+	-	-	-	-	-	-	-	-	-	-
25	9049	IBW9	*14:02		-	-	-	+	-	+	-	-	-	-	-	-	-	-	-
26	9285	WT49	*58:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
27	9191	CH1007	*07:05	*51:01	+	-	-	-	-	-	-	-	+	-	-	-	-	-	-
28	9320	BEL5GB	*44:02	*44:03	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-
29	9050	MOU	*44:03		-	-	-	+	-	-	-	-	-	-	-	-	-	-	-
30	9021	RSH	*42:01		-	-	-	-	-	-	-	+	-	-	-	-	-	-	-
31	9019	DUCAF	*18:01		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-
32	9297	HAG	*41:02		-	-	-	+	-	-	-	+	-	-	-	-	-	-	-
33	9098	MT14B	*40:01		-	-	-	+	-	-	-	-	-	-	-	-	-	-	-
34	9104	DHIF	*38:01		-	-	-	-	+	-	-	-	-	-	-	-	-	-	-
35	9302	SSTO	*44:02		-	-	-	+	-	-	-	-	-	-	-	-	-	-	-
36	9024	KT17	*15:01	*35:01	-	-	-	-	-	-	+	-	-	-	-	-	-	-	-
37	9065	HHKB	*07:02		+	-	-	-	-	-	-	-	-	-	-	-	-	-	-
38	9099	LZL	*15:01		-	-	-	-	-	-	+	-	-	-	-	-	-	-	-
39	9315	CML	*08:01	*27:05	-	+	-	-	+	-	-	+	-	-	-	-	-	-	+
40	9134	WHONP199	*13:02	*46:01	-	-	+	+	-	-	+	-	-	+	-	-	-	-	-
41	9055	H0301	*14:02		-	-	-	+	-	+	-	-	-	-	-	-	-	-	-
42	9066	TAB089	*46:01		-	-	-	-	-	-	+	-	-	-	-	-	-	-	-
43	9076	T7526	*46:01		-	-	-	-	-	-	+	-	-	-	-	-	-	-	-
44	9057	TEM	*38:01		-	-	-	-	+	-	-	-	-	-	-	-	-	-	-
45	9239	SHJO	*42:01	*50:01	-	-	-	+	-	-	-	+	-	-	-	-	-	-	-
46	9013	SCHU	*07:02		+	-	-	-	-	-	-	-	-	-	-	-	-	-	-
47	9045	TUBO	*51:01		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-
48	9303	TER-ND	*35:01	*44:03	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

CELL LINE VALIDATION SHEET																			
HLA-B low resolution SSP kit																			
					Well														
						37	38	39	40	41	42	43	44	45	46	47	48	49	50
				Prod. No.:		201204217	201204218	201208719	201204220	201204221	201204222	201204223	201208724	201204225	201204227	201204228	201204229	201204231	201204232
	IHWC cell line		B*																
1	9001	SA	*07:02			-	-	-	+	-	-	-	+	-	-	+	-	-	-
2	9280	LK707	*52:01	*73:01		-	-	-	-	-	-	-	+	-	-	+	-	-	-
3	9011	E4181324	*52:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-
4	9275	GU373	*15:10	*53:01		+	-	+	-	-	-	-	-	-	-	-	-	-	-
5	9009	KAS011	*37:01			+	-	-	+	-	-	-	-	-	-	-	+	-	-
6	9353	SM	*39:01	*51:01		-	-	-	-	+	-	+	+	-	-	-	-	-	-
7	9020	QBL	*18:01			-	+	-	-	-	-	-	-	-	-	-	-	-	-
8	9025	DEU	*35:01			+	+	-	-	-	-	-	-	-	-	-	-	-	-
9	9026	YAR	*38:01			-	-	-	-	-	+	+	-	-	-	-	-	-	-
10	9107	LKT3	*54:01			-	-	-	-	-	-	-	-	-	-	+	-	-	+
11	9051	PITOUT	*44:03			+	-	-	-	-	-	-	-	-	-	-	+	-	-
12	9052	DBB	*57:01			+	-	-	-	-	-	-	-	-	-	-	-	-	-
13	9004	JESTHOM	*27:05			-	-	-	-	-	-	-	-	-	-	+	-	-	+
14	9071	OLGA	*15:01	*15:20		+	-	+	-	-	-	-	-	-	-	-	-	-	-
15	9075	DKB	*40:01			-	-	-	-	-	-	-	+	+	-	-	-	-	-
16	9037	SWEIG007	*40:02			-	-	-	-	-	-	-	-	-	-	-	-	-	-
17	9282	CTM3953540	*08:01	*55:01		-	-	-	-	-	-	-	+	-	+	+	-	-	+
18	9257	32367	*14:01	*56:01		-	-	-	-	-	-	+	+	-	-	+	-	-	+
19	9038	BM16	*18:01			-	+	-	-	-	-	-	-	-	-	-	-	-	-
20	9059	SLE005	*40:01			-	-	-	-	-	-	-	+	+	-	-	-	-	-
21	9064	AMALA	*15:01			-	-	+	-	-	-	-	-	-	-	-	-	-	-
22	9056	KOSE	*35:03			+	+	-	-	-	-	-	-	-	-	-	-	-	-
23	9124	IHL	*40:02	*56:02		-	-	-	-	-	-	-	-	-	-	+	-	-	+
24	9035	JBUSH	*38:01			-	-	-	-	+	+	+	-	-	-	-	-	-	-
25	9049	IBW9	*14:02			-	-	-	-	-	-	+	+	-	-	-	-	-	-
26	9285	WT49	*58:01			+	-	-	-	-	-	-	-	-	-	-	-	-	-
27	9191	CH1007	*07:05	*51:01		-	-	-	-	-	-	-	+	-	-	+	-	-	-
28	9320	BEL5GB	*44:02	*44:03		+	-	-	-	-	-	-	-	-	-	-	+	-	-
29	9050	MOU	*44:03			+	-	-	-	-	-	-	-	-	-	-	+	-	-
30	9021	RSH	*42:01			-	-	-	-	-	-	-	+	-	+	+	-	-	-
31	9019	DUCAF	*18:01			-	+	-	-	-	-	-	-	-	-	-	-	-	-
32	9297	HAG	*41:02			-	-	-	-	-	-	-	+	-	+	-	-	-	-
33	9098	MT14B	*40:01			-	-	-	-	-	-	-	+	+	-	-	-	-	-
34	9104	DHIF	*38:01			-	-	-	-	+	+	+	-	-	-	-	-	-	-
35	9302	SSTO	*44:02			+	-	-	-	-	-	-	-	-	-	-	+	-	-
36	9024	KT17	*15:01	*35:01		+	+	+	-	-	-	-	-	-	-	-	-	-	-
37	9065	HHKB	*07:02			-	-	-	-	-	-	-	+	-	-	+	-	-	-
38	9099	LZL	*15:01			-	-	+	-	-	-	-	-	-	-	-	-	-	-
39	9315	CML	*08:01	*27:05		-	-	-	-	-	-	-	+	-	+	+	-	-	+
40	9134	WHONP199	*13:02	*46:01		-	-	-	-	-	-	-	-	-	-	+	-	+	-
41	9055	H0301	*14:02			-	-	-	-	-	-	+	+	-	-	-	-	-	-
42	9066	TAB089	*46:01			-	-	-	-	-	-	-	-	-	-	+	-	+	-
43	9076	T7526	*46:01			-	-	-	-	-	-	-	-	-	-	+	-	+	-
44	9057	TEM	*38:01			-	-	-	-	+	+	+	-	-	-	-	-	-	-
45	9239	SHJO	*42:01	*50:01		-	-	-	-	-	-	-	+	-	+	+	-	-	-
46	9013	SCHU	*07:02			-	-	-	-	-	-	-	+	-	-	+	-	-	-
47	9045	TUBO	*51:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-
48	9303	TER-ND	*35:01	*44:03		+	+	-	-	-	-	-	-	-	-	-	+	-	-

101.708-24 – including *Taq* polymerase, IFU-01
 101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

CELL LINE VALIDATION SHEET																			
HLA-B low resolution SSP kit																			
					Well														
						51	52	53	54	55	56	57	58	59	60	61	62	63	64
					Prod. No.:	201204233	201204234	201204235	201204236	201204237	201204238	201204239	201204240	201204241	201204242	201204244	201204245	201204246	201204247
	IHWC cell line		B*																
1	9001	SA	*07:02			-	-	-	-	-	-	-	-	-	-	-	-	-	+
2	9280	LK707	*52:01	*73:01		-	-	+	+	-	+	-	-	-	-	-	-	+	-
3	9011	E4181324	*52:01			-	-	+	+	-	+	-	-	-	-	-	+	+	+
4	9275	GU373	*15:10	*53:01		-	-	-	+	-	-	-	-	-	-	+	-	+	+
5	9009	KAS011	*37:01			-	-	-	-	-	+	-	-	-	-	-	-	+	-
6	9353	SM	*39:01	*51:01		-	-	+	+	-	-	-	-	-	-	-	+	+	+
7	9020	QBL	*18:01			-	-	-	-	-	-	-	-	-	-	+	-	-	+
8	9025	DEU	*35:01			-	-	-	-	-	-	-	-	-	-	+	-	-	+
9	9026	YAR	*38:01			-	-	-	-	-	-	-	-	-	-	-	-	+	-
10	9107	LKT3	*54:01			-	-	-	-	-	-	+	+	-	-	-	-	-	+
11	9051	PITOUT	*44:03			-	-	-	-	-	+	-	-	-	-	-	-	+	-
12	9052	DBB	*57:01			-	-	-	+	-	-	-	-	-	+	-	-	+	-
13	9004	JESTHOM	*27:05			-	-	-	-	-	+	-	-	-	-	-	-	+	-
14	9071	OLGA	*15:01	*15:20		-	-	-	-	-	+	-	-	-	-	-	+	-	+
15	9075	DKB	*40:01			-	-	-	-	+	+	-	-	-	-	-	-	-	+
16	9037	SWEIG007	*40:02			-	-	-	-	+	+	-	-	-	-	-	-	-	+
17	9282	CTM3953540	*08:01	*55:01		-	-	+	-	-	-	-	+	-	-	-	-	-	+
18	9257	32367	*14:01	*56:01		-	-	-	-	-	-	-	-	+	-	-	-	-	+
19	9038	BM16	*18:01			-	-	-	-	-	-	-	-	-	-	+	-	-	+
20	9059	SLE005	*40:01			-	-	-	-	+	+	-	-	-	-	-	-	-	+
21	9064	AMALA	*15:01			-	-	-	-	-	+	-	-	-	-	-	+	-	+
22	9056	KOSE	*35:03			-	-	-	-	-	-	-	-	-	-	+	-	-	+
23	9124	IHL	*40:02	*56:02		-	-	-	-	+	+	-	-	+	-	-	-	-	+
24	9035	JBUSH	*38:01			-	-	-	-	-	-	-	-	-	-	-	-	+	-
25	9049	IBW9	*14:02			-	-	-	-	-	-	-	-	-	-	-	-	-	+
26	9285	WT49	*58:01			-	-	-	+	-	-	-	-	-	-	-	-	+	-
27	9191	CH1007	*07:05	*51:01		-	-	+	+	-	-	-	-	-	-	+	-	+	+
28	9320	BEL5GB	*44:02	*44:03		-	-	-	-	-	+	+	-	-	-	-	-	+	-
29	9050	MOU	*44:03			-	-	-	-	-	+	-	-	-	-	-	-	+	-
30	9021	RSH	*42:01			-	-	-	-	-	-	-	-	-	-	-	-	-	+
31	9019	DUCAF	*18:01			-	-	-	-	-	-	-	-	-	-	+	-	-	+
32	9297	HAG	*41:02			-	-	-	-	+	+	-	-	-	-	-	-	-	+
33	9098	MT14B	*40:01			-	-	-	-	+	+	-	-	-	-	-	-	-	+
34	9104	DHIF	*38:01			-	-	-	-	-	-	-	-	-	-	-	-	+	-
35	9302	SSTO	*44:02			-	-	-	-	-	+	+	-	-	-	-	-	+	-
36	9024	KT17	*15:01	*35:01		-	-	-	-	-	+	-	-	-	-	+	+	-	+
37	9065	HHKB	*07:02			-	-	-	-	-	-	-	-	-	-	-	-	-	+
38	9099	LZL	*15:01			-	-	-	-	-	+	-	-	-	-	-	+	-	+
39	9315	CML	*08:01	*27:05		-	-	-	-	-	+	-	-	-	-	-	-	+	+
40	9134	WHONP199	*13:02	*46:01		-	-	-	-	-	+	-	-	-	-	-	+	+	+
41	9055	H0301	*14:02			-	-	-	-	-	-	-	-	-	-	-	-	-	+
42	9066	TAB089	*46:01			-	-	-	-	-	-	-	-	-	-	-	+	-	+
43	9076	T7526	*46:01			-	-	-	-	-	-	-	-	-	-	-	+	-	+
44	9057	TEM	*38:01			-	-	-	-	-	-	-	-	-	-	-	-	+	-
45	9239	SHJO	*42:01	*50:01		-	-	+	-	+	+	-	-	+	-	-	-	-	+
46	9013	SCHU	*07:02			-	-	-	-	-	-	-	-	-	-	-	-	-	+
47	9045	TUBO	*51:01			-	-	+	+	-	-	-	-	-	-	+	-	+	-
48	9303	TER-ND	*35:01	*44:03		-	-	-	-	-	+	-	-	-	-	+	-	+	+

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

CELL LINE VALIDATION SHEET																				
DR low resolution primer set																				
					Well															
					65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80
				Prod. No.:	201297201	201297202	201297203	201297204	201297205	201297206	201297207	201297208	201297209	201297210	201297211	201297212	201297213	201297214	201205515	201297216
	IHC cell line		DRB1																	
1	9001	SA	*01:01		+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
2	9280	LK707	*15:02	*04:05	-	-	+	-	-	-	-	+	-	-	-	-	-	-	-	-
3	9011	E4181324	*15:02		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
4	9275	GU373	*03:01		-	-	-	-	+	+	-	-	-	-	-	-	-	-	-	-
5	9009	KAS011	*16:01		-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-
6	9353	SM	*04:07	*08:03	-	-	-	-	-	-	-	+	-	+	-	-	-	-	-	-
7	9020	QBL	*03:01		-	-	-	-	+	+	-	-	-	-	-	-	-	-	-	-
8	9025	DEU	*04:01		-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	-
9	9026	YAR	*04:02		-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	-
10	9107	LKT3	*04:05		-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	-
11	9051	PITOUT	*07:01		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-
12	9052	DBB	*07:01		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-
13	9004	JESTHOM	*01:01		+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
14	9071	OLGA	*08:02		-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	+
15	9075	DKB	*09:01		-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-
16	9037	SWEIG007	*11:01		-	-	-	-	-	-	-	-	-	-	-	-	+	-	+	+
17	9282	CTM3953540	*03:01	*13:01	-	-	-	-	+	+	-	-	-	-	-	-	-	-	+	+
18	9257	32367	*09:01	*11:01	-	-	-	-	-	-	-	-	-	-	+	-	+	-	+	+
19	9038	BM16	*12:01		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	-
20	9059	SLE005	*13:02		-	-	-	-	-	+	-	-	-	-	-	-	-	-	+	+
21	9064	AMALA	*14:02		-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-
22	9056	KOSE	*13:02	*14:01	-	-	-	-	-	+	-	-	-	-	-	-	-	-	+	+
23	9124	IHL	*08:03	*14:14	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-
24	9035	JBUSH	*11:01		-	-	-	-	-	-	-	-	-	-	-	-	+	-	+	+
25	9049	IBW9	*07:01		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-
26	9285	WT49	*03:01		-	-	-	-	+	+	-	-	-	-	-	-	-	-	-	-
27	9191	CH1007	*04:05	*10:01	-	-	-	-	-	-	-	+	-	-	-	+	-	-	-	-
28	9320	BEL5GB	*04:16	*07:01	-	-	-	-	-	-	-	+	+	-	-	-	-	-	-	-
29	9050	MOU	*07:01		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-
30	9021	RSH	*03:02		-	-	-	-	+	-	+	-	-	-	-	-	-	-	-	-
31	9019	DUCAF	*03:01		-	-	-	-	+	+	-	-	-	-	-	-	-	-	-	-
32	9297	HAG	*13:03		-	-	-	-	-	+	-	-	-	-	-	-	-	-	+	-
33	9098	MT14B	*04:04		-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	-
34	9104	DHIF	*11:01		-	-	-	-	-	-	-	-	-	-	-	-	+	-	+	+
35	9302	SSTO	*04:03		-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	-
36	9024	KT17	*04:03	*04:06	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	-
37	9065	HHKB	*13:01		-	-	-	-	-	+	-	-	-	-	-	-	-	-	+	+
38	9099	LZL	*14:02		-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-
39	9315	CML	*03:01	*04:01	-	-	-	-	+	+	-	+	-	-	-	-	-	-	-	-
40	9134	WHONP199	*07:01	*09:01	-	-	-	-	-	-	-	-	+	-	+	-	-	-	-	-
41	9055	H0301	*13:02		-	-	-	-	-	+	-	-	-	-	-	-	-	-	+	+
42	9066	TAB089	*08:03		-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-
43	9076	T7526	*09:01		-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-
44	9057	TEM	*14:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
45	9239	SHJO	*07:01		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-
46	9013	SCHU	*15:01		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
47	9045	TUBO	*11:04	*12:01	-	-	-	-	-	-	-	-	-	-	-	-	+	+	+	+
48	9303	TER-ND	*01:03		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-

Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot-specific information

CELL LINE VALIDATION SHEET												
DR low resolution primer set												
					Prod. No.:	Well ¹						
						81	82	83	84	85	86	87
						201297217	201208818	201208819	201297220	201208829	201208830	201208831
	IHCW cell line		DRB1									
1	9001	SA	*01:01			-	-	-	-	-	-	-
2	9280	LK707	*15:02	*04:05		-	-	-	-	-	+	+
3	9011	E4181324	*15:02			-	-	-	-	-	-	+
4	9275	GU373	*03:01			+	-	-	-	+	-	-
5	9009	KAS011	*16:01			-	-	-	-	-	-	+
6	9353	SM	*04:07	*08:03		-	-	-	-	-	+	-
7	9020	QBL	*03:01			+	-	-	-	+	-	-
8	9025	DEU	*04:01			-	-	-	-	-	+	-
9	9026	YAR	*04:02			-	-	-	-	-	+	-
10	9107	LKT3	*04:05			-	-	-	-	-	+	-
11	9051	PITOUT	*07:01			-	-	-	-	-	+	-
12	9052	DBB	*07:01			-	+	-	-	-	+	-
13	9004	JESTHOM	*01:01			-	-	-	-	-	-	-
14	9071	OLGA	*08:02			-	-	-	-	-	-	-
15	9075	DKB	*09:01			-	-	+	-	-	+	-
16	9037	SWEIG007	*11:01			-	-	-	-	+	-	-
17	9282	CTM3953540	*03:01	*13:01		+	-	-	-	+	-	-
18	9257	32367	*09:01	*11:01		-	-	+	-	+	+	-
19	9038	BM16	*12:01			-	-	-	-	+	-	-
20	9059	SLE005	*13:02			+	-	-	-	+	-	-
21	9064	AMALA	*14:02			+	-	+	-	+	-	-
22	9056	KOSE	*13:02	*14:01		+	+	+	+	+	-	-
23	9124	IHL	*08:03	*14:14		+	-	+	+	+	-	-
24	9035	JBUSH	*11:01			-	-	-	-	+	-	-
25	9049	IBW9	*07:01			-	-	-	-	-	+	-
26	9285	WT49	*03:01			+	-	-	-	+	-	-
27	9191	CH1007	*04:05	*10:01		-	-	-	-	-	+	-
28	9320	BEL5GB	*04:16	*07:01		-	-	-	-	-	+	-
29	9050	MOU	*07:01			-	-	-	-	-	+	-
30	9021	RSH	*03:02			+	-	-	-	+	-	-
31	9019	DUCAF	*03:01			+	-	-	-	+	-	-
32	9297	HAG	*13:03			+	-	-	-	+	-	-
33	9098	MT14B	*04:04			-	-	-	-	-	+	-
34	9104	DHIF	*11:01			-	-	-	-	+	-	-
35	9302	SSTO	*04:03			-	-	-	-	-	+	-
36	9024	KT17	*04:03	*04:06		-	-	-	-	-	+	-
37	9065	HHKB	*13:01			+	-	-	-	+	-	-
38	9099	LZL	*14:02			+	-	+	-	+	-	-
39	9315	CML	*03:01	*04:01		+	-	-	-	+	+	-
40	9134	WHONP199	*07:01	*09:01		-	-	+	-	-	+	-
41	9055	H0301	*13:02			+	-	-	-	+	-	-
42	9066	TAB089	*08:03			-	-	-	-	-	-	-
43	9076	T7526	*09:01			-	-	+	-	-	+	-
44	9057	TEM	*14:01			-	+	+	+	+	-	-
45	9239	SHJO	*07:01			-	-	-	-	+	+	-
46	9013	SCHU	*15:01			-	-	-	-	-	-	+</

¹The DRB4*01:03:01:02N allele is amplified by primer mix 82 in the DBB/9052 cell line.

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

CELL LINE VALIDATION SHEET														
DQ low resolution primer set														
					Well									
					88	89	90	91	92	93	94	95		
				Production No.	201194801	201194802	201194803	201194804	201194805	201194806	201194807	201208908		
	IHWC cell line		DQB1											
1	9001	SA	*05:01		+	-	-	-	-	-	-	-	-	-
2	9280	LK707	*06:01	*02:02	-	+	+	-	-	-	-	-	-	-
3	9011	E4181324	*06:01		-	+	-	-	-	-	-	-	-	-
4	9275	GU373	*02:01		-	-	+	-	-	-	-	-	-	-
5	9009	KAS011	*05:02		+	-	-	-	-	-	-	-	-	-
6	9353	SM	*03:02	*06:01	-	+	-	-	+	-	+	-	-	-
7	9020	QBL	*02:01		-	-	+	-	-	-	-	-	-	-
8	9025	DEU	*03:01		-	-	-	+	-	-	+	-	-	-
9	9026	YAR	*03:02		-	-	-	-	+	-	+	-	-	-
10	9107	LKT3	*04:01		-	-	-	-	-	-	-	-	+	-
11	9051	PITOUT	*02:02		-	-	+	-	-	-	-	-	-	-
12	9052	DBB	*03:03		-	-	-	-	-	+	+	-	-	-
13	9004	JESTHOM	*05:01		+	-	-	-	-	-	-	-	-	-
14	9071	OLGA	*04:02		-	-	-	-	-	-	-	-	+	-
15	9075	DKB	*03:03		-	-	-	-	-	+	+	-	-	-
16	9037	SWEIG007	*03:01		-	-	-	+	-	-	+	-	-	-
17	9282	CTM 3953540	*02:01	*06:03	-	+	+	-	-	-	-	-	-	-
18	9257	32367	*06:02	*02:02	-	+	+	-	-	-	-	-	-	-
19	9038	BM16	*03:01		-	-	-	+	-	-	+	-	-	-
20	9059	SLE005	*06:04		-	+	-	-	-	-	-	-	-	-
21	9064	AMALA	*03:01		-	-	-	+	-	-	+	-	-	-
22	9056	KOSE	*05:03	*06:04	+	+	-	-	-	-	-	-	-	-
23	9124	IHL	*05:03	*06:01	+	+	-	-	-	-	-	-	-	-
24	9035	JBUSH	*03:01		-	-	-	+	-	-	+	-	-	-
25	9049	IBW9	*02:02		-	-	+	-	-	-	-	-	-	-
26	9285	WT49	*02:01		-	-	+	-	-	-	-	-	-	-
27	9191	CH1007	*04:01	*05:01	+	-	-	-	-	-	-	-	+	-
28	9320	BEL5GB	*02:02	*03:01	-	-	+	+	-	-	+	-	-	-
29	9050	MOU	*02:02		-	-	+	-	-	-	-	-	-	-
30	9021	RSH	*04:02		-	-	-	-	-	-	-	-	+	-
31	9019	DUCAF	*02:01		-	-	+	-	-	-	-	-	-	-
32	9297	HAG	*03:01		-	-	-	+	-	-	+	-	-	-
33	9098	MT14B	*03:02		-	-	-	-	+	-	+	-	-	-
34	9104	DHIF	*03:01		-	-	-	+	-	-	+	-	-	-
35	9302	SSTO	*03:05		-	-	-	-	+	-	+	-	-	-
36	9024	KT17	*03:02		-	-	-	-	+	-	+	-	-	-
37	9065	HHKB	*06:03		-	+	-	-	-	-	-	-	-	-
38	9099	LZL	*03:01		-	-	-	+	-	-	+	-	-	-
39	9315	CML	*02:01	*03:01	-	-	+	+	-	-	+	-	-	-
40	9134	WHONP199	*02:02	*03:03	-	-	+	-	-	+	+	-	-	-
41	9055	H0301	*06:09		-	+	-	-	-	-	-	-	-	-
42	9066	TAB089	*06:01		-	+	-	-	-	-	-	-	-	-
43	9076	T7526	*03:03		-	-	-	-	-	+	+	-	-	-
44	9057	TEM	*05:03		+	-	-	-	-	-	-	-	-	-
45	9239	SHJO	*02:02		-	-	+	-	-	-	-	-	-	-
46	9013	SCHU	*06:02		-	+	-	-	-	-	-	-	-	-
47	9045	TUBO	*03:01		-	-	-	+	+	-	+	-	-	-
48	9303	TER-ND	*05:01		+	-	-	-	-	-	-	-	-	-

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

CERTIFICATE OF ANALYSIS**Olerup SSP® HLA-A-B-DR-DQ SSP Combi Tray**

Product number: 101.708-24 – including *Taq* pol.
 101.708-24u – without *Taq* pol.
 Lot number: 43R
 Expiry date: 2015-March-01
 Number of tests: 24 tests
 Number of wells per test: 95 +1

Well specifications:

Well No.	Production No.	Well No.	Production No.	Well No.	Production No.
1	2012-071-01	9	2012-041-10	17	2012-019-18
2	2012-019-02	10	2012-019-11	18	2012-019-19
3	2012-071-03	11	2012-019-12	19	2012-019-20
4	2012-019-04	12	2012-019-13	20	2012-019-21
5	2012-019-06	13	2012-041-14	21	2012-019-23
6	2012-019-07	14	2012-019-15		
7	2012-019-08	15	2012-019-16		
8	2012-019-09	16	2012-019-17		

The specificity of each primer solution of the kit has been tested against 48 well characterized IHWC cell line DNAs.

Additional 5'- and/or 3'-primers in primer solutions 1, 3 to 7, 9, 11 and 13 to 21 were tested by separately adding additional 3'-primers or 5'-primers. One or two primers in primer solutions 2, 3, 7, 8, 9, 10, 14, 17 and 18 were not possible to test.

Well No.	Production No.	Well No.	Production No.	Well No.	Production No.
22	2012-042-48	38	2012-042-18	54	2012-042-36
23	2012-042-02	39	2012-087-19	55	2012-042-37
24	2012-042-03	40	2012-042-20	56	2012-042-38
25	2012-042-04	41	2012-042-21	57	2012-042-39
26	2012-042-05	42	2012-042-22	58	2012-042-40
27	2012-042-06	43	2012-042-23	59	2012-042-41
28	2012-042-07	44	2012-087-24	60	2012-042-42
29	2012-042-08	45	2012-042-25	61	2012-042-44
30	2012-072-09	46	2012-042-27	62	2012-042-45
31	2012-042-10	47	2012-042-28	63	2012-042-46
32	2012-042-11	48	2012-042-29	64	2012-042-47
33	2012-042-13	49	2012-042-31		
34	2012-042-14	50	2012-042-32		
35	2012-042-15	51	2012-042-33		
36	2012-042-16	52	2012-042-34		
37	2012-042-17	53	2012-042-35		

The specificity of each primer solution of the kit has been tested against 48 well characterized IHWC cell line DNAs.

Additional 5'-primers and/or 3'-primers in primer solutions 22 to 24, 27, 33 to 35, 39, 41, 45, 48, 49, 53 and 58 were tested by separately adding additional 3'-primers or 5'-primers. One or two primers in primer solutions 49, 56 and 60 were not possible to test.

101.708-24 – including *Taq* polymerase, IFU-01101.708-24u – without *Taq* polymerase, IFU-02Visit www.olerup-ssp.com for

“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

Well No.	Production No.	Well No.	Production No.	Well No.	Production No.
65	2012-972-01	73	2012-972-09	81	2012-972-17
66	2012-972-02	74	2012-972-10	82	2012-088-18
67	2012-972-03	75	2012-972-11	83	2012-088-19
68	2012-972-04	76	2012-972-12	84	2012-972-20
69	2012-972-05	77	2012-972-13	85	2012-088-29
70	2012-972-06	78	2012-972-14	86	2012-088-30
71	2012-972-07	79	2012-055-15	87	2012-088-31
72	2012-972-08	80	2012-972-16		

The specificity of each primer solution of the kit has been tested against 48 well characterized IHWC cell line DNAs.

Additional 5'-primers and/or 3'-primers in primer solutions 65, 67, 68, 70, 73 to 76, 79, 80, 82 and 84 were tested by separately adding additional 3'-primers, respectively additional 5'-primers.

One, two or three of the primers in primer solutions 65, 67, 68, 72 to 74, 77, 79, 80 and 86 were not possible to test.

Well No.	Production No.
88	2011-948-01
89	2011-948-02
90	2011-948-03
91	2011-948-04
92	2011-948-05
93	2011-948-06
94	2011-948-07
95	2011-089-08

The specificity of each primer solution of the kit has been tested against 48 well characterized IHWC cell line DNAs.

One additional 5'-primer in primer solution 89 was tested by separately adding one additional 3'-primer.

The negative control primer pairs, **Production No. 2012-002-01**, can detect contamination with PCR products diluted 10^{-7} .

Results: No false positive or false negative amplifications were obtained.

Date of approval: 2012-November-29

Approved by

Production Quality Control

101.708-24 – including *Taq* polymerase, IFU-01
101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

Declaration of Conformity

Product name: *Olerup* SSP® HLA-A-B-DR-DQ SSP Combi Tray
Product number: 101.708-24/24u
Lot number: 43R

Intended use: HLA-A, HLA-B, HLA-DR and HLA-DQ low resolution histo-
compatibility testing

Manufacturer: *Olerup* SSP AB
Franzengatan 5
SE-112 51 Stockholm, Sweden
Phone: +46-8-717 88 27
Fax: +46-8-717 88 18

We, *Olerup* SSP AB, hereby declare that this product, to which this Declaration of Conformity relates is in conformity with the following Standard(s) and other normative document(s) ISO 9001:2008 and ISO 13485:2012, following the provisions of the 98/79/EC Directive on *in vitro* diagnostic medical devices, Annex II List B, conformity assessed using Annex IV, as transposed into the national laws of the Member States of the European Union.

The Technical Documentation File is maintained at *Olerup* SSP AB, Franzengatan 5, SE-112 51 Stockholm, Sweden.

Notified Body: Lloyd's Register Quality Assurance Limited, Hiramford, Middlemarch Office Village, Siskin Drive, Coventry CV3 4FJ, United Kingdom. (Notified Body number: 0088.)

Stockholm, Sweden
2012-November-29

Ann-Cathrin Jareman
Head of QA and Regulatory Affairs

101.708-24 – including *Taq* polymerase, IFU-01
101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

101.708-24 – including *Taq* polymerase, IFU-01
101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

101.708-24 – including *Taq* polymerase, IFU-01
101.708-24u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: **43R**

Lot-specific information

ADDRESSES:

Manufacturer:

Olerup SSP AB, Franzengatan 5, SE-112 51 Stockholm, Sweden.

Tel: +46-8-717 88 27

Fax: +46-8-717 88 18

E-mail: info-ssp@olerup.com

Web page: <http://www.olerup-ssp.com>

Distributed by:

Olerup GmbH, Löwengasse 47 / 6, AT-1030 Vienna, Austria.

Tel: +43-1-710 15 00

Fax: +43-1-710 15 00 10

E-mail: support-at@olerup.com

Web page: <http://www.olerup.com>

Olerup Inc., 901 S. Bolmar St., Suite R, West Chester, PA 19382

Tel: 1-877-OLERUP1

Fax: 610-344-7989

E-mail: info.us@olerup.com

Web page: <http://www.olerup.com>

For information on *Olerup* SSP distributors worldwide, contact **Olerup GmbH**.