

101.311-24/06 – including *Taq* pol., IFU-01
101.311-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: **09N**

Lot-specific information
Olerup SSP[®] DPB1

Product number: 101.311-24/04 – including *Taq* pol.
101.311-24u/04u – without *Taq* pol.
Lot number: 09N
Expiry date: 2014-August-01
Number of tests: 24 tests – Product No. 101.311-24/24u
4 tests – Product No. 101.311-04/04u
Number of wells per test: 48
Storage - pre-aliquoted primers: dark at -20°C
- PCR Master Mix: -20°C
- Adhesive PCR seals RT
- Product Insert RT

This Product Description is only valid for Lot No. 09N.

CHANGES COMPARED TO THE PREVIOUS OLERUP SSP[®] DPB1 LOT

The Lot-specific information for DPB1 including and without *Taq* polymerase is now described in one common Product Insert.

The DPB1 specificity and interpretation tables have been updated for the DPB1 alleles described since the previous Olerup SSP[®] DPB1 lot was made (**Lot No. 09M**).

The primers of the wells detailed below have been added, exchanged or modified.

Well	5'-primer	3'-primer	rationale
3	Modified	-	Improved specificity of primer pair.
41	-	Added	Primer added for the DPB1*134:01 allele.

Change in revision R01 compared to R00:

1. Primer mix 42 may have a tendency of giving rise to nonspecific amplifications.

Change in revision R02 compared to R01:

1. The DPB1*31:01 allele is not amplified by primer mix 4. This has been corrected in the Specificity and Interpretation Tables.

101.311-24/06 – including *Taq* pol., IFU-01
101.311-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: **09N**

Lot-specific information

101.311-24/06 – including *Taq* pol., IFU-01
101.311-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: **09N**

Lot-specific information

PRODUCT DESCRIPTION

DPB1 SSP typing

CONTENT

The primer set contains 5'- and 3'-primers for identifying the DPB1*01:01 to DPB1*136:01 alleles.

PLATE LAYOUT

Each test consists of 48 PCR reactions in a 48 well cut PCR plate.

1	2	3	4	5	6	7	8
9	10	11	12	13	14	15	16
17	18	19	20	21	22	23	24
25	26	27	28	29	30	31	32
33	34	35	36	37	38	39	40
41	42	43	44	45	46	47	48

The 48 well cut PCR plate is marked with 'DPB1' in silver/gray ink.

Well No. 1 is marked with the Lot No. '09N'.

A faint row of numbers is seen between wells 1 and 2 or wells 7 and 8 of the PCR trays. These stem from the manufacture of the trays, and should be disregarded.

The PCR plates are heat-sealed with a PCR-compatible foil.

Please note: When removing each 48 well PCR plate, make sure that the remaining plates stay sealed. Use a scalpel or a similar instrument to carefully cut the foil between the plates.

INTERPRETATION

Only DPB1 alleles will be amplified by the DPB1 typing kit. Thus, the interpretation of DPB1 typings is not influenced by the DPB2 gene or other HLA class II genes.

UNIQUELY IDENTIFIED ALLELES

All the phenotypically different DPB1 alleles, i.e. **DPB1*01:01 to DPB1*136:01**, recognized by the HLA Nomenclature Committee in January 2012¹ will give rise to unique amplification patterns by the primers in the DPB1 typing kit.

The DPB1 typing kit cannot distinguish the following silent mutations: DPB1*01:01:01-01:01:03, the 02:01:02, 02:01:04-02:01:05 and 02:01:07, 03:01:01-03:01:02, 04:01:01:01-04:01:03, 04:02:01:01-04:02:01:02, 05:01:01-05:01:02, 11:01:01-11:01:02, 20:01:01-20:01:02, and the DPB1*41:01:01-41:01:02 alleles.

¹DPB1 alleles listed on the IMGT/HLA web page 2012-January-12, release 3.7.0, www.ebi.ac.uk/imgt/hla.

101.311-24/06 – including *Taq* pol., IFU-01
101.311-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: 09N

Lot-specific information

RESOLUTION IN HOMO- AND HETEROZYGOTES

A total of 151 alleles generate 135 amplification patterns that can be combined in 9180 homozygous and heterozygous combinations. 716 of these genotypes do not give rise to unique amplification patterns. The different lengths of the specific PCR products were not considered in these calculations.

++++++	++-+---	+--+---+	---+---+	-----	--+-----	*19:01, *26:01:01 = *19:01, *26:01:02
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*26:01:01, *106:01 = *26:01:02, *106:01
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*77:01, *127:01 = *94:01, *127:01
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*08:01, *26:01:01 = *08:01, *26:01:02
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*51:01, *118:01 = *120:01N, *122:01
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*61:01N, *65:01 = *65:01, *92:01
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*65:01, *127:01 = *68:01, *127:01 = *84:01, *127:01
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*40:01, *127:01 = *53:01, *127:01
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*26:01:02, *65:01 = *26:01:02, *84:01
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*26:01:01, *65:01 = *26:01:01, *84:01
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*01:01:01, *65:01 = *01:01:01, *84:01
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*65:01, *118:01 = *84:01, *118:01
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*21:01, *127:01 = *36:01, *127:01
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*61:01N, *127:01 = *92:01, *127:01
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*01:01:01, *111:01 = *03:01:01, *118:01
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*01:01:01, *61:01N = *01:01:01, *92:01
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*61:01N, *118:01 = *92:01, *118:01
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*01:01:01, *35:01:02 = *01:01:01, *132:01
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*89:01, *127:01 = *90:01, *127:01 = *122:01, *127:01 = *127:01, *127:01
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*89:01, *122:01 = *90:01, *122:01
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*01:01:01, *13:01 = *13:01, *90:01
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*01:01:01, *26:01:01 = *01:01:01, *26:01:02
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*01:01:01, *01:01:01 = *01:01:01, *89:01 = *01:01:01, *90:01
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*89:01, *118:01 = *90:01, *118:01 = *118:01, *118:01
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*61:01N, *90:01 = *90:01, *92:01
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*89:01, *90:01 = *90:01, *90:01
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*26:01:01, *123:01 = *26:01:02, *123:01
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*02:01:02, *26:01:01 = *02:01:02, *26:01:02
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*16:01, *26:01:01 = *16:01, *26:01:02
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*02:02, *26:01:01 = *02:02, *26:01:02
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*22:01, *26:01:01 = *22:01, *26:01:02
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*04:02:01:01, *26:01:01 = *04:02:01:01, *26:01:02
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*26:01:01, *105:01 = *26:01:02, *105:01
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*23:01, *26:01:01 = *23:01, *26:01:02
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*05:01:01, *26:01:01 = *05:01:01, *26:01:02
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*26:01:01, *33:01 = *26:01:02, *33:01
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*24:01, *26:01:01 = *24:01, *26:01:02
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*26:01:01, *121:01 = *26:01:02, *121:01
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*04:01:01:01, *26:01:01 = *04:01:01:01, *26:01:02
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*26:01:01, *134:01 = *26:01:02, *134:01
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*26:01:01, *39:01 = *26:01:02, *39:01
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*06:01, *26:01:01 = *06:01, *26:01:02
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*17:01, *26:01:01 = *17:01, *26:01:02
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*21:01, *26:01:01 = *21:01, *26:01:02
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*26:01:01, *30:01 = *26:01:02, *30:01
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*20:01:01, *26:01:01 = *20:01:01, *26:01:02
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*18:01, *26:01:01 = *18:01, *26:01:02
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*26:01:01, *55:01 = *26:01:02, *55:01
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*26:01:01, *44:01 = *26:01:02, *44:01
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*09:01, *26:01:01 = *09:01, *26:01:02
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*10:01, *26:01:01 = *10:01, *26:01:02
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*26:01:01, *114:01 = *26:01:02, *114:01
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*14:01, *26:01:01 = *14:01, *26:01:02
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*03:01:01, *26:01:01 = *03:01:01, *26:01:02
++++++	++-+---	+--+---+	---+---+	-----	--+-----	*26:01:01, *124:01 = *26:01:02, *124:01

101.311-24/06 – including *Taq* pol., IFU-01
 101.311-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: 09N

Lot-specific information

+++++--	-----+	+++++--	---+---+	-----+	---+---+	*26:01:02, *61:01N = *26:01:02, *92:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*26:01:01, *61:01N = *26:01:01, *92:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*26:01:02, *35:01:02 = *26:01:02, *132:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*26:01:01, *35:01:02 = *26:01:01, *132:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*25:01, *26:01:01 = *25:01, *26:01:02
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*26:01:01, *28:01 = *26:01:02, *28:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*11:01:01, *26:01:01 = *11:01:01, *26:01:02
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*13:01, *26:01:01 = *13:01, *26:01:02 = *13:01, *27:01 =
						*13:01, *89:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*26:01:01, *133:01 = *26:01:02, *133:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*27:01, *107:01 = *89:01, *107:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*27:01, *133:01 = *89:01, *133:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*15:01, *26:01:01 = *15:01, *26:01:02
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*26:01:02, *85:01 = *26:01:02, *125:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*26:01:02, *26:01:02 = *26:01:02, *27:01 = *26:01:02, *89:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*26:01:01, *85:01 = *26:01:01, *125:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*26:01:01, *26:01:01 = *26:01:01, *26:01:02 = *26:01:01,
						*27:01 = *26:01:01, *89:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*35:01:02, *56:01 = *85:01, *92:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*27:01, *92:01 = *56:01, *132:01 = *61:01N, *125:01 = *92:01,
						*125:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*56:01, *61:01N = *56:01, *92:01 = *89:01, *92:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*25:01, *125:01 = *56:01, *79:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*27:01, *56:01 = *56:01, *85:01 = *56:01, *125:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*56:01, *56:01 = *56:01, *89:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*27:01, *110:01 = *89:01, *110:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*03:01:01, *04:01:01:01 = *04:01:01:01, *124:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*27:01, *37:01 = *37:01, *89:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*27:01, *61:01N = *61:01N, *85:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*27:01, *132:01 = *89:01, *132:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*27:01, *79:01 = *79:01, *89:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*03:01:01, *133:01 = *13:01, *124:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*77:01, *122:01 = *94:01, *122:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*68:01, *122:01 = *84:01, *122:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*61:01N, *122:01 = *92:01, *122:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*68:01, *92:01 = *84:01, *92:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*77:01, *92:01 = *92:01, *94:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*03:01:01, *76:01 = *14:01, *92:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*61:01N, *76:01 = *76:01, *92:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*61:01N, *92:01 = *92:01, *92:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*19:01, *124:01 = *104:01, *106:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*03:01:01, *19:01 = *19:01, *104:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*03:01:01, *106:01 = *106:01, *124:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*50:01, *77:01 = *50:01, *94:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*68:01, *78:01 = *78:01, *84:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*03:01:01, *68:01 = *03:01:01, *84:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*68:01, *104:01 = *84:01, *104:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*68:01, *111:01 = *84:01, *111:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*68:01, *124:01 = *84:01, *124:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*61:01N, *68:01 = *61:01N, *84:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*68:01, *132:01 = *84:01, *132:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*25:01, *68:01 = *25:01, *84:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*68:01, *79:01 = *79:01, *84:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*50:01, *68:01 = *50:01, *84:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*35:01:02, *50:01 = *50:01, *132:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*02:01:02, *03:01:01 = *02:01:02, *124:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*02:02, *03:01:01 = *02:02, *124:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*77:01, *78:01 = *78:01, *94:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*03:01:01, *77:01 = *03:01:01, *94:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*03:01:01, *105:01 = *105:01, *124:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*77:01, *104:01 = *94:01, *104:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*77:01, *111:01 = *94:01, *111:01
+++++--	-----+	+++++--	---+---+	-----+	---+---+	*77:01, *124:01 = *94:01, *124:01

101.311-24/06 – including *Taq* pol., IFU-01
101.311-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: 09N

Lot-specific information

++-----+	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*61:01N,*77:01=*61:01N,*94:01
++-----+	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*77:01,*132:01=*94:01,*132:01
++-----+	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*25:01,*77:01=*25:01,*94:01
++-----+	+--+-----+	+--+-----	-----+-----	-----+-----	+--+-----+	*77:01,*79:01=*79:01,*94:01
++-----+	+--+-----+	+++++-----	+-----+-----	-----+-----	+--+-----+	*03:01:01,*05:01:01=*05:01:01,*104:01
++-----+	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*25:01,*112:01=*52:01,*116:01
++-----+	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*17:01,*124:01=*124:01,*131:01
++-----+	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*06:01,*132:01=*20:01:01,*88:01
++-----+	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*37:01,*87:01=*52:01,*93:01
++-----+	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*21:01,*114:01=*36:01,*114:01
++-----+	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*25:01,*119:01=*79:01,*119:01
++-----+	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*29:01,*114:01=*44:01,*114:01
++-----+	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*25:01,*88:01=*29:01,*79:01
++-----+	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*70:01,*114:01=*114:01,*124:01
++-----+	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*25:01,*35:01:02=*25:01,*132:01
++-----+	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*03:01:01,*35:01:02=*03:01:01,*132:01
++-----+	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*35:01:02,*124:01=*124:01,*132:01
++-----+	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*03:01:01,*104:01=*104:01,*104:01
++-----+	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*35:01:02,*132:01=*132:01,*132:01
+++++-----	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*04:01:01:01,*13:01=*04:01:01:01,*133:01
+++++-----	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*13:01,*106:01=*106:01,*133:01
+++++-----	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*02:01:02,*13:01=*02:01:02,*133:01
+++++-----	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*02:02,*13:01=*02:02,*133:01
+++++-----	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*13:01,*105:01=*105:01,*133:01
+++++-----	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*17:01,*133:01=*131:01,*133:01
+++++-----	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*30:01,*107:01=*30:01,*133:01
+++++-----	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*13:01,*107:01=*107:01,*107:01
+++++-----	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*44:01,*55:01=*44:01,*58:01
+++++-----	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*29:01,*44:01=*44:01,*44:01
+++++-----	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*27:01,*76:01=*76:01,*85:01
+++++-----	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*27:01,*125:01=*89:01,*125:01
+++++-----	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*17:01,*27:01=*17:01,*89:01
+++++-----	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*27:01,*131:01=*89:01,*131:01
+++++-----	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*27:01,*93:01=*89:01,*93:01
+++++-----	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*27:01,*36:01=*36:01,*89:01
+++++-----	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*21:01,*27:01=*21:01,*89:01
+++++-----	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*27:01,*98:01=*89:01,*98:01
+++++-----	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*27:01,*58:01=*58:01,*89:01
+++++-----	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*27:01,*119:01=*89:01,*119:01
+++++-----	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*27:01,*66:01=*66:01,*85:01
+++++-----	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*09:01,*27:01=*09:01,*89:01
+++++-----	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*10:01,*27:01=*10:01,*89:01
+++++-----	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*27:01,*54:01=*54:01,*89:01
+++++-----	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*27:01,*35:01:01=*35:01:01,*89:01
+++++-----	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*85:01,*85:01=*85:01,*89:01
+++++-----	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*27:01,*27:01=*27:01,*89:01
+++++-----	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*11:01:01,*20:01:01=*11:01:01,*130:01
+++++-----	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*68:01,*130:01=*84:01,*130:01
+++++-----	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*20:01:01,*68:01=*20:01:01,*84:01
+++++-----	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*18:01,*68:01=*18:01,*84:01
+++++-----	+--+-----+	+++++-----	-----+-----	-----+-----	+--+-----+	*77:01,*130:01=*94:01,*130:01
++++						

101.311-24/06 – including *Taq* pol., IFU-01
101.311-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: 09N

Lot-specific information

---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*38:01, *128:01 = *97:01, *120:01N
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*77:01, *99:01 = *94:01, *99:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*60:01, *99:01 = *83:01, *99:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*04:01:01:01, *04:02:01:01 = *105:01, *126:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*60:01, *126:01 = *83:01, *126:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*04:01:01:01, *77:01 = *105:01, *134:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*04:01:01:01, *82:01 = *105:01, *128:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*04:01:01:01, *60:01 = *04:01:01:01, *83:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*39:01, *60:01 = *39:01, *83:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*60:01, *121:01 = *83:01, *121:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*60:01, *134:01 = *83:01, *134:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*60:01, *128:01 = *83:01, *128:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*60:01, *120:01N = *83:01, *120:01N
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*60:01, *96:01 = *83:01, *96:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*31:01, *77:01 = *77:01, *102:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*40:01, *60:01 = *40:01, *83:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*34:01, *121:01 = *34:01, *134:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*31:01, *94:01 = *94:01, *102:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*17:01, *39:01 = *39:01, *131:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*40:01, *99:01 = *53:01, *99:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*51:01, *99:01 = *99:01, *99:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*40:01, *51:01 = *53:01, *126:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*31:01, *134:01 = *102:01, *134:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*31:01, *34:01 = *34:01, *102:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*31:01, *31:01 = *31:01, *102:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*17:01, *33:01 = *33:01, *131:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*33:01, *60:01 = *33:01, *83:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*23:01, *33:01 = *33:01, *71:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*68:01, *76:01 = *76:01, *84:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*51:01, *68:01 = *51:01, *84:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*49:01, *68:01 = *49:01, *84:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*28:01, *68:01 = *28:01, *84:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*53:01, *68:01 = *53:01, *84:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*76:01, *77:01 = *76:01, *94:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*46:01, *49:01 = *46:01, *51:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*49:01, *115:01 = *51:01, *115:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*49:01, *86:01 = *51:01, *86:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*24:01, *60:01 = *24:01, *83:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*49:01, *108:01 = *51:01, *108:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*49:01, *80:01 = *51:01, *80:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*49:01, *77:01 = *49:01, *94:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*51:01, *77:01 = *51:01, *94:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*49:01, *60:01 = *49:01, *83:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*51:01, *60:01 = *51:01, *83:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*28:01, *77:01 = *28:01, *94:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*53:01, *77:01 = *53:01, *94:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*53:01, *60:01 = *53:01, *83:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*17:01, *28:01 = *28:01, *131:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*08:01, *82:01 = *32:01, *68:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*08:01, *60:01 = *08:01, *83:01 = *41:01:01, *68:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*08:01, *46:01 = *46:01, *84:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*08:01, *123:01 = *84:01, *123:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*02:01:02, *08:01 = *02:01:02, *84:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*08:01, *32:01 = *32:01, *84:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*08:01, *115:01 = *84:01, *115:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*08:01, *41:01:01 = *41:01:01, *84:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*19:01, *60:01 = *19:01, *83:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*60:01, *106:01 = *83:01, *106:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*68:01, *71:01 = *75:01, *103:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*08:01, *48:01 = *48:01, *84:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*02:02, *84:01 = *02:02, *106:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*19:01, *19:01 = *19:01, *84:01
---+---+ +---+---+ ---+---+ +---+---+ ---+---+ +---+---+	*57:01, *77:01 = *57:01, *94:01

101.311-24/06 – including *Taq* pol., IFU-01
101.311-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: 09N

Lot-specific information

-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*68:01, *108:01 = *84:01, *108:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*68:01, *80:01 = *80:01, *84:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*68:01, *129:01 = *84:01, *129:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*59:01, *68:01 = *59:01, *84:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*04:02:01:01, *68:01 = *04:02:01:01, *84:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*68:01, *77:01 = *68:01, *94:01 = *77:01, *84:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*68:01, *82:01 = *82:01, *84:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*60:01, *68:01 = *60:01, *84:01 = *68:01, *83:01 = *83:01, *84:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*68:01, *105:01 = *84:01, *105:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*73:01, *80:01 = *75:01, *108:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*73:01, *77:01 = *73:01, *94:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*75:01, *77:01 = *75:01, *94:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*60:01, *75:01 = *75:01, *83:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*68:01, *116:01 = *84:01, *116:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*68:01, *91:01 = *84:01, *91:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*68:01, *98:01 = *84:01, *98:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*08:01, *81:01 = *81:01, *84:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*08:01, *08:01 = *08:01, *84:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*68:01, *119:01 = *84:01, *119:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*68:01, *73:01 = *73:01, *84:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*68:01, *75:01 = *75:01, *84:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*68:01, *70:01 = *70:01, *84:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*14:01, *68:01 = *14:01, *84:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*35:01:01, *68:01 = *35:01:01, *84:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*35:01:02, *68:01 = *35:01:02, *84:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*45:01, *68:01 = *45:01, *84:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*57:01, *68:01 = *57:01, *84:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*68:01, *68:01 = *68:01, *84:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*48:01, *60:01 = *48:01, *83:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*60:01, *81:01 = *81:01, *83:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*16:01, *60:01 = *16:01, *83:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*46:01, *60:01 = *46:01, *83:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*60:01, *123:01 = *83:01, *123:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*02:01:02, *60:01 = *02:01:02, *83:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*32:01, *60:01 = *32:01, *83:01 = *41:01:01, *82:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*60:01, *115:01 = *83:01, *115:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*06:01, *46:01 = *06:01, *115:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*46:01, *64:01N = *64:01N, *115:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*17:01, *123:01 = *123:01, *131:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*02:01:02, *17:01 = *02:01:02, *131:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*17:01, *46:01 = *17:01, *115:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*46:01, *131:01 = *115:01, *131:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*46:01, *86:01 = *86:01, *115:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*46:01, *46:01 = *46:01, *115:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*02:01:02, *23:01 = *02:01:02, *71:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*60:01, *101:01 = *83:01, *101:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*02:02, *60:01 = *02:02, *83:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*02:02, *23:01 = *02:02, *71:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*22:01, *60:01 = *22:01, *83:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*47:01, *60:01 = *47:01, *83:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*60:01, *103:01 = *83:01, *103:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*60:01, *71:01 = *71:01, *83:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*02:02, *17:01 = *02:02, *131:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*17:01, *23:01 = *23:01, *131:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*60:01, *100:01 = *83:01, *100:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*05:01:01, *60:01 = *05:01:01, *83:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*60:01, *97:01 = *83:01, *97:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*38:01, *60:01 = *38:01, *83:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*60:01, *112:01 = *83:01, *112:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*77:01, *116:01 = *94:01, *116:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*60:01, *116:01 = *83:01, *116:01
-----+--+ +-----+ -----+ -+-----+ -----++- -+-----+	*70:01, *77:01 = *70:01, *94:01

101.311-24/06 – including *Taq* pol., IFU-01
101.311-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: 09N

Lot-specific information

-----+	+-----+	+++++---	-----+	-----+	-----+	*77:01, *119:01 = *94:01, *119:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*14:01, *77:01 = *14:01, *94:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*77:01, *91:01 = *91:01, *94:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*35:01:01, *77:01 = *35:01:01, *94:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*35:01:02, *77:01 = *35:01:02, *94:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*77:01, *98:01 = *94:01, *98:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*45:01, *77:01 = *45:01, *94:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*60:01, *63:01 = *63:01, *83:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*60:01, *62:01 = *62:01, *83:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*77:01, *108:01 = *94:01, *108:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*77:01, *80:01 = *80:01, *94:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*60:01, *80:01 = *80:01, *83:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*77:01, *129:01 = *94:01, *129:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*59:01, *77:01 = *59:01, *94:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*04:02:01:01, *77:01 = *04:02:01:01, *94:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*04:02:01:01, *60:01 = *04:02:01:01, *83:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*23:01, *60:01 = *23:01, *83:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*77:01, *82:01 = *82:01, *94:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*60:01, *77:01 = *60:01, *94:01 = *77:01, *83:01 = *83:01, *94:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*77:01, *105:01 = *94:01, *105:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*77:01, *77:01 = *77:01, *94:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*60:01, *82:01 = *82:01, *83:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*60:01, *105:01 = *83:01, *105:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*60:01, *60:01 = *60:01, *83:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*54:01, *109:01 = *54:01, *113:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*30:01, *109:01 = *30:01, *113:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*09:01, *109:01 = *09:01, *113:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*17:01, *109:01 = *17:01, *113:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*109:01, *131:01 = *113:01, *131:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*10:01, *109:01 = *10:01, *113:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*58:01, *109:01 = *58:01, *113:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*55:01, *109:01 = *55:01, *113:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*86:01, *86:01 = *86:01, *109:01 = *86:01, *113:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*109:01, *113:01 = *113:01, *113:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*91:01, *119:01 = *98:01, *119:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*17:01, *30:01 = *30:01, *131:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*09:01, *91:01 = *14:01, *131:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*17:01, *35:01:01 = *17:01, *35:01:02
-----+	+-----+	+++++---	-----+	-----+	-----+	*17:01, *55:01 = *55:01, *131:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*17:01, *17:01 = *17:01, *131:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*35:01:01, *119:01 = *119:01, *119:01
-----+	+-----+	+++++---	-----+	-----+	-----+	*09:01, *35:01:01 = *09:01, *35:01:02
-----+	+-----+	+++++---	-----+	-----+	-----+	*14:01, *35:01:01 = *14:01, *35:01:02
-----+	+-----+	+++++---	-----+	-----+	-----+	*35:01:01, *35:01:01 = *35:01:01, *35:01:02

```
*01:01:01 = *01:01:01-01:01:03
*02:01:02 = *02:01:02-02:01:07
*03:01:01 = *03:01:01-03:01:02
*04:01:01 = *04:01:01-04:01:02
*04:02:01:01 = 04:02:01:01-04:02:01:02
*05:01:01 = *05:01:01-05:01:02
*11:01:01 = *11:01:01-11:01:02
*20:01:01 = *20:01:01-20:01:02
*41:01:01 = *41:01:01-41:01:02
```

101.311-24/06 – including *Taq* pol., IFU-01
 101.311-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **09N**

Lot-specific information

SPECIFICITY TABLE

DPB1 SSP typing

Specificities and sizes of the PCR products of the 48 primer mixes used for DPB1 SSP typing³

Primer Mix	Size of spec. PCR product ¹	Size of control band ²	Amplified DPB1 alleles
1	245 bp	515 bp	*01:01:01-01:01:03, 03:01:01-03:01:02, 13:01, 25:01, 26:01:01-26:01:02, 29:01, 37:01, 44:01, 50:01, 52:01, 56:01, 61:01N, 78:01, 79:01, 88:01, 90:01, 92:01, 104:01, 107:01, 110:01, 111:01, 114:01, 118:01, 122:01, 124:01, 127:01, 132:01, 133:01
2	220 bp	515 bp	*01:01:01-01:01:03, 03:01:01-03:01:02, 18:01, 20:01:01-20:01:02, 25:01, 26:01:01-26:01:02, 27:01, 36:01, 50:01, 52:01, 56:01, 61:01N, 78:01, 79:01, 85:01, 87:01, 89:01, 90:01, 92:01, 104:01, 111:01, 114:01, 118:01, 122:01, 124:01, 127:01, 130:01, 132:01
3^{4,5,6}	100 bp	430 bp	*01:01:01-01:01:03, 04:01:01-04:01:03, 11:01:01-11:01:02, 13:01, 15:01, 26:01:01-26:01:02, 27:01, 31:01, 33:01, 39:01, 40:01, 56:01, 65:01, 66:01, 72:01, 74:01, 85:01, 89:01, 90:01, 96:01, 99:01, 102:01, 107:01, 110:01, 117:01, 118:01, 120:01N, 121:01, 125:01, 126:01, 127:01, 128:01, 133:01, 134:01
4	140 bp	430 bp	*01:01:01-01:01:03, 04:01:01-04:01:03, 24:01, 26:01:01-26:01:02, 27:01, 28:01, 39:01, 40:01, 49:01, 51:01, 53:01, 56:01, 65:01, 66:01, 72:01, 76:01, 85:01, 89:01, 90:01, 92:01, 96:01, 99:01, 102:01, 118:01, 120:01N, 121:01, 122:01, 125:01, 126:01, 127:01, 128:01, 134:01
5	235 bp	515 bp	*01:01:01-01:01:03, 08:01, 19:01, 50:01, 57:01, 65:01, 68:01, 73:01, 75:01, 84:01, 90:01, 106:01, 118:01, 122:01, 127:01
6	165 bp	430 bp	*01:01:01-01:01:03, 13:01, 26:01:01-26:01:02, 65:01, 107:01, 118:01, 122:01, 127:01, 133:01
7	160 bp	430 bp	*01:01:01-01:01:03, 13:01, 26:01:01-26:01:02, 56:01, 65:01, 76:01, 90:01, 92:01, 107:01, 110:01, 118:01, 122:01, 125:01, 127:01, 133:01
8	210 bp	515 bp	*02:01:02-02:02, 04:01:01-04:02:01:02, 05:01:01-05:01:02, 08:01, 16:01, 19:01, 22:01, 23:01, 24:01, 32:01, 33:01, 38:01, 39:01, 40:01, 46:01, 47:01, 48:01, 49:01, 51:01, 53:01, 60:01, 62:01, 63:01, 65:01, 68:01, 71:01, 75:01, 77:01, 80:01, 81:01, 82:01, 84:01, 94:01, 96:01, 97:01, 99:01, 100:01, 101:01, 103:01, 105:01,

101.311-24/06 – including *Taq* pol., IFU-01
 101.311-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **09N**

Lot-specific information

			106:01, 112:01, 115:01, 116:01, 120:01N, 121:01, 123:01, 126:01, 128:01, 134:01
9	165 bp, 160 bp	430 bp	*02:01:02-02:01:05, 02:01:06 ^w , 02:01:07-02:02, 04:02:01:01-04:02:01:02, 05:01:01-05:01:02, 06:01, 16:01, 17:01, 18:01, 19:01, 20:01:01-20:01:02, 21:01, 22:01, 23:01, 30:01, 32:01, 34:01, 36:01, 38:01, 41:01:01-41:01:02, 46:01, 47:01, 48:01, 55:01, 58:01, 59:01, 60:01, 62:01, 63:01, 64:01N, 69:01, 71:01, 77:01, 80:01, 82:01, 83:01, 86:01, 87:01, 91:01, 93:01, 94:01, 95:01, 97:01, 98:01, 100:01, 101:01, 103:01, 105:01, 106:01, 108:01, 109:01, 112:01, 113:01, 115:01, 116:01, 123:01, 129:01, 130:01, 131:01, 136:01
10	220 bp	515 bp	*02:01:02-02:02, 08:01, 16:01, 19:01, 22:01, 32:01, 33:01, 41:01:01-41:01:02, 46:01, 47:01, 48:01, 71:01, 81:01, 95:01, 101:01, 103:01, 106:01, 115:01, 123:01
11	190 bp	515 bp	*02:01:02-02:01:07, 04:01:01:01-04:02:01:02, 23:01, 24:01, 32:01, 33:01, 41:01:01-41:01:02, 46:01, 47:01, 51:01, 59:01, 60:01, 66:01, 71:01, 72:01, 73:01, 75:01, 77:01, 80:01, 81:01, 82:01, 83:01, 86:01, 94:01, 99:01, 105:01, 108:01, 109:01, 113:01, 115:01, 119:01, 120:01N, 121:01, 123:01, 126:01, 128:01, 129:01, 134:01
12⁴	80 bp	515 bp	*02:01:02, 02:01:03 ^w , 02:01:04-02:01:07, 06:01, 08:01, 09:01, 10:01, 16:01, 17:01, 29:01, 32:01, 37:01, 41:01:01-41:01:02, 44:01, 46:01, 48:01, 64:01N, 81:01, 86:01, 88:01, 93:01, 109:01, 113:01, 115:01, 123:01, 131:01
13	180 bp	430 bp	*02:02, 05:01:01-05:01:02, 19:01, 22:01, 24:01, 38:01, 47:01, 97:01, 100:01, 101:01, 106:01
14^{4,8}	105 bp, 175 bp	430 bp	*02:02, 05:01:01-05:01:02, 21:01, 22:01, 30:01, 36:01, 38:01, 54:01, 97:01, 100:01, 101:01, 114:01
15	200 bp	430 bp	*02:02, 48:01, 95:01, 99:01, 100:01, 101:01, 112:01, 116:01
16^{4,6,9}	80 bp, 165 bp, 210 bp	430 bp	*03:01:01-03:01:02, 04:02:01:01-04:02:01:02, 14:01, 18:01, 20:01:01-20:01:02, 25:01, 28:01, 35:01:01-35:01:02, 45:01, 49:01, 50:01, 51:01, 53:01, 57:01, 59:01, 60:01, 61:01N, 68:01, 70:01, 73:01, 75:01, 76:01, 77:01, 78:01, 79:01, 80:01, 82:01, 83:01, 91:01, 92:01, 98:01, 99:01, 104:01, 105:01, 108:01, 111:01, 116:01, 119:01, 122:01, 124:01, 127:01, 129:01, 130:01, 132:01
17	260 bp	515 bp	*03:01:01-03:01:02, 06:01, 09:01, 10:01, 11:01:01-11:01:02, 13:01, 14:01, 17:01, 20:01:01-20:01:02,

101.311-24/06 – including *Taq* pol., IFU-01
 101.311-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **09N**

Lot-specific information

			21:01, 25:01, 26:01:01-26:01:02, 27:01, 29:01, 30:01, 35:01:01-35:01:02, 36:01, 37:01, 44:01, 45:01, 52:01, 54:01, 55:01, 56:01, 58:01, 61:01N, 64:01N, 67:01, 69:01, 70:01, 76:01, 78:01, 79:01, 85:01, 87:01, 88:01, 91:01, 92:01, 93:01, 98:01, 104:01, 107:01, 110:01, 111:01, 114:01, 119:01, 124:01, 125:01, 130:01, 131:01, 132:01, 133:01, 136:01
18	180 bp	430 bp	*03:01:01-03:01:02, 06:01, 09:01, 14:01, 17:01, 20:01:01-20:01:02, 29:01, 35:01:01-35:01:02, 44:01, 61:01N, 64:01N, 69:01, 70:01, 76:01, 78:01, 86:01, 88:01, 91:01, 92:01, 98:01, 104:01, 111:01, 119:01, 124:01, 130:01, 131:01, 132:01
19	200 bp	430 bp	*03:01:01-03:01:02, 06:01, 11:01:01-11:01:02, 14:01, 20:01:01-20:01:02, 25:01, 29:01, 44:01, 45:01, 52:01, 56:01, 61:01N, 64:01N, 67:01, 69:01, 70:01, 74:01, 76:01, 78:01, 87:01, 91:01, 92:01, 104:01, 111:01, 114:01, 119:01, 124:01, 130:01, 136:01
20	235 bp	515 bp	*03:01:01-03:01:02, 09:01, 10:01, 13:01, 14:01, 25:01, 26:01:01-26:01:02, 29:01, 35:01:01-35:01:02, 37:01, 44:01, 45:01, 52:01, 54:01, 56:01, 61:01N, 67:01, 70:01, 76:01, 78:01, 79:01, 88:01, 92:01, 104:01, 107:01, 110:01, 111:01, 114:01, 119:01, 124:01, 125:01, 132:01, 133:01
21⁴	100 bp	430 bp	*03:01:01-03:01:02, 09:01, 14:01, 29:01, 35:01:01-35:01:02, 44:01, 50:01, 57:01, 61:01N, 70:01, 76:01, 78:01, 88:01, 92:01, 104:01, 111:01, 119:01, 124:01, 132:01
22	185 bp	430 bp	*04:01:01:01-04:01:03, 24:01, 33:01, 39:01, 49:01, 51:01, 66:01, 72:01, 81:01, 96:01, 99:01, 117:01, 120:01N, 121:01, 126:01, 128:01, 134:01
23¹⁰	145 bp, 210 bp	430 bp	*05:01:01-05:01:02, 34:01, 36:01, 38:01, 62:01, 63:01, 97:01, 100:01, 112:01, 114:01, 116:01, 127:01
24¹¹	195 bp, 265 bp	515 bp	*05:01:01-05:01:02, 08:01, 16:01, 19:01, 22:01, 31:01, 38:01, 57:01, 63:01, 65:01, 68:01, 84:01, 97:01, 101:01, 102:01, 103:01, 106:01
25¹²	190 bp, 265 bp	515 bp	*05:01:01-05:01:02, 21:01, 22:01, 36:01, 38:01, 44:01, 58:01, 63:01, 70:01, 97:01, 114:01
26^{4,6,13}	75 bp, 205 bp	515 bp	*06:01, 09:01, 17:01, 29:01, 44:01, 46:01, 64:01N, 86:01, 88:01, 101:01, 123:01, 131:01
27⁴	95 bp	430 bp	*06:01, 17:01, 20:01:01-20:01:02, 46:01, 64:01N, 69:01, 80:01, 86:01, 91:01, 98:01, 108:01, 130:01, 131:01
28¹⁴	200 bp, 260 bp	430 bp	*09:01, 10:01, 13:01, 17:01, 21:01, 26:01:01-26:01:02, 27:01, 30:01, 35:01:01-35:01:02, 36:01, 37:01, 54:01, 55:01, 58:01, 66:01, 79:01, 85:01, 86:01, 88:01, 93:01,

101.311-24/06 – including *Taq* pol., IFU-01
 101.311-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **09N**

Lot-specific information

			98:01, 107:01, 109:01, 110:01, 113:01, 117:01, 119:01, 125:01, 129:01, 131:01, 132:01, 133:01
29	220 bp	515 bp	*09:01, 10:01, 17:01, 30:01, 54:01, 55:01, 58:01, 86:01, 113:01, 131:01
30	185 bp	430 bp	*09:01, 14:01, 17:01, 35:01:01, 76:01, 86:01, 91:01, 98:01, 119:01, 131:01
31^{4,7}	120 bp	430 bp	*01:01:01-01:01:03, 05:01:01-05:01:02, 08:01, 10:01, 11:01:01-11:01:02, 13:01, 16:01, 19:01, 21:01, 22:01, 25:01, 26:01:01-26:01:02, 27:01, 30:01, 31:01, 36:01, 37:01, 38:01, 45:01, 52:01, 54:01, 55:01, 56:01, 58:01, 63:01, 65:01, 67:01, 68:01, 79:01, 84:01, 85:01, 87:01, 89:01, 90:01, 93:01, 97:01, 102:01, 103:01, 106:01, 107:01, 110:01, 114:01, 118:01, 122:01, 125:01, 127:01, 133:01, 136:01
32^{4,7,15}	75 bp, 195 bp	430 bp	*11:01:01-11:01:02, 78:01, 130:01
33¹⁶	140 bp, 205 bp	430 bp	*13:01, 33:01, 81:01, 107:01, 110:01, 117:01, 123:01, 133:01
34	210 bp	430 bp	*14:01, 45:01, 67:01, 76:01, 91:01, 119:01
35	195 bp	430 bp	*15:01, 74:01
36	260 bp	515 bp	*15:01, 18:01, 28:01, 34:01, 40:01, 53:01, 62:01
37¹⁷	185 bp, 255 bp	430 bp	*04:02:01:01-04:02:01:02, 08:01 [?] , 15:01, 21:01 [?] , 22:01 [?] , 24:01 [?] , 25:01 [?] , 27:01 [?] , 28:01, 29:01 [?] , 31:01 [?] , 32:01 [?] , 34:01 [?] , 35:01:01 [?] , 36:01 [?] , 37:01 [?] , 38:01 [?] , 40:01, 41:01:01 [?] -41:01:02 [?] , 44:01 [?] , 46:01 [?] , 47:01 [?] , 48:01 [?] , 49:01 [?] , 51:01 [?] , 52:01 [?] , 53:01, 54:01 [?] , 56:01 [?] , 57:01 [?] , 58:01 [?] , 59:01, 60:01 [?] , 61:01N [?] , 62:01 [?] , 63:01 [?] , 64:01N [?] , 65:01 [?] , 66:01 [?] , 67:01 [?] , 68:01 [?] , 69:01 [?] , 70:01 [?] , 71:01 [?] , 72:01 [?] , 73:01 [?] , 74:01, 75:01 [?] , 76:01 [?] , 77:01 [?] , 78:01 [?] , 79:01 [?] , 80:01 [?] , 81:01 [?] , 82:01 [?] , 83:01 [?] , 84:01 [?] , 86:01 [?] , 87:01 [?] , 89:01 [?] , 90:01 [?] , 91:01 [?] , 92:01 [?] , 93:01 [?] , 94:01 [?] , 95:01 [?] , 96:01 [?] , 97:01 [?] , 98:01 [?] , 99:01 [?] , 100:01 [?] , 101:01 [?] , 102:01 [?] , 103:01 [?] , 108:01 [?] , 109:01 [?] , 110:01 [?] , 111:01 [?] , 112:01 [?] , 113:01 [?] , 114:01 [?] , 115:01 [?] , 116:01 [?] , 117:01 [?] , 118:01 [?] , 119:01 [?] , 120:01N [?] , 122:01 [?] , 125:01 [?] , 126:01, 127:01 [?] , 128:01 [?] , 129:01 [?] , 130:01 [?] , 136:01 [?]
38⁴	100 bp	515 bp	*23:01, 34:01, 52:01, 55:01, 58:01, 62:01, 63:01, 67:01, 71:01, 87:01, 95:01, 103:01, 112:01, 121:01
39	210 bp	515 bp	*28:01, 31:01, 34:01, 57:01, 59:01, 72:01, 73:01, 95:01, 102:01, 108:01
40^{4,5,18}	95 bp, 180 bp	430 bp	*06:01 [?] , 08:01 [?] , 10:01 [?] , 11:01:01 [?] -11:01:02 [?] , 16:01 [?] , 17:01 [?] , 18:01 [?] , 20:01:01 [?] -20:01:02 [?] , 21:01 [?] , 22:01 [?] , 23:01 [?] , 24:01 [?] , 25:01 [?] , 26:01:01 [?] -26:01:02 [?] , 27:01 [?] , 28:01 [?] , 29:01 [?] , 30:01, 31:01 [?] , 32:01 [?] , 33:01 [?] , 34:01 [?] ,

101.311-24/06 – including *Taq* pol., IFU-01
 101.311-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **09N**

Lot-specific information

35:01:01 [?] -35:01:02 [?] , 36:01 [?] , 37:01 [?] , 38:01 [?] , 39:01 [?] , 40:01 [?] , 41:01:01 [?] -41:01:02 [?] , 44:01 [?] , 45:01 [?] , 46:01 [?] , 47:01 [?] , 48:01 [?] , 49:01 [?] , 50:01 [?] , 51:01 [?] , 52:01 [?] , 53:01 [?] , 54:01, 55:01 [?] , 56:01 [?] , 57:01 [?] , 58:01 [?] , 59:01 [?] , 60:01 [?] , 61:01N [?] , 62:01 [?] , 63:01 [?] , 64:01N [?] , 65:01 [?] , 66:01 [?] , 67:01 [?] , 68:01 [?] , 69:01 [?] , 70:01 [?] , 71:01 [?] , 72:01 [?] , 73:01 [?] , 74:01 [?] , 75:01 [?] , 76:01 [?] , 77:01 [?] , 78:01 [?] , 79:01 [?] , 80:01 [?] , 81:01 [?] , 82:01 [?] , 83:01 [?] , 84:01 [?] , 85:01 [?] , 86:01 [?] , 87:01 [?] , 88:01 [?] , 89:01 [?] , 90:01 [?] , 91:01 [?] , 92:01 [?] , 93:01 [?] , 94:01 [?] , 95:01 [?] , 96:01 [?] , 97:01 [?] , 98:01 [?] , 99:01 [?] , 100:01 [?] , 101:01 [?] , 102:01 [?] , 103:01 [?] , 107:01, 108:01 [?] , 109:01 [?] , 110:01 [?] , 111:01 [?] , 112:01 [?] , 113:01 [?] , 114:01 [?] , 115:01 [?] , 116:01 [?] , 117:01 [?] , 118:01 [?] , 119:01 [?] , 120:01N [?] , 121:01 [?] , 122:01 [?] , 123:01 [?] , 124:01 [?] , 125:01 [?] , 126:01 [?] , 127:01 [?] , 128:01 [?] , 129:01 [?] , 130:01 [?] , 131:01 [?] , 132:01 [?] , 133:01 [?] , 134:01 [?] , 136:01 [?]			
41^{4,19}	75 bp, 225 bp, 250 bp	515 bp	*31:01, 34:01, 77:01, 94:01, 134:01
42^{7,20}	140 bp, 185 bp, 245 bp, 285 bp	430 bp	*32:01, 82:01, 97:01, 128:01
43^{4,5,21}	65 bp, 190 bp	430 bp	*05:01:01-05:01:02, 06:01 [?] , 08:01 [?] , 10:01 [?] , 11:01:01 [?] - 11:01:02 [?] , 16:01 [?] , 17:01 [?] , 18:01 [?] , 19:01, 20:01:01 [?] - 20:01:02 [?] , 21:01 [?] , 22:01 [?] , 23:01 [?] , 24:01 [?] , 25:01 [?] , 26:01:01 [?] -26:01:02 [?] , 27:01 [?] , 28:01 [?] , 29:01 [?] , 30:01 [?] , 31:01 [?] , 32:01 [?] , 33:01 [?] , 34:01, 35:01:01 [?] -35:01:02 [?] , 36:01 [?] , 37:01 [?] , 38:01 [?] , 39:01 [?] , 40:01 [?] , 41:01:01 [?] - 41:01:02 [?] , 44:01 [?] , 45:01 [?] , 46:01 [?] , 47:01 [?] , 48:01 [?] , 49:01 [?] , 50:01 [?] , 51:01 [?] , 52:01 [?] , 53:01 [?] , 54:01 [?] , 55:01 [?] , 56:01 [?] , 57:01 [?] , 58:01 [?] , 59:01 [?] , 60:01 [?] , 61:01N [?] , 62:01, 63:01 [?] , 64:01N [?] , 65:01 [?] , 66:01 [?] , 67:01 [?] , 68:01 [?] , 69:01 [?] , 70:01 [?] , 71:01 [?] , 72:01 [?] , 73:01 [?] , 74:01 [?] , 75:01 [?] , 76:01 [?] , 77:01 [?] , 78:01 [?] , 79:01 [?] , 80:01 [?] , 81:01 [?] , 82:01 [?] , 83:01 [?] , 84:01 [?] , 85:01 [?] , 86:01 [?] , 87:01 [?] , 88:01 [?] , 89:01 [?] , 90:01 [?] , 91:01 [?] , 92:01 [?] , 93:01 [?] , 94:01 [?] , 95:01 [?] , 96:01 [?] , 97:01 [?] , 98:01 [?] , 99:01 [?] , 100:01 [?] , 101:01 [?] , 102:01 [?] , 103:01 [?] , 104:01, 108:01 [?] , 109:01 [?] , 110:01 [?] , 111:01 [?] , 112:01 [?] , 113:01 [?] , 114:01 [?] , 115:01 [?] , 116:01 [?] , 117:01 [?] , 118:01 [?] , 119:01 [?] , 120:01N [?] , 121:01 [?] , 122:01 [?] , 123:01 [?] , 124:01 [?] , 125:01 [?] , 126:01 [?] , 127:01 [?] , 128:01 [?] , 129:01 [?] , 130:01 [?] , 131:01 [?] , 132:01 [?] , 133:01 [?] , 134:01 [?] , 136:01 [?]
44^{4,22}	110 bp, 185 bp, 245 bp, 270 bp	430 bp	*38:01, 64:01N, 111:01, 118:01, 120:01N

101.311-24/06 – including *Taq* pol., IFU-01
 101.311-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **09N**

Lot-specific information

45 ^{4,23}	120 bp, 190 bp, 245 bp	430 bp	*39:01, 46:01, 49:01, 80:01, 86:01, 96:01, 108:01, 115:01, 117:01
46 ^{4,24}	70 bp, 210 bp	430 bp	*41:01:01-41:01:02, 60:01, 69:01, 83:01, 96:01, 136:01
47 ^{5,25}	165 bp, 255 bp	430 bp	*02:01:02-02:02, 04:01:01:01-04:01:03, 08:01 [?] , 17:01, 21:01 [?] , 22:01 [?] , 23:01, 24:01 [?] , 25:01 [?] , 26:01:01 [?] , 27:01 [?] , 28:01, 29:01 [?] , 30:01, 31:01 [?] , 32:01 [?] , 33:01, 34:01 [?] , 35:01:01 [?] , 36:01 [?] , 37:01 [?] , 38:01 [?] , 39:01, 40:01 [?] , 41:01:01 [?] -41:01:02 [?] , 44:01, 46:01 [?] , 47:01 [?] , 48:01 [?] , 49:01 [?] , 51:01 [?] , 52:01 [?] , 53:01 [?] , 54:01 [?] , 55:01, 56:01 [?] , 57:01 [?] , 58:01 [?] , 60:01 [?] , 61:01N [?] , 62:01 [?] , 63:01 [?] , 64:01N [?] , 65:01 [?] , 66:01 [?] , 67:01 [?] , 68:01 [?] , 69:01 [?] , 70:01 [?] , 71:01 [?] , 72:01 [?] , 73:01 [?] , 74:01 [?] , 75:01 [?] , 76:01 [?] , 77:01 [?] , 78:01 [?] , 79:01 [?] , 80:01 [?] , 81:01 [?] , 82:01 [?] , 83:01 [?] , 84:01 [?] , 86:01 [?] , 87:01 [?] , 89:01 [?] , 90:01 [?] , 91:01 [?] , 92:01 [?] , 93:01 [?] , 94:01 [?] , 95:01 [?] , 96:01 [?] , 97:01 [?] , 98:01 [?] , 99:01 [?] , 100:01 [?] , 101:01 [?] , 102:01 [?] , 103:01 [?] , 105:01, 106:01, 108:01 [?] , 109:01 [?] , 110:01 [?] , 111:01 [?] , 112:01 [?] , 113:01 [?] , 114:01, 115:01 [?] , 116:01 [?] , 117:01 [?] , 118:01 [?] , 119:01 [?] , 120:01N [?] , 121:01, 122:01 [?] , 123:01, 124:01, 125:01 [?] , 127:01 [?] , 128:01 [?] , 129:01 [?] , 130:01 [?] , 133:01, 134:01, 136:01 [?]
48 ^{4,26}	115 bp, 210 bp, 280 bp	430 bp	*56:01, 61:01N, 66:01, 76:01, 85:01, 92:01, 110:01, 125:01

¹Alleles are assigned by the presence of specific PCR product(s). However, the sizes of the specific PCR products may be helpful in the interpretation of DPB1 SSP typings.

When the primers in a primer mix can give rise to specific PCR products of more than one length this is indicated if the size difference is 20 base pairs or more. Size differences shorter than 20 base pairs are not given. For high resolution SSP kits the respective lengths of the specific PCR product(s) of the alleles amplified by these primer mixes are given.

Nonspecific amplifications, i.e. a ladder or a smear of bands, may sometimes be seen. GC-rich primers have a higher tendency of giving rise to nonspecific amplifications than other primers.

PCR fragments longer than the control bands may sometimes be observed. Such bands should be disregarded and do not influence the interpretation of the SSP typings.

PCR fragments migrating faster than the control bands, but slower than a 400 bp fragment may be seen in some gel read-outs. Such bands can be disregarded and do not influence the interpretation of the SSP typings.

Some primers may give rise to primer oligomer artifacts. Sometimes this phenomenon is an inherent feature of the primer pair(s) of a primer mix. More often it is due to other factors such as too low amount of DNA in the PCR reactions, taking too long time in setting up the PCR reactions, working at elevated room temperature or using thermal cyclers that are not pre-heated.

²The internal positive control primer pairs amplify segments of the human growth hormone gene. The two different control primer pairs give rise to either an internal positive control band of 430 base pairs, for most wells, or a band of 515 base pairs, for some wells.

Well number 1 contains the primer pair giving rise to the longer, 515 bp, internal positive control band in order to help in the correct orientation of the DPB1 typing.

In addition, wells number 2, 5, 8, 10 to 12, 17, 20, 24 to 26, 29, 36, 38, 39 and 41 contain the primer pair giving rise to the longer, 515 bp, internal positive control band in order to allow kit

101.311-24/06 – including *Taq* pol., IFU-01
101.311-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: 09N**Lot-specific information**

identification.

In the presence of a specific amplification the intensity of the control band often decreases.

³The 1st, 3rd and 4th exon of many DPB1 alleles is not sequenced and thus the primer pairs in primer mixes 37, 40, 43 and 47 may amplify other DPB1 alleles. We assume that unknown sequences in these exons and in the introns are conserved within loci and within allelic groups.

⁴Specific PCR fragments shorter than 125 base pairs are less intense and not as sharp as longer specific bands.

⁵Primer mixes 3, 40, 43 and 47 may yield less specific PCR products than the other DPB1 primer mixes.

⁶Primer mixes 3, 16 and 26 may give rise to primer dimer formation.

⁷Primer mixes 31, 32 and 42 may give rise to non-specific amplifications.

⁸Primer mix 14: Specific PCR fragment of 105 bp in the DPB1*02:02, 05:01:01-05:01:02, 22:01, 38:01, 97:01, 100:01 and 101:01 alleles. Specific PCR fragment of 175 bp in the DPB1*30:01 and 54:01 alleles. Specific PCR fragments of 105 and 175 bp in the DPB1*21:01, 36:01 and 114:01 alleles.

⁹Primer mix 16: Specific PCR fragment of 80 bp in the DPB1*03:01:01-03:01:02, 04:02:01:01-04:02:01:02, 14:01, 18:01, 20:01:01-20:01:02, 25:01, 28:01, 35:01:01-35:01:02, 45:01, 49:01, 50:01, 51:01, 53:01, 57:01, 59:01, 60:01, 61:01N, 68:01, 70:01, 73:01, 75:01, 76:01, 77:01, 78:01, 79:01, 80:01, 82:01, 83:01, 91:01, 92:01, 98:01, 104:01, 105:01, 108:01, 111:01, 116:01, 119:01, 122:01, 124:01, 129:01, 130:01 and 132:01 alleles. Specific PCR fragment of 165 bp in the DPB1*99:01 allele. Specific PCR fragment of 210 bp in the DPB1*127:01 allele.

¹⁰Primer mix 23: Specific PCR fragment of 145 bp in the DPB1*05:01:01-05:01:02, 34:01, 36:01, 38:01, 62:01, 63:01, 97:01, 100:01, 112:01, 114:01 and 116:01 alleles. Specific PCR fragment of 210 bp in the DPB1*127:01 allele.

¹¹Primer mix 24: Specific PCR fragment of 195 bp in the DPB1*101:01 allele. Specific PCR fragment of 265 bp in the DPB1*05:01:01-05:01:02, 08:01, 16:01, 19:01, 22:01, 31:01, 38:01, 57:01, 63:01, 65:01, 68:01, 84:01, 97:01, 102:01, 103:01 and 106:01 alleles.

¹²Primer mix 25: Specific PCR fragment of 190 bp in the DPB1*05:01:01-05:01:02, 21:01, 22:01, 36:01, 38:01, 44:01, 58:01, 63:01, 97:01 and 114:01 alleles. Specific PCR fragment of 265 bp in the DPB1*70:01 allele.

¹³Primer mix 26: Specific PCR fragment of 75 bp in the DPB1*06:01, 09:01, 17:01, 29:01, 44:01, 46:01, 64:01N, 86:01, 88:01, 101:01 and 131:01 alleles. Specific PCR fragment of 205 bp in the DPB1*123:01 allele.

¹⁴Primer mix 28: Specific PCR fragment of 200 bp in the DPB1*09:01, 10:01, 13:01, 17:01, 21:01, 26:01:01-26:01:02, 27:01, 30:01, 35:01:01-35:01:02, 36:01, 37:01, 54:01, 55:01, 58:01, 66:01, 79:01, 85:01, 86:01, 88:01, 93:01, 98:01, 107:01, 109:01, 110:01, 113:01, 117:01, 125:01, 129:01, 131:01, 132:01 and 133:01 alleles. Specific PCR fragment of 260 bp in the DPB1*119:01 allele.

¹⁵Primer mix 32: Specific PCR fragment of 75 bp in the DPB1*78:01 allele. Specific PCR fragment of 195 bp in the DBB1*11:01:01-11:01:02 and 130:01 alleles.

¹⁶Primer mix 33: Specific PCR fragment of 140 bp in the DPB1*13:01, 33:01, 81:01, 107:01, 110:01, 117:01 and 133:01 alleles. Specific PCR fragment of 205 bp in the DBB1*123:01 allele.

¹⁷Primer mix 37: Specific PCR fragment of 185 bp in the DPB1*15:01, 28:01, 40:01, 53:01 and 74:01 alleles. Specific PCR fragment of 255 bp in the DPB1*04:02:01:01-04:02:01:02, 59:01 and 126:01 alleles. The 3rd exon of many DPB1 alleles is not sequenced and thus the primer pair in primer mix 37 matching the DPB1*04:02 and 59:01 sequence motif in the 3rd exon may amplify other DPB1 alleles.

¹⁸Primer mix 40: Specific PCR fragment of 95 bp in the DPB1*107:01 allele. Specific PCR fragment of 180 bp in the DPB1*30:01 and 54:01 alleles.

The 1st exon of most DPB1 alleles is not sequenced and thus the primer pair in primer mix 40 matching the DPB1*107:01 sequence motif in the 1st exon may amplify other DPB1 alleles.

¹⁹Primer mix 41: Specific PCR fragment of 75 bp in the DPB1*77:01 allele. Specific PCR fragment of 250 bp in the DPB1*134:01 allele. Specific PCR fragment of 225 and 250 bp in the DPB1*31:01, 34:01 and 94:01 alleles.

²⁰Primer mix 42: Specific PCR fragment of 140 bp in the DPB1*97:01 allele. Specific PCR fragment of 185 bp in the DPB1*32:01 allele. Specific PCR fragment of 245 bp in the DPB1*128:01 allele. Specific PCR fragment of 285 bp in the DPB1*82:01 allele.

101.311-24/06 – including *Taq* pol., IFU-01
101.311-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: 09N**Lot-specific information**

²¹Primer mix 43: Specific PCR fragment of 65 bp in the DPB1*05:01:01, 05:01:02, 19:01 and 104:01 alleles. Specific PCR fragment of 190 bp in the DPB1*34:01 and 62:01 alleles. The 4th exon of most DPB1 alleles is not sequenced and thus the primer pair in primer mix 43 matching the DPB1*05:01:01, 05:01:02, 19:01 and 104:01 sequence motif in the 4th exon may amplify other DPB1 alleles.

²²Primer mix 44: Specific PCR fragment of 185 bp in the DPB1*111:01 and 120:01N alleles. Specific PCR fragment of 245 bp in the DPB1*38:01 allele. Specific PCR fragment of 270 bp in the DPB1*64:01N allele. Specific PCR fragment of 110 and 185 bp in the DPB1*118:01 allele.

²³Primer mix 45: Specific PCR fragment of 120 bp in the DPB1*46:01, 80:01, 86:01 and 108:01 alleles. Specific PCR fragment of 190 bp in the DPB1*39:01, 49:01, 96:01 and 117:01 alleles. Specific PCR fragment of 245 bp in the DPB1*115:01 allele.

²⁴Primer mix 46: Specific PCR fragment of 70 bp in the DPB1 *69:01 and 136:01 alleles. Specific PCR fragment of 210 bp in the DPB1*96:01 allele. Specific PCR fragment of 70 and 210 bp in the DPB1*41:01:01-41:01:02, 60:01 and 83:01 alleles.

²⁵Primer mix 47: Specific PCR fragment of 165 bp in the DPB1*44:01 and 114:01 alleles. Specific PCR fragment of 255 bp in the DPB1*02:01:02-02:02, 04:01:01-04:01:02, 17:01, 23:01, 28:01, 30:01, 33:01, 39:01, 55:01, 105:01, 106:01, 121:01, 123:01, 124:01, 133:01 and 134:01 alleles. The 3rd exon of many DPB1 alleles is not sequenced and thus the primer pair in primer mix 47 matching the DPB1*02:01:02-02:02, 04:01:01-04:01:02, 17:01, 23:01, 28:01, 30:01, 33:01, 39:01, 55:01, 105:01, 106:01, 121:01, 123:01, 124:01, 133:01 and 134:01 sequence motif in the 3rd exon may amplify other DPB1 alleles.

²⁶Primer mix 48: Specific PCR fragment of 115 bp in the DPB1*56:01, 66:01, 76:01, 92:01, 110:01 and 125:01 alleles. Specific PCR fragment of 210 bp in the DPB1*61:01N allele. Specific PCR fragment of 280 bp in the DPB1*85:01 allele.

‘w’, might be weakly amplified.

‘?’ , nucleotide sequence information not available for the primer matching sequence.

101.311-24/06 – including *Taq* pol., IFU-01
 101.311-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: 09N

Lot-specific information

INTERPRETATION TABLE																								
DPB1 SSP typing																								
Amplification patterns of the DPB1*01:01 to DPB1*136:01 alleles																								
	Well																							
	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
Length of spec. PCR product(s)	245	220	100	140	235	165	160	210	165	220	185	80	180	105	200	80	260	180	200	235	100	185	145	195
														175		165							210	265
																210								
Length of int. pos. control ¹	515	515	430	430	515	430	430	515	430	515	515	515	430	430	430	430	515	430	430	515	430	430	430	515
5'-primer(s) ²	5'-TgT ³ 9(112)	5'-TgT ³ 9(112)	5'-CgC ³ 36(194)	5'-CgC ³ 36(194)	5'-ggg ³ 11(119)	5'-gTA ³ 35(191)	5'-CgC ³ 36(194)	5'-TTT ³ 9(113)	5'-AgC ³ 35(190)	5'-TTT ³ 9(113)	5'-gTT ³ 35(191)	5'-gAT ³ 55(252)	5'-TTT ³ 9(113)	5'-gTT ³ 11(119)	5'-gAg ³ 28(169)	5'-gAT ³ 12(121)	5'-gTT ³ 11(119)	5'-gTT ³ 11(119)	5'-gTT ³ 11(119)	5'-gTT ³ 11(119)	5'-gAC ³ 57(258)	5'-CgC ³ 36(194)	5'-gAT ³ 12(121)	5'-TTT ³ 9(113)
			5'-CgC ³ 36(194)	5'-CgC ³ 36(194)	5'-ggg ³ 11(119)	5'-gTA ³ 35(191)	5'-CgC ³ 36(194)		5'-CgT ³ 36(194)					5'-gTT ³ 11(119)	5'-AgC ³ 35(190)	5'-gAg ³ 28(169)	5'-gTT ³ 11(119)	5'-gTT ³ 11(119)	5'-gTT ³ 11(119)	5'-gTT ³ 11(119)		5'-CgC ³ 36(194)	5'-AgC ³ 35(190)	
														5'-AgC ³ 35(190)		5'-TgA ³ 56(254)								
3'-primer(s) ³	5'-CAT ³ 76(315)	5'-CTT ³ 69(292)	5'-CAG ³ 55(251)	5'-CTT ³ 69(292)	5'-CAT ³ 76(315)	5'-CAT ³ 76(315)	5'-CAT ³ 76(315)	5'-gAT ³ 65(280)	5'-TAT ³ 76(313)	5'-CTC ³ 69(292)	5'-gTT ³ 84(337)	5'-CTC ³ 69(292)	5'-gCC ³ 55(252)	5'-gCC ³ 55(252)	5'-CgC ³ 84(338)	5'-CTT ³ 69(292)	5'-CgT ³ 84(338)	5'-TAG ³ 57(258)	5'-gAg ³ 65(280)	5'-CAT ³ 76(315)	5'-CAT ³ 76(315)	5'-CgC ³ 84(338)	5'-CTT ³ 69(292)	5'-gTg ³ 60(265)
				5'-CTT ³ 69(292)				5'-ggT ³ 65(281)	5'-CAT ³ 76(313)		5'-CgC ³ 84(338)											5'-CTT ³ 69(292)	5'-CgT ³ 84(338)	
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24



101.311-24/06 – including *Taq* pol., IFU-01
101.311-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: 09N

Lot-specific information

[illegible]

101.311-24/06 – including *Taq* pol., IFU-01
 101.311-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **09N**

Lot-specific information

Length of spec. PCR product(s)	245	220	100	140	235	165	160	210	165	220	185	80	180	105	200	80	260	180	200	235	100	185	145	195
														175		165							210	265
																210								
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
DPB1 allele ⁴																								
*01:01:01-01:01:03	1	2	3	4	5	6	7																	
*02:01:02, 02:01:04								8	9	10	11	12												
02:01:05, 02:01:07																								
*02:01:03								8	9	10	11	w												
*02:01:06								8	w	10	11	12												
*02:02								8	9	10			13	14	15									
*03:01:01-03:01:02	1	2														16	17	18	19	20	21			
*04:01:01:01-04:01:03			3	4				8			11											22		
*04:02:01:01-04:02:01:02								8	9		11					16								
*05:01:01-05:01:02								8	9				13	14									23	24
*06:01									9			12					17	18	19					
*08:01					5			8		10		12												24
*09:01												12					17	18		20	21			
*10:01												12					17			20				
*11:01:01-11:01:02			3														17		19					
*13:01	1		3			6	7										17			20				
*14:01																16	17	18	19	20	21			
*15:01			3																					
*16:01								8	9	10		12												24
*17:01									9			12					17	18						
*18:01		2							9							16								
*19:01					5			8	9	10			13											24
*20:01:01-20:01:02		2							9							16	17	18	19					
*21:01									9				14				17							
*22:01								8	9	10			13	14										24
*23:01								8	9		11													
*24:01				4				8			11		13									22		
*25:01	1	2														16	17		19	20				
*26:01:01	1	2	3	4		6	7										17			20				
*26:01:02	1	2	3	4		6	7										17			20				
*27:01		2	3	4													17							
*28:01				4												16								
*29:01	1											12					17	18	19	20	21			
*30:01									9				14				17							
*31:01			3																					24
*32:01								8	9	10	11	12												
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24

101.311-24/06 – including *Taq* pol., IFU-01
 101.311-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **09N**

Lot-specific information

190	75	95	200	220	185	120	75	140	210	195	260	185	100	210	95	75	140	65	110	120	70	165	115	Length of spec. PCR product(s)
265	205		260				200	205				255			180	225	185	190	185	190	205	255	210	
																250	245		245	245			280	
																	285		270					
25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	Well No.
																								DPB1 allele ⁴
						31																		*01:01:01-01:01:03
																						47		*02:01:02, 02:01:04 02:01:05, 02:01:07
																						47		*02:01:03
																						47		*02:01:06
																						47		*02:02
																								*03:01:01-03:01:02
																						47		*04:01:01-01- 04:01:03
												37												*04:02:01-01- 04:02:01:02
25						31												43						*05:01:01-05:01:02
	26	27													?			?						*06:01
						31						?			?			?				?		*08:01
	26		28	29	30																			*09:01
			28	29		31									?			?						*10:01
						31	32								?			?						*11:01:01-11:01:02
			28			31		33																*13:01
					30				34															*14:01
									35	36	37													*15:01
						31									?			?						*16:01
	26	27	28	29	30										?			?				47		*17:01
										36					?			?						*18:01
						31												43						*19:01
		27													?			?						*20:01:01-20:01:02
25			28			31						?			?			?				?		*21:01
25						31						?			?			?				?		*22:01
												38			?			?				47		*23:01
												?			?			?				?		*24:01
						31						?			?			?				?		*25:01
			28			31									?			?				?		*26:01:01
			28			31									?			?						*26:01:02
			28			31						?			?			?				?		*27:01
									36	37	39	?			?			?				47		*28:01
	26											?			?			?				?		*29:01
			28	29		31									40			?				47		*30:01
						31						?	39	?	41			?				?		*31:01
												?			?			42	?			?		*32:01
25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	Well No.

101.311-24/06 – including *Taq* pol., IFU-01
 101.311-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **09N**

Lot-specific information

Length of spec. PCR product(s)	245	220	100	140	235	165	160	210	165	220	185	80	180	105	200	80	260	180	200	235	100	185	145	195
														175		165							210	265
																210								
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
*33:01			3					8		10	11											22		
*34:01									9														23	
*35:01:01																16	17	18		20	21			
*35:01:02																16	17	18		20	21			
*36:01		2							9					14			17						23	
*37:01	1											12					17			20				
*38:01								8	9				13	14									23	24
*39:01			3	4				8														22		
*40:01			3	4				8																
*41:01:01-41:01:02									9	10	11	12												
*44:01	1											12					17	18	19	20	21			
*45:01																16	17		19	20				
*46:01								8	9	10	11	12												
*47:01								8	9	10	11		13											
*48:01								8	9	10		12			15									
*49:01				4				8								16						22		
*50:01	1	2			5											16					21			
*51:01				4				8			11					16						22		
*52:01	1	2															17		19	20				
*53:01				4				8								16								
*54:01													14				17			20				
*55:01									9								17							
*56:01	1	2	3	4			7										17		19	20				
*57:01					5											16					21		24	
*58:01									9								17							
*59:01									9		11					16								
*60:01								8	9		11					16								
*61:01N	1	2														16	17	18	19	20	21			
*62:01								8	9														23	
*63:01								8	9														23	24
*64:01N									9			12					17	18	19					
*65:01			3	4	5	6	7	8																24
*66:01			3	4							11											22		
*67:01																	17		19	20				
*68:01					5			8								16								24
*69:01									9								17	18	19					
*70:01																16	17	18	19	20	21			
*71:01								8	9	10	11													
*72:01			3	4							11											22		
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24

101.311-24/06 – including *Taq* pol., IFU-01
 101.311-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **09N**

Lot-specific information

190	75	95	200	220	185	120	75	140	210	195	260	185	100	210	95	75	140	65	110	120	70	165	115	Length of spec. PCR product(s)
265	205		260				200	205				255			180		225	185	190	185	190	205	255	210
																	250	245		245	245			280
																		285		270				
25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	Well No.
								33							?			?				47		*33:01
											36	?	38	39	?	41		43				?		*34:01
			28		30							?			?			?				?		*35:01:01
			28												?			?						*35:01:02
25			28			31						?			?			?				?		*36:01
			28			31						?			?			?				?		*37:01
25						31						?			?			?	44			?		*38:01
															?			?		45		47		*39:01
											36	37			?			?				?		*40:01
												?			?			?			46	?		*41:01:01-41:01:02
25	26											?			?			?				47		*44:01
						31				34					?			?						*45:01
	26	27										?			?			?		45		?		*46:01
												?			?			?				?		*47:01
												?			?			?				?		*48:01
												?			?			?		45		?		*49:01
															?			?						*50:01
												?			?			?				?		*51:01
						31						?	38		?			?				?		*52:01
											36	37			?			?				?		*53:01
			28	29		31						?			40			?				?		*54:01
			28	29		31							38		?			?				47		*55:01
						31						?			?			?				?	48	*56:01
												?		39	?			?				?		*57:01
25			28	29		31						?	38		?			?				?		*58:01
												37		39	?			?						*59:01
												?			?			?			46	?		*60:01
												?			?			?				?	48	*61:01N
											36	?	38		?			43				?		*62:01
25						31						?	38		?			?				?		*63:01
	26	27										?			?			?	44			?		*64:01N
						31						?			?			?				?		*65:01
			28									?			?			?				?	48	*66:01
						31				34		?	38		?			?				?		*67:01
						31						?			?			?				?		*68:01
		27										?			?			?			46	?		*69:01
25												?			?			?				?		*70:01
												?	38		?			?				?		*71:01
												?		39	?			?				?		*72:01
25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	Well No.

101.311-24/06 – including *Taq* pol., IFU-01
 101.311-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **09N**

Lot-specific information

Length of spec. PCR product(s)	245	220	100	140	235	165	160	210	165	220	185	80	180	105	200	80	260	180	200	235	100	185	145	195
														175		165							210	265
																210								
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
*73:01					5						11					16								
*74:01			3															19						
*75:01					5			8			11					16								
*76:01				4			7									16	17	18	19	20	21			
*77:01								8	9		11					16								
*78:01	1	2														16	17	18	19	20	21			
*79:01	1	2														16	17			20				
*80:01								8	9		11					16								
*81:01								8		10	11	12										22		
*82:01								8	9		11					16								
*83:01									9		11					16								
*84:01					5			8																24
*85:01		2	3	4													17							
*86:01									9		11	12						18						
*87:01		2							9								17		19					
*88:01	1											12					17	18		20	21			
*89:01		2	3	4																				
*90:01	1	2	3	4	5		7																	
*91:01									9							16	17	18	19					
*92:01	1	2		4			7									16	17	18	19	20	21			
*93:01									9			12					17							
*94:01								8	9		11													
*95:01									9	10					15									
*96:01			3	4				8														22		
*97:01								8	9				13	14									23	24
*98:01									9							16	17	18						
*99:01			3	4				8			11				15	16						22		
*100:01								8	9				13	14	15								23	
*101:01								8	9	10			13	14	15									24
*102:01			3	4																				24
*103:01								8	9	10														24
*104:01	1	2														16	17	18	19	20	21			
*105:01								8	9		11					16								
*106:01					5			8	9	10			13											24
*107:01	1		3			6	7										17			20				
*108:01									9		11					16								
*109:01									9		11	12												
*110:01	1		3				7										17			20				
*111:01	1	2														16	17	18	19	20	21			
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24

101.311-24/06 – including *Taq* pol., IFU-01
 101.311-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **09N**

Lot-specific information

190	75	95	200	220	185	120	75	140	210	195	260	185	100	210	95	75	140	65	110	120	70	165	115	Length of spec. PCR product(s)
265	205		260				200	205				255			180	225	185	190	185	190	205	255	210	
																250	245		245	245			280	
																	285		270					
25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	Well No.
												?		39	?			?				?		*73:01
										35		37			?			?				?		*74:01
												?			?			?				?		*75:01
					30				34			?			?			?				?	48	*76:01
												?			?	41		?				?		*77:01
							32					?			?			?				?		*78:01
			28			31						?			?			?				?		*79:01
		27										?			?			?		45		?		*80:01
								33				?			?			?				?		*81:01
												?			?		42	?				?		*82:01
												?			?			?			46	?		*83:01
						31						?			?			?				?		*84:01
			28			31						?			?			?					48	*85:01
	26	27	28	29	30							?			?			?		45		?		*86:01
						31						?	38		?			?				?		*87:01
	26		28												?			?						*88:01
						31						?			?			?				?		*89:01
						31						?			?			?				?		*90:01
		27			30				34			?			?			?				?		*91:01
												?			?			?				?	48	*92:01
			28			31						?			?			?				?		*93:01
												?			?	41		?				?		*94:01
												?	38	39	?			?				?		*95:01
												?			?			?		45	46	?		*96:01
25						31						?			?		42	?				?		*97:01
		27	28		30							?			?			?				?		*98:01
												?			?			?				?		*99:01
												?			?			?				?		*100:01
	26											?			?			?				?		*101:01
						31						?		39	?			?				?		*102:01
						31						?	38		?			?				?		*103:01
																		43						*104:01
																						47		*105:01
						31																47		*106:01
			28			31	33								40									*107:01
		27										?		39	?			?		45		?		*108:01
			28									?			?			?				?		*109:01
			28			31	33					?			?			?				?	48	*110:01
												?			?			?	44			?		*111:01
25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	Well No.

101.311-24/06 – including *Taq* pol., IFU-01
 101.311-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **09N**

Lot-specific information

Length of spec.	245	220	100	140	235	165	160	210	165	220	185	80	180	105	200	80	260	180	200	235	100	185	145	195
PCR product(s)														175		165							210	265
																210								
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
*112:01								8	9						15								23	
*113:01									9		11	12												
*114:01	1	2												14			17		19	20			23	
*115:01								8	9	10	11	12												
*116:01								8	9						15	16							23	
*117:01			3																			22		
*118:01	1	2	3	4	5	6	7																	
*119:01											11				16		17	18	19	20	21			
*120:01N			3	4				8			11											22		
*121:01			3	4				8			11											22		
*122:01	1	2		4	5	6	7								16									
*123:01								8	9	10	11	12												
*124:01	1	2													16		17	18	19	20	21			
*125:01			3	4			7										17			20				
*126:01			3	4				8			11											22		
*127:01	1	2	3	4	5	6	7								16								23	
*128:01			3	4				8			11											22		
*129:01									9		11				16									
*130:01		2							9						16		17	18	19					
*131:01									9		12						17	18						
*132:01	1	2													16		17	18		20	21			
*133:01	1		3			6	7										17			20				
*134:01			3	4				8			11											22		
*136:01									9								17		19					
DPB1 allele ⁴																								
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24

¹The internal positive control primer pairs amplify segments of the human growth hormone gene. The two different control primer pairs give rise to either an internal positive control band of 430 base pairs, for most wells, or a band of 515 base pairs, for some wells.

Well number 1 contains the primer pair giving rise to the longer, 515 bp, internal positive control band in order to help in the correct orientation of the DPB1 typing.

In addition, wells number 2, 5, 8, 10 to 12, 17, 20, 24 to 26, 29, 36, 38, 39 and 41 contain the primer pair giving rise to the longer, 515 bp, internal positive control band in order to allow kit identification.

In the presence of a specific amplification the intensity of the control band often decreases.

²The codon, and in parenthesis the nucleotide, in the 1st, 2nd, 3rd or 4th exon or in the 1st intron, matching the specificity-determining 3'-end of the primer is given. Codon and nucleotide numbering as on the on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

³The codon, and in parenthesis the nucleotide, in the 1st, 2nd, 3rd or 4th exon, matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Codon and nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

101.311-24/06 – including *Taq* pol., IFU-01
101.311-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: **09N**

Lot-specific information

Left primer				Left specific information																Length of spec. PCR product(s)				
190	75	95	200	220	185	120	75	140	210	195	260	185	100	210	95	75	140	65	110	120	70	165	115	
265	205		260				200	205				255			180	225	185	190	185	190	205	255	210	
																250	245		245	245			280	
																	285		270					
25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	Well No.
												?	38		?			?				?		*112:01
			28	29								?			?			?				?		*113:01
25						31						?			?			?				47		*114:01
												?			?			?		45		?		*115:01
												?			?			?				?		*116:01
			28					33				?			?			?		45		?		*117:01
						31						?			?			?	44			?		*118:01
			28	30				34				?			?			?				?		*119:01
												?			?			?	44			?		*120:01N
													38		?			?				47		*121:01
						31						?			?			?				?		*122:01
26								33							?			?				47		*123:01
															?			?				47		*124:01
			28			31						?			?			?				?	48	*125:01
												37			?			?						*126:01
						31						?			?			?				?		*127:01
												?			?		42	?				?		*128:01
			28									?			?			?				?		*129:01
		27					32					?			?			?				?		*130:01
26	27	28	29	30											?			?						*131:01
			28												?			?						*132:01
			28			31		33							?			?				47		*133:01
															?	41		?				47		*134:01
						31						?			?			?			46	?		*136:01
																								DPB1 allele ⁴
25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	Well No.

⁴The sequence of the DPB1*02011 allele has been shown to be identical to DPB1*02:01:02.

DPB1*0701 has never been assigned.

DPB1*1201 has never been assigned.

The sequence of the DPB1*4201 allele has been shown to be identical to DPB1*31:01.

The sequence of the DPB1*4301 allele has been shown to be identical to DPB1*28:01.

⁵Primer mix 14: Specific PCR fragment of 105 bp in the DPB1*02:02, 05:01:01-05:01:02, 22:01, 38:01, 97:01, 100:01 and 101:01 alleles. Specific PCR fragment of 175 bp in the DPB1*30:01 and 54:01 alleles. Specific PCR fragments of 105 and 175 bp in the DPB1*21:01, 36:01 and 114:01 alleles.

Primer mix 16: Specific PCR fragment of 80 bp in the DPB1*03:01:01-03:01:02, 04:02:01:01-04:02:01:02, 14:01, 18:01, 20:01:01-20:01:02, 25:01, 28:01, 35:01:01-35:01:02, 45:01, 49:01, 50:01, 51:01, 53:01, 57:01, 59:01, 60:01, 61:01N, 68:01, 70:01, 73:01, 75:01, 76:01, 77:01, 78:01, 79:01, 80:01, 82:01, 83:01, 91:01, 92:01, 98:01, 104:01, 105:01, 108:01, 111:01, 116:01, 119:01, 122:01, 124:01, 129:01, 130:01 and 132:01 alleles. Specific PCR fragment of 165 bp in the DPB1*99:01 allele. Specific PCR fragment of 210 bp in the DPB1*127:01 allele.

101.311-24/06 – including *Taq* pol., IFU-01
101.311-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: 09N**Lot-specific information**

Primer mix 23: Specific PCR fragment of 145 bp in the DPB1*05:01:01-05:01:02, 34:01, 36:01, 38:01, 62:01, 63:01, 97:01, 100:01, 112:01, 114:01 and 116:01 alleles. Specific PCR fragment of 210 bp in the DPB1*127:01 allele.

Primer mix 24: Specific PCR fragment of 195 bp in the DPB1*101:01 allele. Specific PCR fragment of 265 bp in the DPB1*05:01:01-05:01:02, 08:01, 16:01, 19:01, 22:01, 31:01, 38:01, 57:01, 63:01, 65:01, 68:01, 84:01, 97:01, 102:01, 103:01 and 106:01 alleles.

Primer mix 25: Specific PCR fragment of 190 bp in the DPB1*05:01:01-05:01:02, 21:01, 22:01, 36:01, 38:01, 44:01, 58:01, 63:01, 97:01 and 114:01 alleles. Specific PCR fragment of 265 bp in the DPB1*70:01 allele.

Primer mix 26: Specific PCR fragment of 75 bp in the DPB1*06:01, 09:01, 17:01, 29:01, 44:01, 46:01, 64:01N, 86:01, 88:01, 101:01 and 131:01 alleles. Specific PCR fragment of 205 bp in the DPB1*123:01 allele.

Primer mix 28: Specific PCR fragment of 200 bp in the DPB1*09:01, 10:01, 13:01, 17:01, 21:01, 26:01:01-26:01:02, 27:01, 30:01, 35:01:01-35:01:02, 36:01, 37:01, 54:01, 55:01, 58:01, 66:01, 79:01, 85:01, 86:01, 88:01, 93:01, 98:01, 107:01, 109:01, 110:01, 113:01, 117:01, 125:01, 129:01, 131:01, 132:01 and 133:01 alleles. Specific PCR fragment of 260 bp in the DPB1*119:01 allele.

Primer mix 32: Specific PCR fragment of 75 bp in the DPB1*78:01 allele. Specific PCR fragment of 195 bp in the DBB1*11:01:01-11:01:02 and 130:01 alleles.

Primer mix 33: Specific PCR fragment of 140 bp in the DPB1*13:01, 33:01, 81:01, 107:01, 110:01, 117:01 and 133:01 alleles. Specific PCR fragment of 205 bp in the DBB1*123:01 allele.

Primer mix 37: Specific PCR fragment of 185 bp in the DPB1*15:01, 28:01, 40:01, 53:01 and 74:01 alleles. Specific PCR fragment of 255 bp in the DPB1*04:02:01:01-04:02:01:02, 59:01 and 126:01 alleles. The 3rd exon of many DPB1 alleles is not sequenced and thus the primer pair in primer mix 37 matching the DPB1*04:02 and 59:01 sequence motif in the 3rd exon may amplify other DPB1 alleles.

Primer mix 40: Specific PCR fragment of 95 bp in the DPB1*107:01 allele. Specific PCR fragment of 180 bp in the DPB1*30:01 and 54:01 alleles.

The 1st exon of most DPB1 alleles is not sequenced and thus the primer pair in primer mix 40 matching the DPB1*107:01 sequence motif in the 1st exon may amplify other DPB1 alleles.

Primer mix 41: Specific PCR fragment of 75 bp in the DPB1*77:01 allele. Specific PCR fragment of 250 bp in the DPB1*134:01 allele. Specific PCR fragment of 225 and 250 bp in the DPB1*31:01, 34:01 and 94:01 alleles.

Primer mix 42: Specific PCR fragment of 140 bp in the DPB1*97:01 allele. Specific PCR fragment of 185 bp in the DPB1*32:01 allele. Specific PCR fragment of 245 bp in the DPB1*128:01 allele. Specific PCR fragment of 285 bp in the DPB1*82:01 allele.

Primer mix 43: Specific PCR fragment of 65 bp in the DPB1*05:01:01, 05:01:02, 19:01 and 104:01 alleles. Specific PCR fragment of 190 bp in the DPB1*34:01 and 62:01 alleles. The 4th exon of most DPB1 alleles is not sequenced and thus the primer pair in primer mix 43 matching the DPB1*05:01:01, 05:01:02, 19:01 and 104:01 sequence motif in the 4th exon may amplify other DPB1 alleles.

Primer mix 44: Specific PCR fragment of 185 bp in the DPB1*111:01 and 120:01N alleles. Specific PCR fragment of 245 bp in the DPB1*38:01 allele. Specific PCR fragment of 270 bp in the DPB1*64:01N allele. Specific PCR fragment of 110 and 185 bp in the DPB1*118:01 allele.

Primer mix 45: Specific PCR fragment of 120 bp in the DPB1*46:01, 80:01, 86:01 and 108:01 alleles. Specific PCR fragment of 190 bp in the DPB1*39:01, 49:01, 96:01 and 117:01 alleles. Specific PCR fragment of 245 bp in the DPB1*115:01 allele.

Primer mix 46: Specific PCR fragment of 70 bp in the DPB1 *69:01 and 136:01 alleles. Specific PCR fragment of 210 bp in the DPB1*96:01 allele. Specific PCR fragment of 70 and 210 bp in the DPB1*41:01:01-41:01:02, 60:01 and 83:01 alleles.

Primer mix 47: Specific PCR fragment of 165 bp in the DPB1*44:01 and 114:01 alleles. Specific PCR fragment of 255 bp in the DPB1*02:01:02-02:02:02, 04:01:01-04:01:02, 17:01, 23:01, 28:01, 30:01, 33:01, 39:01, 55:01, 105:01, 106:01, 121:01, 123:01, 124:01, 133:01 and 134:01 alleles. The 3rd exon of many DPB1 alleles is not sequenced and thus the primer pair in primer mix 47 matching the DPB1*02:01:02-02:02:02, 04:01:01-04:01:02, 17:01, 23:01, 28:01, 30:01, 33:01, 39:01, 55:01, 105:01, 106:01, 121:01, 123:01, 124:01, 133:01 and 134:01 sequence motif in the 3rd exon may amplify other DPB1 alleles.

Primer mix 48: Specific PCR fragment of 115 bp in the DPB1*56:01, 66:01, 76:01, 92:01, 110:01 and 125:01 alleles. Specific PCR fragment of 210 bp in the DPB1*61:01N allele. Specific PCR fragment of 280 bp in the DPB1*85:01 allele.

‘w’, might be weakly amplified.

‘?’, nucleotide sequence information not available for the primer matching sequence.

101.311-24/06 – including *Taq* pol., IFU-01
 101.311-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **09N**

Lot-specific information

CELL LINE VALIDATION SHEET															
DPB1 SSP kit															
					Well										
					1	2	3	4	5	6	7	8	9	10	11
					12	13	14	15	16						
				Lot No.:	201184901	201184902	201294303	201294304	201184905	201184906	201294307	201184908	201184909	201184910	201184911
					201184912	201184913	201184914	201184915	201184916						
IHC cell line					DPB1										
1	9001	SA	*04:02		-	-	-	-	-	-	-	+	+	-	+
2	9280	LK707	*04:01		-	-	+	+	-	-	-	+	-	-	+
3	9011	E4181324	*02:01	*04:01	-	-	+	+	-	-	-	+	+	+	+
4	9275	GU373	*04:02	*17:01	-	-	-	-	-	-	-	+	+	-	+
5	9009	KAS011	*04:01	*14:01	-	-	+	+	-	-	-	+	-	-	+
6	9353	SM	*05:01		-	-	-	-	-	-	-	+	+	-	-
7	9020	QBL	*02:02		-	-	-	-	-	-	-	+	+	-	+
8	9025	DEU	*04:01		-	-	+	+	-	-	-	+	-	-	-
9	9026	YAR	*04:01		-	-	+	+	-	-	-	+	-	-	-
10	9107	LKT3	*05:01		-	-	-	-	-	-	-	+	+	-	-
11	9051	PITOUT	*04:01		-	-	+	+	-	-	-	+	-	-	-
12	9052	DBB	*04:01		-	-	+	+	-	-	-	+	-	-	-
13	9004	JESTHOM	*04:01		-	-	+	+	-	-	-	+	-	-	-
14	9071	OLGA	*03:01	*04:02	+	+	-	-	-	-	-	+	+	-	+
15	9075	DKB	*04:01		-	-	+	+	-	-	-	+	-	-	-
16	9037	SWEIG007	*04:02		-	-	-	-	-	-	-	+	-	+	-
17	9282	CTM3953540	*01:01	*02:01	+	+	+	+	+	+	+	+	+	+	+
18	9257	32367	*04:02	*18:01	-	+	-	-	-	-	-	+	+	-	+
19	9038	BM16	*02:01		-	-	-	-	-	-	-	+	+	+	+
20	9059	SLE005	*03:01		+	+	-	-	-	-	-	-	-	-	+
21	9064	AMALA	*04:02		-	-	-	-	-	-	-	+	+	-	+
22	9056	KOSE	*02:01	*13:01	+	-	+	-	-	+	+	+	+	+	+
23	9124	IHL	*02:01		-	-	-	-	-	-	-	+	+	+	+
24	9035	JBUSH	*04:01		-	-	+	+	-	-	-	+	-	-	-
25	9049	IBW9	*01:01		+	+	+	+	+	+	+	-	-	-	-
26	9285	WT49	*04:01		-	-	+	+	-	-	-	+	-	-	-
27	9191	CH1007	*03:01	*13:01	+	+	+	-	-	+	+	-	-	-	+
28	9320	BEL5GB	*02:01	*04:01	-	-	+	+	-	-	-	+	+	+	+
29	9050	MOU	*02:01		-	-	-	-	-	-	-	+	+	+	+
30	9021	RSH	*01:01	*105:01	+	+	+	+	+	+	+	+	+	+	+
31	9019	DUCAF	*02:02		-	-	-	-	-	-	-	+	+	-	+
32	9297	HAG	*02:01	*03:01	+	+	-	-	-	-	-	+	+	+	+
33	9098	MT14B	*04:02		-	-	-	-	-	-	-	+	-	+	+
34	9104	DHIF	*04:01		-	-	+	+	-	-	-	+	-	-	-
35	9302	SSTO	*04:01		-	-	+	+	-	-	-	+	-	-	-
36	9024	KT17	*05:01		-	-	-	-	-	-	-	+	-	-	+
37	9065	HHKB	*04:01		-	-	+	+	-	-	-	+	-	-	-
38	9099	LZL	*04:02		-	-	-	-	-	-	-	+	+	-	+
39	9315	CML	*04:01	*14:01	-	-	+	+	-	-	-	+	-	-	+
40	9134	WHONP199	*05:01		-	-	-	-	-	-	-	+	-	-	+
41	9055	H0301	*05:01		-	-	-	-	-	-	-	+	+	-	-
42	9066	TAB089	*02:02		-	-	-	-	-	-	-	+	+	-	+
43	9076	T7526	*13:01		+	-	+	-	-	+	+	-	-	-	-
44	9057	TEM	*04:01		-	-	+	+	-	-	-	+	-	-	-
45	9239	SHJO	*03:01	*04:02	+	+	-	-	-	-	-	+	+	-	+
46	9013	SCHU	*04:02		-	-	-	-	-	-	-	+	-	+	+
47	9045	TUBO	*02:01	*04:02	-	-	-	-	-	-	-	+	+	+	+
48	9303	TER-ND	*02:01	*04:01	-	-	+	+	-	-	-	+	+	+	+

101.311-24/06 – including *Taq* pol., IFU-01
 101.311-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **09N**

Lot-specific information

CELL LINE VALIDATION SHEET																					
DPB1 SSP kit																					
					Well																
					17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	
				Lot No.:	201184917	201184918	201184919	201184920	201184921	201294322	201184923	201184924	201184925	201184926	201184927	201184928	201184929	201184930	201294331	201184932	
	IHWC cell line		DPB1																		
1	9001	SA	*04:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
2	9280	LK707	*04:01		-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	
3	9011	E4181324	*02:01	*04:01	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	
4	9275	GU373	*04:02	*17:01	+	+	-	-	-	-	-	-	-	+	+	+	+	+	-	-	
5	9009	KAS011	*04:01	*14:01	+	+	+	+	+	+	-	-	-	-	-	-	-	+	-	-	
6	9353	SM	*05:01		-	-	-	-	-	-	+	+	+	-	-	-	-	-	+	-	
7	9020	QBL	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
8	9025	DEU	*04:01		-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	
9	9026	YAR	*04:01		-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	
10	9107	LKT3	*05:01		-	-	-	-	-	-	+	+	+	-	-	-	-	-	+	-	
11	9051	PITOUT	*04:01		-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	
12	9052	DBB	*04:01		-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	
13	9004	JESTHOM	*04:01		-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	
14	9071	OLGA	*03:01	*04:02	+	+	+	+	+	-	-	-	-	-	-	-	-	-	-	-	
15	9075	DKB	*04:01		-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	
16	9037	SWEIG007	*04:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
17	9282	CTM3953540	*01:01	*02:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
18	9257	32367	*04:02	*18:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
19	9038	BM16	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
20	9059	SLE005	*03:01		+	+	+	+	+	-	-	-	-	-	-	-	-	-	-	-	
21	9064	AMALA	*04:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
22	9056	KOSE	*02:01	*13:01	+	-	-	+	-	-	-	-	-	-	-	+	-	-	+	-	
23	9124	IHL	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
24	9035	JBUSH	*04:01		-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	
25	9049	IBW9	*01:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
26	9285	WT49	*04:01		-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	
27	9191	CH1007	*03:01	*13:01	+	+	+	+	+	-	-	-	-	-	-	+	-	-	+	-	
28	9320	BEL5GB	*02:01	*04:01	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	
29	9050	MOU	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
30	9021	RSH	*01:01	*105:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
31	9019	DUCAF	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
32	9297	HAG	*02:01	*03:01	+	+	+	+	+	-	-	-	-	-	-	-	-	-	-	-	
33	9098	MT14B	*04:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
34	9104	DHIF	*04:01		-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	
35	9302	SSTO	*04:01		-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	
36	9024	KT17	*05:01		-	-	-	-	-	-	+	+	+	-	-	-	-	-	+	-	
37	9065	HHKB	*04:01		-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	
38	9099	LZL	*04:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
39	9315	CML	*04:01	*14:01	+	+	+	+	+	+	-	-	-	-	-	-	-	+	-	-	
40	9134	WHONP199	*05:01		-	-	-	-	-	-	+	+	+	-	-	-	-	-	+	-	
41	9055	H0301	*05:01		-	-	-	-	-	-	+	+	+	-	-	-	-	-	+	-	
42	9066	TAB089	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
43	9076	T7526	*13:01		+	-	-	+	-	-	-	-	-	-	-	+	-	-	+	-	
44	9057	TEM	*04:01		-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	
45	9239	SHJO	*03:01	*04:02	+	+	+	+	+	-	-	-	-	-	-	-	-	-	-	-	
46	9013	SCHU	*04:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
47	9045	TUBO	*02:01	*04:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
48	9303	TER-ND	*02:01	*04:01	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot-specific information

CELL LINE VALIDATION SHEET																						
DPB1 SSP kit																						
					Well																	
						33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	
				Lot No.:		201294333	201184934	201184935	201184936	201294337	201184938	201184939	201184940	201294341	201184942	201184943	201294344	201184945	201184946	201184947	201184948	
	IHWC cell line		DPB1																			
1	9001	SA	*04:02		-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	
2	9280	LK707	*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
3	9011	E4181324	*02:01	*04:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
4	9275	GU373	*04:02	*17:01	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	+	-	
5	9009	KAS011	*04:01	*14:01	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
6	9353	SM	*05:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	
7	9020	QBL	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
8	9025	DEU	*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
9	9026	YAR	*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
10	9107	LKT3	*05:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	
11	9051	PITOUT	*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
12	9052	DBB	*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
13	9004	JESTHOM	*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
14	9071	OLGA	*03:01	*04:02	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	
15	9075	DKB	*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
16	9037	SWEIG007	*04:02		-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	
17	9282	CTM3953540	*01:01	*02:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
18	9257	32367	*04:02	*18:01	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	+	-	
19	9038	BM16	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
20	9059	SLE005	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
21	9064	AMALA	*04:02		-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	
22	9056	KOSE	*02:01	*13:01	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
23	9124	IHL	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
24	9035	JBUSH	*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
25	9049	IBW9	*01:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
26	9285	WT49	*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
27	9191	CH1007	*03:01	*13:01	+	-	-	-	-	-	-	-	+	-	-	+	-	-	-	-	-	
28	9320	BEL5GB	*02:01	*04:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
29	9050	MOU	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
30	9021	RSH	*01:01	*105:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
31	9019	DUCAF	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
32	9297	HAG	*02:01	*03:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
33	9098	MT14B	*04:02		-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	
34	9104	DHIF	*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
35	9302	SSTO	*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
36	9024	KT17	*05:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	
37	9065	HHKB	*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
38	9099	LZL	*04:02		-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	
39	9315	CML	*04:01	*14:01	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
40	9134	WHONP199	*05:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	
41	9055	H0301	*05:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	
42	9066	TAB089	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
43	9076	T7526	*13:01		+	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	
44	9057	TEM	*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
45	9239	SHJO	*03:01	*04:02	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	+	-	
46	9013	SCHU	*04:02		-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	
47	9045	TUBO	*02:01	*04:02	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	+	-	
48	9303	TER-ND	*02:01	*04:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	

101.311-24/06 – including *Taq* pol., IFU-01
 101.311-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **09N**

Lot-specific information

CERTIFICATE OF ANALYSIS

Olerup SSP® DPB1 SSP

Product number: 101.311-24/04 – including *Taq* pol.
 101.311-24u/04u – without *Taq* pol.
 Lot number: 09N
 Expiry date: 2014-August-01
 Number of tests: 24 tests – Product No. 101.311-24/24u
 4 tests – Product No. 101.311-04/04u
 Number of wells per test: 48

Well specifications:

Well No.	Production No.	Well No.	Production No.	Well No.	Production No.
1	2011-849-01	17	2011-849-17	33	2012-943-33
2	2011-849-02	18	2011-849-18	34	2011-849-34
3	2012-943-03	19	2011-849-19	35	2011-849-35
4	2012-943-04	20	2011-849-20	36	2011-849-36
5	2011-849-05	21	2011-849-21	37	2012-943-37
6	2011-849-06	22	2012-943-22	38	2011-849-38
7	2012-943-07	23	2011-849-23	39	2011-849-39
8	2011-849-08	24	2011-849-24	40	2011-849-40
9	2011-849-09	25	2011-849-25	41	2012-943-41
10	2011-849-10	26	2011-849-26	42	2011-849-42
11	2011-849-11	27	2011-849-27	43	2011-849-43
12	2011-849-12	28	2011-849-28	44	2012-943-44
13	2011-849-13	29	2011-849-29	45	2011-849-45
14	2011-849-14	30	2011-849-30	46	2011-849-46
15	2011-849-15	31	2012-943-31	47	2011-849-47
16	2011-849-16	32	2011-849-32	48	2011-849-48

The specificity of each primer solution of the kit has been tested against 48 well characterized IHWC cell line DNAs.

No DNAs carrying the alleles to be amplified by primer solutions 41, 42, 44 and 46 were available. In primer solutions 41, 42 and 46 it was only possible to test the 5'-primers, the 3'-primers were not possible to test. In primer solution 44 it was only possible to test the 3'-primers, the 5'-primers were not possible to test. In primer solutions 15, 16, 23, 25, 26, 32, 33 and 38 one or two 5'-primers were not possible to test. In primer solutions 8, 11, 24, 28 and 48 one 3'-primer was not possible to test. Additional primers in primer solutions 14, 38, 40, 43, 47 and 48 were tested by separately adding 5'-primers or 3'-primers.

101.311-24/06 – including *Taq* pol., IFU-01
101.311-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: **09N**

Lot-specific information

Results: No false positive or false negative amplifications were obtained.

Date of approval: 2012-February-15

Approved by:

Production Quality Control

101.311-24/06 – including *Taq* pol., IFU-01
101.311-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: **09N**

Lot-specific information

Declaration of Conformity

Product name: Olerup SSP® DPB1
Product number: 101.311-24/04, -24u/04u
Lot number: 09N

Intended use: DPB1 high resolution histocompatibility testing

Manufacturer: Olerup SSP AB
Franzengatan 5
SE-112 51 Stockholm, Sweden
Phone: +46-8-717 88 27
Fax: +46-8-717 88 18

We, Olerup SSP AB, hereby declare that this product, to which this Declaration of Conformity relates is in conformity with the following Standard(s) and other normative document(s) ISO 9001:2008 and ISO 13485:2003, following the provisions of the 98/79/EC Directive on *in vitro* diagnostic medical devices, Annex III, as transposed into the national laws of the Member States of the European Union.

The Technical Documentation File is maintained at Olerup SSP AB, Franzengatan 5, SE-112 51 Stockholm, Sweden.

The Authorized Representative located within the Community is: Olerup SSP AB.

Stockholm, Sweden
2012-February-15

Ann-Cathrin Jareman
Head of QA and Regulatory Affairs

101.311-24/06 – including *Taq* pol., IFU-01
101.311-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: **09N**

Lot-specific information

101.311-24/06 – including *Taq* pol., IFU-01
101.311-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: **09N**

Lot-specific information

ADDRESSES:

Manufacturer:

Olerup SSP AB, Franzengatan 5, SE-112 51 Stockholm, Sweden.

Tel: +46-8-717 88 27

Fax: +46-8-717 88 18

E-mail: info-ssp@olerup.com

Web page: <http://www.olerup-ssp.com>

Distributed by:

Olerup GmbH, Löwengasse 47 / 6, AT-1030 Vienna, Austria.

Tel: +43-1-710 15 00

Fax: +43-1-710 15 00 10

E-mail: support-at@olerup.com

Web page: <http://www.olerup.com>

Olerup Inc., 901 S. Bolmar St., Suite R, West Chester, PA 19382

Tel: 1-877-OLERUP1

Fax: 610-344-7989

E-mail: info.us@olerup.com

Web page: <http://www.olerup.com>

For information on *Olerup* SSP distributors worldwide, contact **Olerup GmbH**.