

101.125-24/06 – including *Taq* pol., IFU-01.  
101.125-24u/06u – without *Taq* pol., IFU-02.

Visit <https://labproducts.caredx.com> for  
“Instructions for Use” (IFU)

Lot No.: **7G1**

Lot-specific information

## **Olerup SSP® DRB1\*15**

|                                         |                                                                                        |
|-----------------------------------------|----------------------------------------------------------------------------------------|
| <b>Product number:</b>                  | 101.125-24/06 – including <i>Taq</i> pol.<br>101.125-24u/06u – without <i>Taq</i> pol. |
| <b>Lot number:</b>                      | 7G1                                                                                    |
| <b>Expiry date:</b>                     | 2022-09-01                                                                             |
| <b>Number of tests:</b>                 | 24 test – Product No. 101.125-24/24u<br>6 tests – Product No. 101.125-06/06u           |
| <b>Number of wells per test:</b>        | 31+1                                                                                   |
| <b>Storage - pre-aliquoted primers:</b> | dark at -20°C                                                                          |
| - PCR Master Mix:                       | -20°C                                                                                  |
| - Adhesive PCR seals                    | RT                                                                                     |
| - Product Insert                        | RT                                                                                     |

### **This Product Description is only valid for Lot No. 7G1**

Complete product documentation consists of generic Instructions for Use (IFU), lot specific Product Insert, Worksheet and Certificate.

### **CHANGES COMPARED TO THE PREVIOUS OLERUP SSP® DRB1\*15 LOT (2F2)**

The DRB1\*15 kit is updated for new alleles to enable separation of:

- Null and Alternatively expressed alleles
- The product documentation has been updated for new alleles of IMGT 3.32.0

The format of the Worksheet has been changed.

The DRB1\*15 specificity and interpretation tables have been updated for the DRB1 alleles described since the previous *Olerup SSP®* DRB1\*15 lot was made (**Lot No. 2F2**).

101.125-24/06 – including *Taq* pol., IFU-01.  
101.125-24u/06u – without *Taq* pol., IFU-02.

Visit <https://labproducts.caredx.com> for  
“Instructions for Use” (IFU)

Lot No.: **7G1**

**Lot-specific information**

The primers of the wells detailed below have been exchanged, added or modified compared to the previous lot.

| Well | 5'-primer             | 3'-primer | rationale                                                                              |
|------|-----------------------|-----------|----------------------------------------------------------------------------------------|
| 3    | Added                 | -         | 5'-primer added for the DRB1*15:01:33 allele.                                          |
| 16   | Modified              | -         | 5'-primer modified for improved HLA-specific amplification.                            |
| 22   | Removed,<br>Exchanged | -         | Excess 5'-primer removed, 5'-primer exchanged for improved HLA-specific amplification. |
| 28   | Modified              | -         | 5'-primer modified for improved HLA-specific amplification.                            |

Change in revision R01 compared to R00:

1. Primer mix 7 does not amplify the DRB1\*15:56 allele. Thus, this lot of the DRB1\*15 subtyping kit cannot distinguish the DRB1\*15:56 and DRB1\*15:01:01:01-15:01:35, 15:51-15:53, 15:62, 15:71-15:72, 15:76, 15:79, 15:81-15:83, 15:86-15:87, 15:90, 15:92, 15:95, 15:97-15:98, 15:106-15:109, 15:111, 15:114, 15:116, 15:121, 15:124, 15:127-15:128, 15:132-15:133, 15:135, 15:139, 15:141, 15:144-15:146 and 15:150 alleles. This has been corrected in the Specificity and Interpretation Tables.

Changes in revision R02 compared to R01:

1. The expiration date has been altered due to extension of shelf-life.

101.125-24/06 – including *Taq* pol., IFU-01.  
101.125-24u/06u – without *Taq* pol., IFU-02.

Visit <https://labproducts.caredx.com> for  
“Instructions for Use” (IFU)

Lot No.: **7G1**

**Lot-specific information**

Well **32** contains Negative Control primer pairs, that will amplify more than 95% of the Olerup SSP<sup>®</sup> HLA Class I, DRB, DQB1, DPB1 and DQA1 amplicons as well as all the amplicons generated by the control primer pairs matching the human growth hormone gene.

HLA-specific PCR product sizes range from 75 to 200 base pairs.  
The PCR product generated by the positive control primer pair is 430 base pairs.

| Length of PCR product        | 105                  | 200                     | 105                  | 80                   | 75                   | 80                   | 85                   |
|------------------------------|----------------------|-------------------------|----------------------|----------------------|----------------------|----------------------|----------------------|
| <b>5'-primer<sup>1</sup></b> | <b>164</b>           | <b>340</b>              | <b>440</b>           | <b>45</b>            | <b>45</b>            | <b>43</b>            | <b>36</b>            |
|                              | 5'-CAC <sup>3'</sup> | 5'-Agg <sup>3'</sup>    | 5'-TTA <sup>3'</sup> | 5'-Tgg <sup>3'</sup> | 5'-Tgg <sup>3'</sup> | 5'-Tgg <sup>3'</sup> | 5'-TAC <sup>3'</sup> |
|                              |                      |                         |                      |                      |                      |                      | <b>36</b>            |
|                              |                      |                         |                      |                      |                      |                      | 5'-TAT <sup>3'</sup> |
| <b>3'-primer<sup>2</sup></b> | <b>231</b>           | <b>2<sup>nd</sup> I</b> | <b>507</b>           | <b>59</b>            | <b>58</b>            | <b>57</b>            | <b>47</b>            |
|                              | 5'-TgC <sup>3'</sup> | 5'-AAA <sup>3'</sup>    | 5'-TTg <sup>3'</sup> | 5'-CTC <sup>3'</sup> | 5'-ggC <sup>3'</sup> | 5'-CTC <sup>3'</sup> | 5'-ACA <sup>3'</sup> |
|                              |                      |                         |                      |                      |                      |                      | <b>48</b>            |
|                              |                      |                         |                      |                      |                      |                      | 5'-gCA <sup>3'</sup> |
|                              |                      |                         |                      |                      |                      |                      | <b>48</b>            |
|                              |                      |                         |                      |                      |                      |                      | 5'-gCC <sup>3'</sup> |
|                              |                      |                         |                      |                      |                      |                      | <b>52</b>            |
|                              |                      |                         |                      |                      |                      |                      | 5'-TgT <sup>3'</sup> |
| <b>A*</b>                    | +                    | +                       | +                    |                      |                      |                      |                      |
| <b>B*</b>                    | +                    | +                       | +                    |                      |                      |                      |                      |
| <b>C*</b>                    | +                    | +                       | +                    |                      |                      |                      |                      |
| <b>DRB1</b>                  |                      |                         |                      | +                    | +                    |                      |                      |
| <b>DRB3</b>                  |                      |                         |                      | +                    | +                    |                      |                      |
| <b>DRB5</b>                  |                      |                         |                      | +                    |                      |                      |                      |
| <b>DQB1</b>                  |                      |                         |                      |                      | +                    |                      |                      |
| <b>DPB1</b>                  |                      |                         |                      |                      |                      | +                    |                      |
| <b>DQA1</b>                  |                      |                         |                      |                      |                      |                      | +                    |

<sup>1</sup>The nucleotide position for HLA class I genes and the codon for HLA class II genes, in the 2<sup>nd</sup> or 3<sup>rd</sup> exon, matching the specificity-determining 3'-end of the primer is given. Nucleotide and codon numbering as on the [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla) web site. The sequence of the 3 terminal nucleotides of the primer is given.

<sup>2</sup>The nucleotide position for HLA class I genes and the codon for HLA class II genes, in the 2<sup>nd</sup> or 3<sup>rd</sup> exon or the 2<sup>nd</sup> intron, matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Nucleotide and codon numbering as on the [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla) web site. The sequence of the 3 terminal nucleotides of the primer is given.

101.125-24/06 – including *Taq* pol., IFU-01.  
101.125-24u/06u – without *Taq* pol., IFU-02.

Visit <https://labproducts.caredx.com> for  
“Instructions for Use” (IFU)

Lot No.: **7G1**

Lot-specific information

## PRODUCT DESCRIPTION

### DRB1\*15 SSP subtyping

#### CONTENT

The primer set contains 5'- and 3'-primers for identifying the DRB1\*15:01 to DRB1\*15:150 alleles.

#### PLATE LAYOUT

Each test consists of 32 PCR reactions in a 32 well cut PCR plate.

|           |           |           |           |           |           |           |           |
|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
| <b>1</b>  | <b>2</b>  | <b>3</b>  | <b>4</b>  | <b>5</b>  | <b>6</b>  | <b>7</b>  | <b>8</b>  |
| <b>9</b>  | <b>10</b> | <b>11</b> | <b>12</b> | <b>13</b> | <b>14</b> | <b>15</b> | <b>16</b> |
| <b>17</b> | <b>18</b> | <b>19</b> | <b>20</b> | <b>21</b> | <b>22</b> | <b>23</b> | <b>24</b> |
| <b>25</b> | <b>26</b> | <b>27</b> | <b>28</b> | <b>29</b> | <b>30</b> | <b>31</b> | <b>NC</b> |

The 32 well cut PCR plate is marked with 'DRB1\*15' in silver/gray ink.

Well No. 1 is marked with the Lot No. '7G1'.

Wells 1 to 31 – DRB1\*15 high resolution primers.

Well 32 – Negative Control (NC).

A faint row of numbers is seen between wells 1 and 2 or wells 7 and 8 of the PCR trays. These stem from the manufacture of the trays, and should be disregarded.

The PCR plates are covered with a PCR-compatible foil.

**Please note:** When removing each 32 well PCR plate, make sure that the remaining plates stay covered. Use a scalpel or a similar instrument to carefully cut the foil between the plates.

#### INTERPRETATION

Due to the sharing of sequence motifs between DRB1 alleles, non-DRB1\*15 alleles will be amplified by other primer mixes. For further details see Specificity Table.

#### UNIQUELY IDENTIFIED ALLELES

All the DRB1\*15 alleles, i.e. **DRB1\*15:01 to DRB1\*15:150**, recognized by the HLA Nomenclature Committee in April 2018<sup>1, 2</sup> will be amplified by the primers in the DRB1\*15 subtyping kit<sup>3, 4</sup>.

The DRB1\*15 kit enables separation of the confirmed DRB1\*15 alleles as listed in the IMGT/HLA database 3.26.0. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources. Current allele confirmation status for DRB1\*15 alleles is listed below.

The DRB1\*15 kit also enables identification of null and alternatively expressed alleles.

101.125-24/06 – including *Taq* pol., IFU-01.  
101.125-24u/06u – without *Taq* pol., IFU-02.

Visit <https://labproducts.caredx.com> for  
“Instructions for Use” (IFU)

Lot No.: **7G1**

**Lot-specific information**

The following DRB1\*15 alleles can be distinguished by the different sizes of the HLA-specific PCR product:

| Alleles                                      | Primer mix |
|----------------------------------------------|------------|
| DRB1*15:06:01-15:06:03, 15:77,15:110, 15:125 | 9          |
| DRB1*15:16, 15:35                            | 28         |
| DRB1*15:30, 15:104:01-15:104:02              | 19         |
| DRB1*15:42, 15:64                            | 27         |
| DRB1*15:75, 15:134N                          | 26         |

<sup>1</sup>DRB1 alleles listed on the IMGT/HLA web page 2018-April-16, release 3.32.0, [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla).

<sup>2</sup>Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

<sup>3</sup>The DRB1\*15:70, 15:89 and 15:91 and the DRB1\*16:33 and 16:36 alleles give rise to identical amplification patterns with the DRB1\*15 high resolution kit. These alleles can be distinguished by the DR low resolution and/or DRB1\*16 kits.

<sup>4</sup>This lot of the DRB1\*15 subtyping kit cannot distinguish the DRB1\*15:56 and DRB1\*15:01:01:01-15:01:35, 15:51-15:53, 15:62, 15:71-15:72, 15:76, 15:79, 15:81-15:83, 15:86-15:87, 15:90, 15:92, 15:95, 15:97-15:98, 15:106-15:109, 15:111, 15:114, 15:116, 15:121, 15:124, 15:127-15:128, 15:132-15:133, 15:135, 15:139, 15:141, 15:144-15:146 and 15:150 alleles.

101.125-24/06 – including *Taq* pol., IFU-01.  
101.125-24u/06u – without *Taq* pol., IFU-02.

Visit <https://labproducts.caredx.com> for  
“Instructions for Use” (IFU)

Lot No.: **7G1**

Lot-specific information

## ALLELE CONFIRMATION STATUS

| Allele                  | Status <sup>1</sup> | Allele                  | Status <sup>1</sup> | Allele            | Status <sup>1</sup> | Allele              | Status <sup>1</sup> |
|-------------------------|---------------------|-------------------------|---------------------|-------------------|---------------------|---------------------|---------------------|
| DRB1*15:01:01:01        | Confirmed           | DRB1*15:02:16           | Unconfirmed         | DRB1*15:42        | Unconfirmed         | DRB1*15:91          | Unconfirmed         |
| DRB1*15:01:01:02        | Unconfirmed         | DRB1*15:02:17           | Unconfirmed         | DRB1*15:43        | Unconfirmed         | DRB1*15:92          | Unconfirmed         |
| DRB1*15:01:01:03        | Unconfirmed         | DRB1*15:02:18           | Unconfirmed         | DRB1*15:44        | Unconfirmed         | <b>DRB1*15:93</b>   | <b>Confirmed</b>    |
| DRB1*15:01:01:04        | Unconfirmed         | <b>DRB1*15:03:01:01</b> | <b>Confirmed</b>    | DRB1*15:45        | Unconfirmed         | <b>DRB1*15:94</b>   | <b>Confirmed</b>    |
| <b>DRB1*15:01:02</b>    | <b>Confirmed</b>    | <b>DRB1*15:03:01:02</b> | <b>Confirmed</b>    | DRB1*15:46        | Unconfirmed         | DRB1*15:95          | Unconfirmed         |
| <b>DRB1*15:01:03</b>    | <b>Confirmed</b>    | <b>DRB1*15:03:02</b>    | <b>Confirmed</b>    | DRB1*15:47        | Unconfirmed         | DRB1*15:96          | Unconfirmed         |
| <b>DRB1*15:01:04</b>    | <b>Confirmed</b>    | <b>DRB1*15:04</b>       | <b>Confirmed</b>    | DRB1*15:48        | Unconfirmed         | DRB1*15:97          | Unconfirmed         |
| DRB1*15:01:05           | Unconfirmed         | <b>DRB1*15:05</b>       | <b>Confirmed</b>    | DRB1*15:49        | Unconfirmed         | DRB1*15:98          | Unconfirmed         |
| <b>DRB1*15:01:06</b>    | <b>Confirmed</b>    | <b>DRB1*15:06:01</b>    | <b>Confirmed</b>    | DRB1*15:50N       | Unconfirmed         | DRB1*15:99          | Unconfirmed         |
| <b>DRB1*15:01:07</b>    | <b>Confirmed</b>    | DRB1*15:06:02           | Unconfirmed         | DRB1*15:51        | Unconfirmed         | DRB1*15:100         | Unconfirmed         |
| <b>DRB1*15:01:08</b>    | <b>Confirmed</b>    | DRB1*15:06:03           | Unconfirmed         | DRB1*15:52        | Unconfirmed         | DRB1*15:101         | Unconfirmed         |
| DRB1*15:01:09           | Unconfirmed         | <b>DRB1*15:07:01</b>    | <b>Confirmed</b>    | DRB1*15:53        | Unconfirmed         | <b>DRB1*15:102</b>  | <b>Confirmed</b>    |
| DRB1*15:01:10           | Unconfirmed         | DRB1*15:07:02           | Unconfirmed         | DRB1*15:54        | Unconfirmed         | DRB1*15:103         | Unconfirmed         |
| <b>DRB1*15:01:11</b>    | <b>Confirmed</b>    | DRB1*15:07:03           | Unconfirmed         | <b>DRB1*15:55</b> | <b>Confirmed</b>    | DRB1*15:104         | Unconfirmed         |
| <b>DRB1*15:01:12</b>    | <b>Confirmed</b>    | DRB1*15:08              | Unconfirmed         | <b>DRB1*15:56</b> | <b>Confirmed</b>    | DRB1*15:105         | Unconfirmed         |
| <b>DRB1*15:01:13</b>    | <b>Confirmed</b>    | <b>DRB1*15:09</b>       | <b>Confirmed</b>    | <b>DRB1*15:57</b> | <b>Confirmed</b>    | DRB1*15:106         | Unconfirmed         |
| DRB1*15:01:14           | Unconfirmed         | <b>DRB1*15:10</b>       | <b>Confirmed</b>    | DRB1*15:58        | Unconfirmed         | DRB1*15:107         | Unconfirmed         |
| DRB1*15:01:15           | Unconfirmed         | <b>DRB1*15:11</b>       | <b>Confirmed</b>    | DRB1*15:59        | Unconfirmed         | DRB1*15:108         | Unconfirmed         |
| DRB1*15:01:16           | Unconfirmed         | DRB1*15:12              | Unconfirmed         | DRB1*15:60        | Unconfirmed         | DRB1*15:109         | Unconfirmed         |
| DRB1*15:01:17           | Unconfirmed         | DRB1*15:13              | Unconfirmed         | <b>DRB1*15:61</b> | <b>Confirmed</b>    | DRB1*15:110         | Unconfirmed         |
| DRB1*15:01:18           | Unconfirmed         | <b>DRB1*15:14</b>       | <b>Confirmed</b>    | DRB1*15:62        | Unconfirmed         | DRB1*15:111         | Unconfirmed         |
| DRB1*15:01:19           | Unconfirmed         | <b>DRB1*15:15:01</b>    | <b>Confirmed</b>    | DRB1*15:63        | Unconfirmed         | <b>DRB1*15:112</b>  | <b>Confirmed</b>    |
| DRB1*15:01:20           | Unconfirmed         | DRB1*15:15:02           | Unconfirmed         | <b>DRB1*15:64</b> | <b>Confirmed</b>    | <b>DRB1*15:113N</b> | <b>Confirmed</b>    |
| DRB1*15:01:21           | Unconfirmed         | DRB1*15:16              | Unconfirmed         | DRB1*15:65        | Unconfirmed         | DRB1*15:114         | Unconfirmed         |
| DRB1*15:01:22           | Unconfirmed         | DRB1*15:17N             | Unconfirmed         | DRB1*15:66:01     | Unconfirmed         | DRB1*15:115N        | Unconfirmed         |
| DRB1*15:01:23           | Unconfirmed         | <b>DRB1*15:18</b>       | <b>Confirmed</b>    | DRB1*15:66:02     | Unconfirmed         | DRB1*15:116         | Unconfirmed         |
| DRB1*15:01:24           | Unconfirmed         | DRB1*15:19              | Unconfirmed         | <b>DRB1*15:67</b> | <b>Confirmed</b>    | DRB1*15:117         | Unconfirmed         |
| <b>DRB1*15:01:25</b>    | <b>Confirmed</b>    | <b>DRB1*15:20</b>       | <b>Confirmed</b>    | <b>DRB1*15:68</b> | <b>Confirmed</b>    | <b>DRB1*15:118</b>  | <b>Confirmed</b>    |
| DRB1*15:01:26           | Unconfirmed         | DRB1*15:21              | Unconfirmed         | DRB1*15:69        | Unconfirmed         | DRB1*15:119         | Unconfirmed         |
| DRB1*15:01:27           | Unconfirmed         | <b>DRB1*15:22</b>       | <b>Confirmed</b>    | DRB1*15:70        | Unconfirmed         | DRB1*15:120         | Unconfirmed         |
| DRB1*15:01:28           | Unconfirmed         | <b>DRB1*15:23</b>       | <b>Confirmed</b>    | DRB1*15:71        | Unconfirmed         | DRB1*15:121         | Unconfirmed         |
| DRB1*15:01:29           | Unconfirmed         | <b>DRB1*15:24</b>       | <b>Confirmed</b>    | DRB1*15:72        | Unconfirmed         | DRB1*15:122         | Unconfirmed         |
| DRB1*15:01:30           | Unconfirmed         | DRB1*15:25              | Unconfirmed         | DRB1*15:73        | Unconfirmed         | <b>DRB1*15:123</b>  | <b>Confirmed</b>    |
| DRB1*15:01:31           | Unconfirmed         | DRB1*15:26              | Unconfirmed         | <b>DRB1*15:74</b> | <b>Confirmed</b>    | DRB1*15:124         | Unconfirmed         |
| <b>DRB1*15:02:01:01</b> | <b>Confirmed</b>    | <b>DRB1*15:27</b>       | <b>Confirmed</b>    | <b>DRB1*15:75</b> | <b>Confirmed</b>    | DRB1*15:125         | Unconfirmed         |
| DRB1*15:02:01:02        | Unconfirmed         | DRB1*15:28              | Unconfirmed         | DRB1*15:76        | Unconfirmed         | DRB1*15:126         | Unconfirmed         |
| <b>DRB1*15:02:02</b>    | <b>Confirmed</b>    | DRB1*15:29              | Unconfirmed         | DRB1*15:77        | Unconfirmed         | DRB1*15:127         | Unconfirmed         |
| DRB1*15:02:03           | Unconfirmed         | DRB1*15:30              | Unconfirmed         | <b>DRB1*15:78</b> | <b>Confirmed</b>    | DRB1*15:128         | Unconfirmed         |
| <b>DRB1*15:02:04</b>    | <b>Confirmed</b>    | DRB1*15:31              | Unconfirmed         | DRB1*15:79        | Unconfirmed         | DRB1*15:129N        | Unconfirmed         |
| DRB1*15:02:05           | Unconfirmed         | DRB1*15:32              | Unconfirmed         | DRB1*15:80N       | Unconfirmed         | DRB1*15:130         | Unconfirmed         |
| DRB1*15:02:06           | Unconfirmed         | <b>DRB1*15:33</b>       | <b>Confirmed</b>    | DRB1*15:81        | Unconfirmed         | DRB1*15:131         | Unconfirmed         |
| DRB1*15:02:07           | Unconfirmed         | DRB1*15:34              | Unconfirmed         | DRB1*15:82        | Unconfirmed         | DRB1*15:132         | Unconfirmed         |
| DRB1*15:02:08           | Unconfirmed         | <b>DRB1*15:35</b>       | <b>Confirmed</b>    | DRB1*15:83        | Unconfirmed         | DRB1*15:133         | Unconfirmed         |
| DRB1*15:02:09           | Unconfirmed         | <b>DRB1*15:36</b>       | <b>Confirmed</b>    | DRB1*15:84        | Unconfirmed         | DRB1*15:134N        | Unconfirmed         |
| DRB1*15:02:10           | Unconfirmed         | <b>DRB1*15:37:01</b>    | <b>Confirmed</b>    | DRB1*15:85        | Unconfirmed         | DRB1*15:135         | Unconfirmed         |
| DRB1*15:02:11           | Unconfirmed         | <b>DRB1*15:37:02</b>    | <b>Confirmed</b>    | DRB1*15:86        | Unconfirmed         | <b>DRB1*15:136</b>  | <b>Confirmed</b>    |
| DRB1*15:02:12           | Unconfirmed         | <b>DRB1*15:38</b>       | <b>Confirmed</b>    | DRB1*15:87        | Unconfirmed         | DRB1*15:137N        | Unconfirmed         |
| DRB1*15:02:13           | Unconfirmed         | DRB1*15:39              | Unconfirmed         | DRB1*15:88        | Unconfirmed         | DRB1*15:138N        | Unconfirmed         |
| DRB1*15:02:14           | Unconfirmed         | <b>DRB1*15:40</b>       | <b>Confirmed</b>    | DRB1*15:89        | Unconfirmed         | DRB1*15:139         | Unconfirmed         |
| <b>DRB1*15:02:15</b>    | <b>Confirmed</b>    | DRB1*15:41              | Unconfirmed         | DRB1*15:90        | Unconfirmed         |                     |                     |

<sup>1</sup>Allele status “confirmed” or “unconfirmed” as listed on the IMGT/HLA web page 2016-October-14, release 3.26.0, [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla).

## RESOLUTION IN HOMO- AND HETEROZYGOES

Results file with resolution in DRB1\*15 homo- and heterozygotes is available upon request.

101.125-24/06 – including *Taq* pol., IFU-01.  
101.125-24u/06u – without *Taq* pol., IFU-02.

Visit <https://labproducts.caredx.com> for  
“Instructions for Use” (IFU)

Lot No.: **7G1**

Lot-specific information

## SPECIFICITY TABLE

### DRB1\*15 SSP subtyping

Specificities and sizes of the PCR products of the 31+1 primer mixes used for  
DRB1\*15 SSP subtyping

| Primer Mix | Size of spec. PCR product <sup>1</sup> | Size of control band <sup>2</sup> | Amplified DRB1*15 alleles <sup>3</sup>                                                                                                                                                                                                                                                                                                                                                                                                                           | Other amplified DRB1 alleles                                                                                          |
|------------|----------------------------------------|-----------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------|
| 1          | 140 bp                                 | 515 bp                            | *15:01:01:01-15:02:09, 15:02:11-15:06:03, 15:08-15:10, 15:12-15:27, 15:29-15:33, 15:35-15:47, 15:49-15:58, 15:60-15:68, 15:70-15:87, 15:89-15:95, 15:97-15:114, 15:116-15:129N, 15:131-15:142, 15:144-15:147, 15:149-15:150                                                                                                                                                                                                                                      | *16:09:01-16:10, 16:33, 16:36-16:37                                                                                   |
| 2          | 200 bp                                 | 515 bp                            | *15:01:01:01-15:03:02, 15:06:01-15:09, 15:11-15:14, 15:16, 15:18 <sup>w</sup> , 15:19-15:20, 15:22, 15:24, 15:26, 15:29-15:30, 15:32-15:33, 15:35-15:48, 15:51-15:53, 15:55-15:62, 15:64-15:65, 15:67-15:68, 15:71-15:72, 15:74-15:79, 15:81-15:83, 15:85-15:87, 15:90, 15:92-15:99, 15:101-15:102, 15:104:01-15:114, 15:116, 15:118-15:119, 15:121, 15:123-15:128, 15:130-15:141, 15:143-15:146, 15:149-15:150                                                  |                                                                                                                       |
| 3          | 150 bp                                 | 430 bp                            | *15:01:01:01-15:02:18, 15:06:01-15:09, 15:11-15:14, 15:16, 15:18 <sup>w</sup> , 15:19-15:20, 15:22, 15:24, 15:26, 15:28-15:30, 15:32-15:33, 15:35-15:48, 15:51-15:53, 15:55-15:56, 15:59-15:62, 15:64-15:65, 15:67-15:69, 15:71-15:72, 15:74-15:77, 15:79, 15:81-15:83, 15:85-15:87, 15:90, 15:92-15:93, 15:95-15:99, 15:101-15:102, 15:104:01-15:114, 15:116, 15:118-15:119, 15:121, 15:123-15:128, 15:131-15:136, 15:138N-15:141, 15:143-15:146, 15:149-15:150 | *04:73:01-04:73:02, 04:105:01-04:105:02, 11:82, 13:09, 13:129, 14:37                                                  |
| 4          | 260 bp                                 | 430 bp                            | *15:01:01:01-15:01:35, 15:03:01:01-15:07:03, 15:09-15:10, 15:12-15:13, 15:16-15:18, 15:20-15:25, 15:32-15:33, 15:35-15:37:02, 15:40-15:43, 15:45-15:46, 15:48-15:57, 15:59, 15:61-15:62, 15:64-15:67, 15:70-15:73, 15:75-15:77, 15:79, 15:81-15:84, 15:86-15:98, 15:100, 15:102, 15:106-15:114, 15:116-15:117, 15:120-15:121, 15:123-15:125, 15:127-15:130, 15:132-15:139, 15:141-15:147, 15:150                                                                 | *16:15, 16:23, 16:33-16:34, 16:36, 16:38:01-16:38:02                                                                  |
| 5          | 260 bp                                 | 430 bp                            | *15:02:01:01-15:02:18, 15:08, 15:11, 15:14-15:15:03, 15:19, 15:26-15:27, 15:29-15:31, 15:34, 15:38-15:39, 15:44, 15:47, 15:58, 15:60, 15:63, 15:68, 15:78, 15:80N, 15:99, 15:101, 15:103-15:105:02, 15:115N, 15:118-15:119, 15:122, 15:126, 15:131, 15:140, 15:149                                                                                                                                                                                               | *16:01:01-16:01:06, 16:01:08-16:05:02, 16:07-16:11, 16:13N-16:14, 16:16-16:22, 16:24-16:32, 16:35, 16:37, 16:39-16:47 |

101.125-24/06 – including *Taq* pol., IFU-01.  
101.125-24u/06u – without *Taq* pol., IFU-02.

Visit <https://labproducts.caredx.com> for  
“Instructions for Use” (IFU)

Lot No.: **7G1**

Lot-specific information

|                         |        |               |                                                                                                         |                                                                                                               |
|-------------------------|--------|---------------|---------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------|
| <b>6</b>                | 150 bp | 430 bp        | *15:03:01:01-15:03:02, 15:57, 15:78, 15:94,<br>15:130, 15:137N                                          |                                                                                                               |
|                         | 180 bp |               | *15:14                                                                                                  |                                                                                                               |
| <b>7</b>                | 135 bp | 430 bp        | *15:07:01-15:07:03, 15:11, 15:34, 15:88,<br>15:96, 15:115N                                              | *16:01:01-16:05:02,<br>16:07-16:08, 16:11,<br>16:13N-16:32,<br>16:34-16:35,<br>16:38:01-16:43,<br>16:45-16:47 |
| <b>8</b>                | 200 bp | 430 bp        | *15:05, 15:23, 15:25, 15:27, 15:31, 15:34,<br>15:54, 15:63, 15:66:01-15:66:02, 15:73,<br>15:100, 15:120 |                                                                                                               |
|                         | 245 bp |               | *15:36                                                                                                  |                                                                                                               |
| <b>9<sup>4</sup></b>    | 90 bp  | 430 bp        | *15:06:01-15:06:03, 15:77, 15:94, 15:125                                                                |                                                                                                               |
|                         | 165 bp |               | *15:138N                                                                                                |                                                                                                               |
|                         | 215 bp |               | *15:19, 15:69, 15:110, 15:113N                                                                          |                                                                                                               |
| <b>10</b>               | 170 bp | 430 bp        | *15:12                                                                                                  |                                                                                                               |
|                         | 200 bp |               | *15:04, 15:15:01-15:15:03, 15:88, 15:142,<br>15:147                                                     |                                                                                                               |
| <b>11</b>               | 155 bp | 430 bp        | *15:47, 15:65                                                                                           |                                                                                                               |
|                         | 200 bp |               | *15:49                                                                                                  |                                                                                                               |
| <b>12<sup>4</sup></b>   | 95 bp  | 430 bp        | *15:09, 15:48                                                                                           | *04:162, 11:88,<br>13:177, 16:44                                                                              |
|                         | 220 bp |               | *15:123                                                                                                 |                                                                                                               |
| <b>13<sup>4</sup></b>   | 110 bp | 430 bp        | *15:32                                                                                                  |                                                                                                               |
|                         | 215 bp |               | *15:10, 15:84, 15:103, 15:117, 15:120                                                                   |                                                                                                               |
| <b>14</b>               | 200 bp | 430 bp        | *15:13, 15:17N                                                                                          |                                                                                                               |
| <b>15</b>               | 150 bp | <b>515 bp</b> | *15:68, 15:136                                                                                          |                                                                                                               |
|                         | 205 bp |               | *15:10, 15:21, 15:27, 15:34, 15:54,<br>15:66:01-15:66:02, 15:122                                        | *16:05:01-16:05:02,<br>16:07                                                                                  |
| <b>16<sup>4</sup></b>   | 110 bp | 430 bp        | *15:102                                                                                                 |                                                                                                               |
|                         | 210 bp |               | *15:08, 15:129N                                                                                         |                                                                                                               |
|                         | 260 bp |               | *15:74                                                                                                  |                                                                                                               |
| <b>17</b>               | 160 bp | 430 bp        | *15:18                                                                                                  |                                                                                                               |
|                         | 220 bp |               | *15:123                                                                                                 |                                                                                                               |
| <b>18</b>               | 160 bp | <b>515 bp</b> | *15:115N                                                                                                |                                                                                                               |
|                         | 205 bp |               | *15:20                                                                                                  | *04:226                                                                                                       |
| <b>19<sup>6</sup></b>   | 175 bp | 430 bp        | *15:30                                                                                                  |                                                                                                               |
|                         | 225 bp |               | *15:21, 15:25, 15:37:01-15:37:02, 15:57,<br>15:100, 15:104:01-15:104:02                                 | *16:04:01-16:04:02,<br>16:18, 16:46                                                                           |
| <b>20</b>               | 165 bp | 430 bp        | *15:22                                                                                                  |                                                                                                               |
|                         | 200 bp |               | *15:25, 15:27, 15:34, 15:54, 15:66:01-<br>15:66:02, 15:100, 15:112, 15:120                              |                                                                                                               |
| <b>21</b>               | 165 bp | 430 bp        | *15:03:01:01-15:03:02, 15:23, 15:78, 15:94,<br>15:130, 15:137N, 15:142                                  |                                                                                                               |
| <b>22<sup>5</sup></b>   | 170 bp | 430 bp        | *15:24, 15:93, 15:138N                                                                                  |                                                                                                               |
|                         | 200 bp |               | *15:28, 15:40                                                                                           |                                                                                                               |
| <b>23<sup>4</sup></b>   | 125 bp | 430 bp        | *15:55                                                                                                  |                                                                                                               |
|                         | 170 bp |               | *15:29                                                                                                  |                                                                                                               |
|                         | 220 bp |               | *15:25                                                                                                  |                                                                                                               |
| <b>24</b>               | 140 bp | 430 bp        | *15:61                                                                                                  |                                                                                                               |
|                         | 200 bp |               | *15:26, 15:40, 15:43                                                                                    |                                                                                                               |
| <b>25<sup>4,6</sup></b> | 70 bp  | 430 bp        | *15:39                                                                                                  |                                                                                                               |

101.125-24/06 – including *Taq* pol., IFU-01.  
101.125-24u/06u – without *Taq* pol., IFU-02.

Visit <https://labproducts.caredx.com> for  
“Instructions for Use” (IFU)

Lot No.: **7G1**

Lot-specific information

|                   |        |        |                                |                               |
|-------------------|--------|--------|--------------------------------|-------------------------------|
|                   | 135 bp |        | *15:61, 15:67                  |                               |
|                   | 195 bp |        | *15:33, 15:113N                |                               |
| 26 <sup>4,6</sup> | 75 bp  | 430 bp | *15:38, 15:63                  |                               |
|                   | 125 bp |        | *15:75                         |                               |
|                   | 195 bp |        | *15:134N                       |                               |
| 27 <sup>4,6</sup> | 110 bp | 430 bp | *15:64                         |                               |
|                   | 135 bp |        | *15:67                         |                               |
|                   | 175 bp |        | *15:14, 15:42, 15:93           |                               |
| 28 <sup>6</sup>   | 135 bp | 430 bp | *15:16, 15:118                 |                               |
|                   | 255 bp |        | *15:35, 15:96                  | *16:34, 16:38:01-<br>16:38:02 |
| 29 <sup>4</sup>   | 100 bp | 430 bp | *15:137N                       |                               |
|                   | 200 bp |        | *15:41, 15:50N, 15:80N, 15:112 |                               |
| 30                | 165 bp | 430 bp | *15:45                         |                               |
|                   | 215 bp |        | *15:44, 15:129N                |                               |
| 31 <sup>4</sup>   | 85 bp  | 515bp  | *15:46                         |                               |
| 32 <sup>7</sup>   | -      | -      | <b>Negative Control</b>        |                               |

<sup>1</sup>Alleles are assigned by the presence of specific PCR product(s). However, the sizes of the specific PCR products may be helpful in the interpretation of DRB1\*15 SSP subtypings.

When the primers in a primer mix can give rise to HLA-specific PCR products of more than one length this is indicated if the size difference is more than 20 base pairs. Size differences of 20 base pairs or less are not given. For high resolution SSP kits, the alleles listed are specified according to amplicon length.

Nonspecific amplifications, i.e. a ladder or a smear of bands, may sometimes be seen. GC-rich primers have a higher tendency of giving rise to nonspecific amplifications than other primers.

PCR fragments longer than the control bands may sometimes be observed. Such bands should be disregarded and do not influence the interpretation of the SSP typings.

PCR fragments migrating faster than the control bands, but slower than a 400 bp fragment may be seen in some gel read-outs. Such bands can be disregarded and do not influence the interpretation of the SSP typings. Some primers may give rise to primer oligomer artifacts. Sometimes this phenomenon is an inherent feature of the primer pair(s) of a primer mix. More often it is due to other factors such as too low amount of DNA in the PCR reactions, taking too long time in setting up the PCR reactions, working at elevated room temperature or using thermal cyclers that are not pre-heated.

<sup>2</sup>The internal positive control primer pairs amplify segments of the human growth hormone gene. The internal positive control bands are 430 or 515 base pairs respectively, well distribution as outlined in the table. Well number 1 contains the longer, 515 bp, internal positive control band. The well distribution of the internal controls can help in orientation of the kit on gel photo, as well as allow for kit identification. In the presence of a specific amplification the intensity of the control band often decreases.

<sup>3</sup>For several DRB1 alleles 1st and/or 3rd exon(s) and beyond, as well as intron nucleotide sequences, are not available. In these instances it is not known whether some of the primers of the SSP sets are completely matched with the target sequences or not. Assumption is made that unknown sequences in these regions are conserved within allelic groups.

<sup>4</sup>HLA-specific PCR products shorter than 125 base pairs have a lower intensity and are less sharp than longer PCR products.

<sup>5</sup>Primer mix 22 may have a tendency of unspecific amplification.

<sup>6</sup>Primer mixes 19, 25, 26, 27 and 28 have a tendency to giving rise to primer oligomer formation.

<sup>7</sup>Primer mix 32 contains a negative control, which will amplify more than 95% of HLA amplicons as well as the amplicons generated by the control primer pairs matching the human growth hormone gene. HLA-specific PCR product sizes range from 75 to 200 base pairs and the PCR product generated by the HGH positive control primer pair is 430 base pairs.

‘w’, may be weakly amplified.

101.125-24/06 – including *Taq* pol., IFU-01.  
101.125-24u/06u – without *Taq* pol., IFU-02.

Visit <https://labproducts.caredx.com> for  
“Instructions for Use” (IFU)

Lot No.: **7G1**

Lot-specific information

## PRIMER SPECIFICATION

| Well No.                                 | 1                    | 2                    | 3                    | 4                    | 5                    | 6                    | 7                    | 8                    | 9                    | 10                   | 11                   | 12                   |
|------------------------------------------|----------------------|----------------------|----------------------|----------------------|----------------------|----------------------|----------------------|----------------------|----------------------|----------------------|----------------------|----------------------|
| Length of spec. PCR product              | 140                  | 200                  | 150                  | 260                  | 260                  | 150                  | 135                  | 200                  | 90                   | 170                  | 155                  | 95                   |
|                                          |                      |                      |                      |                      |                      | 180                  |                      | 245                  | 165                  | 200                  | 200                  | 220                  |
|                                          |                      |                      |                      |                      |                      |                      |                      |                      | 215                  |                      |                      |                      |
| Length of int. pos. control <sup>1</sup> | 515                  | 515                  | 430                  | 430                  | 430                  | 430                  | 430                  | 430                  | 430                  | 430                  | 430                  | 430                  |
| 5'-primer(s) <sup>2</sup>                | 13(126)<br>5'-Agg 3' | 13(126)<br>5'-Agg 3' | 28(171)<br>5'-gAT 3' | 13(126)<br>5'-Agg 3' | 13(126)<br>5'-Agg 3' | 20(146)<br>5'-TgC 3' | 13(126)<br>5'-Agg 3' | 13(126)<br>5'-Agg 3' | 6(103)<br>5'-CAG 3'  | 13(126)<br>5'-Agg 3' | 13(126)<br>5'-Agg 3' | 30(175)<br>5'-gAT 3' |
|                                          | 13(126)<br>5'-AgA 3' | 13(126)<br>5'-AgA 3' | 29(174)<br>5'-Agg 3' | 13(126)<br>5'-AgA 3' |                      | 30(175)<br>5'-gAC 3' |                      |                      | 10(115)<br>5'-ggT 3' |                      |                      |                      |
|                                          |                      |                      | 30(175)<br>5'-gAT 3' |                      |                      |                      |                      |                      | 12(121)<br>5'-CTg 3' |                      |                      |                      |
|                                          |                      |                      |                      |                      |                      |                      |                      |                      | 25(160)<br>5'-TgC 3' |                      |                      |                      |
|                                          |                      |                      |                      |                      |                      |                      |                      |                      | 50(236)<br>5'-ggg 3' |                      |                      |                      |
|                                          |                      |                      |                      |                      |                      |                      |                      |                      | 50(236)<br>5'-ggC 3' |                      |                      |                      |
| 3'-primer(s) <sup>3</sup>                | 47(227)<br>5'-ggA 3' | 67(286)<br>5'-gAT 3' | 67(286)<br>5'-gAT 3' | 86(344)<br>5'-CCA 3' | 86(344)<br>5'-CAC 3' | 67(286)<br>5'-gAT 3' | 44(218)<br>5'-CCC 3' | 67(286)<br>5'-gAg 3' | 67(286)<br>5'-gAT 3' | 57(256)<br>5'-gCT 3' | 51(239)<br>5'-CCC 3' | 48(230)<br>5'-CCT 3' |
|                                          |                      |                      |                      |                      |                      |                      | 47(227)<br>5'-ggT 3' | 67(286)<br>5'-gAg 3' |                      | 67(286)<br>5'-gAA 3' | 67(287)<br>5'-ggg 3' | 90(356)<br>5'-CTA 3' |
|                                          |                      |                      |                      |                      |                      |                      |                      | 81(329)<br>5'-TgC 3' |                      |                      |                      |                      |
| Well No.                                 | 1                    | 2                    | 3                    | 4                    | 5                    | 6                    | 7                    | 8                    | 9                    | 10                   | 11                   | 12                   |

| Well No.                                 | 13                   | 14                   | 15                   | 16                   | 17                   | 18                   | 19                   | 20                   | 21                   | 22                   | 23                   | 24                   |
|------------------------------------------|----------------------|----------------------|----------------------|----------------------|----------------------|----------------------|----------------------|----------------------|----------------------|----------------------|----------------------|----------------------|
| Length of spec. PCR product              | 110                  | 200                  | 150                  | 110                  | 160                  | 160                  | 175                  | 165                  | 165                  | 170                  | 125                  | 140                  |
|                                          | 215                  |                      | 205                  | 210                  | 220                  | 205                  | 225                  | 200                  |                      | 200                  | 170                  | 200                  |
|                                          |                      |                      |                      | 260                  |                      |                      |                      |                      |                      |                      | 220                  |                      |
| Length of int. pos. control <sup>1</sup> | 430                  | 430                  | 515                  | 430                  | 430                  | 515                  | 430                  | 430                  | 430                  | 430                  | 430                  | 430                  |
| 5'-primer(s) <sup>2</sup>                | 13(126)<br>5'-Agg 3' | 13(126)<br>5'-Agg 3' | 13(126)<br>5'-Agg 3' | 13(126)<br>5'-Agg 3' | 30(175)<br>5'-gAT 3' | 30(175)<br>5'-gAT 3' | 13(126)<br>5'-Agg 3' | 13(126)<br>5'-Agg 3' | 30(175)<br>5'-gAC 3' | 13(125)<br>5'-gAA 3' | 13(126)<br>5'-Agg 3' | 15(131)<br>5'-gTA 3' |
|                                          |                      |                      |                      |                      |                      |                      |                      |                      |                      | 15(131)<br>5'-gTA 3' |                      | 19(143)<br>5'-CAg 3' |
|                                          |                      |                      |                      |                      |                      |                      |                      |                      |                      | 23(154)<br>5'-AgT 3' |                      | 21(148)<br>5'-ggT 3' |
|                                          |                      |                      |                      |                      |                      |                      |                      |                      |                      | 25(160)<br>5'-TgC 3' |                      | 38(199)<br>5'-CCA 3' |
|                                          |                      |                      |                      |                      |                      |                      |                      |                      |                      | 25(161)<br>5'-gCA 3' |                      |                      |
| 3'-primer(s) <sup>3</sup>                | 37(196)<br>5'-ggT 3' | 64(277)<br>5'-Tg 3'  | 49(233)<br>5'-CCT 3' | 37(196)<br>5'-gAg 3' | 70(295)<br>5'-CTC 3' | 69(293)<br>5'-CT 3'  | 57(258)<br>5'-gCC 3' | 55(251)<br>5'-gCA 3' | 71(298)<br>5'-CgC 3' | 67(286)<br>5'-gAT 3' | 41(208)<br>5'-gTg 3' | 71(298)<br>5'-CgC 3' |
|                                          | 71(299)<br>5'-gCT 3' | 69(294)<br>5'-TgT 3' | 67(286)<br>5'-gAT 3' | 66(283)<br>5'-gTT 3' | 90(356)<br>5'-CTA 3' | 84(337)<br>5'-CCT 3' | 74(307)<br>5'-CAG 3' | 65(282)<br>5'-TCA 3' |                      |                      | 57(257)<br>5'-CgA 3' |                      |
|                                          |                      |                      | 67(286)<br>5'-gAT 3' | 72(301)<br>5'-Cg 3'  |                      |                      | 77(317)<br>5'-AgT 3' | 67(286)<br>5'-gAg 3' |                      |                      | 73(305)<br>5'-ggC 3' |                      |
|                                          |                      |                      | 70(296)<br>5'-TCT 3' | 87(347)<br>5'-TCC 3' |                      |                      | 77(317)<br>5'-AAT 3' |                      |                      |                      |                      |                      |
| Well No.                                 | 13                   | 14                   | 15                   | 16                   | 17                   | 18                   | 19                   | 20                   | 21                   | 22                   | 23                   | 24                   |

101.125-24/06 – including *Taq* pol., IFU-01.  
101.125-24u/06u – without *Taq* pol., IFU-02.

Visit <https://labproducts.caredx.com> for  
“Instructions for Use” (IFU)

Lot No.: **7G1**

Lot-specific information

| Well No.                                    | 25                   | 26                    | 27                   | 28                   | 29                   | 30                   | 31                   |
|---------------------------------------------|----------------------|-----------------------|----------------------|----------------------|----------------------|----------------------|----------------------|
| Length of spec.                             | 70                   | 75                    | 110                  | 135                  | 100                  | 165                  | 85                   |
| PCR product                                 | 135                  | 125                   | 135                  | 255                  | 200                  | 215                  |                      |
|                                             | 195                  | 195                   | 175                  |                      |                      |                      |                      |
| Length of int.<br>pos. control <sup>1</sup> | 430                  | 430                   | 430                  | 430                  | 430                  | 430                  | 515                  |
| 5'-primer(s) <sup>2</sup>                   | 10(115)<br>5'-ggT 3' | 19(144)<br>5'-...T 3' | 20(146)<br>5'-TgC 3' | 13(126)<br>5'-Agg 3' | 13(126)<br>5'-Agg 3' | 13(126)<br>5'-Agg 3' | 72(301)<br>5'-Cgg 3' |
|                                             | 15(130)<br>5'-AgA 3' | 43(216)<br>5'-gAg 3'  | 23(154)<br>5'-AgT 3' |                      |                      |                      |                      |
|                                             | 36(194)<br>5'-ggg 3' | 60(266)<br>5'-gTC 3'  | 23(155)<br>5'-gCA 3' |                      |                      |                      |                      |
|                                             | 38(199)<br>5'-CCA 3' |                       | 36(194)<br>5'-ggg 3' |                      |                      |                      |                      |
|                                             | 58(259)<br>5'-ACA 3' |                       | 44(217)<br>5'-ACA 3' |                      |                      |                      |                      |
| 3'-primer(s) <sup>3</sup>                   | 67(286)<br>5'-gAT 3' | 71(298)<br>5'-CgC 3'  | 67(286)<br>5'-gAT 3' | 45(220)<br>5'-CCT 3' | 34(187)<br>5'-CTA 3' | 55(250)<br>5'-CCA 3' | 86(344)<br>5'-CCA 3' |
|                                             |                      |                       |                      | 47(226)<br>5'-gAg 3' | 62(273)<br>5'-CTC 3' | 72(301)<br>5'-Cg 3'  |                      |
|                                             |                      |                       |                      | 85(341)<br>5'-CAg 3' | 65(282)<br>5'-TCA 3' | 73(305)<br>5'-CgA 3' |                      |
|                                             |                      |                       |                      |                      | 70(295)<br>5'-Tg 3'  |                      |                      |
| Well No.                                    | 25                   | 26                    | 27                   | 28                   | 29                   | 30                   | 31                   |

<sup>1</sup>The internal positive control primer pairs amplify segments of the human growth hormone gene. The internal positive control bands are 430 or 515 base pairs respectively, well distribution as outlined in the table. Well number 1 contains the longer, 515 bp, internal positive control band. The well distribution of the internal controls can help in orientation of the kit on gel photo, as well as allow for kit identification. In the presence of a specific amplification the intensity of the control band often decreases.

<sup>2</sup>The nucleotide position matching the specificity-determining 3'-end of the primer is given. Nucleotide numbering as on the [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla) web site. The sequence of the 3 terminal nucleotides of the primer is given.

<sup>3</sup>The nucleotide position matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Nucleotide numbering as on the [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla) web site. The sequence of the 3 terminal nucleotides of the primer is given.

101.125-24/06 – including *Taq* pol., IFU-01.  
101.125-24u/06u – without *Taq* pol., IFU-02.

Visit <https://labproducts.caredx.com> for  
“Instructions for Use” (IFU)

Lot No.: **7G1**

Lot-specific information

| CELL LINE VALIDATION SHEET             |                             |            |        |            |           |           |           |           |           |           |           |           |           |           |           |           |           |           |           |           |
|----------------------------------------|-----------------------------|------------|--------|------------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
| DRB1*15 SSP subtyping kit <sup>2</sup> |                             |            |        |            |           |           |           |           |           |           |           |           |           |           |           |           |           |           |           |           |
|                                        |                             |            |        |            | Well      |           |           |           |           |           |           |           |           |           |           |           |           |           |           |           |
|                                        |                             |            |        |            | 1         | 2         | 3         | 4         | 5         | 6         | 7         | 8         | 9         | 10        | 11        | 12        | 13        | 14        | 15        | 16        |
|                                        |                             |            |        | Prod. No.: | 201780101 | 201780102 | 201895703 | 201780104 | 201780105 | 201436306 | 201780107 | 201436308 | 201780109 | 201436310 | 201436332 | 201780112 | 201780113 | 201780114 | 201780115 | 201895716 |
|                                        | IHCW cell line <sup>1</sup> |            | DRB1   |            |           |           |           |           |           |           |           |           |           |           |           |           |           |           |           |           |
| 1                                      | 9001                        | SA         | *01:01 |            | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 2                                      | 9280                        | LK707      | *15:02 | *04:05     | +         | +         | +         | -         | +         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 3                                      | 9011                        | E4181324   | *15:02 |            | +         | +         | +         | -         | +         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 4                                      | 9275                        | GU373      | *03:01 |            | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 5                                      | 9009                        | KAS011     | *16:01 |            | -         | -         | -         | -         | +         | -         | +         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 6                                      | 9353                        | SM         | *04:07 | *08:03     | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 7                                      | 9020                        | QBL        | *03:01 |            | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 8                                      | 9025                        | DEU        | *04:01 |            | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 9                                      | 9026                        | YAR        | *04:02 |            | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 10                                     | 9107                        | LKT3       | *04:05 |            | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 11                                     | 9051                        | PITOUT     | *07:01 |            | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 12                                     | 9052                        | DBB        | *07:01 |            | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 13                                     | 9004                        | JESTHOM    | *01:01 |            | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 14                                     | 9071                        | OLGA       | *08:02 |            | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 15                                     | 9075                        | DKB        | *09:01 |            | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 16                                     | 9037                        | SWEIG007   | *11:01 |            | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 17                                     | 9282                        | CTM3953540 | *03:01 | *13:01     | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 18                                     | 9257                        | 32367      | *09:01 | *11:01     | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 19                                     | 9038                        | BM16       | *12:01 |            | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 20                                     | 9059                        | SLE005     | *13:02 |            | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 21                                     | 9064                        | AMALA      | *14:02 |            | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 22                                     | 9056                        | KOSE       | *13:02 | *14:01     | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 23                                     | 9124                        | IHL        | *08:03 | *14:14     | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 24                                     | 9035                        | JBUSH      | *11:01 |            | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 25                                     | 9049                        | IBW9       | *07:01 |            | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 26                                     | 9285                        | WT49       | *03:01 |            | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 27                                     | 9191                        | CH1007     | *04:05 | *10:01     | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 28                                     | 9320                        | BEL5GB     | *04:16 | *07:01     | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 29                                     | 9050                        | MOU        | *07:01 |            | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 30                                     | 9021                        | RSH        | *03:02 |            | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 31                                     | 9019                        | DUCAF      | *03:01 |            | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 32                                     | 9297                        | HAG        | *13:03 |            | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 33                                     | 9098                        | MT14B      | *04:04 |            | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 34                                     | 9104                        | DHIF       | *11:01 |            | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 35                                     | 9302                        | SSTO       | *04:03 |            | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 36                                     | 9024                        | KT17       | *04:03 | *04:06     | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 37                                     | 9065                        | HHKB       | *13:01 |            | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 38                                     | 9099                        | LZL        | *14:02 |            | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 39                                     | 9315                        | CML        | *03:01 | *04:01     | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 40                                     | 9134                        | WHONP199   | *07:01 | *09:01     | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 41                                     | 9055                        | H0301      | *13:02 |            | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 42                                     | 9066                        | TAB089     | *08:03 |            | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 43                                     | 9076                        | T7526      | *09:01 |            | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 44                                     | 9057                        | TEM        | *14:01 |            | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 45                                     | 9239                        | SHJO       | *07:01 |            | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 46                                     | 9013                        | SCHU       | *15:01 |            | +         | +         | +         | +         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 47                                     | 9045                        | TUBO       | *11:04 | *12:01     | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |
| 48                                     | 9303                        | TER-ND     | *01:03 |            | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         |

101.125-24/06 – including *Taq* pol., IFU-01.  
101.125-24u/06u – without *Taq* pol., IFU-02.

Visit <https://labproducts.caredx.com> for  
“Instructions for Use” (IFU)

Lot No.: **7G1**

Lot-specific information

| CELL LINE VALIDATION SHEET             |      |            |        |        |  |  |  |  |  |                            |           |           |           |
|----------------------------------------|------|------------|--------|--------|--|--|--|--|--|----------------------------|-----------|-----------|-----------|
| DRB1*15 SSP subtyping kit <sup>2</sup> |      |            |        |        |  |  |  |  |  |                            |           |           |           |
|                                        |      |            |        |        |  |  |  |  |  | Well                       |           |           |           |
|                                        |      |            |        |        |  |  |  |  |  | 17                         | 18        | 19        | 20        |
|                                        |      |            |        |        |  |  |  |  |  | 21                         | 22        | 23        | 24        |
|                                        |      |            |        |        |  |  |  |  |  | 25                         | 26        | 27        | 28        |
|                                        |      |            |        |        |  |  |  |  |  | 29                         | 30        | 31        | 31        |
|                                        |      |            |        |        |  |  |  |  |  | 201780117                  | 201780118 | 201436319 | 201558220 |
|                                        |      |            |        |        |  |  |  |  |  | 201436321                  | 201895722 | 201436323 | 201436324 |
|                                        |      |            |        |        |  |  |  |  |  | 201558225                  | 201780126 | 201558227 | 201895728 |
|                                        |      |            |        |        |  |  |  |  |  | 201780129                  | 201780130 | 201780131 | 201780131 |
|                                        |      |            |        |        |  |  |  |  |  | Prod. No.:                 |           |           |           |
|                                        |      |            |        |        |  |  |  |  |  | IHW cell line <sup>1</sup> |           |           |           |
|                                        |      |            |        |        |  |  |  |  |  | DRB1                       |           |           |           |
| 1                                      | 9001 | SA         | *01:01 |        |  |  |  |  |  | -                          | -         | -         | -         |
| 2                                      | 9280 | LK707      | *15:02 | *04:05 |  |  |  |  |  | -                          | -         | -         | -         |
| 3                                      | 9011 | E4181324   | *15:02 |        |  |  |  |  |  | -                          | -         | -         | -         |
| 4                                      | 9275 | GU373      | *03:01 |        |  |  |  |  |  | -                          | -         | -         | -         |
| 5                                      | 9009 | KAS011     | *16:01 |        |  |  |  |  |  | -                          | -         | -         | -         |
| 6                                      | 9353 | SM         | *04:07 | *08:03 |  |  |  |  |  | -                          | -         | -         | -         |
| 7                                      | 9020 | QBL        | *03:01 |        |  |  |  |  |  | -                          | -         | -         | -         |
| 8                                      | 9025 | DEU        | *04:01 |        |  |  |  |  |  | -                          | -         | -         | -         |
| 9                                      | 9026 | YAR        | *04:02 |        |  |  |  |  |  | -                          | -         | -         | -         |
| 10                                     | 9107 | LKT3       | *04:05 |        |  |  |  |  |  | -                          | -         | -         | -         |
| 11                                     | 9051 | PITOUT     | *07:01 |        |  |  |  |  |  | -                          | -         | -         | -         |
| 12                                     | 9052 | DBB        | *07:01 |        |  |  |  |  |  | -                          | -         | -         | -         |
| 13                                     | 9004 | JESTHOM    | *01:01 |        |  |  |  |  |  | -                          | -         | -         | -         |
| 14                                     | 9071 | OLGA       | *08:02 |        |  |  |  |  |  | -                          | -         | -         | -         |
| 15                                     | 9075 | DKB        | *09:01 |        |  |  |  |  |  | -                          | -         | -         | -         |
| 16                                     | 9037 | SWEIG007   | *11:01 |        |  |  |  |  |  | -                          | -         | -         | -         |
| 17                                     | 9282 | CTM3953540 | *03:01 | *13:01 |  |  |  |  |  | -                          | -         | -         | -         |
| 18                                     | 9257 | 32367      | *09:01 | *11:01 |  |  |  |  |  | -                          | -         | -         | -         |
| 19                                     | 9038 | BM16       | *12:01 |        |  |  |  |  |  | -                          | -         | -         | -         |
| 20                                     | 9059 | SLE005     | *13:02 |        |  |  |  |  |  | -                          | -         | -         | -         |
| 21                                     | 9064 | AMALA      | *14:02 |        |  |  |  |  |  | -                          | -         | -         | -         |
| 22                                     | 9056 | KOSE       | *13:02 | *14:01 |  |  |  |  |  | -                          | -         | -         | -         |
| 23                                     | 9124 | IHL        | *08:03 | *14:14 |  |  |  |  |  | -                          | -         | -         | -         |
| 24                                     | 9035 | JBUSH      | *11:01 |        |  |  |  |  |  | -                          | -         | -         | -         |
| 25                                     | 9049 | IBW9       | *07:01 |        |  |  |  |  |  | -                          | -         | -         | -         |
| 26                                     | 9285 | WT49       | *03:01 |        |  |  |  |  |  | -                          | -         | -         | -         |
| 27                                     | 9191 | CH1007     | *04:05 | *10:01 |  |  |  |  |  | -                          | -         | -         | -         |
| 28                                     | 9320 | BEL5GB     | *04:16 | *07:01 |  |  |  |  |  | -                          | -         | -         | -         |
| 29                                     | 9050 | MOU        | *07:01 |        |  |  |  |  |  | -                          | -         | -         | -         |
| 30                                     | 9021 | RSH        | *03:02 |        |  |  |  |  |  | -                          | -         | -         | -         |
| 31                                     | 9019 | DUCAF      | *03:01 |        |  |  |  |  |  | -                          | -         | -         | -         |
| 32                                     | 9297 | HAG        | *13:03 |        |  |  |  |  |  | -                          | -         | -         | -         |
| 33                                     | 9098 | MT14B      | *04:04 |        |  |  |  |  |  | -                          | -         | -         | -         |
| 34                                     | 9104 | DHIF       | *11:01 |        |  |  |  |  |  | -                          | -         | -         | -         |
| 35                                     | 9302 | SSTO       | *04:03 |        |  |  |  |  |  | -                          | -         | -         | -         |
| 36                                     | 9024 | KT17       | *04:03 | *04:06 |  |  |  |  |  | -                          | -         | -         | -         |
| 37                                     | 9065 | HHKB       | *13:01 |        |  |  |  |  |  | -                          | -         | -         | -         |
| 38                                     | 9099 | LZL        | *14:02 |        |  |  |  |  |  | -                          | -         | -         | -         |
| 39                                     | 9315 | CML        | *03:01 | *04:01 |  |  |  |  |  | -                          | -         | -         | -         |
| 40                                     | 9134 | WHONP199   | *07:01 | *09:01 |  |  |  |  |  | -                          | -         | -         | -         |
| 41                                     | 9055 | H0301      | *13:02 |        |  |  |  |  |  | -                          | -         | -         | -         |
| 42                                     | 9066 | TAB089     | *08:03 |        |  |  |  |  |  | -                          | -         | -         | -         |
| 43                                     | 9076 | T7526      | *09:01 |        |  |  |  |  |  | -                          | -         | -         | -         |
| 44                                     | 9057 | TEM        | *14:01 |        |  |  |  |  |  | -                          | -         | -         | -         |
| 45                                     | 9239 | SHJO       | *07:01 |        |  |  |  |  |  | -                          | -         | -         | -         |
| 46                                     | 9013 | SCHU       | *15:01 |        |  |  |  |  |  | -                          | -         | -         | -         |
| 47                                     | 9045 | TUBO       | *11:04 | *12:01 |  |  |  |  |  | -                          | -         | -         | -         |
| 48                                     | 9303 | TER-ND     | *01:03 |        |  |  |  |  |  | -                          | -         | -         | -         |

101.125-24/06 – including *Taq* pol., IFU-01.  
101.125-24u/06u – without *Taq* pol., IFU-02.

Visit <https://labproducts.caredx.com> for  
“Instructions for Use” (IFU)

**Lot No.: 7G1**

**Lot-specific information**

<sup>1</sup>The provided cell line HLA specificities are retrieved from the <http://www.ihwg.org/hla> web site. The specificity of an individual cell line may thus be subject to change.

<sup>2</sup>The specificity of each primer solution in the kit has been tested against 48 well characterized cell line DNAs and where applicable, additional cell line DNAs.

No DNAs carrying the alleles to be amplified by primer solutions 8, 9, 11 to 20 and 22 to 31 were available. The specificities of the primers in primer solutions 8, 11 to 13, 15, 16, 19, 20, 22, 23 and 28 were tested by separately adding one or three 5'-primers, respectively one 3'-primer. In primer solutions 14, 17, 18, 29 and 30 it was only possible to test the 5'-primers, the 3'-primers were not possible to test. In primer solutions 9, 24 to 27 and 31 it was only possible to test the 3'-primers, the 5'-primers were not possible to test. In primer solutions 3, 6 and 22 one or more of the 5'-primers were not possible to test, and in primer solutions 2, 3, 7, 8, 11 to 13, 15, 16, 19, 20, 23 and 28 one, two or three of the 3'-primers were not possible to test. Additional primers in primer solution 10 were tested by separately adding one 5'-primer.

101.125-24/06 – including *Taq* pol., IFU-01.  
101.125-24u/06u – without *Taq* pol., IFU-02.

Visit <https://labproducts.caredx.com> for  
“Instructions for Use” (IFU)

Lot No.: **7G1**

Lot-specific information

101.125-24/06 – including *Taq* pol., IFU-01.  
101.125-24u/06u – without *Taq* pol., IFU-02.

Visit <https://labproducts.caredx.com> for  
“Instructions for Use” (IFU)

Lot No.: **7G1**

Lot-specific information

## ADDRESSES:

### Manufacturer:

**CareDx AB**, Franzengatan 5, SE-112 51 Stockholm, Sweden.

**Tel:** +46-8-508 939 00

**Fax:** +46-8-717 88 18

**E-mail:** [orders-se@caredx.com](mailto:orders-se@caredx.com)

**Web page:** <https://labproducts.caredx.com/>

### Distributed by:

**CareDx GmbH**, Löwengasse 47 / 6, AT-1030 Vienna, Austria.

**Tel:** +43-1-710 15 00

**Fax:** +43-1-710 15 00 10

**E-mail:** [orders-at@caredx.com](mailto:orders-at@caredx.com)

**Web page:** <https://labproducts.caredx.com/>

**CareDx Lab Solutions Inc.**, 901 S. Bolmar St., Suite R, West Chester, PA 19382

**Tel:** 1-877-653-78171

**Fax:** 610-344-7989

**E-mail:** [orders-us@caredx.com](mailto:orders-us@caredx.com)

**Web page:** <https://labproducts.caredx.com/>

For information on CareDx distributors worldwide, contact **CareDx GmbH**.