

101.516-24/04 – including *Taq* polymerase, IFU-01
101.516-24u/04u – without *Taq* polymerase, IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: 3F7

Lot-specific information
***Olerup SSP[®]* HLA-B*15**

Product number:	101.516-24/04 – including <i>Taq</i> pol. 101.516-24u/04u – without <i>Taq</i> pol.
Lot number:	3F7
Expiry date:	2020-02-01
Number of tests:	24 tests – Product No. 101.516-24/24u 4 tests – Product No. 101.516-04/04u
Number of wells per test:	95+1
Storage - pre-aliquoted primers:	dark at -20°C
- PCR Master Mix:	-20°C

This Product Description is only valid for Lot No. 3F7.

Complete product documentation consists of generic Instructions for Use (IFU), lot specific Product Insert, Worksheet and Certificate.

CHANGES COMPARED TO THE PREVIOUS OLERUP SSP[®]
HLA-B*15 Lot (0D5)

The HLA-B*15 kit is updated for new alleles to enable separation of:

- Confirmed¹ alleles as listed in the IMGT/HLA database
- Polymorphisms in exons outside of the region encoding the peptide binding domain
- Null and Alternatively expressed alleles

The format of the Worksheet has been changed.

¹As described in section Uniquely Identified Alleles.

The HLA-B*15 specificity and interpretation tables have been updated for the HLA-B alleles described since the previous *Olerup SSP[®]* HLA-B*15 lot was made (**Lot No. 0D5**). The kit design is based on IMGT/HLA database 3.25.0.

101.516-24/04 – including *Taq* polymerase, IFU-01
101.516-24u/04u – without *Taq* polymerase, IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: 3F7

Lot-specific information

The primers of the wells detailed below have been exchanged, added or modified compared to the previous lot.

Well	5'-primer	3'-primer	rationale
1	-	Added	3'-primer added for the B*15:370.
24	-	Added, modified	Modified 3'-primer added for increased yield.
30	Added	-	5'-primer added for the B*15:375N allele.
42	Modified, added	-	5'-primer modified for improved HLA-specific amplification, 5'-primer added for the B*15:17:05 allele.
46	Added	-	5'-primer added for the B*15:17:05 allele.
47	Added	-	5'-primer added for the B*15:02:11 allele.
57	Added	-	5'-primer added for the B*15:369 allele.
59	Added	-	5'-primer added from well 79 for improved HLA-specific amplification.
63	Added	Added	Primer pair added for the B*15:315 allele.
64	Added	-	5'-primer added from well 79 for improved HLA-specific amplification.
65	Added	-	5'-primer added for the B*15:375N allele.
66	-	Added, modified	Modified 3'-primer added for increased yield. Strength of control band has been optimized.
72	Added	-	5'-primer added for the B*15:380N allele.
75	Added	-	5'-primer added for the B*15:367 allele.
77	Modified, added	-	Modified 5'-primer added for increased yield.
78	Added	-	5'-primer added for the B*15:369 allele.
79	Moved	Removed	5'-primer moved to wells 59 and 64 for improved HLA-specific amplification, 3'-primer removed.
80	Added	-	5'-primer added for the B*15:380N allele.
83	Exchanged	Exchanged	Primer pair exchanged for improved HLA-specific amplification.
88	-	Added	3'-primer added for the B*15:377Q allele.
89	-	Added	3'-primer added for the B*15:377Q allele.

101.516-24/04 – including *Taq* polymerase, IFU-01
101.516-24u/04u – without *Taq* polymerase, IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **3F7**

Lot-specific information

Well **96** contains Negative Control primer pairs, that will amplify more than 95% of the Olerup SSP[®] HLA Class I, DRB, DQB1, DPB1 and DQA1 amplicons as well as all the amplicons generated by the control primer pairs matching the human growth hormone gene.

HLA-specific PCR product sizes range from 75 to 200 base pairs.
The PCR product generated by the positive control primer pair is 430 base pairs.

Length of PCR product	105	200	105	80	75	80	85
5'-primer¹	164	340	440	45	45	43	36
	5'-CAC ^{3'}	5'-Agg ^{3'}	5'-TTA ^{3'}	5'-Tgg ^{3'}	5'-Tgg ^{3'}	5'-Tgg ^{3'}	5'-TAC ^{3'}
							36
							5'-TAT ^{3'}
3'-primer²	231	2nd I	507	59	58	57	47
	5'-TgC ^{3'}	5'-AAA ^{3'}	5'-TTg ^{3'}	5'-CTC ^{3'}	5'-ggC ^{3'}	5'-CTC ^{3'}	5'-ACA ^{3'}
							48
							5'-gCA ^{3'}
							48
							5'-gCC ^{3'}
							52
							5'-TgT ^{3'}
A*	+	+	+				
B*	+	+	+				
C*	+	+	+				
DRB1				+	+		
DRB3				+	+		
DRB5				+			
DQB1					+		
DPB1						+	
DQA1							+

¹The nucleotide position for HLA class I genes and the codon for HLA class II genes, in the 2nd or 3rd exon, matching the specificity-determining 3'-end of the primer is given. Nucleotide and codon numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

²The nucleotide position for HLA class I genes and the codon for HLA class II genes, in the 2nd or 3rd exon or the 2nd intron, matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Nucleotide and codon numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

101.516-24/04 – including *Taq* polymerase, IFU-01
101.516-24u/04u – without *Taq* polymerase, IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **3F7**

Lot-specific information

PRODUCT DESCRIPTION

HLA-B*15 SSP subtyping

CONTENT

The primer set contains 5'- and 3'-primers for identifying the B*15:01 to B*15:398 alleles.

PLATE LAYOUT

Each test consists of 96 PCR reactions in a 96 well cut PCR plate.

1	2	3	4	5	6	7	8
9	10	11	12	13	14	15	16
17	18	19	20	21	22	23	24
25	26	27	28	29	30	31	32
33	34	35	36	37	38	39	40
41	42	43	44	45	46	47	48
49	50	51	52	53	54	55	56
57	58	59	60	61	62	63	64
65	66	67	68	69	70	71	72
73	74	75	76	77	78	79	80
81	82	83	84	85	86	87	88
89	90	91	92	93	94	95	NC

The 96 well cut PCR plate is marked with 'HLA-B*15' in silver/gray ink.

Well No. 1 is marked with the Lot No. '3F7'.

Wells 1 to 95 – HLA-B*15 high resolution primers.

Well 96 – Negative Control (NC).

A faint row of numbers is seen between wells 1 and 2 or wells 7 and 8 of the PCR trays. These stem from the manufacture of the trays, and should be disregarded. The PCR plates are covered with a PCR-compatible foil.

INTERPRETATION

Due to the sharing of sequence motifs between HLA-B alleles, non-HLA-B*15 alleles will be amplified by some primer mixes. For further details see Specificity Table.

UNIQUELY IDENTIFIED ALLELES

All the HLA-B*15 alleles, i.e. **B*15:01 to B*15:398**, recognized by the HLA Nomenclature Committee in July 2016¹ will be amplified by the primers in the HLA-B*15 subtyping kit².

The HLA-B*15 kit enables separation of the confirmed HLA-B*15 alleles as listed in the IMGT/HLA database. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources. Current allele confirmation status for HLA-B*15 alleles is listed below.

The HLA-B*15 kit also enables identification of polymorphisms in exons outside of the region encoding the peptide binding domain and of null and alternatively expressed alleles.

101.516-24/04 – including *Taq* polymerase, IFU-01
101.516-24u/04u – without *Taq* polymerase, IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **3F7**

Lot-specific information

The following HLA-B*15 alleles can be distinguished by the different sizes of the HLA-specific PCR product:

Alleles	Primer mix	Alleles	Primer mix
B*15:92, 15:104	52	B*15:169, 15:225	87
B*15:120, 15:258N	93	B*15:206, 15:343	50
B*15:142, 15:147	78	B*15:208, 15:230	70
B*15:166, 15:193	86	B*15:301, 15:302N	64

¹HLA-B alleles listed on the IMGT/HLA web page 2016-July-14, release 3.25.0, www.ebi.ac.uk/imgt/hla.

²Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

101.516-24/04 – including *Taq* polymerase, IFU-01
 101.516-24u/04u – without *Taq* polymerase, IFU-02

Visit www.olerup.com for
 “Instructions for Use” (IFU)

Lot No.: **3F7**

Lot-specific information

ALLELE CONFIRMATION STATUS

Allele	Status ¹	Allele	Status ¹	Allele	Status ¹	Allele	Status ¹
B*15:01:01:01	Confirmed	B*15:02:09	Unconfirmed	B*15:18:01:01	Confirmed	B*15:54	Confirmed
B*15:01:01:02N	Confirmed	B*15:02:10	Unconfirmed	B*15:18:01:02	Unconfirmed	B*15:55	Confirmed
B*15:01:01:03	Unconfirmed	B*15:02:11	Unconfirmed	B*15:18:02	Confirmed	B*15:56	Confirmed
B*15:01:01:04	Confirmed	B*15:03:01:01	Confirmed	B*15:18:03	Confirmed	B*15:57	Unconfirmed
B*15:01:01:05	Confirmed	B*15:03:01:02	Confirmed	B*15:18:04	Confirmed	B*15:58	Confirmed
B*15:01:02	Confirmed	B*15:03:02	Confirmed	B*15:18:05	Unconfirmed	B*15:60	Unconfirmed
B*15:01:03	Unconfirmed	B*15:03:03	Confirmed	B*15:18:06	Unconfirmed	B*15:61	Confirmed
B*15:01:04	Unconfirmed	B*15:03:04	Confirmed	B*15:19	Confirmed	B*15:62	Confirmed
B*15:01:06	Confirmed	B*15:03:05	Confirmed	B*15:20	Confirmed	B*15:63	Confirmed
B*15:01:07	Unconfirmed	B*15:03:06	Unconfirmed	B*15:21	Confirmed	B*15:64:01	Confirmed
B*15:01:08	Unconfirmed	B*15:04:01:01	Confirmed	B*15:23	Confirmed	B*15:64:02	Confirmed
B*15:01:09	Confirmed	B*15:04:01:02	Unconfirmed	B*15:24:01	Confirmed	B*15:65	Unconfirmed
B*15:01:10	Confirmed	B*15:04:02	Confirmed	B*15:24:02	Unconfirmed	B*15:66	Unconfirmed
B*15:01:11	Confirmed	B*15:04:03	Unconfirmed	B*15:25:01	Confirmed	B*15:67	Confirmed
B*15:01:12	Unconfirmed	B*15:04:04	Confirmed	B*15:25:02	Confirmed	B*15:68	Confirmed
B*15:01:13	Confirmed	B*15:05:01	Confirmed	B*15:25:03	Unconfirmed	B*15:69	Confirmed
B*15:01:14	Confirmed	B*15:05:02	Confirmed	B*15:26N	Confirmed	B*15:70	Confirmed
B*15:01:15	Confirmed	B*15:06	Confirmed	B*15:27:01	Confirmed	B*15:71	Confirmed
B*15:01:16	Confirmed	B*15:07:01	Confirmed	B*15:27:02	Unconfirmed	B*15:72	Unconfirmed
B*15:01:17	Unconfirmed	B*15:07:02	Unconfirmed	B*15:27:03	Unconfirmed	B*15:73	Confirmed
B*15:01:18	Unconfirmed	B*15:07:03	Confirmed	B*15:28	Unconfirmed	B*15:74	Confirmed
B*15:01:19	Unconfirmed	B*15:08:01	Confirmed	B*15:29	Confirmed	B*15:75	Confirmed
B*15:01:20	Unconfirmed	B*15:08:02	Unconfirmed	B*15:30	Confirmed	B*15:76	Unconfirmed
B*15:01:21	Unconfirmed	B*15:09	Confirmed	B*15:31	Confirmed	B*15:77	Unconfirmed
B*15:01:22	Unconfirmed	B*15:10:01	Confirmed	B*15:32:01	Confirmed	B*15:78:01	Confirmed
B*15:01:23	Unconfirmed	B*15:10:02	Confirmed	B*15:32:02	Unconfirmed	B*15:78:02	Unconfirmed
B*15:01:24	Unconfirmed	B*15:10:03	Unconfirmed	B*15:33	Confirmed	B*15:78:03	Unconfirmed
B*15:01:25	Unconfirmed	B*15:10:04	Unconfirmed	B*15:34	Confirmed	B*15:79N	Unconfirmed
B*15:01:26	Unconfirmed	B*15:11:01	Confirmed	B*15:35	Confirmed	B*15:80	Unconfirmed
B*15:01:27	Unconfirmed	B*15:11:02	Confirmed	B*15:36	Unconfirmed	B*15:81	Unconfirmed
B*15:01:28	Unconfirmed	B*15:11:03	Unconfirmed	B*15:37	Confirmed	B*15:82	Unconfirmed
B*15:01:29	Unconfirmed	B*15:11:04	Unconfirmed	B*15:38:01	Confirmed	B*15:83	Unconfirmed
B*15:01:30	Unconfirmed	B*15:11:05	Unconfirmed	B*15:38:02	Confirmed	B*15:84	Unconfirmed
B*15:01:31	Confirmed	B*15:11:06	Unconfirmed	B*15:39:01	Confirmed	B*15:85	Unconfirmed
B*15:01:32	Unconfirmed	B*15:11:07	Unconfirmed	B*15:39:02	Unconfirmed	B*15:86	Confirmed
B*15:01:33	Unconfirmed	B*15:12	Confirmed	B*15:40:01	Confirmed	B*15:87	Unconfirmed
B*15:01:34	Unconfirmed	B*15:13:01	Confirmed	B*15:40:02	Unconfirmed	B*15:88	Unconfirmed
B*15:01:35	Unconfirmed	B*15:13:02	Unconfirmed	B*15:42	Confirmed	B*15:89	Confirmed
B*15:01:36	Confirmed	B*15:14	Confirmed	B*15:43	Unconfirmed	B*15:90	Unconfirmed
B*15:01:37	Unconfirmed	B*15:15	Confirmed	B*15:44	Unconfirmed	B*15:91	Unconfirmed
B*15:01:38	Unconfirmed	B*15:16:01:01	Confirmed	B*15:45	Confirmed	B*15:92	Unconfirmed
B*15:01:39	Unconfirmed	B*15:16:01:02	Confirmed	B*15:46	Confirmed	B*15:93	Unconfirmed
B*15:02:01	Confirmed	B*15:16:02	Confirmed	B*15:47:01	Confirmed	B*15:94N	Unconfirmed
B*15:02:02	Unconfirmed	B*15:16:03	Unconfirmed	B*15:47:02	Confirmed	B*15:95	Unconfirmed
B*15:02:03	Unconfirmed	B*15:17:01:01	Confirmed	B*15:48	Confirmed	B*15:96	Unconfirmed
B*15:02:04	Confirmed	B*15:17:01:02	Unconfirmed	B*15:49	Unconfirmed	B*15:97	Unconfirmed
B*15:02:05	Unconfirmed	B*15:17:02	Confirmed	B*15:50	Confirmed	B*15:98	Unconfirmed
B*15:02:06	Unconfirmed	B*15:17:03	Unconfirmed	B*15:51	Unconfirmed	B*15:99	Unconfirmed
B*15:02:07	Unconfirmed	B*15:17:04	Unconfirmed	B*15:52	Confirmed	B*15:101	Unconfirmed
B*15:02:08	Unconfirmed	B*15:17:05	Unconfirmed	B*15:53	Confirmed	B*15:102	Unconfirmed

101.516-24/04 – including *Taq* polymerase, IFU-01
101.516-24u/04u – without *Taq* polymerase, IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **3F7**

Lot-specific information

Allele	Status ¹	Allele	Status ¹	Allele	Status ¹	Allele	Status ¹
B*15:103	Unconfirmed	B*15:154	Confirmed	B*15:202	Unconfirmed	B*15:251	Unconfirmed
B*15:104	Unconfirmed	B*15:155	Confirmed	B*15:203	Unconfirmed	B*15:252	Unconfirmed
B*15:105	Unconfirmed	B*15:156	Unconfirmed	B*15:204	Confirmed	B*15:253	Unconfirmed
B*15:106	Confirmed	B*15:157	Confirmed	B*15:205	Unconfirmed	B*15:254	Unconfirmed
B*15:107	Unconfirmed	B*15:158	Confirmed	B*15:206	Confirmed	B*15:255	Unconfirmed
B*15:108	Confirmed	B*15:159	Confirmed	B*15:207	Confirmed	B*15:256	Unconfirmed
B*15:109	Unconfirmed	B*15:160	Confirmed	B*15:208	Confirmed	B*15:257	Unconfirmed
B*15:110	Unconfirmed	B*15:161	Confirmed	B*15:209N	Unconfirmed	B*15:258N	Unconfirmed
B*15:111N	Unconfirmed	B*15:162	Confirmed	B*15:210	Unconfirmed	B*15:259	Unconfirmed
B*15:112	Unconfirmed	B*15:163	Unconfirmed	B*15:211	Confirmed	B*15:260	Unconfirmed
B*15:113	Confirmed	B*15:164	Unconfirmed	B*15:212	Unconfirmed	B*15:261	Unconfirmed
B*15:114	Unconfirmed	B*15:165	Unconfirmed	B*15:213	Unconfirmed	B*15:262N	Unconfirmed
B*15:115	Unconfirmed	B*15:166	Unconfirmed	B*15:214	Unconfirmed	B*15:263	Confirmed
B*15:116	Unconfirmed	B*15:167	Confirmed	B*15:215	Unconfirmed	B*15:264	Unconfirmed
B*15:117	Unconfirmed	B*15:168	Unconfirmed	B*15:216	Unconfirmed	B*15:265	Unconfirmed
B*15:118	Unconfirmed	B*15:169	Unconfirmed	B*15:217	Unconfirmed	B*15:266	Confirmed
B*15:119	Unconfirmed	B*15:170	Confirmed	B*15:218Q	Unconfirmed	B*15:267	Unconfirmed
B*15:120	Unconfirmed	B*15:171	Confirmed	B*15:219	Unconfirmed	B*15:268	Confirmed
B*15:121	Unconfirmed	B*15:172	Confirmed	B*15:220	Confirmed	B*15:269	Unconfirmed
B*15:122	Unconfirmed	B*15:173	Confirmed	B*15:221	Unconfirmed	B*15:270	Unconfirmed
B*15:123	Unconfirmed	B*15:174	Unconfirmed	B*15:222	Unconfirmed	B*15:271	Confirmed
B*15:124	Unconfirmed	B*15:175	Unconfirmed	B*15:223	Unconfirmed	B*15:272N	Confirmed
B*15:125	Confirmed	B*15:176	Confirmed	B*15:224	Unconfirmed	B*15:273	Unconfirmed
B*15:126	Unconfirmed	B*15:177	Unconfirmed	B*15:225	Unconfirmed	B*15:274	Unconfirmed
B*15:127	Unconfirmed	B*15:178	Unconfirmed	B*15:226N	Unconfirmed	B*15:275:01	Unconfirmed
B*15:128	Unconfirmed	B*15:179:01	Unconfirmed	B*15:227	Unconfirmed	B*15:275:02	Unconfirmed
B*15:129	Unconfirmed	B*15:179:02	Unconfirmed	B*15:228	Confirmed	B*15:276	Unconfirmed
B*15:131	Unconfirmed	B*15:180	Confirmed	B*15:229	Confirmed	B*15:277	Confirmed
B*15:132	Confirmed	B*15:181N	Confirmed	B*15:230	Unconfirmed	B*15:278	Confirmed
B*15:133	Unconfirmed	B*15:182N	Unconfirmed	B*15:231	Unconfirmed	B*15:279	Unconfirmed
B*15:134	Unconfirmed	B*15:183	Confirmed	B*15:232	Confirmed	B*15:280	Unconfirmed
B*15:135	Confirmed	B*15:184	Unconfirmed	B*15:233	Unconfirmed	B*15:281	Unconfirmed
B*15:136	Unconfirmed	B*15:185	Unconfirmed	B*15:234	Unconfirmed	B*15:282	Confirmed
B*15:137	Unconfirmed	B*15:186	Unconfirmed	B*15:235	Confirmed	B*15:283	Unconfirmed
B*15:138	Unconfirmed	B*15:187	Unconfirmed	B*15:236	Confirmed	B*15:284	Unconfirmed
B*15:139	Unconfirmed	B*15:188	Unconfirmed	B*15:237	Confirmed	B*15:285	Confirmed
B*15:140	Unconfirmed	B*15:189	Confirmed	B*15:238	Unconfirmed	B*15:286	Unconfirmed
B*15:141	Confirmed	B*15:190N	Confirmed	B*15:239	Confirmed	B*15:287	Unconfirmed
B*15:142	Confirmed	B*15:191	Unconfirmed	B*15:240	Confirmed	B*15:288	Unconfirmed
B*15:143	Unconfirmed	B*15:192	Unconfirmed	B*15:241	Confirmed	B*15:289	Unconfirmed
B*15:144	Unconfirmed	B*15:193	Unconfirmed	B*15:242:01	Unconfirmed	B*15:290	Confirmed
B*15:145	Unconfirmed	B*15:194	Unconfirmed	B*15:242:02	Unconfirmed	B*15:291	Unconfirmed
B*15:146	Confirmed	B*15:195	Unconfirmed	B*15:243	Unconfirmed	B*15:292	Unconfirmed
B*15:147	Confirmed	B*15:196	Unconfirmed	B*15:244	Unconfirmed	B*15:293	Unconfirmed
B*15:148	Unconfirmed	B*15:197:01	Unconfirmed	B*15:245Q	Unconfirmed	B*15:294N	Unconfirmed
B*15:149N	Unconfirmed	B*15:197:02	Unconfirmed	B*15:246N	Unconfirmed	B*15:295	Unconfirmed
B*15:150	Confirmed	B*15:198	Unconfirmed	B*15:247	Unconfirmed	B*15:296	Unconfirmed
B*15:151	Confirmed	B*15:199	Unconfirmed	B*15:248	Confirmed	B*15:297	Unconfirmed
B*15:152	Unconfirmed	B*15:200	Unconfirmed	B*15:249	Unconfirmed	B*15:298	Unconfirmed
B*15:153	Confirmed	B*15:201	Unconfirmed	B*15:250	Unconfirmed	B*15:299	Unconfirmed

101.516-24/04 – including *Taq* polymerase, IFU-01
101.516-24u/04u – without *Taq* polymerase, IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **3F7**

Lot-specific information

Allele	Status ¹	Allele	Status ¹	Allele	Status ¹	Allele	Status ¹
B*15:300	Confirmed	B*15:330	Unconfirmed	B*15:360	Unconfirmed	B*15:390	Unconfirmed
B*15:301	Confirmed	B*15:331	Unconfirmed	B*15:361	Unconfirmed	B*15:391	Unconfirmed
B*15:302N	Unconfirmed	B*15:332	Confirmed	B*15:362	Unconfirmed	B*15:392	Unconfirmed
B*15:303	Unconfirmed	B*15:333	Unconfirmed	B*15:363:01	Unconfirmed	B*15:393	Unconfirmed
B*15:304N	Unconfirmed	B*15:334	Unconfirmed	B*15:364	Unconfirmed	B*15:394	Unconfirmed
B*15:305	Unconfirmed	B*15:335	Confirmed	B*15:365	Unconfirmed	B*15:395	Unconfirmed
B*15:306	Unconfirmed	B*15:336	Unconfirmed	B*15:366	Unconfirmed	B*15:396	Unconfirmed
B*15:307	Unconfirmed	B*15:337	Unconfirmed	B*15:367	Unconfirmed	B*15:397	Unconfirmed
B*15:308	Unconfirmed	B*15:338	Unconfirmed	B*15:368	Unconfirmed	B*15:398	Unconfirmed
B*15:309	Unconfirmed	B*15:339	Unconfirmed	B*15:369	Confirmed		
B*15:310	Unconfirmed	B*15:340	Unconfirmed	B*15:370	Unconfirmed		
B*15:311	Unconfirmed	B*15:341	Unconfirmed	B*15:371	Unconfirmed		
B*15:312	Unconfirmed	B*15:342	Confirmed	B*15:372	Unconfirmed		
B*15:313	Unconfirmed	B*15:343	Unconfirmed	B*15:373	Unconfirmed		
B*15:314	Unconfirmed	B*15:344	Unconfirmed	B*15:374	Unconfirmed		
B*15:315	Confirmed	B*15:345	Unconfirmed	B*15:375N	Unconfirmed		
B*15:316	Unconfirmed	B*15:346	Unconfirmed	B*15:376	Unconfirmed		
B*15:317	Unconfirmed	B*15:347	Confirmed	B*15:377Q	Unconfirmed		
B*15:318	Unconfirmed	B*15:348	Unconfirmed	B*15:378	Confirmed		
B*15:319	Confirmed	B*15:349	Unconfirmed	B*15:379	Confirmed		
B*15:320	Unconfirmed	B*15:350	Unconfirmed	B*15:380N	Confirmed		
B*15:321Q	Unconfirmed	B*15:351	Unconfirmed	B*15:381	Unconfirmed		
B*15:322	Unconfirmed	B*15:352	Unconfirmed	B*15:382	Unconfirmed		
B*15:323	Unconfirmed	B*15:353	Unconfirmed	B*15:383	Unconfirmed		
B*15:324	Unconfirmed	B*15:354	Unconfirmed	B*15:384	Unconfirmed		
B*15:325	Unconfirmed	B*15:355	Unconfirmed	B*15:385	Unconfirmed		
B*15:326	Unconfirmed	B*15:356	Unconfirmed	B*15:386	Unconfirmed		
B*15:327	Unconfirmed	B*15:357	Unconfirmed	B*15:387	Unconfirmed		
B*15:328	Unconfirmed	B*15:358	Unconfirmed	B*15:388	Unconfirmed		
B*15:329	Unconfirmed	B*15:359	Unconfirmed	B*15:389	Unconfirmed		

¹Allele status “confirmed” or “unconfirmed” as listed on the IMGT/HLA web page 2016-July-14, release 3.25.0, www.ebi.ac.uk/imgt/hla.

101.516-24/04 – including *Taq* polymerase, IFU-01
101.516-24u/04u – without *Taq* polymerase, IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **3F7**

Lot-specific information
SPECIFICITY TABLE

HLA-B*15 SSP subtyping

Specificities and sizes of the PCR products of the 95+1 primer mixes used for
HLA-B*15 SSP subtyping

Primer Mix	Size of spec. PCR product ¹	Size of control band ²	Amplified HLA-B*15 alleles ³	Other amplified Class I alleles
1	400 bp	800 bp	*15:01:01:01-15:01:04, 15:01:06-15:07:03, 15:12-15:15, 15:19-15:20, 15:24:01-15:28, 15:30-15:36, 15:38:01-15:40:02, 15:42-15:43, 15:45-15:50, 15:53-15:55, 15:57-15:58, 15:60-15:65, 15:69-15:71, 15:73, 15:75, 15:77-15:79N, 15:81-15:89, 15:91-15:92, 15:94N, 15:96, 15:98, 15:102-15:107, 15:109-15:113, 15:116-15:118, 15:120-15:123, 15:125-15:129, 15:131-15:132, 15:135-15:138, 15:140-15:142, 15:144-15:147, 15:150-15:151, 15:154-15:160, 15:163-15:167, 15:169-15:175, 15:178-15:179:02, 15:181N-15:185, 15:187-15:188, 15:190N, 15:192-15:195, 15:199, 15:201-15:207, 15:210-15:214, 15:217-15:220, 15:223-15:225, 15:227-15:228, 15:231-15:235, 15:237, 15:239-15:242:02, 15:244-15:251, 15:253, 15:256-15:262N, 15:264-15:267, 15:269-15:272N, 15:274, 15:276-15:283, 15:285-15:289, 15:291, 15:295-15:302N, 15:304N, 15:308-15:310, 15:315-15:322, 15:324-15:328, 15:330-15:334, 15:336, 15:339-15:350, 15:352-15:353, 15:355, 15:357-15:360, 15:363-15:366, 15:368-15:375N, 15:377Q-15:379, 15:381, 15:383-15:387, 15:389-15:393, 15:395, 15:397-15:398	*40:12, C*16:85
2	290 bp	1070 bp	*15:01:01:01-15:01:03, 15:01:06-15:02:11, 15:03:03, 15:04:01:01-15:08:02, 15:10:02, 15:11:01-15:11:03, 15:11:05-15:15, 15:17:01:01-15:17:05, 15:19-15:21, 15:24:01-15:28, 15:30-15:36, 15:38:01-15:40:02, 15:42-15:46, 15:48, 15:50, 15:55-15:58, 15:60, 15:63, 15:65-15:66, 15:70-15:71, 15:73, 15:75-15:79N, 15:81-15:89, 15:92, 15:94N, 15:96-15:97, 15:101-15:102, 15:104-15:107, 15:109-15:113, 15:116-15:118, 15:120-15:122, 15:125-15:126, 15:128-15:129, 15:135-15:150, 15:152, 15:154-15:155, 15:157, 15:159-15:160, 15:162-15:172, 15:174-15:175, 15:177-15:179:02, 15:181N-15:185,	*46:01:01-46:68, C*01:73, C*05:01:28, C*07:02:30, C*08:16:02

101.516-24/04 – including *Taq* polymerase, IFU-01
 101.516-24u/04u – without *Taq* polymerase, IFU-02

Visit www.olerup.com for
 “Instructions for Use” (IFU)

Lot No.: **3F7**

Lot-specific information

			15:187-15:196, 15:199, 15:201-15:209N, 15:211, 15:213-15:219, 15:223-15:225, 15:227-15:228, 15:230-15:234, 15:236- 15:237, 15:239-15:241, 15:244-15:251, 15:256-15:262N, 15:264-15:265, 15:267- 15:273, 15:276-15:280, 15:283-15:289, 15:291, 15:295-15:305, 15:308-15:310, 15:315-15:322, 15:324-15:328, 15:330- 15:334, 15:339-15:346, 15:348-15:350, 15:352-15:353, 15:355-15:361, 15:363- 15:368, 15:370-15:375N, 15:377Q- 15:378, 15:381, 15:383-15:387, 15:390- 15:396, 15:398	
3⁴	120 bp	1070 bp	*15:01:01:01-15:01:04, 15:01:06- 15:01:16, 15:01:17 ^w , 15:01:18-15:01:30, 15:01:32-15:01:39, 15:04:01:01- 15:04:04, 15:07:01-15:08:02, 15:11:01- 15:12, 15:14-15:15, 15:19, 15:24:01- 15:24:02, 15:26N-15:28, 15:30, 15:32:01-15:32:02, 15:34-15:35, 15:38:01-15:38:02, 15:43, 15:45-15:46, 15:50, 15:53-15:54, 15:56-15:58, 15:60, 15:63, 15:66, 15:68, 15:70-15:71, 15:73, 15:75-15:77, 15:79N, 15:81-15:82, 15:85, 15:87, 15:90, 15:92, 15:94N, 15:96-15:97, 15:101-15:102, 15:104- 15:105, 15:109-15:111N, 15:113, 15:117-15:118, 15:120, 15:122, 15:125- 15:126, 15:128-15:129, 15:135, 15:137, 15:140, 15:142-15:149N, 15:152, 15:154, 15:157, 15:159-15:160, 15:163- 15:167, 15:169, 15:171-15:172, 15:174- 15:175, 15:178, 15:180-15:184, 15:187, 15:189-15:193, 15:201-15:203, 15:205- 15:207, 15:209N, 15:211-15:212, 15:215, 15:217, 15:225, 15:227-15:228, 15:231-15:234, 15:236, 15:239, 15:241, 15:244-15:247, 15:251, 15:256-15:262N, 15:264, 15:267, 15:269-15:270, 15:272N, 15:276-15:280, 15:284-15:286, 15:296, 15:298-15:299, 15:303-15:305, 15:309, 15:315-15:318, 15:320-15:321Q, 15:324, 15:326-15:327, 15:331, 15:333- 15:334, 15:336, 15:340, 15:342-15:344, 15:346-15:347, 15:349-15:350, 15:353, 15:355, 15:359, 15:363, 15:365, 15:367- 15:368, 15:370-15:373, 15:375N, 15:377Q, 15:381, 15:385-15:387, 15:390-15:392, 15:398	*07:100, 13:31, 13:41, 18:19, 27:25, 27:115, 35:14:01-35:14:02, 35:43:01- 35:44, 35:62, 35:67, 35:79, 35:86, 35:102, 35:117, 35:135, 35:163, 35:185, 35:213, 35:265, 39:36, 44:146, 45:05, 46:01:01-46:02, 46:04- 46:05, 46:07N, 46:09-46:10, 46:12, 46:14-46:17, 46:20, 46:22-46:24, 46:26-46:32, 46:34-46:42, 46:44, 46:46-46:53, 46:55N-46:57, 46:60- 46:61, 46:63-46:68, 50:47, 51:61:01- 51:61:02, 52:21, 54:06, 54:33, 55:21, 56:03, 56:43, 58:78, A*26:68, A*26:100, A*68:56, C*02:56, C*06:20, C*12:50
4	180 bp 390 bp	800 bp	*15:01:01:02N *15:285, 15:378	
5⁶	385 bp	1070 bp	*15:02:01-15:02:11, 15:08:01-15:11:07, 15:13:01-15:13:02, 15:15, 15:18:01:01- 15:18:06, 15:21, 15:23, 15:29, 15:31, 15:37, 15:44, 15:51-15:52, 15:55, 15:64:01-15:64:02, 15:72, 15:76, 15:80,	C*06:147

101.516-24/04 – including *Taq* polymerase, IFU-01
101.516-24u/04u – without *Taq* polymerase, IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **3F7**

Lot-specific information

			15:88-15:90, 15:93, 15:99, 15:108, 15:112, 15:114-15:115, 15:119, 15:121, 15:124, 15:133-15:134, 15:139, 15:143- 15:144, 15:148, 15:153, 15:161, 15:170, 15:176, 15:186, 15:189, 15:191, 15:194, 15:197:01-15:198, 15:200, 15:209N, 15:213-15:215, 15:221, 15:223, 15:226N, 15:229, 15:238, 15:252, 15:263, 15:265, 15:275:01-15:275:02, 15:283-15:284, 15:290-15:294N, 15:297, 15:301-15:303, 15:305-15:308, 15:311- 15:314, 15:319, 15:323, 15:325, 15:328- 15:330, 15:335, 15:337-15:338, 15:341, 15:345, 15:350-15:351, 15:354, 15:357- 15:358, 15:367, 15:374, 15:378, 15:380N, 15:382, 15:384, 15:388, 15:394	
6	525 bp	1070 bp	*15:02:01-15:02:11, 15:13:01-15:13:02, 15:21, 15:25:01-15:25:03, 15:36, 15:44, 15:62, 15:77, 15:80, 15:85, 15:88-15:89, 15:106, 15:112, 15:121, 15:139, 15:144, 15:154, 15:165, 15:170, 15:194, 15:204, 15:213-15:214, 15:223, 15:240, 15:250, 15:265, 15:271, 15:283, 15:289, 15:291, 15:297, 15:301-15:302N, 15:308, 15:319, 15:325, 15:328, 15:330, 15:341, 15:345, 15:357-15:358, 15:366, 15:374, 15:378, 15:393-15:394	*46:33, 46:67
7	330 bp	1070 bp	*15:03:01:01-15:03:06, 15:09-15:10:04, 15:18:01:01-15:18:06, 15:23, 15:29, 15:37, 15:46-15:47:02, 15:49, 15:51- 15:54, 15:61-15:62, 15:64:01-15:64:02, 15:69, 15:72, 15:74, 15:80, 15:90-15:91, 15:93, 15:98-15:99, 15:103, 15:106, 15:108, 15:114-15:115, 15:119, 15:123- 15:124, 15:127, 15:131-15:134, 15:143, 15:151, 15:153, 15:156, 15:158, 15:161, 15:173, 15:176, 15:186, 15:197:01- 15:198, 15:200, 15:210, 15:212, 15:220- 15:221, 15:226N, 15:229, 15:235, 15:238, 15:242:01-15:243, 15:252- 15:253, 15:255, 15:259, 15:263, 15:266, 15:274-15:275:02, 15:281-15:282, 15:290, 15:292-15:294N, 15:306-15:307, 15:311-15:314, 15:323, 15:329, 15:335, 15:337-15:338, 15:351, 15:354, 15:369, 15:379-15:380N, 15:382, 15:388-15:389, 15:397	*40:12, C*07:46, C*07:356, C*12:184, C*15:99, C*16:85
8 ⁴	105 bp	1070 bp	*15:04:01:01-15:04:04, 15:16:01:01- 15:16:03, 15:67, 15:95, 15:137, 15:155, 15:222, 15:254, 15:293, 15:310, 15:362	*13:04, 13:46, 13:72, 35:37, 35:235, 35:265, 35:276, 39:06:01-39:06:05, 39:33-39:34, 39:50, 39:57, 39:62, 39:64, 39:83, 39:90, 39:102, 39:109, 40:86, 40:159, 40:330, 41:30, 45:02, 45:09, 46:32, 49:04:01-49:05, 50:09, 50:33, 51:13:01-51:13:02, 51:37, 51:63, 51:92:01-51:92:02, 51:97, 51:172, 52:14, 52:39, 54:14, 55:13,

101.516-24/04 – including *Taq* polymerase, IFU-01
101.516-24u/04u – without *Taq* polymerase, IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **3F7**

Lot-specific information

				55:23, 55:32, 55:58, 55:68, 56:22, 56:50, 59:03, 73:01-73:02, C*07:463, C*12:183
9	365 bp	1070 bp	*15:05:01-15:05:02, 15:31, 15:48, 15:51 ^w , 15:52, 15:55, 15:84, 15:86, 15:91, 15:107, 15:114, 15:123-15:124, 15:136, 15:151, 15:155, 15:179:01 ^w -15:179:02 ^w , 15:185-15:186, 15:188, 15:218Q ^w , 15:222, 15:224, 15:235, 15:237, 15:275:01-15:275:02, 15:283, 15:332, 15:348, 15:383, 15:393, 15:396	*46:06, 46:09, 46:11, 46:13:01-46:13:03, 46:18-46:19, 46:21:01-46:21:02, 46:25, 46:33, 46:43
10	300 bp 515 bp	1070 bp	*15:37-15:38:02, 15:185, 15:323, 15:335, 15:364, 15:368 *15:06, 15:27:01-15:27:03, 15:84, 15:109, 15:112, 15:195, 15:327, 15:344, 15:398	
11	140 bp	800 bp	*15:07:01-15:07:03, 15:45, 15:55, 15:68, 15:126, 15:197:02, 15:207, 15:324, 15:329, 15:331, 15:360, 15:383	*07:68:01-07:68:03, 07:100, 07:237, 07:262, 08:71, 14:05, 14:13, 14:53, 18:14, 27:32, 27:138, 35:05:01-35:05:03, 35:22, 35:51, 35:58, 35:66, 35:72, 35:89, 35:97, 35:199, 35:232, 38:19, 39:03, 39:14, 39:24:01-39:24:02, 39:29, 39:37, 39:76, 40:35:01-40:35:02, 40:71, 41:24, 42:10, 42:17, 44:54, 44:106, 44:135, 44:158, 44:184, 44:213, 46:12, 48:12, 48:14, 48:29, 51:64, 53:14, 58:18
12	400 bp	1070 bp	*15:08:01-15:08:02, 15:29, 15:49, 15:56, 15:191, 15:236	C*06:147
13	395 bp	800 bp	*15:09-15:10:04, 15:18:01:01-15:18:06, 15:21, 15:23, 15:37, 15:44, 15:51-15:52, 15:66, 15:72, 15:80, 15:90, 15:93, 15:99, 15:108, 15:114-15:115, 15:119, 15:124, 15:133-15:134, 15:153, 15:161, 15:176, 15:186, 15:189, 15:197:01-15:198, 15:200, 15:221, 15:226N, 15:229, 15:238, 15:243, 15:252, 15:263, 15:275:01-15:275:02, 15:290, 15:292-15:294N, 15:306-15:307, 15:311-15:314, 15:323, 15:329, 15:335, 15:337-15:338, 15:351, 15:354, 15:376, 15:380N, 15:382, 15:388, 15:394 *15:79N	
14	435 bp 465 bp 515 bp	1070 bp	*15:09-15:10:04, 15:30, 15:37, 15:45, 15:48, 15:63, 15:83, 15:90, 15:99, 15:150, 15:188, 15:243, 15:248, 15:252, 15:287, 15:312-15:313, 15:323-15:324, 15:329, 15:338, 15:361 *15:110	
15^{4,6}	105 bp	1070 bp	*15:01:02, 15:09, 15:243, 15:252, 15:277-15:278, 15:361	*07:02:11, 07:05:03, 35:251, 40:01:06, 40:26, 51:06:01, 51:06:03, 51:59, 51:64, 51:135, 51:194, C*02:02:23
16⁴	120 bp	1070 bp	*15:116, 15:124, 15:352	*40:52, 40:59, 46:06, 46:43, 48:23, 57:59

101.516-24/04 – including *Taq* polymerase, IFU-01
101.516-24u/04u – without *Taq* polymerase, IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **3F7**

Lot-specific information

	190 bp		*15:12, 15:19, 15:270, 15:298, 15:304N	
17	435 bp	1070 bp	*15:13:01-15:13:02, 15:16:01:01-15:17:05, 15:23-15:24:02, 15:67, 15:87, 15:95, 15:157, 15:162, 15:168, 15:177, 15:196, 15:208, 15:216, 15:222, 15:230, 15:254, 15:268, 15:273, 15:356, 15:361-15:362, 15:396	
18	195 bp	1070 bp	*15:14, 15:91, 15:131, 15:161	*18:56, 35:45, 35:71, 44:17, 44:43:01-44:43:02, 44:144, 44:230, 45:09, 46:17, 53:22, 58:07, C*02:82 , C*06:165 , C*14:61
19	395 bp	800 bp	*15:16:01:01-15:17:05, 15:67, 15:74, 15:95, 15:162, 15:168, 15:177, 15:196, 15:208, 15:216, 15:222, 15:230, 15:254, 15:268, 15:273, 15:356, 15:361-15:362, 15:396	
20	360 bp 505 bp	1070 bp	*15:105 *15:16:01:01-15:16:03, 15:34, 15:62, 15:67, 15:85-15:86, 15:95, 15:137, 15:194, 15:222, 15:242:02, 15:254, 15:303, 15:362, 15:393	*46:11, 46:18, 46:33, 46:53, 46:67
21	205 bp 425 bp	1070 bp	*15:16:01:01-15:16:03, 15:67, 15:95, 15:222, 15:254, 15:362 *15:43	C*03:158
22	190 bp	800 bp	*15:01:01:01-15:01:04, 15:01:06-15:01:39, 15:03:01:01-15:05:02, 15:07:01-15:11:07, 15:12 ^w , 15:14-15:16:03, 15:18:01:01-15:18:06, 15:19 ^w , 15:20, 15:23-15:24:02, 15:26N-15:30, 15:32:01-15:35, 15:37-15:39:02, 15:43, 15:45-15:46, 15:48, 15:51, 15:53-15:54, 15:56-15:58, 15:60-15:68, 15:70-15:76, 15:78:01-15:82, 15:84-15:85, 15:87, 15:90-15:92, 15:94N-15:99, 15:101-15:105, 15:107-15:111N, 15:113, 15:115-15:116, 15:118-15:120, 15:122-15:123, 15:125-15:127, 15:129, 15:131-15:136, 15:140-15:143, 15:145-15:155, 15:157-15:161, 15:163-15:167, 15:169, 15:171-15:176, 15:178-15:185, 15:187, 15:189-15:195, 15:197:01-15:198, 15:200-15:203, 15:205-15:206, 15:209N-15:210, 15:212, 15:215, 15:217, 15:219-15:223, 15:225-15:229, 15:231-15:232, 15:234-15:237, 15:239, 15:242:01-15:244, 15:246N-15:249, 15:251, 15:253-15:264, 15:266-15:269, 15:270 ^w , 15:272N, 15:274-15:282, 15:284-15:288, 15:290, 15:292-15:296, 15:298 ^w , 15:299-15:300, 15:303, 15:304N ^w , 15:305-15:307, 15:309-15:318, 15:320-15:322, 15:324, 15:326-15:327, 15:329, 15:331-15:338, 15:340, 15:343-15:344, 15:346-15:347, 15:350-15:355, 15:359,	*07:20, 07:24, 07:60, 07:100, 07:131, 07:220, 07:223, 07:235, 08:21, 08:25, 13:03-13:04, 13:48, 13:62, 13:86, 18:36, 18:56-18:58, 18:101, 27:153, 35:01:01:01-35:01:35, 35:01:37-35:02:07, 35:03:01-35:03:05, 35:03:07-35:14:02, 35:16-35:30, 35:32:01-35:32:02, 35:34, 35:36-35:59:02, 35:61:01-35:64:02, 35:67-35:72, 35:74-35:85, 35:89-35:92, 35:94-35:95, 35:97-35:130N, 35:132-35:156, 35:158-35:163, 35:165N-35:187, 35:189, 35:191-35:212, 35:214, 35:216N-35:229, 35:231-35:233, 35:235-35:260, 35:263-35:275, 35:278-35:321, 38:30, 39:32-39:33, 39:74, 40:05, 40:26, 40:28, 40:63, 40:71, 40:92, 40:174, 40:324, 41:17, 44:10, 44:14, 44:17, 44:43:01-44:43:02, 44:45, 44:51, 44:62, 44:77, 44:82, 44:107, 44:123, 44:135, 44:184, 44:213, 44:230, 45:09, 45:14, 46:01:01-46:02, 46:04-46:05, 46:07N-46:17, 46:19-46:20, 46:22-46:28, 46:31-46:50, 46:52-46:68, 48:02:01-48:02:03, 48:25, 49:04:01-49:04:02, 51:01:01:01-51:01:31, 51:01:33-51:01:38, 51:01:40-51:01:49, 51:01:51-51:09:03, 51:11N-51:15, 51:19-51:20, 51:22-51:24:05,

101.516-24/04 – including *Taq* polymerase, IFU-01
101.516-24u/04u – without *Taq* polymerase, IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **3F7**

Lot-specific information

			15:361-15:365, 15:367, 15:369-15:373, 15:375N-15:377Q, 15:379-15:383, 15:385-15:388, 15:390-15:393, 15:397-15:398	51:26-51:30, 51:32-51:33, 51:35, 51:37-51:41N, 51:43-51:46, 51:48-51:80, 51:83-51:84, 51:86-51:92:02, 51:94-51:100, 51:102-51:142, 51:145-51:152, 51:154-51:156, 51:158:01-51:171, 51:173Q-51:196, 51:198-51:208, 52:01:01:01-52:01:10, 52:01:12-52:01:20, 52:01:22-52:07, 52:09-52:25:02, 52:27-52:36, 52:38-52:61, 53:01:01-53:01:06, 53:01:08-53:06, 53:08:01-53:43, 54:03, 54:06, 54:33, 55:08, 55:51, 56:01:01:01-56:09, 56:11, 56:13-56:17, 56:20:01-56:22, 56:24-56:30, 56:33-56:51, 57:01:01-57:01:02, 57:01:04-57:03:02, 57:06-57:12, 57:14:01-57:24, 57:27-57:31, 57:33-57:40, 57:42, 57:44-57:83, 58:01:01:01-58:02:02, 58:04, 58:06-58:19, 58:21-58:26, 58:28:01-58:29, 58:31N-58:34, 58:36, 58:38-58:75, 58:77-58:81, 59:04, 78:01:01-78:06, 78:09, 82:01-82:03, C*03:81, C*03:175, C*03:199, C*03:245, C*03:317, C*15:15
23⁸	230 bp	800 bp	*15:19, 15:304N	*46:24
24	205 bp 305 bp	1070 bp	*15:26N *15:33, 15:72, 15:78:01-15:78:03, 15:107, 15:116, 15:141, 15:150, 15:188, 15:209N, 15:248, 15:287, 15:352	*46:06, 46:45
25	315 bp 385 bp	1070 bp	*15:46, 15:53, 15:106, 15:212 *15:28	C*16:85
26	230 bp	1070 bp	*15:32:01-15:32:02, 15:299	C*03:178, C*07:186, C*07:242, C*07:265
27⁴	115 bp	1070 bp	*15:33, 15:248	*49:22, 51:126, 57:54, C*03:87:01-C*03:87:02, C*05:27, C*05:39, C*08:115
28	240 bp	1070 bp	*15:35, 15:141, 15:290, 15:381	
29	255 bp	1070 bp	*15:37-15:38:02, 15:185, 15:323, 15:335, 15:364, 15:368	*07:218, 07:220, 08:118, 14:05-14:06:01, 14:08:02, 14:13, 14:53, 18:01:01:01-18:01:05, 18:01:07-18:01:24, 18:03-18:07:02, 18:09, 18:12:01-18:15, 18:17N, 18:19-18:20, 18:23N-18:26, 18:28-18:34, 18:36-18:44:01, 18:45-18:51, 18:53-18:55, 18:57, 18:59-18:100, 18:102-18:104, 18:106-18:109, 18:111-18:123, 18:125-18:127, 39:32, 39:43, 39:48, 40:51, 40:113, 40:163, 40:174, 51:06:01, 51:06:03, 51:45, 51:62, 51:64, 51:135, 51:194, 56:31, 67:06
30⁶	320 bp	1070 bp	*15:40:01-15:40:02, 15:47:01-15:47:02, 15:49, 15:52, 15:114, 15:117, 15:124, 15:138, 15:238, 15:241, 15:389	*46:06, 46:25, 46:30

101.516-24/04 – including *Taq* polymerase, IFU-01
101.516-24u/04u – without *Taq* polymerase, IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **3F7**

Lot-specific information

	355 bp		*15:75, 15:312	
	405 bp		*15:375N	
31	320 bp	1070 bp	*15:44, 15:50, 15:69, 15:83, 15:86, 15:93, 15:121, 15:186, 15:188, 15:199, 15:218Q, 15:224, 15:252, 15:265, 15:323, 15:348, 15:368	*46:03, 46:18, 46:21:01-46:21:02, 46:29
	395 bp		*15:114, 15:153	
32	275 bp	1070 bp	*15:196	
	350 bp		*15:190N	
	380 bp		*15:48, 15:108, 15:136, 15:235, 15:249	*46:19
33	200 bp	1070 bp	*15:118	
	255 bp		*15:18:03, 15:42, 15:73, 15:86, 15:224, 15:303	*46:11, 46:18
	340 bp		*15:192, 15:281	
34	160 bp	1070 bp	*15:51, 15:179:01-15:179:02, 15:199, 15:218Q	*08:28, 08:35, 08:37, 08:89, 08:107, 18:56, 18:113, 35:115, 35:260, 37:35, 37:56, 41:24, 41:30, 42:09, 44:17, 44:230, 45:09, 51:97, 53:22, 57:24, C*07:68, C*07:140, C*07:151, C*07:302
35	165 bp	1070 bp	*15:30, 15:58, 15:73, 15:150, 15:303, 15:324	*08:70, 13:31, 13:41, 35:44, 39:36, 40:324, 46:61, 48:19, 51:61:01-51:61:02, 52:21, 54:33, 55:21, 56:43, C*02:68
36 ⁶	360 bp	800 bp	*15:82, 15:260	C*03:186:01
	435 bp		*15:36, 15:89, 15:115, 15:256, 15:339	C*16:85
37	170 bp	1070 bp	*15:47:01-15:47:02, 15:49, 15:65, 15:84, 15:99, 15:221, 15:237, 15:379	*07:02:01-07:02:17, 07:02:19-07:02:25, 07:02:27-07:02:45, 07:02:47-07:03, 07:05:01:01-07:13, 07:15-07:16, 07:18:01-07:18:02, 07:20-07:24, 07:27, 07:29-07:33:03, 07:35-07:42, 07:44N, 07:45 ^w , 07:46-07:47, 07:49N-07:51, 07:52 ^w , 07:53-07:54, 07:56:01-07:65, 07:67N-07:79, 07:81-07:95, 07:97-07:99, 07:101-07:114, 07:116-07:130, 07:132-07:145, 07:148-07:163, 07:165-07:171, 07:173-07:177, 07:179-07:196, 07:198-07:208, 07:210-07:222, 07:224, 07:227-07:234, 07:236-07:259, 07:261-07:270, 07:272N-07:275, 08:20, 08:40, 08:79, 08:156, 13:18, 13:73, 18:13, 18:21, 18:30, 35:08:01-35:08:08, 35:18, 35:45, 35:61:01-35:61:02, 35:66, 35:80, 35:99-35:100, 35:105, 35:142, 35:156, 35:158, 35:172, 35:174, 35:176, 35:187, 35:192, 35:206, 35:225, 35:229, 35:232, 35:262, 35:300-35:301, 37:07, 39:08 ^w , 39:11 ^w , 39:33, 40:15-40:16, 40:23, 40:32, 40:98, 40:128, 40:161, 40:198, 46:03, 46:62, 48:05, 48:08, 51:05, 51:29, 51:54, 51:82, 52:24, 53:11, 53:16,

101.516-24/04 – including *Taq* polymerase, IFU-01
101.516-24u/04u – without *Taq* polymerase, IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **3F7**

Lot-specific information

				53:24, 54:10, 54:20, 55:09, 55:37, 55:52, 56:51, 57:02:01-57:02:02, 57:12-57:13, 57:19, 57:28N, 57:30, 57:42, 57:63, 58:22, C*01:90, C*01:101-01:102, C*01:113, C*01:116, C*04:108, C*05:132, C*15:08
38	225 bp	1070 bp	*15:61, 15:70, 15:76, 15:101, 15:180 ^w , 15:189, 15:236, 15:255, 15:397 ^w	*07:29 ^w , 07:65, 07:72, 07:86, 08:01:01:01 ^w -08:01:34 ^w , 08:04 ^w -08:05 ^w , 08:07 ^w -08:22 ^w , 08:24 ^w -08:53:02 ^w , 08:55 ^w -08:65 ^w , 08:67N ^w -08:71 ^w , 08:73 ^w -08:93 ^w , 08:95 ^w -08:109 ^w , 08:111 ^w -08:116 ^w , 08:118 ^w -08:125 ^w , 08:127 ^w -08:130 ^w , 08:132 ^w , 08:134 ^w -08:138 ^w , 08:140 ^w -08:150 ^w , 08:152 ^w -08:157 ^w , 08:160 ^w -08:165 ^w , 14:11, 14:28 ^w , 27:106 ^w , 35:29:01-35:29:02, 35:69, 35:76, 35:136, 35:186 ^w , 35:198, 35:229, 35:254, 35:294, 35:316, 39:04, 39:63, 40:43 ^w , 40:166 ^w , 41:27 ^w , 44:90, 45:06 ^w , 46:01:01 ^w , 46:01:03 ^w -46:01:22 ^w , 46:03 ^w -46:04 ^w , 46:06 ^w -46:34 ^w , 46:36 ^w -46:49 ^w , 46:51Q ^w -46:68 ^w , 51:20 ^w , 51:208 ^w , 53:23 ^w , 54:01:01-54:01:03, 54:01:05-54:01:07, 54:03-54:11, 54:12 ^w , 54:13-54:22, 54:24-54:29, 54:31-54:35, 55:01:01-55:01:07, 55:01:09-55:02:10, 55:03 ^w , 55:04-55:05, 55:07-55:09, 55:11, 55:12 ^w , 55:13-55:15, 55:16 ^w , 55:17, 55:19-55:24, 55:25 ^w , 55:26-55:30, 55:32-55:48, 55:49 ^w , 55:50-55:68, 55:69 ^w , 55:70-55:75, 55:77-55:78, 56:01:01:01-56:06, 56:08-56:20:02, 56:21 ^w , 56:22, 56:24-56:34, 56:36-56:40, 56:41 ^w , 56:43, 56:45-56:51, 78:01:01-78:01:02, 78:03, 78:07-78:08, 83:01, C*07:441^w, C*16:86^w
39	315 bp 385 bp	800 bp	*15:60, 15:119 *15:95, 15:285, 15:378	C*07:213
40⁵	220 bp 425 bp	1070 bp	*15:71, 15:175, 15:381 *15:13:01-15:13:02, 15:16:01:01-15:17:05, 15:23-15:24:02, 15:67, 15:95, 15:157, 15:162, 15:168, 15:177, 15:196, 15:208, 15:216, 15:222, 15:230, 15:254, 15:268, 15:273, 15:356, 15:361-15:362, 15:396	C*06:77
41	465 bp	1070 bp	*15:01:01:01-15:01:04, 15:01:06-15:08:02, 15:11:01-15:16:03, 15:18:01:01-15:19, 15:21, 15:23-15:29, 15:31-15:36, 15:38:01-15:40:02, 15:43-15:44, 15:46-15:47:02, 15:49-15:57, 15:60-15:62, 15:64:01-15:72,	*35:43:01, 35:185, 46:01:01-46:10, 46:12-46:17, 46:19-46:60, 46:62-46:68, 56:03

101.516-24/04 – including *Taq* polymerase, IFU-01
101.516-24u/04u – without *Taq* polymerase, IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **3F7**

Lot-specific information

			15:74-15:76, 15:78:01-15:82, 15:84-15:85, 15:87-15:89, 15:91-15:98, 15:101-15:129, 15:131-15:132, 15:134-15:136, 15:138-15:149N, 15:151-15:161, 15:163-15:167, 15:169-15:176, 15:178-15:187, 15:189-15:195, 15:197:01-15:207, 15:209N-15:215, 15:217-15:223, 15:225-15:227, 15:228?, 15:231-15:232, 15:234-15:242:02, 15:244-15:247, 15:249-15:251, 15:254-15:268, 15:270-15:272N, 15:274-15:286, 15:288-15:302N, 15:304N-15:311, 15:314-15:322, 15:325-15:328, 15:330-15:333, 15:335-15:337, 15:340-15:341, 15:343-15:355, 15:357-15:360, 15:362-15:382, 15:384-15:395, 15:397-15:398	
42	350 bp	1070 bp	*15:02:01-15:03:06, 15:05:01-15:06, 15:09-15:10:04, 15:13:01-15:13:02, 15:16:01:01-15:18:06, 15:21, 15:23, 15:25:01-15:25:03, 15:29, 15:31, 15:36-15:37, 15:39:01-15:40:02, 15:44, 15:48, 15:52, 15:55, 15:61-15:62, 15:64:01-15:64:02, 15:67, 15:69, 15:72, 15:74, 15:80, 15:86, 15:88-15:91, 15:93, 15:95, 15:98, 15:103, 15:106-15:108, 15:112, 15:114-15:115, 15:119, 15:121, 15:123-15:124, 15:127, 15:131-15:134, 15:136, 15:138-15:139, 15:151, 15:153, 15:155-15:156, 15:158, 15:161-15:162, 15:170, 15:173, 15:176-15:177, 15:185-15:186, 15:188, 15:194-15:198, 15:200, 15:204, 15:208, 15:210, 15:213-15:214, 15:216, 15:219-15:220, 15:222-15:224, 15:226N, 15:229-15:230, 15:235, 15:238, 15:240, 15:242:01-15:243, 15:250, 15:252-15:255, 15:263, 15:265-15:266, 15:268, 15:271, 15:273-15:275:02, 15:281-15:283, 15:288-15:294N, 15:297, 15:301-15:302N, 15:306-15:308, 15:310-15:314, 15:319, 15:323, 15:325, 15:328-15:330, 15:332, 15:335, 15:337-15:339, 15:341, 15:345, 15:348, 15:351, 15:354, 15:356-15:358, 15:360-15:362, 15:364, 15:366, 15:369, 15:374, 15:376, 15:378, 15:380N, 15:382, 15:384, 15:388-15:389, 15:393-15:397	*46:06, 46:08, 46:11, 46:13:01-46:13:03, 46:18-46:19, 46:21:01-46:21:02, 46:25-46:26, 46:33, 46:43
	460 bp		*15:73, 15:86, 15:224, 15:303	*46:11, 46:18
43	525 bp	1070 bp	*15:04:01:01-15:04:04, 15:16:01:01-15:16:03, 15:67, 15:83, 15:95, 15:137, 15:155, 15:222, 15:254, 15:293, 15:303, 15:310, 15:362	*46:11, 46:18, 46:32
44	515 bp	1070 bp	*15:01:01:01-15:01:04, 15:01:06-15:01:13, 15:01:15-15:01:39, 15:03:01:01-15:03:06, 15:05:01-15:05:02, 15:07:01-15:11:02, 15:11:04-15:12, 15:14-15:15,	*35:43:01, 35:185, 46:01:01-46:01:18, 46:01:20-46:10, 46:12-46:17, 46:19-46:31, 46:34-46:66, 46:68, 56:03

101.516-24/04 – including *Taq* polymerase, IFU-01
 101.516-24u/04u – without *Taq* polymerase, IFU-02

Visit www.olerup.com for
 “Instructions for Use” (IFU)

Lot No.: **3F7**

Lot-specific information

			15:17:01:01-15:18:01:02, 15:18:03-15:19, 15:23-15:24:02, 15:26N ^w , 15:28-15:31, 15:33-15:35, 15:37-15:40:02, 15:43, 15:45-15:58, 15:60-15:61, 15:63-15:66, 15:68-15:76, 15:78:01-15:79N, 15:81-15:82, 15:87, 15:90-15:94N, 15:96-15:99, 15:101-15:105, 15:107-15:108, 15:111N, 15:113-15:120, 15:122-15:128, 15:131-15:136, 15:138, 15:140-15:143, 15:145-15:153, 15:156-15:169, 15:171-15:183, 15:185-15:193, 15:196-15:203, 15:205-15:212, 15:215-15:221, 15:224-15:227, 15:228 [?] , 15:229-15:239, 15:241-15:242:01, 15:243-15:249, 15:251-15:253, 15:255-15:264, 15:266-15:270, 15:272N-15:282, 15:284-15:288, 15:290, 15:292, 15:294N-15:296, 15:298, 15:300, 15:304N-15:307, 15:309, 15:311-15:318, 15:320-15:324, 15:326, 15:329, 15:331-15:340, 15:343, 15:346-15:356, 15:359-15:361, 15:363-15:365, 15:367-15:373, 15:375N-15:377Q, 15:379-15:382, 15:384-15:392, 15:396-15:397	
45	475 bp	1070 bp	*15:01:01:01-15:01:04, 15:01:06-15:01:39, 15:03:01:01-15:05:02, 15:07:01-15:12, 15:14-15:16:03, 15:18:01:01-15:19, 15:23-15:24:02, 15:26N-15:30, 15:32:01-15:35, 15:37-15:40:02, 15:43, 15:45-15:54, 15:56-15:58, 15:60-15:76, 15:78:01-15:87, 15:90-15:99, 15:101-15:105, 15:107-15:111N, 15:113-15:120, 15:122-15:127, 15:129, 15:131-15:136, 15:138, 15:140-15:143, 15:145-15:155, 15:157-15:161, 15:163-15:167, 15:169, 15:171-15:176, 15:178-15:195, 15:197:01-15:203, 15:205-15:206, 15:209N-15:210, 15:212, 15:215, 15:217-15:227, 15:228 [?] , 15:229, 15:231-15:232, 15:234-15:239, 15:241-15:249, 15:251-15:264, 15:266-15:270, 15:272N, 15:274-15:282, 15:284-15:288, 15:290, 15:292-15:296, 15:298-15:300, 15:303-15:307, 15:309-15:318, 15:320-15:324, 15:326-15:327, 15:329, 15:331-15:338, 15:340, 15:343-15:344, 15:346-15:348, 15:350-15:355, 15:359, 15:361-15:365, 15:367-15:373, 15:375N-15:377Q, 15:379-15:382, 15:385-15:393, 15:397-15:398	*35:43:01, 35:185, 46:01:01-46:68, 56:03
46	360 bp	800 bp	*15:01:01:01-15:01:04, 15:01:06-15:04:04, 15:06-15:19, 15:21, 15:23-15:30, 15:32:01-15:40:02, 15:43-15:47:02, 15:49-15:50, 15:53-15:54, 15:56-15:58, 15:60-15:74, 15:76-15:82, 15:85, 15:87, 15:89-15:90, 15:92-15:99,	*35:43:01, 35:185, 46:01:01-46:05, 46:07N-46:08, 46:10, 46:12, 46:14-46:17, 46:20, 46:22-46:24, 46:26-46:32, 46:34-46:42, 46:44-46:68, 56:03

101.516-24/04 – including *Taq* polymerase, IFU-01
 101.516-24u/04u – without *Taq* polymerase, IFU-02

Visit www.olerup.com for
 “Instructions for Use” (IFU)

Lot No.: **3F7**

Lot-specific information

			15:101-15:104, 15:106, 15:108-15:110, 15:112-15:113, 15:115-15:122, 15:125-15:129, 15:131-15:135, 15:137-15:144, 15:146-15:150, 15:152-15:154, 15:156-15:161, 15:163-15:175, 15:177-15:178, 15:180-15:184, 15:187, 15:189-15:199, 15:201-15:217, 15:219-15:221, 15:223, 15:225-15:227, 15:228 [?] , 15:229-15:234, 15:236, 15:238-15:249, 15:251-15:274, 15:276-15:282, 15:284-15:311, 15:313-15:331, 15:333-15:347, 15:349-15:382, 15:384-15:392, 15:394-15:395, 15:397-15:398	
47	320 bp	1070 bp	*15:01:01:01-15:01:04, 15:01:06-15:01:30, 15:01:32-15:19, 15:21, 15:23-15:39:02, 15:43, 15:45-15:46, 15:48, 15:51, 15:53-15:58, 15:60-15:68, 15:70-15:82, 15:84-15:85, 15:87-15:92, 15:94N-15:99, 15:101-15:113, 15:115-15:116, 15:118-15:120, 15:122-15:123, 15:125-15:129, 15:131-15:137, 15:139-15:185, 15:187, 15:189-15:198, 15:200-15:217, 15:219-15:223, 15:225-15:227, 15:228 [?] , 15:229-15:237, 15:239-15:240, 15:242:01-15:251, 15:253-15:264, 15:266-15:322, 15:324-15:340, 15:342-15:347, 15:349-15:367, 15:369-15:382, 15:384-15:388, 15:390-15:398	*35:43:01, 35:185, 46:01:01-46:02, 46:04-46:05, 46:07N-46:17, 46:19-46:20, 46:22-46:28, 46:31-46:68, 56:03
48	420 bp	1070 bp	*15:01:01:01-15:01:04, 15:01:06-15:01:16, 15:01:17 ^w , 15:01:18-15:19, 15:21, 15:23-15:32:02, 15:34-15:40:02, 15:43-15:58, 15:60-15:71, 15:73-15:77, 15:79N-15:99, 15:101-15:106, 15:108-15:115, 15:117-15:129, 15:131-15:140, 15:142-15:149N, 15:151-15:187, 15:189-15:227, 15:228 [?] , 15:229-15:247, 15:249-15:286, 15:288-15:294N, 15:296-15:351, 15:353-15:382, 15:384-15:398	*35:43:01, 35:185, 46:01:01-46:05, 46:07N-46:44, 46:46-46:53, 46:55N-46:57, 46:59-46:68, 56:03
49	135 bp	1070 bp	*15:57, 15:135	*13:39, 13:67
50 ⁴	85 bp 185 bp	1070 bp	*15:67, 15:343 *15:206	*35:110, 58:77, A*02:598, A*68:56
51	135 bp 175 bp 200 bp	1070 bp	*15:98, 15:109, 15:163 *15:277-15:278 *15:246N	*35:122, 56:29
52 ⁹	200 bp 345 bp	800 bp	*15:104, 15:307, 15:321Q *15:92, 15:213	*51:142
53	210 bp	800 bp	*15:96	C*03:175, C*03:199, C*03:245
54	250 bp 410 bp	800 bp	*15:187 *15:97	C*12:36
55	175 bp 210 bp 255 bp	1070 bp	*15:214 *15:321Q *15:102, 15:284	*07:120, 18:81, 35:250, 40:150, 51:165, 57:29 *51:30, 73:01-73:02

101.516-24/04 – including *Taq* polymerase, IFU-01
101.516-24u/04u – without *Taq* polymerase, IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **3F7**

Lot-specific information

56⁴	70 bp 260 bp	1070 bp	*15:103 *15:227	*27:69, 37:23, 40:213
57	340 bp 480 bp	1070 bp	*15:111N, 15:300 *15:369	
58	175 bp	1070 bp	*15:113, 15:319	*52:17
59^{5,9}	230 bp 370 bp 420 bp	800 bp	*15:191, 15:234 *15:198, 15:270 *15:94N, 15:294N	C*06:33 *35:161
60	160 bp	800 bp	*15:20, 15:228, 15:271, 15:305	*35:01:01:01-35:01:27, 35:01:29-35:30, 35:32:01-35:34, 35:36-35:42:02, 35:44-35:45, 35:47-35:56, 35:58-35:59:02, 35:61:01-35:72, 35:74-35:78, 35:80-35:86, 35:88-35:134N, 35:136-35:184, 35:186-35:278, 35:280-35:322, 48:02:01-48:02:03, 51:01:01:01-51:24:05, 51:26-51:41N, 51:43-51:46, 51:48-51:156, 51:158:01-51:163, 51:165-51:208, 52:01:01:01-52:06:03, 52:08-52:61, 53:01:01-53:43, 56:05:01-56:05:02, 56:21, 56:36, 58:01:01:01-58:02:02, 58:04-58:10N, 58:12-58:19, 58:21-58:29, 58:31N-58:81, 67:05, 78:01:01-78:09, 81:03, 83:01?
61	160 bp 205 bp	1070 bp	*15:175 *15:125	*40:331 *44:22, 44:105, 49:24
62⁴	105 bp	1070 bp	*15:126	*07:175
63	175 bp 280 bp	1070 bp	*15:127 *15:315	
64⁴	90 bp 135 bp 240 bp 365 bp	800 bp	*15:302N *15:270 *15:301 *15:122	C*03:48, C*05:126, C*06:85
65	405 bp 470 bp 520 bp	1070 bp	*15:375N *15:232 *15:129, 15:342, 15:395	*46:64
66	435 bp	1070 bp	*15:02:01-15:02:11, 15:13:01-15:13:02, 15:20-15:21, 15:25:01-15:25:03, 15:36, 15:44, 15:77, 15:85, 15:88-15:89, 15:118, 15:121, 15:139, 15:144, 15:154, 15:194, 15:204, 15:213-15:214, 15:223, 15:240, 15:250, 15:265, 15:271, 15:283, 15:289, 15:291, 15:297, 15:301-15:302N, 15:319, 15:325, 15:328, 15:330, 15:341, 15:345, 15:357-15:358, 15:366, 15:374, 15:378, 15:393-15:394	*13:01:01-13:01:12, 13:06-13:07N, 13:12:01-13:13:02, 13:17, 13:20-13:23, 13:25-13:26:02, 13:28-13:29, 13:36, 13:39, 13:43, 13:50-13:52, 13:57, 13:60-13:61, 13:63N, 13:73, 13:77-13:78, 13:80, 13:83, 13:86-13:87, 13:90, 13:92, 35:46, 35:207, 35:311, 40:137, 44:08, 44:57, 44:60, 46:33, 46:67, 57:01:01-57:10, 57:12-57:15, 57:17-57:44, 57:46-57:50, 57:52-57:66, 57:68, 57:70-57:81, 58:36
67	265 bp 300 bp	800 bp	*15:53, 15:183, 15:212, 15:336, 15:345 *15:132	*46:10, C*16:85
68⁶	215 bp 380 bp	1070 bp	*15:134 *15:247	C*05:03

101.516-24/04 – including *Taq* polymerase, IFU-01
101.516-24u/04u – without *Taq* polymerase, IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **3F7**

Lot-specific information

69	130 bp	1070 bp	*15:57	
70⁴	90 bp 185 bp	1070 bp	*15:208 *15:138, 15:230, 15:241, 15:297	*07:68:01-07:68:03, 07:214, 07:237, 07:262, 18:35, 18:110, 35:66, 40:77, 40:87:01-40:87:02, 40:121, 40:158, 40:222, 40:237, 44:150, 48:12, 48:14, 48:26, 48:29
71	290 bp 315 bp 400 bp	1070 bp	*15:53, 15:212, 15:336, 15:347 *15:205 *15:139	C*01:124, C*16:85
72⁴	110 bp 200 bp 295 bp 390 bp	1070 bp	*15:226N *15:380N *15:140 *15:201	*35:216N, 39:87N *51:166
73	155 bp	1070 bp	*15:09, 15:30 ^w , 15:42, 15:58 ^w , 15:73 ^w , 15:150 ^w , 15:229, 15:243, 15:252, 15:303 ^w , 15:324 ^w , 15:329, 15:361	*07:05:01:01-07:06:02, 07:32, 07:34, 07:40, 07:53, 07:69, 07:78, 07:80, 07:90, 07:97, 07:105, 07:112, 07:123, 07:137-07:138, 07:140, 07:176, 07:182N, 07:201N, 07:206-07:207, 07:210, 07:213-07:214, 07:222, 07:249, 07:258, 07:262, 07:264, 07:269-07:270, 08:20, 08:79, 13:16, 13:20, 13:31 ^w , 13:48, 13:62, 14:13, 14:53, 27:20, 27:24, 27:153, 35:44 ^w , 38:10, 38:32, 39:35, 39:36 ^w , 39:107, 40:05, 40:15-40:16, 40:23, 40:26, 40:28, 40:51, 40:95, 40:148, 40:158, 40:161, 40:174, 40:183, 40:198, 40:309, 42:19, 46:61 ^w , 48:05, 48:15, 51:01:01:01-51:01:18, 51:01:19 ^w , 51:01:20-51:01:30, 51:01:32-51:01:49, 51:01:51-51:02:02, 51:02:04-51:04, 51:06:01-51:07:02, 51:11N-51:14, 51:16, 51:21-51:24:05, 51:26-51:30, 51:32-51:35, 51:38-51:39, 51:41N, 51:43, 51:46, 51:48-51:53, 51:55- 51:60, 51:61:01 ^w -51:61:02 ^w , 51:62- 51:71, 51:74-51:80, 51:82-51:84, 51:86, 51:88-51:89, 51:95-51:96, 51:98N-51:107, 51:109-51:114, 51:116-51:121, 51:123-51:134, 51:136-51:138, 51:140, 51:142- 51:145, 51:147, 51:149N-51:152, 51:154-51:156, 51:158:01-51:171, 51:173Q-51:180, 51:182-51:193, 51:195-51:196, 51:198-51:205, 51:206 ^w , 51:207-51:208, 52:01:01:01- 52:01:10, 52:01:12-52:02:02, 52:04- 52:09, 52:11, 52:13-52:18, 52:20, 52:21 ^w , 52:22-52:25:02, 52:27-52:43, 52:44 ^w , 52:45-52:56, 52:57 ^w , 52:58- 52:61, 54:20, 54:33 ^w , 55:01:01- 55:01:16, 55:03, 55:05, 55:09, 55:11, 55:15, 55:17, 55:21 ^w , 55:24-55:25,

101.516-24/04 – including *Taq* polymerase, IFU-01
101.516-24u/04u – without *Taq* polymerase, IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **3F7**

Lot-specific information

				55:28-55:29, 55:31, 55:33, 55:36, 55:38, 55:40, 55:44-55:45, 55:52-55:55N, 55:59-55:60, 55:64, 55:66, 55:73-55:76, 55:78, 56:05:01-56:06, 56:21, 56:25, 56:37, 57:57, 57:63, 58:08:01-58:08:02, 78:01:01-78:03, 78:05-78:07, 78:09, C*01:90 , C*01:116 , C*02:68^w , C*03:317 , C*15:24 , C*15:128
74 ⁴	75 bp 230 bp	800 bp	*15:145, 15:176 *15:262N	*35:237, 52:28
75	280 bp 315 bp 390 bp	1070 bp	*15:146 *15:367 *15:201	*51:48, C*02:83 , C*08:24 , C*16:90
76	475 bp	1070 bp	*15:02:01-15:02:11, 15:06, 15:13:01-15:13:02, 15:17:01:01-15:17:05, 15:21, 15:25:01-15:25:03, 15:31, 15:36, 15:44, 15:55, 15:77, 15:88-15:89, 15:106, 15:112, 15:121, 15:128, 15:137, 15:139, 15:144, 15:156, 15:162, 15:168, 15:170, 15:177, 15:196, 15:204, 15:208, 15:213-15:214, 15:216, 15:230, 15:233, 15:240, 15:250, 15:265, 15:271, 15:273, 15:283, 15:289, 15:291, 15:297, 15:301-15:302N, 15:308, 15:319, 15:325, 15:328, 15:330, 15:339, 15:341-15:342, 15:345, 15:356-15:358, 15:360, 15:366, 15:374, 15:378, 15:384, 15:394-15:396	
77 ^{6,8}	185 bp 230 bp	1070 bp	*15:174 *15:123, 15:220, 15:388	
78 ⁴	110 bp	1070 bp	*15:142	*07:29, 07:186, 08:01:11, 08:56:01, 08:133, 13:93, 35:218, 35:256, 51:68, 51:176, 57:49, C*07:02:30
	385 bp 475 bp		*15:147 *15:232, 15:369	*46:42 *46:64
79	215 bp 420 bp	800 bp	*15:182N *15:148, 15:330	
80 ^{4,5}	85 bp 200 bp 225 bp 340 bp	1070 bp	*15:149N *15:380N *15:263 *15:320	
81 ⁶	395 bp	1070 bp	*15:11:01-15:11:07, 15:143, 15:148-15:149N, 15:152, 15:209N, 15:215, 15:284, 15:303, 15:305, 15:350, 15:367	
82 ^{4,5}	80 bp 200 bp 245 bp	1070 bp	*15:81 *15:246N *15:277	*53:35, 58:47 *35:222, 40:76, C*05:30 , C*06:125
83 ⁶	255 bp	1070 bp	*15:159, 15:260	*35:302, 40:79, 41:14, 46:52
84 ⁵	155 bp 195 bp 235 bp	1070 bp	*15:167 *15:211 *15:197:01, 15:215	C*03:168 *35:117
85 ⁵	340 bp 395 bp	1070 bp	*15:157 *15:76, 15:101, 15:255	C*07:335 , C*14:62

101.516-24/04 – including *Taq* polymerase, IFU-01
101.516-24u/04u – without *Taq* polymerase, IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **3F7**

Lot-specific information

	420 bp		*15:164	
86	230 bp	1070 bp	*15:166	
	260 bp		*15:272N	
	400 bp		*15:193	
87 ⁴	85 bp	1070 bp	*15:71, 15:175, 15:204, 15:225, 15:344, 15:381	*18:29, 18:98, 35:50, 35:231, 35:280, 37:22, 40:248, 40:282, 41:12, 41:36, 46:14, 49:03, 51:151, 52:51, 54:02, 55:16, 56:35, 57:39, 57:73, 58:41, C*03:278, C*16:86
	295 bp		*15:169	
88	240 bp	800 bp	*15:171	C*14:67
	450 bp		*15:158, 15:240, 15:377Q	C*05:80
89	365 bp	1070 bp	*15:160, 15:282	C*03:84, C*03:281, C*07:275
	420 bp		*15:94N, 15:294N	
	445 bp		*15:377Q	
90 ⁴	95 bp	1070 bp	*15:202, 15:239, 15:376	*13:58, 18:12:01-18:12:02, 35:10, 35:13, 35:16, 35:28, 35:69, 35:80, 35:226, 35:270, 35:318, 37:01:01-37:01:05, 37:01:07-37:05, 37:07, 37:09-37:27, 37:29-37:59, 40:185, 40:209, 40:308, 44:83, 44:134, 49:03, 52:01:01:01-52:01:28, 52:03-52:15, 52:17, 52:19-52:31:02, 52:33-52:61, 53:17:01-53:17:02, 53:28, 53:38, 78:05-78:06
	195 bp		*15:173	*18:63, 39:90, 50:11, A*01:201, A*24:343
91 ^{6,7}	135 bp	1070 bp	*15:184, 15:203	*18:109, 38:62
	325 bp		*15:170, 15:356	
92 ⁴	95 bp	1070 bp	*15:266, 15:332	*13:37, 38:49, 51:75
	155 bp		*15:172, 15:328	
93 ⁹	255 bp	1070 bp	*15:258N, 15:304N	C*01:124, C*16:85
	295 bp		*15:53, 15:177, 15:212, 15:336, 15:347	
	385 bp		*15:120	
94	260 bp	1070 bp	*15:272N	
	285 bp		*15:123	
	325 bp		*15:178	
95 ⁶	160 bp	1070 bp	*15:358	*46:06, 46:25, 46:30
	340 bp		*15:181N	
	395 bp		*15:40:01-15:40:02, 15:47:01-15:47:02, 15:49, 15:52, 15:114, 15:117, 15:124, 15:138, 15:238, 15:241, 15:389	
96 ¹⁰	-	-	Negative Control	

¹Alleles are assigned by the presence of specific PCR product(s). However, the sizes of the specific PCR products may be helpful in the interpretation of HLA-B*15 SSP typings.

When the primers in a primer mix can give rise to HLA-specific PCR products of more than one length this is indicated if the size difference is more than 20 base pairs. Size differences of 20 base pairs or less are not given. For high resolution SSP kits, the alleles listed are specified according to amplicon length.

Nonspecific amplifications, i.e. a ladder or a smear of bands, may sometimes be seen. GC-rich primers have a higher tendency of giving rise to nonspecific amplifications than other primers.

PCR fragments longer than the control bands may sometimes be observed. Such bands should be disregarded and do not influence the interpretation of the SSP typings.

101.516-24/04 – including *Taq* polymerase, IFU-01
101.516-24u/04u – without *Taq* polymerase, IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: 3F7

Lot-specific information

PCR fragments migrating faster than the control bands, but slower than a 400 bp fragment may be seen in some gel read-outs. Such bands can be disregarded and do not influence the interpretation of the SSP typings. Some primers may give rise to primer oligomer artifacts. Sometimes this phenomenon is an inherent feature of the primer pair(s) of a primer mix. More often it is due to other factors such as too low amount of DNA in the PCR reactions, taking too long time in setting up the PCR reactions, working at elevated room temperature or using thermal cyclers that are not pre-heated.

²The internal positive control primer pairs amplify segments of the human growth hormone gene. The internal positive control bands are 1070 or 800 base pairs respectively, well distribution as outlined in the table. Well number 1 contains the shorter, 800 bp, internal positive control band. The well distribution of the internal controls can help in orientation of the kit on gel photo, as well as allow for kit identification. In the presence of a specific amplification the intensity of the control band often decreases.

³For several HLA Class I alleles 1st and/or 4th exon(s) and beyond, as well as intron nucleotide sequences, are not available. In these instances it is not known whether some of the primers of the SSP sets are completely matched with the target sequences or not. Assumption is made that unknown sequences in these regions are conserved within allelic groups.

⁴HLA-specific PCR products shorter than 125 base pairs have a lower intensity and are less sharp than longer PCR products.

⁵Primer mixes 40, 59, 80, 82, 84 and 85 have a tendency to giving rise to primer oligomer formation.

⁶Primer mixes 5, 15, 30, 36, 68, 77, 81, 83, 91 and 95 may have tendencies of unspecific amplifications.

⁷Primer mix 91 may give rise to a long unspecific amplification product of approximately 800 bp. This should be disregarded when interpreting the B*15 typings.

⁸Primer mixes 23 and 77 may give rise to a lower yield of HLA-specific PCR product than the other HLA-B*15 resolution primer mixes.

⁹In primer mixes 52, 59 and 93 the positive control band may be weaker than for other HLA-B*15 primer mixes.

¹⁰Primer mix 96 contains a negative control, which will amplify more than 95% of HLA amplicons as well as the amplicons generated by the control primer pairs matching the human growth hormone gene. HLA-specific PCR product sizes range from 75 to 200 base pairs and the PCR product generated by the HGH positive control primer pair is 430 base pairs.

'w', might be weakly amplified.

'?', nucleotide sequence information not available for the primer matching sequence.

101.516-24/04 – including *Taq* polymerase, IFU-01
101.516-24u/04u – without *Taq* polymerase, IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **3F7**

Lot-specific information

PRIMER SPECIFICATION

Well No.	1	2	3	4	5	6	7	8	9	10	11	12
Length of spec. PCR product	400	290	120	180	385	525	330	105	365	300	140	400
				390						515		
Length of int. pos. control ¹	800	1070	1070	800	1070	1070	1070	1070	1070	1070	800	1070
5'-primer(s) ²	45 5'-ggA 3'	45 5'-ggA 3'	463 5'-TgA 3'	45 5'-ggA 3'	45 5'-ggA 3'	355 5'-TCA 3'	45 5'-ggA 3'	357 5'-Tgg 3'	2 nd I 5'-CAA 3'	360 5'-CAC 3'	363 5'-AgC 3'	45 5'-ggA 3'
										368 5'-gTT 3'		
										583 5'-gAC 3'		
3'-primer(s) ³	272 5'-Tgg 3'	165 5'-Tgg 3'	538 5'-CCA 3'	1 st I 5'-gAC 3'	259 5'-gTT 3'	3 rd I 5'-ATC 3'	206 5'-CCT 3'	419 5'-Cgg 3'	527 5'-CCA 3'	3 rd I 5'-ATC 3'	463 5'-gCT 3'	272 5'-TgA 3'
	283 5'-Tgg 3'		545 5'-AgA 3'	267 5'-ATg 3'				419 5'-CgA 3'	527 5'-CCA 3'			287 5'-TCC 3'
Well No.	1	2	3	4	5	6	7	8	9	10	11	12

Well No.	13	14	15	16	17	18	19	20	21	22	23	24
Length of spec. PCR product	395	465	105	120	435	195	395	360	205	190	230	205
	435	515		190				505	425			305
Length of int. pos. control ¹	800	1070	1070	1070	1070	1070	800	1070	1070	800	800	1070
5'-primer(s) ²	45 5'-ggA 3'	368 5'-gTg 3'	369 5'-TAC 3'	419 5'-gTC 3'	45 5'-ggA 3'	419 5'-gTC 3'	45 5'-ggA 3'	379 5'-ACC 3'	45 5'-ggA 3'	409 5'-ggC 3'	821 5'-gCT 3'	2 nd I 5'-CAA 3'
		419 5'-gTA 3'						523 5'-CCT 3'				
3'-primer(s) ³	272 5'-TgC 3'	3 rd I 5'-ATC 3'	435 5'-TCT 3'	499 5'-ggA 3'	309 5'-ATC 3'	572 5'-gCg 3'	266 5'-TCC 3'	3 rd I 5'-ATC 3'	83 5'-TgA 3'	559 5'-Cag 3'	916 5'-gAT 3'	369 5'-CCT 3'
	311 5'-ggT 3'			570 5'-CCg 3'			278 5'-TgA 3'		301 5'-gTC 3'			463 5'-gCg 3'
												471 5'-gTT 3'
Well No.	13	14	15	16	17	18	19	20	21	22	23	24

Well No.	25	26	27	28	29	30	31	32	33	34	35	36
Length of spec. PCR product	315	230	115	240	255	320	320	275	200	160	165	360
	385					355	395	350	255			435
						405		380	340			
Length of int. pos. control ¹	1070	1070	1070	1070	1070	1070	1070	1070	1070	1070	1070	800
5'-primer(s) ²	45 5'-ggA 3'	368 5'-gTC 3'	485 5'-CAA 3'	362 5'-gAC 3'	369 5'-TAC 3'	479 5'-..C 3'	488 5'-ggg 3'	2 nd I 5'-CAA 3'	2 nd I 5'-CAA 3'	419 5'-gTC 3'	412 5'-ATA 3'	45 5'-ggA 3'
						526 5'-gTC 3'	560 5'-CAC 3'					
						560 5'-CgA 3'						
3'-primer(s) ³	193 5'-CgT 3'	558 5'-Agg 3'	559 5'-CAg 3'	559 5'-CAg 3'	583 5'-gTg 3'	3 rd I 5'-ATC 3'	3 rd I 5'-ATC 3'	437 5'-AAA 3'	362 5'-TCA 3'	538 5'-gTC 3'	538 5'-CCA 3'	238 5'-CCT 3'
	263 5'-gTA 3'							513 5'-TCT 3'	420 5'-gCT 3'			309 5'-gTg 3'
								544 5'-ggT 3'	505 5'-gCT 3'			
Well No.	25	26	27	28	29	30	31	32	33	34	35	36

101.516-24/04 – including *Taq* polymerase, IFU-01
101.516-24u/04u – without *Taq* polymerase, IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **3F7**

Lot-specific information

Well No.	37	38	39	40	41	42	43	44	45	46	47	48
Length of spec. PCR product	170	225	315	220	465	350	525	515	475	360	320	420
			385	425		460						
Length of int. pos. control ¹	1070	1070	800	1070	1070	1070	1070	1070	1070	800	1070	1070
5'-primer(s) ²	409 5'-ggC 3'	106 5'-CCA 3'	45 5'-ggA 3'	45 5'-ggA 3'	419 5'-gTC 3'	420 5'-TTA 3'	357 5'-Tgg 3'	369 5'-TAC 3'	409 5'-ggC 3'	522 5'-gCT 3'	558 5'-ggA 3'	463 5'-TgA 3'
						522 5'-gCT 3'				527 5'-TgA 3'	560 5'-CCT 3'	
						539 5'-gCT 3'						
3'-primer(s) ³	538 5'-CCg 3'	292 5'-gTC 3'	181 5'-gTT 3'	97 5'-gTg 3'	3 rd 5'-ATC 3'	3 rd 5'-ATC 3'	3 rd 5'-ATC 3'	3 rd 5'-ATC 3'	3 rd 5'-ATC 3'	3 rd 5'-ATC 3'	3 rd 5'-ATC 3'	3 rd 5'-ATC 3'
			203 5'-CTg 3'	302 5'-ggT 3'								
			257 5'-CCT 3'									
			267 5'-ATg 3'									
Well No.	37	38	39	40	41	42	43	44	45	46	47	48

Well No.	49	50	51	52	53	54	55	56	57	58	59	60
Length of spec. PCR product	135	85	135	200	210	250	175	70	340	175	230	160
		185	175	345		410	210	260	480		370	
			200				255				420	
Length of int. pos. control ¹	1070	1070	1070	800	800	800	1070	1070	1070	1070	800	800
5'-primer(s) ²	209 5'-ggC 3'	463 5'-TgA 3'	106 5'-CCA 3'	45 5'-ggA 3'	419 5'-gTC 3'	45 5'-ggA 3'	3 rd 5'-Agg 3'	25 5'-CCg 3'	400 5'-TCT 3'	133 5'-CCA 3'	45 5'-ggA 3'	652 5'-CCg 3'
				3 rd 5'-Agg 3'				652 5'-CCA 3'	539 5'-gTg 3'	145 5'-CAA 3'	673 5'-CCA 3'	661 5'-ACg 3'
								540 5'-TgC 3'			679 5'-ggC 3'	
3'-primer(s) ³	302 5'-ggT 3'	506 5'-TgT 3'	200 5'-TCA 3'	221 5'-ATA 3'	589 5'-CTT 3'	127 5'-CTg 3'	704 5'-CCA 3'	56 5'-TCg 3'	3 rd 5'-ATC 3'	272 5'-Tgg 3'	107 5'-ACT 3'	774 5'-ggT 3'
		608 5'-gCC 3'	200 5'-TCA 3'	723 5'-Tgg 3'		288 5'-gTg 3'	735 5'-TCg 3'	872 5'-TCA 3'			295 5'-TCA 3'	
			239 5'-gCT 3'	727 5'-CCA 3'			779 5'-CTA 3'				916 5'-gAT 3'	
			263 5'-..g 3'	735 5'-TCg 3'			787 5'-TCC 3'					
Well No.	49	50	51	52	53	54	55	56	57	58	59	60

101.516-24/04 – including *Taq* polymerase, IFU-01
101.516-24u/04u – without *Taq* polymerase, IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **3F7**

Lot-specific information

Well No.	61	62	63	64	65	66	67	68	69	70	71	72
Length of spec. PCR product	160	105	175	90	405	435	265	215	130	90	290	110
	205		280	135	470		300	380		185	315	200
				240	520						400	295
				365								390
Length of int. pos. control ¹	1070	1070	1070	800	1070	1070	800	1070	1070	1070	1070	1070
5'-primer(s) ²	106 5'-CCg 3'	363 5'-AgC 3'	25 5'-CCg 3'	45 5'-ggA 3'	363 5'-AAT 3'	209 5'-ggC 3'	45 5'-ggA 3'	419 5'-gTC 3'	209 5'-ggC 3'	463 5'-TgA 3'	45 5'-ggA 3'	261 5'-AAC 3'
	151 5'-gCC 3'		419 5'-gTC 3'	679 5'-ggC 3'	410 5'-gCg 3'	209 5'-ggC 3'		668 5'-ggT 3'				662 5'-CCT 3'
				724 5'-g.C 3'	479 5'-..C 3'							754 5'-ACg 3'
3'-primer(s) ³												852 5'-CAA 3'
	272 5'-Tgg 3'	425 5'-CgC 3'	134 5'-AgT 3'	115 5'-ggA 3'	3 rd I 5'-ATC 3'	353 5'-TgA 3'	142 5'-ggT 3'	593 5'-CgC 3'	299 5'-TCA 3'	510 5'-CAA 3'	167 5'-ACA 3'	331 5'-CTA 3'
			554 5'-CCg 3'	242 5'-CCA 3'		353 5'-TgA 3'	177 5'-AAA 3'	916 5'-gAT 3'		605 5'-gCT 3'	189 5'-gCT 3'	916 5'-gAT 3'
				774 5'-ggT 3'		362 5'-TCA 3'					277 5'-ggC 3'	
Well No.	61	62	63	64	65	66	67	68	69	70	71	72

Well No.	73	74	75	76	77	78	79	80	81	82	83	84
Length of spec. PCR product	155	75	280	475	185	110	215	85	395	80	255	155
		230	315		230	385	420	200		200		195
			390			475		225		245		235
								340				
Length of int. pos. control ¹	1070	800	1070	1070	1070	1070	800	1070	1070	1070	1070	1070
5'-primer(s) ²	412 5'-ATA 3'	368 5'-gTA 3'	662 5'-CCT 3'	409 5'-ggT 3'	47 5'-AgT 3'	97 5'-TCg 3'	45 5'-ggA 3'	267 5'-Ag 3'	45 5'-ggA 3'	106 5'-CCA 3'	106 5'-CCA 3'	363 5'-AgT 3'
		523 5'-CCg 3'	736 5'-gCA 3'		463 5'-TgA 3'	400 5'-TCT 3'		709 5'-AgT 3'				404 5'-CCT 3'
			769 5'-Agg 3'			410 5'-gCg 3'		824 5'-TTg 3'				446 5'-CgT 3'
						495 5'-CAC 3'		852 5'-CAA 3'				
3'-primer(s) ³	527 5'-CCT 3'	559 5'-CAG 3'	916 5'-gAT 3'	3 rd I 5'-ATC 3'	106 5'-CAT 3'	165 5'-Tgg 3'	93 5'-AAC 3'	311 5'-ggT 3'	272 5'-TgT 3'	146 5'-CCg 3'	322 5'-gTg 3'	559 5'-CAG 3'
					605 5'-gCA 3'	3 rd I 5'-ATC 3'	296 5'-CTT 3'	916 5'-gAT 3'		263 5'-..g 3'		
										312 5'-AgT 3'		
Well No.	73	74	75	76	77	78	79	80	81	82	83	84

101.516-24/04 – including *Taq* polymerase, IFU-01
101.516-24u/04u – without *Taq* polymerase, IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **3F7**

Lot-specific information

Well No.	85	86	87	88	89	90	91	92	93	94	95
Length of spec.	340	230	85	240	365	95	135	95	255	260	160
PCR product	395	260	295	450	420	195	325	155	295	285	340
	420	400			445				385	325	395
Length of int.	1070	1070	1070	800	1070	1070	1070	1070	1070	1070	1070
pos. control ¹											
5'-primer(s) ²	45	2 nd I	97	45	45	206	45	463	45	2 nd I	2 nd I
	5'-ggA 3'	5'-CAA 3'	5'-TCC 3'	5'-ggA 3'	5'-ggA 3'	5'-gAC 3'	5'-ggA 3'	5'-TgA 3'	5'-ggA 3'	5'-CAA 3'	5'-CAA 3'
			2 nd I			463	365		369	767	3rd I
			5'-CAA 3'			5'-TgA 3'	5'-gAC 3'		5'-TAC 3'	5'-ggg 3'	5'-Agg 3'
							376				
							5'-gCA 3'				
3'-primer(s) ³	218	394	142	118	242	259	203	514	135	426	502
	5'-gCA 3'	5'-gCT 3'	5'-TgC 3'	5'-gCT 3'	5'-CCT 3'	5'-CTC 3'	5'-TCT 3'	5'-CTg 3'	5'-Ag 3'	5'-TCC 3'	5'-CTA 3'
	272	426	457	320	247	616	463	515	167	487	559
	5'-TgT 3'	5'-TCC 3'	5'-gTA 3'	5'-AgC 3'	5'-ATg 3'	5'-CgT 3'	5'-gCT 3'	5'-CCC 3'	5'-ACA 3'	5'-CgT 3'	5'-CTC 3'
	296	565		323	295			578	173	916	685
	5'-CTg 3'	5'-CAT 3'		5'-AgA 3'	5'-TCA 3'			5'-TgT 3'	5'-TCC 3'	5'-gAT 3'	5'-ggT 3'
				337	320				261		
				5'-CTT 3'	5'-AgC 3'				5'-gTA 3'		
									585		
									5'-AgT 3'		
Well No.	85	86	87	88	89	90	91	92	93	94	95

¹The internal positive control primer pairs amplify segments of the human growth hormone gene. The internal positive control bands are 1070 or 800 base pairs respectively, well distribution as outlined in the table. Well number 1 contains the shorter, 800 bp, internal positive control band. The well distribution of the internal controls can help in orientation of the kit on gel photo, as well as allow for kit identification. In the presence of a specific amplification the intensity of the control band often decreases.

²The nucleotide position matching the specificity-determining 3'-end of the primer is given. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

³The nucleotide position matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

101.516-24/04 – including *Taq* polymerase, IFU-01
101.516-24u/04u – without *Taq* polymerase, IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **3F7**

Lot-specific information

[illegible]

101.516-24/04 – including *Taq* polymerase, IFU-01
101.516-24u/04u – without *Taq* polymerase, IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **3F7**

Lot-specific information

CELL LINE VALIDATION SHEET																				
HLA-B*15 SSP subtyping kit ²																				
					Well															
				Prod. No.:	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
					201558717	201558718	201558719	201558720	201558721	201558722	201558723	201781524	201558725	201558726	201558727	201558728	201558729	201781530	201558731	201558732
	IHCW cell line ¹		B*																	
1	9001	SA	*07:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
2	9280	LK707	*52:01	*73:01	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-
3	9011	E4181324	*52:01		-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-
4	9275	GU373	*15:10	*53:01	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-
5	9009	KAS011	*37:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
6	9353	SM	*39:01	*51:01	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-
7	9020	QBL	*18:01		-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-
8	9025	DEU	*35:01		-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-
9	9026	YAR	*38:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
10	9107	LKT3	*54:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
11	9051	PITOUT	*44:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
12	9052	DBB	*57:01		-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-
13	9004	JESTHOM	*27:05		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
14	9071	OLGA	*15:01	*15:20	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-
15	9075	DKB	*40:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
16	9037	SWEIG007	*40:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
17	9282	CTM3953540	*08:01	*55:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
18	9257	32367	*14:01	*56:01	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-
19	9038	BM16	*18:01		-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-
20	9059	SLE005	*40:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
21	9064	AMALA	*15:01		-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-
22	9056	KOSE	*35:03		-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-
23	9124	IHL	*40:02	*56:02	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-
24	9035	JBUSH	*38:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
25	9049	IBW9	*14:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
26	9285	WT49	*58:01		-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-
27	9191	CH1007	*07:05	*51:01	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-
28	9320	BEL5GB	*44:02	*44:03	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
29	9050	MOU	*44:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
30	9021	RSH	*42:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
31	9019	DUCAF	*18:01		-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-
32	9297	HAG	*41:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
33	9098	MT14B	*40:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
34	9104	DHIF	*38:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
35	9302	SSTO	*44:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
36	9024	KT17	*15:01	*35:01	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-
37	9065	HHKB	*07:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
38	9099	LZL	*15:01		-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-
39	9315	CML	*08:01	*27:05	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
40	9134	WHONP199	*13:02	*46:01	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-
41	9055	H0301	*14:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
42	9066	TAB089	*46:01		-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-
43	9076	T7526	*46:01		-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-
44	9057	TEM	*38:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
45	9239	SHJO	*42:01	*50:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
46	9013	SCHU	*07:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
47	9045	TUBO	*51:01		-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-
48	9303	TER-ND	*35:01	*44:03	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-

101.516-24/04 – including *Taq* polymerase, IFU-01
101.516-24u/04u – without *Taq* polymerase, IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **3F7**

Lot-specific information

CELL LINE VALIDATION SHEET																													
HLA-B*15 SSP subtyping kit ²																													
						Well																							
					Prod. No.:	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48								
						201558733	201558734	201558735	201558736	201558737	201558738	201558739	201558740	201558741	201781542	201558743	201558744	201558745	201781546	201781547	201558748								
	IHCW cell line ¹			B*																									
1	9001	SA	*07:02			-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
2	9280	LK707	*52:01	*73:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
3	9011	E4181324	*52:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
4	9275	GU373	*15:10	*53:01		-	-	-	-	-	-	-	-	-	-	+	-	+	+	+	+	+	+	+	+	+	+		
5	9009	KAS011	*37:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
6	9353	SM	*39:01	*51:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
7	9020	QBL	*18:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
8	9025	DEU	*35:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
9	9026	YAR	*38:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
10	9107	LKT3	*54:01			-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
11	9051	PITOUT	*44:03			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
12	9052	DBB	*57:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
13	9004	JESTHOM	*27:05			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
14	9071	OLGA	*15:01	*15:20		-	-	-	-	-	-	-	-	-	+	-	-	+	+	+	+	+	+	+	+	+	+		
15	9075	DKB	*40:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
16	9037	SWEIG007	*40:02			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
17	9282	CTM3953540	*08:01	*55:01		-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
18	9257	32367	*14:01	*56:01		-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
19	9038	BM16	*18:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
20	9059	SLE005	*40:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
21	9064	AMALA	*15:01			-	-	-	-	-	-	-	-	+	-	-	+	+	+	+	+	+	+	+	+	+	+		
22	9056	KOSE	*35:03			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
23	9124	IHL	*40:02	*56:02		-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
24	9035	JBUSH	*38:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
25	9049	IBW9	*14:02			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
26	9285	WT49	*58:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
27	9191	CH1007	*07:05	*51:01		-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
28	9320	BEL5GB	*44:02	*44:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
29	9050	MOU	*44:03			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
30	9021	RSH	*42:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
31	9019	DUCAF	*18:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
32	9297	HAG	*41:02			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
33	9098	MT14B	*40:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
34	9104	DHIF	*38:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
35	9302	SSTO	*44:02			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
36	9024	KT17	*15:01	*35:01		-	-	-	-	-	-	-	-	+	-	-	+	+	+	+	+	+	+	+	+	+	+		
37	9065	HHKB	*07:02			-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
38	9099	LZL	*15:01			-	-	-	-	-	-	-	-	+	-	-	+	+	+	+	+	+	+	+	+	+	+		
39	9315	CML	*08:01	*27:05		-	-	-	-	-	W	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
40	9134	WHONP199	*13:02	*46:01		-	-	-	-	-	+	-	-	+	-	-	+	+	+	+	+	+	+	+	+	+	+		
41	9055	H0301	*14:02			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
42	9066	TAB089	*46:01			-	-	-	-	-	+	-	-	+	-	-	+	+	+	+	+	+	+	+	+	+	+		
43	9076	T7526	*46:01			-	-	-	-	-	+	-	-	+	-	-	+	+	+	+	+	+	+	+	+	+	+		
44	9057	TEM	*38:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
45	9239	SHJO	*42:01	*50:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
46	9013	SCHU	*07:02			-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
47	9045	TUBO	*51:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
48	9303	TER-ND	*35:01	*44:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		

101.516-24/04 – including *Taq* polymerase, IFU-01
101.516-24u/04u – without *Taq* polymerase, IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: 3F7

Lot-specific information

[illegible]

101.516-24/04 – including *Taq* polymerase, IFU-01
101.516-24u/04u – without *Taq* polymerase, IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **3F7**

Lot-specific information

CELL LINE VALIDATION SHEET																			
HLA-B*15 SSP subtyping kit ²																			
				Well															
				65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80
				Prod. No.:	201781565	201781566	201558767	201558768	201558769	201558770	201558771	201781572	201558773	201558774	201781575	201558776	201781577	201781578	201781579
IHCW cell line ¹				B*															
1	9001	SA	*07:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
2	9280	LK707	*52:01	*73:01	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-
3	9011	E4181324	*52:01		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-
4	9275	GU373	*15:10	*53:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
5	9009	KAS011	*37:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
6	9353	SM	*39:01	*51:01	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-
7	9020	QBL	*18:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
8	9025	DEU	*35:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
9	9026	YAR	*38:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
10	9107	LKT3	*54:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
11	9051	PITOUT	*44:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
12	9052	DBB	*57:01		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
13	9004	JESTHOM	*27:05		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
14	9071	OLGA	*15:01	*15:20	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
15	9075	DKB	*40:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
16	9037	SWEIG007	*40:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
17	9282	CTM3953540	*08:01	*55:01	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-
18	9257	32367	*14:01	*56:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
19	9038	BM16	*18:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
20	9059	SLE005	*40:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
21	9064	AMALA	*15:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
22	9056	KOSE	*35:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
23	9124	IHL	*40:02	*56:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
24	9035	JBUSH	*38:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
25	9049	IBW9	*14:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
26	9285	WT49	*58:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
27	9191	CH1007	*07:05	*51:01	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-
28	9320	BEL5GB	*44:02	*44:03	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
29	9050	MOU	*44:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
30	9021	RSH	*42:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
31	9019	DUCAF	*18:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
32	9297	HAG	*41:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
33	9098	MT14B	*40:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
34	9104	DHIF	*38:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
35	9302	SSTO	*44:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
36	9024	KT17	*15:01	*35:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
37	9065	HHKB	*07:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
38	9099	LZL	*15:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
39	9315	CML	*08:01	*27:05	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
40	9134	WHONP199	*13:02	*46:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
41	9055	H0301	*14:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
42	9066	TAB089	*46:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
43	9076	T7526	*46:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
44	9057	TEM	*38:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
45	9239	SHJO	*42:01	*50:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
46	9013	SCHU	*07:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
47	9045	TUBO	*51:01		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-
48	9303	TER-ND	*35:01	*44:03	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-

101.516-24/04 – including *Taq* polymerase, IFU-01
101.516-24u/04u – without *Taq* polymerase, IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **3F7**

Lot-specific information

CELL LINE VALIDATION SHEET																				
HLA-B*15 SSP subtyping kit ²																				
					Well															
					Prod. No.:	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95
					Prod. No.:	201558781	201781582	201781583	201558784	201558785	201558786	201558787	201781588	201781589	201558790	201558791	201558792	201558793	201558794	201558795
	IHWC cell line ¹		B*																	
1	9001	SA	*07:02			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
2	9280	LK707	*52:01	*73:01		-	-	-	-	-	-	-	-	-	+	-	-	-	-	-
3	9011	E4181324	*52:01			-	-	-	-	-	-	-	-	-	+	-	-	-	-	-
4	9275	GU373	*15:10	*53:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
5	9009	KAS011	*37:01			-	-	-	-	-	-	-	-	-	+	-	-	-	-	-
6	9353	SM	*39:01	*51:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
7	9020	QBL	*18:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
8	9025	DEU	*35:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
9	9026	YAR	*38:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
10	9107	LKT3	*54:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
11	9051	PITOUT	*44:03			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
12	9052	DBB	*57:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
13	9004	JESTHOM	*27:05			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
14	9071	OLGA	*15:01	*15:20		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
15	9075	DKB	*40:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
16	9037	SWEIG007	*40:02			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
17	9282	CTM3953540	*08:01	*55:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
18	9257	32367	*14:01	*56:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
19	9038	BM16	*18:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
20	9059	SLE005	*40:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
21	9064	AMALA	*15:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
22	9056	KOSE	*35:03			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
23	9124	IHL	*40:02	*56:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
24	9035	JBUSH	*38:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
25	9049	IBW9	*14:02			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
26	9285	WT49	*58:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
27	9191	CH1007	*07:05	*51:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
28	9320	BEL5GB	*44:02	*44:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
29	9050	MOU	*44:03			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
30	9021	RSH	*42:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
31	9019	DUCAF	*18:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
32	9297	HAG	*41:02			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
33	9098	MT14B	*40:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
34	9104	DHIF	*38:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
35	9302	SSTO	*44:02			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
36	9024	KT17	*15:01	*35:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
37	9065	HHKB	*07:02			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
38	9099	LZL	*15:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
39	9315	CML	*08:01	*27:05		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
40	9134	WHONP199	*13:02	*46:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
41	9055	H0301	*14:02			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
42	9066	TAB089	*46:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
43	9076	T7526	*46:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
44	9057	TEM	*38:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
45	9239	SHJO	*42:01	*50:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
46	9013	SCHU	*07:02			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
47	9045	TUBO	*51:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
48	9303	TER-ND	*35:01	*44:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-

101.516-24/04 – including *Taq* polymerase, IFU-01
101.516-24u/04u – without *Taq* polymerase, IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: 3F7

Lot-specific information

¹The provided cell line HLA specificities are retrieved from the <http://www.ihwg.org/hla> web site. The specificity of an individual cell line may thus be subject to change.

²The specificity of each primer solution in the kit has been tested against 48 well characterized cell line DNAs and where applicable, additional cell line DNAs.

No DNAs carrying the alleles to be amplified by primer solutions 4, 9, 10, 18, 24 to 27, 30 to 34, 36, 39, 49 to 54, 56 to 59, 61 to 65, 67 to 72, 74, 75, 78 to 80, 82 to 89 and 91 to 95 were available. The specificities of the primers in primer solutions 9, 10, 18, 24 to 27, 30, 31, 33, 34, 36, 49, 53, 59, 61, 63, 64, 65, 67, 68, 70 to 72, 78, 82, 85, 87, 91 and 93 to 95 were tested by separately adding one or two additional 5'-primers, respectively one to three additional 3'-primers. In primer solutions 4, 32, 39, 50 to 52, 54, 56, 62, 69, 79, 86, 88, 89 and 92 it was only possible to test the 5'-primers, the 3'-primers were not possible to test. In primer solutions 57, 58, 74, 75, 80 and 84 it was only possible to test the 3'-primers, the 5'-primers were not possible to test. In primer solutions 1, 3, 12, 13, 19, 21, 24, 25, 33, 36, 55, 59, 63, 64, 67, 68, 70 to 72, 77, 82, 85, 87, 90, 91, 93 to 95 one or more 3'-primers were not possible to test. In primer solutions 10, 14, 20, 30, 31, 42, 46, 47, 59 to 61, 64, 65, 72, 77, 78, 91 and 94 one, two or three 5'-primers were not possible to test. Additional primers in primer solutions 16, 40, 42, 66, 77 and 90 were tested by separately adding one additional 5'-primer or one additional 3'-primer.

101.516-24/04 – including *Taq* polymerase, IFU-01
101.516-24u/04u – without *Taq* polymerase, IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **3F7**

Lot-specific information

ADDRESSES:

Manufacturer:

Olerup SSP AB, Franzengatan 5, SE-112 51 Stockholm, Sweden.

Tel: +46-8-717 88 27

Fax: +46-8-717 88 18

E-mail: olerup-se@caredx.com

Web page: <http://www.olerup.com>

Distributed by:

Olerup GmbH, Löwengasse 47 / 6, AT-1030 Vienna, Austria.

Tel: +43-1-710 15 00

Fax: +43-1-710 15 00 10

E-mail: olerup-at@caredx.com

Web page: <http://www.olerup.com>

Olerup Inc., 901 S. Bolmar St., Suite R, West Chester, PA 19382

Tel: 1-877-OLERUP1

Fax: 610-344-7989

E-mail: olerup-us@caredx.com

Web page: <http://www.olerup.com>

For information on *Olerup* distributors worldwide, contact **Olerup GmbH**.