

101.570-06 – including *Taq* polymerase, IFU-01
101.570-06u – without *Taq* polymerase, IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **4L7**

Lot-specific information

Olerup SSP® HLA-B*55

Product number:	101.570-06 – including <i>Taq</i> polymerase 101.570-06u – without <i>Taq</i> polymerase
Lot number:	4L7
Expiry date:	2024-08-01
Number of tests:	6
Number of wells per test:	31+1
Storage - pre-aliquoted primers:	dark at -20°C
- PCR Master Mix:	-20°C
- Adhesive PCR seals	RT
- Product Insert	RT

This Product Description is only valid for Lot No. 4L7.

Complete product documentation consists of generic Instructions for Use (IFU), lot specific Product Insert, Worksheet and Certificate.

CHANGES COMPARED TO THE PREVIOUS OLERUP SSP® HLA-B*55 Lot (9G5)

- The product documentation has been updated for new alleles of IMGT 3.40.0
- The kit resolution focuses on common and well documented (CWD) alleles¹.

The HLA-B*55 specificity and interpretation tables have been updated for the HLA-B alleles described since the previous *Olerup SSP®* HLA-B*55 lot was made (**Lot No. 9G5**).

The primers of the wells detailed below have been exchanged, added or modified compared to the previous lot.

Well	5'-primer	3'-primer	rationale
10	Added	Added	Primer pair added for from mix 30.
15	Added	-	Primer pair added for the B*55:01:07 allele.
30	Moved	-	Primer pair moved to mix 10 for improved HLA-specific amplification.

101.570-06 – including *Taq* polymerase, IFU-01
101.570-06u – without *Taq* polymerase, IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **4L7**

Lot-specific information

Well **32** contains Negative Control primer pairs, that will amplify a majority of the *Olerup SSP[®]* HLA Class I, DRB, DQB1, DPB1 and DQA1 amplicons as well as all the amplicons generated by the control primer pairs matching the human growth hormone gene.

HLA-specific PCR product sizes range from 75 to 200 base pairs.

The PCR product generated by the positive control primer pair is 200 base pairs.

Length of PCR product	105	200	105	80	75	80	85
5'-primer¹	164	340	440	45	45	43	36
	5'-CAC ^{3'}	5'-Agg ^{3'}	5'-TTA ^{3'}	5'-Tgg ^{3'}	5'-Tgg ^{3'}	5'-Tgg ^{3'}	5'-TAC ^{3'}
							36
							5'-TAT ^{3'}
3'-primer²	231	2nd I	507	59	58	57	47
	5'-TgC ^{3'}	5'-AAA ^{3'}	5'-TTg ^{3'}	5'-CTC ^{3'}	5'-ggC ^{3'}	5'-CTC ^{3'}	5'-ACA ^{3'}
							48
							5'-gCA ^{3'}
							48
							5'-gCC ^{3'}
							52
							5'-TgT ^{3'}
A*	+	+	+				
B*	+	+	+				
C*	+	+	+				
DRB1				+	+		
DRB3				+	+		
DRB5				+			
DQB1					+		
DPB1						+	
DQA1							+

¹The nucleotide position for HLA class I genes and the codon for HLA class II genes, in the 2nd or 3rd exon, matching the specificity-determining 3'-end of the primer is given. Nucleotide and codon numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

²The nucleotide position for HLA class I genes and the codon for HLA class II genes, in the 2nd or 3rd exon or the 2nd intron, matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Nucleotide and codon numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

101.570-06 – including *Taq* polymerase, IFU-01
101.570-06u – without *Taq* polymerase, IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **4L7**

Lot-specific information

PRODUCT DESCRIPTION

HLA-B*55 SSP subtyping

CONTENT

The primer set contains 5'- and 3'-primers for identifying the B*55:01 to B*55:107 alleles.

PLATE LAYOUT

Each test consists of 32 PCR reactions in a 32 well PCR plate.

1	2	3	4	5	6	7	8
9	10	11	12	13	14	15	16
17	18	19	20	21	22	23	24
25	26	27	28	29	30	31	NC

The 32 well PCR plate is marked with ‘HLA-B*55’ in silver/gray ink.

Well No. 1 is marked with the Lot No. ‘4L7’.

Wells 1 to 31 – HLA-B*55 high resolution primers.

Well 32 – Negative Control (NC).

A faint row of numbers is seen between wells 1 and 2 or wells 7 and 8 of the PCR trays. These stem from the manufacture of the trays, and should be disregarded.

The PCR plates are covered with a PCR-compatible foil.

Please note: When removing each 32 well PCR plate, make sure that the remaining plates stay covered. Use a scalpel or a similar instrument to carefully cut the foil between the plates.

INTERPRETATION

Due to the sharing of sequence motifs between HLA-B alleles non-HLA-B*55 alleles will be amplified by some primer mixes. For further details see Specificity Table.

UNIQUELY IDENTIFIED ALLELES

All the HLA-B*55 alleles, i.e. **B*55:01 to B*55:107**, recognized by the HLA Nomenclature Committee in April 2020^{1,2} will be amplified by the primers in the HLA-B*55 subtyping kit³.

The HLA-B*55 kit enables separation of the confirmed HLA-B*55 alleles as listed in the IMGT/HLA database 3.26.0. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources. Current allele confirmation status for HLA-B*55 alleles is listed below.

The HLA-B*55 kit also enables identification of many null and alternatively expressed alleles.

101.570-06 – including *Taq* polymerase, IFU-01
101.570-06u – without *Taq* polymerase, IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **4L7**

Lot-specific information

The following HLA-B*55 alleles can be distinguished by the different sizes of the HLA-specific PCR product:

Alleles	Primer mix
B*55:26, 55:62	25

¹HLA-B alleles listed on the IMGT/HLA web page 2020-April-20, release 3.40.0, www.ebi.ac.uk/imgt/hla.

²Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

³The HLA-B*55 primer set cannot separate the following alleles. These alleles can be distinguished by the HLA-B low resolution kit and/or the HLA-B*15, B*51 and B*56 high resolution kit.

Alleles

B*55:78, 15:42, 51:300
B*55:22, 56:60

101.570-06 – including *Taq* polymerase, IFU-01
101.570-06u – without *Taq* polymerase, IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **4L7**

Lot-specific information

ALLELE CONFIRMATION STATUS

Allele	Status ¹	Allele	Status ¹	Allele	Status ¹	Allele	Status ¹
B*55:01:01	Confirmed	B*55:07	Confirmed	B*55:37	Confirmed	B*55:67	Unconfirmed
B*55:01:02	Unconfirmed	B*55:08	Confirmed	B*55:38	Unconfirmed	B*55:68	Unconfirmed
B*55:01:03	Confirmed	B*55:09	Unconfirmed	B*55:39	Unconfirmed	B*55:69	Unconfirmed
B*55:01:04	Confirmed	B*55:10	Confirmed	B*55:40	Unconfirmed	B*55:70	Unconfirmed
B*55:01:05	Confirmed	B*55:11	Confirmed	B*55:41	Unconfirmed	B*55:71	Unconfirmed
B*55:01:06	Confirmed	B*55:12	Confirmed	B*55:42	Unconfirmed	B*55:72	Confirmed
B*55:01:07	Unconfirmed	B*55:13	Unconfirmed	B*55:43	Unconfirmed	B*55:73	Unconfirmed
B*55:01:08	Unconfirmed	B*55:14	Unconfirmed	B*55:44	Unconfirmed	B*55:74	Unconfirmed
B*55:01:09	Unconfirmed	B*55:15	Unconfirmed	B*55:45	Confirmed	B*55:75	Unconfirmed
B*55:01:10	Unconfirmed	B*55:16	Confirmed	B*55:46	Unconfirmed	B*55:76	Unconfirmed
B*55:01:11	Unconfirmed	B*55:17	Confirmed	B*55:47	Unconfirmed	B*55:77	Unconfirmed
B*55:01:12	Unconfirmed	B*55:18	Confirmed	B*55:48	Unconfirmed	B*55:78	Unconfirmed
B*55:01:13	Unconfirmed	B*55:19	Unconfirmed	B*55:49	Unconfirmed	B*55:80	Unconfirmed
B*55:01:14	Unconfirmed	B*55:20	Confirmed	B*55:50	Confirmed		
B*55:01:15	Unconfirmed	B*55:21	Unconfirmed	B*55:51	Unconfirmed		
B*55:01:16	Unconfirmed	B*55:22	Unconfirmed	B*55:52	Unconfirmed		
B*55:02:01:01	Confirmed	B*55:23	Confirmed	B*55:53	Unconfirmed		
B*55:02:01:02	Unconfirmed	B*55:24	Unconfirmed	B*55:54	Confirmed		
B*55:02:02	Unconfirmed	B*55:25	Confirmed	B*55:55N	Unconfirmed		
B*55:02:03	Confirmed	B*55:26	Unconfirmed	B*55:56	Confirmed		
B*55:02:04	Unconfirmed	B*55:27	Unconfirmed	B*55:57	Confirmed		
B*55:02:05	Unconfirmed	B*55:28	Unconfirmed	B*55:58	Unconfirmed		
B*55:02:06	Unconfirmed	B*55:29	Confirmed	B*55:59	Unconfirmed		
B*55:02:07	Confirmed	B*55:30	Unconfirmed	B*55:60	Unconfirmed		
B*55:02:08	Unconfirmed	B*55:31	Confirmed	B*55:61	Unconfirmed		
B*55:02:09	Unconfirmed	B*55:32	Unconfirmed	B*55:62	Unconfirmed		
B*55:02:10	Unconfirmed	B*55:33	Confirmed	B*55:63	Unconfirmed		
B*55:03	Unconfirmed	B*55:34	Confirmed	B*55:64	Unconfirmed		
B*55:04	Confirmed	B*55:35	Unconfirmed	B*55:65	Unconfirmed		
B*55:05	Unconfirmed	B*55:36	Confirmed	B*55:66	Confirmed		

¹Allele status “confirmed” or “unconfirmed” as listed on the IMGT/HLA web page 2016-October-14, release 3.26.0, www.ebi.ac.uk/imgt/hla.

RESOLUTION IN HOMO- AND HETEROZYGOES

Results file with resolution in HLA-B*55 homo- and heterozygotes is available upon request.

101.570-06 – including *Taq* polymerase, IFU-01
101.570-06u – without *Taq* polymerase, IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **4L7**

Lot-specific information

SPECIFICITY TABLE

HLA-B*55 SSP subtyping

Specificities and sizes of the PCR products of the 31+1 primer mixes used for
HLA-B*55 SSP subtyping

Primer Mix	Size of spec. PCR product ¹	Size of control band ²	Amplified HLA-B*55 alleles ³	Other amplified HLA Class I alleles
1	400 bp	800 bp	*55:01:01:01-55:05, 55:07-55:17, 55:19-55:33, 55:35-55:46, 55:48-55:77, 55:79-55:107	*07:65, 07:209, 07:271, 27:157, 35:76, 40:166 [?] , 44:90, 54:01:01:01-54:23, 54:25-54:29, 54:31-54:41, 56:01:01:01-56:13, 56:15-56:16, 56:18-56:22, 56:24-56:32, 56:34-56:41, 56:43, 56:45-56:53, 56:55:01:02-56:66, 56:69-56:73, 81:01:01:01-81:01:05, 81:03 [?] -81:04N [?] , 81:06 [?] -81:08 [?] , 81:09, 82:01:01:01-82:03, 83:01
2⁴	125 bp	1070 bp	*55:01:01:01-55:03, 55:05, 55:07, 55:09-55:13, 55:15-55:16, 55:18-55:19, 55:21-55:37, 55:39-55:48, 55:50, 55:52-55:71, 55:73-55:85, 55:87-55:107	*07:235, 13:01:01:01-13:02:30, 13:02:32-13:04, 13:06-13:09, 13:11-13:20, 13:22:01-13:23, 13:25-13:34, 13:36-13:61, 13:63N-13:64, 13:66-13:70, 13:72-13:78, 13:80-13:107, 13:109-13:141, 14:22, 14:37, 14:55, 15:16:01:01-15:16:04, 15:20, 15:34, 15:42, 15:62, 15:67, 15:85-15:86, 15:95, 15:137, 15:194, 15:222, 15:242:02, 15:254, 15:303, 15:362, 15:393, 15:407-15:408, 15:485-15:486, 15:550, 15:552, 18:22, 18:52, 18:105, 27:14, 35:01:01:01-35:01:42, 35:01:44-35:05:02, 35:06-35:08:13, 35:10-35:11:03, 35:13-35:15:02, 35:19-35:21, 35:23-35:30:01:02, 35:33-35:38, 35:40N-35:42:02, 35:45-35:50, 35:52, 35:54-35:58, 35:60, 35:62-35:63, 35:65Q-35:66, 35:68:01-35:72, 35:74-35:78, 35:80-35:85, 35:89-35:98, 35:100-35:101:02, 35:103-35:112, 35:115-35:116, 35:119-35:124, 35:126-35:134N, 35:136-35:148, 35:150:01-35:163, 35:166-35:173:02N, 35:175-35:178, 35:180-35:184, 35:186-35:195, 35:197-35:198, 35:200, 35:202-35:204, 35:206-35:207, 35:209-35:214, 35:216N-35:229, 35:231-35:232, 35:235-35:241, 35:243-35:250, 35:252-35:259, 35:261-35:264, 35:267-35:275, 35:277-35:316, 35:318-35:320, 35:322-35:330, 35:332-35:345, 35:347-35:371, 35:373-35:408, 35:410-35:413, 35:415-35:424, 35:426-35:470, 35:472-35:473, 35:477-35:478, 40:28, 44:10, 44:15:01:01-44:15:01:02, 44:18, 44:48, 44:256, 44:276, 44:339-44:340, 44:481, 45:01:01:01-45:25, 46:11, 46:18, 46:33, 46:53, 46:67, 48:02:01-48:02:03, 48:25, 49:01:01:01-49:14, 49:16-49:39, 49:41-49:60N, 49:62-49:67, 50:01:01:01-50:01:01:06, 50:01:03-50:01:13, 50:01:15-50:02:01:02, 50:05-50:19, 50:31-50:50, 50:52-50:59, 50:62-50:69, 51:15, 51:56:01-51:56:03, 51:157, 51:158:02, 51:225, 51:300, 53:01:01:01-53:06,

101.570-06 – including *Taq* polymerase, IFU-01
101.570-06u – without *Taq* polymerase, IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **4L7**

Lot-specific information

				53:08:01-53:43:02, 53:45-53:59, 54:01:01:01-54:05N, 54:07-54:08N, 54:10, 54:12-54:14, 54:16-54:36, 54:39-54:41, 56:01:01:01-56:01:03, 56:01:05-56:02:03, 56:07-56:11, 56:13-56:17, 56:19N, 56:22-56:26, 56:28N-56:30, 56:33-56:47, 56:50-56:56, 56:58-56:62, 56:64-56:70, 56:72-56:73, 57:05, 58:01:01:01-58:01:02, 58:01:04-58:02:02, 58:04-58:07, 58:09, 58:11-58:13, 58:15-58:19, 58:21-58:26, 58:28:01-58:29, 58:31N-58:45:02, 58:47-58:58, 58:60, 58:62-58:87, 58:89-58:113, 58:115-58:123, 59:01:01:01-59:01:01:02, 59:03-59:10N, 82:01:01:01-82:03, C*02:56, C*04:180:01-04:180:02, C*05:217, C*06:20, C*07:81, C*07:168, C*07:450, C*08:123, C*08:139, C*15:126, C*16:98, C*16:102, C*16:110
3	150 bp	1070 bp	*55:01:01:01-55:01:24, 55:03, 55:05, 55:09 ^w , 55:11, 55:15, 55:24-55:25, 55:29, 55:31, 55:33, 55:36, 55:38, 55:40, 55:44-55:45, 55:52 ^w , 55:53-55:55N, 55:59-55:60, 55:64, 55:66, 55:73-55:76, 55:78-55:79, 55:84-55:85, 55:87, 55:90-55:95, 55:97N-55:106	*07:78 ^w , 13:16, 13:20, 13:48, 13:128 ^w , 15:42, 51:62, 51:106:01-51:106:02, 51:223, 51:225, 51:300, 52:25:01-52:25:02, 54:20 ^w , 56:25, 56:37, 56:66, C*01:90^w, C*08:207^w, C*15:91, C*15:128
4	235 bp	1070 bp	*55:03, 55:49	*08:15, 15:289, 35:74, 35:186, 40:73, 46:01:01:01, 46:01:02-46:01:09, 46:02-46:04, 46:06-46:40:01, 46:41N-46:43, 46:52, 46:66
5⁴	95 bp	800 bp	*55:50	*08:83, 13:06, 13:53, 13:131, 14:43, 14:50, 18:74N, 18:78, 38:38, 41:37, A*31:147, C*05:110, C*12:120
	235 bp		*55:04, 55:49, 55:86	*07:19, 07:31, 07:34, 07:43, 08:01:01:01-08:05, 08:07-08:08N, 08:10, 08:13-08:15, 08:17-08:20:02, 08:22-08:24, 08:26:01-08:39, 08:41-08:48, 08:50-08:59:02, 08:61-08:78, 08:80-08:83, 08:85-08:88, 08:90-08:106, 08:108-08:124, 08:126-08:128, 08:130-08:155, 08:157-08:164, 08:166-08:170, 08:172-08:174, 08:176-08:180, 08:182-08:183, 08:185, 08:187-08:204, 08:206-08:231, 08:233-08:253, 14:05, 14:13, 14:53, 18:14, 35:87, 37:09, 38:19, 39:03:01:01-39:03:01:03, 39:14:01:01-39:14:01:03, 39:24:01-39:24:02, 39:29, 39:37:01:01-39:37:01:02, 39:76, 39:120, 39:144, 39:152, 40:39, 41:02:01:01-41:02:09, 41:04, 41:10-41:11, 41:13, 41:15, 41:18, 41:23-41:24, 41:27, 41:31, 41:36, 41:38-41:43, 41:45N-41:47, 41:49, 41:51-41:52, 41:55, 41:58, 42:01:01:01-42:02:02, 42:05:01-42:12, 42:14-42:28, 44:166, 48:37, C*06:264, C*07:294, C*07:526:01, C*07:530, C*07:681, C*08:08:01-08:08:02

101.570-06 – including *Taq* polymerase, IFU-01
101.570-06u – without *Taq* polymerase, IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **4L7**

Lot-specific information

6^{4,5}	100 bp	1070 bp	*55:05	*15:390, C*03:349, C*16:152
7	180 bp	1070 bp	*55:08, 55:51	*07:20, 07:24, 07:60, 07:100, 07:131, 07:220, 07:223, 07:358, 08:21, 08:25, 08:175, 13:62, 13:86, 15:09:01-15:10:05, 15:30:01:01-15:30:01:03, 15:37, 15:45, 15:48, 15:63, 15:90, 15:99, 15:150, 15:243, 15:248, 15:287, 15:312-15:313, 15:324, 15:329, 15:338, 15:361, 15:427, 15:435, 15:439, 15:445, 15:455, 15:503, 15:506, 15:530, 15:541, 27:76, 27:153, 35:02:01:01-35:02:16, 35:04:01:01-35:04:03, 35:09:01-35:09:03, 35:12:01:01-35:12:04, 35:18, 35:22:01:01-35:22:01:02, 35:34, 35:39, 35:44, 35:81, 35:83, 35:95, 35:129N, 35:149, 35:154, 35:162, 35:172, 35:182-35:184, 35:199, 35:201, 35:211, 35:220, 35:233, 35:251, 35:258, 35:266, 35:270, 35:273, 35:285, 35:292, 35:309, 35:311, 35:316, 35:321, 35:323, 35:335, 35:339, 35:357, 35:361, 35:366-35:367, 35:372, 35:374, 35:377-35:379, 35:384, 35:387, 35:391, 35:396, 35:403, 35:407, 35:410, 35:419, 35:425, 35:443, 35:463, 35:474, 35:477, 40:05:01:01-40:05:01:03, 40:26, 40:28, 40:63, 40:92, 40:174, 40:324, 40:343, 41:17, 41:62, 44:62, 44:77, 44:82, 44:107, 44:123, 44:184, 44:213, 44:459, 46:74, 46:81, 49:18:01-49:18:02, 50:14, 51:01:01:01-51:01:14, 51:01:16-51:01:31, 51:01:33-51:01:38, 51:01:40-51:01:58, 51:01:60-51:09:03, 51:11N-51:12, 51:18-51:20, 51:22-51:24:05, 51:26-51:30, 51:32-51:33, 51:35, 51:38-51:41N, 51:43-51:44N, 51:46, 51:48-51:61:02, 51:64-51:80, 51:83-51:84, 51:86-51:91, 51:94-51:96, 51:98N-51:100, 51:102-51:105, 51:107-51:142, 51:145-51:152, 51:154-51:156, 51:158:01-51:171, 51:173Q-51:181, 51:183-51:191, 51:193-51:196, 51:198-51:213, 51:215, 51:217-51:219, 51:221-51:222, 51:224, 51:226, 51:228-51:250, 51:252-51:260, 51:262-51:273N, 51:275-51:279, 51:282-51:294, 51:296-51:297, 51:299, 51:301-51:315, 52:01:01:01-52:01:20, 52:01:22-52:07, 52:09:01-52:13, 52:15-52:24, 52:27-52:36, 52:38, 52:40-52:61, 52:64-52:66, 52:68-52:70, 52:72-52:74, 52:76-52:83, 52:85-52:87, 52:89N-52:96N, 53:19, 53:36, 53:42, 56:05:01-56:06, 56:15, 56:21, 57:02:01-57:03:05, 57:05, 57:07, 57:09, 57:12, 57:17, 57:28N, 57:39, 57:42, 57:46, 57:57, 57:63, 57:66, 57:70, 57:80, 57:82, 57:84, 57:94, 57:96, 57:101, 57:124, 58:08:01-58:08:02, 58:28:01-58:28:02, 58:65, 78:01:01:02-78:06, 78:09-78:10, C*03:81, C*03:152, C*03:317, C*03:382, C*03:388
8	230 bp 425 bp	800 bp	*55:66 *55:12	*56:41
9	215 bp 375 bp	1070 bp	*55:25 *55:11, 55:14	*07:225, 07:227

101.570-06 – including *Taq* polymerase, IFU-01
101.570-06u – without *Taq* polymerase, IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **4L7**

Lot-specific information

10^{4,6}	70 bp	1070 bp	*55:07	*07:13, 07:247, 07:270, 08:188, 15:436, 27:198, 39:121, 54:01:02, 59:06, 67:02:01:01-67:02:01:02, C*01:02:48, C*03:02:18, C*03:03:42, C*03:04:50, C*07:01:22, C*07:02:31, C*08:16:02, C*14:02:27, C*15:02:06
	220 bp		*55:17, 55:38, 55:72	*08:09, 08:84, 15:83:01:01-15:83:01:02, 39:06:02:01-39:06:06, 39:34, 39:57, 39:62, 39:64, 39:90, 39:109, 39:129, 39:132, 39:140-39:141, 39:154, 42:04, 51:21, 51:36, 51:101, 51:143, 51:153, 51:197, 51:280, 51:295, 52:37, 52:71, 52:84, 78:07, C*14:92
11⁴	90 bp	1070 bp	*55:48	*07:78, 13:01:01:01-13:02:32, 13:07N-13:09, 13:11, 13:14-13:20, 13:22:01-13:23, 13:25, 13:27-13:34, 13:36-13:45, 13:47, 13:49N-13:50, 13:52, 13:54-13:58, 13:60-13:61, 13:63N-13:70, 13:73-13:78, 13:80-13:85, 13:87-13:107, 13:109-13:119, 13:121-13:126, 13:128-13:130, 13:132-13:139N, 13:141, 40:48, 40:377, 45:10, 49:07, 50:31, 54:26, 56:60, 59:09
	125 bp 180 bp		*55:15 *55:09, 55:22, 55:24	
12	210 bp	1070 bp	*55:10, 55:82	*07:02:01:01-07:02:27, 07:02:29-07:02:81, 07:04:01-07:07:03, 07:09:01-07:09:02, 07:11-07:12, 07:14-07:15, 07:17-07:26, 07:28, 07:30-07:31, 07:33:01-07:36, 07:39-07:49N, 07:51-07:64, 07:66-07:68:03, 07:73-07:82, 07:84, 07:87-07:124, 07:126-07:142, 07:144-07:179, 07:181N-07:182N, 07:185, 07:187-07:203, 07:205-07:206, 07:208-07:209, 07:211-07:218, 07:220-07:227, 07:229:01-07:258, 07:260-07:291, 07:293-07:305, 07:307-07:312, 07:313 ^w , 07:314-07:351N, 07:353-07:367, 07:369-07:382, 07:384-07:399, 07:401, 38:26, 42:01:01:01-42:01:04, 42:04-42:06, 42:08, 42:10, 42:12-42:16, 42:19-42:25, 42:28, 67:01:01-67:01:03, 67:03-67:07, 81:01:01:01-81:04N, 81:06-81:09
13	150 bp	1070 bp	*55:02:01:01-55:02:14, 55:07, 55:10, 55:12, 55:16, 55:18-55:19, 55:22, 55:26, 55:30, 55:34-55:35, 55:37 ^w , 55:39, 55:41-55:43, 55:47-55:48, 55:50, 55:54, 55:61-55:63, 55:65, 55:67, 55:69-55:72, 55:77, 55:80, 55:82-55:83N, 55:88-55:89N, 55:96, 55:107	*13:01:01:01-13:03, 13:06-13:09, 13:11-13:12:02, 13:14-13:15, 13:17, 13:18 ^w , 13:19, 13:22:01-13:23, 13:25, 13:27-13:30, 13:32-13:34, 13:36-13:40, 13:42-13:45, 13:47, 13:49N-13:50, 13:52-13:58, 13:60-13:61, 13:63N-13:70, 13:73 ^w , 13:74-13:78, 13:80-13:85, 13:87-13:88, 13:90-13:107, 13:109, 13:111-13:126, 13:129-13:130, 13:132-13:139N, 13:141, 15:86, 15:224, 15:485, 15:552, 18:191, 35:60, 35:432, 35:446, 40:48, 40:71, 40:377, 44:10, 44:339-44:340, 46:11, 46:18, 49:11, 49:55, 51:15, 51:157, 54:01:01:01-54:03, 54:05N, 54:07-54:08N, 54:10 ^w , 54:12-54:13, 54:17-54:19, 54:21-54:32, 54:34-54:39, 54:41, 56:01:01:01-56:01:04, 56:01:06-56:02:03, 56:04:01:01-56:04:04, 56:07-56:08, 56:10, 56:14, 56:16-56:17, 56:19N-56:20:02, 56:23-56:24, 56:26-56:29, 56:33-56:36, 56:38N-56:42, 56:44-56:49, 56:51 ^w , 56:52-56:58, 56:60-56:62, 56:63 ^w ,

101.570-06 – including *Taq* polymerase, IFU-01
101.570-06u – without *Taq* polymerase, IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **4L7**

Lot-specific information

				56:64-56:65, 56:67-56:73, 58:73, 59:01:01:01-59:01:01:02, 59:04-59:10N
14⁴	95 bp	1070 bp	*55:07	*07:247, 07:270, 54:01:01:01-54:01:06, 54:01:08-54:23, 54:25 ^w , 54:26-54:34, 54:36-54:41
15	245 bp	1070 bp	*55:01:01:01-55:03, 55:05, 55:07, 55:10-55:13, 55:15-55:21, 55:23, 55:25-55:48, 55:50, 55:52-55:85, 55:87-55:107	*08:09, 08:84, 08:254, 13:46, 13:53, 13:131, 14:37, 15:42, 15:83:01:01-15:83:01:02, 15:485, 35:60, 35:276, 39:06:01-39:06:06, 39:34, 39:57, 39:62, 39:64, 39:83, 39:90, 39:109, 39:127, 39:129, 39:132, 39:140-39:141, 39:154, 41:01:01:01-41:01:07, 41:05-41:07, 41:09, 41:12, 41:14, 41:16, 41:20-41:22, 41:25-41:26, 41:28-41:30, 41:32-41:35, 41:37, 41:44, 41:50, 41:53-41:54, 41:57, 41:59-41:61, 41:64, 42:04, 45:08, 46:18, 50:45, 50:49, 51:21, 51:36, 51:101, 51:143, 51:153, 51:157, 51:197, 51:280, 51:295, 51:300, 52:37, 52:71, 52:84, 54:01:01:01-54:02, 54:04-54:05N, 54:07-54:08N, 54:10-54:25, 54:27-54:32, 54:34-54:37, 54:39, 54:41, 56:23, 56:54, 59:01:01:01-59:03, 59:05-59:08, 59:10N, 78:07, C*05:152, C*07:463, C*07:758, C*12:183, C*14:92
16	245 bp	1070 bp		*13:03-13:04, 13:48, 13:62, 13:120, 15:04:01:01-15:04:04, 15:16:01:01-15:16:04, 15:67, 15:95, 15:137, 15:155, 15:222, 15:254, 15:293, 15:303, 15:310, 15:362, 15:408, 15:499, 15:517, 15:531, 15:550, 15:552, 27:153, 35:37, 35:235, 35:265, 35:273, 39:33, 41:17, 44:15:01:01-44:15:01:02, 44:18, 44:20, 44:47, 44:100, 44:197, 44:339, 44:459, 45:01:01:01-45:07, 45:09, 45:11-45:23, 45:24 ^w , 45:25, 46:11, 46:32, 49:01:01:01-49:01:07, 49:01:09-49:06, 49:08-49:14, 49:16-49:46, 49:48-49:60N, 49:62-49:67, 50:01:01:01-50:02:01:02, 50:04-50:11, 50:13-50:16, 50:18-50:20, 50:32-50:44, 50:46-50:48, 50:50, 50:52-50:69, 51:01:01:01-51:01:31, 51:01:33-51:01:38, 51:01:40-51:01:58, 51:01:60-51:03, 51:05, 51:07:01-51:09:03, 51:11N-51:15, 51:17-51:20, 51:22-51:24:05, 51:26-51:30, 51:32-51:33, 51:35, 51:37-51:41N, 51:43-51:44N, 51:48-51:55, 51:57-51:58, 51:60-51:61:02, 51:63:01-51:63:02, 51:65-51:80, 51:83-51:92:02, 51:94-51:100, 51:102-51:134, 51:136-51:138, 51:140-51:142, 51:145-51:147, 51:149N-51:152, 51:154-51:156, 51:158:01-51:193, 51:195-51:196, 51:198-51:212, 51:214-51:219, 51:221-51:226, 51:228-51:279, 51:282-51:292, 51:294, 51:296-51:299, 51:301-51:309, 51:311-51:315, 52:01:01:01-52:01:20, 52:01:22-52:07, 52:09:01-52:36, 52:38-52:66, 52:68-52:70, 52:72-52:83, 52:85-52:87, 52:89N-52:96N, 54:03, 54:33, 54:40, 56:01:01:01-56:01:12, 56:01:14-56:01:15, 56:05:01-56:08, 56:13-56:17, 56:20:01-56:22, 56:24-56:30, 56:33-56:50, 56:52-56:53, 56:55:01:02, 56:58-56:59, 56:61-56:62, 56:64-56:65, 56:67-56:73,

101.570-06 – including *Taq* polymerase, IFU-01
101.570-06u – without *Taq* polymerase, IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **4L7**

Lot-specific information

				58:08:01-58:08:02, 58:122, 59:04, 78:01:01:02-78:06, 78:09-78:10, C*03:464
17⁴	105 bp	1070 bp	*55:17, 55:20, 55:27-55:28, 55:39, 55:81	*07:84, 08:09, 08:84, 08:254, 13:35, 13:59, 13:62, 15:83:01:01-15:83:01:02, 27:14, 27:81, 27:130, 27:153, 35:273, 40:06:01:01-40:06:08, 40:06:10-40:06:17, 40:06:19-40:06:22, 40:44, 40:53, 40:70:01-40:70:02, 40:75, 40:83, 40:93, 40:95-40:96, 40:103, 40:109-40:110, 40:127, 40:131, 40:148, 40:161-40:162, 40:165, 40:167, 40:177, 40:190, 40:230, 40:244, 40:256N, 40:268-40:270, 40:275, 40:294-40:295, 40:298:01-40:298:02, 40:306-40:307, 40:311, 40:313, 40:317-40:318, 40:340-40:342, 40:348-40:350, 40:354, 40:361N-40:364, 40:366, 40:373-40:374, 40:376, 40:387-40:389, 40:408, 40:411, 40:423, 40:425, 40:427, 40:434, 40:445, 40:449, 40:451, 41:01:01:01-41:01:07, 41:05-41:07, 41:09, 41:12, 41:14, 41:16-41:17, 41:20-41:22, 41:25-41:26, 41:28-41:29, 41:32-41:35, 41:37, 41:53-41:54, 41:56Q-41:57, 41:59-41:61, 41:64, 42:04, 44:20, 44:47, 44:100, 44:197, 44:459, 49:18:01-49:18:02, 50:14, 50:20, 51:01:01:01-51:01:47, 51:01:49-51:03, 51:05, 51:07:01-51:12, 51:14, 51:16-51:24:05, 51:26-51:34, 51:36, 51:38-51:41N, 51:43-51:44N, 51:48-51:55, 51:57-51:58, 51:60-51:61:02, 51:65-51:80, 51:82-51:91, 51:93-51:96, 51:98N, 51:100-51:105, 51:107-51:130, 51:132, 51:134, 51:136-51:138, 51:140-51:147, 51:149N-51:156, 51:158:01-51:171, 51:173Q-51:193, 51:195-51:212, 51:214-51:215, 51:217-51:219, 51:221-51:222, 51:224, 51:226, 51:228-51:250, 51:252-51:280, 51:282-51:292, 51:294-51:299, 51:301-51:309, 51:311-51:315, 52:01:01:01-52:13, 52:15-52:24, 52:26-52:38, 52:40-52:74, 52:76-52:87, 52:89N-52:96N, 54:04, 54:11, 54:15, 56:05:01-56:06, 56:15, 56:21, 58:08:01-58:08:02, 59:02, 78:01:01:02-78:10, C*03:464, C*14:92
18	215 bp 310 bp 415 bp	800 bp	*55:25 *55:19 *55:18, 55:31	*13:125, 18:163, 40:10:01:01-40:10:02, 44:344, 54:30, 56:23, 56:42, 56:44, 56:54, 56:67, 59:01:01:01-59:10N
19	205 bp	1070 bp	*55:30, 55:36	*07:234, 15:244, 15:315, 18:126, 39:54, 39:138, 50:56, C*02:131, C*03:70, C*03:127, C*03:179, C*03:274
20	220 bp	1070 bp	*55:20, 55:56	*08:09, 08:84, 08:254, 13:46, 15:83:01:01-15:83:01:02, 41:01:01:01-41:01:07, 41:05-41:07, 41:09, 41:12, 41:14, 41:16-41:17, 41:20-41:22, 41:25-41:26, 41:28-41:30, 41:32-41:35, 41:37, 41:44, 41:50, 41:53-41:54, 41:56Q-41:57, 41:59-41:61, 41:64, 42:04, 44:15:01:01-44:15:01:02, 44:18, 44:20, 44:100, 45:01:01:01-45:25, 51:08:01:01-51:08:04, 51:20, 51:36, 51:44N,

101.570-06 – including *Taq* polymerase, IFU-01
101.570-06u – without *Taq* polymerase, IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **4L7**

Lot-specific information

				51:97, 51:141, 51:153, 51:263, 51:274-51:275, 52:19, 56:13
21	170 bp 270 bp	800 bp	*55:83N *55:10	*07:209, 07:271, 15:506, 18:01:01:01-18:03:02, 18:05:01:01-18:15, 18:18:01:01-18:28, 18:30-18:79, 18:81-18:97, 18:99-18:101, 18:103-18:106Q, 18:108-18:111, 18:113-18:119, 18:121-18:124, 18:126-18:192, 27:75, 35:25, 35:124-35:125, 35:142, 35:162, 35:267, 35:320, 37:01:01:01-37:21, 37:23-37:40, 37:42N-37:91, 40:149 [?] , 44:130, 44:156, 44:221, 44:224, 44:258, 51:116, 51:176, 51:179, 51:218, 52:16, 52:27, 52:77, 53:43:01-53:43:02, 54:21, 56:16, 57:58, 57:80, 58:12, 81:01:01:01-81:01:05, 81:03 [?] -81:04N [?] , 81:05, 81:06 [?] -81:08 [?] , 82:01:01:01-82:03, 83:01
22⁴	95 bp 225 bp	800 bp	*55:23, 55:27, 55:68 *55:09, 55:21, 55:37, 55:52	*07:84, 13:04, 13:35, 13:72, 13:127, 15:04:01:01-15:04:04, 15:16:01:01-15:16:04, 15:67, 15:95, 15:155, 15:222, 15:254, 15:293, 15:310, 15:362, 15:408, 15:499, 15:517, 15:531, 15:550, 35:37, 35:235, 35:265, 40:44, 40:159, 41:21, 44:20, 44:47, 44:100, 44:197, 45:09, 46:32, 49:04:01-49:05:02, 50:33, 51:37, 51:90, 51:92:01-51:92:02, 51:197, 51:216, 51:261, 51:274, 51:298, 52:12, 52:26, 52:63, 52:75, 54:14-54:15, 56:50, 58:122, C*03:464, C*07:463, C*07:758, C*12:183 *07:78, 07:84, 13:18, 13:31, 13:41, 13:128, 13:131, 15:04:01:01-15:04:04, 15:137, 15:303, 15:499, 15:531, 35:265, 39:33, 40:161, 46:32, 49:41, 50:47, 51:05, 51:29, 51:54, 51:61:01-51:61:02, 51:82, 51:197, 51:217, 51:279, 52:21:01-52:21:02, 52:24, 54:10, 54:20, 54:33, 54:40, 56:43, C*05:152, C*12:183, C*14:92
23	165 bp 220 bp	1070 bp	*55:29 *55:16, 55:21	*13:105, 51:174:02 *07:143, 13:31, 13:41, 13:131, 15:04:01:01-15:04:04, 15:137, 15:303, 15:499, 15:531, 27:157, 35:265, 40:166, 42:02:01:01-42:02:02, 42:09, 42:17-42:18, 42:27, 45:06, 46:32, 50:47, 51:61:01-51:61:02, 52:21:01-52:21:02, 54:02, 54:33, 54:40, 56:35, 56:43, C*12:183
24⁴	100 bp	1070 bp	*55:13, 55:23, 55:27, 55:32, 55:58, 55:68	*07:84, 13:04, 13:35, 13:46, 13:72, 13:127, 13:140, 15:04:01:01-15:04:04, 15:16:01:01-15:16:04, 15:67, 15:95, 15:137, 15:155, 15:222, 15:254, 15:293, 15:310, 15:362, 15:408, 15:499, 15:517, 15:531, 15:550, 35:37, 35:235, 35:265, 35:276, 39:06:01-39:06:06, 39:33-39:34, 39:50, 39:57, 39:62, 39:64, 39:83, 39:90, 39:102, 39:109, 39:127, 39:129, 39:132, 39:140-39:141, 39:154, , 40:44, 40:86, 40:159, 40:330, 41:21, 41:30, 41:50, 44:20, 44:47, 44:100, 44:197, 45:02, 45:09, 46:32, 49:04:01-49:05:02, 49:50, 50:09, 50:33, 51:13:01-51:13:02, 51:37, 51:63:01-51:63:02, 51:90, 51:92:01-51:92:02, 51:97, 51:172, 51:197, 51:216, 51:251, 51:261, 51:274, 51:298, 52:12,

101.570-06 – including *Taq* polymerase, IFU-01
101.570-06u – without *Taq* polymerase, IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **4L7**

Lot-specific information

	235 bp		*55:33	52:14:01-52:14:02, 52:26, 52:39, 52:63, 52:75, 54:14-54:15, 56:22, 56:50, 58:122, 59:03, 73:01:01:01-73:02, C*03:464, C*05:152, C*07:463, C*07:758, C*12:183
25⁵	155 bp 280 bp 380 bp	1070 bp	*55:62 *55:26 *55:18, 55:34, 55:40	*13:125, 27:196, 40:10:01:01-40:10:02, 44:344, 54:35, 56:67
26^{4,5}	75 bp 310 bp 400 bp	1070 bp	*55:42 *55:45 *55:34, 55:47	*13:88
27⁴	95 bp 170 bp	1070 bp	*55:46 *55:43	*50:39, 54:16 *51:127
28^{4,5}	95 bp	800 bp	*55:41, 55:44, 55:58	*13:46, 15:137, 49:54, 51:172, 51:223, 73:01:01:01, 73:02
	215 bp		*55:57	
29	160 bp	1070 bp	*55:83N	*08:101, 13:132, 15:434, 41:55, 44:24, 44:181, 51:39, 54:01:01:01-54:01:08, 54:03-54:04, 54:06- 54:38, 54:40-54:41, C*05:30, C*06:125, C*07:441:01-07:441:02, C*16:86, C*16:144, C*16:147
30⁴	100 bp	1070 bp	*55:35, 55:50	*08:83, 13:06, 13:53, 13:131, 14:43, 14:50, 15:76, 15:101, 18:74N, 18:78, 38:38, 41:37, 56:39, A*31:147, C*05:110, C*12:120
31⁴	105 bp	800 bp	*55:13, 55:58	*13:46, 13:140, 15:137, 35:235, 35:276, 39:06:01- 39:06:06, 39:33-39:34, 39:50, 39:57, 39:62, 39:64, 39:83, 39:90, 39:102, 39:109, 39:127, 39:129, 39:132, 39:140-39:141, 39:154, 40:86, 40:330, 41:50, 45:02, 49:50, 50:09, 51:13:01-51:13:02, 51:92:01-51:92:02, 51:172, 52:14:01-52:14:02, 52:75, 56:22, 59:03, 73:01:01:01-73:02, C*05:152 *51:127
	170 bp 235 bp		*55:43 *55:55N	
32⁷			Negative Control	

¹Alleles are assigned by the presence of specific PCR product(s). However, the sizes of the specific PCR products may be helpful in the interpretation of HLA-B*55 SSP typings.

When the primers in a primer mix can give rise to HLA-specific PCR products of more than one length this is indicated if the size difference is more than 20 base pairs. Size differences of 20 base pairs or less are not given. For high resolution SSP kits, the alleles listed are specified according to amplicon length.

Nonspecific amplifications, i.e. a ladder or a smear of bands, may sometimes be seen. GC-rich primers have a higher tendency of giving rise to nonspecific amplifications than other primers.

PCR fragments longer than the control bands may sometimes be observed. Such bands should be disregarded and do not influence the interpretation of the SSP typings.

PCR fragments migrating faster than the control bands, but slower than a 400 bp fragment may be seen in some gel read-outs. Such bands can be disregarded and do not influence the interpretation of the SSP typings.

Some primers may give rise to primer oligomer artifacts. Sometimes this phenomenon is an inherent feature of the primer pair(s) of a primer mix. More often it is due to other factors such as too low amount of DNA in the PCR reactions, taking too long time in setting up the PCR reactions, working at elevated room temperature or using thermal cyclers that are not pre-heated.

101.570-06 – including *Taq* polymerase, IFU-01
101.570-06u – without *Taq* polymerase, IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **4L7**

Lot-specific information

²The internal positive control primer pairs amplify segments of the human growth hormone gene. The internal positive control bands are 1070 or 800 base pairs respectively, well distribution as outlined in the table. Well number 1 contains the shorter, 800 bp, internal positive control band. The well distribution of the internal controls can help in orientation of the kit on gel photo, as well as allow for kit identification. In the presence of a specific amplification the intensity of the control band often decreases.

³For several HLA Class I alleles 1st and/or 4th exon(s) and beyond, as well as intron nucleotide sequences, are not available. In these instances it is not known whether some of the primers of the SSP sets are completely matched with the target sequences or not. Assumption is made that unknown sequences in these regions are conserved within allelic groups.

⁴HLA-specific PCR products shorter than 125 base pairs have a lower intensity and are less sharp than longer PCR products.

⁵Primer mixes 6, 25, 26 and 28 may have tendencies of unspecific amplifications.

⁶Primer mix 10 may give rise to a long fragment of approximately 650 bp in some HLA-B alleles. This band should not be considered in the interpretation of HLA-B*55 typings.

⁷Primer mix 32 contains a negative control, which will amplify a majority of HLA amplicons as well as the amplicons generated by the control primer pairs matching the human growth hormone gene. HLA-specific PCR product sizes range from 75 to 200 base pairs and the PCR product generated by the HGH positive control primer pair is 200 base pairs.

Abbreviations

‘?’ , nucleotide sequence of the primer matching region is not known.

‘w’ , may be weakly amplified.

101.570-06 – including *Taq* polymerase, IFU-01
101.570-06u – without *Taq* polymerase, IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **4L7**

Lot-specific information

PRIMER SPECIFICATION

Well No.	1	2	3	4	5	6	7	8	9	10	11	12
Length of spec. PCR product	400	125	150	235	95	100	180	230	215	70	90	210
					235			425	375	220	125	
											180	
Length of int. pos. control ¹	800	1070	1070	1070	800	1070	1070	800	1070	1070	1070	1070
5'-primer(s) ²	41	379	420	106	363	245	419	48	48	175	186	103
	5'-CTg 3'	5'-ACC 3'	5'-TTA 3'	5'-CCA 3'	5'-AgC 3'	5'-ggC 3'	5'-gTA 3'	5'-gCC 3'	5'-gCC 3'	5'-TgA 3'	5'-AgA 3'	5'-CCT 3'
		379			506				272	379	420	
		5'-ACC 3'			5'-gCT 3'				5'-CTA 3'	5'-ACg 3'	5'-TTA 3'	
3'-primer(s) ³	272	463	527	299	559	302	559	107	94	206	272	272
	5'-TgT 3'	5'-gCT 3'	5'-CCT 3'	5'-TCA 3'	5'-CgT 3'	5'-ggC 3'	5'-CAG 3'	5'-ACT 3'	5'-gAC 3'	5'-CCC 3'	5'-TgT 3'	5'-TgT 3'
	275							302	362	559	469	
	5'-CCA 3'							5'-ggT 3'	5'-TCA 3'	5'-CgT 3'	5'-CCg 3'	
											559	
											5'-CTC 3'	
Well No.	1	2	3	4	5	6	7	8	9	10	11	12

Well No.	13	14	15	16	17	18	19	20	21	22	23	24
Length of spec. PCR product	150	95	245	245	105	215	205	220	170	95	165	100
						310			270	225	220	235
						415						
Length of int. pos. control ¹	1070	1070	1070	1070	1070	800	1070	1070	800	800	1070	1070
5'-primer(s) ²	420	206	357	357	357	48	134	357	41	357	97	357
	5'-TTA 3'	5'-Agg 3'	5'-Tgg 3'	5'-Tgg 3'	5'-Tgg 3'	5'-gCC 3'	5'-CCA 3'	5'-Tgg 3'	5'-CTg 3'	5'-Tgg 3'	5'-TCC 3'	5'-Tgg 3'
			357				141		668		357	
			5'-Tgg 3'				5'-ATg 3'		5'-ggC 3'		5'-Tgg 3'	
3'-primer(s) ³	523	259	559	559	419	94	299	538	142	412	272	412
	5'-ACA 3'	5'-gTT 3'	5'-CgT 3'	5'-CAG 3'	5'-CgT 3'	5'-gAC 3'	5'-TCT 3'	5'-gTC 3'	5'-TgA 3'	5'-gTC 3'	5'-TgT 3'	5'-gTC 3'
	527				424	187			799	539	479	419
	5'-CCA 3'				5'-gTg 3'	5'-gTT 3'			5'-CTA 3'	5'-TCC 3'	5'-CCA 3'	5'-Cgg 3'
						292					538	419
						5'-gTA 3'					5'-CCA 3'	5'-CgA 3'
						292						550
						5'-gTA 3'						5'-CAT 3'
						292						
						5'-gTA 3'						
Well No.	13	14	15	16	17	18	19	20	21	22	23	24

101.570-06 – including *Taq* polymerase, IFU-01
101.570-06u – without *Taq* polymerase, IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **4L7**

Lot-specific information

Well No.	25	26	27	28	29	30	31
Length of spec.	155	75	95	95	160	100	105
PCR product	280	310	170	215			170
	380	400					235
Length of int. pos. control ¹	1070	1070	1070	800	1070	1070	800
5'-primer(s) ²	48	41	357	362	106	209	357
	5'-gCC 3'	5'-CTg 3'	5'-Tgg 3'	5'-gAC 3'	5'-CCA 3'	5'-ggC 3'	5'-Tgg 3'
	3 rd I	431			668	506	
	5'-Agg 3'	5'-CgA 3'			5'-ggC 3'	5'-gCT 3'	
	160	181	412	409	226	272	419
	5'-gTT 3'	5'-gTT 3'	5'-gTC 3'	5'-ATA 3'	5'-CAC 3'	5'-TgT 3'	5'-CgA 3'
3'-primer(s) ³	259	272	487	427	799	559	487
	5'-CTC 3'	5'-Tgg 3'	5'-CgT 3'	5'-gTg 3'	5'-CTA 3'	5'-CgT 3'	5'-CgT 3'
	259	463		535			549
	5'-CTC 3'	5'-gCT 3'		5'-CTC 3'			5'-AgT 3'
	682						
	5'-CCT 3'						
Well No.	25	26	27	28	29	30	31

¹The internal positive control primer pairs amplify segments of the human growth hormone gene. The internal positive control bands are 1070 or 800 base pairs respectively, well distribution as outlined in the table. Well number 1 contains the shorter, 800 bp, internal positive control band. The well distribution of the internal controls can help in orientation of the kit on gel photo, as well as allow for kit identification. In the presence of a specific amplification the intensity of the control band often decreases.

²The nucleotide position matching the specificity-determining 3'-end of the primer is given. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

³The nucleotide position matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

101.570-06 – including *Taq* polymerase, IFU-01
101.570-06u – without *Taq* polymerase, IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **4L7**

Lot-specific information

CELL LINE VALIDATION SHEET																				
HLA-B*55 SSP subtyping kit ²																				
					Well															
					1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
				Prod. No.:	201777501	201777502	201777503	201777504	201777505	201777506	201777507	201777508	201777509	202023210	201777511	201777512	201777513	201777514	202023215	201777516
	IHCW cell line ¹		B*																	
1	9001	SA	*07:02		-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-
2	9280	LK707	*52:01	*73:01	-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	+
3	9011	E4181324	*52:01		-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	+
4	9275	GU373	*15:10	*53:01	-	+	-	-	-	-	+	-	-	-	-	-	-	-	-	-
5	9009	KAS011	*37:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
6	9353	SM	*39:01	*51:01	-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	+
7	9020	QBL	*18:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
8	9025	DEU	*35:01		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-
9	9026	YAR	*38:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
10	9107	LKT3	*54:01		+	+	-	-	-	-	-	-	-	-	-	-	+	+	+	-
11	9051	PITOUT	*44:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
12	9052	DBB	*57:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
13	9025	JESTHOM	*27:05		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
14	9071	OLGA	*15:01	*15:20	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-
15	9075	DKB	*40:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
16	9037	SWEIG007	*40:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
17	9282	CTM3953540	*08:01	*55:01	+	+	+	-	-	-	-	-	-	-	-	-	-	-	+	-
18	9257	32367	*14:01	*56:01	+	+	-	-	-	-	-	-	-	-	-	-	+	-	-	+
19	9038	BM16	*18:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
20	9059	SLE005	*40:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
21	9064	AMALA	*15:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
22	9056	KOSE	*35:03		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-
23	9124	IHL	*40:02	*56:02	+	+	-	-	-	-	-	-	-	-	-	-	+	-	-	-
24	9035	JBUSH	*38:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
25	9049	IBW9	*14:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
26	9285	WT49	*58:01		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-
27	9191	CH1007	*07:05	*51:01	-	-	-	-	-	-	+	-	-	-	-	+	-	-	-	+
28	9320	BEL5GB	*44:02	*44:03	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
29	9050	MOU	*44:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
30	9021	RSH	*42:01		-	-	-	-	+	-	-	-	-	-	-	+	-	-	-	-
31	9019	DUCAF	*18:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
32	9297	HAG	*41:02		-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-
33	9098	MT14B	*40:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
34	9104	DHIF	*38:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
35	9302	SSTO	*44:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
36	9024	KT17	*15:01	*35:01	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-
37	9065	HHKB	*07:02		-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-
38	9099	LZL	*15:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
39	9315	CML	*08:01	*27:05	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-
40	9134	WHONP199	*13:02	*46:01	-	+	-	+	-	-	-	-	-	-	+	-	+	-	-	-
41	9055	H0301	*14:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
42	9066	TAB089	*46:01		-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-
43	9076	T7526	*46:01		-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-
44	9057	TEM	*38:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
45	9239	SHJO	*42:01	*50:01	-	+	-	-	+	-	-	-	-	-	-	+	-	-	-	+
46	9013	SCHU	*07:02		-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-
47	9045	TUBO	*51:01		-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	+
48	9303	TER-ND	*35:01	*44:03	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-

101.570-06 – including *Taq* polymerase, IFU-01

Visit <https://labproducts.caredx.com> for

101.570-06u – without *Taq* polymerase, IFU-02

“Instructions for Use” (IFU)

Lot No.: 4L7

Lot-specific information

CELL LINE VALIDATION SHEET																															
HLA-B*55 SSP subtyping kit ²																															
					Prod. No.:	Well																									
						17	18	19	20	21	22	23	24	25	26	27	28	29	30	31											
						201777517	201777518	201777519	201777520	201897921	201777522	201777523	201777524	201777525	201777526	201777527	201777528	201897929	202023230	201777531											
	IHWC cell line ¹			B*																											
1	9001	SA	*07:02	*73:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-											
2	9280	LK707	*52:01		+	-	-	-	-	-	-	-	+	-	-	-	+	-	-	+											
3	9011	E4181324	*52:01		+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-											
4	9275	GU373	*15:10	*53:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-											
5	9009	KAS011	*37:01	*51:01	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-											
6	9353	SM	*39:01		+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-											
7	9020	QBL	*18:01		-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-											
8	9025	DEU	*35:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-											
9	9026	YAR	*38:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-											
10	9107	LKT3	*54:01		-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-											
11	9051	PITOUT	*44:03	*15:20	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-											
12	9052	DBB	*57:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-											
13	9025	JESTHOM	*27:05		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-											
14	9071	OLGA	*15:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-											
15	9075	DKB	*40:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-											
16	9037	SWEIG007	*40:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-											
17	9282	CTM3953540	*08:01		*55:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-											
18	9257	32367	*14:01		*56:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-											
19	9038	BM16	*18:01			-	-	-	-	+	-	-	-	-	-	-	-	-	-	-											
20	9059	SLE005	*40:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-											
21	9064	AMALA	*15:01	*56:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-												
22	9056	KOSE	*35:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-											
23	9124	IHL	*40:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-											
24	9035	JBUSH	*38:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-											
25	9049	IBW9	*14:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-											
26	9285	WT49	*58:01	*51:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-												
27	9191	CH1007	*07:05		+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-											
28	9320	BEL5GB	*44:02		*44:03	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-											
29	9050	MOU	*44:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-											
30	9021	RSH	*42:01		-	-																									

101.570-06 – including *Taq* polymerase, IFU-01
101.570-06u – without *Taq* polymerase, IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **4L7**

Lot-specific information

¹The provided cell line HLA specificities are retrieved from the <http://www.ihwg.org/hla> web site. The specificity of an individual cell line may thus be subject to change.

²The specificity of each primer solution in the kit has been tested against 48 well characterized cell line DNAs and where applicable, additional cell line DNAs.

No DNAs carrying the alleles to be amplified by primer solutions 6, 8 to 10, 19, 26, 27 and 30 were available. The specificities of the primers in primer solutions 8, 10, 19, 26, 27 and 30 were tested by separately adding one or two additional 5'-primers, respectively one or three additional 3'-primers.

In primer solution 6 it was only possible to test the 3'-primer, the 5'-primer was not possible to test.

In primer solution 9 it was only possible to test the 5'-primers, the 3'-primers were not possible to test.

In primer solutions 2, 11, 15, 19 and 26 one 5'-primer was not possible to test.

In primer solutions 1, 8, 11, 13, 17, 18, 21, 23 to 29 and 31 one or more 3'-primers were not possible to test.

Additional primers in primer solutions 5, 11, 18, 21, 25 and 29 were tested by separately adding additional 3'-primers and/or additional 5'-primers.

101.570-06 – including *Taq* polymerase, IFU-01
101.570-06u – without *Taq* polymerase, IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **4L7**

Lot-specific information

ADDRESSES:

Manufacturer:

CareDx AB, Franzengatan 5, SE-112 51 Stockholm, Sweden.

Tel: +46-8-508 939 00

Fax: +46-8-717 88 18

E-mail: orders-se@caredx.com

Web page: <https://labproducts.caredx.com/>

Distributed by:

CareDx GmbH, Löwengasse 47 / 6, AT-1030 Vienna, Austria.

Tel: +43-1-710 15 00

Fax: +43-1-710 15 00 10

E-mail: orders-at@caredx.com

Web page: <https://labproducts.caredx.com/>

CareDx Lab Solutions Inc., 901 S. Bolmar St., Suite R, West Chester, PA 19382

Tel: 1-877-653-78171

Fax: 610-344-7989

E-mail: orders-us@caredx.com

Web page: <https://labproducts.caredx.com/>

For information on CareDx distributors worldwide, contact **CareDx GmbH**.