

101.512-24/03– including *Taq* pol., IFU-01
101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: **92R**

Lot-specific information

Olerup SSP[®] HLA-B*07

Product number:	101.512-24/03 – including <i>Taq</i> pol. 101.512-24u/03u – without <i>Taq</i> pol.
Lot number:	92R
Expiry date:	2015-September-01
Number of tests:	24 tests – Product No. 101.512-24/24u 3 tests – Product No. 101.512-03/03u
Number of wells per test:	80
Storage - pre-aliquoted primers:	dark at -20°C
- PCR Master Mix:	-20°C
- Adhesive PCR seals	RT
- Product Insert	RT

This Product Description is only valid for Lot No. 92R.

CHANGES COMPARED TO THE PREVIOUS OLERUP SSP[®] HLA-B*07 LOT (32N)

The HLA-B*07 kit is updated for new alleles to enable separation of:

- Confirmed¹ alleles as listed in the IMGT/HLA database
- Polymorphisms in exons outside of the region encoding the peptide binding domain
- Null and Alternatively expressed alleles

One well has been added to the HLA-B*07 kit, well **80**.

¹As described in section Uniquely Identified Alleles.

The HLA-B*07 specificity and interpretation tables have been updated for the HLA-B alleles described since the previous *Olerup SSP[®]* HLA-B*07 lot was made (**Lot No. 32N**).

The primers of the wells detailed below have been exchanged, added or modified compared to the previous lot.

Well	5'-primer	3'-primer	rationale
7	Modified	-	Increased yield of HLA-specific PCR product.
21	Modified		Increased yield of HLA-specific PCR product.

101.512-24/03– including *Taq* pol., IFU-01
 101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: 92R**Lot-specific information**

33	Added	Added	Primer pair added for the B*07:149 allele
36	-	Added	3'-primer added for the B*07:167N allele.
41	Added	Added	Primer pair added for the B*07:149 allele, 5'-primer added for the B*07:156 allele.
50	Modified	Added	Improved specificity of primer pair, 3'-primer added for the B*07:175 allele.
54	-	Added	3'-primer added for the B*07:167N allele.
65	-	Modified	Improved specificity of primer pair.
68	-	Added	3'-primer added for the B*07:148 allele.
73	Added		Primer pair added for the B*07:119 allele, 5'-primer added for the B*07:169 allele.
74	Added		5'-primer added for the B*07:156 allele. .
75	Modified	Added	Improved specificity of primer pair, 3'-primer added for the B*07:175 allele.
77	Added	Added	Primer pair added for the B*07:148 allele.
79	Added	Added	Primer pair added for the B*07:119 allele.
80	New	New	New primer pair for the B*07:161N allele.

Change in revision R01 compared to R00:

1. Spelling mistake in the section of changes to previous lot, has been corrected.

101.512-24/03– including *Taq* pol., IFU-01
101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: **92R**

Lot-specific information

101.512-24/03– including *Taq* pol., IFU-01
 101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **92R**

Lot-specific information

PRODUCT DESCRIPTION

HLA-B*07 SSP subtyping

CONTENT

The primer set contains 5'- and 3'-primers for identifying the B*07:02 to B*07:176 alleles.

PLATE LAYOUT

Each test consists of 80 PCR reactions in a 96 well cut PCR plate. Wells 81 to 96 are empty.

1	2	3	4	5	6	7	8
9	10	11	12	13	14	15	16
17	18	19	20	21	22	23	24
25	26	27	28	29	30	31	32
33	34	35	36	37	38	39	40
41	42	43	44	45	46	47	48
49	50	51	52	53	54	55	56
57	58	59	60	61	62	63	64
65	66	67	68	69	70	71	72
73	74	75	76	77	78	79	80
empty	empty	empty	empty	empty	empty	empty	empty
empty	empty	empty	empty	empty	empty	empty	empty

The 96 well cut PCR plate is marked with 'HLA-B*07' in silver/gray ink.

Well No. 1 is marked with the Lot Number '92R'.

A faint row of numbers is seen between wells 1 and 2 or wells 7 and 8 of the PCR trays. These stem from the manufacture of the trays, and should be disregarded.

The PCR plates are covered with a PCR-compatible foil.

INTERPRETATION

The interpretation of HLA-B*07 SSP subtypings will be influenced by the B*08, the B*13, several B*14, several B*15, the B*18, the B*27, several B*35, the B*37, the B*38, the B*39, most B*40, the B*41, the B*42, several B*44, the B*45, four B*46, the B*48, the B*49, the B*50, several B*51, three B*52, seven B*53, the B*54, the B*55, the B*56, three B*57, four B*58, five B*59, four B*67, two B*73, five B*78, six B*81 and three B*82 alleles when present on the other haplotype.

In addition, the C*02:60 allele will be amplified by primer mix 31, the C*02:23 allele will be amplified by primer mixes 3, 10 and 18, the C*03:05, C*03:25, C*03:27 and C*03:167 alleles will be amplified by primer mixes 14, 26, 31 and 56, the C*03:124, C*04:49, C*05:15 and C*07:262 alleles will be amplified by primer mix 63, the C*03:137 and C*12:82 allele will be amplified by primer mixes 5 and 72, the C*03:143 will be amplified by primer mix 14, the C*04:77 allele will be amplified by primer mixes 4, 10, 18 and 53, the C*07:38:01-07:38:02 alleles will be amplified by

101.512-24/03– including *Taq* pol., IFU-01
101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: 92R**Lot-specific information**

primer mixes 10, 12, 13, 18, 29, 48 and 53, the C*07:102 allele will be weakly amplified by primer mixes 5 and 72, the C*08:08 allele will be amplified by primer mixes 14, 20, 26, 31 and 59 and the C*15:24 allele will be amplified by primer mix 32.

UNIQUELY IDENTIFIED ALLELES

All the HLA-B*07 alleles, i.e. **B*07:02 to B*07:176**, recognized by the HLA Nomenclature Committee in October 2012¹ will be amplified by the primers in the HLA-B*07 subtyping kit².

The HLA-B*07 kit enables separation of the confirmed HLA-B*07 alleles as listed in the IMGT/HLA database. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources. Current allele confirmation status for HLA-B*07 alleles is listed below.

The HLA-B*07 kit also enables identification of polymorphisms in exons outside of the region encoding the peptide binding domain and of null and alternatively expressed alleles.

The B*07:42 and 07:118 alleles may be distinguished by the different sizes of the specific PCR products generated by primer mix 26.

The B*07:49N and 07:169 alleles may be distinguished by the different sizes of the specific PCR products generated by primer mix 73.

The B*07:94 and 07:135 alleles may be distinguished by the different sizes of the specific PCR products generated by primer mix 63.

The B*07:102 and 07:106 alleles can be distinguished by the different sizes of the specific PCR products generated by primer mix 69.

The B*07:121 and 07:126 alleles may be distinguished by the different sizes of the specific PCR products generated by primer mix 75.

The HLA-B*07 subtyping kit cannot distinguish the following silent mutations: the B*07:02:01-07:02:25, 07:02:27 and 07:02:29-07:02:37, the B*07:05:01-07:05:06, the B*07:18:01-07:18:02, the B*07:22:01-07:22:02 the B*07:56:01-B*07:56:02, the B*07:68:01-07:68:03 or the B*07:96:01-07:96:02 alleles.

¹HLA-B alleles listed on the IMGT/HLA web page 2012-October-14, release 3.10.0, www.ebi.ac.uk/imgt/hla.

²The B*07:50 and B*48:08 alleles will give rise to identical amplification patterns. These alleles can e.g. be distinguished by the HLA-B low resolution primer set.

101.512-24/03– including *Taq* pol., IFU-01
101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: **92R**

Lot-specific information

ALLELE CONFIRMATION STATUS

Allele	Status ¹	Allele	Status ¹	Allele	Status ¹	Allele	Status ¹	Allele	Status ¹
B*07:02:01	Confirmed	B*07:11	Unconfirmed	B*07:58	Unconfirmed	B*07:105	Confirmed	B*07:155	Unconfirmed
B*07:02:02	Confirmed	B*07:12	Confirmed	B*07:59	Unconfirmed	B*07:106	Confirmed	B*07:156	Unconfirmed
B*07:02:03	Confirmed	B*07:13	Confirmed	B*07:60	Unconfirmed	B*07:107	Unconfirmed	B*07:157	Unconfirmed
B*07:02:04	Unconfirmed	B*07:14	Confirmed	B*07:61	Unconfirmed	B*07:108	Unconfirmed	B*07:158	Unconfirmed
B*07:02:05	Unconfirmed	B*07:15	Confirmed	B*07:62	Confirmed	B*07:109	Confirmed	B*07:159	Unconfirmed
B*07:02:06	Unconfirmed	B*07:16	Confirmed	B*07:63	Confirmed	B*07:110	Confirmed	B*07:160	Unconfirmed
B*07:02:07	Confirmed	B*07:17	Confirmed	B*07:64	Confirmed	B*07:111N	Confirmed	B*07:161N	Unconfirmed
B*07:02:08	Confirmed	B*07:18:01	Unconfirmed	B*07:65	Unconfirmed	B*07:112	Unconfirmed	B*07:162	Unconfirmed
B*07:02:09	Unconfirmed	B*07:18:02	Confirmed	B*07:66	Unconfirmed	B*07:113	Unconfirmed	B*07:163	Unconfirmed
B*07:02:10	Confirmed	B*07:19	Confirmed	B*07:67N	Unconfirmed	B*07:114	Unconfirmed	B*07:164	Unconfirmed
B*07:02:11	Unconfirmed	B*07:20	Confirmed	B*07:68:01	Confirmed	B*07:115	Unconfirmed	B*07:165	Unconfirmed
B*07:02:12	Unconfirmed	B*07:21	Unconfirmed	B*07:68:02	Confirmed	B*07:116	Unconfirmed	B*07:166	Unconfirmed
B*07:02:13	Confirmed	B*07:22:01	Unconfirmed	B*07:68:03	Confirmed	B*07:117	Unconfirmed	B*07:167N	Unconfirmed
B*07:02:14	Confirmed	B*07:22:02	Confirmed	B*07:69	Confirmed	B*07:118	Unconfirmed	B*07:168	Unconfirmed
B*07:02:15	Confirmed	B*07:23	Unconfirmed	B*07:70	Unconfirmed	B*07:119	Confirmed	B*07:169	Unconfirmed
B*07:02:16	Unconfirmed	B*07:24	Confirmed	B*07:71	Confirmed	B*07:120	Unconfirmed	B*07:170	Unconfirmed
B*07:02:17	Confirmed	B*07:25	Unconfirmed	B*07:72	Confirmed	B*07:121	Unconfirmed	B*07:171	Unconfirmed
B*07:02:18	Unconfirmed	B*07:26	Confirmed	B*07:73	Confirmed	B*07:122	Confirmed	B*07:172	Unconfirmed
B*07:02:19	Unconfirmed	B*07:27	Unconfirmed	B*07:74	Confirmed	B*07:123	Unconfirmed	B*07:173	Unconfirmed
B*07:02:20	Unconfirmed	B*07:28	Unconfirmed	B*07:75	Unconfirmed	B*07:124	Unconfirmed	B*07:174	Unconfirmed
B*07:02:21	Confirmed	B*07:29	Unconfirmed	B*07:76	Confirmed	B*07:125	Confirmed	B*07:175	Confirmed
B*07:02:22	Unconfirmed	B*07:30	Unconfirmed	B*07:77	Confirmed	B*07:126	Confirmed	B*07:176	Unconfirmed
B*07:02:23	Confirmed	B*07:31	Confirmed	B*07:78	Unconfirmed	B*07:127	Unconfirmed		
B*07:02:24	Unconfirmed	B*07:32	Confirmed	B*07:79	Unconfirmed	B*07:128	Unconfirmed		
B*07:02:25	Unconfirmed	B*07:33	Confirmed	B*07:80	Confirmed	B*07:129	Unconfirmed		
B*07:02:26	Unconfirmed	B*07:34	Unconfirmed	B*07:81	Unconfirmed	B*07:130	Unconfirmed		
B*07:02:27	Confirmed	B*07:35	Unconfirmed	B*07:82	Unconfirmed	B*07:131	Unconfirmed		
B*07:02:28	Confirmed	B*07:36	Confirmed	B*07:83	Unconfirmed	B*07:132	Unconfirmed		
B*07:02:29	Unconfirmed	B*07:37	Confirmed	B*07:84	Unconfirmed	B*07:133	Unconfirmed		
B*07:02:30	Unconfirmed	B*07:38	Confirmed	B*07:85	Unconfirmed	B*07:134	Unconfirmed		
B*07:02:31	Unconfirmed	B*07:39	Unconfirmed	B*07:86	Confirmed	B*07:135N	Unconfirmed		
B*07:02:32	Unconfirmed	B*07:40	Unconfirmed	B*07:87	Confirmed	B*07:136	Unconfirmed		
B*07:02:33	Unconfirmed	B*07:41	Unconfirmed	B*07:88	Confirmed	B*07:137	Unconfirmed		
B*07:02:34	Unconfirmed	B*07:42	Confirmed	B*07:89	Unconfirmed	B*07:138	Unconfirmed		
B*07:02:35	Unconfirmed	B*07:43	Confirmed	B*07:90	Unconfirmed	B*07:139	Unconfirmed		
B*07:02:36	Unconfirmed	B*07:44	Unconfirmed	B*07:91	Unconfirmed	B*07:140	Unconfirmed		
B*07:02:37	Unconfirmed	B*07:45	Unconfirmed	B*07:92	Confirmed	B*07:141	Unconfirmed		
B*07:03	Confirmed	B*07:46	Unconfirmed	B*07:93	Unconfirmed	B*07:142	Unconfirmed		
B*07:04	Confirmed	B*07:47	Confirmed	B*07:94	Confirmed	B*07:143	Unconfirmed		
B*07:05:01	Confirmed	B*07:48	Confirmed	B*07:95	Confirmed	B*07:144	Unconfirmed		
B*07:05:02	Unconfirmed	B*07:49N	Unconfirmed	B*07:96:01	Confirmed	B*07:145	Unconfirmed		
B*07:05:03	Unconfirmed	B*07:50	Unconfirmed	B*07:96:02	Unconfirmed	B*07:146	Unconfirmed		
B*07:05:04	Unconfirmed	B*07:51	Confirmed	B*07:97	Unconfirmed	B*07:147	Unconfirmed		
B*07:05:05	Confirmed	B*07:52	Unconfirmed	B*07:98	Confirmed	B*07:148	Confirmed		
B*07:05:06	Unconfirmed	B*07:53	Confirmed	B*07:99	Unconfirmed	B*07:149	Confirmed		
B*07:06	Confirmed	B*07:54	Confirmed	B*07:100	Unconfirmed	B*07:150	Unconfirmed		
B*07:07	Confirmed	B*07:55	Unconfirmed	B*07:101	Unconfirmed	B*07:151	Unconfirmed		
B*07:08	Confirmed	B*07:56:01	Unconfirmed	B*07:102	Confirmed	B*07:152	Unconfirmed		
B*07:09	Confirmed	B*07:56:02	Unconfirmed	B*07:103	Confirmed	B*07:153	Unconfirmed		
B*07:10	Confirmed	B*07:57	Confirmed	B*07:104	Unconfirmed	B*07:154	Unconfirmed		

¹ Allele status “confirmed” or “unconfirmed” as listed on the IMGT/HLA web page 2012-October-14, release 3.10.0, www.ebi.ac.uk/imgt/hla.

101.512-24/03– including *Taq* pol., IFU-01
 101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **92R**

Lot-specific information

SPECIFICITY TABLE

HLA-B*07 SSP subtyping

Specificities and sizes of the PCR products of the 80 primer mixes used for HLA-B*07 SSP subtyping

Primer Mix	Size of spec. PCR product ¹	Size of control band ²	Amplified HLA-B*07 alleles ^{3,4}	Other amplified HLA Class I alleles ^{3,5}
1	395 bp	800 bp	*07:02:01-07:02:27, 07:02:29-07:02:37, 07:04-07:07, 07:09, 07:10 ^w , 07:11-07:12, 07:14-07:15, 07:17- 07:26, 07:28-07:31, 07:33-07:36, 07:39- 07:49N, 07:51-07:64, 07:66-07:68:03, 07:70- 07:82, 07:84, 07:86- 07:124, 07:126-07:176	*14:21, 38:26, 40:166 [?] , 42:01:01- 42:02, 42:04-42:06, 42:08-42:10, 42:12-42:17, 45:06, 67:01:01- 67:01:02, 67:03, 73:01 ^w -73:02 ^w , 81:02
2	285 bp	1070 bp	*07:02:01-07:06, 07:08- 07:11, 07:13, 07:15- 07:17, 07:20-07:32, 07:34-07:39, 07:41- 07:52, 07:54-07:59, 07:61-07:76, 07:79- 07:83, 07:85-07:99, 07:101-07:121, 07:123- 07:136, 07:138, 07:140- 07:161N, 07:164, 07:166-07:176	*08:79, 35:66, 37:07, 40:15-40:16, 40:31-40:32, 40:45, 40:80, 40:98, 42:05:01-42:05:02, 48:01:01- 48:01:03, 48:05-48:12, 48:14-48:16, 48:18-48:20, 48:22, 48:27-48:28, 81:01-81:06
3	235 bp	1070 bp	*07:02:01-07:04, 07:07- 07:13, 07:15-07:16, 07:18:01-07:18:02, 07:20-07:27, 07:29- 07:31, 07:35-07:39, 07:41-07:52, 07:54- 07:59, 07:61-07:68:03, 07:70-07:73, 07:75- 07:77, 07:79, 07:81- 07:89, 07:91-07:95, 07:98-07:99, 07:101- 07:104, 07:106-07:111N, 07:113-07:114, 07:116- 07:121, 07:124-07:136, 07:141-07:153, 07:155-	*15:138, 15:241, 18:35, 35:66, 37:07, 40:31-40:33, 40:52, 40:59, 40:98, 46:06, 48:08, 48:23, C*02:23

101.512-24/03– including *Taq* pol., IFU-01
 101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **92R**

Lot-specific information

			07:164, 07:166-07:167N, 07:169-07:171, 07:173- 07:175	
4⁶	105 bp	800 bp	*07:02:01-07:02:25, 07:02:27-07:03, 07:05:01-07:18:02, 07:20-07:24, 07:27- 07:32, 07:35-07:39, 07:41-07:42, 07:44- 07:47, 07:49N-07:52, 07:54, 07:56:01-07:59, 07:61-07:65, 07:67N- 07:79, 07:81-07:99, 07:101-07:121, 07:123- 07:130, 07:132-07:138, 07:140-07:145, 07:148- 07:163, 07:166-07:176	*08:79, 35:66, 37:07, 40:15-40:16, 40:23, 40:32, 40:98, 40:128, 40:198, 44:150, 48:05, 48:08, C*04:77
5⁷	220 bp	1070 bp	*07:02:01-07:02:37, 07:04-07:07, 07:09- 07:15, 07:17-07:26, 07:28, 07:30-07:31, 07:33-07:36, 07:39- 07:49N, 07:51-07:64, 07:66-07:68:03, 07:73- 07:84, 07:87-07:176	*27:01-27:11, 27:13-27:15, 27:17, 27:19-27:21, 27:24-27:28, 27:30- 27:38, 27:40-27:58, 27:60-27:76, 27:78-27:91, 27:93-27:99, 38:26, 42:01:01-42:02, 42:04-42:06, 42:08- 42:10, 42:12-42:17, 54:02, 55:10, 56:35, 67:01:01-67:03, 73:01-73:02, 81:01-81:06, C*03:137, C*07:102^w, C*12:82
6	215 bp	1070 bp	*07:03, 07:08, 07:16, 07:27, 07:32, 07:38, 07:50, 07:69, 07:85	*08:94, 13:14, 14:01:01-14:01:02, 14:07N-14:08, 14:10, 14:12, 14:14, 14:19, 14:26, 14:32, 15:68, 15:71, 15:175, 18:01:01:01-18:01:12, 18:01:14-18:03, 18:05-18:15, 18:17N-18:24, 18:26-18:67, 18:69, 18:71-18:83, 27:12, 27:16, 27:18, 27:23, 27:29, 27:92, 35:50, 35:84, 35:162, 35:197, 37:01:01-37:37, 38:01:01-38:04, 38:06-38:25, 38:27- 38:32, 38:35-38:42, 39:01:01:01- 39:01:01:03, 39:01:03-39:03, 39:05:01-39:11, 39:13:01-39:20, 39:22-39:51, 39:53-39:62, 39:64- 39:71, 39:73-39:75, 39:77-39:79, 40:02:01-40:06:06, 40:08-40:09, 40:11:01-40:11:02, 40:13, 40:18- 40:20, 40:24, 40:26-40:27:01, 40:28- 40:29, 40:35, 40:37, 40:39-40:40, 40:44, 40:50, 40:56, 40:64, 40:68, 40:70-40:71, 40:74-40:75, 40:78,

101.512-24/03– including *Taq* pol., IFU-01
 101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **92R**

Lot-specific information

				40:82-40:83, 40:85-40:86, 40:89-40:91, 40:93-40:99, 40:104-40:105, 40:107, 40:109, 40:111, 40:115, 40:119-40:120, 40:122, 40:127, 40:131, 40:133Q, 40:142N-40:145, 40:148-40:149, 40:157, 40:159, 40:161-40:162, 40:165, 40:167, 40:169, 40:173-40:174, 40:176-40:177, 40:180-40:181, 40:184, 40:189-40:190, 40:201-40:203, 40:205-40:206, 40:209-40:211, 40:214, 40:220, 40:224-40:225, 42:07, 42:11, 44:103, 44:130-44:131, 48:01:01-48:21, 48:23-48:28, 51:57, 51:103, 51:116, 53:21
7^{6,12}	105 bp, 215 bp	1070 bp	*07:04, 07:25, 07:29, 07:52, 07:146	*40:136, 41:08, 42:05:01-42:05:02
8¹³	160 bp, 235 bp	1070 bp	*07:05:01-07:06, 07:32, 07:34, 07:59, 07:69, 07:78, 07:80, 07:90, 07:97, 07:105, 07:112, 07:123, 07:137-07:138, 07:140, 07:176	*08:79, 40:01:01-40:01:20, 40:01:22-40:01:25, 40:01:27, 40:01:29, 40:07, 40:10:01-40:10:02, 40:12, 40:14:01-40:16, 40:21-40:23, 40:25, 40:30, 40:34, 40:36, 40:38, 40:43, 40:45-40:49, 40:54-40:55, 40:61-40:63, 40:65-40:67, 40:69, 40:72:01-40:74, 40:76-40:77, 40:79-40:81, 40:84, 40:87:01-40:88, 40:100-40:102, 40:106, 40:108, 40:112-40:114, 40:116-40:118N, 40:121, 40:123, 40:125-40:126, 40:128, 40:132, 40:134-40:141, 40:146-40:147, 40:149-40:156, 40:158, 40:160:01-40:160:02, 40:166, 40:168, 40:170-40:171, 40:175, 40:179, 40:182-40:183, 40:185-40:186, 40:188, 40:191-40:199, 40:204, 40:207-40:208, 40:210, 40:212-40:213, 40:215-40:218, 40:221-40:223, 41:08, 42:05:01-42:05:02, 44:31, 48:01:01-48:01:03, 48:03:01-48:03:02, 48:05-48:07, 48:09-48:12, 48:14-48:16, 48:18-48:22, 48:28, 81:01-81:06
9⁸	355 bp	1070 bp	*07:05:01-07:05:06	*35:46, 37:07
10¹⁴	285 bp, 415 bp	800 bp	*07:07, 07:12, 07:14, 07:18:01-07:18:02, 07:47, 07:77, 07:81,	*08:52, 15:138, 15:230, 15:241, 18:35, 38:41, 40:01:01-40:01:17, 40:01:19-40:01:29, 40:07,

101.512-24/03– including *Taq* pol., IFU-01
 101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **92R**

Lot-specific information

			07:137, 07:162-07:163	40:10:01-40:10:02, 40:12, 40:14:01-40:14:03, 40:21-40:23, 40:25, 40:33, 40:36, 40:38, 40:42-40:43, 40:46-40:49, 40:52, 40:54-40:55, 40:59-40:63, 40:65-40:67, 40:69, 40:72:01-40:74, 40:76-40:77, 40:79, 40:81, 40:84, 40:87:01-40:88, 40:100-40:102, 40:106, 40:108, 40:109 [?] , 40:112-40:114, 40:116-40:118N, 40:121, 40:123-40:126, 40:128, 40:132, 40:134-40:141, 40:146-40:147, 40:149-40:156, 40:158, 40:160:01-40:160:02, 40:166, 40:168, 40:170-40:172, 40:175, 40:179, 40:182-40:183, 40:185-40:188, 40:191-40:199, 40:204, 40:207-40:208, 40:210, 40:212-40:213, 40:215-40:218, 40:221-40:223, 41:08, 44:31, 44:150, 46:06, 48:03:01-48:03:02, 48:17-48:18, 48:23, C*02:23, C*04:77, C*07:38:01-07:38:02
11¹⁵	130 bp, 210 bp	1070 bp	*07:08, 07:10, 07:32, 07:38-07:39	*08:94, 18:07:01-18:07:02, 18:29, 18:79, 27:01-27:05:08, 27:05:10-27:11, 27:13-27:15, 27:17, 27:19-27:21, 27:23-27:28, 27:30-27:76, 27:78-27:84, 27:86-27:91, 27:93-27:99, 35:50, 35:84, 35:162, 35:197, 37:21, 38:06-38:07, 39:58, 40:08, 40:13, 40:68, 40:162, 40:220, 48:06, 51:57, 51:103, 51:116, 53:21, 73:01-73:02
12	225 bp	1070 bp	*07:09, 07:11, 07:17, 07:162	*15:138, 15:241, 18:35, 35:66, 40:38, 40:52, 40:59-40:60, 40:158, 40:222, 46:06, 48:14, 48:23, C*07:38:01-07:38:02
13^{6,16}	85 bp, 290 bp	1070 bp	*07:12, 07:14, 07:18:01-07:18:02, 07:31, 07:34, 07:43, 07:137, 07:162-07:163	*40:30, 40:34, 40:59, 40:137, 40:160:01-40:160:02, 40:199, 41:08, 42:05:01-42:05:02, 44:150, 48:17, C*07:38:01-07:38:02
14⁹	215 bp	800 bp	*07:26, 07:34, 07:43, 07:48, 07:80, 07:131, 07:147, 07:164	*08:13, 08:25, 14:05, 14:13, 15:55, 18:14, 18:22, 27:07:01-27:07:03, 27:11, 27:24, 27:32-27:34, 27:43, 35:05:01-35:05:03, 35:16-35:17, 35:22, 35:30-35:31, 35:51, 35:58, 35:72, 35:89, 35:97, 35:113-35:114, 35:125, 35:199, 39:03, 39:14,

101.512-24/03– including *Taq* pol., IFU-01
 101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **92R**

Lot-specific information

				39:24:01-39:24:02, 39:29, 40:02:01-40:03, 40:05, 40:08-40:09, 40:13, 40:18-40:19, 40:24, 40:27:01-40:27:02, 40:29, 40:31, 40:35, 40:37, 40:39-40:40, 40:45, 40:50, 40:56-40:58, 40:71, 40:78, 40:80, 40:82, 40:85, 40:89-40:91, 40:94, 40:97, 40:104-40:105, 40:107, 40:111, 40:115, 40:119, 40:122, 40:133Q, 40:142N-40:145, 40:157, 40:164, 40:169, 40:173-40:174, 40:176, 40:180-40:181, 40:189, 40:200-40:203, 40:205-40:206, 40:209, 40:211, 40:214, 40:219-40:220, 40:224, 40:226, 44:54, 44:110, 44:135, 48:01:01-48:01:03, 48:04, 48:06-48:07, 48:09-48:16, 48:18, 48:20, 48:22, 48:24, 48:27-48:28, 51:64, 51:81, 51:148, 53:14, 55:04, 55:49, 55:51, 58:18, 58:27, 81:01-81:06, C*03:05, C*03:25, C*03:27, C*03:143, C*03:167, C*08:08
15¹⁷	385 bp, 415 bp	1070 bp	*07:11, 07:47, 07:54, 07:57, 07:75, 07:123	*08:17, 08:47, 39:20, 42:15, 45:06
16	230 bp	1070 bp	*07:16, 07:27, 07:37- 07:38, 07:50, 07:75, 07:91	*13:14, 14:14, 15:68, 15:71, 15:175, 18:01:01:01-18:01:01:02, 18:01:03-18:01:12, 18:01:14-18:02, 18:05-18:15, 18:17N-18:24, 18:26-18:67, 18:69, 18:71-18:83, 27:01, 27:18, 27:23, 27:29, 27:44, 27:92, 35:50, 35:84, 35:162, 35:197, 37:01:01-37:36, 38:01:01-38:02:04, 38:04, 38:06-38:12, 38:15-38:25, 38:27-38:32, 38:34N-38:40, 38:42, 39:05:01-39:05:02, 39:07-39:08, 39:11, 39:13:01-39:13:02, 39:20, 39:37, 39:49, 39:55-39:56, 40:02:01-40:06:06, 40:08-40:09, 40:11:01-40:11:02, 40:13, 40:18-40:20, 40:24, 40:26-40:29, 40:35, 40:37, 40:39-40:40, 40:44, 40:56, 40:64, 40:68, 40:70-40:71, 40:74-40:75, 40:78, 40:82-40:83, 40:85-40:86, 40:89-40:91, 40:93-40:99, 40:103-40:105, 40:107, 40:109, 40:111, 40:115, 40:119-40:120, 40:122, 40:127, 40:131, 40:133Q, 40:142N-40:145,

101.512-24/03– including *Taq* pol., IFU-01
 101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **92R**

Lot-specific information

				40:148-40:149, 40:157, 40:159, 40:161-40:162, 40:165, 40:167, 40:169, 40:173-40:174, 40:176-40:177, 40:180-40:181, 40:184, 40:189-40:190, 40:200-40:203, 40:205-40:206, 40:209-40:211, 40:214, 40:220, 40:224-40:226, 42:11, 44:103, 44:130-44:131, 48:01:01-48:21, 48:23-48:28, 51:57, 51:103, 51:116, 53:21, 57:10
17¹⁸	345 bp, 425 bp	1070 bp	*07:11, 07:27, 07:46, 07:57, 07:75	*08:17, 08:47, 38:17, 39:20, 40:188 [?] , 42:15
18	265 bp	1070 bp	*07:12, 07:78, 07:84, 07:162	*35:66, 40:59-40:60, 40:125, 48:17, C*02:23, C*04:77, C*07:38:01- 07:38:02
19¹⁹	205 bp, 250 bp, 295 bp	800 bp	*07:20, 07:35, 07:140, 07:153	*15:14, 15:91, 15:131, 15:161, 18:56, 35:45, 35:71, 35:169, 40:132, 40:178, 44:10, 44:14, 44:17, 44:31, 44:43:01-44:43:02, 44:45, 44:51, 44:54, 44:62, 44:77, 44:82, 44:106- 44:107, 44:123, 44:135, 44:158, 45:09, 45:14, 46:17, 51:23, 51:108, 52:15, 53:22, 57:07, 58:07, 82:01- 82:03
20²⁰	160 bp, 205 bp, 275 bp	1070 bp	*07:21, 07:51, 07:59	C*08:08
21^{6,21}	100 bp, 165 bp	1070 bp	*07:22:01-07:22:02, 07:52	
22	385 bp	1070 bp	*07:23	*40:79, 40:190, 41:14
23⁶	85 bp	1070 bp	*07:20, 07:24, 07:131	*15:230, 40:63
24^{6,7}	125 bp	800 bp	*07:25, 07:39, 07:48	*18:35, 27:36, 35:66, 40:01:01- 40:01:19, 40:01:21-40:01:29, 40:07, 40:10:01-40:10:02, 40:12, 40:14:01- 40:14:03, 40:21-40:22N, 40:25, 40:30-40:31, 40:33-40:34, 40:36, 40:38, 40:42-40:43, 40:45-40:49, 40:52, 40:54-40:55, 40:59-40:63, 40:65-40:66, 40:69, 40:72:01-40:74, 40:76-40:77, 40:79-40:81, 40:84, 40:87:01-40:88, 40:100-40:102, 40:106, 40:108, 40:112-40:114, 40:116-40:118N, 40:121, 40:123- 40:126, 40:128, 40:132, 40:134- 40:141, 40:146-40:147,

101.512-24/03– including *Taq* pol., IFU-01
 101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **92R**

Lot-specific information

				40:149-40:156, 40:160:01-40:160:02, 40:166, 40:168, 40:170-40:172, 40:175, 40:179, 40:182, 40:185-40:188, 40:191-40:197, 40:199, 40:204, 40:207-40:208, 40:210, 40:212-40:213, 40:215-40:218, 40:221-40:223, 41:08, 42:05:01-42:05:02, 44:31, 46:06, 48:01:01-48:01:03, 48:03:01-48:03:02, 48:06-48:07, 48:09-48:12, 48:14, 48:16-48:23, 48:26-48:28, 81:01-81:06
25 ^{10,22}	305 bp, 425 bp	1070 bp	*07:27, 07:117	*38:17, 40:188 [?]
26 ^{6,9,11,23}	95 bp, 205 bp, 235 bp, 290 bp	800 bp	*07:20, 07:24, 07:28, 07:42, 07:51, 07:60, 07:100, 07:118, 07:131	*08:14, 08:21, 08:25, 08:87, 15:07:01-15:07:02, 15:45, 15:55, 15:68, 15:126, 15:207, 27:32, 35:05:01-35:05:03, 35:22, 35:51, 35:58, 35:72, 35:89, 35:97, 35:114, 35:199, 40:05, 40:71, 40:174, 40:180, 44:54, 44:106, 44:135, 44:158, 46:12, 51:64, 51:148, 53:14, 55:51, 58:18, C*03:05, C*03:25, C*03:27, C*03:167, C*08:08
27 ²⁴	165 bp, 215 bp	800 bp	*07:29, 07:41	
28 ^{6,9}	110 bp	1070 bp	*07:30	
29 ^{6,7,9,25}	85 bp, 200 bp	800 bp	*07:31, 07:34, 07:43- 07:44	*40:199, 41:08, 42:05:01-42:05:02, C*07:38:01-07:38:02
30	190 bp	800 bp	*07:05:01-07:06, 07:32, 07:40, 07:53, 07:69, 07:78, 07:80, 07:90, 07:97, 07:105, 07:112, 07:123, 07:137-07:138, 07:140, 07:176	*08:40, 08:79, 13:01:01-13:01:02, 13:01:04-13:02:15, 13:07N-13:09, 13:11, 13:13-13:23, 13:25, 13:27- 13:33, 13:36-13:45, 13:47, 13:49N- 13:50, 13:52, 13:54-13:61, 13:63N, 27:07:01-27:07:03, 27:11, 27:20, 27:24, 27:27, 27:33, 27:43, 35:31, 35:88, 39:47, 39:50, 40:01:01- 40:01:20, 40:01:22-40:01:25, 40:01:27, 40:01:29-40:02:09, 40:02:11-40:02:12, 40:02:14- 40:02:17, 40:04, 40:06:01:01-40:08, 40:10:01-40:16, 40:19, 40:21-40:23, 40:25, 40:27:01-40:27:02, 40:29- 40:30, 40:34-40:38, 40:40, 40:43, 40:45-40:50, 40:53-40:58, 40:61- 40:62, 40:64-40:70, 40:72:01-40:88,

101.512-24/03– including *Taq* pol., IFU-01
 101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **92R**

Lot-specific information

				40:90-40:91, 40:93-40:97, 40:99-40:119, 40:121-40:123, 40:125-40:130, 40:132, 40:134-40:158, 40:160:01-40:162, 40:164-40:171, 40:173, 40:175-40:179, 40:182-40:183, 40:185-40:186, 40:188-40:198, 40:200-40:226, 41:19, 44:31, 48:01:01-48:01:03, 48:03:01-48:07, 48:09-48:16, 48:18-48:22, 48:24, 48:28, 51:10, 51:16, 51:31, 51:34, 51:81-51:82, 51:93, 52:08, 54:26, 55:09, 55:22, 55:24, 58:27, 81:01-81:06
31²⁶	215 bp, 280 bp	800 bp	*07:19, 07:33, 07:53, 07:55, 07:60, 07:100	*08:70, 08:103, 14:05, 14:13, 15:07:01-15:07:02, 15:45, 15:55, 15:68, 15:126, 15:207, 18:14, 27:07:01-27:07:03, 27:11, 27:24, 27:32-27:34, 27:43, 35:05:01-35:05:03, 35:22, 35:31, 35:51, 35:58, 35:72, 35:89, 35:97, 35:114, 35:199, 37:09, 38:19, 39:03, 39:14, 39:24:01-39:24:02, 39:29, 39:37, 39:76, 40:02:01-40:03, 40:05, 40:08-40:09, 40:13, 40:18-40:19, 40:24, 40:27:01-40:27:02, 40:29, 40:35, 40:37, 40:39-40:40, 40:50, 40:56-40:58, 40:71, 40:78, 40:82, 40:85, 40:89-40:91, 40:94, 40:97, 40:104-40:105, 40:107, 40:111, 40:115, 40:119, 40:122, 40:133Q, 40:142N-40:145, 40:157, 40:164, 40:169, 40:173-40:174, 40:176, 40:180-40:181, 40:189, 40:200-40:203, 40:205-40:206, 40:209, 40:211, 40:214, 40:219-40:220, 40:224, 40:226, 42:08, 44:54, 44:106, 44:135, 44:158, 46:12, 48:04, 48:13, 48:19, 48:24, 51:64, 51:81, 51:148, 53:14, 55:04, 55:49, 55:51, 58:18, 58:27, C*02:60, C*03:05, C*03:25, C*03:27, C*03:167, C*08:08
32¹¹	190 bp	1070 bp	*07:34	*08:01:01-08:05, 08:08N-08:10, 08:12:01-08:13, 08:15-08:20, 08:22-08:24, 08:26-08:27, 08:29-08:36, 08:38-08:39, 08:41-08:48, 08:50-08:73, 08:75-08:76, 08:78, 08:80-08:86N, 08:88, 08:90-08:105, 13:06,

101.512-24/03– including *Taq* pol., IFU-01
 101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **92R**

Lot-specific information

				13:53, 14:13, 15:42, 15:83, 15:86, 15:188, 15:224, 15:252, 18:26, 18:61, 35:60, 35:87, 37:12, 38:01:01, 38:01:03-38:01:05, 38:01:07-38:07, 38:09-38:21, 38:23-38:28, 38:31-38:42, 39:01:01:01-39:01:01:03, 39:01:03-39:06:02, 39:08-39:14, 39:16, 39:18-39:20, 39:22-39:24:02, 39:26-39:28, 39:30-39:31, 39:34-39:37, 39:39:01-39:39:02, 39:41-39:42, 39:44-39:46, 39:48, 39:51-39:54, 39:56-39:66, 39:68-39:73, 39:75-39:79, 40:39, 40:51, 40:163, 40:184, 40:199, 41:01-41:03:02, 41:05-41:09, 41:11-41:16, 41:20, 41:22-41:25, 42:01:01-42:02, 42:04-42:08, 42:10-42:15, 42:17, 46:18, 51:21, 51:36, 51:101, 51:143, 54:01:01-54:02, 54:04-54:05N, 54:07-54:08N, 54:10-54:13, 54:17-54:25, 54:27-54:28, 55:01:01-55:05, 55:07, 55:10-55:13, 55:15-55:21, 55:25-55:26, 55:28-55:40, 55:42-55:45, 55:47-55:50, 55:52-55:57, 56:10, 56:12, 56:23, 58:20, 59:01:01:01-59:03, 59:05, 67:01:01-67:03, 78:07, C*15:24
33 ^{6,7,27}	90 bp, 350 bp	1070 bp	*07:61, 07:149	*08:02, 13:36, 15:89, 15:115, 18:09, 38:02:01-38:02:04, 38:08, 38:15, 38:18, 38:23, 38:29, 38:35, 44:12, 51:54, 51:78:01-51:78:02, 53:09, 53:11-53:13, 56:07
34	155 bp	1070 bp	*07:62, 07:90	
35 ⁶	105 bp	1070 bp	*07:63	
36 ²⁸	180 bp, 250 bp	1070 bp	*07:64, 07:167N	
37	210 bp	1070 bp	*07:65	*18:04, 35:08:05, 35:42:02, 35:185, 39:63, 40:01:11, 40:58, 56:05:01, 56:06, 78:01:01-78:01:02, 78:02:02-78:03, 78:07
38 ^{9,29}	190 bp, 215 bp	1070 bp	*07:66, 07:103, 07:113	*08:64
39 ⁸	425 bp	1070 bp	*07:05:01-07:06, 07:34, 07:40, 07:53, 07:80, 07:90, 07:97, 07:105,	*40:166, 42:01:01-42:02, 42:04-42:05:02, 42:07-42:08, 42:10, 42:12-42:15, 42:17, 54:04, 54:11, 55:04,

101.512-24/03– including *Taq* pol., IFU-01
 101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **92R**

Lot-specific information

			07:112, 07:137-07:138, 07:140, 07:176	55:08, 55:17, 55:20, 55:28, 55:49, 55:51, 56:05:01-56:06, 56:12, 56:15, 56:21, 81:01-81:04N, 81:06
40^{9,30}	265 bp, 395 bp	1070 bp	*07:56:01-07:56:02, 07:72, 07:83	*08:26, 08:56:01-08:56:02, 38:28, 39:63, 40:21 [?] , 40:58 [?] , 40:137 [?] , 41:12, 49:03
41^{6,7,31}	90 bp, 415 bp	1070 bp	*07:58, 07:149, 07:156	*08:02, 13:36, 15:89, 15:115, 18:09, 38:02:01-38:02:04, 38:08, 38:15, 38:18, 38:23, 38:29, 38:35, 44:12, 51:54, 51:78:01-51:78:02, 53:09, 53:11-53:13, 56:07
42³²	195 bp, 225 bp	1070 bp	*07:67N, 07:71	
43⁷	140 bp	1070 bp	*07:68:01-07:68:03, 07:100	*08:71, 14:05, 14:13, 15:07:01- 15:07:02, 15:45, 15:55, 15:68, 15:126, 15:207, 18:14, 27:32, 35:05:01-35:05:03, 35:22, 35:51, 35:58, 35:66, 35:72, 35:89, 35:97, 35:199, 38:19, 39:03, 39:14, 39:24:01-39:24:02, 39:29, 39:37, 39:76, 40:35, 40:71, 41:24, 42:10, 42:17, 44:54, 44:106, 44:135, 44:158, 46:12, 48:12, 48:14, 51:64, 53:14, 58:18
44	150 bp	1070 bp	*07:10	*27:75 ^w , 73:01 ^w -73:02 ^w
45⁶	105 bp	1070 bp	*07:27, 07:50	*08:04, 08:17, 08:54, 15:03:01- 15:03:03, 15:47:01-15:47:02, 15:54, 15:61-15:62, 15:64:01-15:64:02, 15:68-15:69, 15:91, 15:98, 15:103, 15:123, 15:127, 15:131-15:132, 15:151, 15:156, 15:158, 15:173, 15:210, 15:220, 15:235, 15:242, 15:253, 15:259, 18:37, 27:18, 27:29, 37:28, 38:03, 39:02:01-39:02:02, 39:08, 39:13:01-39:13:02, 39:23, 39:39:01-39:39:02, 39:49, 40:12, 40:46, 40:93, 42:11, 44:40, 44:44, 44:130, 44:156, 48:01:01-48:05, 48:07-48:28, 52:16, 55:18
46	425 bp	1070 bp	*07:13, 07:15, 07:92	*08:15, 39:27, 40:73, 67:02, 73:01 ^w - 73:02 ^w
47	215 bp	800 bp	*07:04, 07:19, 07:25, 07:146	*08:01:01-08:05, 08:07-08:08N, 08:10-08:11, 08:14-08:15, 08:17- 08:19N, 08:21-08:24, 08:26-08:39, 08:41-08:48, 08:50-08:54,

101.512-24/03– including *Taq* pol., IFU-01
 101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **92R**

Lot-specific information

				08:56:01-08:59, 08:61-08:69, 08:71-08:78, 08:80-08:83, 08:85-08:88, 08:90-08:105, 35:87, 37:09, 41:02:01-41:02:04, 41:04, 41:10-41:11, 41:13, 41:15, 41:18-41:19, 41:23-41:24, 42:01:01-42:02, 42:05:01-42:12, 42:14-42:17, 44:106, 44:158
48 ³³	200 bp, 280 bp	1070 bp	*07:76-07:77	*40:60, C*07:38:01-07:38:02
49	200 bp	1070 bp	*07:82	*44:129
50 ^{6,34}	75 bp, 105 bp, 165 bp	1070 bp	*07:74, 07:87, 07:175	*08:77, 15:126, 48:27
51 ³⁵	345 bp, 470 bp	1070 bp	*07:93, 07:97	*08:39, 38:27
52	290 bp	800 bp	*07:95	
53 ⁹	235 bp	800 bp	*07:14, 07:17, 07:28, 07:96:01-07:96:02, 07:115	*15:230, 27:36, 40:42, 40:60, 40:124:01-40:124:02, 40:172, 44:150, 48:17, 48:26, C*04:77, C*07:38:01-07:38:02
54 ⁷	190 bp	1070 bp	*07:98, 07:167N	
55	260 bp	1070 bp	*07:99	
56 ^{7,9,11}	230 bp	1070 bp	*07:20, 07:24, 07:45, 07:60, 07:100, 07:131	*08:21, 08:25, 15:07:01-15:07:02, 15:45, 15:55, 15:68, 15:126, 15:207, 35:05:01-35:05:03, 35:22, 35:51, 35:58, 35:72, 35:89, 35:97, 35:114, 35:199, 40:05, 40:71, 40:174, 44:54, 44:106, 44:135, 44:158, 46:12, 51:64, 51:148, 53:14, 55:51, 58:18, C*03:05, C*03:25, C*03:27, C*03:167
57 ³⁶	150 bp, 215 bp	1070 bp	*07:70, 07:109	
58 ^{7,37}	245 bp, 440 bp	1070 bp	*07:79, 07:112	
59 ^{6,38}	95 bp, 190 bp	800 bp	*07:89, 07:113, 07:115	*08:64, 08:69, 14:05, 38:19, 39:03, 39:24:01-39:24:02, 39:37, 39:76, 41:18, 42:06, 48:07, C*08:08
60	145 bp	1070 bp	*07:73	*40:166, 44:90, 45:06
61 ³⁹	205 bp, 420 bp	1070 bp	*07:88, 07:105	
62	415 bp	800 bp	*07:16, 07:27, 07:37-	*08:03, 08:23, 08:54, 14:14, 14:18,

101.512-24/03– including *Taq* pol., IFU-01
 101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **92R**

Lot-specific information

			07:38, 07:50, 07:75, 07:91	15:68, 38:01:01-38:02:04, 38:04-38:12, 38:15-38:25, 38:27- 38:40, 38:42, 39:05:01-39:05:02, 39:07-39:08, 39:11, 39:13:01- 39:13:02, 39:20, 39:37, 39:49, 39:55- 39:56, 40:01:01-40:01:29, 40:07, 40:09 [?] , 40:11:01 [?] -40:11:02 [?] , 40:14:01 [?] -40:15 [?] , 40:16, 40:18 [?] - 40:19 [?] , 40:21 [?] -40:22N [?] , 40:23, 40:24 [?] -40:38 [?] , 40:42 [?] , 40:44 [?] -40:49 [?] , 40:51, 40:52 [?] -40:58 [?] , 40:59-40:60, 40:61 [?] -40:65 [?] , 40:67 [?] -40:69 [?] , 40:71 [?] , 40:72:01-40:73, 40:74 [?] -40:76 [?] , 40:77, 40:78 [?] , 40:79, 40:80 [?] -40:88 [?] , 40:90 [?] -40:93 [?] , 40:95 [?] -40:96 [?] , 40:98 [?] - 40:121 [?] , 40:123 [?] , 40:124:01, 40:124:02 [?] -40:137 [?] , 40:138-40:140, 40:141 [?] , 40:145 [?] , 40:146, 40:148 [?] - 40:149 [?] , 40:150, 40:151 [?] , 40:152- 40:153, 40:154 [?] , 40:155N, 40:156 [?] - 40:165 [?] , 40:167 [?] -40:171 [?] , 40:172, 40:173 [?] -40:185 [?] , 40:186, 40:187 [?] - 40:213 [?] , 40:215 [?] -40:226 [?] , 41:01- 41:25, 42:11, 44:18, 45:01-45:05, 45:07-45:14, 48:01:01-48:28, 49:01:01-49:25, 50:01:01, 50:01:02 [?] , 50:01:03-50:02, 50:04-50:19
63 ^{6,40}	65 bp, 230 bp	1070 bp	*07:94, 07:135N	C*03:124, C*04:49, C*05:15, C*07:262
64	210 bp	1070 bp	*07:69, 07:85	*14:01:01-14:01:02, 14:07N-14:08, 14:10, 14:12, 14:14, 14:19, 14:26, 14:32, 18:33, 27:12, 27:16, 27:92, 38:01:01-38:02:04, 38:04, 38:08- 38:25, 38:27-38:32, 38:34N-38:42, 39:01:01:01-39:01:01:03, 39:01:03- 39:01:15, 39:03, 39:05:01-39:07, 39:09, 39:11, 39:14-39:15, 39:18- 39:19:02, 39:22, 39:24:01-39:38Q, 39:40:01N-39:48, 39:50-39:51, 39:53-39:57, 39:59-39:62, 39:64- 39:71, 39:73-39:75, 39:77-39:79, 40:75
65 ⁴¹	385 bp, 450 bp	1070 bp	*07:54, 07:111N, 07:123	*45:06
66	435 bp	1070 bp	*07:36, 07:38, 07:81	*08:03, 08:52, 08:78, 38:01:01- 38:01:07, 38:05-38:07, 38:09-38:14,

101.512-24/03– including *Taq* pol., IFU-01
 101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **92R**

Lot-specific information

				38:16, 38:19-38:22, 38:24-38:28, 38:30-38:34N, 38:36-38:42, 40:19 [?] , 40:109 [?] , 40:117 [?] , 44:18, 48:18, 49:01:01-49:01:04, 49:03-49:25
67⁸	230 bp	800 bp	*07:29, 07:72, 07:86, 07:114	*08:01:01-08:05, 08:07-08:56:02, 08:58-08:93, 08:95-08:105, 14:11, 14:28, 38:33, 39:04, 39:63, 40:01:01- 40:01:29, 40:07, 40:14:01 [?] -40:15 [?] , 40:16, 40:21 [?] -40:22N [?] , 40:23, 40:25 [?] , 40:30 [?] -40:34 [?] , 40:36 [?] , 40:38 [?] , 40:42 [?] -40:43 [?] , 40:45 [?] -40:49 [?] , 40:51, 40:52 [?] -40:55 [?] , 40:58 [?] , 40:59- 40:60, 40:61 [?] -40:63 [?] , 40:65 [?] , 40:66, 40:67 [?] , 40:69 [?] , 40:72:01-40:73, 40:76 [?] , 40:77, 40:79, 40:80 [?] -40:81 [?] , 40:84 [?] , 40:87:01 [?] -40:88 [?] , 40:92 [?] , 40:100 [?] -40:102 [?] , 40:106 [?] , 40:108 [?] , 40:110 [?] , 40:112 [?] -40:114 [?] , 40:116 [?] - 40:118N [?] , 40:121 [?] , 40:123 [?] , 40:124:01, 40:124:02 [?] -40:126 [?] , 40:128 [?] -40:130 [?] , 40:132 [?] , 40:134 [?] - 40:137 [?] , 40:138-40:140, 40:141 [?] , 40:146-40:147, 40:150, 40:151 [?] , 40:152-40:153, 40:154 [?] , 40:155N, 40:156 [?] , 40:158 [?] , 40:160:01 [?] - 40:160:02 [?] , 40:163 [?] , 40:166 [?] , 40:168 [?] , 40:170 [?] -40:171 [?] , 40:172, 40:175 [?] , 40:178 [?] -40:179 [?] , 40:182 [?] - 40:183 [?] , 40:185 [?] , 40:186, 40:187 [?] - 40:188 [?] , 40:191 [?] -40:194 [?] , 40:196 [?] - 40:199 [?] , 40:204 [?] , 40:207 [?] -40:208 [?] , 40:212 [?] -40:213 [?] , 40:215 [?] -40:219 [?] , 40:221 [?] -40:223 [?] , 41:01-41:25, 44:18, 45:01-45:14, 48:22, 49:01:01-49:23, 49:25, 50:01:01, 50:01:02 [?] , 50:01:03-50:02, 50:04-50:19
68⁴²	335 bp, 365 bp, 415 bp	1070 bp	*07:13, 07:15, 07:101, 07:148, 07:160	*08:16, 40:193 [?] , 67:02
69^{7,43}	190 bp, 245 bp, 465 bp	1070 bp	*07:102, 07:106, 07:108	*08:18, 73:01-73:02
70^{6,44}	125 bp, 205 bp	800 bp	*07:104, 07:114	
71⁴⁵	245 bp, 280 bp	1070 bp	*07:107-07:108	

101.512-24/03– including *Taq* pol., IFU-01
 101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **92R**

Lot-specific information

72	165 bp	1070 bp	*07:13, 07:110	*38:11, 54:02, 67:02, 81:06, C*03:137, C*07:102^w, C*12:82
73^{6,46}	85 bp, 240 bp	1070 bp	*07:49N, 07:119, 07:169	
74⁴⁷	345 bp, 425 bp	1070 bp	*07:120, 07:156	*40:150, 57:29
75^{6,46}	105 bp, 250 bp	1070 bp	*07:121, 07:126, 07:175	*15:126
76⁴⁷	170 bp, 285 bp	1070 bp	*07:122, 07:128	
77⁴⁸	335 bp, 405 bp	800 bp	*07:116, 07:148	*40:193 [?]
78^{6,49}	110 bp, 380 bp	1070 bp	*07:125, 07:129	
79^{6,50}	65 bp, 95 bp	1070 bp	*07:119, 07:130	
80	405 bp	1070 bp	*07:161N	

¹Alleles are assigned by the presence of specific PCR product(s). However, the sizes of the specific PCR products may be helpful in the interpretation of HLA-B*07 SSP typings.

When the primers in a primer mix can give rise to HLA-specific PCR products of more than one length this is indicated if the size difference is more than 20 base pairs. Size differences of 20 base pairs or less are not given. For high resolution SSP kits the respective lengths of the HLA-specific PCR product(s) are given for the alleles amplified by these primer mixes.

Nonspecific amplifications, i.e. a ladder or a smear of bands, may sometimes be seen. GC-rich primers have a higher tendency of giving rise to nonspecific amplifications than other primers.

PCR fragments longer than the control bands may sometimes be observed. Such bands should be disregarded and do not influence the interpretation of the SSP typings.

PCR fragments migrating faster than the control bands, but slower than a 400 bp fragment may be seen in some gel read-outs. Such bands can be disregarded and do not influence the interpretation of the SSP typings.

Some primers may give rise to primer oligomer artifacts. Sometimes this phenomenon is an inherent feature of the primer pair(s) of a primer mix. More often it is due to other factors such as too low amount of DNA in the PCR reactions, taking too long time in setting up the PCR reactions, working at elevated room temperature or using thermal cyclers that are not pre-heated.

²The internal positive control primer pairs amplify segments of the human growth hormone gene. The two different control primer pairs give rise to either an internal positive control band of 1070 base pairs, for most wells, or a band of 800 base pairs, for some wells.

Well number 1 contains the primer pair giving rise to the shorter, 800 bp, internal positive control band in order to help in the correct orientation of the HLA-B*07 subtyping.

In addition, wells number 4, 10, 14, 19, 24, 26, 27, 29 to 31, 47, 52, 53, 59, 62, 67, 70 and 77 contain the primer pair giving rise to the shorter, 800 bp, internal positive control band in order to allow kit identification. In the presence of a specific amplification the intensity of the control band often decreases.

³For several HLA B alleles 1st and/or 4th exon(s) and beyond, as well as intron nucleotide sequences, are not available. In these instances it is not known whether some of the primers of the SSP sets are completely matched with the target sequences or not. We assume that unknown sequences in these regions are conserved within allelic groups.

⁴The B*07:42 and 07:118 alleles may be distinguished by the different sizes of the specific PCR products generated by primer mix 26.

The B*07:49N and 07:169 alleles may be distinguished by the different sizes of the specific PCR products generated by primer mix 73.

101.512-24/03– including *Taq* pol., IFU-01
 101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: 92R**Lot-specific information**

The B*07:94 and 07:135 alleles may be distinguished by the different sizes of the specific PCR products generated by primer mix 63.

The B*07:102 and 07:106 alleles can be distinguished by the different sizes of the specific PCR products generated by primer mix 69.

The B*07:121 and 07:126 alleles may be distinguished by the different sizes of the specific PCR products generated by primer mix 75.

The B*07:50 and B*48:08 alleles will give rise to identical amplification patterns. These alleles can e.g. be distinguished by the HLA-B low resolution primer set.

⁵Due to the sharing of sequence motifs between HLA-B alleles non-HLA-B*07 alleles will be amplified by primer mixes 1 to 19, 22 to 26, 29 to 33, 37 to 41, 43 to 51, 53, 56, 59, 60, 62, 64 to 69, 72, 74, 75 and 77. In addition, the C*02:60 allele will be amplified by primer mix 31, the C*02:23 allele will be amplified by primer mixes 3, 10 and 18, the C*03:05, C*03:25, C*03:27 and C*03:167 alleles will be amplified by primer mixes 14, 26, 31 and 56, the C*03:124, C*04:49, C*05:15 and C*07:262 alleles will be amplified by primer mix 63, the C*03:137 and C*12:82 allele will be amplified by primer mixes 5 and 72, the C*03:143 will be amplified by primer mix 14, the C*04:77 allele will be amplified by primer mixes 4, 10, 18 and 53, the C*07:38:01-07:38:02 alleles will be amplified by primer mixes 10, 12, 13, 18, 29, 48 and 53, the C*07:102 allele will be weakly amplified by primer mixes 5 and 72, the C*08:08 allele will be amplified by primer mixes 14, 20, 26, 31 and 59 and the C*15:24 allele will be amplified by primer mix 32.

⁶HLA-specific PCR products shorter than 125 base pairs have a lower intensity and are less sharp than longer PCR products.

⁷Primer mixes 24, 29, 33, 41, 43, 54, 56, 58 and 69 have a tendency to giving rise to primer oligomer formation.

⁸Primer mixes 9, 39 and 67 may give rise to a lower yield of HLA-specific PCR product than the other B*07 primer mixes

⁹Primer mixes 14, 28, 29, 38, 40, 53 and 56 may have tendencies of unspecific amplifications. most pronounced in primer mixes 29 and 53.

¹⁰Primer mix 25 may give rise to unspecific amplification in B*73 alleles.

¹¹Primer mix 32 may give rise to a faint band in many B*14, B*37, C*04, C*05, C*08 and C*18 alleles.

¹²Primer mix 7: Specific PCR fragment of 105 bp in the B*07:04, 07:25, 07:52 and 07:146 and the B*40:136, 41:08 and 42:05:01-42:05:02 alleles. Specific PCR fragment of 215 bp in the B*07:29 allele.

¹³Primer mix 8: Specific PCR fragment of 160 bp in the B*07:59 allele. Specific PCR fragment of 235 bp in the B*07:05:01-07:06, 07:32, 07:34, 07:69, 07:78, 07:80, 07:90, 07:97, 07:105, 07:112, 07:123, 07:137-07:138, 07:140, 07:176 and the B*08:79, 40:01:01-40:01:20, 40:01:22-40:01:25, 40:01:27, 40:01:29, 40:07, 40:10:01-40:10:02, 40:12, 40:14:01-40:16, 40:21-40:23, 40:25, 40:30, 40:34, 40:36, 40:38, 40:43, 40:45-40:49, 40:54-40:55, 40:61-40:63, 40:65-40:67, 40:69, 40:72:01-40:74, 40:76-40:77, 40:79-40:81, 40:84, 40:87:01-40:88, 40:100-40:102, 40:106, 40:108, 40:112-40:114, 40:116-40:118N, 40:121, 40:123, 40:125-40:126, 40:128, 40:132, 40:134-40:141, 40:146-40:147, 40:149-40:156, 40:158, 40:160:01-40:160:02, 40:166, 40:168, 40:170-40:171, 40:175, 40:179, 40:182-40:183, 40:185-40:186, 40:188, 40:191-40:199, 40:204, 40:207-40:208, 40:210, 40:212-40:213, 40:215-40:218, 40:221-40:223, 41:08, 42:05:01-42:05:02, 44:31, 48:01:01-48:01:03, 48:03:01-48:03:02, 48:05-48:07, 48:09-48:12, 48:14-48:16, 48:18-48:22, 48:28 and 81:01-81:06 alleles.

¹⁴Primer mix 10: Specific PCR fragment of 285 bp in the B*07:07, 07:77, *07:12, 07:14, 07:18:01-07:18:02, 07:137, 07:162-07:163 and the B*15:138, 15:230, 15:241, 18:35, 40:01:01-40:01:17, 40:01:19-40:01:29, 40:07, 40:10:01-40:10:02, 40:12, 40:14:01-40:14:03, 40:21-40:23, 40:25, 40:33, 40:36, 40:38, 40:42-40:43, 40:46-40:49, 40:52, 40:54-40:55, 40:59, 40:60-40:63, 40:65-40:67, 40:69, 40:72:01-40:74, 40:76-40:77, 40:79, 40:81, 40:84, 40:87:01-40:88, 40:100-40:102, 40:106, 40:108, 40:112-40:114, 40:116-40:118N, 40:121, 40:123-40:126, 40:128, 40:132, 40:134-40:136, 40:137, 40:138-40:141, 40:146-40:147, 40:149-40:156, 40:158, 40:160:01-40:160:02, 40:166, 40:168, 40:170-40:172, 40:175, 40:179, 40:182-40:183, 40:185-40:188, 40:191-40:199, 40:204, 40:207-40:208, 40:210, 40:212-40:213, 40:215-40:218, 40:221-40:223, 41:08, 44:31, 46:06, 44:150 48:03:01-48:03:02, 48:17 48:23 and in the C*02:23, C*04:77, C*07:38:01-07:38:02 alleles. Specific PCR fragment of 415 bp in the B*07:47, 07:81 and the B*08:52, 38:41, 40:109² and 48:18 alleles

¹⁵Primer mix 11: Specific PCR fragment of 130 bp in B*07:39 allele. Specific PCR fragment of 210 bp in the B*07:08, 07:10 07:32, 07:38 and the B*08:94, 18:07:01-18:07:02, 18:29, 18:79, 27:01-27:05:08, 27:05:10-27:11, 27:13-27:15, 27:17, 27:19-27:21, 27:23, 27:24-27:28, 27:30-27:36, 27:37

101.512-24/03– including *Taq* pol., IFU-01
 101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: 92R**Lot-specific information**

27:38-27:76, 27:78-27:84, 27:86-27:91, 27:93-27:99, 35:50, 35:84, 35:162, 35:197, 37:21, 38:06-38:07, 39:58, 40:08, 40:13, 40:68, 40:162, 40:220, 48:06, 51:57, 51:103, 51:116, 53:21 and 73:01-73:02 alleles

¹⁶Primer mix 13: Specific PCR fragment of 85 bp in the B*07:31, 07:34 and 07:43 and the B*40:199, 41:08 and 42:05:01-42:05:02 and in the C*07:38:01-07:38:02 alleles. Specific PCR fragment of 290 bp in the B*07:12, 07:14, 07:18:01-07:18:02, 07:137 and 07:162-07:163 and the B*40:30, 40:34, 40:59, 40:137, 40:160:01-40:160:02, 44:150 and 48:17 alleles.

¹⁷Primer mix 15: Specific PCR fragment of 385 bp in the B*07:54 and 07:123 and the B*45:06 alleles. Specific PCR fragment of 415 bp in the B *07:11, 07:47, 07:57 and 07:75 and the B*08:17, 08:47, 39:20 and 42:15 alleles.

¹⁸Primer mix 17: Specific PCR fragment of 345 bp in the B*07:46 allele. Specific PCR fragment of 425 bp in the the B*07:11, 07:27, 07:57 and 07:75 and the B*08:17, 08:47, 38:17, 39:20, 40:188² and 42:15 alleles.

¹⁹Primer mix 19: Specific PCR fragment of 205 bp in the B*15:14, 15:91, 15:131, 15:161, 18:56, 35:45, 35:71, 35:169, 40:132, 40:178, 44:10, 44:14, 44:17, 44:31, 44:43:01-44:43:02, 44:45, 44:51, 44:62, 44:77, 44:82, 44:107, 44:123, 45:09, 45:14, 46:17, 51:23, 51:108, 52:15, 53:22, 57:07, 58:07 and 82:01-82:03 alleles. Specific PCR fragment of 250 bp in the B* 44:54, 44:106 and 44:158 alleles. Specific PCR fragment of 205 and 250 bp in the B*07:20, 07:140 and 07:153 and the B*44:135 alleles. Specific PCR fragment of 250 and 295 bp in the B*07:35 allele.

²⁰Primer mix 20: Specific PCR fragment of 160 bp in the B*07:59 allele. Specific PCR fragment of 205 bp in the B*07:51 and the C*08:08 alleles. Specific PCR fragment of 275 bp in the B*07:21 allele.

²¹Primer mix 21: Specific PCR fragment of 100 bp in the B*07:52 allele. Specific PCR fragment of 165 bp in the B*07:22:01-07:22:02 alleles.

²²Primer mix 25: Specific PCR fragment of 305 bp in the B*07:117 allele. Specific PCR fragment of 425 bp in the B*07:27, 38:17 and 40:188² alleles.

²³Primer mix 26: Specific PCR fragment of 95 bp in the B*07:28 and 07:118 and B*08:14, 08:87, 27:32 and 40:180 alleles. Specific PCR fragment of 205 bp in the B*07:51 and C*08:08 alleles. Specific PCR fragment of 235 bp in the B*07:20, 07:24, 07:60, 07:100 and 07:131 and the B*08:21, 08:25, 15:07:01-15:07:02, 15:45, 15:55, 15:68, 15:126, 15:207, 35:05:01-35:05:02, 35:22, 35:51, 35:58, 35:72, 35:89, 35:97, 35:114, 40:05, 40:71, 40:174, 46:12, 51:64, 53:14, 55:51 and 58:18 and in the C*03:05, C*03:25 and C*03:27 alleles. Specific PCR fragment of 290 bp in the B*07:42 allele. Specific PCR fragments of 95 and 235 bp in the B*44:54 and 44:106 44:135 alleles.

²⁴Primer mix 27: Specific PCR fragment of 165 bp in the B*07:41 allele . Specific PCR fragment of 215 bp in the B*07:29 allele.

²⁵Primer mix 29: Specific PCR fragment of 85 bp in the B*07:31, 07:34 and 07:43 and the B*40:199 41:08 and 42:05:01-42:05:02 and in the C*07:38:01-07:38:02 alleles. Specific PCR fragment of 200 bp in the B*07:44 allele.

²⁶Primer mix 31: Specific PCR fragment of 215 bp in the B*07:55 and the B*08:70 and 48:19 and in the C*02:60 alleles.

Specific PCR fragment of 280 bp in the B*07:19, 07:33, 07:53, 07:60 and the B*08:103, 14:05, 14:13, 15:07:01-15:07:02, 15:45, 15:55, 15:68, 15:126, 15:207, 18:14, 27:07:01-27:07:03, 27:11, 27:24, 27:32-27:34, 27:43, 35:05:01-35:05:03, 35:22, 35:31, 35:51, 35:58, 35:72, 35:89, 35:97, 35:114, 35:199, 37:09, 38:19, 39:03, 39:14, 39:24:01-39:24:02, 39:29, 39:37, 39:76, 40:02:01-40:03, 40:05, 40:08-40:09, 40:13, 40:18-40:19, 40:24, 40:27:01-40:27:02, 40:29, 40:35, 40:37, 40:39-40:40, 40:50, 40:56-40:58, 40:71, 40:78, 40:82, 40:85, 40:89-40:91, 40:94, 40:97, 40:104-40:105, 40:107, 40:111, 40:115, 40:119, 40:122, 40:133Q, 40:142N-40:145, 40:157, 40:164, 40:169, 40:173-40:174, 40:176, 40:180-40:181, 40:189, 40:200-40:203, 40:205-40:206, 40:209, 40:211, 40:214, 40:219-40:220, 40:224, 40:226, 42:08, 44:54, 44:106, 44:135, 44:158, 46:12, 48:04, 48:13, 48:24, 51:64, 51:81, 51:148, 53:14, 55:04, 55:49, 55:51, 58:18, 58:27, C*02:60, C*03:05, 03:25, 03:27, 03:167 and in the C*08:08 alleles. Specific PCR fragments of 215 and 280 bp in the B*07:100 and the B*15:07:01-15:07:02, 15:45, 15:68, 15:126, 15:207 and 46:12 alleles.

²⁷Primer mix 33: Specific PCR fragment of 90 bp in the B *07:149 and the B*08:02, 13:36, 15:89, 15:115, 18:09, 38:02:01-38:02:04, 38:08, 38:15, 38:18, 38:23, 38:29, 38:35, 44:12, 51:54, 51:78:01-51:78:02, 53:09, 53:11-53:13, 56:07 alleles. Specific PCR fragment of 350 bp in the B*07:61 allele.

²⁸Primer mix 36: Specific PCR fragment of 180 bp in the B*07:167N. Specific PCR fragment of 180 bp in the B*07:64 allele.

101.512-24/03– including *Taq* pol., IFU-01
 101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: 92R**Lot-specific information**

²⁹Primer mix 38: Specific PCR fragment of 190 bp in the B*07:103 and 07:113 and the B*08:64 alleles. Specific PCR fragment of 215 bp in the B*07:66 allele.

³⁰Primer mix 40: Specific PCR fragment of 265 bp in the B*07:56:01-07:56:02, 07:72 and the B*08:26, 08:56:01-08:56:02, 38:28, 39:63, 40:21[?], 40:58[?], 40:137[?], 41:12 and 49:03 alleles. Specific PCR fragment of 395 bp in the B*07:83 allele.

³¹Primer mix 41: Specific PCR fragment of 90 bp in the B*07:149 and the B*08:02, 13:36, 15:89, 15:115, 18:09, 38:02:01-38:02:04, 38:08, 38:15, 38:18, 38:23, 38:29, 38:35, 44:12, 51:54, 51:78:01-51:78:02, 53:09, 53:11-53:13 and 56:07 alleles. Specific PCR fragment of 415 bp in the B*07:58 and 07:156 alleles

³²Primer mix 42: Specific PCR fragment of 195 bp in the B*07:67N allele. Specific PCR fragment of 225 bp in the B*07:71 allele.

³³Primer mix 48: Specific PCR fragment of 200 bp in the B*07:76 allele. Specific PCR fragment of 280 bp in the B*07:77 and the B*40:60 and C*07:38:01-07:38:02 alleles.

³⁴Primer mix 50: Specific PCR fragment of 75 bp in the B*07:74 and the B*08:77 and 48:27 alleles. Specific PCR fragment of 105 bp in the B*07:175 and the 15:126 alleles. Specific PCR fragment of 165 bp in the B*07:87 allele.

³⁵Primer mix 51: Specific PCR fragment of 345 bp in the B*07:97 allele. Specific PCR fragment of 470 bp in the B*07:93 and the B*08:39 and 38:27 alleles.

³⁶Primer mix 57: Specific PCR fragment of 150 bp in the B*07:109 allele. Specific PCR fragment of 215 bp in the B*07:70 allele.

³⁷Primer mix 58: Specific PCR fragment of 245 bp in the B*07:112 allele. Specific PCR fragment of 440 bp in the B*07:79 allele.

³⁸Primer mix 59: Specific PCR fragment of 95 bp in the B*07:115 and the B*08:69, 14:05, 38:19, 39:03, 39:24:01-39:24:02, 39:37, 39:76, 41:18, 42:06 and 48:07 and in the C*08:08 alleles. Specific PCR fragment of 190 bp in the B*07:89 and 07:113 and the B*08:64 alleles.

³⁹Primer mix 61: Specific PCR fragment of 205 bp in the B*07:88 allele. Specific PCR fragment of 420 bp in the B*07:105 allele.

⁴⁰Primer mix 63: Specific PCR fragment of 65 bp in the B*07:94 and in the C*03:124, C*04:49, C*05:15 and C*07:262 alleles. Specific PCR fragment of 230 bp in the B*07:135N allele.

⁴¹Primer mix 65: Specific PCR fragment of 385 bp in the B*07:54 and 07:123 and the B*45:06 alleles. Specific PCR fragment of 450 bp in the B*07:111N allele.

⁴²Primer mix 68: Specific PCR fragment of 335 bp in the B*07:148 and the B*40:193[?] alleles. Specific PCR fragment of 365 bp in the B*07:101 allele. Specific PCR fragment of 415 bp in the B*07:13, 07:15 and 07:160 and the B*08:16 and 67:02 alleles.

⁴³Primer mix 69: Specific PCR fragment of 190 bp in the B*07:106 allele. Specific PCR fragment of 245 bp in the B*07:108 allele. Specific PCR fragment of 465 bp in the B*07:102 and the B*08:18 and 73:01-73:02 alleles.

⁴⁴Primer mix 70: Specific PCR fragment of 125 bp in the B*07:104 allele. Specific PCR fragment of 205 bp in the B*07:114 allele.

⁴⁵Primer mix 71: Specific PCR fragment of 245 bp in the B*07:108 allele. Specific PCR fragment of 280 bp in the B*07:107 allele.

⁴⁶Primer mix 73: Specific PCR fragment of 85 bp in the B*07:119 and *07:169 alleles. Specific PCR fragment of 240 bp in the B *07:49N allele.

⁴⁷Primer mix 74: Specific PCR fragment of 345 bp in the B*07:120 and the B*40:150 and 57:29 alleles. Specific PCR fragment of 425 bp in the B*07:156 allele.

⁴⁶Primer mix 75: Specific PCR fragment of 105 bp in the B*07:126 and 07:175 and the B*15:126 alleles. Specific PCR fragment of 250 bp in the B*07:121 allele.

⁴⁷Primer mix 76: Specific PCR fragment of 170 bp in the B*07:128 allele. Specific PCR fragment of 285 bp in the B*07:122 allele.

⁴⁸Primer mix 77: Specific PCR fragment of 335 bp in the B*07:148 and the B*40:193[?] alleles. Specific PCR fragment of 405 bp in the B*07:116 allele.

⁴⁹Primer mix 78: Specific PCR fragment of 110 bp in the B*07:125 allele. Specific PCR fragment of 380 bp in the B*07:129 allele.

⁵⁰Primer mix 79: Specific PCR fragment of 65 bp in the B*07:130 allele. Specific PCR fragment of 95 bp in the B*07:119 allele.

‘w’, might be weakly amplified.

‘?’ , nucleotide sequence information not available for the primer matching sequence.

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot-specific information

CE
0088

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot-specific information

For *In Vitro* Diagnostic Use.

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot-specific information

[illegible]

101.512-24/03– including *Taq* pol., IFU-01
101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: 92R

Lot-specific information

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40
*07:156	1	2	3	4	5																																			
*07:161N	1	2	3	4	5																																			
*07:162	1		3	4	5					10		12	13				18																							
*07:167N	1	2	3	4	5																															36				
*07:175	1	2	3	4	5																																			
*08:01:01-08:01:20, 08:05, 08:08N, 08:10, 08:19N, 08:22, 08:24, 08:27, 08:29-08:36, 08:38, 08:41-08:46, 08:48, 08:50-08:51, 08:53:01-08:53:02, 08:58-08:59, 08:61-08:63, 08:65-08:68, 08:72N-08:73, 08:75-08:76, 08:80-08:83, 08:85-08:86N, 08:88, 08:90-08:93, 08:95-08:102, 08:104-08:105																																		32						
*08:02																																		32	33					
*08:03																																		32						
*08:04																																		32						
*08:07, 08:11, 08:28, 08:37, 08:74																																								
*08:09, 08:12:01-08:12:03, 08:20, 08:55, 08:60, 08:84, 39:04																																		32						
*08:13														14																				32						
*08:14, 08:87																										26														
*08:15																																		32						
*08:16																																		32						
*08:17														15			17																32							
*08:18																																		32						
*08:21																										26														
*08:23, 41:02:01-41:02:04, 41:11, 41:13, 41:15, 41:23																																		32						
*08:25														14												26														
*08:26, 08:56:01-08:56:02																																		32					40	
*08:39																																		32						
*08:40																																								
*08:47														15			17																32							
*08:49, 08:89, 14:11, 14:28																																								
*08:52										10																								32						
*08:54																																		32						
*08:57, 35:87																																		32						
*08:64																																		32					38	
*08:69																																		32						
*08:70																																		31	32					
*08:71																																		32						
*08:77																																								
*08:78																																		32						
*08:79		2		4				8																										30						
*08:94						6					11																							32						
*08:103																																		31	32					
*13:01:01-13:01:02, 13:01:04-13:02:15, 13:07N-13:09, 13:11, 13:13, 13:15-13:23, 13:25, 13:27-13:33, 13:37-13:45, 13:47, 13:49N-13:50, 13:52, 13:54-13:61, 13:63N, 35:88, 51:10, 51:16, 51:31, 51:34, 51:82, 51:93, 52:08, 54:26, 55:09, 55:22, 55:24																																		30						
*13:06, 13:53, 15:42, 15:83, 15:86, 15:188, 15:224, 15:252, 35:60, 39:12, 39:52, 39:72, 46:18, 51:21, 51:36, 51:101, 51:143, 54:01:01-54:01:03, 54:05N, 54:07-54:08N, 54:10, 54:12-54:13, 54:17-54:25, 54:27-54:28, 55:01:01-55:03, 55:05, 55:07, 55:11-55:13, 55:15-55:16, 55:19, 55:21, 55:25-55:26, 55:29-55:40, 55:42-55:45, 55:47-55:48, 55:50, 55:52-55:57, 56:10, 56:23, 58:20, 59:01:01-59:03, 59:05, C*15:24																																		32						
*13:14, 40:04, 40:06:01-01-40:06:06, 40:70						6										16																	30							
*13:36																																	30		33					
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40

101.512-24/03– including *Taq* pol., IFU-01
 101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: 92R

Lot-specific information

41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	Well No.			
41																																74						80	*07:156				
																																								*07:161N			
																																								*07:162			
																																								*07:167N			
									50																								75							*07:175			
						47																				67														*08:01:01-08:01:20, 08:05, 08:08N, 08:10, 08:19N, 08:22, 08:24, 08:27, 08:29-08:36, 08:38, 08:41-08:46, 08:48, 08:50-08:51, 08:53:01-08:53:02, 08:58-08:59, 08:61-08:63, 08:65-08:68, 08:72N, 08:73, 08:75-08:76, 08:80-08:83, 08:85-08:86N, 08:88, 08:90-08:93, 08:95-08:102, 08:104-08:105			
41						47																				67															*08:02		
						47															62				66	67														*08:03			
				45		47																				67															*08:04		
						47																				67															*08:07, 08:11, 08:28, 08:37, 08:74		
																										67															*08:09, 08:12:01-08:12:03, 08:20, 08:55, 08:60, 08:84, 39:04		
																										67															*08:13		
						47																				67															*08:14, 08:87		
				46	47																					67															*08:15		
																										67	68															*08:16	
				45		47																				67															*08:17		
						47																				67		69													*08:18		
						47									56											67															*08:21		
						47															62					67															*08:23, 41:02:01-41:02:04, 41:11, 41:13, 41:15, 41:23		
															56											67																*08:25	
						47																				67																*08:26, 08:56:01-08:56:02	
						47			51																	67															*08:39		
																										67																*08:40	
						47																				67																*08:47	
																										67																*08:49, 08:89, 14:11, 14:28	
						47																			66	67															*08:52		
				45		47																				67															*08:54		
						47																				67																*08:57, 35:87	
						47																				67																*08:64	
						47																				67																*08:69	
																										67																*08:70	
				43		47																				67															*08:71		
						47			50																	67															*08:77		
						47																			66	67															*08:78		
																										67																*08:79	
						47																				67																*08:94	
						47																				67																*08:103	
																																											*13:01:01-13:01:02, 13:01:04-13:02:15, 13:07N-13:09, 13:11, 13:13, 13:15-13:23, 13:25, 13:27-13:33, 13:37-13:45, 13:47, 13:49N-13:50, 13:52, 13:54-13:61, 13:63N, 35:88, 51:10, 51:16, 51:31, 51:34, 51:82, 51:93, 52:08, 54:26, 55:09, 55:22, 55:24
																																											*13:06, 13:53, 15:42, 15:83, 15:86, 15:188, 15:224, 15:252, 35:60, 39:12, 39:52, 39:72, 46:18, 51:21, 51:36, 51:101, 51:143, 54:01:01-54:01:03, 54:05N, 54:07-54:08N, 54:10, 54:12-54:13, 54:17-54:25, 54:27-54:28, 55:01:01-55:03, 55:05, 55:07, 55:11-55:13, 55:15-55:16, 55:19, 55:21, 55:25-55:26, 55:29-55:40, 55:42-55:45, 55:47-55:48, 55:50, 55:52-55:57, 56:10, 56:23, 58:20, 59:01:01-59:03, 59:05, C*15:24
																																											*13:14, 40:04, 40:06:01-01-40:06:06, 40:70
41																																											*13:36
41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	Well No.			

101.512-24/03– including *Taq* pol., IFU-01
 101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: 92R

Lot-specific information

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40			
*14:01:01-14:01:02, 14:07N-14:08, 14:10, 14:12, 14:19, 14:26, 14:32, 27:12, 27:16, 39:15, 39:25N, 39:32-39:33, 39:38Q, 39:40:01N-39:40:02N, 39:43, 39:67, 39:74						6																																					
*14:05													14																		31												
*14:13													14																		31	32											
*14:14, 39:07, 39:55						6										16																											
*14:18																																											
*14:21	1																																										
*15:03:01-15:03:03, 15:47:01-15:47:02, 15:54, 15:61-15:62, 15:64:01-15:64:02, 15:69, 15:98, 15:103, 15:123, 15:127, 15:132, 15:151, 15:156, 15:158, 15:173, 15:210, 15:220, 15:235, 15:242, 15:253, 15:259, 44:40, 44:44, 44:156, 52:16																																											
*15:07:01-15:07:02, 15:45, 15:207, 46:12																										26					31												
*15:14, 15:161, 35:45, 35:71, 35:169, 44:10, 44:14, 44:17, 44:43:01-44:43:02, 44:45, 44:51, 44:62, 44:77, 44:82, 44:107, 44:123, 46:17, 51:23, 51:108, 52:15, 53:22, 57:07, 58:07, 82:01-82:03																		19																									
*15:55, 35:05:01-35:05:03, 35:22, 35:51, 35:58, 35:72, 35:89, 35:97, 35:199, 51:64, 53:14, 58:18													14													26					31												
*15:68						6										16										26					31												
*15:71, 15:175, 18:01:01:01-18:01:01:02, 18:01:03-18:01:12, 18:01:14-18:02, 18:05-18:06, 18:08, 18:10-18:13, 18:15, 18:17N-18:21, 18:23N-18:24, 18:27-18:28, 18:30-18:32, 18:34, 18:36, 18:38-18:55, 18:57-18:60, 18:62-18:67, 18:69, 18:71-18:78, 18:80-18:83, 37:01:01-37:06, 37:08, 37:10-37:11, 37:13-37:20, 37:22-37:27, 37:29-37:36, 40:20, 44:103, 44:131						6										16																											
*15:89, 15:115, 44:12, 51:54, 51:78:01-51:78:02, 53:09, 53:11-53:13, 56:07																																											
*15:91, 15:131																		19																									
*15:126																										26					31												
*15:138, 15:241			3							10	12																																
*15:230										10														23																			
*18:01:02, 18:03, 37:37, 39:17						6																																					
*18:04, 35:08:05, 35:42:02, 35:185, 78:01:01-78:01:02, 78:02:02-78:03																																										37	
*18:07:01-18:07:02, 18:29, 18:79, 27:23, 35:50, 35:84, 35:162, 35:197, 37:21, 51:57, 51:103, 51:116, 53:21						6					11					16																											
*18:09						6										16																									33		
*18:14						6								14		16																31											
*18:22						6								14		16																											
*18:26, 18:61, 37:12						6										16																	32										
*18:33, 27:92						6										16																											
*18:35			3			6				10	12					16								24																			
*18:37, 27:18, 27:29, 37:28, 44:130						6										16																											
*18:56						6										16			19																								
*27:01, 27:44						5					11					16																											
*27:02:01-27:05:08, 27:05:10-27:06, 27:08-27:10, 27:13-27:15, 27:17, 27:19, 27:21, 27:25-27:26, 27:28, 27:30-27:31, 27:35, 27:37-27:38, 27:40-27:42, 27:45-27:58, 27:60-27:74, 27:76, 27:78-27:84, 27:86-27:91, 27:93-27:99						5					11																																
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40			

101.512-24/03– including *Taq* pol., IFU-01
 101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **92R**

Lot-specific information

41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	Well No.							
																						64																		*14:01:01-14:01:02, 14:07N-14:08, 14:10, 14:12, 14:19, 14:26, 14:32, 27:12, 27:16, 39:15, 39:25N, 39:32-39:33, 39:38Q, 39:40:01N-39:40:02N, 39:43, 39:67, 39:74							
		43															59																								*14:05						
		43																			62	64																			*14:13						
																					62																					*14:14, 39:07, 39:55					
																																											*14:18				
																																											*14:21				
				45																																							*15:03:01-15:03:03, 15:47:01-15:47:02, 15:54, 15:61-15:62, 15:64:01-15:64:02, 15:69, 15:98, 15:103, 15:123, 15:127, 15:132, 15:151, 15:156, 15:158, 15:173, 15:210, 15:220, 15:235, 15:242, 15:253, 15:259, 44:40, 44:44, 44:156, 52:16				
		43													56																												*15:07:01-15:07:02, 15:45, 15:207, 46:12				
																																												*15:14, 15:161, 35:45, 35:71, 35:169, 44:10, 44:14, 44:17, 44:43:01-44:43:02, 44:45, 44:51, 44:62, 44:77, 44:82, 44:107, 44:123, 46:17, 51:23, 51:108, 52:15, 53:22, 57:07, 58:07, 82:01-82:03			
		43													56																												*15:55, 35:05:01-35:05:03, 35:22, 35:51, 35:58, 35:72, 35:89, 35:97, 35:199, 51:64, 53:14, 58:18				
		43	45											56							62																							*15:68			
																																												*15:71, 15:175, 18:01:01:01-18:01:01:02, 18:01:03-18:01:12, 18:01:14-18:02, 18:05-18:06, 18:08, 18:10-18:13, 18:15, 18:17N-18:21, 18:23N-18:24, 18:27-18:28, 18:30-18:32, 18:34, 18:36, 18:38-18:55, 18:57-18:60, 18:62-18:67, 18:69, 18:71-18:78, 18:80-18:83, 37:01:01-37:06, 37:08, 37:10-37:11, 37:13-37:20, 37:22-37:27, 37:29-37:36, 40:20, 44:103, 44:131			
41																																												*15:89, 15:115, 44:12, 51:54, 51:78:01-51:78:02, 53:09, 53:11-53:13, 56:07			
				45																																										*15:91, 15:131	
		43						50						56																															*15:126		
												53																																		*15:138, 15:241	
																																														*15:230	
																																														*18:01:02, 18:03, 37:37, 39:17	
																																														*18:04, 35:08:05, 35:42:02, 35:185, 78:01:01-78:01:02, 78:02:02-78:03	
																																														*18:07:01-18:07:02, 18:29, 18:79, 27:23, 35:50, 35:84, 35:162, 35:197, 37:21, 51:57, 51:103, 51:116, 53:21	
41																																															*18:09
		43																																													*18:14
																																														*18:22	
																																														*18:26, 18:61, 37:12	
																							64																							*18:33, 27:92	
																																														*18:35	
				45																																										*18:37, 27:18, 27:29, 37:28, 44:130	
																																														*18:56	
																																														*27:01, 27:44	
																																														*27:02:01-27:05:08, 27:05:10-27:06, 27:08-27:10, 27:13-27:15, 27:17, 27:19, 27:21, 27:25-27:26, 27:28, 27:30-27:31, 27:35, 27:37-27:38, 27:40-27:42, 27:45-27:58, 27:60-27:74, 27:76, 27:78-27:84, 27:86-27:91, 27:93-27:99	
41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	Well No.							

101.512-24/03– including *Taq* pol., IFU-01
101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: 92R

Lot-specific information

[illegible]

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot-specific information

CE
0088

101.512-24/03– including *Taq* pol., IFU-01
101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: **92R**

Lot-specific information

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40
*55:08, 56:05:02, 56:15, 56:21																																								39
*55:10					5																												32							
*55:18																																32								
*55:51													14													26					31								39	
*56:05:01, 56:06																																						37	39	
*57:10																16																								
*57:29																																								
*67:01:01-67:01:02, 67:03	1				5																												32							
*67:02					5																												32							
*73:01-73:02	w				5						11																						32							
*78:07																																								
*81:01, 81:03-81:04N		2			5			8						14											24						30								39	
*81:02	1	2			5			8						14											24						30								39	
*81:05		2			5			8						14											24						30									
*81:06		2			5			8						14											24						30								39	
C*02:23			3							10								18																						
C*02:60																																	31							
C*03:124, C*04:49, C*05:15, C*07:262																																								
C*03:137, C*12:82					5																																			
C*04:77				4						10								18																						
C*07:38:01-07:38:02										10		12	13				18														29									
C*07:102					w																																			
C*08:08														14						20							26				31									
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40

¹The internal positive control primer pairs amplify segments of the human growth hormone gene. The two different control primer pairs give rise to either an internal positive control band of 1070 base pairs, for most wells, or a band of 800 base pairs, for some wells.

Well number 1 contains the primer pair giving rise to the shorter, 800 bp, internal positive control band in order to help in the correct orientation of the HLA-B*07 subtyping.

In addition, wells number 4, 10, 14, 19, 24, 26, 27, 29 to 31, 47, 52, 53, 59, 62, 67, 70 and 77 contain the primer pair giving rise to the shorter, 800 bp, internal positive control band in order to allow kit identification.

In the presence of a specific amplification the intensity of the control band often decreases.

²The sequence of B*0701 allele has been shown to be in error.

³HLA-B*07 alleles in bold lettering are listed as confirmed alleles on the IMGT/HLA web page www.ebi.ac.uk/imgt/hla, release 3.10.0, October 2012.

⁴The B*07:50 and B*48:08 alleles will give rise to identical amplification pattern with the HLA-B*07 subtyping kit. These two alleles can be distinguished e.g with the HLA-B low resolution kit.

⁵The B*07:42 and 07:118 alleles may be distinguished by the different sizes of the specific PCR products generated by primer mix 26.

⁶The B*07:49N and 07:169 alleles may be distinguished by the different sizes of the specific PCR products generated by primer mix 73.

⁷The B*07:94 and 07:135 alleles may be distinguished by the different sizes of the specific PCR products generated by primer mix 63.

⁸The B*07:102 and 07:106 alleles can be distinguished by the different sizes of the specific PCR products generated by primer mix 69.

⁹The B*07:121 and 07:126 alleles may be distinguished by the different sizes of the specific PCR products generated by primer mix 75.

¹⁰Primer mix 7: Specific PCR fragment of 105 bp in the B*07:04, 07:25, 07:52 and 07:146 and the B*40:136, 41:08 and 42:05:01-42:05:02 alleles. Specific PCR fragment of 215 bp in the B*07:29 allele.

Primer mix 8: Specific PCR fragment of 160 bp in the B*07:59 allele. Specific PCR fragment of 235 bp in the B*07:05:01-07:06, 07:32, 07:34, 07:69, 07:78, 07:80, 07:90, 07:97, 07:105, 07:112, 07:123, 07:137-07:138, 07:140, 07:176 and the B*08:79, 40:01:01-40:01:20, 40:01:22-40:01:25, 40:01:27, 40:01:29, 40:07, 40:10:01-40:10:02, 40:12, 40:14:01-40:16, 40:21-40:23, 40:25, 40:30, 40:34, 40:36, 40:38, 40:43, 40:45-40:49, 40:54-40:55, 40:61-40:63, 40:65-40:67, 40:69, 40:72:01-40:74, 40:76-40:77, 40:79-40:81, 40:84, 40:87:01-40:88, 40:100-40:102, 40:106, 40:108, 40:112-40:114, 40:116-40:118N, 40:121, 40:123, 40:125-40:126, 40:128, 40:132, 40:134-40:141, 40:146-40:147, 40:149-40:156, 40:158, 40:160:01-40:160:02, 40:166, 40:168, 40:170-40:171, 40:175, 40:179, 40:182-40:183, 40:185-40:186, 40:188, 40:191-40:199, 40:204, 40:207-40:208, 40:210, 40:212-40:213, 40:215-40:218, 40:221-40:223, 41:08, 42:05:01-42:05:02, 44:31, 48:01:01-48:01:03, 48:03:01-48:03:02, 48:05-48:07, 48:09-48:12, 48:14-48:16, 48:18-48:22, 48:28 and 81:01-81:06 alleles.

101.512-24/03– including *Taq* pol., IFU-01
101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: **92R**

Lot-specific information

41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	Well No.
																																								*55:08, 56:05:02, 56:15, 56:21
																																								*55:10
																																								*55:18
																																								*55:51
																																								*56:05:01, 56:06
																																								*57:10
																																								*57:29
																																								*67:01:01-67:01:02, 67:03
																																								*67:02
																																								*73:01-73:02
																																								*78:07
																																								*81:01, 81:03-81:04N
																																								*81:02
																																								*81:05
																																								*81:06
																																								C*02:23
																																								C*02:60
																																								C*03:124, C*04:49, C*05:15,
																																								C*07:262
																																								C*03:137, C*12:82
																																								C*04:77
																																								C*07:38:01-07:38:02
																																								C*07:102
																																								C*08:08
41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	Well No.

Primer mix 11: Specific PCR fragment of 130 bp in B*07:39 allele. Specific PCR fragment of 210 bp in the B*07:08, *07:10 07:32, 07:38 and the B*08:94, 18:07:01-18:07:02, 18:29, 18:79, 27:01-27:05:08, 27:05:10-27:11, 27:13-27:15, 27:17, 27:19-27:21, 27:23, 27:24-27:28, 27:30-27:36, 27:37 27:38-27:76, 27:78-27:84, 27:86-27:91, 27:93-27:99, 35:50, 35:84, 35:162, 35:197, 37:21, 38:06-38:07, 39:58, 40:08, 40:13, 40:68, 40:162, 40:220, 48:06, 51:57, 51:103, 51:116, 53:21 and 73:01-73:02 alleles

Primer mix 13: Specific PCR fragment of 85 bp in the B*07:31, 07:34 and 07:43 and the B*40:199, 41:08 and 42:05:01-42:05:02 and in the C*07:38:01-07:38:02 alleles. Specific PCR fragment of 290 bp in the B*07:12, 07:14, 07:18:01-07:18:02, 07:137 and 07:162-07:163 and the B*40:30, 40:34, 40:59, 40:137, 40:160:01-40:160:02, 44:150 and 48:17 alleles.

Primer mix 15: Specific PCR fragment of 385 bp in the B*07:54 and 07:123 and the B*45:06 alleles. Specific PCR fragment of 415 bp in the B *07:11, 07:47, 07:57 and 07:75 and the B*08:17, 08:47, 39:20 and 42:15 alleles.

Primer mix 17: Specific PCR fragment of 345 bp in the B*07:46 allele. Specific PCR fragment of 425 bp in the the B*07:11, 07:27, 07:57 and 07:75 and the B*08:17, 08:47, 38:17, 39:20, 40:188² and 42:15 alleles.

Primer mix 19: Specific PCR fragment of 205 bp in the B*15:14, 15:91, 15:131, 15:161, 18:56, 35:45, 35:71, 35:169, 40:132, 40:178, 44:10, 44:14, 44:17, 44:31, 44:43:01-44:43:02, 44:45, 44:51, 44:62, 44:77, 44:82, 44:107, 44:123, 44:135, 45:09, 45:14, 46:17, 51:23, 51:108, 52:15, 53:22, 57:07, 58:07 and 82:01-82:03 alleles.

Specific PCR fragment of 250 bp in the B* 44:54, 44:106, 44:135 and 44:158 alleles. Specific PCR fragment of 205 and 250 bp in the B*07:20, 07:140 and 07:153 and the B*44:135 alleles. Specific PCR fragment of 250 and 295 bp in the B*07:35 allele.

Primer mix 20: Specific PCR fragment of 160 bp in the B*07:59 allele. Specific PCR fragment of 205 bp in the B*07:51 and the C*08:08 alleles. Specific PCR fragment of 275 bp in the B*07:21 allele.

Primer mix 21: Specific PCR fragment of 100 bp in the B*07:52 allele. Specific PCR fragment of 165 bp in the B*07:22:01-07:22:02 alleles.

Primer mix 25: Specific PCR fragment of 305 bp in the B*07:117 allele. Specific PCR fragment of 425 bp in the B*07:27, 38:17 and 40:188² alleles.

Primer mix 26: Specific PCR fragment of 95 bp in the B*07:28 and 07:118 and B*08:14, 08:87, 27:32 and 40:180 alleles. Specific PCR fragment of 205 bp in the B*07:51 and C*08:08 alleles. Specific PCR fragment of 235 bp in the B*07:20, 07:24, 07:60, 07:100 and 07:131 and the B*08:21, 08:25, 15:07:01-15:07:02, 15:45, 15:55, 15:68, 15:126, 15:207, 35:05:01-35:05:02, 35:22, 35:51, 35:58, 35:72, 35:89, 35:97, 35:114, 40:05, 40:71, 40:174, 46:12, 51:64, 53:14, 55:51 and 58:18 and in the C*03:05, C*03:25 and C*03:27 alleles.

Specific PCR fragment of 290 bp in the B*07:42 allele. Specific PCR fragments of 95 and 235 bp in the B*44:54 and 44:106 44:135 alleles.

101.512-24/03– including *Taq* pol., IFU-01
101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: 92R**Lot-specific information**

Primer mix 27: Specific PCR fragment of 165 bp in the B*07:41 allele . Specific PCR fragment of 215 bp in the B*07:29 allele.

Primer mix 29: Specific PCR fragment of 85 bp in the B*07:31, 07:34 and 07:43 and the B*40:199 41:08 and 42:05:01-42:05:02 and in the C*07:38:01-07:38:02 alleles. Specific PCR fragment of 200 bp in the B*07:44 allele.

Primer mix 31: Specific PCR fragment of 215 bp in the B*07:55 and the B*08:70 and 48:19 and in the C*02:60 alleles.

Specific PCR fragment of 280 bp in the B*07:19, 07:33, 07:53, 07:60 and the B*08:103, 14:05, 14:13, 15:07:01-15:07:02, 15:45, 15:55, 15:68, 15:126, 15:207, 18:14, 27:07:01-27:07:03, 27:11, 27:24, 27:32-27:34, 27:43, 35:05:01-35:05:03, 35:22, 35:31, 35:51, 35:58, 35:72, 35:89, 35:97, 35:114, 35:199, 37:09, 38:19, 39:03, 39:14, 39:24:01-39:24:02, 39:29, 39:37, 39:76, 40:02:01-40:03, 40:05, 40:08-40:09, 40:13, 40:18-40:19, 40:24, 40:27:01-40:27:02, 40:29, 40:35, 40:37, 40:39-40:40, 40:50, 40:56-40:58, 40:71, 40:78, 40:82, 40:85, 40:89-40:91, 40:94, 40:97, 40:104-40:105, 40:107, 40:111, 40:115, 40:119, 40:122, 40:133Q, 40:142N-40:145, 40:157, 40:164, 40:169, 40:173-40:174, 40:176, 40:180-40:181, 40:189, 40:200-40:203, 40:205-40:206, 40:209, 40:211, 40:214, 40:219-40:220, 40:224, 40:226, 42:08, 44:54, 44:106, 44:135, 44:158, 46:12, 48:04, 48:13, 48:24, 51:64, 51:81, 51:148, 53:14, 55:04, 55:49, 55:51, 58:18, 58:27, C*02:60, C*03:05, 03:25, 03:27, 03:167 and in the C*08:08 alleles. Specific PCR fragments of 215 and 280 bp in the B*07:100 and the B*15:07:01-15:07:02, 15:45, 15:68, 15:126, 15:207 and 46:12 alleles.

Primer mix 33: Specific PCR fragment of 90 bp in the B *07:149 and the B*08:02, 13:36, 15:89, 15:115, 18:09, 38:02:01-38:02:04, 38:08, 38:15, 38:18, 38:23, 38:29, 38:35, 44:12, 51:54, 51:78:01-51:78:02, 53:09, 53:11-53:13, 56:07 alleles. Specific PCR fragment of 350 bp in the B*07:61 allele.

Primer mix 36: Specific PCR fragment of 180 bp in the B*07:167N. Specific PCR fragment of 180 bp in the B*07:64

Primer mix 38: Specific PCR fragment of 190 bp in the B*07:103 and 07:113 and the B*08:64 alleles. Specific PCR fragment of 215 bp in the B*07:66 allele.

Primer mix 40: Specific PCR fragment of 265 bp in the B*07:56:01-07:56:02, 07:72 and the B*08:26, 08:56:01-08:56:02, 38:28, 39:63, 40:21[?], 40:58[?], 40:137[?], 41:12 and 49:03 alleles. Specific PCR fragment of 395 bp in the B*07:83 allele.

Primer mix 41: Specific PCR fragment of 90 bp in the B*07:149 and the B*08:02, 13:36, 15:89, 15:115, 18:09, 38:02:01-38:02:04, 38:08, 38:15, 38:18, 38:23, 38:29, 38:35, 44:12, 51:54, 51:78:01-51:78:02, 53:09, 53:11-53:13 and 56:07 alleles. Specific PCR fragment of 415 bp in the B*07:58 and *07:156 alleles

Primer mix 42: Specific PCR fragment of 195 bp in the B*07:67N allele. Specific PCR fragment of 225 bp in the B*07:71 allele.

Primer mix 48: Specific PCR fragment of 200 bp in the B*07:76 allele. Specific PCR fragment of 280 bp in the B*07:77 and the B*40:60 and C* 07:38:01-07:38:02 alleles.

Primer mix 50: Specific PCR fragment of 75 bp in the B*07:74 and the B*08:77 and 48:27 alleles. Specific PCR fragment of 105 bp in the B*07:175 and the 15:126 alleles. Specific PCR fragment of 165 bp in the B*07:87 allele.

Primer mix 51: Specific PCR fragment of 345 bp in the B*07:97 allele. Specific PCR fragment of 470 bp in the B*07:93 and the B*08:39 and 38:27 alleles.

Primer mix 57: Specific PCR fragment of 150 bp in the B*07:109 allele. Specific PCR fragment of 215 bp in the B*07:70 allele.

Primer mix 58: Specific PCR fragment of 245 bp in the B*07:112 allele. Specific PCR fragment of 440 bp in the B*07:79 allele.

Primer mix 59: Specific PCR fragment of 95 bp in the B*07:115 and the B*08:69, 14:05, 38:19, 39:03, 39:24:01-39:24:02, 39:37, 39:76, 41:18, 42:06 and 48:07 and in the C*08:08 alleles. Specific PCR fragment of 190 bp in the B*07:89 and 07:113 and the B*08:64 alleles.

Primer mix 61: Specific PCR fragment of 205 bp in the B*07:88 allele. Specific PCR fragment of 420 bp in the B*07:105 allele.

Primer mix 63: Specific PCR fragment of 65 bp in the B*07:94 and in the C*03:124, C*04:49, C*05:15 and C*07:262 alleles. Specific PCR fragment of 230 bp in the B*07:135N allele.

Primer mix 65: Specific PCR fragment of 385 bp in the B*07:54 and 07:123 and the B*45:06 alleles. Specific PCR fragment of 450 bp in the B*07:111N allele.

Primer mix 68: Specific PCR fragment of 335 bp in the B*07:148 and the B*40:193[?] alleles

101.512-24/03– including *Taq* pol., IFU-01
101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: 92R**Lot-specific information**

Specific PCR fragment of 365 bp in the B*07:101 allele. Specific PCR fragment of 415 bp in the B*07:13, 07:15 and 07:160 and the B*08:16 and 67:02 alleles.

Primer mix 69: Specific PCR fragment of 190 bp in the B*07:106 allele. Specific PCR fragment of 245 bp in the B*07:108 allele. Specific PCR fragment of 465 bp in the B*07:102 and the B*08:18 and 73:01-73:02 alleles.

Primer mix 70: Specific PCR fragment of 125 bp in the B*07:104 allele. Specific PCR fragment of 205 bp in the B*07:114 allele.

Primer mix 71: Specific PCR fragment of 245 bp in the B*07:108 allele. Specific PCR fragment of 280 bp in the B*07:107 allele.

Primer mix 73: Specific PCR fragment of 85 bp in the B*07:119 and *07:169 alleles. Specific PCR fragment of 240 bp in the B*07:49N allele.

Primer mix 74: Specific PCR fragment of 345 bp in the B*07:120 and the B*40:150 and 57:29 alleles. Specific PCR fragment of 425 bp in the B*07:156 allele.

Primer mix 75: Specific PCR fragment of 105 bp in the B*07:126 and 07:175 and the B*15:126 alleles. Specific PCR fragment of 250 bp in the B*07:121 allele.

Primer mix 76: Specific PCR fragment of 170 bp in the B*07:128 allele. Specific PCR fragment of 285 bp in the B*07:122 allele.

Primer mix 77: Specific PCR fragment of 335 bp in the B*07:148 and the B*40:193[?] alleles. Specific PCR fragment of 405 bp in the B*07:116 allele.

Primer mix 78: Specific PCR fragment of 110 bp in the B*07:125 allele. Specific PCR fragment of 380 bp in the B*07:129 allele.

Primer mix 79: Specific PCR fragment of 65 bp in the B*07:130 allele. Specific PCR fragment of 95 bp in the B*07:119 allele.

‘w’, might be weakly amplified.

‘?’, nucleotide sequence information not available for the primer matching sequence.

101.512-24/03– including *Taq* pol., IFU-01
 101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **92R**

Lot-specific information

Primers

Well No.	1	2	3	4	5	6	7	8	9	10	11	12
Length of spec. PCR product(s)	395	285	235	105	220	215	105	160	355	285	130	225
							215	235		415	210	
Length of int. pos. control ¹	800	1070	1070	800	1070	1070	1070	1070	1070	800	1070	1070
5'-primer(s) ²	44	363	412	540	103	103	97	412	693	44	103	419
	5'-ggC 3'	5'-AgC 3'	5'-ATg 3'	5'-Cgg 3'	5'-CCT 3'	5'-CCT 3'	5'-TCg 3'	5'-ATA 3'	5'-ggT 3'	5'-ggC 3'	5'-CCT 3'	5'-gTC 3'
					103	103	540	636		363	103	
					5'-CCT 3'	5'-CCT 3'	5'-gAC 3'	5'-CAC 3'		5'-Agg 3'	5'-CCT 3'	
							545			363		
							5'-AgT 3'			5'-Agg 3'		
3'-primer(s) ³	272	605	605	605	282	277	272	605	916	287	191	605
	5'-TgT 3'	5'-gCT 3'	5'-gCT 3'	5'-gCT 3'	5'-gCC 3'	5'-ggT 3'	5'-TgT 3'	5'-gCT 3'	5'-gAT 3'	5'-TCC 3'	5'-CgA 3'	5'-gCT 3'
							605	754		302	272	
							5'-gCT 3'	5'-AgA 3'		5'-ggC 3'	5'-TgC 3'	
										605	272	
										5'-gCT 3'	5'-TgA 3'	
Well No.	1	2	3	4	5	6	7	8	9	10	11	12

Well No.	13	14	15	16	17	18	19	20	21	22	23	24
Length of spec. PCR product(s)	85	215	385	230	345	265	205	160	100	385	85	125
	290		415		425		250	205	165			
							295	275				
Length of int. pos. control ¹	1070	800	1070	1070	1070	1070	800	1070	1070	1070	1070	800
5'-primer(s) ²	355	363	44	103	44	379	363	25	481	322	560	103
	5'-TCA 3'	5'-AgC 3'	5'-ggC 3'	5'-CCT 3'	5'-ggC 3'	5'-ACC 3'	5'-AgC 3'	5'-CCg 3'	5'-CgA 3'	5'-gCC 3'	5'-CCT 3'	5'-CCT 3'
	560	363		103			409	363	545			103
	5'-CAC 3'	5'-AgC 3'		5'-CCT 3'			5'-ggC 3'	5'-AgC 3'	5'-AgT 3'			5'-CCT 3'
								636				527
								5'-CAC 3'				5'-TgT 3'
3'-primer(s) ³	605	539	259	292	218	605	572	131	605	419	605	191
	5'-gCT 3'	5'-TCA 3'	5'-CTC 3'	5'-gTA 3'	5'-gCA 3'	5'-gCT 3'	5'-gCg 3'	5'-ggA 3'	5'-gCT 3'	5'-CgT 3'	5'-gCT 3'	5'-CgA 3'
			287	292	301		616	527				605
			5'-TCC 3'	5'-gTA 3'	5'-gTC 3'		5'-AgA 3'	5'-CCg 3'				5'-gCT 3'
			302	292	302			754				
			5'-ggT 3'	5'-gTA 3'	5'-ggT 3'			5'-AgA 3'				
Well No.	13	14	15	16	17	18	19	20	21	22	23	24

Well No.	25	26	27	28	29	30	31	32	33	34	35	36
Length of spec. PCR product(s)	305	95	165	110	85	190	215	190	90	155	105	180
	425	205	215		200		280		350			250
		235										
		290										
Length of int. pos. control ¹	1070	800	800	1070	800	800	800	1070	1070	1070	1070	1070
5'-primer(s) ²	44	363	97	215	560	412	363	412	261	483	207	363
	5'-ggC 3'	5'-AgC 3'	5'-TCg 3'	5'-gCC 3'	5'-CAC 3'	5'-ATA 3'	5'-AgC 3'	5'-ATA 3'	5'-AAC 3'	5'-gAA 3'	5'-gAC 3'	5'-AgC 3'
			148		852				756	503		
			5'-TgC 3'		5'-CAg 3'				5'-ACT 3'	5'-CCg 3'		
3'-primer(s) ³	181	418	272	282	605	559	538	559	311	605	272	502
	5'-gTT 3'	5'-gTC 3'	5'-TgT 3'	5'-gCC 3'	5'-gCT 3'	5'-CTC 3'	5'-CCA 3'	5'-CgT 3'	5'-Cgg 3'	5'-gCT 3'	5'-TgT 3'	5'-CTA 3'
	301	527			916		603		970			570
	5'-gTC 3'	5'-CCg 3'			5'-gAC 3'		5'-gTC 3'		5'-TCT 3'			5'-CCg 3'
		559										
		5'-CAg 3'										
		613										
		5'-gCT 3'										
Well No.	25	26	27	28	29	30	31	32	33	34	35	36

101.512-24/03– including *Taq* pol., IFU-01
 101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **92R**

Lot-specific information

Well No.	37	38	39	40	41	42	43	44	45	46	47	48
Length of	210	190	425	265	90	195	140	150	105	425	215	200
spec. PCR		215		395	415	225						280
product(s)												
Length of int.	1070	1070	1070	1070	1070	1070	1070	1070	1070	1070	800	1070
pos. control ¹												
5'-primer(s) ²	141 5'-ATT 3'	363 5'-AgC 3'	272 5'-CTA 3'	44 5'-ggC 3'	261 5'-AAC 3'	89 5'-gAC 3'	363 5'-AgC 3'	165 5'-ACC 3'	206 5'-AgA 3'	44 5'-ggC 3'	363 5'-AgC 3'	113 5'-CCA 3'
					625 5'-CCg 3'	118 5'-C.g 3'						368 5'-gTC 3'
					650 5'-CCT 3'							
3'-primer(s) ³	311 5'-ggT 3'	506 5'-Tgg 3'	412 5'-gTT 3'	142 5'-TgC 3'	311 5'-Cgg 3'	272 5'-TgT 3'	463 5'-gCT 3'	272 5'-TgC 3'	272 5'-Tgg 3'	299 5'-TCA 3'	538 5'-gTC 3'	272 5'-TgT 3'
		516 5'-gCg 3'		272 5'-TgA 3'	916 5'-gAC 3'							605 5'-gCT 3'
		539 5'-TCT 3'										
Well No.	37	38	39	40	41	42	43	44	45	46	47	48

Well No.	49	50	51	52	53	54	55	56	57	58	59	60
Length of	200	75	345	290	235	190	260	230	150	245	95	145
spec. PCR		105	470						215	440	190	
product(s)		165										
Length of int.	1070	1070	1070	800	800	1070	1070	1070	1070	1070	800	1070
pos. control ¹												
5'-primer(s) ²	341 5'-ggg 3'	363 5'-AgC 3'	44 5'-ggC 3'	355 5'-CCT 3'	409 5'-ggT 3'	363 5'-AgC 3'	386 5'-gCA 3'	363 5'-AgC 3'	98 5'-CTg 3'	44 5'-ggC 3'	363 5'-AgC 3'	167 5'-gCT 3'
									161 5'-Cgg 3'			
	2 nd 5'-TCg 3'	398 5'-ggC 3'	221 5'-ACC 3'	605 5'-gCT 3'	605 5'-gCT 3'	502 5'-CTA 3'	605 5'-gCT 3'	544 5'-ggT 3'	272 5'-TgT 3'	121 5'-gCC 3'	419 5'-CgA 3'	272 5'-TgT 3'
3'-primer(s) ³		425 5'-CgC 3'	343 5'-T 3'			513 5'-TCA 3'		559 5'-CAg 3'		317 5'-ggC 3'	506 5'-Tgg 3'	
		488 5'-CCA 3'									514 5'-CTT 3'	
Well No.	49	50	51	52	53	54	55	56	57	58	59	60

Well No.	61	62	63	64	65	66	67	68	69	70	71	72
Length of	205	415	65	210	385	435	230	335	190	125	245	165
spec. PCR	420		230		450			365	245	205	280	
product(s)								415	465			
Length of int.	1070	800	1070	1070	1070	1070	800	1070	1070	800	1070	1070
pos. control ¹												
5'-primer(s) ²	44 5'-ggC 3'	44 5'-ggC 3'	256 5'-ACT 3'	103 5'-CCT 3'	44 5'-ggC 3'	44 5'-ggC 3'	44 5'-ggC 3'	44 5'-ggC 3'	44 5'-ggC 3'	106 5'-CCA 3'	44 5'-ggC 3'	103 5'-CCT 3'
			363 5'-AgC 3'	103 5'-CCT 3'					363 5'-AgC 3'	188 5'-Cgg 3'	363 5'-AgC 3'	103 5'-CCT 3'
	82 5'-ggC 3'	292 5'-gTA 3'	282 5'-GCC 3'	272 5'-TgC 3'	259 5'-CTC 3'	309 5'-ATC 3'	106 5'-CAT 3'	212 5'-gCA 3'	119 5'-ggT 3'	272 5'-TgT 3'	119 5'-ggT 3'	226 5'-CAC 3'
3'-primer(s) ³	293 5'-ggA 3'	292 5'-gTA 3'	553 5'-CTA 3'		327 5'-TTT 3'			241 5'-CgA 3'	341 5'-CgT 3'		603 5'-TTC 3'	
		292 5'-gTA 3'						289 5'-AgC 3'	511 5'-CCg 3'			
Well No.	61	62	63	64	65	66	67	68	69	70	71	72

101.512-24/03– including *Taq* pol., IFU-01
 101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **92R**

Lot-specific information

Well No.	73	74	75	76	77	78	79	80
Length of	85	345	105	170	335	110	65	405
spec. PCR	240	425	250	285	405	380	95	
product(s)								
Length of int.	1070	1070	1070	1070	800	1070	1070	1070
pos. control ¹								
5'-primer(s) ²	244	625	363	363	44	206	244	693
	5'-CgC 3'	5'-CCg 3'	5'-AgC 3'	5'-AgC 3'	5'-ggC 3'	5'-AgA 3'	5'-CgC 3'	5'-ggT 3'
	693	704		880		671	693	
	5'-ggT 3'	5'-TgT 3'		5'-CCT 3'		5'-CAg 3'	5'-ggT 3'	
3'-primer(s) ³	299	916	425	607	212	277	299	962
	5'-TCT 3'	5'-gAC 3'	5'-CgC 3'	5'-CAC 3'	5'-gCA 3'	5'-ggT 3'	5'-TCT 3'	5'-CAA 3'
	728		427	916	278	916	719	
	5'-CCT 3'		5'-gTT 3'	5'-gAC 3'	5'-ggA 3'	5'-gAC 3'	5'-AgA 3'	
	893		571					
	5'-CCT 3'		5'-CCg 3'					
Well No.	73	74	75	76	77	78	79	80

¹The internal positive control primer pairs amplify segments of the human growth hormone gene. The two different control primer pairs give rise to either an internal positive control band of 1070 base pairs, for most wells, or a band of 800 base pairs, for some wells.

Well number 1 contains the primer pair giving rise to the shorter, 800 bp, internal positive control band in order to help in the correct orientation of the HLA-B*07 subtyping.

In addition, wells number 4, 10, 14, 19, 24, 26, 27, 29 to 31, 47, 52, 53, 59, 62, 67, 70 and 77 contain the primer pair giving rise to the shorter, 800 bp, internal positive control band in order to allow kit identification.

In the presence of a specific amplification the intensity of the control band often decreases.

²The nucleotide position, in the 1st, 2nd, 3rd or 4th exon, matching the specificity-determining 3'-end of the primer is given. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

³The nucleotide position, in 2nd, 3rd, 4th or 5th exon or the 2nd intron, matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot-specific information

[illegible]

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot-specific information

CELL LINE VALIDATION SHEET																																	
HLA-B*07 SSP subtyping kit																																	
						Well																											
						17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32												
					Prod. No.:	200967517	200967518	200967519	200967520	201314121	200967522	200967523	200967524	201182125	201297026	201297027	200967528	201182129	200967530	200967531	200967532												
	IHWC cell line		B*																														
1	9001	SA	*07:02			-	-	-	-	-	-	-	-																				
2	9280	LK707	*52:01	*73:01		-	-	-	-	-	-	-	-																				
3	9011	E4181324	*52:01			-	-	-	-	-	-	-	-																				
4	9275	GU373	*15:10	*53:01		-	-	-	-	-	-	-	-																				
5	9009	KAS011	*37:01			-	-	-	-	-	-	-	-																				
6	9353	SM	*39:01	*51:01		-	-	-	-	-	-	-	-																				
7	9020	QBL	*18:01			-	-	-	-	-	-	-	-																				
8	9025	DEU	*35:01			-	-	-	-	-	-	-	-																				
9	9026	YAR	*38:01			-	-	-	-	-	-	-	-																				
10	9107	LKT3	*54:01			-	-	-	-	-	-	-	-																				
11	9051	PITOUT	*44:03			-	-	-	-	-	-	-	-																				
12	9052	DBB	*57:01			-	-	-	-	-	-	-	-																				
13	9025	JESTHOM	*27:05			-	-	-	-	-	-	-	-																				
14	9071	OLGA	*15:01	*15:20		-	-	-	-	-	-	-	-																				
15	9075	DKB	*40:01			-	-	-	-	-	-	-	+							+	-	-											
16	9037	SWEIG007	*40:02			-	-	-	-	-	-	-	-							+	+	-											
17	9282	CTM3953540	*08:01	*55:01		-	-	-	-	-	-	-	-							-	-	+											
18	9257	32367	*14:01	*56:01		-	-	-	-	-	-	-	-							-	-	-											
19	9038	BM16	*18:01			-	-	-	-	-	-	-	-							-	-	-											
20	9059	SLE005	*40:01			-	-	-	-	-	-	-	+							+	-	-											
21	9064	AMALA	*15:01			-	-	-	-	-	-	-	-							-	-	-											
22	9056	KOSE	*35:03			-	-	-	-	-	-	-	-							-	-	-											
23	9124	IHL	*40:02	*56:02		-	-	-	-	-	-	-	-							+	+	-											
24	9035	JBUSH	*38:01			-	-	-	-	-	-	-	-							-	-	+											
25	9049	IBW9	*14:02			-	-	-	-	-	-	-	-							-	-	-											
26	9285	WT49	*58:01			-	-	-	-	-	-	-	-							-	-	-											
27	9191	CH1007	*07:05	*51:01		-	-	-	-	-	-	-	-							+	-	-											
28	9320	BEL5GB	*44:02	*44:03		-	-	-	-	-	-	-	-							-	-	-											
29	9050	MOU	*44:03			-	-	-	-	-	-	-	-							-	-	-											
30	9021	RSH	*42:01			-	-	-	-	-	-	-	-							-	-	+											
31	9019	DUCAF	*18:01			-	-	-	-	-	-	-	-							-	-	-											
32	9297	HAG	*41:02			-	-	-	-	-	-	-	-							-	-	+											
33	9098	MT14B	*40:01			-	-	-	-	-	-	-	+							+	-	-											
34	9104	DHIF	*38:01			-	-	-	-	-	-	-	-							-	-	+											
35	9302	SSTO	*44:02			-	-	-	-	-	-	-	-							-	-	-											
36	9024	KT17	*15:01	*35:01		-	-	-	-	-	-	-	-							-	-	-											
37	9065	HHKB	*07:02			-	-	-	-	-	-	-	-							-	-	-											
38	9099	LZL	*15:01			-	-	-	-	-	-	-	-							-	-	-											
39	9315	CML	*08:01	*27:05		-	-	-	-	-	-	-	-							-	-	+											
40	9134	WHONP199	*13:02	*46:01		-	-	-	-	-	-	-	-							+	-	-											
41	9055	H0301	*14:02			-	-	-	-	-	-	-	-							-	-	-											
42	9066	TAB089	*46:01			-	-	-	-	-	-	-	-							-	-	-											
43	9076	T7526	*46:01			-	-	-	-	-	-	-	-							-	-	-											
44	9057	TEM	*38:01			-	-	-	-	-	-	-	-							-	-	+											
45	9239	SHJO	*42:01	*50:01		-	-	-	-	-	-	-	-							-	-	+											
46	9013	SCHU	*07:02			-	-	-	-	-	-	-	-							-	-	-											
47	9045	TUBO	*51:01			-	-	-	-	-	-	-	-							-	-	-											
48	9303	TER-ND	*35:01	*44:03		-	-	-	-	-	-	-	-							-	-	-											

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot-specific information

CELL LINE VALIDATION SHEET																								
HLA-B*07 SSP subtyping kit																								
					Well																			
					Prod. No.:	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48			
						201314133	200967534	201314135	201314136	200967537	201182138	200967539	200967540	201314141	200967542	200967543	200967544	200967545	200967546	200967547	201297048			
	IHWC cell line		B*																					
1	9001	SA	*07:02			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
2	9280	LK707	*52:01	*73:01		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-		
3	9011	E4181324	*52:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
4	9275	GU373	*15:10	*53:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
5	9009	KAS011	*37:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
6	9353	SM	*39:01	*51:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
7	9020	QBL	*18:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
8	9025	DEU	*35:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
9	9026	YAR	*38:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
10	9107	LKT3	*54:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
11	9051	PITOUT	*44:03			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
12	9052	DBB	*57:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
13	9025	JESTHOM	*27:05			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
14	9071	OLGA	*15:01	*15:20		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
15	9075	DKB	*40:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
16	9037	SWEG007	*40:02			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
17	9282	CTM3953540	*08:01	*55:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	-		
18	9257	32367	*14:01	*56:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
19	9038	BM16	*18:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
20	9059	SLE005	*40:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
21	9064	AMALA	*15:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
22	9056	KOSE	*35:03			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
23	9124	IHL	*40:02	*56:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
24	9035	JBUSH	*38:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
25	9049	IBW9	*14:02			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
26	9285	WT49	*58:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
27	9191	CH1007	*07:05	*51:01		-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-		
28	9320	BEL5GB	*44:02	*44:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
29	9050	MOU	*44:03			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
30	9021	RSH	*42:01			-	-	-	-	-	-	+	-	-	-	-	-	-	-	+	-	-		
31	9019	DUCAF	*18:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
32	9297	HAG	*41:02			-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	-		
33	9098	MT14B	*40:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
34	9104	DHIF	*38:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
35	9302	SSTO	*44:02			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
36	9024	KT17	*15:01	*35:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
37	9065	HHKB	*07:02			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
38	9099	LZL	*15:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
39	9315	CML	*08:01	*27:05		-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	-		
40	9134	WHONP199	*13:02	*46:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
41	9055	H0301	*14:02			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
42	9066	TAB089	*46:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
43	9076	T7526	*46:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
44	9057	TEM	*38:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
45	9239	SHJO	*42:01	*50:01		-	-	-	-	-	-	+	-	-	-	-	-	-	-	+	-	-		
46	9013	SCHU	*07:02			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
47	9045	TUBO	*51:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
48	9303	TER-ND	*35:01	*44:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		

101.512-24/03– including *Taq* pol., IFU-01
 101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **92R**

Lot-specific information

CELL LINE VALIDATION SHEET																				
HLA-B*07 SSP subtyping kit																				
					Well															
					49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64
				Prod. No.:	200967549	201314150	200967551	201182152	200967553	201314154	200967555	201297056	201182157	201297058	201182159	200967560	201182161	200967562	201297063	201297064
	IHWC cell line		B*																	
1	9001	SA	*07:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
2	9280	LK707	*52:01	*73:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
3	9011	E4181324	*52:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
4	9275	GU373	*15:10	*53:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
5	9009	KAS011	*37:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
6	9353	SM	*39:01	*51:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+
7	9020	QBL	*18:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
8	9025	DEU	*35:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
9	9026	YAR	*38:01		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	+
10	9107	LKT3	*54:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
11	9051	PITOUT	*44:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
12	9052	DBB	*57:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
13	9025	JESTHOM	*27:05		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
14	9071	OLGA	*15:01	*15:20	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
15	9075	DKB	*40:01		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	-
16	9037	SWEIG007	*40:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
17	9282	CTM3953540	*08:01	*55:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
18	9257	32367	*14:01	*56:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+
19	9038	BM16	*18:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
20	9059	SLE005	*40:01		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	-
21	9064	AMALA	*15:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
22	9056	KOSE	*35:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
23	9124	IHL	*40:02	*56:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
24	9035	JBUSH	*38:01		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	+
25	9049	IBW9	*14:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
26	9285	WT49	*58:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
27	9191	CH1007	*07:05	*51:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
28	9320	BEL5GB	*44:02	*44:03	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
29	9050	MOU	*44:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
30	9021	RSH	*42:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
31	9019	DUCAF	*18:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
32	9297	HAG	*41:02		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	-
33	9098	MT14B	*40:01		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	-
34	9104	DHIF	*38:01		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	+
35	9302	SSTO	*44:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
36	9024	KT17	*15:01	*35:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
37	9065	HHKB	*07:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
38	9099	LZL	*15:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
39	9315	CML	*08:01	*27:05	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
40	9134	WHONP199	*13:02	*46:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
41	9055	H0301	*14:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
42	9066	TAB089	*46:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
43	9076	T7526	*46:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
44	9057	TEM	*38:01		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	+
45	9239	SHJO	*42:01	*50:01	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	-
46	9013	SCHU	*07:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
47	9045	TUBO	*51:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
48	9303	TER-ND	*35:01	*44:03	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-

101.512-24/03– including *Taq* pol., IFU-01
 101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **92R**

Lot-specific information

CELL LINE VALIDATION SHEET																					
HLA-B*07 SSP subtyping kit																					
					Well																
						65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80
				Prod. No.:		201314165	200967566	201297067	201314168	201297069	201297070	201297071	201297072	201314173	201314174	201314175	201297076	201314177	201314178	201314179	201314180
	IHWC cell line		B*																		
1	9001	SA	*07:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
2	9280	LK707	*52:01	*73:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
3	9011	E4181324	*52:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
4	9275	GU373	*15:10	*53:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
5	9009	KAS011	*37:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
6	9353	SM	*39:01	*51:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
7	9020	QBL	*18:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
8	9025	DEU	*35:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
9	9026	YAR	*38:01		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
10	9107	LKT3	*54:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
11	9051	PITOUT	*44:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
12	9052	DBB	*57:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
13	9025	JESTHOM	*27:05		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
14	9071	OLGA	*15:01	*15:20	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
15	9075	DKB	*40:01		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-
16	9037	SWEIG007	*40:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
17	9282	CTM3953540	*08:01	*55:01	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-
18	9257	32367	*14:01	*56:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
19	9038	BM16	*18:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
20	9059	SLE005	*40:01		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-
21	9064	AMALA	*15:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
22	9056	KOSE	*35:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
23	9124	IHL	*40:02	*56:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
24	9035	JBUSH	*38:01		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
25	9049	IBW9	*14:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
26	9285	WT49	*58:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
27	9191	CH1007	*07:05	*51:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
28	9320	BEL5GB	*44:02	*44:03	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
29	9050	MOU	*44:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
30	9021	RSH	*42:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
31	9019	DUCAF	*18:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
32	9297	HAG	*41:02		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-
33	9098	MT14B	*40:01		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-
34	9104	DHIF	*38:01		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
35	9302	SSTO	*44:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
36	9024	KT17	*15:01	*35:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
37	9065	HHKB	*07:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
38	9099	LZL	*15:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
39	9315	CML	*08:01	*27:05	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-
40	9134	WHONP199	*13:02	*46:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
41	9055	H0301	*14:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
42	9066	TAB089	*46:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
43	9076	T7526	*46:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
44	9057	TEM	*38:01		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
45	9239	SHJO	*42:01	*50:01	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-
46	9013	SCHU	*07:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
47	9045	TUBO	*51:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
48	9303	TER-ND	*35:01	*44:03	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-

101.512-24/03– including *Taq* pol., IFU-01
 101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
 “Instructions for Use” (IFU)

Lot No.: **92R**

Lot-specific information

CERTIFICATE OF ANALYSIS

Olerup SSP® HLA-B*07 SSP

Product number: 101.512-24/03 – including *Taq* pol.
 101.512-24u/03u – without *Taq* pol.
 Lot number: 92R
 Expiry date: 2015-September-01
 Number of tests: 24 tests – Product No. 101.512-24/24u
 3 tests – Product No. 101.512-03/03u
 Number of wells per test: 80

Well specifications:

Well No.	Production No.	Well No.	Production No.	Well No.	Production No.
1	2009-675-01	31	2009-675-31	61	2011-821-61
2	2009-675-02	32	2009-675-32	62	2009-675-62
3	2009-675-03	33	2013-141-33	63	2012-970-63
4	2009-675-04	34	2009-675-34	64	2012-970-64
5	2009-675-05	35	2013-141-35	65	2013-141-65
6	2009-675-06	36	2013-141-36	66	2009-675-66
7	2013-141-07	37	2009-675-37	67	2012-970-67
8	2009-675-08	38	2011-821-38	68	2013-141-68
9	2012-970-09	39	2009-675-39	69	2012-970-69
10	2009-675-10	40	2009-675-40	70	2012-970-70
11	2009-675-11	41	2013-141-41	71	2012-970-71
12	2009-675-12	42	2009-675-42	72	2012-970-72
13	2009-675-13	43	2009-675-43	73	2013-141-73
14	2011-821-14	44	2009-675-44	74	2013-141-74
15	2009-675-15	45	2009-675-45	75	2013-141-75
16	2009-675-16	46	2009-675-46	76	2012-970-76
17	2009-675-17	47	2009-675-47	77	2013-141-77
18	2009-675-18	48	2012-970-48	78	2013-141-78
19	2009-675-19	49	2009-675-49	79	2013-141-79
20	2009-675-20	50	2013-141-50	80	2013-141-80
21	2013-141-21	51	2009-675-51		
22	2009-675-22	52	2011-821-52		
23	2009-675-23	53	2009-675-53		
24	2009-675-24	54	2013-141-54		
25	2011-821-25	55	2009-675-55		
26	2012-970-26	56	2012-970-56		
27	2012-970-27	57	2011-821-57		
28	2009-675-28	58	2012-970-58		
29	2011-821-29	59	2011-821-59		
30	2009-675-30	60	2009-675-60		

The specificity of each primer solution of the kit has been tested against 48 well characterized cell line DNAs.

101.512-24/03– including *Taq* pol., IFU-01
101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: 92R**Lot-specific information**

No DNAs carrying the alleles to be amplified by primer solutions 7, 12, 13, 15, 17, 18, 20 to 23, 25, 27 to 29, 34 to 36, 38, 40, 42, 48 to 55, 57 to 61, 63, 65, 68 and 70 to 80 were available.

The specificities of the primers in primer solutions 7, 12, 13, 15, 17, 18, 20, 21, 23, 25, 27, 29, 36, 38, 40, 48, 50 to 53, 57, 59, 60, 63, 68, 70, 72, 73 and 75 to 79 were tested by separately adding additional 5'-primers, respectively additional 3'-primers. In primer solutions 54, 58, 61, 65, 71 and 80 it was only possible to test the 5'-primer, the 3'-primer were not possible to test. In primer solutions 22, 28, 34, 35, 42, 49, 55 and 74 it was only possible to test the 3'-primers, the 5'-primers were not possible to test. In primer solutions 21, 27, 29, 41, 48, 57, 63, 70, 73, 76, 78 and 79 one or two 5'-primers were not possible to test, and in primer solutions 8, 10, 11, 15, 17, 19, 20, 24 to 26, 33, 36, 38, 40, 50, 51, 56, 59, 63, 68, 69, 73 and 75 to 79 one or two 3'-primers were not possible to test. Additional primers in primer solutions 8, 10, 11, 19, 24, 26, 31, 33, 41 and 69 were tested by separately adding one additional 5'-primer, or one additional 3'-primer.

Results: No false positive or false negative amplifications were obtained.

Date of approval: 2013-April-12

Approved by:

Production Quality Control

101.512-24/03– including *Taq* pol., IFU-01
101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: **92R**

Lot-specific information

Declaration of Conformity

Product name: *Olerup* SSP® HLA-B*07
Product number: 101.512-24/24u, 03/03u
Lot number: 92R

Intended use: HLA-B*07 high resolution histocompatibility testing

Manufacturer: *Olerup* SSP AB
Franzengatan 5
SE-112 51 Stockholm, Sweden
Phone: +46-8-717 88 27
Fax: +46-8-717 88 18

We, *Olerup* SSP AB, hereby declare that this product, to which this Declaration of Conformity relates is in conformity with the following Standard(s) and other normative document(s) ISO 9001:2008 and ISO 13485:2012, following the provisions of the 98/79/EC Directive on *in vitro* diagnostic medical devices, Annex II List B, conformity assessed using Annex IV, as transposed into the national laws of the Member States of the European Union.

The Technical Documentation File is maintained at *Olerup* SSP AB, Franzengatan 5, SE-112 51 Stockholm, Sweden.

Notified Body: Lloyd's Register Quality Assurance Limited, Hiramford, Middlemarch Office Village, Siskin Drive, Coventry CV3 4FJ, United Kingdom.
(Notified Body number: 0088.)

Stockholm, Sweden
2013-April-12

Ann-Cathrin Jareman
Head of QA and Regulatory Affairs

101.512-24/03– including *Taq* pol., IFU-01
101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: **92R**

Lot-specific information

101.512-24/03– including *Taq* pol., IFU-01
101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: **92R**

Lot-specific information

101.512-24/03– including *Taq* pol., IFU-01
101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: **92R**

Lot-specific information

101.512-24/03– including *Taq* pol., IFU-01
101. 512-24u/03u – without *Taq* pol., IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: **92R**

Lot-specific information

ADDRESSES:

Manufacturer:

Olerup SSP AB, Franzengatan 5, SE-112 51 Stockholm, Sweden.

Tel: +46-8-717 88 27

Fax: +46-8-717 88 18

E-mail: info-ssp@olerup.com

Web page: <http://www.olerup-ssp.com>

Distributed by:

Olerup GmbH, Löwengasse 47 / 6, AT-1030 Vienna, Austria.

Tel: +43-1-710 15 00

Fax: +43-1-710 15 00 10

E-mail: support-at@olerup.com

Web page: <http://www.olerup.com>

Olerup Inc., 901 S. Bolmar St., Suite R, West Chester, PA 19382

Tel: 1-877-OLERUP1

Fax: 610-344-7989

E-mail: info.us@olerup.com

Web page: <http://www.olerup.com>

For information on *Olerup* SSP distributors worldwide, contact **Olerup GmbH**.