

101.615-24/04 – including *Taq* polymerase
101.615-24u/04u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **7N8**

Lot-specific Information

Olerup SSP® HLA-C*07

Product number:	101.615-24/04 – including <i>Taq</i> pol. 101.615-24u/04u – without <i>Taq</i> pol.
Lot number:	7N8
Expiry date:	2027-01-01
Number of tests:	24 tests – Product No. 101.615-24/24u 4 tests – Product No. 101.615-04/04u
Number of wells per test:	95+1
Storage - pre-aliquoted primers:	dark, between -15°C and -25°C
- PCR Master Mix:	between -15°C and -25°C
- Adhesive PCR seals	RT

This Product Description is only valid for Lot No. 7N8.

Complete product documentation consists of generic Instructions for Use (IFU), lot specific Product Insert, Worksheet and Certificate.

CHANGES COMPARED TO THE PREVIOUS OLERUP SSP® HLA-C*07 Lot (4N5)

- The product documentation has been updated for new alleles of IMGT 3.50.0.
- The kit resolution focuses on common and well documented (CWD) alleles¹.

The HLA-C*07 primer set, specificity and interpretation tables have been updated for the HLA-C alleles described since the previous *Olerup SSP®* HLA-C*07 lot was made (**Lot No. 4N5**).

The primers of the wells detailed below have been exchanged, added or modified compared to the previous lot.

Well	5'-primer	3'-primer	rationale
95	-	Modified	3'-primer modified for improved specificity of HLA-specific primer pair.

¹S. J. Mack, P. Cano, J. A. Hollenbach et al.
Common and well-documented HLA alleles: 2012 update to
the CWD catalogue. Tissue Antigens, 2013, 81, 194–203



101.615-24/04 – including *Taq* polymerase
101.615-24u/04u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **7N8**

Lot-specific Information

Well **96** contains Negative Control primer pairs, that will amplify the majority of the Olerup SSP® HLA Class I, DRB, DQB1, DPB1 and DQA1 amplicons as well as all the amplicons generated by the control primer pairs matching the human growth hormone gene.

HLA-specific PCR product sizes range from 75 to 200 base pairs.
The PCR product generated by the positive control primer pair is 200 base pairs.

Length of PCR product	105	200	105	80	75	80	85
5'-primer¹	164	340	440	45	45	43	36
	5'-CAC ^{3'}	5'-Agg ^{3'}	5'-TTA ^{3'}	5'-Tgg ^{3'}	5'-Tgg ^{3'}	5'-Tgg ^{3'}	5'-TAC ^{3'}
							36
							5'-TAT ^{3'}
3'-primer²	231	2nd I	507	59	58	57	47
	5'-TgC ^{3'}	5'-AAA ^{3'}	5'-TTg ^{3'}	5'-CTC ^{3'}	5'-ggC ^{3'}	5'-CTC ^{3'}	5'-ACA ^{3'}
							48
							5'-gCA ^{3'}
							48
							5'-gCC ^{3'}
							52
							5'-TgT ^{3'}
A*	+	+	+				
B*	+	+	+				
C*	+	+	+				
DRB1				+	+		
DRB3				+	+		
DRB5				+			
DQB1					+		
DPB1						+	
DQA1							+

¹The nucleotide position for HLA class I genes and the codon for HLA class II genes, in the 2nd or 3rd exon, matching the specificity-determining 3'-end of the primer is given. Nucleotide and codon numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

²The nucleotide position for HLA class I genes and the codon for HLA class II genes, in the 2nd or 3rd exon or the 2nd intron, matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Nucleotide and codon numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.



101.615-24/04 – including *Taq* polymerase
101.615-24u/04u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **7N8**

Lot-specific Information

PRODUCT DESCRIPTION

HLA-C*07 SSP subtyping

CONTENT

The primer set contains 5'- and 3'-primers for identifying the HLA-C*07:01 to HLA-C*07:1042N alleles.

PLATE LAYOUT

Each test consists of 96 PCR reactions in a 96 well cut PCR plate.

1	2	3	4	5	6	7	8
9	10	11	12	13	14	15	16
17	18	19	20	21	22	23	24
25	26	27	28	29	30	31	32
33	34	35	36	37	38	39	40
41	42	43	44	45	46	47	48
49	50	51	52	53	54	55	56
57	58	59	60	61	62	63	64
65	66	67	68	69	70	71	72
73	74	75	76	77	78	79	80
81	82	83	84	85	86	87	88
89	90	91	92	93	94	95	NC

The 96 well cut PCR plate is marked with 'HLA-C*07' in silver/gray ink.

Well No. 1 is marked with the Lot No. '7N8'.

Wells 1 to 95 – HLA-C*07 high resolution primers.

Well 96 – Negative Control (NC).

A faint row of numbers is seen between wells 1 and 2 or wells 7 and 8 of the PCR trays. These stem from the manufacture of the trays and should be disregarded. The PCR plates are heat-sealed with a PCR-compatible foil.

INTERPRETATION

Due to the sharing of sequence motifs between HLA-C alleles, non-HLA-C*07 alleles will be amplified by some primer mixes. For further details see Specificity Table.

UNIQUELY IDENTIFIED ALLELES

All the HLA-C*07 alleles i.e., **C*07:01 to C*07:1042N**, recognized by the HLA Nomenclature Committee in October 2022^{1,2} will be amplified by the primers in the HLA-C*07 SSP kit³.



101.615-24/04 – including *Taq* polymerase
101.615-24u/04u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **7N8**

Lot-specific Information

The HLA-C*07 kit enables separation of the confirmed HLA-C*07 alleles as listed in the IMGT/HLA database 3.26.0. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources. Current allele confirmation status for HLA-C*07 alleles is listed below.

The HLA-C*07 kit also enables identification of many null and alternatively expressed alleles.

The following HLA-C*07 alleles can be distinguished by the different sizes of the HLA-specific PCR product:

Alleles	Primer mix	Alleles	Primer mix
C*07:29:01:01-07:29:02, 07:347N, 07:1003	71	C*07:120, 07:462, 07:878	61
C*07:41, 07:447	1	C*07:128, 07:279	66
C*07:52, 07:337	18	C*07:132, 07:179	55
C*07:58, 07:713Q	86	C*07:138, 07:484N	72
C*07:74, 07:308	39	C*07:147, 07:241	60
C*07:78:01-07:78:02, 07:235Q, 07:559	40	C*07:149, 07:381	68
C*07:83, 07:649	85	C*07:164N, 07:212	58
C*07:88, 07:474	41	C*07:165, 07:180	60
C*07:93:02, 07:207	42	C*07:186, 07:451N	58
C*07:99, 07:183	43	C*07:226Q, 07:573	86
C*07:100, 07:161	49	C*07:227N, 07:254	30
C*07:104N, 07:119, 07:893	92	C*07:245, 07:342:01-07:342:02	45
C*07:117, 07:126	51	C*07:252, 07:679	55
C*07:118, 07:203	52	C*07:457, 07:702N, 07:784	92

¹HLA-C alleles listed on the IMGT/HLA web page 2022-October-12, release 3.50.0, www.ebi.ac.uk/imgt/hla.

²Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

³The following alleles will give rise to identical amplification patterns. These alleles can e.g. be distinguished by the HLA-C low resolution kit.

Alleles

C*07:292, 07:922, 07:937, 07:1031, C*12:181



101.615-24/04 – including *Taq* polymerase
101.615-24u/04u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **7N8**

Lot-specific Information

ALLELE CONFIRMATION STATUS

Allele	Status ¹	Allele	Status ¹	Allele	Status ¹	Allele	Status ¹
C*07:01:01:01	Confirmed	C*07:01:45	Unconfirmed	C*07:02:29	Unconfirmed	C*07:04:01:01	Confirmed
C*07:01:01:02	Unconfirmed	C*07:01:46	Unconfirmed	C*07:02:30	Confirmed	C*07:04:01:02	Confirmed
C*07:01:01:03	Unconfirmed	C*07:01:47	Unconfirmed	C*07:02:31	Unconfirmed	C*07:04:02	Confirmed
C*07:01:01:04	Unconfirmed	C*07:01:48	Unconfirmed	C*07:02:32	Confirmed	C*07:04:03	Unconfirmed
C*07:01:01:05	Unconfirmed	C*07:01:49	Unconfirmed	C*07:02:33	Unconfirmed	C*07:04:04	Unconfirmed
C*07:01:01:06	Unconfirmed	C*07:01:50	Unconfirmed	C*07:02:34	Unconfirmed	C*07:04:05	Unconfirmed
C*07:01:01:07	Unconfirmed	C*07:01:51	Confirmed	C*07:02:35	Unconfirmed	C*07:04:06	Confirmed
C*07:01:02	Confirmed	C*07:01:52	Unconfirmed	C*07:02:36	Unconfirmed	C*07:04:07	Unconfirmed
C*07:01:03	Unconfirmed	C*07:01:53	Unconfirmed	C*07:02:37	Unconfirmed	C*07:04:08	Unconfirmed
C*07:01:04	Unconfirmed	C*07:01:54	Unconfirmed	C*07:02:38	Unconfirmed	C*07:04:09	Unconfirmed
C*07:01:05	Confirmed	C*07:01:55	Unconfirmed	C*07:02:39	Confirmed	C*07:04:10	Unconfirmed
C*07:01:06	Unconfirmed	C*07:01:56	Unconfirmed	C*07:02:40	Unconfirmed	C*07:04:11	Unconfirmed
C*07:01:07	Confirmed	C*07:01:57	Unconfirmed	C*07:02:41	Unconfirmed	C*07:05	Confirmed
C*07:01:08	Confirmed	C*07:01:58	Unconfirmed	C*07:02:42	Unconfirmed	C*07:06	Confirmed
C*07:01:09	Confirmed	C*07:02:01:01	Confirmed	C*07:02:43	Unconfirmed	C*07:07	Confirmed
C*07:01:10	Confirmed	C*07:02:01:02	Unconfirmed	C*07:02:44	Confirmed	C*07:08	Unconfirmed
C*07:01:11	Unconfirmed	C*07:02:01:03	Confirmed	C*07:02:45	Unconfirmed	C*07:09	Unconfirmed
C*07:01:12	Unconfirmed	C*07:02:01:04	Unconfirmed	C*07:02:46	Confirmed	C*07:10	Confirmed
C*07:01:13	Unconfirmed	C*07:02:01:05	Unconfirmed	C*07:02:47	Unconfirmed	C*07:11	Unconfirmed
C*07:01:14	Unconfirmed	C*07:02:01:06	Unconfirmed	C*07:02:48	Unconfirmed	C*07:12	Confirmed
C*07:01:15	Confirmed	C*07:02:01:07	Confirmed	C*07:02:49	Confirmed	C*07:13	Confirmed
C*07:01:16	Unconfirmed	C*07:02:01:08	Unconfirmed	C*07:02:50	Unconfirmed	C*07:14	Unconfirmed
C*07:01:17	Unconfirmed	C*07:02:01:09	Unconfirmed	C*07:02:51	Unconfirmed	C*07:15	Confirmed
C*07:01:18	Confirmed	C*07:02:02	Unconfirmed	C*07:02:52	Confirmed	C*07:16	Confirmed
C*07:01:19	Unconfirmed	C*07:02:03	Unconfirmed	C*07:02:53	Unconfirmed	C*07:17:01	Confirmed
C*07:01:20	Unconfirmed	C*07:02:04	Confirmed	C*07:02:54	Unconfirmed	C*07:17:02	Unconfirmed
C*07:01:21	Unconfirmed	C*07:02:05	Unconfirmed	C*07:02:55	Unconfirmed	C*07:17:03	Unconfirmed
C*07:01:22	Unconfirmed	C*07:02:06	Confirmed	C*07:02:56	Unconfirmed	C*07:18	Confirmed
C*07:01:23	Confirmed	C*07:02:07	Confirmed	C*07:02:57	Unconfirmed	C*07:19	Confirmed
C*07:01:24	Unconfirmed	C*07:02:08	Confirmed	C*07:02:58	Unconfirmed	C*07:20	Confirmed
C*07:01:25	Unconfirmed	C*07:02:09	Unconfirmed	C*07:02:59	Unconfirmed	C*07:21	Unconfirmed
C*07:01:26	Unconfirmed	C*07:02:10	Confirmed	C*07:02:60	Unconfirmed	C*07:22	Confirmed
C*07:01:27	Unconfirmed	C*07:02:11	Confirmed	C*07:02:61	Unconfirmed	C*07:23	Unconfirmed
C*07:01:28	Unconfirmed	C*07:02:12	Unconfirmed	C*07:02:62	Confirmed	C*07:24	Unconfirmed
C*07:01:29	Confirmed	C*07:02:13	Confirmed	C*07:02:63	Confirmed	C*07:25	Unconfirmed
C*07:01:30	Confirmed	C*07:02:14	Unconfirmed	C*07:02:64	Unconfirmed	C*07:26:01	Confirmed
C*07:01:31	Confirmed	C*07:02:15	Unconfirmed	C*07:02:65	Confirmed	C*07:26:02	Unconfirmed
C*07:01:32	Unconfirmed	C*07:02:16	Confirmed	C*07:02:66	Unconfirmed	C*07:27:01	Unconfirmed
C*07:01:33	Unconfirmed	C*07:02:17	Confirmed	C*07:02:67	Unconfirmed	C*07:27:02	Confirmed
C*07:01:34	Unconfirmed	C*07:02:18	Confirmed	C*07:02:68	Unconfirmed	C*07:28	Confirmed
C*07:01:35	Unconfirmed	C*07:02:19	Unconfirmed	C*07:02:69	Unconfirmed	C*07:29:01	Unconfirmed
C*07:01:36	Confirmed	C*07:02:20	Confirmed	C*07:02:70	Unconfirmed	C*07:29:02	Unconfirmed
C*07:01:37	Unconfirmed	C*07:02:21	Unconfirmed	C*07:02:71	Unconfirmed	C*07:30	Unconfirmed
C*07:01:38	Unconfirmed	C*07:02:22	Unconfirmed	C*07:02:72	Unconfirmed	C*07:31:01	Confirmed
C*07:01:39	Unconfirmed	C*07:02:23	Unconfirmed	C*07:02:73	Unconfirmed	C*07:31:02	Unconfirmed
C*07:01:40	Unconfirmed	C*07:02:24	Confirmed	C*07:02:74	Unconfirmed	C*07:32N	Confirmed
C*07:01:41	Unconfirmed	C*07:02:25	Unconfirmed	C*07:02:75	Unconfirmed	C*07:33N	Confirmed
C*07:01:42	Confirmed	C*07:02:26	Unconfirmed	C*07:02:76	Confirmed	C*07:35	Confirmed
C*07:01:43	Unconfirmed	C*07:02:27	Unconfirmed	C*07:02:77	Unconfirmed	C*07:36	Confirmed
C*07:01:44	Confirmed	C*07:02:28	Confirmed	C*07:03	Unconfirmed	C*07:37	Unconfirmed



101.615-24/04 – including *Taq* polymerase
101.615-24u/04u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Allele	Status ¹	Allele	Status ¹	Allele	Status ¹	Allele	Status ¹
C*07:38:01	Unconfirmed	C*07:83	Confirmed	C*07:131:01	Confirmed	C*07:177	Unconfirmed
C*07:38:02	Unconfirmed	C*07:84	Unconfirmed	C*07:131:02	Unconfirmed	C*07:178	Unconfirmed
C*07:39	Confirmed	C*07:85	Confirmed	C*07:132	Confirmed	C*07:179	Confirmed
C*07:40	Confirmed	C*07:86	Unconfirmed	C*07:133	Unconfirmed	C*07:180	Confirmed
C*07:41	Unconfirmed	C*07:87	Unconfirmed	C*07:134	Confirmed	C*07:181	Unconfirmed
C*07:42	Unconfirmed	C*07:88	Confirmed	C*07:135	Unconfirmed	C*07:182	Unconfirmed
C*07:43:01	Confirmed	C*07:89	Confirmed	C*07:136	Confirmed	C*07:183	Unconfirmed
C*07:43:02	Unconfirmed	C*07:90	Confirmed	C*07:137:01	Unconfirmed	C*07:184	Confirmed
C*07:44	Unconfirmed	C*07:91	Confirmed	C*07:137:02	Confirmed	C*07:185	Unconfirmed
C*07:45	Confirmed	C*07:92	Unconfirmed	C*07:138	Unconfirmed	C*07:186	Confirmed
C*07:46	Confirmed	C*07:93	Confirmed	C*07:139	Unconfirmed	C*07:187	Unconfirmed
C*07:47	Unconfirmed	C*07:94	Unconfirmed	C*07:140	Confirmed	C*07:188	Unconfirmed
C*07:48	Confirmed	C*07:95	Confirmed	C*07:141:01	Confirmed	C*07:189	Confirmed
C*07:49	Confirmed	C*07:96:01	Confirmed	C*07:141:02	Confirmed	C*07:190	Unconfirmed
C*07:50	Unconfirmed	C*07:96:02	Unconfirmed	C*07:142	Unconfirmed	C*07:191N	Confirmed
C*07:51	Confirmed	C*07:97	Confirmed	C*07:143	Confirmed	C*07:192	Confirmed
C*07:52	Unconfirmed	C*07:98N	Unconfirmed	C*07:144	Unconfirmed	C*07:193	Unconfirmed
C*07:53	Unconfirmed	C*07:99	Unconfirmed	C*07:145	Confirmed	C*07:194	Confirmed
C*07:54	Confirmed	C*07:100	Confirmed	C*07:146	Unconfirmed	C*07:195	Unconfirmed
C*07:55N	Confirmed	C*07:101	Confirmed	C*07:147	Confirmed	C*07:196	Confirmed
C*07:56:01	Unconfirmed	C*07:102	Confirmed	C*07:148	Confirmed	C*07:197	Unconfirmed
C*07:56:02	Confirmed	C*07:103	Unconfirmed	C*07:149	Confirmed	C*07:198N	Confirmed
C*07:57	Unconfirmed	C*07:104N	Confirmed	C*07:150Q	Unconfirmed	C*07:199:01	Unconfirmed
C*07:58	Unconfirmed	C*07:105	Unconfirmed	C*07:151	Confirmed	C*07:199:02	Confirmed
C*07:59	Confirmed	C*07:106	Confirmed	C*07:152N	Unconfirmed	C*07:200	Unconfirmed
C*07:60	Confirmed	C*07:107	Unconfirmed	C*07:153	Confirmed	C*07:201	Unconfirmed
C*07:61N	Confirmed	C*07:108:01	Confirmed	C*07:154	Confirmed	C*07:202	Confirmed
C*07:62	Unconfirmed	C*07:108:02	Confirmed	C*07:155	Confirmed	C*07:203	Confirmed
C*07:63	Unconfirmed	C*07:109	Confirmed	C*07:156	Unconfirmed	C*07:204:01	Unconfirmed
C*07:64	Unconfirmed	C*07:110	Unconfirmed	C*07:157	Unconfirmed	C*07:204:02	Unconfirmed
C*07:65	Confirmed	C*07:111	Unconfirmed	C*07:158	Unconfirmed	C*07:205	Confirmed
C*07:66	Confirmed	C*07:112	Confirmed	C*07:159	Unconfirmed	C*07:206	Unconfirmed
C*07:67	Confirmed	C*07:113	Unconfirmed	C*07:160	Unconfirmed	C*07:207	Confirmed
C*07:68	Confirmed	C*07:114	Unconfirmed	C*07:161	Unconfirmed	C*07:208	Unconfirmed
C*07:69	Confirmed	C*07:115	Unconfirmed	C*07:162	Unconfirmed	C*07:209	Unconfirmed
C*07:70	Confirmed	C*07:116	Confirmed	C*07:163	Unconfirmed	C*07:210	Unconfirmed
C*07:71	Unconfirmed	C*07:117	Confirmed	C*07:164N	Unconfirmed	C*07:211	Unconfirmed
C*07:72	Confirmed	C*07:118	Confirmed	C*07:165	Confirmed	C*07:212	Confirmed
C*07:73	Confirmed	C*07:119	Confirmed	C*07:166	Unconfirmed	C*07:213	Unconfirmed
C*07:74	Unconfirmed	C*07:120	Confirmed	C*07:167	Unconfirmed	C*07:214	Unconfirmed
C*07:75	Confirmed	C*07:121Q	Confirmed	C*07:168	Unconfirmed	C*07:215	Unconfirmed
C*07:76:01	Confirmed	C*07:122	Confirmed	C*07:169	Unconfirmed	C*07:216	Confirmed
C*07:76:02	Confirmed	C*07:123	Confirmed	C*07:170	Unconfirmed	C*07:217	Unconfirmed
C*07:77	Unconfirmed	C*07:124	Confirmed	C*07:171	Confirmed	C*07:218	Confirmed
C*07:78:01	Unconfirmed	C*07:125	Confirmed	C*07:172:01	Confirmed	C*07:219	Unconfirmed
C*07:78:02	Confirmed	C*07:126	Confirmed	C*07:172:02	Confirmed	C*07:220	Unconfirmed
C*07:79	Confirmed	C*07:127	Unconfirmed	C*07:173	Confirmed	C*07:221	Unconfirmed
C*07:80	Confirmed	C*07:128	Unconfirmed	C*07:174	Unconfirmed	C*07:222	Confirmed
C*07:81	Unconfirmed	C*07:129	Confirmed	C*07:175	Unconfirmed	C*07:223	Unconfirmed
C*07:82	Unconfirmed	C*07:130	Confirmed	C*07:176	Confirmed	C*07:224	Confirmed



101.615-24/04 – including *Taq* polymerase
101.615-24u/04u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **7N8**

Lot-specific Information

Allele	Status ¹	Allele	Status ¹	Allele	Status ¹	Allele	Status ¹
C*07:225	Unconfirmed	C*07:270	Unconfirmed	C*07:319	Unconfirmed	C*07:367	Unconfirmed
C*07:226	Confirmed	C*07:271	Unconfirmed	C*07:320	Unconfirmed	C*07:368:01	Unconfirmed
C*07:227N	Unconfirmed	C*07:272	Unconfirmed	C*07:321	Unconfirmed	C*07:368:02	Unconfirmed
C*07:228	Unconfirmed	C*07:273	Unconfirmed	C*07:322	Confirmed	C*07:369	Unconfirmed
C*07:229	Confirmed	C*07:274	Unconfirmed	C*07:323	Unconfirmed	C*07:370	Unconfirmed
C*07:230	Unconfirmed	C*07:275	Unconfirmed	C*07:324	Unconfirmed	C*07:371	Confirmed
C*07:231	Confirmed	C*07:276	Confirmed	C*07:325	Unconfirmed	C*07:372	Unconfirmed
C*07:232	Confirmed	C*07:277	Unconfirmed	C*07:326	Confirmed	C*07:373	Unconfirmed
C*07:233	Unconfirmed	C*07:278	Unconfirmed	C*07:327	Unconfirmed	C*07:374	Confirmed
C*07:234	Confirmed	C*07:279	Confirmed	C*07:328	Unconfirmed	C*07:375	Unconfirmed
C*07:235Q	Unconfirmed	C*07:280	Unconfirmed	C*07:329N	Unconfirmed	C*07:376	Unconfirmed
C*07:236	Unconfirmed	C*07:281	Unconfirmed	C*07:330:01	Unconfirmed	C*07:377	Unconfirmed
C*07:237	Unconfirmed	C*07:282	Unconfirmed	C*07:330:02	Unconfirmed	C*07:378	Unconfirmed
C*07:238	Confirmed	C*07:283	Unconfirmed	C*07:331	Unconfirmed	C*07:379	Unconfirmed
C*07:239	Unconfirmed	C*07:284	Unconfirmed	C*07:332	Unconfirmed	C*07:380	Unconfirmed
C*07:240	Unconfirmed	C*07:285	Confirmed	C*07:333	Unconfirmed	C*07:381	Confirmed
C*07:241	Confirmed	C*07:286	Unconfirmed	C*07:334	Unconfirmed	C*07:382	Unconfirmed
C*07:242	Confirmed	C*07:287	Unconfirmed	C*07:335	Unconfirmed	C*07:383	Unconfirmed
C*07:243	Confirmed	C*07:288	Confirmed	C*07:336	Unconfirmed	C*07:384	Unconfirmed
C*07:244	Unconfirmed	C*07:289	Unconfirmed	C*07:337	Unconfirmed	C*07:385	Unconfirmed
C*07:245	Unconfirmed	C*07:290	Unconfirmed	C*07:338	Unconfirmed	C*07:386	Confirmed
C*07:246:01	Unconfirmed	C*07:291	Unconfirmed	C*07:339	Confirmed	C*07:387	Unconfirmed
C*07:246:02	Unconfirmed	C*07:292	Unconfirmed	C*07:340	Confirmed	C*07:388	Unconfirmed
C*07:247	Unconfirmed	C*07:293	Confirmed	C*07:341:01	Unconfirmed	C*07:389	Unconfirmed
C*07:248	Confirmed	C*07:294	Unconfirmed	C*07:341:02	Unconfirmed	C*07:390	Unconfirmed
C*07:249	Unconfirmed	C*07:296	Unconfirmed	C*07:342	Confirmed	C*07:391	Unconfirmed
C*07:250	Confirmed	C*07:297	Unconfirmed	C*07:343	Confirmed	C*07:392	Unconfirmed
C*07:251	Unconfirmed	C*07:298	Unconfirmed	C*07:344	Unconfirmed	C*07:393N	Unconfirmed
C*07:252	Confirmed	C*07:299	Unconfirmed	C*07:345	Unconfirmed	C*07:394	Unconfirmed
C*07:253	Unconfirmed	C*07:300	Unconfirmed	C*07:346	Unconfirmed	C*07:395	Unconfirmed
C*07:254	Confirmed	C*07:301	Confirmed	C*07:347N	Unconfirmed	C*07:396	Confirmed
C*07:255	Unconfirmed	C*07:302	Unconfirmed	C*07:348	Unconfirmed	C*07:397	Confirmed
C*07:256	Confirmed	C*07:303	Unconfirmed	C*07:349	Unconfirmed	C*07:398	Unconfirmed
C*07:257:01	Unconfirmed	C*07:304	Unconfirmed	C*07:350N	Unconfirmed	C*07:399	Unconfirmed
C*07:257:02	Unconfirmed	C*07:305	Unconfirmed	C*07:351	Confirmed	C*07:400	Unconfirmed
C*07:257:03	Unconfirmed	C*07:306	Unconfirmed	C*07:352	Confirmed	C*07:401	Unconfirmed
C*07:257:04	Unconfirmed	C*07:307	Confirmed	C*07:353	Unconfirmed	C*07:402	Unconfirmed
C*07:258	Unconfirmed	C*07:308	Unconfirmed	C*07:354	Unconfirmed	C*07:403	Unconfirmed
C*07:259	Confirmed	C*07:309	Unconfirmed	C*07:355	Confirmed	C*07:404	Unconfirmed
C*07:260:01	Unconfirmed	C*07:310	Confirmed	C*07:356	Unconfirmed	C*07:405	Unconfirmed
C*07:260:02	Unconfirmed	C*07:311	Confirmed	C*07:357	Unconfirmed	C*07:406	Unconfirmed
C*07:261	Unconfirmed	C*07:312	Unconfirmed	C*07:358	Unconfirmed	C*07:407	Confirmed
C*07:262	Unconfirmed	C*07:313	Unconfirmed	C*07:359	Confirmed	C*07:408	Unconfirmed
C*07:263	Unconfirmed	C*07:314:01	Confirmed	C*07:360	Confirmed	C*07:409	Confirmed
C*07:264N	Confirmed	C*07:314:02	Unconfirmed	C*07:361	Confirmed	C*07:410	Unconfirmed
C*07:265	Unconfirmed	C*07:314:03	Confirmed	C*07:362	Unconfirmed	C*07:411	Confirmed
C*07:266	Unconfirmed	C*07:315	Unconfirmed	C*07:363	Unconfirmed	C*07:412	Confirmed
C*07:267	Unconfirmed	C*07:316	Unconfirmed	C*07:364	Unconfirmed	C*07:413	Unconfirmed
C*07:268	Unconfirmed	C*07:317	Unconfirmed	C*07:365	Unconfirmed	C*07:414	Unconfirmed
C*07:269	Unconfirmed	C*07:318	Unconfirmed	C*07:366	Unconfirmed	C*07:415	Unconfirmed



101.615-24/04 – including *Taq* polymerase
101.615-24u/04u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **7N8**

Lot-specific Information

Allele	Status ¹	Allele	Status ¹	Allele	Status ¹	Allele	Status ¹
C*07:416	Unconfirmed	C*07:455	Unconfirmed	C*07:494Q	Unconfirmed	C*07:534	Unconfirmed
C*07:417	Unconfirmed	C*07:456	Confirmed	C*07:495	Unconfirmed	C*07:535	Unconfirmed
C*07:418	Unconfirmed	C*07:457	Confirmed	C*07:496	Unconfirmed	C*07:536	Unconfirmed
C*07:419	Unconfirmed	C*07:458	Unconfirmed	C*07:497	Unconfirmed	C*07:537	Unconfirmed
C*07:420	Unconfirmed	C*07:459	Unconfirmed	C*07:498	Unconfirmed	C*07:538	Unconfirmed
C*07:421	Unconfirmed	C*07:460	Unconfirmed	C*07:499	Unconfirmed	C*07:539	Unconfirmed
C*07:422	Unconfirmed	C*07:461	Unconfirmed	C*07:500	Unconfirmed	C*07:540	Unconfirmed
C*07:423	Unconfirmed	C*07:462	Confirmed	C*07:501	Unconfirmed	C*07:541	Unconfirmed
C*07:424	Unconfirmed	C*07:463	Unconfirmed	C*07:502	Unconfirmed	C*07:542	Unconfirmed
C*07:425	Confirmed	C*07:464	Unconfirmed	C*07:503	Unconfirmed	C*07:544	Unconfirmed
C*07:426	Unconfirmed	C*07:465	Unconfirmed	C*07:504	Unconfirmed	C*07:545	Unconfirmed
C*07:427	Unconfirmed	C*07:466	Unconfirmed	C*07:505	Unconfirmed	C*07:546	Unconfirmed
C*07:428	Unconfirmed	C*07:467	Unconfirmed	C*07:506	Unconfirmed	C*07:547	Unconfirmed
C*07:429	Confirmed	C*07:468	Unconfirmed	C*07:507N	Unconfirmed	C*07:548	Unconfirmed
C*07:430	Unconfirmed	C*07:469	Unconfirmed	C*07:508	Unconfirmed	C*07:549	Unconfirmed
C*07:431	Unconfirmed	C*07:470	Unconfirmed	C*07:509	Unconfirmed	C*07:550	Unconfirmed
C*07:432	Unconfirmed	C*07:471	Unconfirmed	C*07:510	Unconfirmed	C*07:551N	Unconfirmed
C*07:433	Unconfirmed	C*07:472	Unconfirmed	C*07:511	Unconfirmed	C*07:552	Unconfirmed
C*07:434	Unconfirmed	C*07:473	Unconfirmed	C*07:512	Confirmed	C*07:553	Unconfirmed
C*07:435	Unconfirmed	C*07:474	Unconfirmed	C*07:513Q	Unconfirmed	C*07:554	Unconfirmed
C*07:436	Unconfirmed	C*07:475	Unconfirmed	C*07:514	Unconfirmed	C*07:555	Unconfirmed
C*07:437N	Unconfirmed	C*07:476N	Confirmed	C*07:515	Unconfirmed	C*07:556	Unconfirmed
C*07:438	Unconfirmed	C*07:477	Unconfirmed	C*07:516	Unconfirmed	C*07:557	Unconfirmed
C*07:439	Confirmed	C*07:478	Unconfirmed	C*07:517	Unconfirmed	C*07:558	Unconfirmed
C*07:440	Unconfirmed	C*07:479	Unconfirmed	C*07:518	Unconfirmed	C*07:559	Unconfirmed
C*07:441	Unconfirmed	C*07:480	Unconfirmed	C*07:519	Unconfirmed	C*07:560	Unconfirmed
C*07:442	Unconfirmed	C*07:481	Unconfirmed	C*07:520	Unconfirmed	C*07:561	Unconfirmed
C*07:443	Unconfirmed	C*07:482	Unconfirmed	C*07:521	Unconfirmed	C*07:562	Unconfirmed
C*07:444	Unconfirmed	C*07:483N	Confirmed	C*07:522	Unconfirmed	C*07:563	Unconfirmed
C*07:445	Confirmed	C*07:484N	Unconfirmed	C*07:523	Unconfirmed		
C*07:446	Unconfirmed	C*07:485	Confirmed	C*07:524	Confirmed		
C*07:447	Unconfirmed	C*07:486	Unconfirmed	C*07:525	Unconfirmed		
C*07:448	Unconfirmed	C*07:487	Unconfirmed	C*07:526	Unconfirmed		
C*07:449	Unconfirmed	C*07:488	Unconfirmed	C*07:527	Unconfirmed		
C*07:450	Confirmed	C*07:489	Unconfirmed	C*07:528	Unconfirmed		
C*07:451N	Confirmed	C*07:490	Unconfirmed	C*07:529	Unconfirmed		
C*07:452N	Unconfirmed	C*07:491:01N	Unconfirmed	C*07:530	Unconfirmed		
C*07:453	Unconfirmed	C*07:491:02N	Unconfirmed	C*07:531	Unconfirmed		
C*07:454:01	Unconfirmed	C*07:492	Unconfirmed	C*07:532	Unconfirmed		
C*07:454:02	Unconfirmed	C*07:493	Unconfirmed	C*07:533	Unconfirmed		

¹Allele status “confirmed” or “unconfirmed” as listed on the IMGT/HLA web page 2016-October-14, release 3.26.0, www.ebi.ac.uk/imgt/hla.

RESOLUTION IN HOMO- AND HETEROZYGOTES

Results file with resolution in HLA-C*07 homo- and heterozygotes is available upon request.



101.615-24/04 – including *Taq* polymerase
101.615-24u/04u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 7N8

Lot-specific Information

SPECIFICITY TABLE

HLA-C*07 SSP subtyping

Specificities and sizes of the PCR products of the 95+1 primer mixes used for HLA-C*07 SSP subtyping

Primer Mix	Size of spec. PCR product ¹	Size of control band ²	Amplified HLA-C*07 alleles ³	Other amplified HLA-Class I alleles
1	245 bp	800 bp	*07:01:01:01-07:01:02:14, 07:01:04-07:01:88, 07:01:90-07:01:101, 07:01:103-07:01:115, 07:02:01:01-07:02:20, 07:02:22-07:02:104:02, 07:02:106:01-07:33N, 07:35-07:40, 07:42-07:165, 07:167-07:294, 07:296-07:347N, 07:349-07:622, 07:624-07:765, 07:767-07:798, 07:800-07:831, 07:833N-07:942, 07:944-07:947, 07:949-07:1018, 07:1020-07:1042N	*12:181
	425 bp		*07:01:01:01-07:01:02:14, 07:01:04-07:01:40, 07:01:42-07:01:69, 07:01:71-07:02:40, 07:02:42-07:02:130, 07:02:131 ^w , 07:02:132-07:06:06, 07:08, 07:10-07:33N, 07:35-07:48, 07:50-07:75, 07:77-07:78:01, 07:79-07:114, 07:116-07:209, 07:211-07:222, 07:224-07:237, 07:239-07:246:02, 07:248-07:294, 07:296-07:314:03, 07:316-07:327, 07:329N-07:356, 07:358-07:360, 07:362-07:376, 07:378-07:402, 07:404-07:405:02, 07:407-07:414, 07:416-07:477, 07:478 ^w , 07:479-07:558:01:02, 07:559 ^w , 07:560-07:597, 07:599-07:655, 07:656 ^w , 07:657-07:695:01:02, 07:697Q-07:720, 07:722-07:726N, 07:728-07:740, 07:742-07:817, 07:819Q-07:884, 07:886N-07:913, 07:914 ^w , 07:915-07:994, 07:996-07:999, 07:1001N-07:1042N	*08:252, 12:349, 14:121
2	290 bp	800 bp	*07:01:01:01-07:01:02:14, 07:01:04, 07:01:06-07:01:35, 07:01:37-07:01:112, 07:01:114-07:01:116, 07:05-07:06:01:12, 07:06:03-07:07, 07:09, 07:14, 07:16, 07:18:01:01-07:18:10, 07:20:01-07:22, 07:24, 07:26:01:01-07:28, 07:30-07:31:02, 07:35-07:36, 07:40, 07:43:01-07:44, 07:52-07:53, 07:55N, 07:58-07:60, 07:65, 07:69-07:71, 07:73:01-07:73:02, 07:77-07:78:02, 07:81-07:83, 07:85-07:86, 07:89, 07:91-07:96:02, 07:98N, 07:103-07:104N, 07:106, 07:109:01-07:109:02, 07:111, 07:113, 07:115-07:116, 07:118-07:120, 07:122, 07:124, 07:128-07:129, 07:131:01-07:132, 07:134, 07:140:01-07:141:03, 07:148, 07:150Q-07:151, 07:153, 07:156, 07:162-07:166, 07:170, 07:173, 07:176, 07:179-07:180, 07:182, 07:184, 07:188-07:191N, 07:196-07:197, 07:200-07:201, 07:203-07:207, 07:210, 07:212, 07:214-07:215, 07:219, 07:222-07:225, 07:227N-07:228, 07:231,	*02:02:28, 03:03:29, 03:04:19, 03:113:01, 05:01:17, 05:09:02, 05:17, 06:02:23, 06:153:02, 08:15:02, 08:72:02, 12:03:10, 12:181, 15:02:33:01-15:02:33:02, 15:65, 16:01:30, 16:02:11, A*02:01:126, A*25:27:02, A*26:01:64, B*15:17:04, B*18:01:48, B*35:01:38, B*35:03:15, B*35:08:04, B*35:15:02, B*39:136, B*58:01:17



For *In Vitro* Diagnostic Use
MA124 v02 SSP PI Template
Date: January 2023, Rev. No: 00

101.615-24/04 – including *Taq* polymerase
101.615-24u/04u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **7N8**

Lot-specific Information

07:235Q-07:237, 07:246:01-07:249, 07:253-
07:257:04, 07:263, 07:266-07:269, 07:276-
07:282, 07:292-07:293, 07:296-07:301,
07:303-07:304, 07:310-07:311, 07:313,
07:315, 07:330:01-07:332, 07:337,
07:343:01:01-07:343:01:02, 07:362, 07:366,
07:369-07:371, 07:377, 07:380, 07:385-
07:387, 07:404, 07:407-07:408, 07:410-
07:412, 07:417-07:419, 07:421-07:422,
07:424, 07:427, 07:430-07:431, 07:435,
07:437N, 07:440-07:445, 07:448, 07:450:01-
07:450:02, 07:453-07:455, 07:458, 07:460-
07:462, 07:468-07:471, 07:475-07:476N,
07:479, 07:481, 07:488, 07:490-07:493,
07:498-07:499, 07:502, 07:506:01-07:506:02,
07:508, 07:512, 07:514:01-07:515, 07:517,
07:521:01, 07:522, 07:527-07:528, 07:531,
07:536-07:538, 07:540, 07:545-07:547,
07:549-07:550, 07:554-07:557, 07:559,
07:561, 07:570-07:571, 07:579, 07:581,
07:588-07:589, 07:591, 07:602-07:604,
07:607, 07:610, 07:613-07:615, 07:617-
07:621, 07:623-07:624, 07:633N-07:636,
07:638, 07:642-07:643, 07:647, 07:649-
07:650:02, 07:652-07:654, 07:657-07:658,
07:660, 07:662-07:663Q, 07:668-07:671,
07:673, 07:677, 07:682:01:01-07:683, 07:685,
07:687, 07:690N-07:691, 07:694, 07:696,
07:699, 07:701, 07:705-07:706, 07:708,
07:712-07:714, 07:718, 07:720-07:721,
07:725, 07:728, 07:732, 07:734-07:735,
07:739, 07:749N-07:750N, 07:756-07:757,
07:760, 07:762, 07:771, 07:774, 07:777-
07:778, 07:781-07:782, 07:787N-07:788,
07:794-07:796N, 07:800, 07:802-07:806,
07:808, 07:811-07:816, 07:819Q-07:821N,
07:823-07:825, 07:827, 07:830, 07:832-
07:834, 07:841, 07:843-07:848, 07:855-
07:857, 07:859-07:860, 07:863N-07:864,
07:869-07:870, 07:872-07:873, 07:875,
07:877-07:879, 07:882, 07:884, 07:887,
07:893-07:894, 07:898, 07:903, 07:909-
07:910, 07:914-07:915, 07:918-07:919,
07:921-07:922, 07:925, 07:928, 07:930-
07:935, 07:937, 07:943-07:945, 07:947,
07:949, 07:952, 07:954N-07:955, 07:957,
07:959, 07:961, 07:967, 07:970N, 07:973,
07:975, 07:977, 07:979-07:981N, 07:983-
07:985:02, 07:990-07:993:02, 07:996,
07:1000-07:1002, 07:1005N-07:1008, 07:1013,
07:1017, 07:1022-07:1023, 07:1025, 07:1030-
07:1031, 07:1034-07:1036, 07:1039N,
07:1041-07:1042N



101.615-24/04 – including *Taq* polymerase
101.615-24u/04u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **7N8**

Lot-specific Information

3⁴	110 bp	800 bp	*07:01:01:01-07:01:02:14, 07:01:04-07:01:24, 07:01:26-07:01:116, 07:06:01:01-07:06:06, 07:09, 07:18:01:01-07:22, 07:24, 07:26:01:01- 07:26:03, 07:28, 07:30, 07:35-07:36, 07:40, 07:44, 07:52-07:53, 07:55N, 07:57-07:59, 07:65, 07:69-07:71, 07:73:01-07:73:02, 07:77- 07:78:02, 07:81-07:83, 07:85-07:86, 07:89, 07:91-07:96:02, 07:98N, 07:103-07:104N, 07:106, 07:108:01-07:108:02, 07:110-07:113, 07:115-07:116, 07:118-07:120, 07:122, 07:124, 07:128-07:129, 07:131:01-07:132, 07:134, 07:140:01-07:141:03, 07:148, 07:150Q-07:151, 07:153, 07:156, 07:162- 07:166, 07:170, 07:173, 07:176, 07:179- 07:180, 07:182, 07:188-07:191N, 07:196- 07:197, 07:200-07:201, 07:203-07:207, 07:210, 07:212, 07:214-07:215, 07:219, 07:222-07:225, 07:227N-07:228, 07:230- 07:231, 07:235Q-07:237, 07:246:01-07:250, 07:253-07:257:04, 07:263, 07:266-07:269, 07:276-07:282, 07:293, 07:296-07:301, 07:303-07:304, 07:310-07:311, 07:313, 07:317, 07:330:01-07:332, 07:337, 07:343:01:01-07:343:01:02, 07:353, 07:362, 07:366, 07:369-07:371, 07:377, 07:380, 07:385-07:387, 07:404, 07:407-07:408, 07:410-07:412, 07:417-07:419, 07:421-07:422, 07:424, 07:427, 07:430-07:431, 07:435, 07:440-07:445, 07:448, 07:450:01-07:450:02, 07:453-07:455, 07:458, 07:460-07:463, 07:468-07:471, 07:475, 07:479, 07:481, 07:488, 07:490-07:493, 07:498^w, 07:499, 07:502, 07:506:01-07:506:02, 07:508, 07:512, 07:514:01-07:515, 07:517, 07:521:01-07:522, 07:526:01-07:528, 07:531, 07:536-07:537, 07:540, 07:545-07:547, 07:549-07:550, 07:554-07:555, 07:557, 07:559, 07:561, 07:568, 07:570-07:571, 07:579, 07:581, 07:587-07:589, 07:591, 07:602-07:604, 07:607, 07:610, 07:613-07:615, 07:617- 07:621, 07:623-07:624, 07:627, 07:633N- 07:636, 07:638-07:639, 07:642-07:643, 07:647, 07:649-07:650:02, 07:652-07:654, 07:657-07:658, 07:660, 07:662-07:663Q, 07:669-07:671, 07:673, 07:677, 07:682:01:01- 07:683, 07:685, 07:687, 07:690N-07:691, 07:694, 07:696, 07:699, 07:701, 07:705- 07:706, 07:708, 07:712-07:714, 07:718, 07:720-07:721, 07:725-07:726N, 07:728, 07:731-07:732, 07:734-07:735, 07:739, 07:746N, 07:749N-07:750N, 07:756-07:758, 07:760, 07:762, 07:771, 07:774, 07:776N- 07:778, 07:781-07:782, 07:786-07:788, 07:794-07:796N, 07:800, 07:802-07:806, 07:808, 07:811-07:816, 07:819Q-07:821N, 07:823-07:825, 07:827, 07:832-07:834,	*04:112, 06:02:13, 06:44, 06:252, 12:287, 15:62, 15:169, 16:64, 18:05
----------------------	---------------	---------------	---	--



101.615-24/04 – including *Taq* polymerase
101.615-24u/04u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **7N8**

Lot-specific Information

			07:841, 07:843-07:848, 07:855-07:857, 07:859-07:860, 07:863N-07:864, 07:869- 07:870, 07:872-07:873, 07:875, 07:877- 07:879, 07:882, 07:884, 07:887, 07:893- 07:894, 07:898, 07:903, 07:909-07:910, 07:915, 07:918-07:921, 07:925, 07:928, 07:930-07:935, 07:943-07:945, 07:947, 07:949, 07:952, 07:954N-07:955, 07:957, 07:959, 07:961, 07:967, 07:970N, 07:973, 07:977, 07:979-07:981N, 07:983-07:985:02, 07:990-07:992, 07:996, 07:1000-07:1002, 07:1005N-07:1008, 07:1017, 07:1022-07:1023, 07:1025, 07:1030, 07:1034-07:1036, 07:1039N, 07:1041-07:1042N	
4	185 bp	1070 bp	*07:02:01:01-07:02:93, 07:02:95-07:02:97, 07:02:99-07:02:126, 07:02:128-07:02:144, 07:10, 07:13:01:01-07:13:01:02, 07:15, 07:19, 07:23, 07:25, 07:29:01:01-07:29:02, 07:32N- 07:33N, 07:37:01:02-07:39, 07:42, 07:46- 07:51, 07:54, 07:56:01-07:56:02, 07:61N- 07:62, 07:64, 07:66-07:68, 07:72, 07:74- 07:76:02, 07:79-07:80, 07:84, 07:87-07:88, 07:90, 07:97, 07:99-07:100, 07:102, 07:105, 07:107, 07:114, 07:117, 07:121Q, 07:123:01- 07:123:02, 07:125:01-07:127:02, 07:130, 07:133:01-07:133:02, 07:135-07:138, 07:143- 07:147, 07:149, 07:152N, 07:154-07:155, 07:157-07:161, 07:167-07:169, 07:172:01- 07:172:02, 07:174-07:175, 07:178, 07:183, 07:185-07:187, 07:192-07:195, 07:198N, 07:202, 07:208-07:209, 07:211, 07:213, 07:216-07:218, 07:220-07:221, 07:226Q, 07:229, 07:232-07:234, 07:238-07:241, 07:243-07:245, 07:251-07:252, 07:259-07:262, 07:264N-07:265, 07:270-07:271, 07:273- 07:275, 07:283-07:291, 07:302, 07:305- 07:306, 07:308-07:309, 07:312, 07:314:01- 07:314:03, 07:316-07:322, 07:325-07:327, 07:333-07:335, 07:339-07:342:02, 07:344- 07:352, 07:356, 07:359-07:360, 07:363, 07:367-07:368:02, 07:372-07:376, 07:379, 07:381-07:384, 07:388-07:393N, 07:397- 07:402, 07:405:01-07:405:02, 07:409, 07:413- 07:416, 07:423, 07:425, 07:429, 07:432Q- 07:434, 07:436, 07:438-07:439, 07:446, 07:449, 07:451N-07:452N, 07:456-07:457, 07:464-07:465, 07:472-07:474, 07:477-07:478, 07:482-07:486, 07:489:01-07:489:02, 07:494Q-07:497, 07:500, 07:503-07:505, 07:507N, 07:509-07:511, 07:513Q, 07:516, 07:518-07:520, 07:524-07:525, 07:529- 07:530:01, 07:532-07:533, 07:539, 07:541- 07:544, 07:548, 07:551N, 07:553, 07:560, 07:564-07:567, 07:572-07:578, 07:580, 07:582Q-07:584, 07:590, 07:592-07:598, 07:601, 07:605-07:606, 07:608,	*17:48, B*40:60, B*40:245:01-40:245:02



101.615-24/04 – including *Taq* polymerase
101.615-24u/04u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **7N8**

Lot-specific Information

			07:611-07:612, 07:616, 07:625, 07:628- 07:632Q, 07:637, 07:639-07:641, 07:644- 07:646, 07:648, 07:659:01:01-07:659:01:02, 07:661, 07:665-07:667:02, 07:675N-07:676, 07:678-07:679, 07:684, 07:686N, 07:688- 07:689, 07:692, 07:695:01:01-07:695:01:02, 07:697Q, 07:700, 07:702N-07:704, 07:707, 07:709, 07:711, 07:715:01-07:717, 07:719, 07:722-07:724, 07:727, 07:729:01:01-07:730, 07:733N, 07:736:01:01-07:738, 07:740-07:741, 07:743N-07:745N, 07:747N-07:748, 07:752N- 07:755, 07:759, 07:761, 07:763-07:770N, 07:772-07:773N, 07:775, 07:779:01:01- 07:779:01:02, 07:783-07:785, 07:789-07:793, 07:798-07:799, 07:801, 07:807N, 07:809- 07:810, 07:817-07:818, 07:822, 07:826, 07:828-07:829, 07:835-07:837, 07:840N, 07:842, 07:849N, 07:851, 07:853-07:854, 07:861-07:862, 07:866-07:867, 07:871, 07:874, 07:880, 07:883, 07:885-07:886N, 07:888-07:891, 07:900, 07:902, 07:904- 07:908, 07:911-07:913, 07:916, 07:923- 07:924, 07:927, 07:929, 07:936N, 07:938- 07:939, 07:941-07:942, 07:946, 07:950, 07:956, 07:960, 07:962-07:966, 07:968- 07:969, 07:971-07:972, 07:976, 07:978, 07:982, 07:987-07:989, 07:994-07:995, 07:997-07:999, 07:1003-07:1004, 07:1009, 07:1011, 07:1014Q-07:1016, 07:1018- 07:1021, 07:1024, 07:1026-07:1029, 07:1032, 07:1037N-07:1038, 07:1040N	
5^{4,6}	95 bp	800 bp	*07:53, 07:216	*01:215, 16:01:35, 16:02:03, A*24:172:02, B*44:347, B*44:527, B*50:76, B*51:55 A*25:76, B*18:63, B*39:90, B*44:486, B*50:11
	170 bp		*07:57	
	245 bp		*07:03	
6⁴	100 bp	1070 bp	*07:04:01:01-07:04:28, 07:11-07:12, 07:63, 07:68, 07:101, 07:139, 07:142, 07:181, 07:199:01-07:199:02, 07:272, 07:302, 07:323- 07:324, 07:328-07:329N, 07:336, 07:338, 07:354-07:355, 07:357, 07:361, 07:364- 07:365, 07:378, 07:394-07:395, 07:403, 07:406, 07:420, 07:426, 07:428, 07:459, 07:466-07:467, 07:480, 07:487, 07:501, 07:534-07:535, 07:562-07:563, 07:569, 07:585-07:586, 07:600:01N-07:600:02N, 07:622, 07:626, 07:651, 07:655-07:656, 07:664, 07:672N, 07:674, 07:693, 07:698, 07:742, 07:751N, 07:780, 07:797N, 07:831, 07:838-07:839N, 07:852, 07:858, 07:868, 07:876, 07:881N, 07:892, 07:895-07:897, 07:901, 07:917, 07:926, 07:941, 07:948,	*05:62, 05:146, 05:163, 08:142, B*44:148



101.615-24/04 – including *Taq* polymerase
101.615-24u/04u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **7N8**

Lot-specific Information

7	190 bp	1070 bp	07:951, 07:953, 07:974Q, 07:1010, 07:1012, 07:1033 *07:05, 07:538	B*27:36, B*40:469, B*48:21, B*48:26
	285 bp	1070 bp	*07:45, 07:589, 07:611	*01:158, 02:57, 04:40, 04:242, 06:123, 06:358, A*02:84, B*15:554, B*38:23, B*39:73
8	505 bp		*07:06:01:01-07:06:06, 07:330:02, 07:343:01:01-07:343:01:02, 07:613, 07:620, 07:706, 07:721, 07:735, 07:739, 07:760, 07:776N, 07:803, 07:827, 07:841, 07:887, 07:932, 07:959, 07:979	
	215 bp 425 bp	800 bp	*07:33N, 07:936N *07:07, 07:09, 07:49, 07:76:01-07:76:02, 07:210, 07:238, 07:247, 07:315, 07:328, 07:403, 07:406, 07:598	*04:397, 04:432, 04:442, 06:17, 18:01:01:01-18:16
10	185 bp	1070 bp	*07:08, 07:108:01-07:108:02	*04:58, 04:160:01- 04:160:02, 04:368, B*40:100
11	285 bp 305 bp	800 bp	*07:24, 07:112, 07:218, 07:396 *07:10, 07:43:01-07:43:02, 07:184, 07:196, 07:367, 07:568	*03:03:29, 03:04:19, 03:113:01, 15:02:33:01- 15:02:33:02, B*35:01:38, B*35:03:15, B*35:08:04, B*35:15:02, B*37:04:02, B*57:01:18, B*57:03:04, B*58:01:17
12	130 bp	1070 bp	*07:12, 07:41, 07:447, 07:892, 07:953	*03:251, 03:314, 03:607, 04:120, 05:01:01:01- 05:01:26, 05:01:28- 05:01:60, 05:01:62- 05:01:72, 05:03-05:08, 05:10-05:16, 05:18:01- 05:22:02, 05:24-05:48N, 05:50-05:51Q, 05:53- 05:61, 05:63-05:94, 05:96-05:109, 05:111- 05:124, 05:126-05:132, 05:135, 05:137-05:142, 05:144-05:150, 05:152- 05:200, 05:202Q-05:205, 05:207-05:236, 05:238- 05:243, 05:246-05:261, 05:263N-05:277, 06:129, 08:01:01:01-08:01:07, 08:01:09-08:06, 08:08:01- 08:14, 08:16:01-08:46, 08:48-08:50, 08:52N- 08:63, 08:65-08:71, 08:73-08:75, 08:77-08:83, 08:85-08:103:01:02, 08:105-08:123, 08:125- 08:156, 08:158-08:159, 08:161N-08:168,



101.615-24/04 – including *Taq* polymerase
101.615-24u/04u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **7N8**

Lot-specific Information

				08:170-08:187, 08:189, 08:191-08:215, 08:217- 08:242, 08:245, 08:247- 08:256, 15:130, 15:229, A*29:10:01-29:10:02, A*29:161, B*14:32, B*15:337, B*15:627, B*18:83, B*44:148, B*51:355, B*58:76
13^{4,6}	85 bp	1070 bp	*07:11, 07:56:01-07:56:02, 07:60, 07:79, 07:109:01-07:109:02, 07:682:01:01-07:682:02, 07:736:01:01-07:736:01:02, 07:739	*03:69, 03:274, 03:461, 04:27, 04:52, 04:405, 06:53:01-06:53:02, 12:58, 12:252, 17:26, A*24:443, A*68:275
14	215 bp	1070 bp	*07:02:01:01-07:02:41, 07:02:42 ^w , 07:02:43- 07:05, 07:08, 07:10-07:15, 07:17:01:01- 07:17:04, 07:23, 07:25, 07:27:01-07:27:02, 07:29:01:01-07:29:02, 07:31:01-07:32N, 07:37:01:02-07:39, 07:41-07:43:02, 07:45- 07:51, 07:54, 07:56:01-07:56:02, 07:62-07:63, 07:66-07:68, 07:72, 07:74-07:76:02, 07:79- 07:80, 07:84, 07:87-07:88, 07:90, 07:97, 07:99-07:102, 07:105, 07:107, 07:114, 07:117, 07:121Q, 07:123:01-07:123:02, 07:125:01- 07:127:02, 07:130, 07:133:01-07:133:02, 07:135-07:138, 07:142-07:147, 07:149, 07:152N, 07:154-07:155, 07:157-07:161, 07:167-07:169, 07:171-07:172:02, 07:174- 07:175, 07:177-07:178, 07:181, 07:183- 07:187, 07:192-07:195, 07:198N-07:199:02, 07:202, 07:208-07:209, 07:211, 07:213, 07:216-07:218, 07:220-07:221, 07:226Q, 07:229, 07:232-07:234, 07:238-07:245, 07:251-07:252, 07:258-07:262, 07:264N- 07:265, 07:270-07:275, 07:283-07:290, 07:291 ^w , 07:294, 07:302, 07:305-07:309, 07:312, 07:315-07:316, 07:318-07:321, 07:323-07:329N, 07:333-07:334, 07:335 ^w , 07:336, 07:338-07:342:02, 07:344-07:350N, 07:352, 07:354-07:361, 07:363-07:365, 07:367-07:368:02, 07:372-07:376, 07:378- 07:379, 07:381-07:384, 07:388-07:390, 07:392-07:401, 07:403, 07:405:01-07:406, 07:409, 07:413-07:416, 07:420, 07:423, 07:425-07:426, 07:428-07:429, 07:432Q- 07:434, 07:436, 07:438-07:439, 07:446- 07:447, 07:451N-07:452N, 07:456-07:457, 07:459, 07:464-07:467, 07:472-07:474, 07:477-07:478, 07:480, 07:482-07:487, 07:489:01-07:489:02, 07:494Q-07:497, 07:500-07:501, 07:503-07:505, 07:507N, 07:509-07:511, 07:513Q, 07:516, 07:519- 07:520, 07:523-07:525, 07:529-07:530:02, 07:532-07:535, 07:538-07:539, 07:541-07:544, 07:548, 07:551N-07:553, 07:556, 07:558:01:01-07:558:01:02, 07:560,	*02:14:01-02:14:02, 02:107, 02:164, 03:39, 03:67, 03:344:01:01- 03:344:01:02, 04:42:01- 04:42:02, 04:220, 05:43, 06:02:01:01-06:02:01:92, 06:02:03-06:02:103, 06:04:01-06:43:02, 06:45- 06:75, 06:78-06:131, 06:133-06:146, 06:148- 06:192, 06:194-06:251, 06:253-06:339, 06:341- 06:353N, 06:355-06:360, 08:37, 12:16:01-12:16:02, 12:147, 12:195:01- 12:195:03, 12:217, 14:147, 16:21, 16:80, 17:20, 18:01:01:01-18:04, 18:06-18:08, 18:10-18:16



101.615-24/04 – including *Taq* polymerase
101.615-24u/04u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **7N8**

Lot-specific Information

			07:562-07:567, 07:569, 07:572-07:577, 07:580, 07:582Q, 07:584-07:586, 07:590, 07:592-07:601, 07:605-07:606, 07:608-07:609, 07:611-07:612, 07:616, 07:622, 07:625- 07:626, 07:628-07:632Q, 07:637, 07:640- 07:641, 07:644-07:646, 07:648, 07:651, 07:655-07:656, 07:659:01:01-07:659:01:02, 07:661, 07:664-07:667:02, 07:672N, 07:674- 07:676, 07:678-07:681, 07:684, 07:686N, 07:688-07:689, 07:692-07:693, 07:695:01:01- 07:695:01:02, 07:697Q-07:698, 07:700, 07:702N-07:704, 07:707, 07:709-07:711, 07:715:01-07:717, 07:719, 07:722, 07:727, 07:729:01:01-07:730, 07:733N, 07:736:01:01- 07:738, 07:740-07:745N, 07:747N-07:748, 07:751N-07:755, 07:759, 07:761, 07:763- 07:770N, 07:772-07:773N, 07:775, 07:779:01:01-07:780, 07:783-07:785, 07:789- 07:793, 07:797N-07:799, 07:801, 07:807N, 07:809-07:810, 07:817-07:818, 07:822, 07:826, 07:828-07:831, 07:835-07:840N, 07:842, 07:849N-07:854, 07:858, 07:861- 07:862, 07:865-07:868, 07:871, 07:874, 07:876, 07:880-07:881N, 07:883, 07:885- 07:886N, 07:888-07:890, 07:892, 07:895- 07:897, 07:899-07:902, 07:904-07:908, 07:911-07:913, 07:916-07:917, 07:923-07:924, 07:926-07:927, 07:929, 07:938-07:942, 07:946, 07:948, 07:950-07:951, 07:953, 07:956, 07:958, 07:960, 07:962, 07:964- 07:966, 07:968-07:969, 07:971-07:972, 07:974Q-07:976, 07:978, 07:982, 07:986:01- 07:988, 07:994-07:995, 07:997-07:999, 07:1003-07:1004, 07:1009-07:1016, 07:1018- 07:1021, 07:1024, 07:1026, 07:1028-07:1029, 07:1032-07:1033, 07:1038, 07:1040N	
15^{4,7}	80 bp	1070 bp	*07:50, 07:615, 07:735, 07:781	*01:217:01:01- 01:217:01:02, 02:175, 03:562, 04:320, 05:261, 16:111, 17:38:01:01- 17:38:01:02, 17:49
	290 bp 460 bp 240 bp		*07:14 *07:89, 07:927	
16⁷		800 bp	*07:13:01:01-07:13:01:02, 07:28, 07:289, 07:536, 07:725, 07:834, 07:857, 07:1028	*04:01:43, 05:01:17, 05:09:02, 05:17, 08:15:02, 08:72:02, 15:02:33:01- 15:02:33:02, B*35:03:15, B*37:04:02, B*49:01:22
17⁶	130 bp	1070 bp	*07:04:01:01-07:04:25, 07:04:27-07:04:28, 07:11, 07:63, 07:68, 07:101, 07:139, 07:142, 07:181, 07:199:01-07:199:02, 07:260:01- 07:260:02, 07:272, 07:302, 07:323-07:324, 07:328-07:329N, 07:336, 07:338, 07:354- 07:355, 07:357-07:358, 07:361, 07:364- 07:365, 07:378, 07:394-07:395, 07:403,	*05:23, 05:62, 05:134, 05:151, 08:07, 08:47, 08:104, 08:188



101.615-24/04 – including *Taq* polymerase
101.615-24u/04u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **7N8**

Lot-specific Information

			07:406, 07:420, 07:426, 07:428, 07:459, 07:466-07:467, 07:480, 07:487, 07:501, 07:523, 07:534-07:535, 07:552, 07:562- 07:563, 07:569, 07:585-07:586, 07:600:01N- 07:600:02N, 07:622, 07:626, 07:651, 07:655- 07:656, 07:664, 07:672N, 07:674, 07:693, 07:698, 07:742, 07:751N, 07:780, 07:797N, 07:831, 07:838-07:839N, 07:852, 07:858, 07:868, 07:876, 07:881N, 07:895-07:897, 07:901, 07:917, 07:926, 07:941, 07:948, 07:951, 07:974Q, 07:1010, 07:1012, 07:1033	
18⁴	120 bp 230 bp	800 bp	*07:52, 07:243 *07:15, 07:337, 07:500	*03:490, 03:548, 05:03 *02:77, B*07:77, B*07:193, B*40:60, B*40:245:01-40:245:02
	265 bp		*07:486	
19⁵	220 bp	800 bp	*07:26:01:01-07:26:03, 07:92, 07:96:01- 07:96:02, 07:314:01-07:314:03, 07:317, 07:351, 07:578, 07:583, 07:723, 07:756, 07:821N, 07:845, 07:870, 07:891, 07:894	*08:252, 12:349, 18:09
	465 bp		*07:07, 07:16, 07:51, 07:177, 07:181, 07:328, 07:367, 07:598, 07:723, 07:914, 07:993:01- 07:993:02	*08:252, 12:349, 14:121
20^{5,6}	145 bp	1070 bp	*07:03, 07:17:01:01-07:17:04, 07:140:01- 07:141:03, 07:304, 07:307, 07:512, 07:599, 07:671, 07:812, 07:850, 07:932, 07:940, 07:986:01-07:986:02	*02:02:28, 03:03:29, 03:04:19, 03:113:01, 04:01:43, 05:01:17, 05:09:02, 05:17, 06:02:23, 06:153:02, 08:15:02, 08:72:02, 12:03:10, 15:02:33:01-15:02:33:02, 16:01:30, 16:02:11, A*25:27:02, A*26:01:64, A*30:01:03, A*31:01:32, A*33:03:09, B*07:02:43, B*07:06:03, B*07:33:03, B*07:139:02, B*08:01:25, B*15:17:04, B*18:01:48, B*35:01:38, B*35:03:15, B*35:08:04, B*35:15:02, B*37:04:02, B*39:09:03, B*49:01:22, B*51:01:63, B*57:01:18, B*57:03:04, B*58:01:17
	170 bp 190 bp		*07:44, 07:48 *07:81, 07:168, 07:450:01-07:450:02, 07:959	B*15:17:04, B*18:01:48, B*35:01:38, B*35:03:15, B*35:08:04, B*35:15:02, B*37:04:02, B*39:09:03, B*49:01:22, B*51:01:63, B*57:01:18, B*57:03:04, B*58:01:17
21⁴	100 bp 225 bp	800 bp	*07:54 *07:06:01:01-07:06:06, 07:18:01:01-07:19, 07:28, 07:71, 07:112, 07:173, 07:330:02, 07:517, 07:607, 07:613, 07:620-07:621, 07:687, 07:706, 07:721, 07:735, 07:739,	



101.615-24/04 – including *Taq* polymerase
101.615-24u/04u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **7N8**

Lot-specific Information

			07:760, 07:776N, 07:781, 07:783, 07:803, 07:805, 07:815, 07:819Q, 07:827, 07:830, 07:841, 07:845, 07:847, 07:887, 07:919, 07:932, 07:952, 07:955, 07:959, 07:977, 07:979, 07:1005N, 07:1035-07:1036	
22	205 bp 250 bp	800 bp	*07:84 *07:20:01-07:20:02, 07:96:01-07:96:02, 07:127:01-07:127:02, 07:263, 07:390, 07:562, 07:578	*02:49, 02:75, 02:115, 03:15, 03:32:01:01- 03:32:01:02, 03:45, 03:60, 03:136, 03:297, 04:03:01:01-04:03:11, 04:06:01-04:06:03, 04:16, 04:80, 04:147, 04:160:01- 04:160:02, 04:171, 04:190, 04:256, 04:286, 04:294, 04:299, 04:335, 04:337, 04:351:01- 04:351:02, 04:357, 04:363, 04:381, 04:383, 04:393, 04:400, 04:402, 04:473, 04:489, 06:70:01- 06:70:02, 06:132:01- 06:132:02, 06:305, 14:116, 15:25, 15:62, 15:169, A*01:110 , A*11:135
23⁶	155 bp 185 bp	1070 bp	*07:23 *07:21	
24⁶	135 bp	1070 bp	*07:25, 07:404	A*02:01:126, A*03:01:46, A*03:154:02, A*11:32:02, A*24:02:62, A*24:347:01, A*24:459:01, A*25:27:02, A*26:01:64, A*68:01:12
25⁴	105 bp 150 bp 235 bp	1070 bp	*07:46, 07:356 *07:30 *07:02:01:17N, 07:311	*03:342, B*56:08^w , B*56:14
26⁴	100 bp 250 bp 300 bp	1070 bp	*07:32N *07:55N *07:31:01-07:31:02, 07:177, 07:514:01- 07:514:02	*06:02:23, 06:153:02, 12:03:10, 16:01:30, 16:02:11
27	130 bp	1070 bp	*07:37:01:02, 07:69, 07:176, 07:307, 07:512, 07:755, 07:850, 07:967, 07:1016	*04:01:43, 04:400 ^w , 05:01:17, 05:09:02, 08:15:02, B*07:02:43 , B*07:06:03, B*07:33:03 , B*07:139:02, B*35:08:04
28	240 bp 175 bp 260 bp 290 bp	800 bp	*07:267, 07:322 *07:44 *07:35, 07:47:01, 07:352, 07:461 *07:396	*04:198 ^w , 12:115
29	255 bp	1070 bp	*07:36, 07:42, 07:355	*04:258, A*11:350
30^{4,6}	60 bp 260 bp	800 bp	*07:75, 07:254, 07:951 *07:227N	*04:66



101.615-24/04 – including *Taq* polymerase
101.615-24u/04u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **7N8**

Lot-specific Information

31	420 bp	1070 bp	*07:39-07:40, 07:177, 07:210 ^w , 07:238 ^w , 07:244, 07:276, 07:328 ^w , 07:335, 07:480, 07:563, 07:975, 07:993:01-07:993:02	*08:252, 14:121
32⁴	115 bp 155 bp 245 bp	1070 bp	*07:134 *07:189 *07:07, 07:09, 07:49, 07:76:01-07:76:02, 07:210, 07:238, 07:247, 07:315, 07:328, 07:403, 07:406, 07:598	*06:29 *02:14:01-02:14:02, 02:107, 02:164, 04:42:01- 04:42:02, 04:220, 05:43, 06:02:01-01-06:02:01:92, 06:02:03-06:02:05, 06:02:07-06:02:11, 06:02:13-06:02:31, 06:02:33-06:02:64, 06:02:66-06:02:103, 06:04:01-06:10, 06:12- 06:51, 06:53:01-06:75, 06:78-06:121, 06:123, 06:125-06:131, 06:133- 06:146, 06:148-06:216, 06:218-06:247, 06:249- 06:251, 06:253-06:279, 06:281N-06:360, 15:23:01-15:23:02, 15:63, 15:138, 15:158, 17:20, 18:01:01-01-18:08, 18:10- 18:16, B*08:17, B*08:47, B*08:295
	285 bp		*07:60, 07:79	*03:39, 03:67, B*08:39
33	205 bp 230 bp	1070 bp	*07:66 *07:82	
34⁵	335 bp	1070 bp	*07:70, 07:503	*06:36
35	320 bp 360 bp	1070 bp	*07:72 *07:71	
36⁵	265 bp	800 bp	*07:20:01-07:20:02, 07:64, 07:73:01-07:73:02, 07:92, 07:96:01-07:96:02, 07:172:01- 07:172:02, 07:390, 07:578, 07:583, 07:723, 07:891	*04:397, 04:432, 04:442, 08:252, 12:349, 14:121, 18:04
37⁶	235 bp 265 bp	1070 bp	*07:121Q *07:38:01-07:38:02, 07:63, 07:68, 07:260:01- 07:260:02, 07:302, 07:941	*05:111, 05:207, 08:38, B*07:77, B*07:193, B*40:60, B*40:245:01- 40:245:02
38⁴	110 bp	1070 bp	*07:256, 07:497	*01:53, 04:379, 06:189, B*07:175, B*15:126, B*46:63
	150 bp		*07:65, 07:717	*04:222, B*07:189, B*39:148
39⁴	100 bp 195 bp 255 bp	800 bp	*07:74 *07:308 *07:59	
40^{4,6}	120 bp 205 bp 440 bp	800 bp	*07:235Q *07:62 *07:07, 07:09, 07:76:01-07:76:02, 07:78:01- 07:78:02, 07:315, 07:328, 07:406, 07:559, 07:598, 07:656, 07:914	B*14:70Q *04:397, 04:432, 04:442, 06:17, 18:01:01-01-18:16



101.615-24/04 – including *Taq* polymerase
101.615-24u/04u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **7N8**

Lot-specific Information

41	140 bp 170 bp 255 bp	1070 bp	*07:91, 07:474 *07:234 *07:88	*04:419, A*01:311, A*25:32
42⁴	100 bp 205 bp 250 bp	1070 bp	*07:15, 07:207, 07:850 *07:77, 07:93:01-07:93:02, 07:905:01:01- 07:905:01:02 *07:35, 07:123:01-07:123:02, 07:173, 07:294, 07:352, 07:725, 07:775, 07:850	*02:02:28, B*07:02:43, B*07:06:03, B*07:33:03, B*07:139:02, B*35:15:02 *15:02:33:01-15:02:33:02, B*07:02:43, B*07:06:03, B*07:33:03, B*07:139:02, B*08:01:25, B*18:01:48, B*35:01:38, B*35:03:15, B*35:08:04, B*35:15:02, B*39:09:03, B*39:136, B*51:01:63, B*57:01:18, B*57:03:04, B*58:01:17
43^{4,5}	105 bp 150 bp 175 bp 220 bp	800 bp	*07:94, 07:183 *07:99, 07:189 *07:80, 07:339, 07:410, 07:662 *07:445	
44⁵	165 bp 195 bp 290 bp	1070 bp	*07:98N *07:65, 07:77, 07:717 *07:396	
45^{5,6}	240 bp 475 bp	1070 bp	*07:245 *07:95, 07:342:01-07:342:02, 07:599	*04:34, 04:122, 04:212, 04:265, 06:153:01- 06:153:02, A*26:112:01, B*18:115, B*35:178, B*35:282, B*35:316, B*38:79, B*57:102, B*73:01:01:01-73:03
46⁵	315 bp 350 bp 420 bp	1070 bp	*07:202 *07:90, 07:757 *07:152N	
47⁵	200 bp 255 bp 420 bp 455 bp	1070 bp	*07:86 *07:222, 07:375 *07:152N *07:97	
48	155 bp 220 bp 255 bp 410 bp	1070 bp	*07:122, 07:373 *07:87 *07:222, 07:375 *07:85, 07:244, 07:276, 07:480	
49⁴	70 bp 155 bp 230 bp	800 bp	*07:133:01-07:133:02, 07:330:01-07:330:02, 07:748 *07:101, 07:148, 07:161, 07:583 *07:100, 07:259, 07:316, 07:338, 07:579	*03:03:29, 03:04:19 *04:129, 05:42, 05:46:01:01-05:46:01:02, 06:67, 08:05, 08:21, 08:137, 12:127, 12:203, 17:05 *01:17, 01:21, 01:23, 01:128, 01:152, 01:157, 12:215



101.615-24/04 – including *Taq* polymerase
101.615-24u/04u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **7N8**

Lot-specific Information

50	345 bp 425 bp	800 bp	*07:310 *07:124, 07:155	
51⁴	65 bp 235 bp 285 bp	1070 bp	*07:126 *07:121Q, 07:409, 07:483N *07:117	
52⁴	65 bp 205 bp 280 bp	1070 bp	*07:203, 07:426 *07:118 *07:106, 07:691	*06:194
53⁴	70 bp 110 bp 155 bp	1070 bp	*07:133:01-07:133:02, 07:330:01-07:330:02, 07:748 *07:140:01-07:140:02, 07:151, 07:154, 07:929 *07:122, 07:373	*03:03:29, 03:04:19 B*08:01:25, B*37:04:02
54^{4,5}	120 bp 200 bp 235 bp 285 bp 340 bp 370 bp 430 bp	1070 bp	*07:153, 07:324, 07:738 *07:446 *07:64, 07:374, 07:402 *07:593N *07:125:01-07:125:02, 07:356, 07:531, 07:1010 *07:129 *07:01:02:01-07:01:02:14, 07:01:27, 07:01:52, 07:01:80-07:01:81, 07:01:83, 07:01:93, 07:01:103, 07:01:112-07:01:113, 07:73:02, 07:89, 07:93:01, 07:141:03, 07:450:02, 07:453, 07:514:02, 07:547, 07:627, 07:668, 07:806, 07:814, 07:856N, 07:877, 07:903, 07:933, 07:985:02, 07:993:01, 07:1006, 07:1042N	B*35:364, B*35:379, B*35:462, B*42:31, B*51:142, B*52:80 B*08:109 *04:397, 04:432, 04:442, 18:12 B*52:103N *04:397, 04:432, 04:442
55	160 bp 210 bp 290 bp 530 bp	1070 bp	*07:179, 07:234, 07:679 *07:252 *07:132 *07:343:01:01-07:343:01:02	*04:353, 06:07 A*11:179
56⁶	210 bp 260 bp	1070 bp	*07:150Q *07:123:01-07:123:02, 07:143, 07:173, 07:294, 07:421, 07:524, 07:725, 07:775, 07:850	*15:02:33:01-15:02:33:02, B*07:02:43, B*07:06:03, B*07:33:03, B*07:139:02, B*08:01:25, B*18:01:48, B*35:01:38, B*35:03:15, B*35:08:04, B*35:15:02, B*39:09:03, B*39:136, B*51:01:63, B*57:01:18, B*57:03:04, B*58:01:17
57⁴	70 bp 175 bp 230 bp 270 bp	800 bp	*07:01:25 *07:130, 07:915 *07:259 *07:421, 07:524	*05:01:17, 05:09:02, 05:17, 08:15:02
58⁴	80 bp 450 bp	1070 bp	*07:137:01-07:137:02, 07:186, 07:212, 07:265, 07:748, 07:868 *07:164N, 07:451N	*03:162, 03:220, 03:434, 05:151, B*35:509, B*38:43, B*44:247, B*49:29, B*50:70, B*51:181, B*52:74



101.615-24/04 – including *Taq* polymerase
101.615-24u/04u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **7N8**

Lot-specific Information

59⁴	110 bp	800 bp	*07:81, 07:145, 07:168, 07:191N, 07:450:01-07:450:02, 07:733N, 07:959	B*15:88, B*15:162, B*15:200, B*39:82, B*44:247, B*44:295:01-44:295:02, B*47:09, B*49:45, B*51:115, B*51:181, B*56:59, B*58:05
60	210 bp 130 bp	1070 bp	*07:485 *07:165, 07:241, 07:265, 07:304, 07:307, 07:512, 07:671, 07:755, 07:812, 07:815, 07:850, 07:940, 07:967, 07:986:01-07:986:02	*02:02:28, 03:03:29, 03:04:19, 03:113:01, 04:01:43, 05:01:17, 05:09:02, 06:02:23, 08:15:02, 12:03:10, 15:02:33:01-15:02:33:02, B*07:02:43, B*07:06:03, B*07:33:03, B*07:139:02, B*15:17:04, B*49:01:22, B*51:01:63
61⁴	300 bp 370 bp 120 bp	1070 bp	*07:147, 07:485 *07:180 *07:81, 07:243, 07:256, 07:462, 07:497, 07:878	*01:53, 04:379, 05:149, 06:189, 16:76, B*07:175, B*07:414, B*15:126, B*40:304, B*46:63
62	245 bp 200 bp 285 bp	1070 bp	*07:120 *07:84, 07:205 *07:106, 07:267, 07:322, 07:691	*06:267N A*26:101 *04:198, A*01:342, A*26:94
63	275 bp 320 bp	800 bp	*07:24, 07:218, 07:301, 07:425 *07:198N, 07:202	
64⁶	135 bp 230 bp	1070 bp	*07:69, 07:1016 *07:192, 07:272, 07:316, 07:326, 07:338, 07:409, 07:579, 07:898, 07:1037N	*04:400 ^w *01:17, 01:21, 01:23, 01:128, 01:152, 01:157, 06:18 ^w , 12:215
65⁶	150 bp 185 bp	1070 bp	*07:136, 07:264N *07:371, 07:494Q	B*18:224 *06:296, 16:16Q, B*15:218Q, B*35:333Q
66⁶	215 bp 150 bp 185 bp	800 bp	*07:551N, 07:632Q *07:136 *07:128, 07:232, 07:513Q, 07:633N	B*18:224 *06:323N, 12:42Q, B*18:153, B*35:461N, B*46:07N, B*46:51Q
67	225 bp 180 bp 215 bp 290 bp	1070 bp	*07:279 *07:285 *07:393N *07:456	
68⁴	75 bp 230 bp 260 bp 290 bp	1070 bp	*07:381 *07:149 *07:133:01, 07:242, 07:748 *07:456	*03:178, 03:300, 03:474
69	135 bp 185 bp 240 bp	1070 bp	*07:170-07:171, 07:397, 07:412 *07:494Q *07:02:01:17N	*16:16Q, B*15:218Q, B*35:333Q



101.615-24/04 – including *Taq* polymerase
101.615-24u/04u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **7N8**

Lot-specific Information

70⁶	215 bp 345 bp 405 bp	1070 bp	*07:131:01-07:131:02, 07:194, 07:368:01-07:368:02 *07:310 *07:61N, 07:1037N	
71⁴	105 bp 120 bp 240 bp	1070 bp	*07:67 *07:347N *07:29:01:01-07:29:02, 07:177, 07:184, 07:248, 07:850, 07:1003, 07:1013	*03:03:29, 03:04:19, 03:113:01, B*07:02:43, B*07:33:03, B*07:139:02
72⁴	280 bp 110 bp	1070 bp	*07:425 *07:01:110, 07:22, 07:67, 07:137:01-07:138, 07:812, 07:847, 07:940	*02:02:28, 06:02:23, 12:03:10, A*02:01:126, A*31:01:32, A*33:03:09, A*68:01:12, B*39:09:03
	245 bp 290 bp		*07:484N *07:250, 07:558:01:01-07:558:01:02	
73	185 bp	1070 bp	*07:107, 07:224, 07:349, 07:513Q	*06:22, 12:42Q, 12:293, 14:119, 16:31, B*15:193, B*35:283, B*46:51Q
	210 bp		*07:458, 07:632Q	
74⁶	210 bp 355 bp	1070 bp	*07:33N, 07:293, 07:519, 07:936N *07:80, 07:662	
75	140 bp 220 bp 265 bp 295 bp	1070 bp	*07:116, 07:329N *07:544 *07:486 *07:708	*12:211 *02:200, 03:492, 03:497, 05:273, 06:349, 12:254, 12:280
76⁴	115 bp 280 bp	1070 bp	*07:437N *07:102	
77⁶	200 bp 240 bp 430 bp	1070 bp	*07:446 *07:533, 07:708 *07:01:02:01-07:01:02:14, 07:01:27, 07:01:52, 07:01:80-07:01:81, 07:01:83, 07:01:93, 07:01:103, 07:01:112-07:01:113, 07:73:02, 07:89, 07:93:01, 07:141:03, 07:450:02, 07:453, 07:514:02, 07:547, 07:627, 07:668, 07:806, 07:814, 07:856N, 07:877, 07:903, 07:933, 07:985:02, 07:993:01, 07:1006, 07:1042N	B*08:109
	485 bp		*07:419	
78	435 bp	1070 bp	*07:01:13, 07:02:128, 07:68, 07:76:02, 07:199:01-07:199:02, 07:941	*02:02:13, 02:02:29, 06:101, 06:127:01:01-06:127:02, 06:136, 06:144, 06:304, 12:02:11, 18:01:01:01-18:04, 18:06-18:16, B*15:434
79	170 bp	1070 bp	*07:600:01N-07:600:02N, 07:603N	
80⁶	155 bp 345 bp	1070 bp	*07:350N *07:593N	
81⁶	150 bp	1070 bp	*07:231	*16:10, B*44:247, B*44:295:01-44:295:02
82⁵	170 bp 220 bp	1070 bp	*07:288 *07:393N, 07:483N	



101.615-24/04 – including *Taq* polymerase
101.615-24u/04u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **7N8**

Lot-specific Information

83⁶	170 bp 215 bp 250 bp	1070 bp	*07:214, 07:429 *07:445, 07:476N *07:361	
84	130 bp 185 bp	1070 bp	*07:386, 07:439, 07:700 *07:107, 07:224, 07:633N	*06:22, 06:323N, 12:293, 14:119, 16:31, B*15:193 , B*35:283 , B*35:461N , B*46:07N
	225 bp		*07:311, 07:551N	
85⁶	200 bp 275 bp	1070 bp	*07:83, 07:458 *07:360, 07:582Q, 07:649	
86	135 bp 285 bp	1070 bp	*07:58, 07:573 *07:226Q, 07:582Q, 07:713Q	
87	175 bp 210 bp 275 bp	1070 bp	*07:339, 07:410 *07:411, 07:507N *07:340, 07:992	*06:352
88	150 bp	1070 bp	*07:452N, 07:491:01N-07:491:02N	*01:111N, 05:259N, 06:301N, 14:47:01N- 14:47:02N, B*15:209N , B*39:139N , B*40:507N , B*46:55N
	205 bp 285 bp		*07:507N *07:407	
89⁴	75 bp 220 bp	1070 bp	*07:229, 07:387 *07:272, 07:326, 07:374, 07:898, 07:1037N	*04:273 ^w , 12:131 *06:18 ^w , 06:120 ^w , 18:12 ^w
90	140 bp	1070 bp	*07:102, 07:351	*03:137, 03:182, 03:459, 03:504, B*07:13 , B*07:247 , B*07:270 , B*27:198 , B*39:121 , B*54:02
91⁶	160 bp	1070 bp	*07:153, 07:324, 07:738	
92⁴	100 bp	1070 bp	*07:104N, 07:457, 07:702N, 07:893	*04:127, 05:84, 05:180N, 08:208N, B*15:272N , B*39:111 , B*40:286N
	170 bp		*07:119, 07:784	*02:99, 05:157, 08:174, B*07:87 , B*08:224 , B*27:55 , B*27:149
93	155 bp	1070 bp	*07:491:01N-07:491:02N	*01:111N, 05:259N, 06:301N, 14:47:01N- 14:47:02N, B*15:209N , B*39:139N , B*40:507N , B*46:55N
	200 bp		*07:453	
94	530 bp	1070 bp	*07:01:01:01-07:01:02:14, 07:01:04-07:05, 07:07-07:17:04, 07:18:04, 07:20:01-07:33N, 07:35-07:294, 07:296-07:330:01, 07:331- 07:342:02, 07:344-07:347N, 07:349-07:516, 07:518-07:606, 07:608-07:612, 07:614-07:619, 07:622-07:686N, 07:688-07:705, 07:707- 07:720, 07:722-07:734, 07:736:01:01-07:738, 07:740-07:759, 07:761-07:775, 07:777-07:780, 07:782, 07:784-07:802, 07:804N, 07:806- 07:814, 07:816-07:818, 07:820N-07:826, 07:828-07:829, 07:831-07:840N, 07:842- 07:844, 07:846:01-07:846:02, 07:848-07:886N,	*12:181



101.615-24/04 – including *Taq* polymerase
101.615-24u/04u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **7N8**

Lot-specific Information

			07:888-07:918, 07:920-07:931, 07:933-07:951, 07:953-07:954N, 07:956-07:958, 07:960- 07:976, 07:978, 07:980-07:1004, 07:1006- 07:1034, 07:1037N-07:1042N
95⁶	245 bp	1070 bp	*07:316, 07:359, 07:533, 07:607, 07:708
96⁸	-	-	Negative Control

¹Alleles are assigned by the presence of specific PCR product(s). However, the sizes of the specific PCR products may be helpful in the interpretation of HLA-C*07 SSP typings.

When the primers in a primer mix can give rise to HLA-specific PCR products of more than one length this is indicated if the size difference is more than 20 base pairs. Size differences of 20 base pairs or less are not given. For high resolution SSP kits, the alleles listed are specified according to amplicon length.

Nonspecific amplifications, i.e. a ladder or a smear of bands, may sometimes be seen. GC-rich primers have a higher tendency of giving rise to nonspecific amplifications than other primers.

PCR fragments longer than the control bands may sometimes be observed. Such bands should be disregarded and do not influence the interpretation of the SSP typings.

PCR fragments migrating faster than the control bands, but slower than a 400 bp fragment may be seen in some gel read-outs. Such bands can be disregarded and do not influence the interpretation of the SSP typings.

Some primers may give rise to primer oligomer artifacts. Sometimes this phenomenon is an inherent feature of the primer pair(s) of a primer mix. More often it is due to other factors such as too low amount of DNA in the PCR reactions, taking too long time in setting up the PCR reactions, working at elevated room temperature or using thermal cyclers that are not pre-heated.

²The internal positive control primer pairs amplify segments of the human growth hormone gene. The internal positive control bands are 1070 or 800 base pairs respectively, well distribution as outlined in the table. Well number 1 contains the shorter, 800 bp, internal positive control band. The well distribution of the internal controls can help in orientation of the kit on gel photo, as well as allow for kit identification. In the presence of a specific amplification the intensity of the control band often decreases.

³For several HLA Class I alleles 1st and/or 4th exon(s) and beyond, as well as intron nucleotide sequences, are not available. In these instances it is not known whether some of the primers of the SSP sets are completely matched with the target sequences or not. Assumption is made that unknown sequences in these regions are conserved within allelic groups.

⁴HLA-specific PCR products shorter than 125 base pairs have a lower intensity and are less sharp than longer PCR products.

⁵Primer mixes 19, 20, 34, 36, 43, 44, 45, 46, 47, 54 and 82 have a tendency to giving rise to primer oligomer formation.

⁶Primer mixes 5, 13, 17, 20, 23, 24, 30, 37, 40, 45, 56, 64, 65, 66, 70, 74, 77, 80, 81, 83, 85, 91 and 95 may have tendencies of unspecific amplifications, most pronounced for primer mix 95.

⁷Primer mixes 15 and 16 may give rise to a lower yield of HLA-specific PCR product than the other C*07 primer mixes.

⁸Primer mix 96 contains a negative control, which will amplify the majority of HLA amplicons as well as the amplicons generated by the control primer pairs matching the human growth hormone gene. HLA-specific PCR product sizes range from 75 to 200 base pairs and the PCR product generated by the HGH positive control primer pair is 200 base pairs.

Abbreviations

w: might be weakly amplified.



101.615-24/04 – including *Taq* polymerase
101.615-24u/04u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 7N8

Lot-specific Information

PRIMER SPECIFICATION

Well No.	1	2	3	4	5	6	7	8	9	10	11	12
Length of spec. PCR product	245	290	110	185	95	100	190	285	215	185	285	130
	425				170			505	425		305	
					245							
Length of int. pos. control ¹	800	800	800	1070	800	1070	1070	1070	800	1070	800	1070
5'-primer(s) ²	47	363	267	368	47	540	363	355	47	368	355	512
	5'-Agg 3'	5'-AAT 3'	5'-CAA 3'	5'-gTC 3'	5'-Agg 3'	5'-gAC 3'	5'-AAT 3'	5'-CCT 3'	5'-Agg 3'	5'-gTT 3'	5'-TCA 3'	5'-gTg 3'
	648	368	270	368	486			992			365	
	5'-CAC 3'	5'-gTA 3'	5'-AAC 3'	5'-gTC 3'	5'-ACC 3'			5'-TAA 3'			5'-gAC 3'	
											385	
											5'-ggT 3'	
3'-primer(s) ³	302	618	341	512	121	601	512	601	90	512	618	601
	5'-ggC 3'	5'-CT 3'	5'-CgT 3'	5'-CCA 3'	5'-ggC 3'	5'-CTT 3'	5'-CCA 3'	5'-CTC 3'	5'-Ag 3'	5'-CCA 3'	5'-CT 3'	5'-CTT 3'
	853			513	539			1016	302			
	5'-CAT 3'			5'-TCT 3'	5'-TCT 3'			5'-CAC 3'	5'-ggT 3'			
					616							
					5'-CgT 3'							
Well No.	1	2	3	4	5	6	7	8	9	10	11	12

Well No.	13	14	15	16	17	18	19	20	21	22	23	24
Length of spec. PCR product	85	215	80	240	130	120	220	145	100	205	155	135
			290			230	465		225	250	185	
			460			265		190				
Length of int. pos. control ¹	1070	1070	1070	800	1070	800	800	1070	800	800	1070	1070
5'-primer(s) ²	289	97	47	419	512	368	47	463	368	134	199	524
	5'-Agg 3'	5'-TCg 3'	5'-Agg 3'	5'-gTT 3'	5'-gTT 3'	5'-gTC 3'	5'-Agg 3'	5'-TgA 3'	5'-gTC 3'	5'-CCA 3'	5'-gTT 3'	5'-CCA 3'
	1049	97	369			3 rd I		486	5 th I	176	504	
	5'-g 3'	5'-TTg 3'	5'-TAC 3'			5'-Cgg 3'		5'-TAC 3'	5'-gTC 3'	5'-gCC 3'	5'-CAC 3'	
			1049					490				
			5'-C 3'					5'-CgA 3'				
								512				
								5'-gTg 3'				
3'-primer(s) ³	343	270	337	618	601	445	97	618	427	341	341	618
	5'-T 3'	5'-TAC 3'	5'-CTT 3'	5'-CT 3'	5'-CTT 3'	5'-ggA 3'	5'-gTA 3'	5'-CT 3'	5'-gTg 3'	5'-CgT 3'	5'-CgT 3'	5'-CT 3'
	343		616			559	341		1043		618	
	5'-g 3'		5'-TgT 3'			5'-CTC 3'	5'-Cgg 3'		5'-CAA 3'		5'-CT 3'	
	1087		1087			668						
	5'-AgC 3'		5'-AgT 3'			5'-TgA 3'						
						779						
						5'-CTA 3'						
						812						
						5'-CCg 3'						
Well No.	13	14	15	16	17	18	19	20	21	22	23	24



101.615-24/04 – including *Taq* polymerase
101.615-24u/04u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **7N8**

Lot-specific Information

Well No.	25	26	27	28	29	30	31	32	33	34	35	36
Length of spec.	105	100	130	175	255	60	420	115	205	335	320	265
PCR product	150	250	240	260		260		155	230		360	
	235	300		290				245				
								285				
Length of int. pos. control ¹	1070	1070	1070	800	1070	800	1070	1070	1070	1070	1070	800
5'-primer(s) ²	206	361	92	368	122	124	47	97	431	341	47	47
	5'-AgA 3'	5'-AgT 3'	5'-gTg 3'	5'-gTC 3'	5'-CCA 3'	5'-gCT 3'	5'-Agg 3'	5'-TCg 3'	5'-CgT 3'	5'-ggA 3'	5'-Agg 3'	5'-Agg 3'
	419	411	193	390	410	322			688			
	5'-gTC 3'	5'-TAg 3'	5'-CCA 3'	5'-gAA 3'	5'-gTT 3'	5'-gCC 3'			5'-CCA 3'			
		560	539	410								
		5'-CAC 3'	5'-gCg 3'	5'-gTg 3'								
				486								
				5'-TAC 3'								
3'-primer(s) ³	270	618	289	618	341	341	289	172	618	387	197	142
	5'-TAC 3'	5'-CT 3'	5'-AgC 3'	5'-CT 3'	5'-CgT 3'	5'-CgT 3'	5'-AgT 3'	5'-CAT 3'	5'-CT 3'	5'-TTg 3'	5'-gAT 3'	5'-TgC 3'
	527		618		618		296	212	853		235	
	5'-CCC 3'		5'-CT 3'		5'-CT 3'		5'-CTT 3'	5'-gCC 3'	5'-CAT 3'		5'-CTg 3'	
	614							302				
	5'-CgA 3'							5'-ggT 3'				
	3 rd I							341				
	5'-TgC 3'							5'-Tgg 3'				
Well No.	25	26	27	28	29	30	31	32	33	34	35	36
Well No.	37	38	39	40	41	42	43	44	45	46	47	48
Length of spec.	235	110	100	120	140	100	105	165	240	315	200	155
PCR product	265	150	195	205	170	205	150	195	475	350	255	220
			255	440	255	250	175	290		420	420	255
							220				455	410
Length of int. pos. control ¹	1070	1070	800	800	1070	1070	800	1070	1070	1070	1070	1070
5'-primer(s) ²	368	355	125	47	125	409	97	368	341	47	47	47
	5'-gTC 3'	5'-CCC 3'	5'-CgT 3'	5'-Agg 3'	5'-CgA 3'	5'-ggC 3'	5'-TCg 3'	5'-gTC 3'	5'-ggA 3'	5'-Agg 3'	5'-Agg 3'	5'-Agg 3'
	388		697	486	211	410	512	463	652		457	439
	5'-CCT 3'		5'-TCg 3'	5'-ACC 3'	5'-Agg 3'	5'-gTg 3'	5'-gTT 3'	5'-TgT 3'	5'-CCA 3'		5'-AgC 3'	5'-ATC 3'
			794		520	445		464				503
			5'-CTC 3'		5'-CgA 3'	5'-TCA 3'		5'-gCT 3'				5'-CCg 3'
						463		493				
						5'-TgT 3'		5'-CTT 3'				
						560						
						5'-CgA 3'						
3'-primer(s) ³	563	425	341	79	341	618	212	618	527	190	131	131
	5'-CgA 3'	5'-CgC 3'	5'-CgT 3'	5'-gTA 3'	5'-CgT 3'	5'-CT 3'	5'-gCC 3'	5'-CT 3'	5'-CCA 3'	5'-ggT 3'	5'-ggT 3'	5'-ggT 3'
	601	464	853	312	618		231		853	226	295	275
	5'-CTT 3'	5'-AgA 3'	5'-CAT 3'	5'-AgT 3'	5'-CT 3'		5'-Tgg 3'		5'-CAT 3'	5'-CAA 3'	5'-TCA 3'	5'-gCC 3'
	603			317			233		295	330	296	
	5'-TTg 3'			5'-CgA 3'			5'-CCC 3'		5'-TCA 3'	5'-TgT 3'	5'-CTT 3'	
				317			277			618	618	
				5'-CAA 3'			5'-gCT 3'			5'-CT 3'	5'-CT 3'	
				563			578					
				5'-Cgg 3'			5'-TgT 3'					
							613					
							5'-gCC 3'					
Well No.	37	38	39	40	41	42	43	44	45	46	47	48



101.615-24/04 – including *Taq* polymerase
101.615-24u/04u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **7N8**

Lot-specific Information

Well No.	49	50	51	52	53	54	55	56	57	58	59	60
Length of spec.	70	345	65	65	70	120	160	210	70	80	110	130
PCR product	155	425	235	205	110	200	210	260	175	450	210	300
	230		285	280	155	235	290		230			370
						285	530		270			
						340						
						370						
						430						
Length of int.	800	800	1070	1070	1070	1070	1070	1070	800	1070	800	1070
pos. control ¹												
5'-primer(s) ²	89	47	368	179	503	47	91	388	270	47	359	47
	5'-gAA 3'	5'-Agg 3'	5'-gTC 3'	5'-gTC 3'	5'-CCg 3'	5'-Agg 3'	5'-ggg 3'	5'-CCT 3'	5'-AAT 3'	5'-Agg 3'	5'-CCT 3'	5'-Agg 3'
	176			317	540	648	173	397	388	512	454	359
	5'-gCA 3'			5'-gCA 3'	5'-gAC 3'	5'-CAC 3'	5'-Cgg 3'	5'-gCg 3'	5'-CCT 3'	5'-gTT 3'	5'-ACT 3'	5'-CCT 3'
	429			458	553		211	409	429		463	527
	5'-gAA 3'			5'-ggT 3'	5'-TgC 3'		5'-Agg 3'	5'-ggC 3'	5'-gAA 3'		5'-TgA 3'	5'-TgA 3'
	440				589		232	450	485		464	
	5'-TTT 3'				5'-TgA 3'		5'-AgA 3'	5'-CTg 3'	5'-CAA 3'		5'-gCA 3'	
	589						992					
	5'-TgA 3'						5'-TAA 3'					
3'-primer(s) ³	289	221	394	341	618	112	341	618	299	327	527	180
	5'-AgC 3'	5'-ACT 3'	5'-gCT 3'	5'-CgT 3'	5'-CT 3'	5'-CCA 3'	5'-CgT 3'	5'-CT 3'	5'-TCA 3'	5'-TTT 3'	5'-CCg 3'	5'-TCC 3'
	289	299	553	3 rd I		218	1043		618	544		246
	5'-AgC 3'	5'-TCT 3'	5'-CTA 3'	5'-CTC 3'		5'-gCT 3'	5'-CAg 3'		5'-CT 3'	5'-ggT 3'		5'-TAA 3'
	618		563			244				559		618
	5'-CT 3'		5'-CgA 3'			5'-CTg 3'				5'-CAg 3'		5'-CT 3'
			569			727						
			5'-ACC 3'			5'-CCA 3'						
			614			809						
			5'-CAg 3'			5'-CAA 3'						
						893						
						5'-CCT 3'						
						912						
						5'-gTA 3'						
Well No.	49	50	51	52	53	54	55	56	57	58	59	60



101.615-24/04 – including *Taq* polymerase
101.615-24u/04u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **7N8**

Lot-specific Information

Well No.	61	62	63	64	65	66	67	68	69	70	71	72
Length of spec.	120	200	275	135	150	150	180	75	135	215	105	110
PCR product	245	285	320	230	185	185	215	230	185	345	120	245
					215	225	290	260	240	405	240	290
								290			280	
Length of int.	1070	1070	800	1070	1070	800	1070	1070	1070	1070	1070	1070
pos. control ¹												
5'-primer(s) ²	355	92	47	89	419	419	368	368	419	47	380	368
	5'-CCC 3'	5'-gTg 3'	5'-Agg 3'	5'-gAA 3'	5'-gTC 3'	5'-gTC 3'	5'-gTC 3'	5'-gTC 3'	5'-gTC 3'	5'-gTC 3'	5'-Agg 3'	5'-CCC 3'
		176	380	113							419	415
		5'-gCC 3'	5'-CCC 3'	5'-CCA 3'							5'-gTA 3'	5'-ACT 3'
		185	385	193							538	538
		5'-CAA 3'	5'-ggT 3'	5'-CCA 3'							5'-A.C 3'	5'-AgT 3'
		458		368							554	544
		5'-ggT 3'		5'-gTC 3'							5'-ggT 3'	5'-gAA 3'
												554
												5'-ggT 3'
3'-primer(s) ³	425	341	146	289	530	530	505	403	506	85	618	618
	5'-CgC 3'	5'-CgT 3'	5'-CCg 3'	5'-AgC 3'	5'-CCA 3'	5'-CCA 3'	5'-gCT 3'	5'-gCC 3'	5'-Tgg 3'	5'-CAA 3'	5'-CT 3'	5'-CT 3'
	445	3 rd I	190	550	535	556	544	557	512	85		
	5'-ggA 3'	5'-CTC 3'	5'-ggT 3'	5'-CAT 3'	5'-CTA 3'	5'-gCT 3'	5'-gC 3'	5'-TgT 3'	5'-CAT 3'	5'-CAg 3'		
	448		202	569	563	564	616	589	563	100		
	5'-CAC 3'		5'-TCA 3'	5'-ACC 3'	5'-CgT 3'	5'-ACC 3'	5'-TgT 3'	5'-CTT 3'	5'-CgT 3'	5'-ggg 3'		
	559		618		565	564		616	3 rd I	221		
	5'-Cgg 3'		5'-CT 3'		5'-CAA 3'	5'-ACT 3'		5'-TgT 3'	5'-TgC 3'	5'-ACT 3'		
					595	604				280		
					5'-g 3'	5'-CgC 3'				5'-CTA 3'		
					596							
					5'-CC 3'							
Well No.	61	62	63	64	65	66	67	68	69	70	71	72

Well No.	73	74	75	76	77	78	79	80	81	82	83	84
Length of spec.	185	210	140	115	200	435	170	155	150	170	170	130
PCR product	210	355	220	280	240			345		220	215	185
			265		430						250	225
			295		485							
Length of int.	1070	1070	1070	1070	1070	1070	1070	1070	1070	1070	1070	1070
pos. control ¹												
5'-primer(s) ²	419	47	512	103	648	213	97	3 rd I	418	368	97	419
	5'-gTC 3'	5'-Agg 3'	5'-gTT 3'	5'-CCT 3'	5'-CAC 3'	5'-CCC 3'	5'-TCg 3'	5'-Cgg 3'	5'-Agg 3'	5'-gTC 3'	5'-TCg 3'	5'-gTC 3'
	1018		3 rd I	268								
	5'-gTg 3'		5'-Cgg 3'	5'-gA 3'								
3'-primer(s) ³	564	85	610	341	809	355	224	704	527	496	226	506
	5'-ACC 3'	5'-CAC 3'	5'-CTA 3'	5'-CgT 3'	5'-CAA 3'	5'-gAg 3'	5'-CCT 3'	5'-CCg 3'	5'-CCg 3'	5'-gAg 3'	5'-CAT 3'	5'-TgT 3'
	565	90	614		845	355	225	893		544	265	508
	5'-CAT 3'	5'-Ag 3'	5'-CgT 3'		5'-ACA 3'	5'-gAg 3'	5'-ACT 3'	5'-CCT 3'		5'-gC 3'	5'-CTA 3'	5'-CTg 3'
	596	233	768		912					553	277	564
	5'-CC 3'	5'-CCC 3'	5'-gTg 3'		5'-gTA 3'					5'-CTA 3'	5'-gCT 3'	5'-ACT 3'
	1055		812		967						308	565
	5'-TgC 3'		5'-CCg 3'		5'-gAT 3'						5'-gCg 3'	5'-CAT 3'
	1076		845									595
	5'-ACA 3'		5'-ACA 3'									5'-g 3'
												614
												5'-CgA 3'
Well No.	73	74	75	76	77	78	79	80	81	82	83	84



101.615-24/04 – including *Taq* polymerase
101.615-24u/04u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **7N8**

Lot-specific Information

Well No.	85	86	87	88	89	90	91	92	93	94	95
Length of spec.	200	135	175	150	75	140	160	100	155	530	245
PCR product	275	285	210	205	220			170	200		
			275	285							
Length of int.	1070	1070	1070	1070	1070	1070	1070	1070	1070	1070	1070
pos. control ¹											
5'-primer(s) ²	373	373	97	97	112	103	732	361	47	992	648
	5'-gCg 3'	5'-gCA 3'	5'-TCg 3'	5'-TCg 3'	5'-CCT 3'	5'-CCT 3'	5'-gAT 3'	5'-AgA 3'	5'-Agg 3'	5'-TAT 3'	5'-CAC 3'
	385	373		355	113				355		
	5'-ggA 3'	5'-gCg 3'		5'-CCC 3'	5'-CCA 3'				5'-CCC 3'		
	458	524			257						
	5'-ggg 3'	5'-CCC 3'			5'-CCA 3'						
	1018										
	5'-gTg 3'										
3'-primer(s) ³	618	618	231	208	289	201	853	422	74	1043	845
	5'-CT 3'	5'-CT 3'	5'-Tgg 3'	5'-CTA 3'	5'-AgC 3'	5'-CTC 3'	5'-CAT 3'	5'-AgA 3'	5'-AgT 3'	5'-CAg 3'	5'-ACA 3'
	1076		259	259				426	470		862
	5'-ACA 3'		5'-CTA 3'	5'-CTA 3'				5'-TCC 3'	5'-TCT 3'		5'-CTT 3'
			278	340				488	471		
			5'-ggT 3'	5'-gTT 3'				5'-CCA 3'	5'-gTT 3'		
			332	470							
			5'-TCC 3'	5'-TCT 3'							
				471							
				5'-gTT 3'							
Well No.	85	86	87	88	89	90	91	92	93	94	95

¹The internal positive control primer pairs amplify segments of the human growth hormone gene. The internal positive control bands are 1070 or 800 base pairs respectively, well distribution as outlined in the table. Well number 1 contains the shorter, 800 bp, internal positive control band. The well distribution of the internal controls can help in orientation of the kit on gel photo, as well as allow for kit identification. In the presence of a specific amplification the intensity of the control band often decreases.

²The nucleotide position matching the specificity-determining 3'-end of the primer is given. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

³The nucleotide position matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.



101.615-24/04 – including *Taq* polymerase
101.615-24u/04u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 7N8

Lot-specific Information

[illegible]

101.615-24/04 – including *Taq* polymerase
101.615-24u/04u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 7N8

Lot-specific Information

[illegible]

101.615-24/04 – including *Taq* polymerase
101.615-24u/04u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 7N8

Lot-specific Information

CELL LINE VALIDATION SHEET																											
HLA-C*07 SSP subtyping kit ²																											
					Well																						
					33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48							
				Prod. No.:	202016533	202016534	202016535	202016536	202016537	202134238	202016539	202134240	202016541	202016542	202016543	202016544	202016545	202016546	202016547	202016548							
	IHC cell line ¹		C*																								
1	9001	SA	*07:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
2	9280	LK707	*07:01	*15:05	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-		
3	9011	E4181324	*12:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
4	9275	GU373	*03:04	*04:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
5	9009	KAS011	*06:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
6	9353	SM	*03:04	*07:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
7	9020	QBL	*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
8	9025	DEU	*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
9	9026	YAR	*12:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
10	9107	LKT3	*01:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
11	9051	PITOUT	*16:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
12	9052	DBB	*06:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
13	9004	JESTHOM	*01:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
14	9071	OLGA	*01:02	*03:04	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
15	9075	DKB	*03:04		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
16	9037	SWEG007	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
17	9282	CTM3953540	*03:03	*07:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
18	9257	32367	*01:02	*07:05	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
19	9038	BM16	*07:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
20	9059	SLE005	*03:04		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
21	9064	AMALA	*03:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
22	9056	KOSE	*12:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
23	9124	IHL	*01:02	*15:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
24	9035	JBUSH	*12:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
25	9049	IBW9	*08:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
26	9285	WT49	*07:18		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
27	9191	CH1007	*07:04	*15:29	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
28	9320	BEL5GB	*05:01	*16:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
29	9050	MOU	*16:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
30	9021	RSH	*17:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
31	9019	DUCAF	*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
32	9297	HAG	*17:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
33	9098	MT14B	*03:04		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
34	9104	DHIF	*12:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
35	9302	SSTO	*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
36	9024	KT17	*03:03	*04:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
37	9065	HHKB	*07:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
38	9099	LZL	*03:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
39	9315	CML	*02:02	*07:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
40	9134	WHONP199	*01:02	*06:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
41	9055	H0301	*08:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
42	9066	TAB089	*01:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
43	9076	T7526	*01:02	*08:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
44	9057	TEM	*12:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
45	9239	SHJO	*06:02	*17:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
46	9013	SCHU	*07:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
47	9045	TUBO	*07:04	*15:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
48	9303	TER-ND	*04:01	*16:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		



101.615-24/04 – including *Taq* polymerase
101.615-24u/04u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 7N8

Lot-specific Information

CELL LINE VALIDATION SHEET																							
HLA-C*07 SSP subtyping kit ²																							
					Well																		
					49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64			
				Prod. No.:	202016549	202016550	202016551	202016552	202016553	202016554	202016555	202016556	202016557	202016558	202134259	202016560	202134261	202016562	202016563	202016564			
	IHCW cell line ¹		C*																				
1	9001	SA	*07:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
2	9280	LK707	*07:01	*15:05	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
3	9011	E4181324	*12:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
4	9275	GU373	*03:04	*04:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
5	9009	KAS011	*06:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
6	9353	SM	*03:04	*07:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
7	9020	QBL	*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
8	9025	DEU	*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
9	9026	YAR	*12:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
10	9107	LKT3	*01:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
11	9051	PITOUT	*16:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
12	9052	DBB	*06:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
13	9004	JESTHOM	*01:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
14	9071	OLGA	*01:02	*03:04	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
15	9075	DKB	*03:04		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
16	9037	SWEIG007	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
17	9282	CTM3953540	*03:03	*07:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
18	9257	32367	*01:02	*07:05	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
19	9038	BM16	*07:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
20	9059	SLE005	*03:04		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
21	9064	AMALA	*03:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
22	9056	KOSE	*12:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
23	9124	IHL	*01:02	*15:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
24	9035	JBUSH	*12:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
25	9049	IBW9	*08:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
26	9285	WT49	*07:18		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
27	9191	CH1007	*07:04	*15:29	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
28	9320	BEL5GB	*05:01	*16:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
29	9050	MOU	*16:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
30	9021	RSH	*17:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
31	9019	DUCAF	*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
32	9297	HAG	*17:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
33	9098	MT14B	*03:04		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
34	9104	DHIF	*12:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
35	9302	SSTO	*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
36	9024	KT17	*03:03	*04:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
37	9065	HHKB	*07:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
38	9099	LZL	*03:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
39	9315	CML	*02:02	*07:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
40	9134	WHONP199	*01:02	*06:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
41	9055	H0301	*08:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
42	9066	TAB089	*01:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
43	9076	T7526	*01:02	*08:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
44	9057	TEM	*12:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
45	9239	SHJO	*06:02	*17:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
46	9013	SCHU	*07:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
47	9045	TUBO	*07:04	*15:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
48	9303	TER-ND	*04:01	*16:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-



101.615-24/04 – including *Taq* polymerase
101.615-24u/04u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 7N8

Lot-specific Information

CELL LINE VALIDATION SHEET																						
HLA-C*07 SSP subtyping kit ²																						
					Well																	
					65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80		
				Prod. No.:	202016565	202016566	202016567	202016568	202016569	202016570	202134271	202016572	202016573	202016574	202134275	202016576	202134277	202016578	202016579	202134280		
	IHCW cell line ¹		C*																			
1	9001	SA	*07:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
2	9280	LK707	*07:01	*15:05	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
3	9011	E4181324	*12:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
4	9275	GU373	*03:04	*04:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
5	9009	KAS011	*06:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
6	9353	SM	*03:04	*07:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
7	9020	QBL	*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
8	9025	DEU	*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
9	9026	YAR	*12:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
10	9107	LKT3	*01:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
11	9051	PITOUT	*16:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
12	9052	DBB	*06:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
13	9004	JESTHOM	*01:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
14	9071	OLGA	*01:02	*03:04	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
15	9075	DKB	*03:04		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
16	9037	SWEIG007	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
17	9282	CTM3953540	*03:03	*07:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
18	9257	32367	*01:02	*07:05	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
19	9038	BM16	*07:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
20	9059	SLE005	*03:04		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
21	9064	AMALA	*03:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
22	9056	KOSE	*12:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
23	9124	IHL	*01:02	*15:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
24	9035	JBUSH	*12:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
25	9049	IBW9	*08:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
26	9285	WT49	*07:18		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
27	9191	CH1007	*07:04	*15:29	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
28	9320	BEL5GB	*05:01	*16:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
29	9050	MOU	*16:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
30	9021	RSH	*17:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
31	9019	DUCAF	*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
32	9297	HAG	*17:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
33	9098	MT14B	*03:04		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
34	9104	DHIF	*12:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
35	9302	SSTO	*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
36	9024	KT17	*03:03	*04:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
37	9065	HHKB	*07:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
38	9099	LZL	*03:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
39	9315	CML	*02:02	*07:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
40	9134	WHONP199	*01:02	*06:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
41	9055	H0301	*08:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
42	9066	TAB089	*01:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
43	9076	T7526	*01:02	*08:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
44	9057	TEM	*12:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
45	9239	SHJO	*06:02	*17:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
46	9013	SCHU	*07:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
47	9045	TUBO	*07:04	*15:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
48	9303	TER-ND	*04:01	*16:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		



For *In Vitro* Diagnostic Use
MA124 v02 SSP PI Template
Date: January 2023, Rev. No: 00

101.615-24/04 – including *Taq* polymerase
101.615-24u/04u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **7N8**

Lot-specific Information

CELL LINE VALIDATION SHEET																				
HLA-C*07 SSP subtyping kit ²																				
						Well														
						81	82	83	84	85	86	87	88	89	90	91	92	93	94	95
					Prod. No.:	202016581	202016582	202016583	202016584	202016585	202016586	202016587	202016588	202016589	202016590	202016591	202016592	202016593	202016594	202237995
	IHCW cell line ¹		C*																	
1	9001	SA	*07:02			-	-	-	-	-	-	-	-	-	-	-	-	-	+	-
2	9280	LK707	*07:01	*15:05		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-
3	9011	E4181324	*12:02			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
4	9275	GU373	*03:04	*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
5	9009	KAS011	*06:02			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
6	9353	SM	*03:04	*07:02		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-
7	9020	QBL	*05:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
8	9025	DEU	*04:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
9	9026	YAR	*12:03			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
10	9107	LKT3	*01:02			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
11	9051	PITOUT	*16:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
12	9052	DBB	*06:02			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
13	9004	JESTHOM	*01:02			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
14	9071	OLGA	*01:02	*03:04		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
15	9075	DKB	*03:04			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
16	9037	SWEIG007	*02:02			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
17	9282	CTM3953540	*03:03	*07:01		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-
18	9257	32367	*01:02	*07:05		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-
19	9038	BM16	*07:01			-	-	-	-	-	-	-	-	-	-	-	-	-	+	-
20	9059	SLE005	*03:04			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
21	9064	AMALA	*03:03			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
22	9056	KOSE	*12:03			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
23	9124	IHL	*01:02	*15:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
24	9035	JBUSH	*12:03			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
25	9049	IBW9	*08:02			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
26	9285	WT49	*07:18			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
27	9191	CH1007	*07:04	*15:29		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-
28	9320	BEL5GB	*05:01	*16:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
29	9050	MOU	*16:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
30	9021	RSH	*17:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
31	9019	DUCAF	*05:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
32	9297	HAG	*17:03			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
33	9098	MT14B	*03:04			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
34	9104	DHIF	*12:03			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
35	9302	SSTO	*05:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
36	9024	KT17	*03:03	*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
37	9065	HHKB	*07:02			-	-	-	-	-	-	-	-	-	-	-	-	-	+	-
38	9099	LZL	*03:03			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
39	9315	CML	*02:02	*07:01		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-
40	9134	WHONP199	*01:02	*06:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
41	9055	H0301	*08:02			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
42	9066	TAB089	*01:02			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
43	9076	T7526	*01:02	*08:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
44	9057	TEM	*12:03			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
45	9239	SHJO	*06:02	*17:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
46	9013	SCHU	*07:02			-	-	-	-	-	-	-	-	-	-	-	-	-	+	-
47	9045	TUBO	*07:04	*15:02		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-
48	9303	TER-ND	*04:01	*16:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-



101.615-24/04 – including *Taq* polymerase
101.615-24u/04u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **7N8**

Lot-specific Information

¹The provided cell line HLA specificities are retrieved from the <http://www.ihwg.org/hla> web site. The specificity of an individual cell line may thus be subject to change.

²The specificity of each primer solution in the kit has been tested against 48 well characterized cell line DNAs and where applicable, additional cell line DNAs.

No DNAs carrying the alleles to be amplified by primer solutions 5, 10, 11, 13, 15, 16, 18, 19, 23 to 31, 34 to 36, 38, 39, 41 to 44, 46 to 93 and 95 were available.

The specificities of the primers in primer solutions 5, 10, 11, 13, 15, 16, 18, 19, 24 to 27, 31, 36, 38, 42 to 44, 47 to 50, 53 to 61, 63, 64, 68, 71, 72, 74, 78, 81, 85 and 89 to 91 were tested by separately adding one or more 5'-primers, and one or more 3'-primers respectively.

In primer solutions 34, 35, 46, 51, 65 to 67, 69, 70, 73, 75, 77, 79, 80, 82 to 84, 87, 88, 92, 93 and 95 it was only possible to test the 5'-primers, the 3'-primers were not possible to be tested.

In primer solutions 23, 28 to 30, 39, 41, 52, 62, 76 and 86 it was only possible to test the 3'-primers, the 5'-primers were not possible to be tested.

In primer solutions 3, 11, 14, 20, 22, 26, 27, 33, 42, 44, 47 to 49, 53, 55 to 57, 59, 60, 63, 64, 71, 72, 85 and 89 one or more 5'-primers were not possible to be tested, and in primer solutions 4, 5, 9, 13, 15, 18, 21, 25, 31, 32, 37, 38, 40, 43, 47, 48, 50, 54, 58, 60, 61, 63, 64, 68, 74 and 85 one or more 3'-primers were not possible to be tested.

In addition, one or more 3'- and/or 5'-primers in primer solutions 2, 8, 20, 21, 33, 37, 40 and 45 were tested by separately adding one 5'-primer and/or one or two 3'-primers.



101.615-24/04 – including *Taq* polymerase
101.615-24u/04u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **7N8**

Lot-specific Information

ADDRESSES:

Manufacturer:

Caredx AB, Franzengatan 5, SE-112 51 Stockholm, Sweden.

Tel: +46-8-508 939 00

Fax: +46-8-717 88 18

E-mail: orders-se@caredx.com

Web page: www.caredx.com

CareDx Lab Solutions Inc., 901 S. Bolmar St., Suite R, West Chester, PA 19382

Tel: 1-877-653-78171

Fax: 610-344-7989

E-mail: orders-us@caredx.com

Web page: www.caredx.com

For information on CareDx distributors worldwide, contact **CareDx AB**.

