

## Olerup SSP<sup>®</sup> DQB1\*02

Product number: 101.213-24 – including *Taq* polymerase  
101.213-24u – without *Taq* polymerase  
Lot number: 63N  
Expiry date: 2014-October-01  
Number of tests: 24  
Number of wells per test: 7

### CHANGES COMPARED TO THE PREVIOUS DQB1\*02 LOT (77K):

Well	5'-primer	3'-primer	rationale
7	New	New	New primer pair for the DQB1*02:06 allele.

THE NUMBER OF WELLS is increased from 6 to 7.

### ALLELE COVERAGE:

DQB1\*02:01 to DQB1\*02:06, i.e. all the currently recognized DQB1\*02 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2012-April-12, release 3.8.0.

### RESOLUTION IN DQB1\*02 HOMO- AND HETEROZYGOTES:

Very good.

### INFLUENCE ON THE INTERPRETATION OF DQB1\*02 SUBTYPINGS BY NON-DQB1\*02 ALLELES:

None of importance.

### MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.