

Olerup SSP® DQB1*03

Product number: 101.214-24/06 – including *Taq* pol.
101.214-24u/06u – without *Taq* pol.
Lot number: 56K
Expiry date: 2013-April-01
Number of tests: 24 test – Product No. 101.214-24
6 tests – Product No. 101.214-06
Number of wells per test: 24

CHANGES COMPARED TO THE PREVIOUS DQB1*03 LOT (22G):

Well	5'-primer	3'-primer	rationale
24	New	New	New primer pair for the DQB1*03:27 allele.

THE NUMBER OF WELLS has been increased from 23 to 24.

ALLELE COVERAGE:

DQB1*03:01 to DQB1*03:27, i.e. all the currently recognized DQB1*03 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2010-July-16, release 3.1.0.

RESOLUTION IN DQB1*03 HOMOZYGOTES:

The DQB1*03:01,03:01, DQB1*03:01,03:02, DQB1*03:02,03:02, DQB1*03:02,03:03 and DQB1*03:03,03:03 genotypes all give rise to unique amplification patterns.

INFLUENCE ON THE INTERPRETATION OF DQB1*03 SUBTYPINGS BY NON-DQB1*03 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions have been received.