

DQB1*05
101.211-24 – including *Taq* polymerase
101.211-24u – without *Taq* polymerase
Lot No.: 42F

Release Note

Page 1 of 1

Olerup SSP[®] DQB1*05

Product number: 101.211-24 – including *Taq* polymerase
101.211-24u – without *Taq* polymerase
Lot number: 42F
Expiry date: 2011-January-01
Number of tests: 24
Number of wells per test: 6

CHANGES COMPARED TO THE PREVIOUS DQB1*05 LOT (09E):

The DQB1*05 primer set is unchanged.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

DQB1*0501 to DQB1*0505, i.e. all the currently recognized DQB1*05 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2009-January-16, release 2.24.0.

RESOLUTION IN DQB1*05 HOMO- AND HETEROZYGOTES:

Very good.

INFLUENCE ON THE INTERPRETATION OF DQB1*05 SUBTYPINGS BY NON-DQB1*05 ALLELES:

None.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.