

101.231-24/04 – including *Taq* polymerase
 101.231-24u/04u – without *Taq* polymerase

Lot No.: **17V**

Olerup SSP® DQA1

Product number: 101.231-24/04 – including *Taq* pol.
 101.231-24u/04u – without *Taq* pol.

Lot number: 17V

Expiry date: 2016-July-01

Number of tests: 24 tests – Product No. 101.231-24/24u
 4 tests – Product No. 101.231-04/04u

Number of wells per test: 32

CHANGES COMPARED TO THE PREVIOUS DQA1 LOT (82R):

Well	5'-primer	3'-primer	rationale
16	Added	Added	Primer pair added from well 31.
29	-	Added	3'-primer added for the DQA1*01:12 allele.
31	Added	Added	Primer pair added from well 32.
32	-	-	Negative Control.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

All the currently recognized DQA1 alleles, i.e. DQA1*01:01 to 01:12, DQA1*02:01, DQA1*03:01 to 03:03, DQA1*04:01 to DQA1*04:04, DQA1*05:01 to 05:11 and DQA1*06:01 to DQA1*06:02, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2013-October-11, release 3.14.0.

RESOLUTION IN DQA1 HOMOZYGOTES:

Good.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS: