

101.231-24/04 – including *Taq* polymerase
 101.231-24u/04u – without *Taq* polymerase

Lot No.: **3F0**

Olerup SSP[®] DQA1

Product number: 101.231-24/04 – including *Taq* pol.
 101.231-24u/04u – without *Taq* pol.
Lot number: 3F0
Expiry date: 2019-12-01
Number of tests: 24 tests – Product No. 101.231-24/24u
 4 tests – Product No. 101.231-04/04u
Number of wells per test: 31+1

CHANGES COMPARED TO THE PREVIOUS DQA1 LOT (9D6):

Well	5'-primer	3'-primer	rationale
19	-	Added	3'-primer added for the DQA1*01:15N allele.
28	Added	Added	Primer pair added for the DQA1*01:14 allele.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

All the currently recognized DQA1 alleles¹, i.e. DQA1*01:01 to 01:15N, DQA1*02:01, DQA1*03:01 to 03:03, DQA1*04:01 to DQA1*04:04, DQA1*05:01 to 05:11 and DQA1*06:01 to 06:02, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2017-January-20, release 3.27.0.

¹Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

RESOLUTION IN DQA1 HOMOZYGOTES:

Good.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.