

101.125-24/06 – including *Taq* polymerase
101.125-24u/06u – without *Taq* polymerase

Lot No.: **09K**

Olerup SSP[®] DRB1*15

Product number: 101.125-24/06 – including *Taq* pol.
101.125-24u/06u – without *Taq* pol.
Lot number: 09K
Expiry date: 2012-April-01
Number of tests: 24 test – Product No. 101.125-24
6 tests – Product No. 101.125-06
Number of wells per test: 30

CHANGES COMPARED TO THE PREVIOUS DRB1*15 LOT (55F):

Well	5'-primer	3'-primer	rationale
8	-	Added	New primer for the DRB1*15:36 allele.
13	-	Added	New primer for the DRB1*15:32 allele.
19	-	Added	New primer for the DRB1*15:37 allele.
22	Added	-	New primer for the DRB1*15:40 allele.
24	Added	-	New primers for the DRB1*15:40 and DRB1*15:43 alleles.
25	New	New	New primer pairs for the DRB1*15:33 and DRB1*15:39 alleles.
26	New	New	New primer pair for the DRB1*15:38 allele.
27	New	New	New primer pair for the DRB1*15:42 allele.
28	New	New	New primer pair for the DRB1*15:35 allele.
29	New	New	New primer pair for the DRB1*15:41 allele.
30	New	New	New primer pair for the DRB1*15:44 allele.

THE NUMBER OF WELLS is increased from 24 to 30.

ALLELE COVERAGE:

DRB1*15:01 to DRB1*15:44, i.e. all the currently recognized DRB1*15 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2010-April-01, release 3.0.0.

RESOLUTION IN DRB1*15 HOMOZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF DRB1*15 SUBTYPINGS BY NON-DRB1*15 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.